

This document was too large to scan as a whole document, therefore it required breaking into smaller sections.

Document number: SD-WM-DP-145

Section 4 of 12

Title: REVISED 60 DAY SAFETY SCREENING +

FERROCYANIDE RESULTS FOR TANK 241-BY-108

ROTARY SAMPLES CORE 98 + CORE 104

Date: 02/02/96 Revision: R001

Originator: BALDWIN JH

Co: WHC

Recipient:

Co:

References: ECN-629001

THIS PAGE INTENTIONALLY
LEFT BLANK.



Westinghouse
Hanford Company

P.O. Box 1970 Richland, WA 99352

WHC-SD-WM-DP-145, REV.1

PART II

WHC-SD-WM-DP-145, REV. 0

ANALYTICAL SERVICES

**60-DAY SAFETY SCREENING AND FERROCYANIDE RESULTS FOR TANK 241-BY-108,
ROTARY SAMPLES, CORE 98 AND CORE 104**

Date Printed: OCTOBER 18, 1995

2-1

INTENTIONALLY
LEFT BLANK.

TABLE OF CONTENTS

Narrative	1
Sample Data Summary	8
Inorganic Analyses	40
Differential Scanning Calorimetry	
DSC Worklist # 2702	41
DSC Worklist # 2867	42
DSC Worklist # 2868	44
DSC Worklist # 2869	46
DSC Worklist # 2870	47
DSC Worklist # 2001	48
DSC Worklist # 2025	54
DSC Worklist # 2026	60
DSC Worklist # 2027	68
DSC Worklist # 2032	76
DSC Worklist # 2038	81
DSC Worklist # 2112	86
DSC Worklist # 2114	92
DSC Worklist # 2115	100
DSC Worklist # 2142	104
DSC Worklist # 2273	112
DSC Worklist # 2317	117
DSC Worklist # 2318	123
DSC Worklist # 2320	131
DSC Worklist # 2321	138
DSC Worklist # 2322	145
DSC Worklist # 2323	152
DSC Worklist # 2324	158
Thermogravimetric Analysis (TGA)	
TGA Worklist # 2002	163
TGA Worklist # 2028	169
TGA Worklist # 2029	173
TGA Worklist # 2030	180
TGA Worklist # 2033	184
TGA Worklist # 2034	192
TGA Worklist # 2035	199
TGA Worklist # 2036	205
TGA Worklist # 2037	210
TGA Worklist # 2155	217
TGA Worklist # 2157	223
TGA Worklist # 2158	227
TGA Worklist # 2159	235

WHC-SD-WM-DP-145, REV. 0
TABLE OF CONTENTS (Continued)

TGA Worklist # 2160	240
TGA Worklist # 2161	248
TGA Worklist # 2200	254
TGA Worklist # 2251	260
TGA Worklist # 2252	264
TGA Worklist # 2269	272
TGA Worklist # 2276	276
TGA Worklist # 2418	281
TGA Worklist # 2419	289
TGA Worklist # 2420	295
TGA Worklist # 2421	301
TGA Worklist # 2422	308
TGA Worklist # 2424	314
TGA Worklist # 2425	320
TGA Worklist # 2426	326
TGA Worklist # 2428	333
TGA Worklist # 2615	337
TGA Worklist # 2685	343

This document consists of pages 1 through 346, plus page 138.1.

WHC-SD-WM-DP-145, REV. 1

WHC-SD-WM-DP-145, REV. 0

NARRATIVE

2-4

THIS PAGE INTENTIONALLY
LEFT BLANK.

60-DAY SAFETY SCREEN RESULTS FOR TANK 241-BY-108
ROTARY SAMPLES, CORE 98 AND CORE 104ANALYTICAL SUMMARY

Core samples from tank 241-BY-108 (BY-108) were received at the 222-S Laboratories and underwent safety screening analyses, consisting of differential scanning calorimetry (DSC), thermogravimetric analysis (TGA), and total alpha activity (*Tank 241-BY-108 Tank Characterization Plan* (TCP)[1]). Results of total organic carbon (TOC) and cyanide (CN) analyses are included for those samples which exceeded DSC action limit.

As required by the *Tank Safety Screening Data Quality Objective* (DQO) [2], a 95% confidence interval was calculated for the sample results exceeding an action limit. The precision requirements of the SAP were satisfied by comparing a one-sided 95% confidence interval of the mean for each sample to the action limit, rather than requiring a relative percent difference between sample and duplicate results of less than 10%. The DSC analysis at the 95% confidence level found the DSC results of three samples exceeded the maximum limit stated in the DQO (Table 4). Notifications by the Chemist and Project Coordinator were made as required.

Before samples were removed from BY-108, an industrial hygiene technician field tested the tank vapors. The technician purged the vapor probe sample tube for five minutes than field measured vapor stream contents using a combustible gas indicator (CGI) and an organic vapor meter (OVM). The technician verbally reported an LEL of 5.0%, an oxygen content of 20 R, and a total organic carbon content of 71.8 ppm using the OVM. Draeger tubes were not used to estimate the ammonia content.

When compared to the decision rules in the DQO, none of the data indicate that the tank should be considered "unsafe." The tank can be considered "safe" once it has been determined that no flammability concern exists in the tank vapor space.

SCOPE

This document serves as the 60-day report deliverable for the tank BY-108 core samples collected on July 27 through August 16, 1995 (Core 98, Segments 1-4 and Core 104 Segments 1-5). The 222-S Laboratories received, extruded, and analyzed each sample in accordance with the TCP. Included in this report are the primary safety screening results obtained from the analyses, and copies of all DSC and TGA raw data scans as requested in the SAP. Also included are the TOC and CN results for those samples that exceeded the DSC action limit. Any additional analyses conducted by the 222-S Laboratories on the BY-108 core samples will be included in a revision to this report. Partial BY-108 core segments from cores 97, 100 and 102 were received from BY-108 but were not analyzed as they duplicated segments of cores 98 and 104. Core 99, a core taken from the same riser as core 98, was sent to PNL's 325 Laboratory for analysis and the results of those analyses will be reported by PNL in a separate document.

SAMPLE RECEIPT, EXTRUSION, AND SUBSAMPLING

The subject core samples reported herein were taken from Tank BY-108. The two samples are identified as Core 98 and Core 104. These core segment samples were received at the 222-S Laboratory between July 28 and August 21, 1995 and extruded between August 1 and August 24, 1995. Table 1 provides the sampling and extrusion report for the two core samples. With the exception of segment 1, of both cores, the recovery was good and each segment was broken into quarter segments and homogenized. Segment 1, for both cores, was treated as a "whole" segment and homogenized. Additional extrusion results are presented in Table 2.

As shown in Table 1, there were some quarter segments missing upon extrusion, that is, there was a gap where that quarter segment would have been expected. Where drainable liquid was obtained, it is believed the liquid originally occupied this gap and accounts for the missing quarter segment. Subsamples for laboratory analyses and archiving were created per the tank TCP.

TABLE 1: SUMMARY DESCRIPTION OF AUGER SAMPLES

Core and Segment* Number	Riser	Sample Total Weight (Grams)	Sample Collection General Description
Core 98 Segment 1(W)	12A	366.8	Extruded approx. 5-6 inches of sample. Sample was wet. Texture of sample was crystalline and tended to "melt" on the sample tray. Sample was light brown to dirty white in color. Collected 170 ml of drainable liquid. Color of liquid was Yellow and opaque.
Core 98 Segment 2 (A),(C), (D)	12A	416.0	Collected 100 ml of drainable liquid, which was yellow brown in color and opaque. Extruded approx. 12 inches of sample. Facies present; Sample was divided into quarter segments. Sample texture ranged from a brown sludge to a material resembling a mixture of brown, beige and dirty white saltcake.
Core 98 Segment 3 (A),(C), (D)	12A	406.5	Collected 60 ml of drainable liquid, which was light brown in color and opaque. Extruded approx. 16.5 inches of sample. Facies present; Sample was divided into quarter segments. Sample texture ranged from a brown sludge to a material resembling a mixture of brown, beige and dirty white saltcake.

WHC-SD-WM-DP-145, REV. 1
WHC-SD-WM-DP-145, REV. 0

Core 98 Segment 4 (A), (B), (C), (D)	12A	490.0	No drainable liquid observed or collected. Extruded approx. 19.0 inches of sample; Sample retained its shape. Facies present; Sample was divided into quarter segments. Sample texture ranged from a brown sludge to a material resembling a mixture of brown sludge with white saltcake.
Core 104 Segment 1(W)	7	131.3	Extruded approx. 6-7 inches of sample. Sample was tan in color, dry, crystalline and granular in shape. No drainable liquids.
Core 104 Segment 2 (A), (B), (C)	7	383.5	No drainable liquid. Extruded approx. 19.0 inches of sample; Upper half segment retained its shape, lower half segment crumbled during extrusion. Facies present in upper quarter segment (crumbly saltcake to sludge). Sample was divided into three quarter segments. Sample texture ranged from a brown sludge to a material resembling a mixture of brown sludge with saltcake.
Core 104 Segment 3 (A), (C), (D)	7	288.8	No drainable liquid. Extruded approx. 17.0 inches of sample; Sample was wet, granular and did not retain its shape. Sample was divided into three quarter segments. Sample texture resembled a brown saltcake.
Core 104 Segment 4 (A), (C), (D)	7	318.3	No drainable liquid. Extruded approx. 18.0 inches of sample. Sample was wet, granular, crumbly and ranged from a yellow to brown crystalline saltcake. Sample was divided into three quarter segments.
Core 104 Segment 5 (A), (B), (C), (D)	7	422.7	No drainable liquid. Extruded approx. 18 inches of sample. Lower half retained its shape and was medium to dark brown, whereas the upper half was dark brown and partially retained its shape. Texture of material resembled a sludge. Sample was divided into four quarter segments.

* (x) represents the quarter segment where "W" represents a whole segment and "A", "B", "C" and "D" represent the quarter segment location with "A" the top of the segment and "D" the bottom of the segment.

ANALYTICAL RESULTS

Analytical results are presented in Table 2.

Differential Scanning Calorimetry (DSC)

DSC analyses were performed under a nitrogen atmosphere using procedure LA-514-113, Rev. B-1 or procedure LA-514-114, Rev. B-0. The results are shown in Table 3 and the raw data scans are attached. The samples were analyzed in duplicate. Any exotherms on the scans would be visible as a rise (Mettler) or a sink (Perkin Elmer) from the baseline established at the beginning and ending of the scan.

Three of the samples exhibited exotherms above the action limit of 481 j/g, therefore, the upper 95% confidence level values calculated for each sample and the results are presented in Table 2.

Thermogravimetric Analysis (TGA)

Weight percent water is calculated from weight loss by TGA. These analyses were performed under a nitrogen atmosphere using procedure LA-560-112, Rev. A-2 or LA-514-114, Rev B-0. The samples and their related "immediate" samples from the un-homogenized extrusions were analyzed in duplicate. "Immediate" samples are samples for TGA analysis that were immediately taken directly from the un-homogenized sample as it was extruded onto the extrusion tray. TGA as well as DSC analyses were performed on homogenized samples for each sample. The results are presented in Table 2, and the raw data scans are attached.

Alpha Total

Analyses for total alpha activity were performed. Samples were prepared by fusion using procedure LA-549-141, Rev. D-0, and analyses were performed using procedure LA-508-101, Rev. D-2. Two fusions were prepared per sample (for duplicate results). Since all of the results were well below the safety screening limit of 41 $\mu\text{Ci}/\text{g}$, reruns to improve standard and spike recoveries or relative percent difference (RPD) were deemed unnecessary. The total alpha results are presented in Table 2.

Total Organic Carbon (TOC)

Three samples which had DSC results above the action limit were submitted for TOC analysis. Analyses were performed in duplicate according to procedure LA-342-100, Rev. C-0. The results are presented in Table 3. The RPD for sample S95T001977 exceeded the required $\pm 10\%$ range. A "triplicate" analysis was performed with a TOC result of 1.39 $\mu\text{g}/\text{g}$ which closely resembles the "duplicate" value.

Cyanide

Three samples which had DSC results above the action limit were submitted for cyanide analysis. Analyses were performed in duplicate according to procedure LA-695-102, Rev. D-0. The results are presented in Table 3.

Responsible Project Coordinator: J. H. Baldwin

REFERENCE

- [1] J. H. Baldwin, "Tank 241-BY-108 Tank Characterization Plan, WHC-SD-WM-TP-275, Rev. OD, Westinghouse Hanford Company, Richland, Washington, July 19, 1995.
- [2] H. Babad, J. W. Hunt, and K. S. Redus, *Tank Safety Screening Data Quality Objective*, WHC-SD-WM-SP-004, Rev. 1, Westinghouse Hanford Company, Richland, Washington, April 27, 1995.

WHC-SD-WM-DP-145, REV.1

WHC-SD-WM-DP-145, REV. 0

Table 4. 95% Upper Confidence Interval Limits
for DSC (Units are in joules/g).

Sample Description	Sample Number	Differential Scanning Calorimetry (DSC) Dry Basis				TGA* Average (%H2O)
		Sample Result (J/g)	Dup. Result (J/g)	Average Result (J/g)	95 % UL	
Quarter Segment B from Core 98 Segment 4	S95T001420	509.1	522.6	558.0	697.7	35.61
Rerun	S95T001420	481.3	409.5	445.4	672.0	35.61
Quarter Segment B from Core 104 Segment 5	S95T001976	593.5	522.6	558.0	781.8	36.66
Quarter Segment C from Core 104 Segment 5	S95T001977	548.7	507.9	528.3	657.1	35.50

TGA = Thermogravimetric Analysis
Quarter Segment B is approximately the second 5 inches of material of the segment. Quarter segment C is approximately the third 5 inches of material from the top of the segment.

WHC-SD-WM-DP-145, REV. 1

WHC-SD-WM-DP-145, REV. 0

SAMPLE DATA SUMMARY

2-11

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 98
SEGMENT #: 1(DL)

Table 2

SEGMENT PORTION: Drainable Liquid

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper									
S95T001373			% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	31.90	34.04	32.97	6.49	n/a	n/a	n/a
S95T001373			DSC Exotherm on Perkin Elmer	Joules/g	None	None	97.54	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001373			DSC Exotherm Dry Calculated	Joules/g Dry	None	None	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a

=> Limit violated
=> Selected Limit

2-12

WHC-SD-WM-DRP-145, REV.0

WHC-SD-WM-DRP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 1(W)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001389		% Water by TGA using Nettler	%	None	None	101.1	n/a	37.72	36.28	37.00	3.89	n/a	n/a	n/a

W Whole Segment: W Whole Segment

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001390		% Water by TGA on Perkin Elmer	%	None	None	96.42	n/a	21.85	19.76	20.80	10.0	n/a	n/a	n/a
S95T001390		DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	100.1	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001390		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001391	F	Alpha of Digested Solid	uCi/g	None	None	108.2	1.34e-03	2.22e-03	3.71e-03	2.97e-03	50.3	76.70	2.90e-03	90.8

=> Limit violated
=> Selected Limit

7-13

WHC-SD-WM-DP-145, REV. L
INTERIM

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Page: 3

CORE NUMBER: 98
SEGMENT #: 2(A)

Table 2

SEGMENT PORTION: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001404			% Water by TGA on Perkin Elmer	%	None	None	100.8	n/a	14.64	14.87	14.75	1.56	n/a	n/a	n/a
S95T001404			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480e+00	95.68	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001404			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480e+00	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001405	F		Alpha of Digested Solid	uCi/g	None	None	94.98	<2.26e-03	2.38e-03	1.86e-03	2.12e-03	24.5	81.30	3.75e-03	117.5

=> Limit violated
=> Selected Limit

2-T-4

WHC-SD-WM-DP-145, REV.0
WHC-SD-WM-DP-145, REV.L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 2(C)

SEGMENT PORTION: C Third Quarter of Segment

Sample#	RA#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
SP5T001399	X	Water by TGA using Mettler	%	None	None	100.7	n/a	43.90	39.74	41.82	9.95	n/a	n/a	n/a
SP5T001399		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	91.10	85.60	88.35	6.23	n/a	n/a	n/a
SP5T001399		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	93.15	n/a	53.00	49.80	51.40	6.23	n/a	n/a	n/a
SP5T001400	F	Alpha of Digested Solid	uCi/g	None	None	94.98	< 2.26e-03	< 2.63e-3	3.98e-03	n/a	n/a	82.80	3.92e-03	358.8

=> Limit violated

=> Selected Limit

ST-15

WHC-SD-WM-DP-145, REV. 0

WHC-SD-WM-DP-145, REV. 1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 2(D)

SEGMENT PORTION: D Bottom Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper										
S95T001393			% Water by TGA using Mettler	%	None	None	98.84	n/a	38.92	36.38	37.65	6.75	n/a	n/a	n/a	n/a
S95T001396			% Water by TGA on Perkin Elmer	%	None	None	100.8	n/a	8.560	9.360	8.960	8.93	n/a	n/a	n/a	n/a
S95T001396			% Water by TGA using Mettler	%	None	None	101.7	n/a	40.98	44.12	42.55	7.38	n/a	n/a	n/a	n/a
S95T001396			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	98.01	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a	n/a
S95T001396			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a	n/a
S95T001398	F		Alpha of Digested Solid	uCi/g	None	None	108.2	1.34e-03	2.67e-02	2.82e-02	2.75e-02	5.46	87.40	6.96e-03	39.6	

=> Limit violated

=> Selected Limit

The TGA - Water by Perkin Elmer represents weight loss at the first inflection point.

The TGA - Water by Mettler - represents total weight loss.

216

WHC-SD-WM-DP-145, REV.0

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 2(DL)

SEGMENT PORTION: Drainable Liquid

Sample#	R/A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001427		% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	27.11	25.03	26.07	7.98	n/a	n/a	n/a
S95T001427		DSC Exotherm on Perkin Elmer	Joules/g	None	None	98.45	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001427		DSC Exotherm Dry Calculated	Joules/g Dry	None	None	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a

=> Limit violated
=> Selected Limit

2-17

WHC-SD-WM-DP-145, REV. A

INTERIM

BY-108 Analytical Summary Cores 98 and 104

BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 3(A)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper	Standard %									
S95T001380			% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	15.74	15.92	15.83	1.01	n/a	n/a	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper	Standard %									
S95T001431			% Water by TGA using Mettler	%	None	None	101.3	n/a	43.77	44.47	44.12	1.59	n/a	n/a	n/a	n/a
S95T001431			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	323.2	323.0	323.1	0.06	n/a	n/a	n/a	n/a
S95T001431			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	106.5	n/a	180.6	180.5	180.6	0.06	n/a	n/a	n/a	n/a
S95T001480	F		Alpha of Digested Solid	uCi/g	None	None	103.2	<2.74e-03	7.63e-03	3.93e-03	5.78e-03	64.0	n/a	6.85e-03	85.5	

=> Limit violated
=> Selected Limit

2-18

WHC-SD-WM-DP-145, REV.1

INTERIM

8Y-108 Analytical Summary Cores 98 and 104
8Y-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 3(C)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001381		% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	15.18	16.82	16.00	0.10	n/a	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001432		% Water by TGA using Mettler	%	None	None	100.9	n/a	41.31	38.53	39.92	6.96	n/a	n/a	n/a
S95T001432		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	114.3	116.8	115.5	2.16	n/a	n/a	n/a
S95T001432		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	105.1	n/a	68.70	70.20	69.45	2.16	n/a	n/a	n/a
S95T001481	F	Alpha of Digested Solid	uCi/g	None	None	103.2	<2.74e-03	<4.66e-3	<5.93e-3	n/a	n/a	n/a	7.32e-03	500.0

=> Limit violated
=> Selected Limit

2-19

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 98
SEGMENT #: 3(D)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper								
S95T001382		% Water by TGA using Mettler	%	None	None	100.8	n/a	39.16	39.44	39.30	0.71	n/a	n/a

D Bottom Quarter of Segment: D Bottom Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
				Lower	Upper									
S95T001433		% Water by TGA using Mettler	%	None	None	101.8	n/a	34.52	38.56	36.54	11.1	n/a	n/a	
S95T001433		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	235.7	226.1	230.9	4.16	n/a	n/a	
S95T001433		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	106.9	n/a	149.6	143.5	146.6	4.16	n/a	n/a	
S95T001482	F	Alpha of Digested Solid	uCi/g	None	None	100.0	<2.24e-02	<2.06e-2	<1.73e-2	n/a	n/a	96.90	2.50e-02	245.9

=> Limit violated
=> Selected Limit

2-20

WHC-SD-WM-DP-145, REV. 1
WHC-SD-WM-DP-145, REV. 0

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 3(DL)

SEGMENT PORTION: Drainable Liquid

Sample#	R A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001430		% Water by TGA using Mettler	%	None	None	101.4	n/a	38.58	39.47	39.02	2.28	n/a	n/a	n/a
S95T001430		DSC Exotherm Dry Calculated	Joules/g Dry	None	None	n/a	n/a	121.7	114.2	118.0	6.36	n/a	n/a	n/a
S95T001430		DSC Exotherm using Mettler	Joules/g	None	None	95.96	n/a	74.20	69.60	71.90	6.40	n/a	n/a	n/a



=> Limit violated
=> Selected Limit

R-21

INTERIM

WHC-SD-WM-DP-145, REV. D

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 98
SEGMENT #: 4(A)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper									
S95T001407			% Water by TGA on Perkin Elmer	%	None	None	96.40	n/a	16.31	17.58	16.95	7.49	n/a	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper									
S95T001419			% Water by TGA using Mettler	%	None	None	97.31	n/a	28.88	36.23	32.55	22.6	n/a	n/a	n/a
S95T001419 1			% Water by TGA using Mettler	%	None	None	101.2	n/a	36.64	36.69	36.66	0.14	n/a	n/a	n/a
S95T001419			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480e-03	n/a	n/a	295.1	390.8	343.0	27.9	n/a	n/a	n/a
S95T001419 1			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480e-03	101.6	n/a	199.0	263.6	231.3	27.9	n/a	n/a	n/a
S95T001419			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480e-03	97.01	n/a	292.9	162.2	227.5	57.4	n/a	n/a	n/a
S95T001423	F		Alpha of Digested Solid	uCi/g	None	None	107.9	<5.29e-03	<1.45e-2	<7.36e-3	n/a	n/a	91.20	1.30e-02	500.0

=> Limit violated
=> Selected Limit

2-22

WHC-SD-WM-DP-145, REV. 0

WHC-SD-WM-DP-145, REV. 1

INTERIM

BY-108 Analytical Summary Cores 98 and 104

BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 4(B)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001408			% Water by TGA using Mettler	%	None	None	101.3	n/a	37.28	36.03	36.66	3.41	n/a	n/a

B Second Quarter of Segment: B Second Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001420			% Water by TGA using Mettler	%	None	None	97.31	n/a	35.41	35.81	35.61	1.12	n/a	n/a	
S95T001420 1			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	481.3	409.5	445.4	16.1	n/a	n/a	
S95T001420			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	509.1	438.1	473.6	15.0	n/a	n/a	
S95T001420 1			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	101.6	n/a	309.9	263.7	286.8	16.1	n/a	n/a	
S95T001420			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	97.01	n/a	327.8	282.1	305.0	15.0	n/a	n/a	
S95T001424	F		Alpha of Digested Solid	uCi/g	None	None	107.9	<5.29e-03	<7.47e-03	6.10e-02	n/a	n/a	92.20	1.06e-02	234.9

=> Limit violated
=> Selected Limit

2-23

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 4(C)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001409			% Water by TGA using Mettler	%	None	None	101.3	n/a	37.47	37.91	37.69	1.17	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001421			% Water by TGA using Mettler	%	None	None	101.5	n/a	38.31	39.82	39.06	3.87	n/a	n/a	
S95T001421			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	221.9	226.4	224.2	2.01	n/a	n/a	
S95T001421			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	102.3	n/a	135.4	138.2	136.8	2.05	n/a	n/a	
S95T001425	F		Alpha of Digested Solid	uCi/g	None	None	99.28	<4.22e-02	2.41e-01	3.13e-01	2.77e-01	26.0	111.2	9.66e-02	37.6

=> Limit violated
=> Selected Limit

2-24

WHC-SD-WM-DP-145, REV. D
WHC-SD-WM-DP-145, REV. 1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 98
SEGMENT #: 4(D)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001410		% Water by TGA on Perkin Elmer	%	None	None	100.6	n/a	19.95	20.93	20.44	4.79	n/a	n/a	n/a
S95T001410		% Water by TGA using Mettler	%	None	None	101.7	n/a	35.32	39.45	37.39	11.0	n/a	n/a	n/a

D Bottom Quarter of Segment: D Bottom Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001422		% Water by TGA using Mettler	%	None	None	101.5	n/a	36.49	36.40	36.45	0.25	n/a	n/a	n/a
S95T001422		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	680.0	n/a	n/a	98.51	107.6	103.1	8.82	n/a	n/a	n/a
S95T001422		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	102.3	n/a	62.60	68.40	65.50	8.85	n/a	n/a	n/a
S95T001426	F	Alpha of Digested Solid	uci/g	None	None	91.22	<9.64e-02	3.52e-01	3.68e-01	3.60e-01	4.44	97.80	2.44e-01	56.8

=> Limit violated

=> Selected Limit

The TGA - Water by Perkin-Elmer represents weight loss at the first inflection point.

The TGA - Water by Mettler - represents total weight loss.

7-25

WHC-SD-WM-DP-145, REV.1
WHC-SD-WM-DP-145, REV.2

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 1(W)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001586			% Water by TGA on Perkin Elmer	%	None	None	100.4	n/a	15.74	15.02	15.38	4.68	n/a	n/a

W Whole Segment: W Whole Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001925			% Water by TGA on Perkin Elmer	%	None	None	101.7	n/a	24.90	23.78	24.34	4.60	n/a	n/a	
S95T001925			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	97.79	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T001925			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T002001	F		Alpha of Digested Solid	uCi/g	None	None	125.4	<4.59e-04	<6.35E-04	<1.51E-3	n/a	n/a	4.78e-03	7.44e-04	1.96E+02

=> Limit violated
=> Selected Limit

2-26

WHC-SD-WM-DP-145, REV. D

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 2(A)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001587			% Water by TGA on Perkin Elmer	%	None	None	100.4	n/a	10.41	12.77	11.59	20.4	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001966			% Water by TGA using Mettler	%	None	None	101.3	n/a	20.84	19.25	20.05	7.93	n/a	n/a	
S95T001966			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T001966			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	106.9	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T002002	F		Alpha of Digested Solid	uCi/g	None	None	125.4	<4.59e-04	3.66e-03	3.18e-03	3.42e-03	14.0	7.48e-02	6.88e-04	2.78E+01

=> Limit violated
=> Selected Limit

2-27

WHC-SD-WM-DP-145, REV.0

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 2(B)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001588			% Water by TGA using Mettler	%	None	None	101.2	n/a	29.88	30.68	30.28	2.64	n/a	n/a	n/a

B Second Quarter of Segment: B Second Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001967			% Water by TGA using Mettler	%	None	None	101.2	n/a	33.54	29.90	31.72	11.5	n/a	n/a	n/a
S95T001967			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001967			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	106.9	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002003	F		Alpha of Digested Solid	uCi/g	None	None	107.4	<2.79e-03	4.74e-03	5.12e-03	4.93e-03	7.71	93.30	3.80e-03	1.30E+02

=> Limit violated
=> Selected Limit

2-28

WHC-SD-WM-DP-145, REV. L

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 104
SEGMENT #: 2(C)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001589			% Water by TGA using Mettler	%	None	None	101.2	n/a	22.95	23.67	23.31	3.09	n/a	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001968			% Water by TGA on Perkin Elmer	%	None	None	101.3	n/a	15.40	9.880	12.64	43.7	n/a	n/a	n/a
S95T001968			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	99.26	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001968			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002004	F		Alpha of Digested Solid	uCi/g	None	None	107.4	<2.79e-03	6.20e-03	6.17e-03	6.18e-03	0.49	91.74	3.83e-03	5.99E+01

=> Limit violated
=> Selected Limit

2-29

WHC-SD-WM-DP-145, REV. D

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 104
SEGMENT #: 3(A)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper								
S95T001861		% Water by TGA using Mettler	%	None	None	101.4	n/a	41.30	41.77	41.53	1.13	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
				Lower	Upper									
S95T001969		% Water by TGA on Perkin Elmer	%	None	None	101.2	n/a	7.150	7.790	7.470	8.57	n/a	n/a	
S95T001969		DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480E-03	99.16	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T001969		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480E-03	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T002005	F	Alpha of Digested Solid	uCi/g	None	None	104.7	<2.95E-03	<3.30E-03	1.83E-03	n/a	n/a	103.1	4.28E-03	1.87E+02

=> Limit violated
=> Selected Limit

2-30

WHC-SD-WM-DP-145, REV. O

WHC-SD-WM-DP-145, REV. 1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 3(C)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001862	X	Water by TGA using Mettler	%	None	None	101.4	n/a	26.93	25.93	26.43	3.78	n/a	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper	Standard %								
S95T001970	X	Water by TGA using Mettler	%	None	None	97.07	n/a	11.12	7.680	9.400	36.6	n/a	n/a	n/a
S95T001970		DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	99.79	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001970		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002006	F	Alpha of Digested Solid	uCi/g	None	None	110.1	6.54e-03	4.05e-03	<4.69E-3	n/a	n/a	91.74	9.66e-03	1.40E+02

=> Limit violated
=> Selected Limit

2-31

WHC-SD-WM-DP-145, REV. D

INTERIM

BY-108 Analytical Summary Cores 98 and 104

BY-108 (R)

Table 2

CORE NUMBER: 104
SEGMENT #: 3(D)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001863			% Water by TGA on Perkin Elmer	%	None	None	101.7	n/a	22.20	15.60	18.90	34.9	n/a	n/a
S95T001863			% Water by TGA using Mettler	%	None	None	n/a	n/a	32.06	n/a	n/a	n/a	n/a	n/a

D Bottom Quarter of Segment: D Bottom Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001971			% Water by TGA using Mettler	%	None	None	97.07	n/a	8.070	7.820	7.945	3.15	n/a	n/a	
S95T001971			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	99.79	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T001971			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	
S95T002007	F		Alpha of Digested Solid	uci/g	None	None	110.1	6.54e-03	1.86e-02	1.82e-02	1.84e-02	2.17	90.18	1.01e-02	5.90E+01

=> Limit violated
=> Selected Limit

The TGA - Water by Perkin-Elmer represents weight loss at the first inflection point.

The TGA - Water by Mettler - represents total weight loss.

2-32

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 4(A)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001904			% Water by TGA on Perkin Elmer	%	None	None	100.3	n/a	32.34	34.59	33.47	6.72	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001972			% Water by TGA using Mettler	%	None	None	100.7	n/a	31.01	33.09	32.05	6.49	n/a	n/a	
S95T001972			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	60.93	0.00e+00	30.46	200	n/a	n/a	
S95T001972			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	105.1	n/a	41.40	0.00e+00	20.70	200	n/a	n/a	
S95T002008	F		Alpha of Digested Solid	uCi/g	None	None	123.1	<1.91e-03	2.14e-03	<2.83E-3	n/a	n/a	78.57	3.85e-03	1.13E+02

=> Limit violated
=> Selected Limit

2-33

WHC-SD-WM-DP-145, REV. D

WHC-SD-WM-DP-145, REV. D

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 4(C)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper								
S95T001905			% Water by TGA using Mettler	%	None	None	100.8	n/a	43.26	45.28	44.27	4.56	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%	
					Lower	Upper									
S95T001973			% Water by TGA using Mettler	%	None	None	100.7	n/a	41.10	41.13	41.12	0.07	n/a	n/a	
S95T001973			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	180.2	194.5	187.3	7.63	n/a	n/a	
S95T001973			DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	105.1	n/a	106.1	114.5	110.3	7.62	n/a	n/a	
S95T002009	F		Alpha of Digested Solid	uCi/g	None	None	123.1	<1.91e-03	<4.66E-03	<3.50E-3	n/a	n/a	100.0	7.55e-03	5.00E+02

=> Limit violated
=> Selected Limit

2-34

WHC-SD-WM-DP-145, REV.0

WHC-SD-WM-DP-145, REV.1

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 4(D)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001906		% Water by TGA using Mettler	%	None	None	100.8	n/a	36.22	33.78	35.00	6.97	n/a	n/a	n/a

D Bottom Quarter of Segment: D Bottom Quarter of Segment

Sample#	R A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
				Lower	Upper									
S95T001974		% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	9.770	8.800	9.285	10.4	n/a	n/a	n/a
S95T001974		DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	99.86	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001974		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002010	F	Alpha of Digested Solid	uCi/g	None	None	115.4	5.50e-04	7.44e-04	<9.72E-4	n/a	n/a	91.96	9.16e-04	1.10E+02

=> Limit violated
=> Selected Limit

2-35

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)
Table 2

CORE NUMBER: 104
SEGMENT #: 5(A)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001915			% Water by TGA on Perkin Elmer	%	None	None	101.7	n/a	12.81	13.13	12.97	2.47	n/a	n/a	n/a

A Top Quarter of Segment: A Top Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001975			% Water by TGA on Perkin Elmer	%	None	None	101.1	n/a	9.370	8.980	9.175	4.25	n/a	n/a	n/a
S95T001975			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	99.86	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001975			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002011	F		Alpha of Digested Solid	uCi/g	None	None	115.4	5.50e-04	1.77e-02	1.08e-02	1.43e-02	48.4	69.20	9.11e-04	1.22E+01

=> Limit violated
=> Selected Limit

2-36

WHC-SD-WM-DP-145, REV. 0

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 2

CORE NUMBER: 104
SEGMENT #: 5(B)

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001916	1	% Water by TGA using Mettler	%	None	None	n/a	n/a	36.23	n/a	n/a	n/a	n/a	n/a	n/a	n/a
S95T001916		% Water by TGA using Mettler	%	None	None	100.7	n/a	35.23	37.56	36.39	15.3	n/a	n/a	n/a	n/a

B Second Quarter of Segment: B Second Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001976		% Water by TGA using Mettler	%	None	None	101.1	n/a	37.69	35.62	36.66	5.65	n/a	n/a	n/a	n/a
S95T001976		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	593.5	552.6	573.0	7.14	n/a	n/a	n/a	n/a
S95T001976		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	93.50	n/a	375.9	350.0	362.9	7.14	n/a	n/a	n/a	n/a
S95T002012	F	Alpha of Digested Solid	uCi/g	None	None	95.95	<7.63e-02	1.79e-01	1.81e-01	1.80e-01	1.11	100.0	1.03e-01	5.18E+01	

=> Limit violated
=> Selected Limit

6.37

WHC-SD-WM-DP-145, REV. D

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 5(C)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001917	1	% Water by TGA using Mettler	%		None	None	n/a	n/a		n/a	n/a	n/a	n/a	n/a	n/a
S95T001917		% Water by TGA using Mettler	%		None	None	100.7	n/a	33.04	26.36	29.70	22.5	n/a	n/a	n/a

C Third Quarter of Segment: C Third Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001977		% Water by TGA using Mettler	%		None	None	101.1	n/a	35.50	35.50	35.50	0.00	n/a	n/a	n/a
S95T001977		DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	480.0	n/a	n/a	548.7	507.9	528.3	7.72	n/a	n/a	n/a
S95T001977		DSC Exotherm using Mettler	Joules/g	-1.0e+03	480.0	93.50	n/a	n/a	353.9	327.6	340.8	7.72	n/a	n/a	n/a
S95T002013	F	Alpha of Digested Solid	uCi/g	None	None	95.95	<7.63e-02	2.66e-01	3.06e-01	2.86e-01	14.0	96.88	9.82e-02	3.55E+01	

=> Limit violated
=> Selected Limit

2-38

WHC-SD-WM-DP-145, REV. L

INTERIM

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 5(D)

Table 2

SEGMENT PORTION: Immediate Sampling (to check moisture loss)

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001918			% Water by TGA using Mettler	%	None	None	100.4	n/a	29.31	35.21	32.26	18.3	n/a	n/a	n/a

D Bottom Quarter of Segment: D Bottom Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits			Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper	Standard %								
S95T001978			% Water by TGA on Perkin Elmer	%	None	None	100.3	n/a	7.080	7.490	7.285	5.63	n/a	n/a	n/a
S95T001978			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+03	480.0	100.7	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T001978			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+03	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T002014	F		Alpha of Digested Solid	uCi/g	None	None	124.3	<9.89e-02	3.90e-01	3.74e-01	3.82e-01	4.19	138.8	1.82e-01	4.41E+01

=> Limit violated
=> Selected Limit

2-39

WHC-SD-WM-DP-145, REV. D

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 98
SEGMENT #: 4(B)TOC & CN

Table 3

SEGMENT PORTION: B Second Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001420			TOC by Persulfate/Coulometry	ug/g	-1.0e+03	3.0e+04	99.67	21.50	2.09e+04	2.07e+04	2.08e+04	0.96	n/a	40.00	n/a
S95T001420			Cyanide by Microdist. & Spec.	ug/g	-1.0e+03	3.9e+04	109.5	1.35e-01	2.01e+03	1.72e+03	1.86e+03	16.5	108.8	182.0	n/a

=> Limit violated
=====
=> Selected Limit

2-40

WHC-SD-WM-DP-145, REV. 0
145, REV. 1

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

Table 3

CORE NUMBER: 104

SEGMENT #: 5(B)TOC & CN

SEGMENT PORTION: B Second Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count	Err%
					Lower	Upper										
S95T001976			TOC by Persulfate/Coulometry	ug/g	-1.0e+03	3.0e+04	90.67	25.80	1.72e+04	1.34e+04	1.53e+04	24.8	86.80	40.00	n/a	
S95T001976			Cyanide by Microdist. & Spec.	ug/g	-1.0e+03	3.9e+04	107.3	1.29e-01	99.90	94.70	97.30	5.34	116.9	5.420	n/a	

=> Limit violated

=> Selected Limit

WHC-SD-WM-DP-145, REV. 0

WHC-SD-WM-DP-145, REV. 1

2-41

BY-108 Analytical Summary Cores 98 and 104
BY-108 (R)

CORE NUMBER: 104
SEGMENT #: 5(C)TOC & CN

Table 3

SEGMENT_PORTION: C Third Quarter of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T001977			TOC by Persulfate/Coulometry	ug/g	-1.0e+03	3.0e+04	90.67	25.80	1.36e+04	1.32e+04	1.34e+04	2.99	n/a	40.00	n/a
S95T001977			Cyanide by Microdist. & Spec.	ug/g	-1.0e+03	3.9e+04	109.2	1.08e-01	125.0	117.0	121.0	6.61	81.50	3.010	n/a



=> Limit violated
=> Selected Limit

WHC SD WM-DP-145, REV. 0
L-42

WHC-SD-WM-DP-145, REV.1

WHC-SD-WM-DP-145, REV. 0

INORGANIC ANALYSES

2-43

THIS PAGE INTENDED ONLY
FOR THE USE OF
THE LEFT BLANK


LABCORE Data Entry Template for Worklist#**2702**Analyst: BDV Instrument: DSC01 Book # Method: LA-514-113 Rev/Mod

Worklist Comment: Dry DSC for BY-108. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
95000118	BY-108 (R)	1 SAMPLE	S95T001976	0 DSC-02	SOLID	<u>N/A</u>	<u>593.5</u>	<u> </u>	Joules/g Dry
95000118	BY-108 (R)	2 DUP	S95T001976	0 DSC-02	SOLID	<u>593.5</u>	<u>552.6</u>	<u>N/A</u>	Joules/g Dry
95000118	BY-108 (R)	3 SAMPLE	S95T001977	0 DSC-02	SOLID	<u>N/A</u>	<u>548.7</u>	<u> </u>	Joules/g Dry
95000118	BY-108 (R)	4 DUP	S95T001977	0 DSC-02	SOLID	<u>548.7</u>	<u>507.9</u>	<u>N/A</u>	Joules/g Dry

Data entered + verified by

Final page for worklist #

2702


Analyst Signature Date 10-2-95

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-44

LABCORE Data Entry Template for Worklist#

2867

Analyst: BDV

Instrument: DSC01

Book # Method: LA-514-113 Rev/Mod

Worklist Comment: Calculated dry DSC for BY-108. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
95000118	BY-108 (R)	1 SAMPLE	S95T001925	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	2 DUP	S95T001925	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000104	BY-108 (R)	3 SAMPLE	S95T001431	0 DSC-02	SOLID	N/A	323.2		Joules/g Dry
95000104	BY-108 (R)	4 DUP	S95T001431	0 DSC-02	SOLID	323.2	323.0	N/A	Joules/g Dry
95000104	BY-108 (R)	5 SAMPLE	S95T001432	0 DSC-02	SOLID	N/A	114.3		Joules/g Dry
95000104	BY-108 (R)	6 DUP	S95T001432	0 DSC-02	SOLID	114.3	116.8	N/A	Joules/g Dry
95000104	BY-108 (R)	7 SAMPLE	S95T001433	0 DSC-02	SOLID	N/A	235.7		Joules/g Dry
95000104	BY-108 (R)	8 DUP	S95T001433	0 DSC-02	SOLID	235.7	226.1	N/A	Joules/g Dry
95000104	BY-108 (R)	9 SAMPLE	S95T001421	0 DSC-02	SOLID	N/A	221.9		Joules/g Dry
95000104	BY-108 (R)	10 DUP	S95T001421	0 DSC-02	SOLID	221.9	226.4	N/A	Joules/g Dry
95000104	BY-108 (R)	11 SAMPLE	S95T001422	0 DSC-02	SOLID	N/A	98.51		Joules/g Dry
95000104	BY-108 (R)	12 DUP	S95T001422	0 DSC-02	SOLID	98.51	107.6	N/A	Joules/g Dry
95000104	BY-108 (R)	13 SAMPLE	S95T001390	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000104	BY-108 (R)	14 DUP	S95T001390	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000104	BY-108 (R)	15 SAMPLE	S95T001396	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000104	BY-108 (R)	16 DUP	S95T001396	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000104	BY-108 (R)	17 SAMPLE	S95T001399	0 DSC-02	SOLID	N/A	91.10		Joules/g Dry
95000104	BY-108 (R)	18 DUP	S95T001399	0 DSC-02	SOLID	91.10	85.60	N/A	Joules/g Dry

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2867**

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
-------	---------	--------	---------	--------------------	--------	--------	-------	----	------

Final page for worklist # 2867

Blandina Valenzuela 10-12-95
 Analyst Signature Date

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
 R = Replicate Number, A = Aliquot Code.

2-46

LABCORE Data Entry Template for Worklist#

2868

Analyst: BDV Instrument: DSC01 Book # Method: LA-514-113 Rev/Mod

Worklist Comment: Calculated dry DSC for BY-108. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
95000104	BY-108 (R)	1 SAMPLE	S95T001404	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000104	BY-108 (R)	2 DUP	S95T001404	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	3 SAMPLE	S95T001966	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	4 DUP	S95T001966	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	5 SAMPLE	S95T001967	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	6 DUP	S95T001967	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	7 SAMPLE	S95T001968	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	8 DUP	S95T001968	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	9 SAMPLE	S95T001969	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	10 DUP	S95T001969	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	11 SAMPLE	S95T001970	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	12 DUP	S95T001970	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	13 SAMPLE	S95T001971	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	14 DUP	S95T001971	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	15 SAMPLE	S95T001972	0	DSC-02	SOLID	N/A	60.93		Joules/g Dry
95000118	BY-108 (R)	16 DUP	S95T001972	0	DSC-02	SOLID	60.93	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	17 TRIP	S95T001972	0	DSC-02	SOLID	60.93	59.60	N/A	Joules/g Dry
95000118	BY-108 (R)	18 SAMPLE	S95T001973	0	DSC-02	SOLID	N/A	180.2		Joules/g Dry

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2868

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
95000118	BY-108 (R)	19 DUP	S95T001973	0	DSC-02	SOLID	180.2	194.5	N/A Joules/g Dry

Final page for worklist # **2868**

Blandina Valenzuela 10-12-95
Analyst Signature Date

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
R = Replicate Number, A = Aliquot Code.

2 - 48

LABCORE Data Entry Template for Worklist#

2869

Analyst: BV Instrument: DSC01 Book # _____

Method: LA-514-113 Rev/Mod _____

Worklist Comment: Calculated dry DSC for BY-108. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
95000118	BY-108 (R)	1 SAMPLE	S95T001974	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	2 DUP	S95T001974	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	3 SAMPLE	S95T001975	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	4 DUP	S95T001975	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000118	BY-108 (R)	5 SAMPLE	S95T001978	0 DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000118	BY-108 (R)	6 DUP	S95T001978	0 DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry

Final page for worklist # 2869

Blandina Valenzuela 10-12-95
Analyst Signature Date

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2870**Analyst: BDV

Instrument: DSC01

Book # Method: LA-514-113 Rev/Mod

Worklist Comment: Calculated dry DSC for BY-108. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
95000101	BY-108 (R)	1 SAMPLE	S95T001365	0	DSC-02	LIQUID	N/A	Ø		Joules/g Dry
95000101	BY-108 (R)	2 DUP	S95T001365	0	DSC-02	LIQUID	Ø	Ø	N/A	Joules/g Dry
95000104	BY-108 (R)	3 SAMPLE	S95T001427	0	DSC-02	LIQUID	N/A	Ø		Joules/g Dry
95000104	BY-108 (R)	4 DUP	S95T001427	0	DSC-02	LIQUID	Ø	Ø	N/A	Joules/g Dry
95000104	BY-108 (R)	5 SAMPLE	S95T001430	0	DSC-02	LIQUID	N/A	121.7		Joules/g Dry
95000104	BY-108 (R)	6 DUP	S95T001430	0	DSC-02	LIQUID	121.7	114.2	N/A	Joules/g Dry
95000104	BY-108 (R)	7 SAMPLE	S95T001373	0	DSC-02	LIQUID	N/A	Ø		Joules/g Dry
95000104	BY-108 (R)	8 DUP	S95T001373	0	DSC-02	LIQUID	Ø	Ø	N/A	Joules/g Dry

Final page for worklist # **2870**Blandina Valenquele 10-12-95

Analyst Signature

Date

Analyst Signature

Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2001

Analyst: JASInstrument: DSC0 1
10-16-95 BDVBook # 12N/YAMethod: LA-514-113 Rev/Mod B-C-B-1

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID	<u>28.45</u>	<u>27.6</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001419 0		DSC-01	SOLID	<u>N/A</u>	<u>292.9</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001419 0		DSC-01	SOLID	<u>292.9</u>	<u>162.2</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001420 0		DSC-01	SOLID	<u>N/A</u>	<u>327.8</u>		Joules/g
95000104	BY-108 (R)	5 DUP	S95T001420 0		DSC-01	SOLID	<u>327.8</u>	<u>282.1</u>	<u>N/A</u>	Joules/g

Final page for worklist # 2001

Jas 8-16-95
 Analyst Signature Date

LJ 8-22-95
 Analyst Signature Date

S95T001419 produced an endothermic region at 125.3°C with a delta H of 801.7 J/g. Because of high RPDs the samples will be rerun.

Data Entry Comments: S95T001420 produced an endothermic region at 130.4°C with a delta H of 822.1 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 252 TO 256.

BEST AVAILABLE COPY


DSC STD 12N14A

6.560 mg

Rate: 10.0 °C/min

File: 00017.001 DSC METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001419 SAM N2

32.400 mg

Rate: 10.0 °C/min

File: 00019.001 DSC METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-1/92, REV. 1

BEST AVAILABLE COPY


S95T001419 DUP N2

33.600 mg

Rate: 10.0 °C/min

File: 00022.001 DSC METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-142, REV. L


BEST AVAILABLE COPY

S95T001420 SAM N2

35.740 mg

Rate: 10.0 °C/min

File: 00024.001 DSC METTLER 16-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. L

BEST AVAILABLE COPY


S95T001420 DUP N2

26.820 mg

Rate: 10.0 °C/min

File: 00026.001 DSC METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-145, REV. L

LABCORE Data Entry Template for Worklist#

2025

Analyst: RJ McCown Instrument: DSC0 1 Book # 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>28.9</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001419 1	DSC-01	SOLID	<u>N/A</u>	<u>199.0</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001419 1	DSC-01	SOLID	<u>199.0</u>	<u>263.6</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001420 1	DSC-01	SOLID	<u>N/A</u>	<u>309.9</u>		Joules/g
95000104	BY-108 (R)	5 DUP	S95T001420 1	DSC-01	SOLID	<u>309.9</u>	<u>263.7</u>	<u>N/A</u>	Joules/g

Final page for worklist # 2025

RJ McCown 9/6/95
Analyst Signature DateJ. M. Lye 9/7/95
Analyst Signature Date

Verified 9/7/95 J. M. Lye

Data Entry Comments:

S95T001419 sample has an exotherm of 853.5 J/g at 129°C; duplicate has 770.5 J/g at 126°C. JNET
S95T001420 sample has an endotherm of 712.2 J/g at 122°; duplicate has 823.9 J/g at 119°C. JNET

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

DSC STD 12N14A

6.580 mg

Rate: 10.0 °C/min


File: 00065.001

Ident: 0.0

DSC METTLER 06-Sep-95

222-S Laboratory

RJMcGown
9/16/95



BEST AVAILABLE COPY

S95T001419 SAM N2

34.981 mg


Rate: 10.0 °C/min

File: 00067.001 DSC METTLER 06-Sep-95

Ident: 0.0

222-S Laboratory

R. McLellan
9/6/95



S95T001419 DUP N2

25.079 mg

Rate: 10.0 °C/min

File: 00069.001 DSC METTLER 06-Sep-95

Ident: 0.0

222-S Laboratory

RJMcClown
9/6/95

exo>

09 -e

20 . mW

Integration
Delta H 19324 mJ
770.5 J/g
Peak 126.7 °C
-45.0 mW

Integration
Delta H 6612 mJ
263.6 J/g
Peak 287.4 °C
7.6 mW

100. 200. 300. 400. °C

S95T001420 SAM N2


34.051 mg


Rate: 10.0 °C/min

File: 00071.001 DSC METTLER 06-Sep-95

Ident: 0.0 222-S Laboratory

RJMcClown
9/6/95





LABCORE Data Entry Template for Worklist#

2026

Analyst: JDS Instrument: DSC0 3 Book # 12N14-AMethod: LA-514-114 Rev/Mod B-D

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-03	SOLID	<u>28.49</u> <small>805 8-24-95</small>	<u>28.49</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001390 0	DSC-03	SOLID	<u>N/A</u>	<u>X</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001390 0	DSC-03	SOLID	<u>28.45</u>	<u>X</u>	<u>N/A</u>	Joules/g
		4 STD		DSC-03	SOLID	<u>28.16</u> <small>8-24-95 BDV</small>	<u>28.11</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	5 SAMPLE	S95T001396 0	DSC-03	SOLID	<u>N/A</u>	<u>X</u>		Joules/g
95000104	BY-108 (R)	6 DUP	S95T001396 0	DSC-03	SOLID	<u>X</u>	<u>X</u>	<u>N/A</u>	Joules/g

Final page for worklist # 2026

See attached for signatures

Analyst Signature Date 8-24-95

Analyst Signature Date

Verified by Blandina Valenzuela 8-28-95

S95T001390 produced three endothermic regions at 119.6 °C with a delta H of 610.5 J/g; second at 187.4 °C with a delta H of 89.2 J/g and third at 276.1 °C with a delta H of 134.6 J/g.

Data Entry Comments:


S95T001396 produced one large endothermic exothermic at 146.66 °C with a delta H of 1049.0 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2026**Analyst: JDS Instrument: DSC0 _____ Book # 12NNNMethod: LA-514-113 Rev/Mod B-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID		N/A	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001390 0		DSC-01	SOLID	N/A		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001390 0		DSC-01	SOLID		N/A	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001396 0		DSC-01	SOLID	N/A		Joules/g
95000104	BY-108 (R)	5 DUP	S95T001396 0		DSC-01	SOLID		N/A	Joules/g

Final page for worklist # **2026**8-23-95

Analyst Signature

Date



Date

Other instrument
was used

8-24-95BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

5-64

Curve 1: DSC


File info: IND081001 Fri Aug 18 09:28:07 1995

Sample Weight: 6.480 mg

12N14A Indium at 10C/min

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES TO

BEST AVAILABLE COPY



N₂, EXOTHERM DOWN

TEMP1: 150.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 170.0 °C

JD SPELLMAN

PERKIN ELMER

222-S Lab

Fri Aug 18 09:29:31 1995


Curve 1: DSC

File info: SAM081802 Fri Aug 18 12:08:21 1995

Sample Weight: 23.670 mg

S95T001390 SAM N2

BEST AVAILABLE COPY



exotherm down, N2 purge gas
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 500.0 °C

Temperature (°C)

JD SPELLMAN
PERKIN ELMER
222-S Lab

Fri Aug 18 12:52:03 1995

WHC-SD-WM-DR. 145 REV. L

Curve 1: DSC


File info: SAM081803 Fri Aug 18 13:58:49 1995

Sample Weight: 44.260 mg

S95T001390 DUP

BEST AVAILABLE COPY

49-c



WHC-SD-WM-DP-45, REV.1

exotherm down, N2 purge gas
TEMP1: 25.0 C TIME1: 0.0 min RATE1: 10.0 C/min

Temperature (°C)

JD SPELLMAN
PERKIN ELMER
222-S Lab
Fri Aug 18 14:02:07 1995


Curve 1: DSC

File info: IND082301 Wed Aug 23 09:07:10 1995

Sample Weight: 6.480 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY



N2, EXOTHERM DOWN

TEMP1: 150.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 170.0 °C

JD SPELLMAN
PERKIN ELMER
222-S Lab
Wed Aug 23 10:56:59 1995


Curve 1: DSC

File info: SAM082301 Wed Aug 23 12:41:39 1995

Sample Weight: 39.800 mg

S95T001396 SAM

BEST AVAILABLE COPY



exotherm down, N2 purge gas


TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 500.0 °C

JD SPELLMAN
PERKIN ELMER
222-S Lab
Wed Aug 23 13:43:14 1995

Curve 1: DSC
File info: SAM082302 Wed Aug 23 14:35:05 1995
Sample Weight: 17.900 mg
S95T001396 SAM

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

JD SPELLMAN
PERKIN ELMER
222-S Lab
Wed Aug 23 14:48:15 1995

WHC-SD-WM-DP-145, REV. C

LABCORE Data Entry Template for Worklist#

2027

Analyst: JDSInstrument: DSC0 1Book # I2N14AMethod: LA-514-113 Rev/Mod B-1LA-514-114 Rev/Mod B-0 bdv

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>26.5</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001399 0	DSC-01	SOLID	<u>N/A</u>	<u>53.0</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001399 0	DSC-01	SOLID	<u>53.0</u>	<u>49.8</u>	<u>N/A</u>	Joules/g
		4 STD		DSC-03	SOLID	<u>28.45</u>	<u>27.22</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	5 SAMPLE	S95T001404 0	DSC-03	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000104	BY-108 (R)	6 DUP	S95T001404 0	DSC-03	SOLID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u>	Joules/g

Final page for worklist #

2027

See attached for signatures
Analyst Signature Date 8-31-95R. Jones 8-31-95
Analyst Signature DateVerified by Blandina Valenzuela
8-31-95

S95T001399 produced one endothermic region at 133.5°C with a delta H of 1265.7 J/g.

Data Entry Comments: S95T001404 produced ^{two} ~~one~~ endothermic region one at 134.36°C
^{8-31-95 bdv}
 with a delta H of 823.3 J/g and second at 279.9°C with a delta H of 14.5 J/g.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
 R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2027**Analyst: Jds Instrument: DSC0 Book # 12N14A

Method: LA-514-113 Rev/Mod

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID			N/A	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001399 0	DSC-01	SOLID	N/A			Joules/g
95000104	BY-108 (R)	3 DUP	S95T001399 0	DSC-01	SOLID			N/A	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001404 0	DSC-01	SOLID	N/A			Joules/g
95000104	BY-108 (R)	5 DUP	S95T001404 0	DSC-01	SOLID			N/A	Joules/g

Final page for worklist # 2027Jds

8.29.95

Analyst Signature

Date

Analyst Signature

Date

Other instrument was
used.

8-31-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 273 TO 278.

BEST AVAILABLE COPY

DSC STD 12N14A

6.560 mg

Rate: 10.0 °C/min


File: 00028.001

Ident: 0.0

DSC METTLER

18-Aug-95

222-S Laboratory



WHC-SD-WM-DP-145, REV 1

BEST AVAILABLE COPY


S95T001399 SAM N2

23.800 mg

Rate: 10.0 °C/min

File: 00030.001 DSC METTLER 18-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR. 145, REV. 1

BEST AVAILABLE COPY

S95T001399 DUP N2


30.900 mg

Rate: 10.0 °C/min

File: 00032.001 DSC METTLER 18-Aug-95

Ident: 0.0

222-S Laboratory



WHC-SD-WM-DR. 145, REV. 1


Curve 1: DSC

File info: IND082801 Mon Aug 28 11:14:02 1995

Sample Weight: 6.480 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

N2, EXOTHERM DOWN

TEMP1: 190.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 170.0 °C

Temperature (°C)

JD SPELLMAN


PERKIN ELMER

222-S Lab

Mon Aug 28 11:14:53 1995


Curve 1: DSC
File info: qsav2 Mon Aug 28 14:04:18 1995
Sample Weight: 42.390 mg
S95T001404 SAM

BEST AVAILABLE COPY



Curve 1: DSC
File info: SAM082802 Mon Aug 28 15:14:22 1995
Sample Weight: 34.400 mg
S95T001404 DUP

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN ELMER
222-S Lab
Mon Aug 28 15:18:02 1995

WHC-SD-WM-DR-145, REV. L

LABCORE Data Entry Template for Worklist# 2032Analyst: JDS Instrument: DSC0 3 Book # 12N14-AMethod: LA-514-114 Rev/Mod B-0Worklist Comment: Please run BY-108 ^{DSC} TGAs under N2. bdv
8-31-95 BDV

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-03	LIQUID	<u>28.45</u>	<u>27.75</u>	<u>N/A</u> Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001373 0		DSC-03	LIQUID	<u>N/A</u>	<u>Ø</u>	<u>N/A</u> Joules/g
95000104	BY-108 (R)	3 DUP	S95T001373 0		DSC-03	LIQUID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u> Joules/g

Final page for worklist # 2032See attached for signaturesAnalyst Signature Date 8-31-95L. Jones8-31-95

Analyst Signature

Date

Verified by Blandina Valenzuela
BDV
8-31-95Data Entry Comments: Sample produced one endothermic region at 140.9 °C
with a delta H of 1376.1 J/gUnits shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2032**Analyst: Jds Instrument: DSC0 Book # 12 N 14AMethod: LA-514-113 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGA's under N2. bdv

DSC ~~SMT~~ 8/19/95

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	LIQUID		N/A	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001373 0		DSC-01	LIQUID	N/A		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001373 0		DSC-01	LIQUID		N/A	Joules/g

Final page for worklist # 2032

Jdh Smt 8-29-95 _____ Analyst Signature Date Analyst Signature Date

Other instrument was
used

8-31-95BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-80

Curve 1: DSC


File info: ind082901 Tue Aug 29 08:04:09 1995

Sample Weight: 6.480 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-81 TO 2-83.



N2, EXOTHERM DOWN

TEMP1: 150.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 170.0 °C

Temperature (°C)

JD SPELLMAN

PERKIN ELMER

222-S Lab

Tue Aug 29 14:25:32 1995

WHC-SD-WM-DP-145, REV. 1


Curve 1: DSC

File info: SAM082903 Tue Aug 29 14:00:56 1995

Sample Weight: 39.240 mg

S95T001373 SAM

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. C

exotherm down, N₂ purge gas
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

JD SPELLMAN
PERKIN ELMER
222-S Lab
Tue Aug 29 14:07:09 1995


Curve 1: DSC

File info: SAM082904 Tue Aug 29 15:10:34 1995

Sample Weight: 32.460 mg

S95T001373 SAM

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP1: 30.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 300.0 °C

Temperature (°C)

JD SPELLMAN
PERKIN ELMER
222-S Lab
Tue Aug 29 15:14:56 1995

WHC-SD-WM-DP-145, REV. 1

LABCORE Data Entry Template for Worklist#

2038

Analyst: JDS Instrument: DSC0 3 Book # 12N14A

Method: LA-514-114 Rev/Mod B-D

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-03	LIQUID	<u>28.45</u>	<u>28.01</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001427 0	DSC-03	LIQUID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001427 0	DSC-03	LIQUID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u>	Joules/g

Final page for worklist # **2038**

See attached for signatures

Analyst Signature Date 8-31-95

R. Jones 8-31-95
Analyst Signature Date

Verified by Blandina Valenzuela
BDV
8-31-95

Data Entry Comments: Sample produced one endothermic region at 123.8°C with a delta H of 17.70.2 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2038

Analyst: Jds Instrument: DSC0 Book # 12N14A

Method: LA-514-113 Rev/Mod B-0

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	LIQUID			N/A	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001427 0	DSC-01	LIQUID	N/A			Joules/g
95000104	BY-108 (R)	3 DUP	S95T001427 0	DSC-01	LIQUID			N/A	Joules/g

Final page for worklist # **2038**

Jds 8-29-95
Analyst Signature Date

Analyst Signature Date

Other instrument was
used.

8-31-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: DSC


File info: ind082901 Tue Aug 29 08:04:09 1995

Sample Weight: 6.480 mg

12Ni4A Indium at 10C/min

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 286 TO 288.

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

N2, EXOTHERM DOWN

TEMP1 150.0 C TIME1: 0.0 min RATE1: 10.0 C/min

TEMP2 170.0 C

Temperature (°C)

JO SPELLMAN
PERKIN ELMER
222-S Lab
Tue Aug 29 08:12:14 1995


Curve 1: DSC

File info: SAM082901 Tue Aug 29 09:09:07 1995

Sample Weight: 18.560 mg

S95T001427 SAM

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. A

exotherm down, N₂ purge gas

TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 200.0 °C

JD SPELLMAN
PERKIN ELMER
222-S Lab
Tue Aug 29 09:13:54 1995

Curve 1: DSC


File info: SAM082902 Tue Aug 29 10:06:53 1995

Sample Weight: 24.560 mg

S95T001427 SAM DSC
8/24/95

BEST AVAILABLE COPY

J-88



WHC-SD-VIM-DP-145, REV. 1

exotherm down, N₂ purge gas
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PEAKIN ELMER
222-S Lab
Tue Aug 29 10:59:20 1995

LABCORE Data Entry Template for Worklist#

2112

Analyst: SMF Instrument: DSC0 1 Book # 12N14AMethod: LA-514-113 Rev/Mod B-1

Worklist Comment: Please run BY-108 samples under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID	<u>28.45</u>	<u>29.1</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001421	0	DSC-01	SOLID	<u>N/A</u>	<u>80.2</u>	<u>135.4</u>	Joules/g
95000104	BY-108 (R)	3 DUP	S95T001421	0	DSC-01	SOLID	<u>135.4</u>	<u>138.2</u>	<u>8-30-95 BDV</u>	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001422	0	DSC-01	SOLID	<u>80.2</u>	<u>102.4</u>	<u>8-30-95 N/A</u>	Joules/g
95000104	BY-108 (R)	5 DUP	S95T001422	0	DSC-01	SOLID	<u>62.6</u>	<u>62.6</u>	<u>8-30-95 BDV</u>	Joules/g
						SOLID	<u>62.6</u>	<u>68.4</u>	<u>N/A</u>	Joules/g

Final page for worklist # 2112Susie M. Fulton 8-29-95

Analyst Signature Date

Dany Hammel 8-30-95

Analyst Signature Date

Verified by Blandina Valenzuela

8-31-95

S95T001421 sample produced one large endotherm at 119.3°C with a delta H of 658.1 J/g. The duplicate exotherm result is the sum of two exotherms: 589.6 J/g.

Data Entry Comments:

S95T001422 sample produced one large endotherm at 129.6°C with a delta 798.9 J/g.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 290 TO 294.

BEST AVAILABLE COPY


DSC STD 12N14A

6.560 mg

Rate: 10.0 °C/min

File: 00092.001 DSC METTLER 29-Aug-95

Ident: 0.0 222-S Laboratory




BEST AVAILABLE COPY

S95T001421 N2
38.320 mg

Rate: 10.0 °C/min

File: 00094.001 DSC METTLER 29-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-245, REV. L

BEST AVAILABLE COPY


S95T001421 DUP N2

35.477 mg

Rate: 10.0 °C/min

File: 00096.001 DSC METTLER 29-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001422 N2

25.677 mg

Rate: 10.0 °C/min

File: 00098.001 DSC METTLER 29-Aug-95

Ident: 0.0 222-S Laboratory



26-93

WHC-SD-WM-DP-1425, REV. L

BEST AVAILABLE COPY


S95T001422 DUP N2

33.529 mg

Rate: 10.0 °C/min

File: 00100.001 DSC METTLER 29-Aug-95

Ident: 0.0 222-S Laboratory



WDC-SD-WMDP-145, REV. L

LABCORE Data Entry Template for Worklist#**2114**Analyst: JKS Instrument: DSC0 1 Book # 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S	TYPE	SAMPLE#	R	A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1	STD				DSC-01	SOLID			N/A	Joules/g
95000104	BY-108 (R)	2	SAMPLE	S95T001431	0		DSC-01	SOLID	N/A			Joules/g
95000104	BY-108 (R)	3	DUP	S95T001431	0		DSC-01	SOLID			N/A	Joules/g
95000104	BY-108 (R)	4	SAMPLE	S95T001432	0		DSC-01	SOLID	N/A			Joules/g
95000104	BY-108 (R)	5	DUP	S95T001432	0		DSC-01	SOLID			N/A	Joules/g

Final page for worklist # **2114**JKS9-5-95

Analyst Signature

Date

Analyst Signature

Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-95

LABCORE Data Entry Template for Worklist#**2114**Analyst: JDS Instrument: DSC0 1 Book #: 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>30.3</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001431 0	DSC-01	SOLID	<u>N/A</u>	<u>180.6</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001431 0	DSC-01	SOLID	<u>180.6</u>	<u>180.5</u>	<u>N/A</u>	Joules/g
		4 STD		DSC-01	SOLID	<u>28.45</u>	<u>29.9</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	5 SAMPLE	S95T001432 0	DSC-01	SOLID	<u>N/A</u>	<u>68.7</u>		Joules/g
95000104	BY-108 (R)	6 DUP	S95T001432 0	DSC-01	SOLID	<u>68.7</u>	<u>70.2</u>	<u>N/A</u>	Joules/g

Final page for worklist # **2114**

See Attached for Signatures
 Analyst Signature Date
 Entered Verified 9/15/95 John Fye

Analyst Signature	Date
-------------------	------

S95T001431 produced a large endothermic region at 123.3°C
 with a delta H of 976.5 J/g

Data Entry Comments: S95T001432 produced a large endothermic region at
 119.3°C with a delta H of 943.3 J/g.

3inf 9/16/95

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 291 TO 292.

BEST AVAILABLE COPY

DSC STD 12N14A

6.560 mg

Rate: 10.0 °C/min


File: 00001.001

DSC METTLER

31-Aug-95

Ident: 0.0

222-S Laboratory



WHC-SD-WM-DP. 145, REV. 1

BEST AVAILABLE COPY


S95T001431 SAM N2

36.800 mg

Rate: 10.0 °C/min

File: 00003.001 DSC METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001431 DUP N2

18.650 mg

Rate: 10.0 °C/min

File: 00005.001 DSC METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-11-M-DP-144, REV. 1

BEST AVAILABLE COPY

DSC STD 12N14A


6.580 mg

Rate: 10.0 °C/min

File: 00040.001 DSC METTLER 04-Sep-95

Ident: 0.0

222-S Laboratory




WHC-SD-WM-DP-145, REV. L

Th. S. 9-4-95

BEST AVAILABLE COPY

WHC-SD-W/M-DP-145, REV. 1

S95T001432 SAM N2
42.150 mg Rate: 10.0 °C/min



BEST AVAILABLE COPY

S95T001432 DUP N2


26.700 mg

Rate: 10.0 °C/min

File: 00044.001 DSC METTLER 04-Sep-95

Ident: 0.0

222-S Laboratory



WHC-SD-WN-DR-145, REV. L

LABCORE Data Entry Template for Worklist#

2115

Analyst: SJS

Instrument: DSC0 1

Book # 12N14A

Method: LA-514-113 Rev/Mod B-1

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>30.4</u>	<u>N/A</u>	Joules/g
95000104	BY-108 (R)	2 SAMPLE	S95T001433 0	DSC-01	SOLID	<u>N/A</u>	<u>149.6</u>		Joules/g
95000104	BY-108 (R)	3 DUP	S95T001433 0	DSC-01	SOLID	<u>149.6</u>	<u>143.5</u>	<u>N/A</u>	Joules/g

Final page for worklist #

2115

TJL SJS
Analyst Signature

8-73 8-30-95
Date

Dawn Hammett
Analyst Signature

8/30/95
Date

Verified by Blandina Valenzuela
(8-31-95)

Data Entry Comments: Sample produced one endothermic region at 131.1°C
with a delta H of 972.5°C. J/g
8/31/95
BV

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-104 TO 2-106.

BEST AVAILABLE COPY


DSC STD 12N14A

6.560 mg

Rate: 10.0 °C/min

File: 00102.001 DSC METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. L

BEST AVAILABLE COPY


S95T001433 SAM N2

17.260 mg

Rate: 10.0 °C/min

File: 00104.001 DSC METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001433 DUP N2

17.810 mg

Rate: 10.0 °C/min

File: 00106.001 DSC METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

LABCORE Data Entry Template for Worklist#

2142Analyst: JDS Instrument: DSC0 1 Book #: 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	LIQUID	<u>28.45</u>	<u>28.1</u>	<u>9.12-95</u>	
							<u>27.3</u>	<u>27.3</u>	<u>BDV</u>	<u>N/A</u>
95000101	BY-108 (R)	2 SAMPLE	S95T001365	0	DSC-01	LIQUID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000101	BY-108 (R)	3 DUP	S95T001365	0	DSC-01	LIQUID	<u>Ø</u>	<u>Ø</u>		N/A Joules/g
		4 STD			DSC-01	LIQUID	<u>28.45</u>	<u>27.3</u>		<u>N/A</u> Joules/g
95000104	BY-108 (R)	5 SAMPLE	S95T001430	0	DSC-01	LIQUID	<u>N/A</u>	<u>74.2</u>		Joules/g
95000104	BY-108 (R)	6 DUP	S95T001430	0	DSC-01	LIQUID	<u>74.2</u>	<u>69.6</u>		<u>N/A</u> Joules/g

Final page for worklist # **2142**See attached for signatures
Analyst Signature Date 9-12-95
BDV
Analyst Signature 9-13-95
DateVerified 9/13/95. J.M. Lye

S95T001365 produced one endotherm at 99.3°C with a delta H of 1295.15 J/g.

Data Entry Comments: S95T001430 produced one endotherm at 127.3°C with a delta H of 1366.0 J/g.

LABCORE Data Entry Template for Worklist#

2142

Analyst: Sds Instrument: DSC0 C-O Book #: 12N14A

Method: LA-514-113 Rev/Mod

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	LIQUID			N/A	Joules/g
95000101	BY-108 (R)	2 SAMPLE	S95T001365 0	DSC-01	LIQUID	N/A			Joules/g
95000101	BY-108 (R)	3 DUP	S95T001365 0	DSC-01	LIQUID			N/A	Joules/g
95000104	BY-108 (R)	4 SAMPLE	S95T001430 0	DSC-01	LIQUID	N/A			Joules/g
95000104	BY-108 (R)	5 DUP	S95T001430 0	DSC-01	LIQUID			N/A	Joules/g

Final page for worklist #

2142

Jah Sds

7-8-95

Analyst Signature

Date

Analyst Signature

Date

Data Entry Comments:

Sample S95T001365 has not been released
 from the Hot Cells yet. 7-8-95

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

J 108

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2109 TO 2114.

BEST AVAILABLE COPY


DSC STD 12N14A

6.630 mg

Rate: 10.0 °C/min

File: 00001.001 DSC METTLER 11-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001365 N2

18.420 mg

Rate: 10.0 °C/min

File: 00003.001 DSC METTLER 11-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001365 DUP N2

30.250 mg

Rate: 10.0 °C/min

File: 00004.001 DSC METTLER 11-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

DSC STD 12N14A

6.580 mg

Rate: 10.0 °C/min

File: 00086.001 DSC METTLER 07-Sep-95

Ident: 0.0 222-S Laboratory

exo >

5. mW

2-112

120. 140. 160. 180. °C

Integration
Delta H 180 mJ
27.3 J/g
Peak 159.0 °C
-12.1 mW

JH Spel 9-7-95

WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001430 SAM N2

21.500 mg

Rate: 10.0 °C/min

File: 00088.001 DSC METTLER 07-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001430 DUP N2

32.800 mg

Rate: 10.0 °C/min

File: 00090.001 DSC METTLER 07-Sep-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. L

LABCORE Data Entry Template for Worklist#**2273**Analyst: SMF Instrument: DSC01/3 Book # 12N14AMethod: LA-514-114 Rev/Mod C-O 9-11-95 BDV

Worklist Comment: BY-108 DSC, please run under N2. JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-03	SOLID	<u>28.45</u>	<u>27.82</u>	<u>N/A</u> Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001925 0		DSC-03	SOLID	<u>N/A</u>	<u>Ø</u>	<u>Joules/g</u>
95000118	BY-108 (R)	3 DUP	S95T001925 0		DSC-03	SOLID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u> Joules/g

Final page for worklist # 2273

See attached for signatures
 Analyst Signature Date 9-11-95 BDV

J. Jan 9-14-95
 Analyst Signature Date

Verified 9/14/95. Jan M. Faye

Data Entry Comments: Sample produced three endothermic regions; one at 122.2°C with a delta H of 648.4 J/g, second at 284.2°C with a delta H of 7.2 J/g, and third at 319.7°C with a delta H of 65.2 J/g.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2273**Analyst: SMF Instrument: DSC01 Book # 12 N14AMethod: LA-514-113 Rev/Mod LA-514-114 1C-C
SUF 9-8-95

Worklist Comment: BY-108 DSC, please run under N2. JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID		N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001925	0	DSC-01	SOLID	N/A		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001925	0	DSC-01	SOLID		N/A	Joules/g

Final page for worklist # 2273Smulton 9-8-95

Analyst Signature Date

Analyst Signature Date

Other instruments were used
 9-11-95
 BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-116

Curve 1: DSC


File info: IND090801 Fri Sep 8 00:23:54 1995

Sample Weight: 6.500 mg

12N14A Indium at 10C/min

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 217 TO 219.

BEST AVAILABLE COPY



N2, EXOTHERM DOWN

TEMP: 150.0 °C TIMES: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON

PERKIN ELMER

222-S Lab

Fri Sep 8 03:28:32 1995

Curve 1: DSC


File info: SAM090802 Fri Sep 8 04:03:19 1995

Sample Weight: 21.690 mg

S95T001925 at 10C/min

BEST AVAILABLE COPY

2118



WHC-SD-WM-DR-145, REV. 1

exotherm down, N2 purge gas
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN ELMER
222-S Lab
Fri Sep 8 05:25:06 1995

Curve 1: DSC


File info: SAM0908003 Fri Sep 8 06:19:42 1995

Sample Weight: 14.760 mg

S95T001925 DUP at 10C/min

BEST AVAILABLE COPY

LTT-2-119



WHC-SD-WM-DR-145, REV. L

exotherm down, N2 purge gas
TEMP: 25.0 °C TIMES: 0.0 min RATE: 10.0 °C/min
TIME: 500.0 s

Temperature (°C)

SM FULTON
PERKIN ELMER
222-S Lab
Fri Sep 8 06:32:16 1995

LABCORE Data Entry Template for Worklist#

2317

Analyst: RJInstrument: DSC0 1Book # 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID	<u>28.45</u>	<u>30.4</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001966	0	DSC-01	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001966	0	DSC-01	SOLID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001967	0	DSC-01	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000118	BY-108 (R)	5 DUP	S95T001967	0	DSC-01	SOLID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u>	Joules/g

Final page for worklist #

2317

Rae Ann Green 9-29-95

Analyst Signature

Date

Dany Hammitt 10-5-95

Analyst Signature

Date

Verified by Blandina Valenzuela

10-5-95

S95T001966 produced two endotherms, one at 135.9°C with a delta H of 470.5 J/g and second at 274°C with a delta H of 44.2 J/g.

Data Entry Comments:

S95T001967 produced two endotherms, one at 108.0°C with a delta H of 596.0 J/g and second at 308.8°C with a delta H of 62.7 J/g.

10-2-95

bow

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-120

BEST AVAILABLE COPY


DSC STD 12N14A N2

8.630 mg

Rate: 10.0 °C/min

File: 00008.001 DSC METTLER 29-Sep-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


DSC S95T001966 SAM

49.470 mg

Rate: 10.0 °C/min

File: 00010.001 DSC METTLER 29-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

DSC S95T001966 DUP N2


29.200 mg

Rate: 10.0 °C/min

File: 00012.001 DSC METTLER 29-Sep-95

Ident: 0.0

222-S Laboratory



BEST AVAILABLE COPY

DSC S95T001967 SAM N2


19.840 mg

Rate: 10.0 °C/min

File: 00014.001 DSC METTLER 29-Sep-95

Ident: 0.0

222-S Laboratory



Rae Ann Greer

400. 9/29/95

°C

2-124

WHC-SD-WM-DP-145, REV. A

BEST AVAILABLE COPY


DSC S95T001967 DUP N2

13.264 mg

Rate: 10.0 °C/min

File: 00015.001 DSC METTLER 29-Sep-95

Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#

2318

Analyst: PJM Instrument: DSC0 3 Book # 12N14AMethod: LA-514-114 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-03	SOLID	<u>28.45</u>	<u>28.24</u>	N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001968 0	DSC-03	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001968 0	DSC-03	SOLID	<u>Ø</u>	<u>Ø</u>	N/A	Joules/g
		4 STD		DSC-03	SOLID	<u>28.45</u>	<u>28.21</u>	N/A	Joules/g
95000118	BY-108 (R)	5 SAMPLE	S95T001969 0	DSC-03	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000118	BY-108 (R)	6 DUP	S95T001969 0	DSC-03	SOLID	<u>Ø</u>	<u>Ø</u>	N/A	Joules/g

Final page for worklist # 2318See attached for signatures

Analyst Signature Date

Blandina Valenzuela 10-5-95

Analyst Signature Date

Verified by Blandina Valenzuela

(10-5-95)

S95T001968 produced two endotherms: one at 110.54°C with a delta H of 741.5 J/g and second at 219.3°C with a delta H 15.1 J/g.

Data Entry Comments: S95T001969 produced one endotherm at 134.8 °C with a delta H of 812.2 J/g.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-126

LABCORE Data Entry Template for Worklist#

2318

Analyst: RJMcClown
114

Instrument: DSC0

Book # 12N14AMethod: LA-514-113 Rev/Mod C-O
Atm 9/30/95

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001968 0	DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	3 DUP	S95T001968 0	DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001969 0	DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	5 DUP	S95T001969 0	DSC-01	SOLID			N/A	Joules/g

Final page for worklist #

2318

RJMcClown 9/30/95
Analyst Signature Date

Analyst Signature Date

Other instrument
was used.

10-5-95

BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-127


Curve 1: DSC

File info: IND092802 Thu Sep 28 16:37:00 1995

Sample Weight: 8.670 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY



N2, EXOTHERM DOWN

TEMP1: 158.0 G TIME1: 0.0 min RATE1: 10.0 C/min

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab

Thu Sep 28 16:42:20 1995

Curve 1: DSC


File info: SAM092808 Thu Sep 28 19:40:00 1995

Sample Weight: 19.760 mg

-6967001831 *Recd*
5957001968 9/28/95

BEST AVAILABLE COPY

Q-129



exotherm down, N2 purge gas
TEMP: 35.0 S TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

P.J. MCCOWN
PERKIN ELMER
222-S Lab
Thu Sep 28 19:56:09 1995

WHC-SD-WM-DP-145, REV.1

Curve 1: DSC


File info: SAM092809 Thu Sep 28 21:04:33 1995

Sample Weight: 43.400 mg

695T001831-DUP *Run*
595T001968 DUP

595T001968 DUP

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP1: 50.0 °C TEMP2: 0.0 °C RATE1: 10.0 °C/min
RATE2: 500.0 °C

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Thu Sep 28 21:38:09 1995

WHC-SD-WM-DP-145, REV 1

Curve 1: DSC


File info: IND093001 Sat Sep 30 18:05:39 1995

Sample Weight: 0.670 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY

T-131



WHC-SD-WM-DP-145, REV. 1

N₂, EXOTHERM DOWN

TEMP: 150.0 G TIMES: 0.0 min RATE: 10.0 C/min

Temperature (°C)

PJ MCCOWN

PERKIN ELMER

222-S Lab

Sat Sep 30 18:17:43 1995


Curve 1: DSC

File info: SAM093001 Sat Sep 30 19:29:16 1995

Sample Weight: 33.860 mg

S96T001969 SAM

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP: 35.0 °C TIMES: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sat Sep 30 19:55:44 1995


Curve 1: DSC

File info: SAM093002 Sat Sep 30 20:51:36 1995

Sample Weight: 29.360 mg

S95T001969 DUP

BEST AVAILABLE COPY



exotherm down, N₂ purge gas
TEMP: 35.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN

PERKIN ELMER

222-S Lab

Sat Sep 30 22:26:05 1995

WHC-SD-WM-DP-145, REV. 1

LABCORE Data Entry Template for Worklist#

2320

Analyst: ADP Instrument: DSC0 3 Book # 12 N14AMethod: LA-514-114 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-03	SOLID	<u>28.45</u>	<u>28.39</u>	N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	\$95T001970	0	DSC-03	N/A	Ø		Joules/g
95000118	BY-108 (R)	3 DUP	\$95T001970	0	DSC-03	SOLID	Ø	N/A	Joules/g
95000118	BY-108 (R)	4 SAMPLE	\$95T001971	0	DSC-03	N/A	Ø		Joules/g
95000118	BY-108 (R)	5 DUP	\$95T001971	0	DSC-03	SOLID	Ø	N/A	Joules/g

Final page for worklist #

2320

See attached for signatures

Analyst Signature Date

Dave Hammatt 10-5-95

Analyst Signature Date

Verified by Blandina
Valenzuela
10-5-95

\$95T001970 produced one endotherm at 102.7°C with a delta H
473.2 J/g.

Data Entry Comments: \$95T001971 produced one endotherm at 133.3°C with a delta H
of 611.4 J/g.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2320

Analyst: ADP Instrument: DSC0 Book # 12N14A

Method: LA-514-113 Rev/Mod LO

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001970 0	DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	3 DUP	S95T001970 0	DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001971 0	DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	5 DUP	S95T001971 0	DSC-01	SOLID			N/A	Joules/g

Final page for worklist # **2320**

Anthony Parikh 9-29-95

Analyst Signature

Date

Analyst Signature

Date

Other instrument
was used.

10-5-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: DSC

File info: IN0092902 Fri Sep 29 01:32:29 1995


Sample Weight: 8.670 mg

12N14A Indium at 10C/min

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2136 TO 2140

BEST AVAILABLE COPY

2-136



WHC-SD-WM-DP-145, REV. 1

Anthony Purinton 9-29-95

N2, EXOTHERM DOWN
TEMP: 150:8 8 TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)


AD PURINTON
PERKIN ELMER

222-S Lab

Fri Sep 29 01:38:14 1995

Curve 1: DSC
File info: SAM092901 Fri Sep 29 02:35:58 1995
Sample Weight: 14.380 mg
S95T001970

BEST AVAILABLE COPY



Curve 1: DSC


File info: SAM092902 Fri Sep 29 03:28:30 1995

Sample Weight: 18.910 mg

S95T001970 DUP

BEST AVAILABLE COPY

2CT-138



WHC-SD-VIM-DP-146, REV.1

exotherm down, N₂ purge gas
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min
TEMP: 400.0 °C

Temperature (°C)

AD PURINTON
PERKIN ELMER
222-S Lab
Fri Sep 29 03:34:58 1995


Curve 1: DSC

File info: SAM092903 Fri Sep 29 04:29:04 1995

Sample Weight: 27.510 mg

S95T001971 SAM

BEST AVAILABLE COPY



WHC-CD-WM-DP-146, REV. 1


exotherm down, N₂ purge gas
TEMP: 25.0 °C TIMES: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

AD PURINTON
PERKIN ELMER
222-S Lab
Fri Sep 29 04:30:16 1995

Curve 1: DSC
File info: SAM092904 Fri Sep 29 05:23:04 1995
Sample Weight: 27.450 mg
S95T001971 DUP

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

exotherm down, N₂ purge gas
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 400.0 °C

Temperature (°C)

AD PURINTON
PERKIN ELMER
222-S Lab
Fri Sep 29 05:28:46 1995

LABCORE Data Entry Template for Worklist#**2321**Analyst: DCDInstrument: DSC0 1Book # 12N14AMethod: LA-514-113 Rev/Mod C-0

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>29.9</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001972 0	DSC-01	SOLID	<u>N/A</u>	<u>41.4</u>		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001972 0	DSC-01	SOLID	<u>41.4</u>	<u>0</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	4 TRIP	S95T001972 0	DSC-01	SOLID	<u>41.4</u>	<u>40.5</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	5 SAMPLE	S95T001973 0	DSC-01	SOLID	<u>N/A</u>	<u>106.1</u>		Joules/g
95000118	BY-108 (R)	6 DUP	S95T001973 0	DSC-01	SOLID	<u>106.1</u>	<u>114.5</u>	<u>N/A</u>	Joules/g

Final page for worklist # **2321**

See attached for signatures
 Analyst Signature DCD Date 10-2-95
 Verified by Blandina Valenzuela
10-4-95

DCD
 Analyst Signature _____ Date _____

S95T001972 produced one endotherm at 131.9°C with a delta H of 945.6 J/g. The sample was run in triplicate due to high RPD's, the sample and Data Entry Comments: triplicate results were close.

S95T001973 produced one endotherm at 108.3°C with a delta H of 972.8 J/g, this peak has two distinct peaks within it.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2 141

LABCORE Data Entry Template for Worklist#

2321

Analyst: DCD Instrument: DSC0 _____ Book # 12 N14 A

Method: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S	TYPE	SAMPLE#	R	A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
			1 STD				DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	2	SAMPLE	S95T001972	0		DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	3	DUP	S95T001972	0		DSC-01	SOLID			N/A	Joules/g
95000118	BY-108 (R)	4	SAMPLE	S95T001973	0		DSC-01	SOLID	N/A			Joules/g
95000118	BY-108 (R)	5	DUP	S95T001973	0		DSC-01	SOLID			N/A	Joules/g

Final page for worklist # **2321**


Analyst Signature

Date

10-1-95

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-142

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-143 TO 2-148.


BEST AVAILABLE COPY

DSC STD 12N14A

6.670 mg

Rate: 10.0 °C/min

File: 00037.001 DSC METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV.1

2-143

David C. Dunham 10-1 5

BEST AVAILABLE COPY

S95T001972 SAM N2

30.696 mg

Rate: 10.0 °C/min


File: 00039.001

Ident: 0.0

DSC METTLER

01-Oct-95

222-S Laboratory



2-144

WHC-SC-NM-DR-145, REV. 1

BEST AVAILABLE COPY


S95T001972 DUP N2

28.062 mg

Rate: 10.0 °C/min

File: 00041.001 DSC METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory



J - 145


W4C-SD-V1A-DP-145, REV. I

BEST AVAILABLE COPY

S95T001972 TRIP N2
35.257 mg

Rate: 10.0 °C/min

File: 00043.001 DSC METTLER 01-oct-95
Ident: 0.0 222-S Laboratory



•C
400.
300.
200.
100.

BEST AVAILABLE COPY


S95T001973 SAM N2

10.203 mg

Rate: 10.0 °C/min

File: 00044.001 DSC METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory



WMC-CD-VIM-DP. 145, REV. L

BEST AVAILABLE COPY


S95T001973 DUP N2

15.217 mg

Rate: 10.0 °C/min

File: 00046.001 DSC METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory



F-148

WHC-SD-VIM-DP-145, REV. 1

LABCORE Data Entry Template for Worklist#**2322**Analyst: PJMInstrument: DSC0 3Book # 12N14AMethod: LA-514-114 Rev/Mod GO

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-03	SOLID	<u>28.45</u>	<u>28.41</u>	<u>N/A</u> Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001974	0	DSC-03	SOLID	<u>N/A</u>	<u>0</u>	<u>Joules/g</u>
95000118	BY-108 (R)	3 DUP	S95T001974	0	DSC-03	SOLID	<u>0</u>	<u>0</u>	<u>N/A</u> Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001975	0	DSC-03	SOLID	<u>N/A</u>	<u>0</u>	<u>Joules/g</u>
95000118	BY-108 (R)	5 DUP	S95T001975	0	DSC-03	SOLID	<u>0</u>	<u>0</u>	<u>N/A</u> Joules/g

Final page for worklist # **2322**See attached for signatures

Analyst Signature

Date

10-1-95

Dan Hamm 10-4-95

Analyst Signature

Date

Verified by Blandina
Valenzuela
10-5-95

10/2
S95T001974 produced one endotherm at 112.7°C with a delta H of 742.1 J/g

Data Entry Comments:

S95T001975 produced one endotherm at ~~112.7°C~~ ^{146.6°C} with a delta H of 718.5 J/g
DH/10/95

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#**2322**Analyst: RJMcClain

Instrument: DSC0

Book # 12N14AMethod: LA-514-113 Rev/Mod C-O
Rev 10/1/95

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID		N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001974	0	DSC-01	SOLID	N/A		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001974	0	DSC-01	SOLID		N/A	Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001975	0	DSC-01	SOLID	N/A		Joules/g
95000118	BY-108 (R)	5 DUP	S95T001975	0	DSC-01	SOLID		N/A	Joules/g

Final page for worklist #

2322RJMcClain
Analyst Signature10/1/95
Date

Analyst Signature _____ Date _____

Other instrument

was used

10-5-95

BWN

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: DSC

File info: IND100101 Sun Oct 1 00:07:50 1995

Sample Weight: 8.670 mg


12N14A Indium at 10C/min

BEST AVAILABLE COPY

PJ McCown
10/1/95

TCT-6

WHC-CD-WM-DP. 146, REV. 1



N2, EXOTHERM DOWN

TEMP: 150.0 °C TIME: 0.0 min RATE: 10.0 °C/min

TOP: 170.0 °C

Temperature (°C)

PJ MCCOWN
PEAKIN ELMER
222-S Lab

Sun Oct 1 00:33:38 1995

Curve 1: DSC

File info: SAM100103 Sun Oct 1 04:52:42 1995


Sample Weight: 22.710 mg

S95T001974 SAM

BEST AVAILABLE COPY

PJ McCown
10/1/95

WHC-SD-WM-DP-145, REV. 1



exotherm down, N₂ purge gas
TEMP1: 37.0 8 TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 200.0 8

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sun Oct 1 04:54:44 1995


Curve 1: DSC

File info: SAM100104 Sun Oct 1 06:03:33 1995

Sample Weight: 18.540 mg

S95T001974 DUP

BEST AVAILABLE COPY



WHC-CD-WM-DP-145, REV. 1

C-153

exotherm down, N₂ purge gas
TEMP: 35:8 8 TIME: 0.0 min RATE: 10.0 °C/min

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sun Oct 1 06:19:03 1995

Curve 1: DSC

File info: SAM100101 Sun Oct 1 01:44:17 1995


Sample Weight: 28.500 mg

Rm
10/1/95 S95T004969 SAM
S95T001975

BEST AVAILABLE COPY

PJ McCown
10/1/95

WHD-324/M-DP-146, REV. 1




exotherm down, N2 purge gas
TEMP: 35.0 °C TIME: 0.0 min RATE: 10.0 °C/min
TEMP: 200.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sun Oct 1 01:45:40 1995

Curve 1: DSC
File info: SAM100102 Sun Oct 1 02:42:55 1995
Sample Weight: 25.280 mg
S95T001975 DUP

BEST AVAILABLE COPY



WHD-CD-VHM-DP-146, REV. L

exotherm down, N₂ purge gas
TEMP: 25.8 8 TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sun Oct 1 03:38:00 1995

LABCORE Data Entry Template for Worklist#

2323

Analyst: RJMcClureInstrument: DSC0 1Book # 12N14AMethod: LA-514-113 Rev/Mod C-O

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-01	SOLID	<u>28.45</u>	<u>26.6</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001976 0	DSC-01	SOLID	<u>N/A</u>	<u>375.9</u>		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001976 0	DSC-01	SOLID	<u>375.9</u>	<u>350.0</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	4 SAMPLE	S95T001977 0	DSC-01	SOLID	<u>N/A</u>	<u>353.9</u>		Joules/g
95000118	BY-108 (R)	5 DUP	S95T001977 0	DSC-01	SOLID	<u>353.9</u>	<u>327.6</u>	<u>N/A</u>	Joules/g

Final page for worklist # 2323

RJMcClure 10/1/95
 Analyst Signature Date
 Verified by Blandina Valenzuela 10-4-95

Tom D. McThomas 10/2/95
 Analyst Signature Date

S95T001976 produced one endotherm at 118.2°C with a delta H of
 888.8 J/g.

Data Entry Comments: S95T001977 produced one endotherm at 134.0°C with a delta H
 of 3932.3 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2157 TO 2161.

BEST AVAILABLE COPY

DSC STD 12N14A

6.670 mg

Rate: 10.0 °C/min


File: 00027.001 DSC METTLER 01-Oct-95

Ident: 0.0

222-S Laboratory

RJMcClure
10/1/95

2157



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001976 SAM N2

15.896 mg

Rate: 10.0 °C/min

File: 00033.001 DSC METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory



85T-5

WITTMAYER / 145, REV. 1

BEST AVAILABLE COPY

S95T001976 DUP N2

25.883 mg

Rate: 10.0 °C/min

File: 00035.001 DSC METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

exo >

END

20. mW


Integration
Delta H 9060 mJ
350.0 J/g
Peak 362.3 °C
11.8 mW

Integration
Delta H 22067 mJ
852.6 J/g
Peak 136.9 °C
-49.2 mW


100. 200. 300. 400. °C

WWS-CO-V4M-DP-145 REV/

BEST AVAILABLE COPY



BEST AVAILABLE COPY



LABCORE Data Entry Template for Worklist#

2324

Analyst: PJM Instrument: DSC0 3 Book # 12N14A

Method: LA-514-114 Rev/Mod C-D

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		DSC-03	SOLID	<u>28.45</u>	<u>28.65</u>	<u>N/A</u>	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001978 0	DSC-03	SOLID	<u>N/A</u>	<u>Ø</u>		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001978 0	DSC-03	SOLID	<u>Ø</u>	<u>Ø</u>	<u>N/A</u>	Joules/g

Final page for worklist # **2324**

See attached for signatures

Analyst Signature Date 10/1/95

Tracy Johnson 10-04-95
Analyst Signature Date

Verified by Blandina Valenzuela 10-4-95

Data Entry Comments: Sample produced one endotherm at 133.9°C with
a delta H of 786.9 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-162

LABCORE Data Entry Template for Worklist#

2324

Analyst: RJ McCown

Instrument: DSC0

Book # 12N14A

Method: LA-514-113 Rev/Mod C-O
Run 10/1/95

Worklist Comment: Please run BY-108 DSCs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID		N/A	Joules/g
95000118	BY-108 (R)	2 SAMPLE	S95T001978 0		DSC-01	SOLID	N/A		Joules/g
95000118	BY-108 (R)	3 DUP	S95T001978 0		DSC-01	SOLID		N/A	Joules/g

Final page for worklist #

2324

RJ McCown
Analyst Signature

10/1/95
Date

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.


Curve 1: DSC

File info: IND100102 Sun Oct 1 16:18:22 1995

Sample Weight: 8.670 mg

12N14A Indium at 10C/min

BEST AVAILABLE COPY




Curve 1: DSC

File info: SAM100105 Sun Oct 1 18:00:08 1995

Sample Weight: 26.540 mg

S95T001978 SAM

2-165



WHC-SAM-DP-445, REV. 1

exotherm down, N2 purge gas
TEMP: 250.0 g TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN
PERKIN ELMER
222-S Lab
Sun Oct 1 18:03:59 1995


Curve 1: DSC

File info: SAM100106 Sun Oct 1 19:12:38 1995

Sample Weight: 27.160 mg

S95T001978 DUP

BEST AVAILABLE COPY



exotherm down. N2 purge gas
TEMP: 35.0 8 TIME: 0.0 min RATE: 10.0 °C/min
TEMP: 200.0 8

Temperature (°C)

PJ MCCOWN
PEAKIN ELMER
222-S Lab
Sun Oct 1 19:57:44 1995

WHC-SD-WM-DP-142, REV. 1

LABCORE Data Entry Template for Worklist#**2002**Analyst: JDSInstrument: TGA0 1Book #: 65N84Method: LA-560-112 Rev/Mod A.2

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>58.13</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001419 0	TGA-01	SOLID	<u>N/A</u>	<u>28.88</u>		%
95000104	BY-108 (R)	3 DUP	S95T001419 0	TGA-01	SOLID	<u>28.88</u>	<u>36.23</u>	<u>N/A</u>	%
95000104	BY-108 (R)	4 SAMPLE	S95T001420 0	TGA-01	SOLID	<u>N/A</u>	<u>35.41</u>		%
95000104	BY-108 (R)	5 DUP	S95T001420 0	TGA-01	SOLID	<u>35.41</u>	<u>35.81</u>	<u>N/A</u>	%

Final page for worklist # **2002**

Jah Sd 8-16-95
 Analyst Signature Date

L. J. Jones 8-22-95
 Analyst Signature Date

Verified by Blandina Valenzuela
8-28-95

S95T001419 produced a second weight loss step ^{BOV} ⁸⁻¹⁷⁻⁹⁵ of 6.22% at approximately 335°C
 Sample will be rerun due to high RPD's.

Data Entry Comments: S95T001420 produced a second weight loss step of 9.09% at
 approximately 320°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
 R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-168 TO 2-172.

BEST AVAILABLE COPY

TGA STD 65N8A


46.824 mg

Rate: 10.0 °C/min

File: 00018.001 TG METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height-27.22 mg
-58.13 %
ResiC. 19.61 mg
41.87 %
Dpeak 112.5 °C



WHC-SD-VMM-DP-142, REV. 1

BEST AVAILABLE COPY


S95T001419 SAM N2

32.486 mg

Rate: 10.0 °C/min

File: 00020.001 TG METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-VIM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001419 DUP N2

40.469 mg

Rate: 10.0 °C/min

File: 00023.001 TG METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-1425, REV. A

BEST AVAILABLE COPY


S95T001420 SAM N2

46.602 mg

Rate: 10.0 °C/min

File: 00025.001 TG METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001420 DUP N2

48.181 mg

Rate: 10.0 °C/min

File: 00027.001 TG METTLER 16-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. J

LABCORE Data Entry Template for Worklist#

2028

Analyst: RJM Instrument: TGA0 Book #: 65NB-AMethod: LA-560-112 Rev/Mod A-2

Worklist Comment: Please run BY-108 TGAs under N2. bdb

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	59.74	60.44	N/A %
95000104	BY-108 (R)	2 SAMPLE	S95T001419 1		TGA-01	SOLID	N/A	36.64	
95000104	BY-108 (R)	3 DUP	S95T001419 1		TGA-01	SOLID	36.64	36.69	N/A %

Final page for worklist # 2028

RJM 8/22/95
 Analyst Signature Date

RJM 8-22-95
 Analyst Signature Date

Verified by Blandina Valenzuela
 8-28-95

Data Entry Comments: The sample and duplicate were run to an upper temp. limit of 250°C, since ⁸⁻²⁸⁻⁹⁵ we are concerned with the RPDS of the water weight loss the chemist did not run to upper temperature limit of 500°C due to ALARA concerns.
 Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2173

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-174 TO 2-176.


BEST AVAILABLE COPY

TGASTD 65N8A

17.683 mg

Rate: 10.0 °C/min

File: 00046.001 TG METTLER 22-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-145, REV. 1

BEST AVAILABLE COPY


S95T001419 N2

31.806 mg

Rate: 10.0 °C/min

File: 00047.001 TG METTLER 22-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV 1


BEST AVAILABLE COPY

S95T001419DUP N2

33.507 mg

Rate: 10.0 °C/min

File: 00048.001 TG METTLER 22-Aug-95
Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#

2029

Analyst: SMF Instrument: TGA0 3 Book # 65N8-AMethod: LA-514-114 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.41</u>	N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001380 0	TGA-03	SOLID	<u>15.74</u>	<u>15.74</u>	GmF 9/5/95	%
95000104	BY-108 (R)	3 DUP	S95T001380 0	TGA-03	SOLID	<u>38.17</u>	<u>38.17</u>	15.92	GmF 9/5/95
95000104	BY-108 (R)	4 SAMPLE	S95T001381 0	TGA-03	SOLID	<u>38.17</u>	<u>39.61</u>	N/A	%
95000104	BY-108 (R)	5 DUP	S95T001381 0	TGA-03	SOLID	<u>15.74</u>	<u>15.18</u>	GmF 9/5/95	%
					SOLID	<u>15.74</u>	<u>16.82</u>	GmF 9/5/95	%
					SOLID	<u>42.35</u>	<u>43.14</u>	N/A	%

Final page for worklist # 2029

See attached for signatures

Analyst Signature

Date

8-21-95

BDV

J. Jones

9-23-95

Analyst Signature

Date

Corrected &Verified 9/5/95 J. M. Lye

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2029

Analyst: SMF Instrument: TGA0 Book # 165N87Method: LA-560-112 Rev/Mod BcD

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID			N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001380	0	TGA-01	SOLID	N/A		%
95000104	BY-108 (R)	3 DUP	S95T001380	0	TGA-01	SOLID		N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001381	0	TGA-01	SOLID	N/A		%
95000104	BY-108 (R)	5 DUP	S95T001381	0	TGA-01	SOLID		N/A	%

Final page for worklist # 2029

Susie M. Fulton 8-19-95

Analyst Signature Date

Analyst Signature Date

Other instrument was
used.

8-21-95

BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-178

Curve 1: TGA


File info: TER081901 Sat Aug 19 08:57:31 1995

Sample Weight: 14.556 mg

65N8-A Terlg

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-179 TO 2-183.

BEST AVAILABLE COPY



N2
TEMP1: 30.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 250.0 C


Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Sat Aug 19 09:01:03 1995

WHC-SD-WM-DP-145, REV. J

Curve 1: TGA
File info: SAM081901 Sat Aug 19 11:10:26 1995
Sample Weight: 14.827 mg
S95T001380 10C/min

BEST AVAILABLE COPY



N2
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)


SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 10:26:52 1995

WHC-SD-WM-DP-145, REV. 1

Curve 1: TGA
File info: SAM081902 Sat Aug 19 12:24:43 1995
Sample Weight: 18.051 mg
S95T001380DUP 10C/min

BEST AVAILABLE COPY

C-181



WHD-SD-WM-DP-146, REV. 1

N2
TEMP1: 20.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min


Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 10:39:16 1995

Curve 1: TGA
File info: SAM081903 Sat Aug 19 14:10:05 1995
Sample Weight: 17.864 mg
S95T001381 10C/min

BEST AVAILABLE COPY

2-182



WHC-SD-WM-DP-145, REV. J

N2
TEMP: 500.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 13:43:26 1995

Curve 1: TGA


File info: SAM081904 Sat Aug 19 15:25:27 1995

Sample Weight: 15.711 mg

S95T001381DUP 10C/min

BEST AVAILABLE COPY

2-183



N2

TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

TEMP: 500.0 °C

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 13:50:07 1995

WHC-SD-WM-DP-145, REV 1

LABCORE Data Entry Template for Worklist#**2030**Analyst: SMF Instrument: TGA0 1 Book # 65N8AMethod: LA-560-112 Rev/Mod A-2

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	8.21-95 SOLID BDV	<u>59.74</u> <u>60.</u>	<u>60.23</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001382	0	TGA-01	SOLID	<u>N/A</u>	<u>39.16</u>		%
95000104	BY-108 (R)	3 DUP	S95T001382	0	TGA-01	SOLID	<u>39.16</u>	<u>39.44</u>	<u>N/A</u>	%

Final page for worklist # **2030**

Susie M. Fulton 8/19/95
 Analyst Signature Date

L.J.Jones 8-21-95
 Analyst Signature Date

Verified by Blandina Valenzuela
 8-21-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-184

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2185 TO 2187.


BEST AVAILABLE COPY

TGASTD 65N8A

19.766 mg

Rate: 10.0 °C/min

File: 00034.001 TG METTLER 19-Aug-95
Ident: 0.0 222-S Laboratory



Laurie M. Fulton 8/19/95

BEST AVAILABLE COPY


S95T001382 N2

29.085 mg

Rate: 10.0 °C/min

File: 00035.001 TG METTLER 19-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001382DUP N2

26.614 mg

Rate: 10.0 °C/min

File: 00036.001 TG METTLER 19-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-145, REV 1

LABCORE Data Entry Template for Worklist#

2033

Analyst: JDS Instrument: TGA0 3 Book # 65N8AMethod: LA-514-114 Rev/Mod B-D

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S	TYPE	SAMPLE#	R	A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
			1 STD				TGA-03	SOLID	<u>59.74</u>	<u>57.60</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2	SAMPLE	S95T001390	0		TGA-03	SOLID	<u>N/A</u>	<u>21.85</u>		%
95000104	BY-108 (R)	3	DUP	S95T001390	0		TGA-03	SOLID	<u>21.85</u>	<u>19.76</u>	<u>N/A</u>	%
		4	STD				TGA-03	SOLID	<u>59.74</u>	<u>57.56</u>	^{8-24-95 BDV} <u>N/A</u>	%
95000104	BY-108 (R)	5	SAMPLE	S95T001407	0		TGA-03	SOLID	<u>N/A</u>	<u>16.31</u>		%
95000104	BY-108 (R)	6	DUP	S95T001407	0		TGA-03	SOLID	<u>16.31</u>	<u>17.58</u>	<u>N/A</u>	%

Final page for worklist # 2033See attached for signaturesAnalyst Signature Date 8-24-95Dany HanusaDate 8-24-95Verified by Blandina Valenzuela BDV
8-28-95S95T001390 produced two additional weight loss steps of 16.77%
and 13.45% before 250°CData Entry Comments: S95T001407 produced one additional weight loss step of
22.45% before 250°C

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-188

LABCORE Data Entry Template for Worklist#

2033Analyst: Jds

Instrument: TGA0

Book # 65N8AMethod: LA-560-112 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID			N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001390 0	TGA-01	SOLID	N/A			%
95000104	BY-108 (R)	3 DUP	S95T001390 0	TGA-01	SOLID			N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001407 0	TGA-01	SOLID	N/A			%
95000104	BY-108 (R)	5 DUP	S95T001407 0	TGA-01	SOLID			N/A	%

Final page for worklist #

2033Jds8-23-95

Date

Analyst Signature

Date

Other instrument

was used

8-24-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: TGA

File info: TER081801 Fri Aug 18 09:43:05 1995


Sample Weight: 16.603 mg

65NB-A TGA

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-190 TO 2-195.

BEST AVAILABLE COPY

WHC-SD-WM-DP-145, REV. 1



N2
TEMP1: 35.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 250.0 C

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Fri Aug 18 10:25:43 1995


Curve 1: TGA

File info: SAM081801 Fri Aug 18 11:21:19 1995

Sample Weight: 15.720 mg

S95T001390 SAM N2

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

N2

TEMP1: 50.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 24 11:50:05 1995


Curve 1: TGA

File info: SAM081802 Fri Aug 18 13:49:37 1995

Sample Weight: 20.453 mg

S95T001390 DUP

BEST AVAILABLE COPY



N2
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 24 12:14:28 1995

Curve 1: TGA TER 082301


File info: IN0082301 Wed Aug 23 08:54:21 1995

Sample Weight: 17.856 mg

65N8-A Terliq

BEST AVAILABLE COPY

297.193



WHC-SD-WM-DP-145, REV. 1

N2
TEMP1: 98.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 200.0 C


Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 23 11:02:17 1995

Curve 1: TGA
File info: SAM082301 Wed Aug 23 13:03:23 1995
Sample Weight: 19.261 mg
S95T001407 SAM

BEST AVAILABLE COPY

2-194



WHC-SD-W/M-DP. 145, REV. L

N2 10C/MIN
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 24 12:56:52 1995

Curve 1: TGA

File info: SAM082302 Wed Aug 23 14:44:21 1995


Sample Weight: 17.046 mg

S95T001407 DUP

BEST AVAILABLE COPY

J-195

WHC-SD-WM-DR-145, REV. 1



N2 10C/MIN

TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN

PERKIN-ELMER

7 Series Thermal Analysis System

Thu Aug 24 12:51:51 1995

LABCORE Data Entry Template for Worklist#

2034

Analyst: SMF Instrument: TGA0 1 Book # 65N8AMethod: LA-560-112 Rev/Mod A-Z

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.51</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001408	0	TGA-01	SOLID	<u>N/A</u>	<u>37.28</u>		
95000104	BY-108 (R)	3 DUP	S95T001408	0	TGA-01	SOLID	<u>37.28</u>	<u>36.093</u>	<u>N/A</u>	%
95000104	BY-108 (R)	4 SAMPLE	S95T001409	0	TGA-01	SOLID	<u>N/A</u>	<u>37.47</u>		%
95000104	BY-108 (R)	5 DUP	S95T001409	0	TGA-01	SOLID	<u>37.47</u>	<u>37.91</u>	<u>N/A</u>	%

Final page for worklist # 2034

Smulter8-20-95

Analyst Signature

Date

Verified for Blandina Valenzuela8-21-95bdvLJ Jones8-21-95

Analyst Signature

Date

S95T001408 produced a second weight loss step of 5.26%
 at approximately 340°C

Data Entry Comments:

S95T001409 produced a8-20-95
BDV

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
 R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-197 TO 2-201.


BEST AVAILABLE COPY

TGA STD 65N8A

13.971 mg

Rate: 10.0 °C/min

File: 00039.001 TG METTLER 20-Aug-95
Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001408 N2

23.501 mg

Rate: 10.0 °C/min

File: 00040.001 TG METTLER 20-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

S95T001408DUP N2

25.439 mg


Rate: 10.0 °C/min

File: 00041.001 TG METTLER 20-Aug-95
Ident: 0.0 222-S Laboratory

Step Analysis
Height -9.17 mg
-36.03 %
ResiC. 16.27 mg
63.97 %
Dpeak 75.0 °C

Step Analysis
Height -1.19 mg
-4.66 %
ResiC. 14.64 mg
57.57 %

66T-C



WHC-SD-WM-DP. /45, REV. 1

BEST AVAILABLE COPY

S95T001409 N2

26.766 mg


Rate: 10.0 °C/min

File: 00042.001 TG METTLER 20-Aug-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height-10.03 mg
-37.47 %
Resid. 16.74 mg
62.53 %
Dpeak 113.0 °C

003-7



NHC-SD-WM-DP. 145, REV. L

BEST AVAILABLE COPY


S95T001409DUP N2

31.227 mg

Rate: 10.0 °C/min

File: 00043.001 TG METTLER 20-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. A

WHC-SD-WM-DP-145, REV.1

THE FOLLOWING ANALYSES ARE INCLUDED IN THE DATA PACKAGE BUT THE RESULTS HAVE
NOT BEEN REPORTED IN THE FINAL SUMMARY REPORTS.

THIS PAGE INTENTIONALLY
LEFT BLANK

LABCORE Data Entry Template for Worklist#

2035

Analyst: RDM Instrument: TGA0 A-2 Book # 65N8-AMethod: LA-560-112 Rev/Mod A-2

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.73</u>	N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001410 0	TGA-01	SOLID	<u>N/A</u>	<u>35.32</u>		%
95000104	BY-108 (R)	3 DUP	S95T001410 0	TGA-01	SOLID	<u>35.32</u>	<u>39.45</u>	N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001396 0	TGA-01	SOLID	<u>N/A</u>	<u>40.98</u>		%
95000104	BY-108 (R)	5 DUP	S95T001396 0	TGA-01	SOLID	<u>40.98</u>	<u>44.12</u>	N/A	%

Final page for worklist # 2035

RDM

8/23/95

Analyst Signature

LJ

8-24-95

Analyst Signature

Verified & put out
 for review
Jean M. Faye
 9/4/95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-203

BEST AVAILABLE COPY

TGASTD 65N8A


17.169 mg

Rate: 10.0 °C/min

File: 00049.001 TG METTLER 23-Aug-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height-10.43 mg
-60.73 %
ResiC. 6.74 mg
39.27 %
Dpeak 80.8 °C



WHC-SD-WM-DP-145, REV. 1

R. D. Dugay


BEST AVAILABLE COPY

S95T001410 N2

41.448 mg

Rate: 10.0 °C/min

File: 00052.001 TG METTLER 23-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-1452, REV. 1

R. Johnson


BEST AVAILABLE COPY

S95T001410DUP N2

47.732 mg

Rate: 10.0 °C/min

File: 00053.001 TG METTLER 23-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

R. D. Bryan


BEST AVAILABLE COPY

S95T001396 N2

51.338 mg


Rate: 10.0 °C/min

File: 00050.001 TG METTLER 23-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-145, REV. 1

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 201 TO 293.

BEST AVAILABLE COPY

TGA STD 65N8A

30.452 mg


Rate: 10.0 °C/min

File: 00029.001 TG METTLER 18-Aug-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height-18.32 mg
-60.16 %
ResiC. 12.12 mg
39.79 %
Dpeak 105.0 °C

X-211



WHC-SD-WM-DP-1445, REV. 1


BEST AVAILABLE COPY

S95T001399 SAM N2

29.504 mg

Rate: 10.0 °C/min

File: 00031.001 TG METTLER 18-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-1452, REV. 1

BEST AVAILABLE COPY

S95T001399 DUP N2

35.937 mg


Rate: 10.0 °C/min

File: 00033.001 TG METTLER 18-Aug-95

Ident: 0.0 222-S Laboratory

2-213

Step Analysis
Height-14.28 mg
-39.74 %
ResiC. 21.66 mg
60.26 %
Dpeak 137.0 °C



WHC-SD-WM-DP. 1/95 REV. L

LABCORE Data Entry Template for Worklist#

2037

Analyst: JDSInstrument: TGA0 3Book # 65N8AMethod: LA-514-114 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	LIQUID	<u>59.74</u>	<u>60.39</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001373 0	TGA-03	LIQUID	<u>N/A</u>	<u>31.9</u>		%
95000104	BY-108 (R)	3 DUP	S95T001373 0	TGA-03	LIQUID	<u>31.9</u>	<u>34.01</u>	<u>N/A</u>	%
95000104	BY-108 (R)	4 SAMPLE	S95T001427 0	TGA-03	LIQUID	<u>N/A</u>	<u>27.11</u>		%
95000104	BY-108 (R)	5 DUP	S95T001427 0	TGA-03	LIQUID	<u>27.11</u>	<u>25.03</u>	<u>N/A</u>	%

Final page for worklist # **2037**See attached for signatures

Analyst Signature

Date

8-31-95

Analyst Signature

Date

Entered & Verified Jawn Lye
9/6/95

Data Entry Comments:

S95T001373 SAM has second wt loss step of 20.91% finishing by 250°C;
DUP has 18.7%. S95T001427 SAM has a small step of 25.62% ending
by 250°C; duplicate has 27.14%.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-214

LABCORE Data Entry Template for Worklist#**2037**Analyst: JDS Instrument: TGA0 Book #: 65N8AMethod: LA-560-112 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	LIQUID			N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001427 0	TGA-01	LIQUID	N/A			%
95000104	BY-108 (R)	3 DUP	S95T001427 0	TGA-01	LIQUID			N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001373 0	TGA-01	LIQUID	N/A			%
95000104	BY-108 (R)	5 DUP	S95T001373 0	TGA-01	LIQUID			N/A	%

Final page for worklist # 2037JDS

8-29-95

Analyst Signature

Date

Analyst Signature

Date

(Other instrument was
used.)8-31-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-215


Curve 1: TGA

File info: ter082901 Tue Aug 29 08:24:03 1995

Sample Weight: 26.716 mg

65N8-A Terliq

BEST AVAILABLE COPY



N2 10C/MIN

TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 300.0 °C

Temperature (°C)

JD SPELLMAN

PERKIN-ELMER

7 Series Thermal Analysis System

Tue Aug 29 08:24:24 1995

WHC-SD-WM-DP-145, REV. 1


Curve 1: TGA

File info: SAM082903 Tue Aug 29 13:56:54 1995

Sample Weight: 12.024 mg

S95T001373 SAM

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1


Curve 1: TGA

File info: SAN082904 Tue Aug 29 15:05:31 1995

Sample Weight: 14.372 mg

S95T001373 DUP

BEST AVAILABLE COPY



N2 10C/MIN

TEMP1: 30.0 °C TEMP2: 300.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 31 09:34:19 1995

WHC-SD-WM-DP-145, REV. 1


Curve 1: TGA

File info: SAM082901 Tue Aug 29 09:27:54 1995

Sample Weight: 12.247 mg

S95T001427 SAM

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1


N2 10C/MIN
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 31 10:13:03 1995

Curve 1: TGA
File info: SAM082902 Tue Aug 29 11:55:45 1995
Sample Weight: 15.783 mg
S95T001427 DUP

BEST AVAILABLE COPY



WHCSD-WMDP-145, REV. L

N2 10C/MIN
TEMP1: 35.0 C TIME1: 0.0 min RATE1: 10.0 C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 31 10:00:03 1995

LABCORE Data Entry Template for Worklist#

2155Analyst: SMF Instrument: TGA0 BV 3/1 Book # 65N8-AMethod: LA-560-112 Rev/Mod A-2

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.65</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001421	0	TGA-01	SOLID	<u>N/A</u>	<u>38.31</u>		%
95000104	BY-108 (R)	3 DUP	S95T001421	0	TGA-01	SOLID	<u>38.31</u>	<u>39.82</u>	<u>N/A</u>	%
95000104	BY-108 (R)	4 SAMPLE	S95T001422	0	TGA-01	SOLID	<u>N/A</u>	<u>36.49</u>		%
95000104	BY-108 (R)	5 DUP	S95T001422	0	TGA-01	SOLID	<u>36.49</u>	<u>36.40</u>	<u>N/A</u>	%

Final page for worklist # **2155**Susie M. Dalton Date

8-29-95

Dan Hamm Date

8-30-95

Analyst Signature

Date

Analyst Signature

Date

Verified by Blandina Valenzuela
8-31-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-221

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 222 TO 226.


BEST AVAILABLE COPY

TGA STD 65N8A

23.768 mg

Rate: 10.0 °C/min

File: 00093.001 TG METTLER 30-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR-145, REV 1


Susie M. Fulton 8-29-95

BEST AVAILABLE COPY

S95T001421 N2
28.560 mg

Rate: 10.0 °C/min

File: 00095.001 TG METTLER 30-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY


S95T001421 DUP N2

35.205 mg

Rate: 10.0 °C/min

File: 00097.001 TG METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

S95T001422 N2

9.544 mg


Rate: 10.0 °C/min

File: 00099.001 TG METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height -3.48 mg
-36.49 %
ResiC. 6.06 mg
63.51 %

2225



BEST AVAILABLE COPY


S95T001422 DUP N2

36.802 mg

Rate: 10.0 °C/min

File: 00101.001 TG METTLER 30-Aug-95
Ident: 0.0 222-S Laboratory

2226



WHC-SD-WM-DP. 145, REV. 1

LABCORE Data Entry Template for Worklist#**2157**Analyst: Jds Instrument: TGA0 1 Book #: 65 N 8AMethod: LA-560-112 Rev/Mod A2

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.82</u>	<u>N/A</u> %
95000104	BY-108 (R)	2 SAMPLE	S95T001433 0		TGA-01	SOLID	<u>N/A</u>	<u>34.52</u>	
95000104	BY-108 (R)	3 DUP	S95T001433 0		TGA-01	SOLID	<u>34.52</u>	<u>38.56</u>	<u>N/A</u> %

Final page for worklist #**2157**Jds

8-30-95

Analyst Signature

Date

D. Hammatt

8-30-95

Analyst Signature

Date

Verified by Blandina Valenzuela

8-31-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-227

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 7-228 TO 7-230.


BEST AVAILABLE COPY

TGA STD 65N8A

25.498 mg

Rate: 10.0 °C/min

File: 00103.001 TG METTLER 30-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. /

BEST AVAILABLE COPY


S95T001433 SAM N2

19.487 mg

Rate: 10.0 °C/min

File: 00105.001 TG METTLER 30-Aug-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR. 145, REV. 1


BEST AVAILABLE COPY

S95T001433 DUP N2

19.489 mg

Rate: 10.0 °C/min

File: 00107.001 TG METTLER 30-Aug-95
Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#

2158

Analyst: JDS Instrument: TGA0 1 Book # 65N8AMethod: LA-560-112 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.37</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001389	0	TGA-01	SOLID	<u>N/A</u>	<u>37.72</u>		%
95000104	BY-108 (R)	3 DUP	S95T001389	0	TGA-01	SOLID	<u>37.72</u>	<u>36.28</u>	<u>N/A</u>	%
		4 STD			TGA-01	SOLID	<u>59.05</u> <small>74 6W 9-A-95</small>	<u>59.05</u>	<u>N/A</u>	%
95000104	BY-108 (R)	5 SAMPLE	S95T001393	0	TGA-01	SOLID	<u>N/A</u>	<u>38.92</u>		%
95000104	BY-108 (R)	6 DUP	S95T001393	0	TGA-01	SOLID	<u>38.92</u>	<u>36.38</u>	<u>N/A</u>	%

Final page for worklist # 2158

See attached

Analyst Signature

Date

Re Jones

9-11-95

Analyst Signature

Date

Verified 9/11/95 Jan M. Lye

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-231

LABCORE Data Entry Template for Worklist#**2158**Analyst: JDS Instrument: TGA0 Book # 65 NPAMethod: LA-560-112 Rev/Mod B-D

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID		N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001389 0		TGA-01	SOLID	N/A		%
95000104	BY-108 (R)	3 DUP	S95T001389 0		TGA-01	SOLID		N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001393 0		TGA-01	SOLID	N/A		%
95000104	BY-108 (R)	5 DUP	S95T001393 0		TGA-01	SOLID		N/A	%

Final page for worklist # **2158**Jdh Sjd

4-8-95

Analyst Signature

Date

Analyst Signature

Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2232

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-233 TO 2-238.


BEST AVAILABLE COPY

TGA STD 65N8A

24.084 mg

Rate: 10.0 °C/min

File: 00087.001 TG METTLER 07-Sep-95
Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

S95T001389 SAM N2

23.233 mg


Rate: 10.0 °C/min

File: 00089.001 TG METTLER 07-Sep-95

Ident: 0.0 222-S Laboratory

2-234

Step Analysis
Height -8.76 mg
-37.72 %
ResiC. 14.43 mg
62.13 %
Dpeak 131.0 °C



WHC-SD-WM-DP- 1/15, REV. 1

BEST AVAILABLE COPY


S95T001389 DUP N2

19.006 mg

Rate: 10.0 °C/min

File: 00091.001 TG METTLER 07-Sep-95

Ident: 0.0 222-S Laboratory




BEST AVAILABLE COPY

TGA STD 65N8A

29.063 mg

Rate: 10.0 °C/min

File: 00105.001 TG METTLER 08-Sep-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP- /45, REV. 1

BEST AVAILABLE COPY


S95T001393 SAM N2

30.459 mg

Rate: 10.0 °C/min

File: 00107.001 TG METTLER 08-Sep-95

Ident: 0.0 222-S Laboratory



WHC SD-WM-DP-145, REV. L

BEST AVAILABLE COPY


S95T001393 DUP N2

24.051 mg

Rate: 10.0 °C/min

File: 00109.001 TG METTLER 08-Sep-95
Ident: 0.0 222-S Laboratory

2-238



WHC-SD-WM-DP. /45, REV. 1

LABCORE Data Entry Template for Worklist#

2159

Analyst: JDS Instrument: TGA0 1 Book # 6SN8A

Method: LA-560-112 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.28</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001432 0	TGA-01	SOLID	<u>N/A</u>	<u>41.31</u>		%
95000104	BY-108 (R)	3 DUP	S95T001432 0	TGA-01	SOLID	<u>41.31</u>	<u>38.53</u>	<u>N/A</u>	%

Final page for worklist # **2159**

See Attached for Signature
Analyst Signature Date JMF 9/7/95

Verified by Blandina Valenzuela
9-28-95

Lis Jones 9-26-95
Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-239

LABCORE Data Entry Template for Worklist#

2159

Analyst: JLS Instrument: TGA0 / Book # 65N8A

Method: LA-560-112 Rev/Mod B-C

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.28</u>	N/A %
95000104	BY-108 (R)	2 SAMPLE	S95T001432	0	TGA-01	SOLID	<u>N/A</u>		
95000104	BY-108 (R)	3 DUP	S95T001432	0	TGA-01	SOLID			N/A %
95000104	BY-108 (R)	4 SAMPLE	S95T001431	0	TGA-01	SOLID	<u>N/A</u>		
95000104	BY-108 (R)	5 DUP	S95T001431	0	TGA-01	SOLID			N/A %

Final page for worklist # **2159**

Jah Spt 9-5-95

Analyst Signature

Date

Analyst Signature

Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-240

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 20241 TO 20243.

BEST AVAILABLE COPY


TGA STD 65N8A

26.668 mg

Rate: 10.0 °C/min

File: 00002.001 TG METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001431 SAM N2

35.044 mg

Rate: 10.0 °C/min

File: 00004.001 TG METTLER 31-Aug-95
Ident: 0.0 222-S Laboratory

2-242



WHC-SD-WM-DP-245, REV. 1


BEST AVAILABLE COPY

S95T001431 DUP N2

17.221 mg

Rate: 10.0 °C/min

File: 00006.001 TG METTLER 31-Aug-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145 REV C

LABCORE Data Entry Template for Worklist#**2160**Analyst: SMFInstrument: TGA0 3Book # 65NB-AMethod: LA-514-114 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.00</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001586 0	TGA-03	SOLID	<u>N/A</u>	<u>15.74</u>		x
95000118	BY-108 (R)	3 DUP	S95T001586 0	TGA-03	SOLID	<u>15.74</u>	<u>15.02</u>	<u>N/A</u>	x
95000118	BY-108 (R)	4 SAMPLE	S95T001587 0	TGA-03	SOLID	<u>N/A</u>	<u>10.41</u>		x
95000118	BY-108 (R)	5 DUP	S95T001587 0	TGA-03	SOLID	<u>10.41</u>	<u>12.77</u>	<u>N/A</u>	x
95000118	BY-108 (R)	6 TRIP	S95T001587 0	TGA-03	SOLID	<u>10.41</u>	<u>13.00</u>	<u>N/A</u>	x

Final page for worklist #

2160See attached for signatures

Analyst Signature

Date

8-31-95

LSJ

8-31-95

Date

Verified by Blandina Valenzuela
^{BDV}
8-31-95

Data Entry Comments: S95T001587 produced a second weight loss step of 25.41% at approximately 175°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-244

LABCORE Data Entry Template for Worklist#

2160

Analyst: SMF Instrument: TGA0 Book # 65N8-A

Method: LA-560-112 Rev/Mod LA-514-114/B0
SMF 8-30-95

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001586	0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	3 DUP	S95T001586	0	TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	4 SAMPLE	S95T001587	0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	5 DUP	S95T001587	0	TGA-01	SOLID			N/A	%

Final page for worklist # **2160**

Smulfor 8-30-95
Analyst Signature Date

Analyst Signature Date

Other instrument was
used

8-31-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-245

Curve 1: TGA


File info: TER083001 Wed Aug 30 16: 42: 54 1995

Sample Weight: 20.662 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2246 TO 2251.

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. 1

N2 10C/MIN
TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 300.0 °C

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 16: 46: 35 1995


Curve 1: TGA

File info: SAM083001 Wed Aug 30 17:46:11 1995

Sample Weight: 21.617 mg

S95T001586

BEST AVAILABLE COPY



N2 10C/MIN
TEMP: 30.0 C TZME: 0.0 min RATE: 10.0 C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Aug 30 18:16:32 1995

WHC-SD-WM-DP / 145, REV. 1


Curve 1: TGA

File info: SAM083002 Wed Aug 30 19:09:47 1995

Sample Weight: 16.559 mg

S95T001586 DUP

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV 1

N2 10C/MIN

TEMP1: 35.0 C TIMES: 0.0 min RATE1: 10.0 C/min

Temperature (°C)

SM FULTON

PERKIN-ELMER

7 Series Thermal Analysis System
Wed Aug 30 19:22:54 1995


Curve 1: TGA

File info: SAM083003 Wed Aug 30 20:18:32 1995

Sample Weight: 29.059 mg

S95T001587

BEST AVAILABLE COPY



N2 10C/MIN
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 31 10:42:08 1995

WHC-SD-WM-DR-145, REV. L


Curve 1: TGA

File info: SAM083004 Wed Aug 30 22:10:07 1995

Sample Weight: 20.160 mg

S95T001587 DUP

BEST AVAILABLE COPY



WHC-SD-WM-DR-145, REV. 1

N2 10C/MIN

TEMP1: 35.0 C TIME1: 0.0 min RATE1: 10.0 C/min

TEMP2: 500.0 C

Temperature (°C)

SM FULTON
PERKIN-ELMER

7 Series Thermal Analysis System
Thu Aug 31 10:56:08 1995

Curve 1: TGA


File info: SAM083005 Wed Aug 30 23:16:10 1995

Sample Weight: 18.168 mg

S95T001587 TRIP

BEST AVAILABLE COPY

L-251



WHC-SD-WM-DP-145, REV. L

N2 100/C/MIN
TEMP: 35.0 C TIMES: 0.0 min RATE: 10.0 C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Aug 31 11:01:42 1995


LABCORE Data Entry Template for Worklist#**2161**Analyst: SMF Instrument: TGA0 1 Book # 65N8-AMethod: LA-560-112 Rev/Mod A-2

Worklist Comment: Please run BY-108 TGAs under N2. bdv


GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.47</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001588 0	TGA-01	SOLID	<u>N/A</u>	<u>29.88</u>		%
95000118	BY-108 (R)	3 DUP	S95T001588 0	TGA-01	SOLID	<u>29.88</u>	<u>30.68</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001589 0	TGA-01	SOLID	<u>N/A</u>	<u>22.95</u>		%
95000118	BY-108 (R)	5 DUP	S95T001589 0	TGA-01	SOLID	<u>22.95</u>	<u>23.67</u>	<u>N/A</u>	%

Final page for worklist # **2161**


Susie M. Dalton 8-30-95
Analyst Signature Date



L.J. Jones 8-31-95
Analyst Signature Date



Verified 9/1
9/1/95 J.M. Saxe

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-252

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2053 TO 2057.


BEST AVAILABLE COPY

TGA STD 65N8A

13.534 mg

Rate: 10.0 °C/min

File: 00108.001 TG METTLER 31-Aug-95
Ident: 0.0 222-S Laboratory




BEST AVAILABLE COPY

S95T001588 N2

23.835 mg

Rate: 10.0 °C/min

File: 00109.001 TG METTLER 31-Aug-95
Ident: 0.0 222-S Laboratory




BEST AVAILABLE COPY

S95T001588 DUP N2

29.475 mg

Rate: 10.0 °C/min

File: 00110.001 TG METTLER 31-Aug-95
Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001589 N2

35.722 mg

Rate: 10.0 °C/min

File: 00111.001 TG METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001589 DUP N2

37.429 mg

Rate: 10.0 °C/min

File: 00112.001 TG METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



THIS PAGE INTENTIONALLY
LEFT BLANK

BEST AVAILABLE COPY


S95T001589 DUP N2

37.429 mg

Rate: 10.0 °C/min

File: 00112.001 TG METTLER 31-Aug-95

Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#**2200**Analyst: RDM Instrument: TGA0 1 Book # 65N8-AMethod: LA-560-112 Rev/Mod B-D

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	LIQUID	<u>59.74</u>	<u>60.60</u>	<u>N/A</u>	%
95000101	BY-108 (R)	2 SAMPLE	S95T001365 0	TGA-01	LIQUID	<u>N/A</u>	<u>100.0</u>		%
95000101	BY-108 (R)	3 DUP	S95T001365 0	TGA-01	LIQUID	<u>100.0</u>	<u>100.0</u>	<u>N/A</u>	%
95000104	BY-108 (R)	4 SAMPLE	S95T001430 0	TGA-01	LIQUID	<u>N/A</u>	<u>38.58</u>		%
95000104	BY-108 (R)	5 DUP	S95T001430 0	TGA-01	LIQUID	<u>38.58</u>	<u>39.47</u>	<u>N/A</u>	%

Final page for worklist # **2200**RDM9/20/95

Analyst Signature

Date

L. Jones9-20-95

Analyst Signature

Date

9-21-95

BDV

Verified by Blandina Valenzuela9-21-95

Data Entry Comments:

S95T001430 produced a second weight loss step of 13.66% at approximately 143°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2039 TO 2063.


TGA STD 65N8-A

16.887 mg

Rate: 10.0 °C/min

File: 00093.001 TG METTLER 20-Sep-95

Ident: 0.0 222-S Laboratory



Ramsey analysis


S95T001365 N2

9.803 mg

Rate: 10.0 °C/min

File: 00096.001 TG METTLER 20-Sep-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-145, REV. 1

WHC-SD-WM-DP-145, REV. 0

R. D. Brown 9/20/95

S95T001365DUP N2


9.895 mg

Rate: 10.0 °C/min

File: 00097.001 TG METTLER 20-Sep-95

Ident: 0.0 222-S Laboratory

2-261



WHC-SD-WM-DP-145, REV. 0
WHC-SD-WM-DP-145, REV. 1

Ramanahalli

S95T001430 N2


18.741 mg

Rate: 10.0 °C/min

File: 00094.001 TG METTLER 20-Sep-95

Ident: 0.0 222-S Laboratory

292-262



WHC-SD-WM-DP-145, REV. 0


S95T001430DUP N2

16.479 mg

Rate: 10.0 °C/min

File: 00095.001 TG METTLER 20-Sep-95

Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#**2251**Analyst: JDS Instrument: DSC01 DSC03 Book # 65N8AMethod: LA-514-114 Rev/Mod B-O

Worklist Comment: BY-108 TGA. JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.228</u>	N/A %
95000104	BY-108 (R)	2 SAMPLE	S95T001404 0		TGA-03	SOLID	<u>N/A</u>	<u>14.64</u>	
95000104	BY-108 (R)	3 DUP	S95T001404 0		TGA-03	SOLID	<u>14.64</u>	<u>14.87</u>	N/A %

Final page for worklist # 2251

Analyst Signature


Date 8-28-95

Analyst Signature

Date



Data Entry Comments:



Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-264

Curve 1: TGA


File info: TER082801 Mon Aug 28 11:34:25 1995

Sample Weight: 26.486 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-265 TO 2-267.

BEST AVAILABLE COPY



WHC-SD-WM-DP. /45, REV. 1

N2 10C/MIN
TEMP1: 26.0 °C TEMP2: 300.0 °C


TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Aug 28 11:40:09 1995

Curve 1: TGA
File info: SAM082801 Mon Aug 28 13:53:23 1995
Sample Weight: 25.609 mg
S95T001404 SAM

BEST AVAILABLE COPY



WHC-SD-WM-DP. /45, REV. 1

N2 10C/MIN
TEMP1: 500.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER

7 Series Thermal Analysis System
Thu Aug 31 09:08:00 1995


Curve 1: TGA

File info: SAM082802 Mon Aug 28 15:09:45 1995

Sample Weight: 18.061 mg

S95T001404 DUP

BEST AVAILABLE COPY



N2 100/MIN
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

JD SPELLMAN

PERKIN-ELMER

7 Series Thermal Analysis System
Thu Aug 31 09:14:14 1995

WHC-SD:V1M-DP-145, REV. L

LABCORE Data Entry Template for Worklist#

2252

Analyst: JDS Instrument: TGA01 3
 Method: LA-514-114 Rev/Mod C-O 9-11-95
8v

Worklist Comment: BY-108 rerun tga.JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.24</u>	<u>N/A</u>	%
95000104	BY-108 (R)	2 SAMPLE	S95T001396 0	TGA-03	SOLID	<u>N/A</u>	<u>8.56</u>		%
95000104	BY-108 (R)	3 DUP	S95T001396 0	TGA-03	SOLID	<u>8.56</u>	<u>9.36</u>	<u>N/A</u>	%
		4 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.09</u>	<u>N/A</u>	%
95000104	BY-108 (R)	5 SAMPLE	S95T001410 0	TGA-03	SOLID	<u>N/A</u>	<u>19.95</u>		%
95000104	BY-108 (R)	6 DUP	S95T001410 0	TGA-03	SOLID	<u>19.95</u>	<u>20.93</u>	<u>N/A</u>	%

Final page for worklist # **2252**

See attached for signatures
 Analyst Signature Date 9-11-95

Dan J. Foye 9-13-95
 Analyst Signature Date

Verified 9/14/95 JMF

S95T001396 produced a second weight loss step of 25.08% at approximately 175°C.

Data Entry Comments: S95T001410 produced a second weight loss step of 13.7% at approximately 150°C

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2252

Analyst: Jds Instrument: TGA01 Book # 6SN8A

Method: LA-560-112 Rev/Mod C-O

Worklist Comment: BY-108 rerun tga.JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID			N/A	%
95000104	BY-108 (R)	2 SAMPLE	S95T001396 1	TGA-01	SOLID	N/A			%
95000104	BY-108 (R)	3 DUP	S95T001396 1	TGA-01	SOLID			N/A	%
95000104	BY-108 (R)	4 SAMPLE	S95T001410 1	TGA-01	SOLID	N/A			%
95000104	BY-108 (R)	5 DUP	S95T001410 1	TGA-01	SOLID			N/A	%

Final page for worklist # 2252

Jds

9-8-95

Analyst Signature

Date

Analyst Signature

Date

Other instrument was
used.

9-11-95

BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-269

Curve 1: TGA

File info: TER090802 Fri Sep 8 08:42:53 1995


Sample Weight: 25.202 mg

65N8-A Terlig

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-270 TO 2-276.

BEST AVAILABLE COPY

2-270



WHC-SD-WM-DP-145, REV. 1


N2 10C/MIN
TEMP: 25.0 C TIME: 0.0 min RATE: 10.0 C/min
TEMP: 200.0 C

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Fri Sep 8 08:48:57 1995

Curve 1: TGA
File info: SAM090803 Fri Sep 8 09:51:48 1995
Sample Weight: 26.305 mg
S95T001396 SAM AT 10C/MIN

BEST AVAILABLE COPY



N2
TEMP1: 35.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 500.0 C

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 11 13:35:39 1995

WHC-SD-WM-DP-145, REV.1

Curve A: TGA


TAKE STATUS: SAN-S20040 - UNIT: STEP 40, 100°C, 1000°C

sample weight: 23.363 g

SDS: 0.01395, DLS: 0.00115

BEST AVAILABLE COPY

26-272
Hazardous Waste



TAKE STATUS: SAN-S20040 - UNIT: STEP 40, 100°C, 1000°C
sample weight: 23.363 g
SDS: 0.01395, DLS: 0.00115

Temperature (°C)

TAKE STATUS: SAN-S20040 - UNIT: STEP 40, 100°C, 1000°C
sample weight: 23.363 g
SDS: 0.01395, DLS: 0.00115
Mon Sep 11 13:41:29 1995

WHC-SD-WM-DR-145, REV. 1


Curve 1: TGA

File info: TER090702 Thu Sep 7 08:34:36 1995

Sample Weight: 27.355 mg

65N8-A Terliq

BEST AVAILABLE COPY



N2 10C/MIN
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Thu Sep 7 08:52:30 1995

Curve 1: TGA


File info: SAM090703 Thu Sep 7 10:19:48 1995

Sample Weight: 15.276 mg

S95T004791 DUP at 10C/min

1410 SAM

BEST AVAILABLE COPY




WHC-SD-WM-DP-145, REV.L

N2
TEMP1: 35.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 500.0 C
JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 11 13:47:31 1995

Curve 1: TGA
File info: SAM090704 Thu Sep 7 12:07:13 1995
Sample Weight: 22.855 mg
S95T001410 DUP AT 10C/MIN

BEST AVAILABLE COPY



WHC-SD-WM-DP-145, REV. A

N2
TEMP1: 35.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 500.0 °C

Temperature (°C)

JD SPELLMAN
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 11 13:52:51 1995

LABCORE Data Entry Template for Worklist#

2269

Analyst: JDS Instrument: TGA01 1 Book # 65N8A

Method: LA-560-112 Rev/Mod B-0

Worklist Comment: BY-108 TGA.JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.52</u>	<u>N/A</u> %
95000104	BY-108 (R)	2 SAMPLE	S95T001431 0		TGA-01	SOLID	<u>N/A</u>	<u>43.77</u>	<u></u> %
95000104	BY-108 (R)	3 DUP	S95T001431 0		TGA-01	SOLID	<u>43.77</u>	<u>44.47</u>	<u>N/A</u> %

Final page for worklist # **2269**

JDS 9-4-95
Analyst Signature Date

Jerry Hammontree 9-13-95
Analyst Signature Date

Verified 9/14/95 JDS

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-276

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-277 TO 2-279.

BEST AVAILABLE COPY


TGA STD 65N8A

22.966 mg

Rate: 10.0 °C/min

File: 00041.001 TG METTLER 04-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

S95T001432 SAM N2


28.348 mg

Rate: 10.0 °C/min

File: 00043.001 TG METTLER 04-Sep-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height-12.41 mg
-43.77 %
ResiC. 15.93 mg
56.19 %
Dpeak 93.0 °C



WHC-SD-WM-DP-145, REV. 1

BEST AVAILABLE COPY

S95T001432 DUP N2


15.786 mg

Rate: 10.0 °C/min

File: 00045.001 TG METTLER 04-Sep-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height ~7.02 mg
-44.47 %
ResiC. 8.77 mg
55.53 %
Dpeak 77.0 °C



LABCORE Data Entry Template for Worklist#

2276

Analyst: SMF

Instrument: TGA01 3

Book # 65N8A

Method: LA-514-114 Rev/Mod C-0

9-11-95

BDV

Worklist Comment: BY-108 TGA, please run under N2. JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.76</u>	<u>N/A</u> %
95000118	BY-108 (R)	2 SAMPLE	S95T001925 0		TGA-03	SOLID	<u>N/A</u>	<u>24.90</u>	<u>N/A</u> %
95000118	BY-108 (R)	3 DUP	S95T001925 0		TGA-03	SOLID	<u>24.90</u>	<u>23.78</u>	<u>N/A</u> %

Final page for worklist # **2276**

See attached for signatures
Analyst Signature Date 9-11-95

J. Jones 9-14-95
Analyst Signature Date

Verified 9/14/95 John M. Euge

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Q-280

LABCORE Data Entry Template for Worklist#

2276

Analyst: SMF Instrument: TGA01 Book # 65N8A

Method: LA-560-112 Rev/Mod
SMF 9-8-95 LA-514-114/CO

Worklist Comment: BY-108 TGA, please run under N2. JMF

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID		N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001925	0	TGA-01	SOLID	N/A		%
95000118	BY-108 (R)	3 DUP	S95T001925	0	TGA-01	SOLID		N/A	%

Final page for worklist # 2276

Smfulton/9-8-95

Analyst Signature Date

Analyst Signature Date

Other instrument was used.

9-11-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-281

Curve 1: TGA

File info: TER090801 Fri Sep 8 00:45:23 1995


Sample Weight: 16.699 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 7-282 TO 7-284.

BEST AVAILABLE COPY

WHO-SD-WM-DP-145, REV. 1



N2 10C/MIN

TEMP1: 25.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP2: 200.0 °C

Temperature (°C)

SM FULTON


PERKIN-ELMER

7 Series Thermal Analysis System

Fri Sep 8 02:04:38 1995

Curve 1: TGA
File info: SAM090801 Fri Sep 8 04:07:31 1995
Sample Weight: 17.787 mg
S95T001925 at 10C/min

BEST AVAILABLE COPY



N2
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Fri Sep 8 05:13:14 1995


Curve 1: TGA

File info: SAM090802 Fri Sep 8 06:22:09 1995

Sample Weight: 14.494 mg

S95T001925 DUP at 10C/min

2-284



WHC-SJ-WM-DP-145, REV. 1

N2
TEMP1: 25.0 C TIME1: 0.0 min RATE1: 10.0 C/min
TEMP2: 500.0 C

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Fri Sep 8 06:27:00 1995

LABCORE Data Entry Template for Worklist#

2418

Analyst: PJM Instrument: TGA0 3 Book # 65N8A

Method: LA-514-114 Rev/Mod C-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.54</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001968	0	TGA-03	SOLID	<u>N/A</u>	<u>15.4</u>		%
95000118	BY-108 (R)	3 DUP	S95T001968	0	TGA-03	SOLID	<u>15.4</u>	<u>9.88</u>	<u>N/A</u>	%
		4 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.44</u>	<u>N/A</u>	%
95000118	BY-108 (R)	5 SAMPLE	S95T001969	0	TGA-03	SOLID	<u>N/A</u>	<u>7.15</u>		%
95000118	BY-108 (R)	6 DUP	S95T001969	0	TGA-03	SOLID	<u>7.15</u>	<u>7.79</u>	<u>N/A</u>	%

Final page for worklist # 2418

See attached for signatures
Analyst Signature Date 9-30-95


Analyst Signature Date

Verified by Blandina Valenzuela
10-4-95

S95T001968 produced a second weight loss step of 20.6% at approximately 185°C

Data Entry Comments: S95T001969 produced a second weight loss step of 32.48% at approximately 190°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2418

Analyst: RL McCown

Instrument: TGA0

Book # 65N8-AMethod: LA-560-112 Rev/Mod C-O
Rev 9/30/95

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001968 0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	3 DUP	S95T001968 0	TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	4 SAMPLE	S95T001969 0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	5 DUP	S95T001969 0	TGA-01	SOLID			N/A	%

Final page for worklist # 2418RL McCown 9/30/95
Analyst Signature Date

Analyst Signature Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: TGA


File info: TER092802 Thu Sep 28 16:56:52 1995

Sample Weight: 17.614 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-287 TO 2-292

BEST AVAILABLE COPY



N2 10C/MIN
TEMP: 288.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN

PERKIN-ELMER

7 Series Thermal Analysis System
Thu Sep 28 17:28:38 1995

WHC-SJ-WM-DP-145, REV. 1

Curve 1: TGA

File info: SAM092803 Thu Sep 28 19:49:14 1995


Sample Weight: 23.963 mg

S95T001891 1968

10-2-95
BDV

Sample #
entered wrong

BEST AVAILABLE COPY



WHC-SJ-WM-DP-145, REV.1

N2
TEMP1: 35.0 S TEMP2: 500.0 S TIME1: 0.0 min RATE1: 10.0 °C/min

TEMP1: 35.0 S TEMP2: 500.0 S TIME1: 0.0 min RATE1: 10.0 °C/min

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Tue Oct 3 15:15:08 1995

Curve 1: TGA


File info: SAM092804 Thu Sep 28 21:20:07 1995

Sample Weight: 29.920 mg

10-2-45 REV
S95T001834 DUP
1968

Sample #
entered
incorrectly

BEST AVAILABLE COPY



VHCC-SDA-V:M.DP. /45/, REV. 1

N2
TEMP: 25.0 G TIME: 0.0 min RATE: 10.0 C/min

Temperature (°C)

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Tue Oct 3 15:20:15 1995

Curve 1: TGA


File info: TER093001 Sat Sep 30 18:26:00 1995

Sample Weight: 20.754 mg

65NB-A Terlid

BEST AVAILABLE COPY

R-290



WHC-SDN/M-DP-1452, REV. 1


N2 10C/MIN
TEMP: 200.0 8 TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

P.J. MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Sat Sep 30 18:28:05 1995

Curve 1: TGA
File info: SAM093001 Sat Sep 30 19:33:09 1995
Sample Weight: 25.996 mg
S95T001969 SAM

BEST AVAILABLE COPY



N2
TEMP1: 25.0 S TIME1: 0.0 min RATE1: 10.0 °C/min
TEMP2: 500.0 S

Temperature (°C)

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Tue Oct 3 15:02:01 1995

WHD-C-SD-V10-DP-145, REV. L

Curve 1: TGA


File info: SAM093002 Sat Sep 30 20:55:47 1995

Sample Weight: 21.756 mg

S95T001969 DUP

BEST AVAILABLE COPY

262-2



WIC-CD-WIN-DP-145, REV. 1

N2
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Tue Oct 3 15:07:33 1995

LABCORE Data Entry Template for Worklist#

2419

Analyst: ADP Instrument: TGA0 1 Book #: 65N8AMethod: LA-560-112 Rev/Mod BD

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>57.99</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001970 0	TGA-01	SOLID	<u>N/A</u>	<u>11.12</u>		%
95000118	BY-108 (R)	3 DUP	S95T001970 0	TGA-01	SOLID	<u>11.12</u>	<u>7.68</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001971 0	TGA-01	SOLID	<u>N/A</u>	<u>8.07</u>		%
95000118	BY-108 (R)	5 DUP	S95T001971 0	TGA-01	SOLID	<u>8.07</u>	<u>7.82</u>	<u>N/A</u>	%

Final page for worklist # 2419

 9-29-95
Analyst Signature Date 10-11-95
Analyst Signature DateVerified by Blandina Valenzuela
10-11-95S95T001970 produced a second weight loss step of 23.68% at approximately 14150°C
10/15/95Data Entry Comments: S95T001971 produced a second weight loss step of 29.22%
at approximately 150°C

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-294 TO 2-298.

BEST AVAILABLE COPY

TGA STD 65N8A N2


20.012 mg

Rate: 10.0 °C/min

File: 00103.001 TG METTLER 29-Sep-95
Ident: 0.0 222-S Laboratory

Step Analysis
Height-11.60 mg
-57.99 %
ResiC. 8.41 mg
42.01 %
Dpeak 89.2°C

2-294



Anthony Parente 9-29-95

WHD-SD-114-DP-LH, REV. 1

BEST AVAILABLE COPY

S95T001970 SAM N2


25.383 mg

Rate: 10.0 °C/min

File: 00104.001 TG METTLER 29-Sep-95

Ident: 0.0 222-S Laboratory

Step Analysis
Height -2.82 mg
-11.12 %
ResiC. 22.54 mg
88.80 %



THIS DOCUMENT IS AN OFFICIAL RECORD OF THE U.S. GOVERNMENT.


BEST AVAILABLE COPY

S95T001970 DUP N2

26.589 mg

Rate: 10.0 °C/min

File: 00105.001 TG METTLER 29-Sep-95
Ident: 0.0 222-S Laboratory



WHS-C-114-DP-1453, REV. 1
T. KEY, 1/95


BEST AVAILABLE COPY

S95T001971 SAM N2

16.826 mg

Rate: 10.0 °C/min

File: 00106.001 TG METTLER 29-Sep-95
Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY

S95T001971 DUP N2

27.763 mg

Rate: 10.0 °C/min

Step Analysis

Height -2.17 mg

-7.82 %


ResiC. 25.56 mg

92.08 %

File: 00107.001 TG METTLER 29-Sep-95

Ident: 0.0 222-S Laboratory

2.298



W.H.C. 222-S-LDP-145 REV. 1

LABCORE Data Entry Template for Worklist#**2420**Analyst: DJD Instrument: TGA0 1 Book #: 65N8-AMethod: LA-560-112 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv


GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.14</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001972	0	TGA-01	SOLID	<u>N/A</u>	<u>31.01</u>		%
95000118	BY-108 (R)	3 DUP	S95T001972	0	TGA-01	SOLID	<u>31.01</u>	<u>33.09</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001973	0	TGA-01	SOLID	<u>N/A</u>	<u>41.10</u>		%
95000118	BY-108 (R)	5 DUP	S95T001973	0	TGA-01	SOLID	<u>41.10</u>	<u>41.13</u>	<u>N/A</u>	%

Final page for worklist # 2420 Analyst Signature

10-1-95 Date

 Analyst Signature

10/02/95 Date

Verified by Blandina

Valenzuela
10-4-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-299

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2300 TO 2304

BEST AVAILABLE COPY

TGA STD 65N8A


18.840 mg

Rate: 10.0 °C/min

File: 00038.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory

Step Analysis
Height-11.33 mg
-60.14 %
ResiC. 7.51 mg
39.86 %
Dpeak 89.2 °C

00038



WMC-STD-NM-DP-145, REV. L

David C. Dunbar 10 - 95

BEST AVAILABLE COPY


S95T001972 SAM N2

27.616 mg

Rate: 10.0 °C/min

File: 00040.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-146, REV L

BEST AVAILABLE COPY


S95T001972 DUP N2

31.428 mg

Rate: 10.0 °C/min

File: 00042.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory

2-302



WIGGLED DP. /45.EV. 1


BEST AVAILABLE COPY

S95T001973 SAM N2

12.808 mg

Rate: 10.0 °C/min

File: 00045.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory



WHD-SD-VIM-DP-145, REV. 1


BEST AVAILABLE COPY

S95T001973 DUP N2

13.002 mg

Rate: 10.0 °C/min

File: 00047.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory



WHD-SD-WM-DP-145, REV. 1

LABCORE Data Entry Template for Worklist#

2421

Analyst: PJM Instrument: TGA0 3 Book #: 65N8AMethod: LA-514-114 Rev/Mod C-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.37</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001974	0	TGA-03	SOLID	<u>N/A</u>	<u>9.77</u>	%
95000118	BY-108 (R)	3 DUP	S95T001974	0	TGA-03	SOLID	<u>9.77</u>	<u>8.80</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001975	0	TGA-03	SOLID	<u>N/A</u>	<u>9.37</u>	%
95000118	BY-108 (R)	5 DUP	S95T001975	0	TGA-03	SOLID	<u>9.37</u>	<u>8.98</u>	%

Final page for worklist # 2421

See attached for signatures

Analyst Signature Date

Paul Bennett 10-5-95

Analyst Signature Date

Verified by Blandine
Valenzuela
 10-5-95

S95T001974 produced an additional weight loss step of 31.08 wt %
 at approximately 190°C.

Data Entry Comments: S95T001975 produced an additional weight loss step of 26.84
 weight %. at approximately 190°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,
 R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2421

Analyst: BJMcClain Instrument: TGA0 Book # 65N8A

Method: LA-500-114 Rev/Mod C-O
Rm 10/1/95

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID		N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001974 0		TGA-01	SOLID	N/A		%
95000118	BY-108 (R)	3 DUP	S95T001974 0		TGA-01	SOLID		N/A	%
95000118	BY-108 (R)	4 SAMPLE	S95T001975 0		TGA-01	SOLID	N/A		%
95000118	BY-108 (R)	5 DUP	S95T001975 0		TGA-01	SOLID		N/A	%

Final page for worklist # 2421

BJMcClain 10/1/95
Analyst Signature Date

Analyst Signature Date

Other instrument

was used

10-5-95

BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

Curve 1: TGA

File info: TER100101 Sun Oct 1 00:35:12 1995


Sample Weight: 14.638 mg

65NB-A Terliq

BEST AVAILABLE COPY

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 230/231 234

RJ McCown
10/1/95



N2 10C/MIN
TEMP: 300.0 °C TIME: 0.0 min RATE: 10.0 °C/min


Temperature (°C)

PJ MCCOWN
PEAKIN-ELMER
7 Series Thermal Analysis System
Sun Oct 1 00:43:03 1995

WHC-CDW-M-DP-145, REV. L

Curve 1: TGA
File info: SAM100103 Sun Oct 1 04:56:20 1995
Sample Weight: 19.263 mg
S95T001974 SAM

BEST AVAILABLE COPY



WIC-CDD-VN-DR-145, REV 1

N2
TEMP1: 35.0 8 TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)


PJ MCCOWN

PERKIN-ELMER

7 Series Thermal Analysis System
Wed Oct 4 10:06:14 1995

Curve 1: TGA
File info: SAM100104 Sun Oct 1 06:08:22 1995
Sample Weight: 24.140 mg
S95T001974 DUP

BEST AVAILABLE COPY



WHC-CD-WM-DP-145, REV. I

N2
TEMP: 25.0 G TIME: 0.0 min RATE: 10.0 °/min

Temperature (°C)

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Oct 4 10:23:51 1995


Curve 1: TGA

File info: SAM100101 Sun Oct 1 01:47:33 1995

Sample Weight: 21.038 mg

S95T001975 SAM

BEST AVAILABLE COPY



N2
TEMP: 300.0 8 TIME: 0.0 min RATE: 10.0 °C/min

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Oct 4 10:29:10 1995

WHC-CD-WM-DP-145, REV. 1


Curve 1: TGA

File info: SAM100102 Sun Oct 1 02:58:28 1995

Sample Weight: 17.844 mg

S95T001975 DUP

BEST AVAILABLE COPY



WIGGINS, D.P. / 145, REV. 1

N2
TEMP: 25.0 8 TIME: 0.0 min RATE: 10.0 °/min

Temperature (°C)

PJ MCCOWN
PERKIN-ELMER
7 Series Thermal Analysis System
Wed Oct 4 10:34:36 1995

LABCORE Data Entry Template for Worklist#**2422**Analyst: RMcClownInstrument: TGA0 1Book # 6518AMethod: LA-560-112 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.38</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001976	0	TGA-01	SOLID	<u>N/A</u>	<u>37.69</u>	%
95000118	BY-108 (R)	3 DUP	S95T001976	0	TGA-01	SOLID	<u>37.69</u>	<u>35.62</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001977	0	TGA-01	SOLID	<u>N/A</u>	<u>35.50</u>	%
95000118	BY-108 (R)	5 DUP	S95T001977	0	TGA-01	SOLID	<u>35.50</u>	<u>35.50</u>	%

Final page for worklist # **2422**RMcClown 10/1/95
Analyst Signature DateJLJ 10/12/95
Analyst Signature DateVerified by Blandina
Valenzuela 10-4-95

Data Entry Comments: S95T001976 produced a second weight loss step of 4.58% at approximately 400°C.

S95T001977 produced a second weight loss step of 5.37% at approximately 390°C

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

BEST AVAILABLE COPY

TGA STD 65N8A


22.239 mg

Rate: 10.0 °C/min

File: 00028.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJ McCown
10/1/95



TM-6-244-OP-146.FEV.1

BEST AVAILABLE COPY

S95T001976 SAM N2


28.645 mg

Rate: 10.0 °C/min

File: 00034.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

R. McClowen
10/1/95



BEST AVAILABLE COPY

S95T001976 DUP N2


37.174 mg

Rate: 10.0 °C/min


File: 00036.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJ McLean
10/1/95



BEST AVAILABLE COPY



BEST AVAILABLE COPY


S95T001977 DUP N2
Run 10/1/95

45.867 mg

Rate: 10.0 °C/min

File: 00032.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory

DJ McClain
10/1/95



LABCORE Data Entry Template for Worklist#

2424

Analyst: RJMcClain

Instrument: TGA0 1

Book # 65NBA

Method: LA-560-112 Rev/Mod B-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.20</u>	N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001905	0	TGA-01	N/A	<u>43.26</u>		%
95000118	BY-108 (R)	3 DUP	S95T001905	0	TGA-01	SOLID	<u>43.26</u>	<u>45.28</u>	N/A %
95000118	BY-108 (R)	4 SAMPLE	S95T001906	0	TGA-01	SOLID	<u>N/A</u>	<u>36.22</u>	
95000118	BY-108 (R)	5 DUP	S95T001906	0	TGA-01	SOLID	<u>36.22</u>	<u>33.78</u>	
						<u>37.3</u>	<u>33.22</u>	N/A	%
						10-2-95 PN			

Final page for worklist # **2424**

RJMcClain 10/1/95
Analyst Signature Date

Altick
Analyst Signature Date

Verified by Blandina Valenzuela
10-4-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

J-318

BEST AVAILABLE COPY

TGA STD 65N8A


19.829 mg

Rate: 10.0 °C/min

File: 00048.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJMcClown
10/1/95



BEST AVAILABLE COPY

S95T001905 SAM N2


26.156 mg

Rate: 10.0 °C/min

File: 00049.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJMcConaughay
10/1/95



WHC-SD-WM-DR-145, REV. 1

BEST AVAILABLE COPY

S95T001905 DUP N2


25.652 mg

Rate: 10.0 °C/min

File: 00050.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJMcClure
10/1/95



WHC-SD-WM-DP- 145, REV. 1

BEST AVAILABLE COPY


S95T001906 SAM N2

26.528 mg

Rate: 10.0 °C/min

File: 00051.001 TG METTLER 01-Oct-95
Ident: 0.0 222-S Laboratory

RHMcClure
10/1/95



WHC-SD-WM-DP-145, REV. L

BEST AVAILABLE COPY

S95T001906 DUP N2


30.065 mg

Rate: 10.0 °C/min

File: 00052.001 TG METTLER 01-Oct-95

Ident: 0.0 222-S Laboratory

RJ McCown
10/1/95



LABCORE Data Entry Template for Worklist#**2425**Analyst: SWF Instrument: TGA0 1 Book # 65N8AMethod: LA-560-112 Rev/Mod B-D

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.58</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001861 0	TGA-01	SOLID	<u>N/A</u>	<u>41.30</u>		%
95000118	BY-108 (R)	3 DUP	S95T001861 0	TGA-01	SOLID	<u>41.30</u>	<u>41.77</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001862 0	TGA-01	SOLID	<u>N/A</u>	<u>26.93</u>		%
95000118	BY-108 (R)	5 DUP	S95T001862 0	TGA-01	SOLID	<u>26.93</u>	<u>25.93</u>	<u>N/A</u>	%

Final page for worklist # 2425

Susie M. Fulton 9-24-95
 Analyst Signature Date
2330

R. Jones 9-25-95
 Analyst Signature Date

Verified by Blandine Valenzuela
 (9-26-95)

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2-325 TO 2-329.

BEST AVAILABLE COPY

TGA STD 65N8A

22.225 mg

Rate: 10.0 °C/min


File: 00018.001

TG METTLER

24-Sep-95

Ident: 0.0

222-S Laboratory



50.

100.

150.

Susie M. Fulton 9-29-95

2-325

W.H.C.-SD-W.M.-DP-145 REV. 1


S95T001861 N2

19.353 mg

Rate: 10.0 °C/min

File: 00020.001 TG METTLER 24-Sep-95

Ident: 0.0 222-S Laboratory




S95T001861 DUP N2

34.720 mg

Rate: 10.0 °C/min

File: 00022.001 TG METTLER 24-Sep-95

Ident: 0.0 222-S Laboratory




S95T001862 N2

52.658 mg

Rate: 10.0 °C/min

File: 00024.001 TG METTLER 24-Sep-95
Ident: 0.0 222-S Laboratory



WHC-SD-WM-DRP14/5, REV. 1


S95T001862 DUP N2

25.266 mg

Rate: 10.0 °C/min

File: 00026.001 TG METTLER 24-Sep-95

Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#**2426**Analyst: SMF Instrument: TGA0 3 Book # 65N8AMethod: LA-514-114 Rev/Mod C-O

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-03	SOLID	<u>59.74</u>	<u>60.78</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001863 0	TGA-03	SOLID	<u>N/A</u>	<u>22.2</u>		%
95000118	BY-108 (R)	3 DUP	S95T001863 0	TGA-03	SOLID	<u>22.2</u>	<u>15.6</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001915 0	TGA-03	SOLID	<u>N/A</u>	<u>13.13</u>	<u>7.25-45</u>	<u>BDV</u> %
95000118	BY-108 (R)	5 DUP	S95T001915 0	TGA-03	SOLID	<u>12.81</u>	<u>13.13</u>	<u>12.81</u>	<u>N/A</u> %

Final page for worklist # **2426**

See attached for signatures
 Analyst Signature Date 9-25-95 BDV

L.J. 9-25-95
 Analyst Signature Date

Verified by Blandina
 Valenzuela 9-26-95

S95T001863 produced one ⁹⁻²⁵⁻⁹⁵ other weight loss step of 18.09% at approximately 160°C

Data Entry Comments: S95T001915 produced one other weight loss step of 21.35% at approximately 160°C

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

LABCORE Data Entry Template for Worklist#

2426

Analyst: SME Instrument: TGA0 _____ Book # 65N8D 9-25-95
BDVMethod: LA 560 112 Rev/Med LA-514-114/CDSMF 9-24-95
Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	2 SAMPLE	S95T001863 0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	3 DUP	S95T001863 0	TGA-01	SOLID			N/A	%
95000118	BY-108 (R)	4 SAMPLE	S95T001915 0	TGA-01	SOLID	N/A			%
95000118	BY-108 (R)	5 DUP	S95T001915 0	TGA-01	SOLID			N/A	%

Final page for worklist # 2426

Susie M. Dalton 9-24-95
Analyst Signature Date 2330

Analyst Signature Date

Other instrument
was used.9-25-95
BDV

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-331

Curve 1: TGA


File info: TER092401 Sun Sep 24 16:19:30 1995

Sample Weight: 25.017 mg

65NB-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2332 TO 2336

BEST AVAILABLE COPY



N2 10C/MIN
TEMP1: 25.0 C TIME1: 0.0 min RATE1: 10.0 C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Sun Sep 24 16:21:10 1995

WHC-SD-WM-DP. 145, REV. A


Curve 1: TGA

File info: SAM092401 Sun Sep 24 17:21:55 1995

Sample Weight: 15.861 mg

S95T001863, at 10C/min

2000-333



N2

TEMP1: 300.0 C TIMES: 0.0 min RATE1: 10.0 C/min

TEMP2: 500.0 C

Temperature (°C)

SM FULTON

PERKIN-ELMER

7 Series Thermal Analysis System

Mon Sep 25 09:58:40 1995


Curve 1: TGA

File info: SAM092402 Sun Sep 24 18:49:01 1995

Sample Weight: 21.728 mg

S95T001863 DUP, at 10C/min

2-334



WHC-SD-WM-DP- 145, REV. 1

N2
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)


SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 25 10:03:44 1995

Curve 1: TGA

File info: SAM092403 Sun Sep 24 19:59:28 1995

Sample Weight: 27.899 mg

S95T001915, at 10C/min



N2
TEMP1: 25.0 8 TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 25 09:42:45 1995


WHC-SD-WM-DP-145, REV. 1

Curve 1: TGA

File info: SAM092404 Sun Sep 24 21:45:20 1995

Sample Weight: 19.110 mg

S95T001915 DUP. at 10C/min



WHC-SD-WM-DP-145, REV. 1

N2
TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON
PERKIN-ELMER
7 Series Thermal Analysis System
Mon Sep 25 09:48:53 1995

LABCORE Data Entry Template for Worklist#**2428**Analyst: ADP Instrument: TGA0 1 Book # 65N8AMethod: LA-560-112 Rev/Mod BO

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>59.98</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001918 0	TGA-01	SOLID	<u>N/A</u>	<u>29.31</u>		%
95000118	BY-108 (R)	3 DUP	S95T001918 0	TGA-01	SOLID	<u>29.31</u>	<u>35.21</u>	<u>N/A</u>	%

Final page for worklist #**2428**Ortly Parnts 9-19-95
Analyst Signature DateJean 9-20-95
Analyst Signature DateVerified by Blandina Valenzuela
9-21-95Data Entry Comments: 9/20/95
BDV

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-337

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2338 TO 2340.

BEST AVAILABLE COPY

TGA STD 65N8-A N2

20.024 mg


Rate: 10.0 °C/min

File: 00090.001 TG METTLER 19-Sep-95

Ident: 0.0 222-S Laboratory

Step Analysis
Heignt: 12.01 mg
-69.98 %
Resid.: 8.01 mg
40.02 %
Dpeak = 85.0 °C

2338



WHC-SD-WM-DP-245, REV. L


S95T001918 SAM N2

10.392 mg

Rate: 10.0 °C/min

File: 00091.001 TG METTLER 19-Sep-95

Ident: 0.0 222-S Laboratory




S95T001918 DUP N2

18.728 mg

Rate: 10.0 °C/min

File: 00092.001 TG METTLER 19-Sep-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DR. 145 REV. 1

LABCORE Data Entry Template for Worklist#

2615

Analyst: RDM Instrument: TGA0 1 Book # 65N8-AMethod: LA-560-112 Rev/Mod B-0

Worklist Comment: Please run BY-108 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD		TGA-01	SOLID	<u>59.74</u>	<u>60.53</u>	<u>N/A</u>	%
95000118	BY-108 (R)	2 SAMPLE	S95T001863 0	TGA-01	SOLID	<u>N/A</u>	<u>32.06</u>		%
95000118	BY-108 (R)	3 DUP	S95T001863 0	TGA-01	SOLID	<u>32.06</u>	<u>35.91</u>	<u>N/A</u>	%
95000118	BY-108 (R)	4 SAMPLE	S95T001916 1	TGA-01	SOLID	<u>N/A</u>	<u>36.23</u>		%
95000118	BY-108 (R)	5 DUP	S95T001916 1	TGA-01	SOLID	<u>36.23</u>	<u>27.02</u>	<u>N/A</u>	%

Final page for worklist # 2615

RDM 10/7/95
 Analyst Signature Date

RJ 10-10-95
 Analyst Signature Date

Verified by Blandina Valenzuela
10-11-95

Data Entry Comments: S95T001916 produced a second weight loss of 7.81%
at 380°C.

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-341

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 232 TO 236.


BEST AVAILABLE COPY

TGA STD 65N8-A

17.159 mg

Rate: 10.0 °C/min

File: 00027.001 TG METTLER 07-Oct-95
Ident: 0.0 222-S Laboratory



WHC-SD-WW-DP. 145, REV. 1

R. M. 10/7/95

BEST AVAILABLE COPY


S95T001863 N2

63.603 mg

Rate: 10.0 °C/min

File: 00030.001 TG METTLER 07-Oct-95

Ident: 0.0 222-S Laboratory



WHC-SD-WM-DP-245, REV.1


S95T001863DUP N2

54.215 mg

Rate: 10.0 °C/min

File: 00031.001 TG METTLER 07-Oct-95

Ident: 0.0 222-S Laboratory



W.H.C.-S/N-W.M.-D.P. /HES REV 1


BEST AVAILABLE COPY

S95T001916 N2

38.036 mg

Rate: 10.0 °C/min

File: 00028.001 TG METTLER 07-Oct-95
Ident: 0.0 222-S Laboratory




BEST AVAILABLE COPY

S95T001916DUP N2

27.437 mg

Rate: 10.0 °C/min

File: 00029.001 TG METTLER 07-Oct-95
Ident: 0.0 222-S Laboratory



LABCORE Data Entry Template for Worklist#**2685**Analyst: RJMInstrument: TGA0 1Book # FSN-BPA 8/30/95
65NB-AMethod: LA-560-112 Rev/Mod B-D

Worklist Comment: Please run BY-108 TGAs under N2. gjh

GROUP	PROJECT	S TYPE	SAMPLE#	R A -----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.49</u>	<u>N/A</u> %
95000118	BY-108 (R)	2 SAMPLE	S95T001966 0		TGA-01	SOLID	<u>N/A</u>	<u>20.84</u>	<u></u> %
95000118	BY-108 (R)	3 DUP	S95T001966 0		TGA-01	SOLID	<u>20.84</u>	<u>19.25</u>	<u>N/A</u> %

Final page for worklist #**2685**RJM

1/30/95

RJM

9/30/95

Analyst Signature

Date

Analyst Signature

Date

Verified by Blandina Valenzuela

10-4-95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

2-347

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 2348 TO 2352.


BEST AVAILABLE COPY

TGA STD 65N8-A

16.865 mg

Rate: 10.0 °C/min

File: 00022.001 TG METTLER 30-Sep-95
Ident: 0.0 222-S Laboratory



R. D. Miller 9/30/95

BEST AVAILABLE COPY


S95T001966 N2

27.742 mg

Rate: 10.0 °C/min

File: 00025.001 TG METTLER 30-Sep-95

Ident: 0.0 222-S Laboratory



BEST AVAILABLE COPY


S95T001966 DUP N2

34.288 mg

Rate: 10.0 °C/min

File: 00026.001 TG METTLER 30-Sep-95

Ident: 0.0 222-S Laboratory



WMC-SD-WM-DP-445, REV. 1