



Westinghouse
Hanford Company

P.O. Box 1970 Richland, WA 99352

012192 21
C 98-5



222-S/RCRA ANALYTICAL LABORATORIES

3
-
0
PROJECT: SINGLE-SHELL TANK WASTE

5
-
TANK: 241-U-110

2
-
CORE: 5

1
SEGMENT: 4

6
CUSTOMER ID. NUMBER: 89-041

BEST AVAILABLE COPY

REPORT REVISION: 2

DATE PRINTED: AUGUST 27, 1990

Table of Contents

Signature Form.....	ii
Introduction.....	1
Sampling and Custody Data.....	4
Sample Data Summary.....	7
Physical Test Results.....	9
Extrusion of Segment.....	10
Appearance of Solid.....	10
Homogenization.....	10
Single Shell Tank Segment - Subsamples.....	11
VOA Sampling.....	11
Photograph of Tank 241-U-110 Core 5 Segment 4.....	12
Particle Size Analysis.....	13
Undigested Sample Analysis Results.....	18
Untreated Sample Results Summary.....	19
pH Analysis.....	20
Percent Water.....	21
KOH Fusion Results.....	23
Fusion Analysis Results Summary.....	24
Fusion Analysis Results of Laboratory Solution.....	24.1
Fusion Dissolution I.....	25
Fusion Dissolution II.....	26
Total Alpha and Beta Analysis.....	27
GEA Analysis	35
Uranium Analysis.....	70
Water Digestion Test Analysis.....	71
Water Digestion Results Summary.....	72
Water Digestion Analysis of Laboratory Solutions.....	73
Water Digestion I.....	74
Water Digestion II.....	74.1
Ion Chromatographic Analysis.....	75
Total Organic Carbon Analysis - Not Acidified.....	103
Total Organic Carbon Analysis - Acidified.....	111
Acid Digestion Test Results.....	118
Acid Digestion Data Summary.....	119
Acid Digestion Analysis.....	120
Inductively Coupled Plasma Spectrographic Analysis.....	121
ICP Raw Data Summary.....	122

I have reviewed this report and certify that the package is in compliance with Quality Assurance Project Plan - WHC-SD-CP-QAPP-002. I found it to be a true and accurate accounting both technically and for completeness of the laboratory analyses performed on this sample.

Shirley A. Cervantes
Shirley A. Cervantes
Data Coordinator

Date August 13, 1990

Cary M. Seidel
Cary M. Seidel
Unit Manager

Date August 13, 1990

Larry H. Taylor
Larry H. Taylor
Laboratory Q.A. Officer

Date August 30, 1990

INTRODUCTION

6
6
9
5
7
2
1
1
9

INTRODUCTION

Westinghouse Hanford Company Analytical Laboratories are supporting the characterization efforts of the single shell tanks. The characterization of tank 241-U-110 was performed under Phase 1A and 1B of the "Waste Characterization Plan for the Hanford Site Single-Shelled Tanks" (WHC-EP-0210).

Tank 241-U-110 has a 500,000 gallon capacity. Construction was completed in 1944. The tank received first cycle waste, REDOX high-level waste, coating waste, and laboratory waste until 1975. Between July 7, 1975 and February 2, 1976, P-10 pumps were installed, and 41,700 gallons of liquid waste were pumped from the tank. Tank 241-U-110 still contains an estimated 195,000 gallons of waste.

The Analytical Laboratories performs all analytical analysis to the specifications of the "Quality Assurance Project Plan", WHC-SD-CP-QAPP-002. In accordance with WHC-SD-CP-QAPP-002 the following laboratory policies are being followed. Spikes are performed on either the undissolved sample, or the sample after dissolution as directed by the chemist. If the spike addition is found to be less than 20% of an analyte concentration, the spike recovery is not reported due to errors introduced by the precision of the sample analysis. The concentration of spike additions will be re-evaluated before the start of phase 1C. Two spiking routines are being used during phase 1A and 1B. For the following analyses, Ion Chromatography, Inductively Coupled Plasma, Mercury Hydride, Total Organic Carbon, and Carbonate analyses the solid sample is spiked independently from the sample digestion. Any non-homogeneity of the sample could adversely affect the spike recoveries. For the radio-isotopic analysis and other analyses not specified above the spikes were performed by spiking an aliquot of sample after digestion.

The laboratory does not report sample results from batch analyses that are questionable. The results from questionable batches are discarded and the analysis is repeated. Sample cards (laboratory travelers) for the repeated analysis are reissued for analysis after they have been stamped "rerun". Laboratory travelers are issued using a computerized routine according to a "sample point". This sample point label (segment-n) on the laboratory travelers and on the GEA analysis reports has no relationship to the sampling activities or the sample identification. All results in this data package relate only to the sample identified as segment 4 from core 5 taken from tank 241-U-110.

The organic analysis of this sample will be performed by Pacific Northwest Laboratories (PNL). Due to instrument and procedure problems, PNL has been unable to separate organics from the normal paraffin hydrocarbon present in the samples. The results from the organic analysis will be provided when available.

Samples analyzed for Total Organic Carbon between November 1, 1989 and February 22, 1990 were not acidified. The results from these analyses include total organic carbon, carbonate, and dissolved carbon dioxide from the air. The validity of these analyses are subject to interpretation. The Total Organic Carbon procedure was corrected and these analyses were repeated and

results are reported in this package.

All sample results reported here by weight are reported as the "wet weight" of the sample. Some samples did noticeably lose moisture during the process of aliquoting and weighing the sample for digestion. The percent moisture was determined at the earliest opportunity to minimize any errors introduced by the loss of moisture. Drying samples before analysis resulted in radiation exposure increases of about a factor of ten. In order to reduce and control radiation exposure to laboratory personnel, the samples were not dried before aliquouting and digestion. This policy may result in some laboratory results being biased high.

This report is formatted into sections corresponding to the type of dissolutions performed prior to analysis. A brief summary of analytical results is reported, followed by calibration data and an analysis batch report. Any notable observations regarding an analysis are noted on the batch report for that analysis. Copies of laboratory travelers can be found in Appendix A.

6 1 4 0 6 5 1 2 1 1 6

SAMPLING AND CUSTODY DATA

9 1 1 2 3 5 7 0 5 7 0

CHAIN-OF-CUSTODY RECORD FOR CORE SAMPLING

(1) Shipment Number 5-022-89 (2) Sample Number 89-041 (3) Supervisor Brian Z Hall
 (4) Tank 241-U-110 (5) Riser 19 (6) Segment 4 (7) Cask Serial Number 1009C

Radiation Survey Data:	(8) FIELD	(20) LABORATORY	(9) Shipment Description:
Over Top Dose Rate	<u>.5 mR/hr</u>	<u>NA</u>	A. Work Package Number <u>2W-89-009571W</u>
Side Dose Rate	<u>7 mR/hr</u>	<u>10 mR/hr</u>	B. Cask Seal Number <u>For Future Use</u>
Bottom Dose Rate	<u>2.5 mR/hr</u>	<u>3.5 mR/hr</u>	C. Sampler Number Used <u>48</u>
Smearable Contamination	<u>1 Dst</u> (alpha)	<u>< Dst</u> (alpha)	D. Date and Time Sampler Unseated <u>11-7-89 2:26 PM</u>
	<u>1 Dst</u> (beta-gamma)	<u>< Dst</u> (beta-gamma)	E. Expected Liquid Content <u>0-20%</u>
	RPT <u>B. Hall</u> (Signature)	RPT <u>Robert</u> (Signature)	F. Expected Solid Content <u>80-100%</u>
			G. Dose Rate Through Drill String <u>100 mR/hr</u>
			H. Expected Sample Length <u>19"</u>

(10) INFORMATION (Include statement of laboratory tests to be performed.*)

Characterize segment according to WHC-EP-0210

*Reference laboratory work request, if available.

Comments: Core #5

(11) POINT OF ORIGIN <u>241-U-110</u> <u>200W Area</u>	(12) SENDER NAME <u>Brian L. Hall</u> SENDER SIGNATURE <u>Brian Z Hall</u>	(13) DATE AND TIME RELEASED <u>11-7-89</u> <u>10:00 AM</u>	(14) DESTINATION <u>222S Lab</u> <u>200W Area</u>	(15) RECIPIENT NAME <u>Vida Boyle</u> RECIPIENT SIGNATURE <u>Vida Boyle</u>	(17) DATE AND TIME RECEIVED <u>11/7/89</u> <u>1040</u>
--	---	--	---	--	--

(15) Seal Intact Upon Release? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	(18) Seal Intact Upon Receipt? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	(19) Seal Data Consistent with this Record? Shipment No. <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No			Sample No. <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
---	---	--	--	--	---

**Single Shell Tank Waste Characterization
Summary of Core Sample**

**Phase
I-A**

Tank ID:	241-U-110
Riser ID:	19
Core ID:	5

Date Sampling Initiated:	9/19/89
Date Sampling Completed:	11/7/89

Segment 1	Lab Serial No. F0003	Segment 8	Lab Serial No.
	Customer ID. No. 89-038		Customer ID. No.
	Last Segment? No		Last Segment?
Segment 2	Lab Serial No. F0027	Segment 9	Lab Serial No.
	Customer ID. No. 89-039		Customer ID. No.
	Last Segment? No		Last Segment?
Segment 3	Lab Serial No. F5001	Segment 10	Lab Serial No.
	Customer ID. No. 89-040		Customer ID. No.
	Last Segment? No		Last Segment?
Segment 4	Lab Serial No. F5033	Segment 11	Lab Serial No.
	Customer ID. No. 89-041		Customer ID. No.
	Last Segment? Yes		Last Segment?
Segment 5	Lab Serial No.	Segment 12	Lab Serial No.
	Customer ID. No.		Customer ID. No.
	Last Segment?		Last Segment?
Segment 6	Lab Serial No.	Segment 13	Lab Serial No.
	Customer ID. No.		Customer ID. No.
	Last Segment?		Last Segment?
Segment 7	Lab Serial No.	Segment 14	Lab Serial No.
	Customer ID. No.		Customer ID. No.
	Last Segment?		Last Segment?

Interim

Rev. B 3/27/90
SST-1

Prepared by:

Signature

H. S. Rich

Printed Name

Date: 4/2/90

Verified by:

Signature

C. M. Seidel

Printed Name

Date: 4/2/90

Approved by:

Signature

L. H. Taylor

Printed Name

Date: 5/9/90

SAMPLE DATA SUMMARY

SUMMARY DATA REPORT

Tank 241-U-110
 Core 5
 Segment 4
 Customer Id. 89-041

Untreated	Sample	Acid Digestion Results			
		Sample	Duplicate	Sample	Duplicate
pH	12.71	12.92	Aluminum	34730 ug/g	37858 ug/g
Percent Water	38.80%	39.10%	Barium	22 ug/g	16 ug/g
			Bismuth	20942 ug/g	19470 ug/g
			Boron	LT	LT
			Cadmium	LT	LT
			Calcium	155 ug/g	152 ug/g
			Chromium	446 ug/g	429 ug/g
			Copper	27 ug/g	15 ug/g
			Iron	7927 ug/g	8145 ug/g
			Lanthanum	LT	LT
			Lead	131 ug/g	LT
			Lithium	4 ug/g	LT
			Magnesium	103 ug/g	116 ug/g
			Manganese	1993 ug/g	2048 ug/g
			Molybdenum	20 ug/g	11 ug/g
			Nickel	59 ug/g	44 ug/g
			Phosphorous	9243 ug/g	9534 ug/g
			Potassium	198 ug/g	LT
			Silver	LT	LT
			Sodium	65496 ug/g	65413 ug/g
			Strontium	243 ug/g	245 ug/g
			Tantalum	63 ug/g	13 ug/g
			Tin	31 ug/g	17 ug/g
			Titanium	11 ug/g	LT
			Zinc	57 ug/g	LT
			Zirconium	44 ug/g	LT

LT: Less Than

Total Organic Carbon/Carbonate

4710 ug/g 5030 ug/g

Samples were not acidified before analysis

Total Organic Carbon (Second analysis)

859 ug/g 1100 ug/g

Samples were acidified before analysis

4
0
0
5
5
2
1
9

PHYSICAL TEST RESULTS

**Single Shell Tank
Extrusion of Segment -- Physical Tests**

**Phase
I-A**

Lab Segment Serial No.: F5033

Customer ID: 89-041

Analyst: Rich L. Weiss

Date Extruded: 11-07-89

Drainable Liquid Liquid Submitted for Segment Analysis? -- No

Gross <10 ml	Tare	Net
Serial	Date/Time _____ / _____	Estimated
Specific	Calculated	

Appearance of Liquid:

N/A

Dimensions of Segment

Complete Segment Obtained? No	Length: 16 inches	Calculated Volume: 12.57 cubic in
Remarks None		

Appearance of Solid:

Bottom 4" very soft, runny, light to medium brown; Middle, firmer, medium brown; Top 5" rubbery firm (very cohesive), black in color. Small hard chunks (possibly crystals) of approximately 1/8" in diameter.

Penetrometer

3.8

lbs/sq in

Remarks: None

Homogenization

Procedure: T038A-00712 Revision: F	Quantity of Material	not recorded	grams
Date Homogenized: 11-09-89	Time Homogenized:	5.0	Minutes
Operator: John R. Smith			

Laboratory Notebook Reference

WHC-N-313-4

Notebook No.

3

Page No.

Prepared by: R. L. Rich H.S. Rich Date: 04-02-90
 Signature Printed Name

Verified by: C.M. Seidel C.M. Seidel Date: 04-02-90
 Signature Printed Name

Approved by: L. H. Taylor L.H. Taylor Date: 5/9/90
 Signature Printed Name

Interim

3/27/90

Rev C

SST-3

15
14
13
12
11
10
9
8
7
6
5
4
3
2
1
0

10

**Single Shell Tank
Segment -- Subsamples**

**Phase
1A**

Customer ID: 89-041

Lab Segment Serial No. F5033

Volatile Organic Analysis

VOA Sample	Laboratory Serial Number: 89-041-30	Date Sampled: 11-07-89
------------	-------------------------------------	------------------------

Sample shipped to PNL for analysis

Particle Size Distribution Analysis

Particle Size Sample	Laboratory Serial Number: F5033	Date Sampled: 11-07-89
----------------------	---------------------------------	------------------------

Sample analysis performed at 222-S

Homogenized Solids

Undigested Solids Analysis

Laboratory Serial Number for Sample: F5033	Date Sampled: 11/09/89
--	------------------------

Laboratory Serial Number of Duplicate Sample: F5034

Fusion Analysis of Solids

Laboratory Serial Number for Sample: F5038	Date Sampled: 11/09/89
--	------------------------

Laboratory Serial Number of Duplicate Sample: F5039

Laboratory Serial Number of Spiked Sample: N/A
--

Acid Digestion Analysis of Solids

Laboratory Serial Number for Sample: F5048	Date Sampled: 11/09/89
--	------------------------

Laboratory Serial Number of Duplicate Sample: F5049

Laboratory Serial Number of Spiked Sample: F5050
--

Water Digestion Analysis of Solids

Laboratory Serial Number for Sample: F5043	Date Sampled: 11/09/89
--	------------------------

Laboratory Serial Number of Duplicate Sample: F5044

Laboratory Serial Number of Spiked Sample: F5045
--

Laboratory Notebook Reference


WHC-W-313-4

Notebook No.

3

Page No.


Prepared by:



H. S. Rich
Printed Name

Date: 04/04/90


Verified by:



C. M. Seidel
Printed Name


Date: 04-06-90

Approved by:



L. H. Taylor
Printed Name

Date: 5/9/90



TANK 241-U-110, CORE 5, SEGMENT 4

9 1 1 2 1 5 3 0 6 2 7

Brinkmann
Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS
VIA BRINKMANN 2010
STATISTICS

SAMPLE NAME : SST,B000033,F5033,H2O,DBB/SBK
FILE NAME : A:\SST\F5033.002

DATE	:	15/11/1989	ACQ. RANGE	: 0.5-150	COUNTS	: 86935
TIME	:	11:40	ACQ. MODE	: SAMPLE	S.N.F.	: 0.85
CONFIG.	:	1 (0.7 S1)	ACQ. TIME	: 410 SEC	S.D.U.	: 4559
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.8E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 8.9E-03 %

MEAN Diameter S.D.

Number, Length	:	1.17 μ m	1.14 μ m
Number, Area	:	1.63 μ m	1.23 μ m
Number, Volume	:	2.92 μ m	2.09 μ m
Length, Area	:	2.28 μ m	4.01 μ m
Length, Volume	:	4.61 μ m	4.64 μ m
Area, Volume	:	9.35 μ m	14.65 μ m
Volume, Moment	:	32.30 μ m	25.66 μ m

	MEDIAN Diameter	MODE	CONFIDENCE	
Number	:	0.88 μ m	0.75 μ m	100.00%
Area	:	3.94 μ m	4.75 μ m	86.06%
Volume	:	27.98 μ m	37.84 μ m	99.36%

particle diameter < 150 μ m 100 11/15/89

SAMPLE NAME : SST,B000033,F5033,H2O,DBB/SBK
FILE NAME : A:\SST\F5033.002

DATE	:	15/11/1989	ACQ. RANGE	:	0.5-150	COUNTS	:	86935
TIME	:	11:40	ACQ. MODE	:	SAMPLE	S.N.F.	:	0.85
CONFIG.	:	1 (0.7 S1)	ACQ. TIME	:	410 SEC	S.D.U.	:	4559
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	:	4	CONCENTR.:	:	6.8E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	:	95.00%(V)	SOLIDS	:	0.9E-03 %

PROBABILITY NUMBER DENSITY GRAPH

Name: SST,B000033,F5033,H2O,DBB/SBK

6.8E+06 #/ml(99.9%)

Mode at 0.75 μ m


<< SCALE RANGE (μ m): 0 - 10 >>

Local Median : 0.80 μ m

Local Mean(n1): 1.14 μ m

Local S.D.(n1): 0.08 μ m

Local Conf(n1): 100.00 %



SAMPLE NAME : SST,B000033,F5033,H2O,DBB/SDK
FILE NAME : A:\SST\F5033.002


DATE	:	15/11/1989	: ACQ. RANGE	:	0.5-150	: COUNTS	:	86935
TIME	:	11:40	: ACQ. MODE	:	SAMPLE	: S.N.F.	:	0.85
CONFIG.	:	1 (0.7 S1)	: ACQ. TIME	:	410 SEC	: S.D.U.	:	4559
CELL TYPE	:	MAGNETIC (3)	: SAMPLE SIZE	:	4	: CONCENTR.:	:	6.8E+06 #/ml
SAMPLE TYPE	:	REGULAR	: REQ. CONF.	:	95.00%(V)	: SOLIDS	:	8.9E-03 %

PROBABILITY NUMBER DISTRIBUTION GRAPH

Name: SST,B000033,F5033,H2O,DBB/SDK
6.8E+06 #/ml(99.9%)

Local Median : 0.88 μ m
Local Mean(nl): 1.14 μ m
Local S.D.(nl): 0.88 μ m
Local Conf(nl): 100.00 %

<< SCALE RANGE (μ m): 0 - 10 >>




SAMPLE NAME : SST,B000033,F5033,H2O,DBB/SBK
FILE NAME : A:\SST\F5033.002

DATE	:	15/11/1989	ACQ. RANGE	:	0.5-150	COUNTS	:	86935
TIME	:	11:40	ACQ. MODE	:	SAMPLE	S.N.F.	:	0.85
CONFIG.	:	1 (0.7 S1)	ACQ. TIME	:	410 SEC	S.D.U.	:	4559
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	:	4	CONCENTR.	:	6.8E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	:	95.00%(V)	SOLIDS	:	8.9E-03 %

PROBABILITY VOLUME DENSITY GRAPH

Name: SST,B000033,F5033,H2O,DBB/SBK
8.9E-05 cc/ml(100.0%) Mean(v_m): 2.92 μ m Median : 27.98 μ m
Mode at 5.00 μ m S.D.(v_m): 2.09 μ m Mean(v_m): 32.30 μ m
<< SCALE RANGE (μ m): ADJUSTED >> S.D.(v_m): 25.66 μ m
Conf(v_m): 99.36 %




SAMPLE NAME : SST,B000033,F5033,H2O,DBB/SBK
FILE NAME : A:\SST\F5033.002

DATE	:	15/11/1989	ACQ. RANGE	:	0.5-150	COUNTS	:	86935
TIME	:	11:40	ACQ. MODE	:	SAMPLE	S.N.F.	:	0.05
CONFIG.	:	1 (0.7 S1)	ACQ. TIME	:	410 SEC	S.D.U.	:	4559
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	:	4	CONCENTR.:	:	6.8E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	:	95.00%(V)	SOLIDS	:	8.9E-03 %

PROBABILITY VOLUME DISTRIBUTION GRAPH

Name: SST,B000033,F5033,H2O,DBB/SBK
0.9E-05 cc/ml(100.0%) Median : 27.98 μ m
Mean(μ v): 2.92 μ m Mean(μ m): 32.30 μ m
S.D.(μ v): 2.09 μ m S.D.(μ m): 25.66 μ m
(μ m): ADJUSTED Conf(μ m): 99.36 %



UNDIGESTED SAMPLE ANALYSIS

9 1 1 2 3 5 9 0 6 7 4

61

Single Shell Tank Project

Untreated Sample Results

4/03/90

Tank: 241-U-110
 Core: 5
 Segment: 4
 Customer ID: 89-041

	Check Standard	Blank	Sample	Sample Duplicate	Check Standard
Laboratory ID:	F5032	F5053	F5033	F5034	F5036
pH	100.30%	5.67	12.71	12.92	100.00%
%Water	99.10%	8.8 mg	38.80%	39.10%	95.00%

Approved by: H.S. Rich H.S. Rich

Date: 04/01/90

Verified by: C.M. Seidel C.M. Seidel

Date: 5/2/90

Approved by: L.H. Taylor L.H. Taylor

Date: 5/9/90

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

pH of the solid sample

Instrument	AL10653
Procedure / Rev	LA-212-103/A-3
Technologist	6C269 M. Franz
Date	11/22/89
Temperature	24.7
Starting Time	13:00
Ending Time	15:30
Chemist	R. E. Brandt

	Description	Lab. Id.
1	Initial Check Standard	F5032
2	Blank	F5053
3	Sample of 89-041	F5033
4	Duplicate of 89-041	F5034
5	Final Check Standard	F5036
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume of Standard
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Standard	72C11/5.0 ml			5.0 ml

Interim

Rev. E 4/04/90

Prepared by: R. E. Brandt
Signature

H. S. Rich
Printed Name

Date: 04/06/90

Verified by: C. M. Seidel
Signature

C. M. Seidel
Printed Name

Date: 04/06/90

Approved by: L. H. Taylor
Signature

L. H. Taylor
Printed Name

Date: 5/9/90

SST-102

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

% Water in Sample 89-041

Instrument	N/A
Procedure / Rev	LA-564-101/D-1
Technologist	6B598 R.D. Hale
Date	11/28/89
Temperature	120 C
Starting Time	10:00 11/27/89
Ending Time	10:00 11/28/89
Chemist	R.E. Brandt

	Description	Lab. Id.
1	Initial Check Standard	F5032
2	Reagent Blank	F5053
3	Sample 89-041	F5033
4	Duplicate of 89-041	F5034
5	Final Check Standard	F5036
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	11C11AG/1.0g.			1.0 gram

Prepared by: R.D. Hale
Signature

H. S. Rich
Printed Name

Date: May 4, 1990

Verified by: C. M. Seidel
Signature

C. M. Seidel
Printed Name

Date: May 4, 1990

Approved by: L.H. Taylor
Signature

L.H. Taylor
Printed Name

Date: 5/9/90

Interim

Rev.E

4/04/90

SST-102

This page was intentionally left blank.

7
7
0
6
0
9
5
8
2
1
1
2
1
9

KOH FUSION ANALYSIS

Single Shell Tank Project

Fusion Analysis
Laboratory Results of Solids
Units are Sample Wet Weight

Tank: 241-U-110
 Core: 5
 Segment: 4
 Customer ID: 89-041

Laboratory ID:	Check Standard F5037	Blank F5052	Sample F5038	Sample Duplicate F5039	Spike of Sample F5040	Check Standard F5041
Fusion Digestion (First Digestion)			3.03 g/L	3.4 g/L		
Total Alpha	95.10%	<2.42E-04 uci/L	3.63E-01 uci/g	4.38E-01 uci/g	*	86.40%
Total Beta	95.40%	<5.58E-04 uci/L	1.08E+03 uci/g	1.14E+03 uci/g	*	94.80%
Uranium	108.30%	0.359 ug/L	2.16E+03 ug/g	2.89E+03 ug/g	113.30%	98.60%
Fusion Digestion (Second Digestion)			1.92 g/L	1.73 g/L		
GEA						
Cs-137	102.80%	<6.32E-02 uci/L	4.84E+01 uci/g	4.29E+01 uci/g	106.00%	102.10%

* Ratio of Standard To Spike Insufficient To Calculate Spike Recovery.

9 1 1 2 3 5 9 0 6 1 0

24.1

Single Shell Tank Project

Fusion Analysis
Sample Results on Laboratory Digestions

Tank: 241-U-110
 Core: 5
 Segment: 4
 Customer ID: 89-041

Laboratory ID:	Check	Blank	Sample	Sample	Spike of	Check
	Standard F5037	F5052	F5038	Duplicate F5039	Sample F5040	Standard F5041
Fusion Digestion (First Digestion)			3.03 g/L	3.4 g/L		
Total Alpha	95.10%	<2.42E-04 uci/L	1.10E+00 uci/L	1.49E+00 uci/L	*	86.40%
Total Beta	95.40%	<5.58E-04 uci/L	3.28E+03 uci/L	3.88E+03 uci/L	*	94.80%
Uranium	108.30%	0.359 ug/L	6.55E-03 ug/L	9.82E-03 ug/L	113.30%	98.60%
Fusion Digestion (Second Digestion)			1.92 g/L	1.73 g/L		
GEA						
Cs-137	102.80%	<6.32E-02 uci/L	9.29E+01 uci/L	7.42E+01 uci/L	106.00%	102.10%

* Ratio of Standard to Sample was insufficient to calculate the spike recovery.

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Fusion Dissolution

Instrument	N/A
Procedure / Rev	LA-549-141/A-0
Technologist	6B598 R. D. Hale
Date	11/22/89
Temperature	450 C
Starting Time	10:30
Ending Time	2:30
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Reagent Blank	F5052
2	Sample 89-041	F5038
3	Duplicate 89-041	F5039
4		
5		
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim	Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
	N/A				

Rev E 4/04/90	Prepared by: <u>S. A. Catlow</u> Signature	H. S. Rich Printed Name	Date: 04/06/90
SST-102	Verified by: <u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 04/06/90
	Approved by: <u>L. H. Taylor</u> Signature	L. H. Taylor Printed Name	Date: 04/06/90

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Second Fusion Dissolution

Instrument	N/A
Procedure / Rev	LA-549-141/A-0
Technologist	6B598 R. D. Hale
Date	4-30-90
Temperature	450 C
Starting Time	08:30
Ending Time	10:00
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Reagent Blank	F5052
2	Sample 89-041	F5038
3	Duplicate of 89-041	F5039
4		
5		
6		
7		
8		
9		
10		
11		



	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
N/A				

Interim

Rev E
4/04/90

SST-102

Prepared by:		H. S. Rich Printed Name	Date: May 4, 1990
Verified by:		C. M. Seidel Printed Name	Date: May 4, 1990
Approved by:		L. H. Taylor Printed Name	Date: 5/9/90

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	WA57277
Procedure / Rev	LA-548-101/A-2
Technologist	6C269 M. Franz
Date	11/29/89
Temperature	N/A
Starting Time	07:30
Ending Time	14:30
Chemist	S. A. Catlow

Total Alpha + Total Beta
 Fusion Dissolution
 Detector 16
 Samples were prepared in batch, but
 counted randomly.

	Description	Lab. Id.
1	Initial Check Standard	F5037
2	Reagent Blank	F5052
3	Sample 89-041	F5038
4	Duplicate of 89-041	F5039
5	Spike of 89-041	F5040
6	Final Check Standard	F5041
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	83B44/10 mL			10 mL
Spike	83B44/10 mL	Sample/5 μ L		10.005 mL

Interim

4/04/90

Rev E

SST-102

Prepared by: H. S. Rich
Signature

H. S. Rich
Printed Name

Date: May 4, 1990

Verified by: C. M. Seidel
Signature

C. M. Seidel
Printed Name

Date: May 4, 1990

Approved by: L. H. Taylor
Signature

L. H. Taylor
Printed Name

Date: 5/9/90

Phase
I-ASingle Shell Tank
Calibration RecordAnalyte: Am²⁴¹; Co⁶⁰

Procedure LQ-508-002

Revision: A-0

Instrument: Detector #16

Property Number: WA63999

Technologist: R. A. Jones

Payroll Number: 65801

Date: 2-27-88

Calibration Standard ID: 36B40A4; 36B40B4; 36B40C4

Analyte Concentration: N/A

Type of Calibration: Efficiency

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3			
4	SEE ATTACHED SHEETS		
5			
6			
7			
8			
9			
10			

Interim

Comments:

Rev. (Draft) 1/18/89

Prepared by:

Signature

H. S. Rich

Printed Name

Date: May 4, 1990

Verified by:

Signature

C. M. Seidel

Printed Name

Date: May 4, 1990

SST-103

Approved by:

Signature

L. H. Taylor

Printed Name

Date: May 9, 1990

CALIBRATION SHEET FOR ALPHA/BETA SYSTEMS: USING PROCEDURE LQ-508-002

DETECTOR No. 16

RADIOMUCLIDE: Co-60
 HALF LIFE: 1925
 COUNT TIME: 5
 CPM BKG: 8
 CPM 1" BKG: 8

TIME ZERO DATE (HD): 15883
 DATE COUNTED 5" (HD) 15954
 DATE COUNTED 2" (HD) 16298
 DATE COUNTED 1" (HD) 15917

CALIBRATED BY: RA JONES HD 0 = 09/25/44

STANDARD ID	SIZE	DATE	TIME	COUNTS @ 0 DEG.	COUNTS @ 90 DEG.	COUNTS @ 180 DEG.	COUNTS @ 270 DEG.
32B40A8	1	04/24/88	1200	113109	111016	110133	111447
32B40B7	1	04/24/88	1200	204475	206864	207083	206310
32B40C7	1	04/24/88	1200	300283	298335	294777	296162
32B40A4	2	05/31/88	1200	90965	89136	88796	88938
32B40B3	2	05/31/88	1200	168475	166598	166421	168782
32B40C4	2	05/31/88	1200	254942	252194	253139	256756
32B40A5	5	06/17/88	1200	81186	80053	81855	80106
32B40B6	5	06/17/88	1200	158490	161830	157312	163104
32B40C5	5	06/17/88	1200	228479	233088	233428	231241

STANDARD ID	SIZE	STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
32B40A8	1"	69550	22272	1.01	22547	0.3242
32B40B7	1"	134700	41224	1.01	41731	0.3098
32B40C7	1"	201000	59465	1.01	60197	0.2995
AVERAGE, 1" =		0.3112 +/- 0.95%	0.0243	7.81 %	ON	04/24/88
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
32B40A3	2"	71060	17873	1.16	20753	0.2921
32B40B3	2"	135500	33495	1.16	38893	0.2870
32B40C4	2"	200800	50833	1.16	59025	0.2940
AVERAGE, 2" =		0.2910 +/- 0.95%	0.0070	2.41 %	ON	05/10/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
32B40A5	5"	70160	16141	1.03	16559	0.2360
32B40B6	5"	135700	32018	1.03	32847	0.2421
32B40C5	5"	201900	46293	1.03	47492	0.2352
AVERAGE, 5" =		0.2378 +/- 0.95%	0.0073	3.08 %	ON	05/31/88
NEW EFFS FOR DET		16 Co-60	1" =	0.3112	2" =	0.2910
				5" =		0.2378

CALIBRATION SHEET FOR ALPHA/BETA SYSTEMS: USING PROCEDURE LQ-508-002

DETECTOR No. 16 5" STD TIME ZERO DATE (HD) : 13367
 RADIONUCLIDE: Am-241 TIME ZERO DATE (HD) : 15883
 HALF LIFE: 157800 DATE COUNTED 5" (HD) : 15860
 COUNT TIME: 5 DATE COUNTED 2" (HD) : 16298
 CPM BKG: 0.5 DATE COUNTED 1" (HD) : 15917
 CPM 1" BKG: 0.5 CALIBRATED BY: RA JONES HD 0 = 09/25/44

STANDARD ID	SIZE	DATE	TIME	COUNTS @ 0 DEG.	COUNTS @ 90 DEG.	COUNTS @ 180 DEG.	COUNTS @ 270 DEG.
36B40A7	1	04/25/88	1200	88249	88650	88197	88524
36B40A8	1	04/25/88	1200	87657	86696	87396	87208
36B40B7	1	04/25/88	1200	153949	153455	153519	153879
36B40B8	1	04/25/88	1200	159185	159612	159046	158672
36B40C7	1	04/25/88	1200	224092	223920	224469	225263
36B40C8	1	04/25/88	1200	219897	221702	221080	220950
36B40A4	2	05/10/89	1200	66063	66703	66300	66314
36B40B4	2	05/10/89	1200	122898	122457	122769	122929
36B40C4	2	05/10/89	1200	178443	179133	178587	177410
43B16N	5	02/27/88	1200	67302	66084	69151	70103
43B16O	5	02/27/88	1200	112955	117083	116239	116860
43B16Q	5	02/27/88	1200	185670	190640	192034	187717

7 0 6 1 7 0 1 2 9 6

STANDARD ID	SIZE	STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
36B40A7	1"	62710	17681	1.00	17683	0.2820
36B40A8	1"	60570	17447	1.00	17450	0.2881
36B40B7	1"	109900	30740	1.00	30744	0.2797
36B40B8	1"	109800	31825	1.00	31830	0.2899
36B40C7	1"	159700	44887	1.00	44893	0.2811
36B40C8	1"	159800	44181	1.00	44188	0.2765
AVERAGE, 1" =		0.2829 +/- @95%	0.0100	3.54 %	ON	04/24/88
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
36B40A4	2"	61210	13268	1.00	13292	0.2172
36B40B4	2"	111700	24552	1.00	24596	0.2202
36B40C4	2"	160600	35678	1.00	35743	0.2226
AVERAGE, 2" =		0.2200 +/- @95%	0.0053	2.41 %	ON	05/10/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
43B16N	5"	62090	13632	1.01	13782	0.2220
43B16O	5"	106400	23157	1.01	23412	0.2200
43B16Q	5"	181900	37803	1.01	38219	0.2101
AVERAGE, 5" =		0.2174 +/- @95%	0.0125	5.74 %	ON	02/27/88
NEW EFFS FOR DET		16 Am-241	1" =	0.2829	2" =	0.2200
				5" =		0.2174

ALPHA BETA EFF. DETECTOR #16

ISOTOPE	1"	2"	5"	DATE1"	DATE2"	DATE5"
AM241	.2861	.2211	.2107	022290	022290	022290
C060	.2991	.2877	.2502	022290	022290	022290
SK85	.	.				
Y83	.	.				
SR89	.	.				
SR90	.	.				
Y90	.	.				
SRY90	.	.				
TC99	.	.				
RU103	.	.				
RURH106	.	.				
S8125	.	.				
CS134	.	.				
CS137	.4316	.		032489		
CEPR144	.	.				
PM147	.2117	.		032489		
SM151	.	.				
SN113	.	.				
ZRNB95	.	.				

BIA

ALPHA BETA EFF. DETECTOR #16

ISOTOPE	1"	2"	5"	DATE1"	DATE2"	DATES"
AM241	.2829	.2200	.2174	042688	051289	022988
C060	.3112	.2910	.2378	060188	051289	061688
SR85	-	-	-			
Y88	-	-	-			
SR89	-	-	-			
SR90	-	-	-			
Y90	-	-	-			
SRY90	-	-	-			
TC99	-	-	-			
RU103	-	-	-			
RURH106	-	-	-			
SB125	-	-	-			
CS134	-	-	-			
CS137	.4316	-	-	032489		
CEPR144	-	-	-			
PM147	.2117	-	-	032489		
SM151	-	-	-			
SN113	-	-	-			
ZRN895	-	-	-			

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	401934
Procedure / Rev	LA-548-121/C-2
Technologist	6B598 R. D. Hale
Date	5/2/90
Temperature	N/A
Starting Time	13:00
Ending Time	14:00
Chemist	S. A. Catlow

GEA Analysis

Second Fusion Dissolution

Samples are prepared in batch, but counted randomly.

	Description	Lab. Id.
1	Initial Check Standard	F5037
2	Reagent Blank	F5052
3	Sample of 89-041	F5038
4	Duplicate of 89-041	F5039
5	Spike 89-041	F5040
6	Ending Check Standard	F5041
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	122B44/500uL			22 mL
Spike	122B44/500uL	Sample/100uL		22 mL

Prepared by:	<u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 8-10-90
Verified by:	<u>Cary M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 8-10-90
Approved by:	<u>R. H. Pugh</u> Signature	K. H. Taylor Printed Name	Date: 8-10-90

Interim

4/04/90

Rev. E

SST-102

**Single Shell Tank
Calibration Record**

**Phase
I-A**

Analyte: Mixed Isotope Standards

Procedure	LQ-508-003	Revision:	A-0
Instrument:	GEA Detector #3	Property Number:	WA-77228
Technologist:	J. L. Anderson	Payroll Number:	61413

Date: 02/07/89

Calibration Standard ID: 56B40D1

Analyte Concentration: N/A

Type of Calibration: Gamma Energy Analysis (Efficiency)

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3	See Attached		
4			
5			
6			
7			
8			
9			
10			

Comments:

Interim

Rev. (Draft) 1/18/89

Prepared by: H. S. Rich H. S. Rich Printed Name Date: May 4, 1990

Verified by: C. M. Seidel C. M. Seidel Printed Name Date: May 4, 1990

Approved by: L. H. Taylor L. H. Taylor Printed Name Date: May 9, 1990

SST-103

DETECTOR: 3
 GEOMETRY CODE: 41
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 1
 CALIBRATION DATE: 2-Jul-89
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL
 STANDARD ID: 56B40 D1

ENERGY (KEV)	EFFICIENCY (COUNTS/GAMMA)
59.536	2.833765E-02
88.032	2.881764E-02
122.0614	2.756557E-02
165.853	2.270614E-02
279.1967	
391.668	1.285730E-02
513.99	
661.65	7.841011E-03
898.021	5.779292E-03
1173.237	4.773005E-03
1332.501	4.278530E-03
1836.129	3.371238E-03

EQUATION 0-165 KEV

$$\begin{aligned} \text{LOG(EFF)} = & -1.113845E+01 \\ & + 3.484260E+00 * \text{LOG(ENERGY)} \\ & + -3.990659E-01 * \text{LOG(ENERGY)}^2 \end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned} \text{LOG(EFF)} = & -2.052334E+01 \\ & + 9.121738E+00 * \text{LOG(ENERGY)} \\ & + -1.553578E+00 * \text{LOG(ENERGY)}^2 \\ & + 8.018036E-02 * \text{LOG(ENERGY)}^3 \end{aligned}$$

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

DETECTOR: 3
 GEOMETRY CODE: 42
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 2
 CALIBRATION DATE: 2-Jul-89
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL
 STANDARD ID: 56B40 D1

ENERGY (KEV)	EFFICIENCY (COUNTS/GAMMA)
59.536	7.455306E-03
88.032	7.462748E-03
122.0614	7.578302E-03
165.853	6.965814E-03
279.1967	
391.668	3.596591E-03
513.99	
661.65	2.318396E-03
898.021	1.824191E-03
1173.237	1.461179E-03
1332.501	1.321243E-03
1836.129	1.011332E-03

EQUATION 0-165 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -6.838496\text{E+00} \\ & + 8.819509\text{E-01} * \text{LOG(ENERGY)} \\ & + -9.970528\text{E-02} * \text{LOG(ENERGY)}^2\end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned}\text{LOG(EFF)} = & 3.082260\text{E-01} \\ & + -1.410839\text{E+00} * \text{LOG(ENERGY)} \\ & + 1.042898\text{E-01} * \text{LOG(ENERGY)}^2 \\ & + -5.874725\text{E-03} * \text{LOG(ENERGY)}^3\end{aligned}$$

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

DETECTOR: 3
 GEOMETRY CODE: 43
 GEOMETRY DESCRIPTION: 22 ML LIQUID, POS 3
 CALIBRATION DATE: 2-Jul-89
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL
 STANDARD ID: 56B40 D1

ENERGY (KEV)	EFFICIENCY (COUNTS/GAMMA)
59.536	2.020462E-03
88.032	1.924344E-03
122.0614	2.027231E-03
165.853	1.712371E-03
279.1967	
391.668	1.056509E-03
513.99	
661.65	7.115743E-04
898.021	5.243928E-04
1173.237	4.551585E-04
1332.501	4.223636E-04
1836.129	3.139091E-04

EQUATION 0-165 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -5.300788\text{E+00} \\ & + -3.550643\text{E-01} * \text{LOG(ENERGY)} \\ & + 3.272635\text{E-02} * \text{LOG(ENERGY)}^2\end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -9.815549\text{E+00} \\ & + 2.402920\text{E+00} * \text{LOG(ENERGY)} \\ & + -4.428877\text{E-01} * \text{LOG(ENERGY)}^2 \\ & + 2.059131\text{E-02} * \text{LOG(ENERGY)}^3\end{aligned}$$

*
*
* G A M M A S P E C T R U M A N A L Y S I S *
*
* *

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:32:12

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0
DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41
SPECTRUM SIZE: 4096 CHANNELS
ORDER OF SMOOTHING FUNCTION: 5
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK
PEAK CONFIDENCE FACTOR: 85.0%
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV
ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

LTD CALCULATION PERFORMED

MULTIPLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3615

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5037-6530

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 5.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 16:49:10

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3008. SECONDS

DEAD TIME: 0.27 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89

EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:32:12

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1C	59.59	30.38	1.34	1659.	387.	14.9	
2C	63.96	32.56	1.34	1598.	2579.	4.8	CE-144
3C	73.00	37.08	1.34	1471.	616.	9.9	I-129
4C	1126.58	563.34	1.45	689.	428.	9.2	CS-134, EU-152
5C	1138.93	569.51	1.45	731.	851.	7.5	CS-134, BI-207
6	1209.70	604.88	1.60	778.	5902.	1.5	CS-134
7	1323.58	661.79	1.69	569.	11089.	1.0	CS-137
8C	1591.91	795.90	1.66	533.	3952.	2.4	CS-134
9C	1604.02	801.96	1.66	497.	357.	12.0	CS-134
10	2346.60	1173.25	1.98	405.	5529.	1.5	CO-60
11	2665.04	1332.54	1.92	115.	5058.	1.5	CO-60
12	2730.71	1365.39	1.82	45.	120.	13.4	CS-134
13	2921.99	1461.09	2.01	32.	677.	4.1	K-40

ERROR QUOTATION AT 1.00 SIGMA
 PEAK CONFIDENCE LEVEL AT 85.0%

M - MULTIPLET ANALYSIS CONVERGED NORMALLY

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:32:12

SAMPLE: F5037-6530

DATA COLLECTED ON 2-MAY-90 AT 16:49:10

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

R A D I O N U C L I D E A N A L Y S I S R E P O R T

NUCLIDE ACTIVITY CONCENTRATION IN $\mu\text{Ci}/\text{L}$
DECAY

MEASURED ERROR CORRECTED ERROR

AC-228	LLD<8.95E-01		LLD<8.95E-01	
AG-108M	LLD<1.81E-01		LLD<1.81E-01	
AG-110M	LLD<8.65E-01		LLD<8.65E-01	
AM-241	LLD<1.77E-01		LLD<1.77E-01	
AM-243	LLD<9.61E-02		LLD<9.61E-02	
AR-41	LLD<1.83E-01		LLD<1.83E-01	
AU-198	LLD<1.65E-01		LLD<1.65E-01	
BA-133	LLD<2.15E-01		LLD<2.15E-01	
BA-139	LLD<4.61E-01		LLD<4.61E-01	
BA-140	LLD<7.33E-01		LLD<7.33E-01	
BA-141	LLD<5.00E-01		LLD<5.00E-01	
BE-7	LLD<1.76E+00		LLD<1.76E+00	
BI-207	LLD<1.93E-01		LLD<1.93E-01	
BI-212	LLD<2.76E+00		LLD<2.76E+00	
BI-214	LLD<1.32E+00		LLD<1.32E+00	
CD-109	LLD<1.92E+00		LLD<1.92E+00	
CE-139	LLD<1.04E-01		LLD<1.04E-01	
CE-141	LLD<1.66E-01		LLD<1.66E-01	
CEPR144	LLD<1.42E+00		LLD<1.42E+00	
CO-56	LLD<1.96E-01		LLD<1.96E-01	
CO-57	LLD<8.91E-02		LLD<8.91E-02	
CO-58	LLD<2.00E-01		LLD<2.00E-01	
CA-60	2.13E+01	+3.39E-01	2.13E+01	+3.39E-01
CR-51	LLD<1.25E+00		LLD<1.25E+00	
CS-134	1.25E+01	+3.15E-01	1.25E+01	+3.15E-01
CS-136	LLD<1.98E-01		LLD<1.98E-01	
CS-137	2.99E+01	+3.70E-01	2.99E+01	+3.70E-01
CS-138	LLD<2.61E-01		LLD<2.61E-01	
EU-152	LLD<9.68E-01		LLD<9.68E-01	
EU-154	LLD<4.34E-01		LLD<4.34E-01	
EU-155	LLD<3.17E-01		LLD<3.17E-01	
FE-59	LLD<5.07E-01		LLD<5.07E-01	
HF-181	LLD<2.18E-01		LLD<2.18E-01	
HG-203	LLD<1.40E-01		LLD<1.40E-01	
I-131	LLD<1.78E-01		LLD<1.78E-01	
I-132	LLD<3.32E-01		LLD<3.32E-01	
I-133	LLD<1.91E-01		LLD<1.91E-01	
I-134	LLD<3.07E-01		LLD<3.07E-01	
I-135	LLD<4.89E-01		LLD<4.89E-01	
K-40	2.85E+01	+1.20E+00	2.85E+01	+1.20E+00
KR-85	LLD<4.11E+01		LLD<4.11E+01	
KR-85M	LLD<1.27E-01		LLD<1.27E-01	
KR-87	LLD<3.96E-01		LLD<3.96E-01	
KR-89	LLD<5.47E+00		LLD<5.47E+00	
LA-140	LLD<1.15E-01		LLD<1.15E-01	
LA-142	LLD<4.30E-01		LLD<4.30E-01	
MN-54	LLD<2.18E-01		LLD<2.18E-01	

MN-56	LLD<2.22E-01	LLD<2.22E-01
NA-22	LLD<1.68E-01	LLD<1.68E-01
NA-24	LLD<1.20E-01	LLD<1.20E-01
NB-94	LLD<1.84E-01	LLD<1.84E-01
NB-95	LLD<1.88E-01	LLD<1.88E-01
NB-97	LLD<1.05E+00	LLD<1.05E+00
NP-238	LLD<9.51E-01	LLD<9.51E-01
NP-239	LLD<7.97E-01	LLD<7.97E-01
PA-233	LLD<3.29E-01	LLD<3.29E-01
PA-234M	LLD<4.32E+01	LLD<4.32E+01
PB-210	LLD<4.53E+00	LLD<4.53E+00
PB-212	LLD<2.70E-01	LLD<2.70E-01
PB-214	LLD<3.89E-01	LLD<3.89E-01
PO-210	LLD<1.56E+04	LLD<1.56E+04
PO-214	LLD<6.00E+03	LLD<6.00E+03
PO-216	LLD<1.18E+04	LLD<1.18E+04
PU-239	LLD<1.22E+03	LLD<1.22E+03
PU-241	LLD<4.42E+04	LLD<4.42E+04
RA-224	LLD<2.85E+00	LLD<2.85E+00
RA-226	LLD<2.97E+00	LLD<2.97E+00
RB-88	LLD<1.23E+00	LLD<1.23E+00
RB-89	LLD<1.17E+00	LLD<1.17E+00
RN-220	LLD<1.59E+02	LLD<1.59E+02
RU-103	LLD<1.85E-01	LLD<1.85E-01
RWRH106	LLD<3.52E+00	LLD<3.52E+00
SB-124	LLD<3.79E-01	LLD<3.79E-01
SB-125	LLD<1.30E+00	LLD<1.30E+00
SC-46	LLD<2.49E-01	LLD<2.49E-01
SE-75	LLD<1.94E-01	LLD<1.94E-01
SN-113	LLD<2.40E-01	LLD<2.40E-01
SR-85	LLD<1.80E-01	LLD<1.80E-01
SR-91	LLD<2.96E-01	LLD<2.96E-01
SR-92	LLD<1.79E-01	LLD<1.79E-01
TA-182	LLD<7.30E-01	LLD<7.30E-01
TC-99M	LLD<9.31E-02	LLD<9.31E-02
TE-123M	LLD<1.00E-01	LLD<1.00E-01
TE-125M	LLD<2.50E+01	LLD<2.50E+01
TE-132	LLD<1.29E-01	LLD<1.29E-01
TH-228	LLD<5.15E+00	LLD<5.15E+00
TL-208	LLD<1.99E-01	LLD<1.99E-01
U-235	LLD<1.80E-01	LLD<1.80E-01
U-237	LLD<5.22E-01	LLD<5.22E-01
W-187	LLD<6.08E-01	LLD<6.08E-01
XE-131M	LLD<4.39E+00	LLD<4.39E+00
XE-133	LLD<1.72E-01	LLD<1.72E-01
XE-133M	LLD<1.10E+00	LLD<1.10E+00
XE-135	LLD<1.27E-01	LLD<1.27E-01
XE-138	LLD<9.73E-01	LLD<9.73E-01
Y-88	LLD<1.17E-01	LLD<1.17E-01
Y-91	LLD<6.81E+01	LLD<6.81E+01
Y-91M	LLD<2.24E-01	LLD<2.24E-01
ZN-65	LLD<5.35E-01	LLD<5.35E-01
ZR-95	LLD<3.28E-01	LLD<3.28E-01
ZR-97	LLD<1.78E-01	LLD<1.78E-01
<hr/>		
TOTAL	9.23E+01 +-1.34E+00	9.23E+01 +-1.34E+00

EBAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 2.15E-09 UC/LI

TOTAL MEASURED ACTIVITY = 9.23E+01 (+-1.34E+00) UC/LI

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
59.59	30.38	387.	14.9	6.34E+00
63.96	32.56	2579.	4.8	4.01E+01
73.00	37.08	616.	9.9	8.81E+00
1126.58	563.34	428.	9.2	1.57E+01
1138.93	569.51	851.	7.5	3.15E+01
1604.02	801.96	357.	12.0	1.81E+01
2730.71	1365.39	120.	13.4	9.55E+00

*
*
* G A M M A S P E C T R U M A N A L Y S I S *
*
* *

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:36:32

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0

DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41

SPECTRUM SIZE: 4096 CHANNELS

ORDER OF SMOOTHING FUNCTION: 5

NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK

PEAK CONFIDENCE FACTOR: 85.0%

IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV

ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

LLD CALCULATION PERFORMED

MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3616

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5052

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E+00

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 17:59:04

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3003. SECONDS

DEAD TIME: 0.10 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89

EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:36:32

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	704.00	352.21	1.37	218.	119.	20.3	PB-214
2	1021.91	511.04	1.73	145.	129.	16.5	RN-222, I-133, TL-208, NA-22, ZN-65, RH-106
3	1166.98	583.53	1.92	119.	110.	17.6	TL-208
4	1218.83	609.44	1.57	95.	161.	12.2	BI-214, RU-103
5	1323.96	661.98	1.40	102.	89.	20.4	CS-137
6	1822.73	911.29	1.66	72.	82.	20.1	
7	2921.84	1461.02	2.09	16.	657.	4.0	K-40

ERROR QUOTATION AT 1.00 SIGMA
PEAK CONFIDENCE LEVEL AT 85.0%

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:36:32

SAMPLE: F5052

DATA COLLECTED ON 2-MAY-90 AT 17:59:04
DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

RADIONUCLIDE ANALYSIS REPORT

NUCLIDE ACTIVITY CONCENTRATION IN uCi/LI

DECAY

MEASURED ERROR CORRECTED ERROR

AC-228	LLD<1.50E-01		LLD<1.50E-01	
AG-108M	LLD<3.49E-02		LLD<3.49E-02	
AG-110M	LLD<5.72E-02		LLD<5.72E-02	
AM-241	LLD<5.04E-02		LLD<5.04E-02	
AM-243	LLD<3.04E-02		LLD<3.04E-02	
AR-41	LLD<6.76E-02		LLD<6.76E-02	
AU-198	LLD<3.17E-02		LLD<3.17E-02	
BA-133	LLD<5.48E-02		LLD<5.48E-02	
BA-139	LLD<1.47E-01		LLD<1.47E-01	
BA-140	LLD<1.45E-01		LLD<1.45E-01	
BA-141	LLD<1.42E-01		LLD<1.42E-01	
BE-7	LLD<3.45E-01		LLD<3.45E-01	
BI-207	LLD<3.61E-02		LLD<3.61E-02	
BI-212	LLD<6.39E-01		LLD<6.39E-01	
BI-214	3.83E-01	+-4.68E-02	3.83E-01	+-4.68E-02
CD-109	LLD<6.11E-01		LLD<6.11E-01	
CE-139	LLD<3.32E-02		LLD<3.32E-02	
CE-141	LLD<5.49E-02		LLD<5.49E-02	
CEPR144	LLD<4.61E-01		LLD<4.61E-01	
CO-56	LLD<4.42E-02		LLD<4.42E-02	
CO-57	LLD<2.94E-02		LLD<2.94E-02	
CO-58	LLD<4.79E-02		LLD<4.79E-02	
CO-60	LLD<5.12E-02		LLD<5.12E-02	
CR-51	LLD<2.96E-01		LLD<2.96E-01	
CS-134	LLD<5.57E-02		LLD<5.57E-02	
CS-136	LLD<4.60E-02		LLD<4.60E-02	
CS-137	1.20E-01	+-2.46E-02	1.20E-01	+-2.46E-02
CS-138	LLD<1.06E-01		LLD<1.06E-01	
EU-152	LLD<2.75E-01		LLD<2.75E-01	
EU-154	LLD<1.52E-01		LLD<1.52E-01	
EU-155	LLD<1.06E-01		LLD<1.06E-01	
FE-59	LLD<1.05E-01		LLD<1.05E-01	
HF-181	LLD<3.92E-02		LLD<3.92E-02	
HG-203	LLD<3.60E-02		LLD<3.60E-02	
I-131	LLD<3.68E-02		LLD<3.68E-02	
I-132	LLD<4.25E-02		LLD<4.25E-02	
I-133	LLD<3.60E-02		LLD<3.60E-02	
I-134	LLD<6.60E-02		LLD<6.60E-02	
I-135	LLD<2.35E-01		LLD<2.35E-01	
K-40	1.39E+01	+-5.74E-01	1.39E+01	+-5.74E-01
KR-85	LLD<1.05E+01		LLD<1.05E+01	
KR-85M	LLD<4.11E-02		LLD<4.11E-02	
KR-87	LLD<8.08E-02		LLD<8.08E-02	
KR-89	LLD<1.46E+00		LLD<1.46E+00	
LA-140	LLD<4.30E-02		LLD<4.30E-02	
LA-142	LLD<8.92E-02		LLD<8.92E-02	
MN-54	LLD<4.60E-02		LLD<4.60E-02	

MN-56	LLD<4.98E-02	LLD<4.98E-02		
NA-22	LLD<5.90E-02	LLD<5.90E-02		
NA-24	LLD<5.61E-02	LLD<5.61E-02		
NB-94	LLD<4.21E-02	LLD<4.21E-02		
NB-95	LLD<4.19E-02	LLD<4.19E-02		
NB-97	LLD<6.94E-02	LLD<6.94E-02		
NP-238	LLD<1.67E-01	LLD<1.67E-01		
NP-239	LLD<2.15E-01	LLD<2.15E-01		
PA-233	LLD<8.11E-02	LLD<8.11E-02		
PA-234M	LLD<7.26E+00	LLD<7.26E+00		
PB-210	LLD<9.15E-01	LLD<9.15E-01		
PB-212	LLD<7.15E-02	LLD<7.15E-02		
PB-214	LLD<8.71E-02	LLD<8.71E-02		
PO-210	LLD<3.90E+03	LLD<3.90E+03		
PO-214	LLD<4.46E+02	LLD<4.46E+02		
PO-216	LLD<2.34E+03	LLD<2.34E+03		
PU-239	LLD<4.06E+02	LLD<4.06E+02		
PU-241	LLD<1.40E+04	LLD<1.40E+04		
RA-224	LLD<7.74E-01	LLD<7.74E-01		
RA-226	LLD<8.18E-01	LLD<8.18E-01		
RB-88	LLD<4.00E-01	LLD<4.00E-01		
RB-89	LLD<2.03E-01	LLD<2.03E-01		
RN-220	LLD<3.23E+01	LLD<3.23E+01		
RU-103	LLD<3.71E-02	LLD<3.71E-02		
RURH106	LLD<7.30E-01	LLD<7.30E-01		
SB-124	LLD<3.75E-02	LLD<3.75E-02		
SB-125	LLD<4.06E-01	LLD<4.06E-01		
SC-46	LLD<5.69E-02	LLD<5.69E-02		
SE-75	LLD<4.96E-02	LLD<4.96E-02		
SN-113	LLD<4.72E-02	LLD<4.72E-02		
SR-85	LLD<4.59E-02	LLD<4.59E-02		
SR-91	LLD<6.24E-02	LLD<6.24E-02		
SR-92	LLD<7.00E-02	LLD<7.00E-02		
TA-182	LLD<1.69E-01	LLD<1.69E-01		
TC-99M	LLD<3.10E-02	LLD<3.10E-02		
TE-123M	LLD<3.25E-02	LLD<3.25E-02		
TE-125M	LLD<8.51E+00	LLD<8.51E+00		
TE-132	LLD<3.47E-02	LLD<3.47E-02		
TH-228	LLD<1.67E+00	LLD<1.67E+00		
TL-208	1.35E-01	+-2.37E-02	1.35E-01	+-2.37E-02
U-235	LLD<5.01E-02	LLD<5.01E-02		
U-237	LLD<1.53E-01	LLD<1.53E-01		
W-187	LLD<1.39E-01	LLD<1.39E-01		
XE-131M	LLD<1.41E+00	LLD<1.41E+00		
XE-133	LLD<5.40E-02	LLD<5.40E-02		
XE-133M	LLD<2.94E-01	LLD<2.94E-01		
XE-135	LLD<3.43E-02	LLD<3.43E-02		
XE-138	LLD<2.48E-01	LLD<2.48E-01		
Y-88	LLD<3.80E-02	LLD<3.80E-02		
Y-91	LLD<2.46E+01	LLD<2.46E+01		
Y-91M	LLD<4.72E-02	LLD<4.72E-02		
ZN-65	LLD<1.36E-01	LLD<1.36E-01		
ZR-95	LLD<7.71E-02	LLD<7.71E-02		
ZR-97	LLD<4.36E-02	LLD<4.36E-02		
<hr/>				
TOTAL	1.45E+01	+-5.77E-01	1.45E+01	+-5.77E-01

EBAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 3.64E-13 UC/LI

TOTAL MEASURED ACTIVITY = 1.45E+01 (+-5.77E-01) UC/LI

% TECH. SPEC. = ***** (+****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
704.00	352.21	119.	20.3	2.88E+00
1021.91	511.04	129.	16.5	4.33E+00
1822.73	911.29	82.	20.1	4.61E+00

*
*
* G A M M A S P E C T R U M A N A L Y S I S *
*
* *

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:40:00

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0

DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41

SPECTRUM SIZE: 4096 CHANNELS

ORDER OF SMOOTHING FUNCTION: 5

NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK

PEAK CONFIDENCE FACTOR: 85.0%

IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV

ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

LTD CALCULATION PERFORMED

MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3617

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5038

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 18:59:52

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3004. SECONDS

DEAD TIME: 0.13 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89

EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:40:00

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1C	64.04	32.60	1.27	2136.	1267.	8.9	CE-144
2C	72.97	37.06	1.27	1830.	292.	17.1	I-129
3	1218.88	609.46	1.17	159.	161.	14.1	BI-214, RU-103
4	1323.57	661.78	1.65	179.	6970.	1.2	CS-137
5	2921.82	1461.01	1.79	27.	631.	4.2	K-40

ERROR QUOTATION AT 1.00 SIGMA
PEAK CONFIDENCE LEVEL AT 85.0%

C - MULTIPLET ANALYSIS CONVERGED NORMALLY

6
5
4
3
2
1
0
-1
-2
-3

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:40:00

SAMPLE: F5038

DATA COLLECTED ON 2-MAY-90 AT 18:59:52
DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

RADIOMUCLIDE ANALYSIS REPORT

NUCLIDE ACTIVITY CONCENTRATION IN uCi/LI
DECAY

MEASURED ERROR CORRECTED ERROR

AC-228	LLD<1.69E+00		LLD<1.69E+00	
AG-108M	LLD<5.52E-01		LLD<5.52E-01	
AG-110M	LLD<3.39E+00		LLD<3.39E+00	
AM-241	LLD<8.83E-01		LLD<8.83E-01	
AM-243	LLD<4.56E-01		LLD<4.56E-01	
AR-41	LLD<5.57E-01		LLD<5.57E-01	
AU-198	LLD<5.21E-01		LLD<5.21E-01	
BA-133	LLD<6.94E-01		LLD<6.94E-01	
BA-139	LLD<1.82E+00		LLD<1.82E+00	
BA-140	LLD<1.90E+00		LLD<1.90E+00	
BA-141	LLD<1.87E+00		LLD<1.87E+00	
Br-7	LLD<5.26E+00		LLD<5.26E+00	
BI-207	LLD<4.81E-01		LLD<4.81E-01	
BI-212	LLD<6.64E+00		LLD<6.64E+00	
BI-214	3.83E+00	+-5.41E-01	3.83E+00	+-5.41E-01
CD-109	LLD<8.47E+00		LLD<8.47E+00	
CE-139	LLD<4.12E-01		LLD<4.12E-01	
CE-141	LLD<6.69E-01		LLD<6.69E-01	
CEPR144	LLD<5.79E+00		LLD<5.79E+00	
CO-56	LLD<4.50E-01		LLD<4.50E-01	
Co-57	LLD<3.66E-01		LLD<3.66E-01	
CO-58	LLD<4.15E-01		LLD<4.15E-01	
Cr-60	LLD<5.51E-01		LLD<5.51E-01	
Cr-51	LLD<3.83E+00		LLD<3.83E+00	
CS-134	LLD<5.43E-01		LLD<5.43E-01	
CS-136	LLD<4.09E-01		LLD<4.09E-01	
CS-137	9.40E+01	+-1.34E+00	9.40E+01	+-1.34E+00
CS-138	LLD<9.58E-01		LLD<9.58E-01	
EU-152	LLD<2.75E+00		LLD<2.75E+00	
EU-154	LLD<1.75E+00		LLD<1.75E+00	
EU-155	LLD<1.41E+00		LLD<1.41E+00	
FE-59	LLD<1.05E+00		LLD<1.05E+00	
HF-181	LLD<6.55E-01		LLD<6.55E-01	
HG-203	LLD<4.72E-01		LLD<4.72E-01	
I-131	LLD<5.36E-01		LLD<5.36E-01	
I-132	LLD<1.18E+00		LLD<1.18E+00	
I-133	LLD<5.21E-01		LLD<5.21E-01	
I-134	LLD<7.03E-01		LLD<7.03E-01	
I-135	LLD<2.09E+00		LLD<2.09E+00	
K-40	1.33E+02	+-5.77E+00	1.33E+02	+-5.77E+00
KR-85	LLD<1.33E+02		LLD<1.33E+02	
KR-85M	LLD<5.15E-01		LLD<5.15E-01	
KR-87	LLD<1.22E+00		LLD<1.22E+00	
KR-89	LLD<1.91E+01		LLD<1.91E+01	
LA-140	LLD<5.08E-01		LLD<5.08E-01	
LA-142	LLD<1.18E+00		LLD<1.18E+00	
MN-54	LLD<4.52E-01		LLD<4.52E-01	

MN-56	LLD<5.07E-01	LLD<5.07E-01
NA-22	LLD<6.32E-01	LLD<6.32E-01
NA-24	LLD<5.04E-01	LLD<5.04E-01
NB-94	LLD<4.21E-01	LLD<4.21E-01
NB-95	LLD<4.63E-01	LLD<4.63E-01
NB-97	LLD<4.11E+00	LLD<4.11E+00
NP-238	LLD<1.76E+00	LLD<1.76E+00
NP-239	LLD<2.66E+00	LLD<2.66E+00
PA-233	LLD<1.08E+00	LLD<1.08E+00
PA-234M	LLD<8.92E+01	LLD<8.92E+01
PB-210	LLD<1.53E+01	LLD<1.53E+01
PB-212	LLD<9.61E-01	LLD<9.61E-01
PB-214	LLD<1.21E+00	LLD<1.21E+00
PO-210	LLD<3.67E+04	LLD<3.67E+04
PO-214	LLD<4.75E+03	LLD<4.75E+03
PO-216	LLD<2.61E+04	LLD<2.61E+04
PU-239	LLD<5.12E+03	LLD<5.12E+03
PU-241	LLD<1.71E+05	LLD<1.71E+05
RA-224	LLD<1.01E+01	LLD<1.01E+01
RA-226	LLD<1.12E+01	LLD<1.12E+01
RB-88	LLD<4.00E+00	LLD<4.00E+00
RB-89	LLD<2.53E+00	LLD<2.53E+00
RN-220	LLD<4.11E+02	LLD<4.11E+02
RU-103	LLD<5.46E-01	LLD<5.46E-01
RURH106	LLD<8.90E+00	LLD<8.90E+00
SB-124	LLD<5.15E-01	LLD<5.15E-01
SB-125	LLD<4.87E+00	LLD<4.87E+00
SC-46	LLD<5.84E-01	LLD<5.84E-01
SE-75	LLD<6.28E-01	LLD<6.28E-01
SN-113	LLD<7.58E-01	LLD<7.58E-01
SR-85	LLD<5.82E-01	LLD<5.82E-01
SR-91	LLD<8.35E-01	LLD<8.35E-01
SR-92	LLD<6.68E-01	LLD<6.68E-01
TA-182	LLD<1.83E+00	LLD<1.83E+00
TC-99M	LLD<3.80E-01	LLD<3.80E-01
TE-123M	LLD<4.03E-01	LLD<4.03E-01
TE-125M	LLD<1.11E+02	LLD<1.11E+02
TE-132	LLD<4.48E-01	LLD<4.48E-01
TH-228	LLD<2.39E+01	LLD<2.39E+01
TL-208	LLD<5.23E-01	LLD<5.23E-01
U-235	LLD<6.85E-01	LLD<6.85E-01
U-237	LLD<1.83E+00	LLD<1.83E+00
W-187	LLD<1.56E+00	LLD<1.56E+00
XE-131M	LLD<1.76E+01	LLD<1.76E+01
XE-133	LLD<8.16E-01	LLD<8.16E-01
XE-133M	LLD<3.89E+00	LLD<3.89E+00
XE-135	LLD<4.38E-01	LLD<4.38E-01
XE-138	LLD<3.24E+00	LLD<3.24E+00
Y-88	LLD<3.80E-01	LLD<3.80E-01
Y-91	LLD<1.88E+02	LLD<1.88E+02
Y-91M	LLD<6.31E-01	LLD<6.31E-01
ZN-65	LLD<1.25E+00	LLD<1.25E+00
ZR-95	LLD<7.90E-01	LLD<7.90E-01
ZR-97	LLD<4.05E-01	LLD<4.05E-01

TOTAL 2.31E+02 +-5.95E+00 2.31E+02 +-5.95E+00

E BAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 2.81E-08 UC/LI

TOTAL MEASURED ACTIVITY = 2.31E+02 (+-5.95E+00) UC/LI

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
64.04	32.60	1267.	8.9	1.97E+01
72.97	37.06	292.	17.1	4.18E+00

*
* G A M M A S P E C T R U M A N A L Y S I S *
*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:44:48

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0
DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41
SPECTRUM SIZE: 4096 CHANNELS
ORDER OF SMOOTHING FUNCTION: 5
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK
PEAK CONFIDENCE FACTOR: 85.0%
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV
ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

LLD CALCULATION PERFORMED
MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3618
ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5039
GEOMETRY DESCRIPTION:
SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E-01
STANDARD SIZE: 1.0000E+00 EA
ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 20:02:57

COLLECT LIVE TIME: 3000. SECONDS
REAL TIME: 3004. SECONDS
DEAD TIME: 0.13 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89
EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:44:48

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	64.17	32.67	1.30	1797.	1154.	5.8	CE-144
2	1219.12	609.58	1.22	154.	150.	14.8	BI-214, RU-103
3	1323.58	661.79	1.74	132.	5581.	1.4	CS-137
4	1822.90	911.38	2.13	63.	95.	17.1	
5	2921.93	1461.07	2.33	14.	669.	4.0	K-40
6	3528.93	1764.85	1.56	8.	56.	16.4	BI-214

ERROR QUOTATION AT 1.00 SIGMA
 PEAK CONFIDENCE LEVEL AT 85.0%

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:44:48

SAMPLE: F5039

DATA COLLECTED ON 2-MAY-90 AT 20:02:57

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

R A D I O N U C L I D E A N A L Y S I S R E P O R T

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI DECAY		
	MEASURED	ERROR	CORRECTED

AC-228	LLD<1.39E+00		LLD<1.39E+00	
AG-108M	LLD<5.26E-01		LLD<5.26E-01	
AG-110M	LLD<3.01E+00		LLD<3.01E+00	
AM-241	LLD<8.97E-01		LLD<8.97E-01	
AM-243	LLD<4.58E-01		LLD<4.58E-01	
AR-41	LLD<6.04E-01		LLD<6.04E-01	
AU-198	LLD<4.93E-01		LLD<4.93E-01	
BA-133	LLD<6.83E-01		LLD<6.83E-01	
BA-139	LLD<1.79E+00		LLD<1.79E+00	
BA-140	LLD<1.91E+00		LLD<1.91E+00	
BA-141	LLD<1.81E+00		LLD<1.81E+00	
BE-7	LLD<5.20E+00		LLD<5.20E+00	
BI-207	LLD<4.31E-01		LLD<4.31E-01	
BI-212	LLD<6.59E+00		LLD<6.59E+00	
BI-214	3.57E+00	+5.30E-01	3.57E+00	+5.30E-01
CD-109	LLD<8.47E+00		LLD<8.47E+00	
CE-139	LLD<4.04E-01		LLD<4.04E-01	
CE-141	LLD<6.72E-01		LLD<6.72E-01	
CEPR144	LLD<5.64E+00		LLD<5.64E+00	
CO-56	LLD<4.42E-01		LLD<4.42E-01	
CO-57	LLD<3.60E-01		LLD<3.60E-01	
CO-58	LLD<3.85E-01		LLD<3.85E-01	
CO-60	LLD<5.12E-01		LLD<5.12E-01	
CR-51	LLD<3.76E+00		LLD<3.76E+00	
CS-134	LLD<5.10E-01		LLD<5.10E-01	
CS-136	LLD<3.83E-01		LLD<3.83E-01	
CS-137	7.53E+01	+1.17E+00	7.53E+01	+1.17E+00
CS-138	LLD<1.10E+00		LLD<1.10E+00	
EU-152	LLD<2.75E+00		LLD<2.75E+00	
EU-154	LLD<1.52E+00		LLD<1.52E+00	
EU-155	LLD<1.41E+00		LLD<1.41E+00	
FE-59	LLD<1.07E+00		LLD<1.07E+00	
HF-181	LLD<6.35E-01		LLD<6.35E-01	
HG-203	LLD<4.82E-01		LLD<4.82E-01	
I-131	LLD<5.46E-01		LLD<5.46E-01	
I-132	LLD<1.11E+00		LLD<1.11E+00	
I-133	LLD<4.72E-01		LLD<4.72E-01	
I-134	LLD<6.60E-01		LLD<6.60E-01	
I-135	LLD<2.09E+00		LLD<2.09E+00	
K-40	1.41E+02	+5.76E+00	1.41E+02	+5.76E+00
KR-85	LLD<1.22E+02		LLD<1.22E+02	
KR-85M	LLD<5.09E-01		LLD<5.09E-01	
KR-87	LLD<1.17E+00		LLD<1.17E+00	
KR-89	LLD<1.87E+01		LLD<1.87E+01	
LA-140	LLD<6.35E-01		LLD<6.35E-01	
LA-142	LLD<1.07E+00		LLD<1.07E+00	
MN-54	LLD<5.33E-01		LLD<5.33E-01	

MN-56	LLD<4.98E-01	LLD<4.98E-01
NA-22	LLD<5.45E-01	LLD<5.45E-01
NA-24	LLD<5.04E-01	LLD<5.04E-01
NB-94	LLD<4.36E-01	LLD<4.36E-01
NB-95	LLD<4.70E-01	LLD<4.70E-01
NB-97	LLD<3.65E+00	LLD<3.65E+00
NP-238	LLD<1.57E+00	LLD<1.57E+00
NP-239	LLD<2.56E+00	LLD<2.56E+00
PA-233	LLD<1.07E+00	LLD<1.07E+00
PA-234M	LLD<8.40E+01	LLD<8.40E+01
PB-210	LLD<1.43E+01	LLD<1.43E+01
PB-212	LLD<9.31E-01	LLD<9.31E-01
PB-214	LLD<1.19E+00	LLD<1.19E+00
PO-210	LLD<3.52E+04	LLD<3.52E+04
PO-214	LLD<4.11E+03	LLD<4.11E+03
PO-216	LLD<2.38E+04	LLD<2.38E+04
PU-239	LLD<4.97E+03	LLD<4.97E+03
PU-241	LLD<1.71E+05	LLD<1.71E+05
RA-224	LLD<9.91E+00	LLD<9.91E+00
RA-226	LLD<1.08E+01	LLD<1.08E+01
RB-88	LLD<5.19E+00	LLD<5.19E+00
RB-89	LLD<2.37E+00	LLD<2.37E+00
RN-220	LLD<4.03E+02	LLD<4.03E+02
RU-103	LLD<5.18E-01	LLD<5.18E-01
RURH106	LLD<8.48E+00	LLD<8.48E+00
SB-124	LLD<4.42E-01	LLD<4.42E-01
SB-125	LLD<4.91E+00	LLD<4.91E+00
SC-46	LLD<5.54E-01	LLD<5.54E-01
SE-75	LLD<6.25E-01	LLD<6.25E-01
SN-113	LLD<7.08E-01	LLD<7.08E-01
SR-85	LLD<5.36E-01	LLD<5.36E-01
SR-91	LLD<8.07E-01	LLD<8.07E-01
SR-92	LLD<7.00E-01	LLD<7.00E-01
TA-182	LLD<1.69E+00	LLD<1.69E+00
TC-99M	LLD<3.74E-01	LLD<3.74E-01
TE-123M	LLD<3.78E-01	LLD<3.78E-01
TE-125M	LLD<1.10E+02	LLD<1.10E+02
TE-132	LLD<4.29E-01	LLD<4.29E-01
TH-228	LLD<2.42E+01	LLD<2.42E+01
TL-208	LLD<4.92E-01	LLD<4.92E-01
U-235	LLD<6.56E-01	LLD<6.56E-01
U-237	LLD<1.77E+00	LLD<1.77E+00
W-187	LLD<1.60E+00	LLD<1.60E+00
XE-131M	LLD<1.69E+01	LLD<1.69E+01
XE-133	LLD<7.93E-01	LLD<7.93E-01
XE-133M	LLD<3.81E+00	LLD<3.81E+00
XE-135	LLD<4.27E-01	LLD<4.27E-01
XE-138	LLD<3.15E+00	LLD<3.15E+00
Y-88	LLD<4.92E-01	LLD<4.92E-01
Y-91	LLD<1.96E+02	LLD<1.96E+02
Y-91M	LLD<6.10E-01	LLD<6.10E-01
ZN-65	LLD<1.34E+00	LLD<1.34E+00
ZR-95	LLD<8.39E-01	LLD<8.39E-01
ZR-97	LLD<3.91E-01	LLD<3.91E-01

TOTAL 2.20E+02 +-5.90E+00 2.20E+02 +-5.90E+00

EBAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 3.33E-08 UC/LI

TOTAL MEASURED ACTIVITY = 2.20E+02 (+-5.90E+00) UC/LI

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
64.17	32.67	1154.	5.8	1.79E+01
1822.90	911.38	95.	17.1	5.36E+00
3528.93	1764.85	56.	16.4	5.36E+00

*
*
* G A M M A S P E C T R U M A N A L Y S I S *
*
* *

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:48:52

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0
DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41
SPECTRUM SIZE: 4096 CHANNELS
ORDER OF SMOOTHING FUNCTION: 5
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK
PEAK CONFIDENCE FACTOR: 85.0%
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV
ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

LLD CALCULATION PERFORMED

MULTIPLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3619

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5040

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 21:09:01

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3011. SECONDS

DEAD TIME: 0.37 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89

EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:48:52

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1C	59.12	30.15	1.39	3364.	655.	14.5	
2C	63.96	32.56	1.39	3218.	4583.	4.0	CE-144
3C	72.71	36.93	1.39	2952.	980.	9.9	TE-125M, I-129
4C	1127.03	563.57	1.52	806.	504.	8.4	CS-134, EU-152
5C	1139.02	569.56	1.52	746.	902.	7.2	CS-134, BI-207
6	1209.73	604.89	1.56	807.	5620.	1.5	CS-134
7	1323.63	661.82	1.64	655.	18292.	0.8	CS-137
8C	1591.96	795.93	1.65	513.	3898.	2.3	CS-134
9C	1604.25	802.07	1.65	467.	411.	7.6	CS-134
10	2346.74	1173.33	1.92	375.	5356.	1.5	CO-60
11	2665.32	1332.68	1.98	94.	5056.	1.4	CO-60
12	2729.84	1364.96	2.00	40.	115.	13.5	CS-134
13	2922.05	1461.13	2.16	36.	654.	4.2	K-40

ERROR QUOTATION AT 1.00 SIGMA
 PEAK CONFIDENCE LEVEL AT 85.0%

C - MULTIPLET ANALYSIS CONVERGED NORMALLY

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:48:52

SAMPLE: F5040

DATA COLLECTED ON 2-MAY-90 AT 21:09:01

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

RADIONUCLIDE ANALYSIS REPORT

NUCLIDE ACTIVITY CONCENTRATION IN uCi/LI

DECAY

MEASURED ERROR CORRECTED ERROR

AC-228	LLD<4.41E+00		LLD<4.41E+00	
AG-108M	LLD<9.99E-01		LLD<9.99E-01	
AG-110M	LLD<5.48E+00		LLD<5.48E+00	
AM-241	LLD<1.14E+00		LLD<1.14E+00	
AM-243	LLD<5.95E-01		LLD<5.95E-01	
AR-41	LLD<9.25E-01		LLD<9.25E-01	
AU-198	LLD<9.08E-01		LLD<9.08E-01	
BA-133	LLD<1.18E+00		LLD<1.18E+00	
BA-139	LLD<2.57E+00		LLD<2.57E+00	
BA-140	LLD<3.95E+00		LLD<3.95E+00	
BA-141	LLD<2.81E+00		LLD<2.81E+00	
BE-7	LLD<9.44E+00		LLD<9.44E+00	
BI-207	LLD<9.69E-01		LLD<9.69E-01	
Bi-212	LLD<1.47E+01		LLD<1.47E+01	
Bi-214	LLD<6.53E+00		LLD<6.53E+00	
CD-109	LLD<1.11E+01		LLD<1.11E+01	
CE-139	LLD<5.83E-01		LLD<5.83E-01	
CE-141	LLD<9.23E-01		LLD<9.23E-01	
CEPR144	LLD<7.85E+00		LLD<7.85E+00	
CO-56	LLD<1.02E+00		LLD<1.02E+00	
CO-57	LLD<4.96E-01		LLD<4.96E-01	
CO-58	LLD<1.00E+00		LLD<1.00E+00	
CO-60	1.06E+02	+-1.69E+00	1.06E+02	+-1.69E+00
CR-51	LLD<6.68E+00		LLD<6.68E+00	
CS-134	6.19E+01	+-1.47E+00	6.19E+01	+-1.47E+00
CS-136	LLD<1.04E+00		LLD<1.04E+00	
CS-137	2.47E+02	+-2.60E+00	2.47E+02	+-2.60E+00
CS-138	LLD<1.10E+00		LLD<1.10E+00	
EU-152	LLD<4.64E+00		LLD<4.64E+00	
EU-154	LLD<2.04E+00		LLD<2.04E+00	
EU-155	LLD<1.86E+00		LLD<1.86E+00	
FE-59	LLD<2.50E+00		LLD<2.50E+00	
HF-181	LLD<1.21E+00		LLD<1.21E+00	
HG-203	LLD<7.72E-01		LLD<7.72E-01	
I-131	LLD<9.80E-01		LLD<9.80E-01	
I-132	LLD<2.05E+00		LLD<2.05E+00	
I-133	LLD<1.04E+00		LLD<1.04E+00	
I-134	LLD<1.53E+00		LLD<1.53E+00	
I-135	LLD<2.81E+00		LLD<2.81E+00	
K-40	1.38E+02	+-5.97E+00	1.38E+02	+-5.97E+00
KR-85	LLD<2.12E+02		LLD<2.12E+02	
KR-85M	LLD<7.10E-01		LLD<7.10E-01	
KR-87	LLD<2.14E+00		LLD<2.14E+00	
KR-89	LLD<3.00E+01		LLD<3.00E+01	
LA-140	LLD<6.16E-01		LLD<6.16E-01	
LA-142	LLD<2.29E+00		LLD<2.29E+00	
MN-54	LLD<1.05E+00		LLD<1.05E+00	

MN-56	LLD<1.15E+00	LLD<1.15E+00
NA-22	LLD<7.85E-01	LLD<7.85E-01
NA-24	LLD<5.61E-01	LLD<5.61E-01
NB-94	LLD<9.13E-01	LLD<9.13E-01
NB-95	LLD<9.64E-01	LLD<9.64E-01
NB-97	LLD<6.64E+00	LLD<6.64E+00
NP-238	LLD<4.88E+00	LLD<4.88E+00
NP-239	LLD<4.46E+00	LLD<4.46E+00
PA-233	LLD<1.83E+00	LLD<1.83E+00
PA-234M	LLD<2.24E+02	LLD<2.24E+02
PB-210	LLD<2.56E+01	LLD<2.56E+01
PB-212	LLD<1.48E+00	LLD<1.48E+00
PB-214	LLD<2.14E+00	LLD<2.14E+00
PO-210	LLD<7.90E+04	LLD<7.90E+04
PO-214	LLD<2.98E+04	LLD<2.98E+04
PO-216	LLD<5.66E+04	LLD<5.66E+04
PU-239	LLD<6.70E+03	LLD<6.70E+03
PU-241	LLD<2.44E+05	LLD<2.44E+05
RA-224	LLD<1.57E+01	LLD<1.57E+01
RA-226	LLD<1.62E+01	LLD<1.62E+01
RB-88	LLD<5.20E+00	LLD<5.20E+00
RB-89	LLD<5.74E+00	LLD<5.74E+00
RN-220	LLD<8.43E+02	LLD<8.43E+02
RU-103	LLD<9.98E-01	LLD<9.98E-01
RURH106	LLD<1.79E+01	LLD<1.79E+01
SB-124	LLD<1.88E+00	LLD<1.88E+00
SB-125	LLD<7.19E+00	LLD<7.19E+00
SC-46	LLD<1.30E+00	LLD<1.30E+00
SE-75	LLD<1.06E+00	LLD<1.06E+00
SN-113	LLD<1.31E+00	LLD<1.31E+00
SR-85	LLD<9.30E-01	LLD<9.30E-01
SR-91	LLD<1.60E+00	LLD<1.60E+00
SR-92	LLD<7.32E-01	LLD<7.32E-01
TA-182	LLD<3.70E+00	LLD<3.70E+00
TC-99M	LLD<5.18E-01	LLD<5.18E-01
TE-123M	LLD<5.54E-01	LLD<5.54E-01
TE-125M	LLD<1.46E+02	LLD<1.46E+02
TE-132	LLD<6.82E-01	LLD<6.82E-01
TH-228	LLD<3.15E+01	LLD<3.15E+01
TL-208	LLD<1.06E+00	LLD<1.06E+00
U-235	LLD<9.87E-01	LLD<9.87E-01
U-237	LLD<2.77E+00	LLD<2.77E+00
W-187	LLD<3.11E+00	LLD<3.11E+00
XE-131M	LLD<2.44E+01	LLD<2.44E+01
XE-133	LLD<1.04E+00	LLD<1.04E+00
XE-133M	LLD<6.09E+00	LLD<6.09E+00
XE-135	LLD<7.04E-01	LLD<7.04E-01
XE-138	LLD<5.28E+00	LLD<5.28E+00
Y-88	LLD<4.92E-01	LLD<4.92E-01
Y-91	LLD<3.23E+02	LLD<3.23E+02
Y-91M	LLD<1.21E+00	LLD<1.21E+00
ZN-65	LLD<2.69E+00	LLD<2.69E+00
ZR-95	LLD<1.74E+00	LLD<1.74E+00
ZR-97	LLD<9.05E-01	LLD<9.05E-01

TOTAL 5.53E+02 +-6.89E+00 5.53E+02 +-6.89E+00

E BAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 2.49E-09 UC/LI

TOTAL MEASURED ACTIVITY = 5.53E+02 (+-6.89E+00) UC/LI

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
59.12	30.15	655.	14.5	1.08E+01
63.96	32.56	4583.	4.0	7.13E+01
72.71	36.93	980.	9.9	1.41E+01
1127.03	563.57	504.	8.4	1.85E+01
1139.02	569.56	902.	7.2	3.34E+01
1604.25	802.07	411.	7.6	2.07E+01
2729.84	1364.96	115.	13.5	9.09E+00

*
* * GAMMA SPECTRUM ANALYSIS *
* *

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:53:26

A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 3.0
DETECTOR NUMBER: 3 / GEOMETRY NUMBER: 41
SPECTRUM SIZE: 4096 CHANNELS
ORDER OF SMOOTHING FUNCTION: 5
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK
PEAK CONFIDENCE FACTOR: 85.0%
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV
ERROR QUOTATION: 1.00 SIGMA UNCERTAINTY

~~LLD~~ CALCULATION PERFORMED

~~MULTIPLLET~~ ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD3621

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F5041

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 5.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 2-MAY-90 AT 22:21:39

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3009. SECONDS

DEAD TIME: 0.30 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-OCT-89

EFFICIENCY CALIBRATION PERFORMED 31-JUL-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:53:26

P E A K A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1C	59.59	30.38	1.38	1712.	423.	14.2	
2C	64.05	32.61	1.38	1627.	2771.	4.5	CE-144
3C	72.76	36.95	1.38	1459.	712.	8.9	TE-125M, I-129
4C	1126.84	563.47	1.47	671.	493.	6.9	CS-134, EU-152
5C	1138.95	569.52	1.47	674.	905.	6.4	CS-134, BI-207
6C	1209.75	604.90	1.49	651.	5777.	1.8	CS-134
7C	1218.81	609.43	1.49	591.	148.	22.6	BI-214, RU-103
8	1323.65	661.83	1.61	532.	11013.	1.0	CS-137
9C	1592.02	795.96	1.66	514.	3990.	2.2	CS-134
10C	1604.31	802.10	1.66	460.	402.	5.3	CS-134
11	2346.65	1173.28	1.98	375.	5600.	1.5	CO-60
12	2665.29	1332.67	2.18	104.	5175.	1.4	CO-60
13	2731.09	1365.59	2.08	49.	107.	15.0	CS-134
14	2921.95	1461.07	1.87	47.	638.	4.4	K-40

ERROR QUOTATION AT 1.00 SIGMA
 PEAK CONFIDENCE LEVEL AT 85.0%

C -- MULTIPLET ANALYSIS CONVERGED NORMALLY

24
 --
 --
 65

222-S COUNTING ROOM WESTINGHOUSE HANFORD

04-MAY-90 12:53:26

SAMPLE: F5041

DATA COLLECTED ON 2-MAY-90 AT 22:21:39

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

RADIONUCLIDE ANALYSIS REPORT

NUCLIDE ACTIVITY CONCENTRATION IN uCi/LI
DECAY

MEASURED ERROR CORRECTED ERROR

AC-228	LLD<9.05E-01		LLD<9.05E-01	
AG-108M	LLD<1.79E-01		LLD<1.79E-01	
AG-110M	LLD<8.55E-01		LLD<8.55E-01	
AM-241	LLD<1.73E-01		LLD<1.73E-01	
AM-243	LLD<9.41E-02		LLD<9.41E-02	
AR-41	LLD<1.87E-01		LLD<1.87E-01	
AU-198	LLD<1.70E-01		LLD<1.70E-01	
BA-133	LLD<2.18E-01		LLD<2.18E-01	
BA-139	LLD<4.64E-01		LLD<4.64E-01	
BA-140	LLD<7.22E-01		LLD<7.22E-01	
BA-141	LLD<5.00E-01		LLD<5.00E-01	
BE-7	LLD<1.76E+00		LLD<1.76E+00	
BI-207	LLD<1.85E-01		LLD<1.85E-01	
BI-212	LLD<2.85E+00		LLD<2.85E+00	
BI-214	7.02E-01	+1.59E-01	7.02E-01	+1.59E-01
CD-109	LLD<1.93E+00		LLD<1.93E+00	
CE-139	LLD<1.05E-01		LLD<1.05E-01	
CE-141	LLD<1.70E-01		LLD<1.70E-01	
CEPR144	LLD<1.43E+00		LLD<1.43E+00	
CO-56	LLD<2.05E-01		LLD<2.05E-01	
CO-57	LLD<8.78E-02		LLD<8.78E-02	
CO-58	LLD<2.03E-01		LLD<2.03E-01	
CU-60	2.18E+01	+3.43E-01	2.18E+01	+3.43E-01
CR-51	LLD<1.21E+00		LLD<1.21E+00	
CS-134	1.27E+01	+2.87E-01	1.27E+01	+2.87E-01
CS-136	LLD<2.00E-01		LLD<2.00E-01	
CS-137	2.97E+01	+3.67E-01	2.97E+01	+3.67E-01
CS-138	LLD<2.06E-01		LLD<2.06E-01	
EU-152	LLD<8.62E-01		LLD<8.62E-01	
EU-154	LLD<3.80E-01		LLD<3.80E-01	
EU-155	LLD<3.20E-01		LLD<3.20E-01	
FE-59	LLD<4.93E-01		LLD<4.93E-01	
HF-181	LLD<2.24E-01		LLD<2.24E-01	
HG-203	LLD<1.44E-01		LLD<1.44E-01	
I-131	LLD<1.74E-01		LLD<1.74E-01	
I-132	LLD<3.45E-01		LLD<3.45E-01	
I-133	LLD<1.85E-01		LLD<1.85E-01	
I-134	LLD<3.09E-01		LLD<3.09E-01	
I-135	LLD<5.30E-01		LLD<5.30E-01	
K-40	2.69E+01	+1.21E+00	2.69E+01	+1.21E+00
KR-85	LLD<3.95E+01		LLD<3.95E+01	
KR-85M	LLD<1.27E-01		LLD<1.27E-01	
KR-87	LLD<3.85E-01		LLD<3.85E-01	
KR-89	LLD<5.57E+00		LLD<5.57E+00	
LA-140	LLD<1.54E-01		LLD<1.54E-01	
LA-142	LLD<4.18E-01		LLD<4.18E-01	
MN-54	LLD<2.11E-01		LLD<2.11E-01	

MN-56	LLD<2.32E-01	LLD<2.32E-01
NA-22	LLD<1.50E-01	LLD<1.50E-01
NA-24	LLD<1.27E-01	LLD<1.27E-01
NB-94	LLD<1.74E-01	LLD<1.74E-01
NB-95	LLD<1.83E-01	LLD<1.83E-01
NB-97	LLD<1.04E+00	LLD<1.04E+00
NP-238	LLD<9.40E-01	LLD<9.40E-01
NP-239	LLD<8.24E-01	LLD<8.24E-01
PA-233	LLD<3.33E-01	LLD<3.33E-01
PA-234M	LLD<4.46E+01	LLD<4.46E+01
PB-210	LLD<4.57E+00	LLD<4.57E+00
PB-212	LLD<2.71E-01	LLD<2.71E-01
PB-214	LLD<3.98E-01	LLD<3.98E-01
PO-210	LLD<1.60E+04	LLD<1.60E+04
PO-214	LLD<5.94E+03	LLD<5.94E+03
PO-216	LLD<1.16E+04	LLD<1.16E+04
PU-239	LLD<1.21E+03	LLD<1.21E+03
PU-241	LLD<4.31E+04	LLD<4.31E+04
RA-224	LLD<2.85E+00	LLD<2.85E+00
RA-226	LLD<2.89E+00	LLD<2.89E+00
RB-88	LLD<1.34E+00	LLD<1.34E+00
RB-89	LLD<1.12E+00	LLD<1.12E+00
RN-220	LLD<1.63E+02	LLD<1.63E+02
RU-103	LLD<1.88E-01	LLD<1.88E-01
QURH106	LLD<3.47E+00	LLD<3.47E+00
SB-124	LLD<3.74E-01	LLD<3.74E-01
SB-125	LLD<1.30E+00	LLD<1.30E+00
SC-46	LLD<2.65E-01	LLD<2.65E-01
SE-75	LLD<1.95E-01	LLD<1.95E-01
SN-113	LLD<2.42E-01	LLD<2.42E-01
SR-85	LLD<1.73E-01	LLD<1.73E-01
SR-91	LLD<3.00E-01	LLD<3.00E-01
SR-92	LLD<1.55E-01	LLD<1.55E-01
TA-182	LLD<7.55E-01	LLD<7.55E-01
TC-99M	LLD<9.45E-02	LLD<9.45E-02
TE-123M	LLD<9.92E-02	LLD<9.92E-02
TE-125M	LLD<2.55E+01	LLD<2.55E+01
TE-132	LLD<1.25E-01	LLD<1.25E-01
TM-228	LLD<5.21E+00	LLD<5.21E+00
TL-208	LLD<2.11E-01	LLD<2.11E-01
U-235	LLD<1.78E-01	LLD<1.78E-01
U-237	LLD<5.14E-01	LLD<5.14E-01
W-187	LLD<5.65E-01	LLD<5.65E-01
XE-131M	LLD<4.46E+00	LLD<4.46E+00
XE-133	LLD<1.68E-01	LLD<1.68E-01
XE-133M	LLD<1.14E+00	LLD<1.14E+00
XE-135	LLD<1.29E-01	LLD<1.29E-01
XE-138	LLD<9.81E-01	LLD<9.81E-01
Y-88	LLD<1.27E-01	LLD<1.27E-01
Y-91	LLD<6.30E+01	LLD<6.30E+01
Y-91M	LLD<2.27E-01	LLD<2.27E-01
ZN-65	LLD<5.23E-01	LLD<5.23E-01
ZR-95	LLD<3.35E-01	LLD<3.35E-01
ZR-97	LLD<1.81E-01	LLD<1.81E-01
<hr/>		
TOTAL	9.18E+01 +-1.35E+00	9.18E+01 +-1.35E+00

EBAR = ***** MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 2.10E-09 UC/LI

TOTAL MEASURED ACTIVITY = 9.18E+01 (+-1.35E+00) UC/LI

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
59.59	30.38	423.	14.2	6.93E+00
64.05	32.61	2771.	4.5	4.31E+01
72.76	36.95	712.	8.9	1.02E+01
1126.84	563.47	493.	6.9	1.81E+01
1138.95	569.52	905.	6.4	3.35E+01
1604.31	802.10	402.	5.3	2.03E+01
2731.09	1365.59	107.	15.0	8.50E+00

% TECH. SPEC. = ***** (+-****)

ERROR QUOTATION AT 1.00 SIGMA
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
59.59	30.38	423.	14.2	6.93E+00
64.05	32.61	2771.	4.5	4.31E+01
72.76	36.95	712.	8.9	1.02E+01
1126.84	563.47	493.	6.9	1.81E+01
1138.95	569.52	905.	6.4	3.35E+01
1604.31	802.10	402.	5.3	2.03E+01
2731.09	1365.59	107.	15.0	8.50E+00

50

40

30

20

10

0

-10

-20

-30

-40

-50

-60

-70

-80

-90

-100

-110

-120

-130

-140

-150

-160

-170

-180

-190

-200

-210

-220

-230

-240

-250

-260

-270

-280

-290

-300

-310

-320

-330

-340

-350

-360

-370

-380

-390

-400

-410

-420

-430

-440

-450

-460

-470

-480

-490

-500

-510

-520

-530

-540

-550

-560

-570

-580

-590

-600

-610

-620

-630

-640

-650

-660

-670

-680

-690

-700

-710

-720

-730

-740

-750

-760

-770

-780

-790

-800

-810

-820

-830

-840

-850

-860

-870

-880

-890

-900

-910

-920

-930

-940

-950

-960

-970

-980

-990

-1000

-1010

-1020

-1030

-1040

-1050

-1060

-1070

-1080

-1090

-1100

-1110

-1120

-1130

-1140

-1150

-1160

-1170

-1180

-1190

-1200

-1210

-1220

-1230

-1240

-1250

-1260

-1270

-1280

-1290

-1300

-1310

-1320

-1330

-1340

-1350

-1360

-1370

-1380

-1390

-1400

-1410

-1420

-1430

-1440

-1450

-1460

-1470

-1480

-1490

-1500

-1510

-1520

-1530

-1540

-1550

-1560

-1570

-1580

-1590

-1600

-1610

-1620

-1630

-1640

-1650

-1660

-1670

-1680

-1690

-1700

-1710

-1720

-1730

-1740

-1750

-1760

-1770

-1780

-1790

-1800

-1810

-1820

-1830

-1840

-1850

-1860

-1870

-1880

-1890

-1900

-1910

-1920

-1930

-1940

-1950

-1960

-1970

-1980

-1990

-2000

-2010

-2020

-2030

-2040

-2050

-2060

-2070

-2080

-2090

-2100

-2110

-2120

-2130

-2140

-2150

-2160

-2170

-2180

-2190

-2200

-2210

-2220

-2230

-2240

-2250

-2260

-2270

-2280

-2290

-2300

-2310

-2320

-2330

-2340

-2350

-2360

-2370

-2380

-2390

-2400

-2410

-2420

-2430

-2440

-2450

-2460

-2470

-2480

-2490

-2500

-2510

-2520

-2530

-2540

-2550

-2560

-2570

-2580

-2590

-2600

-2610

-2620

-2630

-2640

-2650

-2660

-2670

-2680

-2690

-2700

-2710

-2720

-2730

-2740

-2750

-2760

-2770

-2780

-2790

-2800

-2810

-2820

-2830

-2840

-2850

-2860

-2870

-2880

-2890

-2900

-2910

-2920

-2930

-2940

-2950

-2960

-2970

-2980

-2990

-3000

-3010

-3020

-3030

-3040

-3050

-3060

-3070

-3080

-30

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041



Uranium Analysis on the Fusion Dissolution

Instrument	WA77344
Procedure / Rev	LA-925-106/A-2
Technologist	6C269
Date	11/28/89
Temperature	N/A
Starting Time	08:00
Ending Time	14:30
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Initial Check Standard	F5037
2	Reagent Blank	F5052
3	Sample 89-041	F5038
4	Duplicate of 89-041	F5039
5	Spike of 89-041	F5040
6	Final Check Standard	F5041
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	58B38/1 ul			5.7 ml
Spike	58B38/1 ul	Sample/10 ul		5.8 ml

Prepared by:		H. S. Rich Printed Name	Date: May 4, 1990
Verified by:		C. M. Seidel Printed Name	Date: May 4, 1990
Approved by:		L. H. Taylor Printed Name	Date: May 9, 1990

Interim

Rev. E 4-04-90

SST-102

WATER DIGESTION TEST ANALYSIS

9 1 1 2 1 5 9 0 6 1 7

9 1 1 2 2 5 9 0 . 6 7 8

72

Single Shell Tank Project

WATER DIGESTION
Laboratory Results of Solids
Units are Sample Wet Weight

Tank 241-U-110

Core 5

Segment 4

Customer ID: 89-041

Laboratory Segment Serial No.: F5033

Laboratory ID:	Check Standard F5042	Blank F5054	Sample F5043	Sample Duplicate F5044	Spike of Sample F5045	Check Standard F5046
----------------	-------------------------	----------------	-----------------	---------------------------	--------------------------	-------------------------

Water Digestion I			9.59 g/L	8.97 g/L	8.73 g/L	
-------------------	--	--	----------	----------	----------	--

Ion Chromatograph

Fluoride	101.60%	<0.1 ppm	2.31E+04 ug/g	2.04E+04 ug/g	N.C.	87.70%
Chloride	98.40%	<0.1 ppm	<1.05E+03 ug/g	<1.12E+03 ug/g	115.90%	102.00%
Nitrate	98.00%	<1.0 ppm	6.49E+04 ug/g	6.05E+04 ug/g	94.25%	95.40%
Phosphate	94.60%	<1.0 ppm	4.97E+04 ug/g	3.96E+04 ug/g	67.60%	95.40%
Sulfate	92.70%	<1.0 ppm	5.60E+03 ug/g	5.11E+03 ug/g	89.80%	96.00%
Total Organic Carbon/ Carbonate (Not Acidified)	98.50%	2.25E-02 g/L	4.70E+03 ug/g	5.03E+03 ug/g	95.20%	102.30%

Water Digestion II			9.6 g/L	12.03 g/L	10.43 g/L	
--------------------	--	--	---------	-----------	-----------	--

Total Organic Carbon/ Carbonate (Acidified)	98.30%	3.5 ug	8.59E+02 ug/g	1.10E+03 ug/g	100.90%	97.30%
--	--------	--------	---------------	---------------	---------	--------

N.C. = Not Calculated

Single Shell Tank Project

Water Digestion
Sample Results on Laboratory Digestion

Tank 241-U-110

Core 5

Segment 4

Customer ID: 89-041 Laboratory Segment Serial No.: F5033

	Check Standard	Blank	Sample	Sample Duplicate	Spike of Sample	Check Standard
Laborator ID:	F5042	F5054	F5043	F5044	F5045	F5046

Ion Chromatograph

Fluoride	101.60%	<0.1 ppm	222 ppm	183 ppm	N.C.	87.70%
Chloride	98.40%	<0.1 ppm	<10.1 ppm	<10.1 ppm	115.90%	102.00%
Nitrate	98.00%	<1.0 ppm	622 ppm	543 ppm	94.25%	95.40%
Phosphate	94.60%	<1.0 ppm	477 ppm	355 ppm	67.60%	95.40%
Sulfate	92.70%	<1.0 ppm	53.7 ppm	45.8 ppm	89.80%	96.00%

Total Organic Carbon/ Carbonate (Not Acidified)	98.50%	2.25E-02 g/L	4.51E-02 g/L	4.51E-02 g/L	95.20%	102.30%
--	--------	--------------	--------------	--------------	--------	---------

Water Digestion II			9.6 g/L	12.03 g/L	10.43 g/L	
--------------------	--	--	---------	-----------	-----------	--

Total Organic Carbon/ Carbonate (Acidified)	98.30%	3.5 ug	8.25E-03 g/L	1.32E-02 g/L	100.90%	97.30%
--	--------	--------	--------------	--------------	---------	--------

N.C.= Not Calculated

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	N/A
Procedure / Rev	LA-504-101/A-2
Technologist	80725
Date	12/04/89
Temperature	25 C
Starting Time	12/01/89 13:30 hours
Ending Time	12/04/89 08:45 hours
Chemist	H. S. Rich

Water Digestion of sample 89-041.
 Note: Sample is not spiked prior to digestion. This procedure provides a sample to be spiked later with the appropriate elements.

	Description	Lab. Id.
1	Reagent Blank	F5054
2	Sample 89-041	F5043
3	Duplicate of 89-041	F5044
4	Spike of 89-041	F5045
5		
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
N/A				
Spike (see note)				

Prepared by:	<u>H. S. Rich</u> Signature	H. S. Rich Printed Name	Date: May 4, 1990
Verified by:	<u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: May 4, 1990
Approved by:	<u>L. H. Taylor</u> Signature	L. H. Taylor Printed Name	Date: May 9, 1990

Analytical Batch

Lab Segment Serial No.: F5026

Customer ID.: 89-041

Instrument	N/A
Procedure / Rev	LA-504-101 / A-2
Technologist	E. Colvin
Date	06/29/90
Temperature	Not Reported.
Starting Time	06/27/90: 08:00
Ending Time	06/29/90; 11:30
Chemist	H. Rich

Water Digestion - Second Leach.

NOTE: Sample is not spiked prior to digestion. This procedure provides a sample to be spiked later with the appropriate elements.

	Description	Lab. Id.
1	Blank	89-041
2	Sample	89-040
3	Duplicate	89-040
4	Spike	89-040
5	Sample	89-041
6	Duplicate	89-041
7	Spike	89-041
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
Standard Type				
N/A				

Prepared by:	<u>S. A. Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 8-09-90
Verified by:	<u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 8-09-90
Approved by:	<u>L. H. Taylor</u> Signature	L. H. Taylor Printed Name	Date: 8-30-90

Interim

Rev E 4/04/90

SST-102

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	WB24721
Procedure / Rev	LA-533-105/A-3
Technologist	6B107 N. E. Wright
Date	12/05/89
Temperature	26 C
Starting Time	10:30
Ending Time	14:30
Chemist	H. S. Rich

Ion Chromatograph Analysis

Water Digestion

Note: Ion chromatograph was not calibrated for nitrite. Nitrite results are not reported.

	Description	Lab. Id.
1	Initial Check Standard	F5042
2	Reagent Blank	F5054
3	Sample 89-041	F5043
4	Duplicate of 89-041	F5044
5	Spike of 89-041	F5045
6	Ending Check Standard	F5046
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	6C11HC/100uL			10.1 mL
Spike	40C9-A/50 uL	Sample/10 uL		5.06 mL

Interim

4-04-90

Rev.E

SST-102

Prepared by: H. S. Rich H. S. Rich Date: May 4, 1990

Printed Name

Verified by: C. M. Seidel C. M. Seidel Date: May 4, 1990

Printed Name

Approved by: L. H. Taylor L. H. Taylor Date: May 9, 1990

Printed Name

**Single Shell Tank
Calibration Record**

**Phase
I-A**

Analyte: Ion Chromatograph

Procedure LA-533-105

Revision: A-3

Instrument: Dionex 4000

Property Number: WB24721

Technologist: Nora Wright

Payroll Number: 6B107

Date: Dec. 1, 1989

Calibration Standard ID: Book number 35C9-60 issued 11-22-89

Analyte Concentration: F=49.6; Cl=61.0; NO₃=500.5; PO₄=500.6; SO₄=500.5 (in ppm)

Type of Calibration: Quadratic least squares

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3			
4	See Attached	Calibration Sheets.	
5			
6			
7			
8			
9			
10			

Interim

Comments:

Rev. (Draft) 1/18/89

Prepared by:

Signature

H. S. Rich

Printed Name

Date: 5-08-90

Verified by:

Signature

C. M. Seidel

Printed Name

Date: 5-08-90

Approved by:

Signature

L. H. Taylor

Printed Name

Date: 5/9/90

SST-105

Detector Parameters

Number of Detectors.....	1
Detector 1 Type.....	CDM-1

Report Options

Run Time (minutes).....	10.00
Detector 1 real time plot scale.....	20.00
Print Report.....	Yes
Print Replot.....	Yes
AutoScale Replot to Highest Peak.....	Yes
Print Retention Times on Chromatogram.....	Yes
List Peaks Not Found in this run.....	No
Report Unknowns found in run.....	Yes
Record Raw Data.....	Yes
Raw Data File Name: c:\dx\data\89120102.d09	
Record Result Data.....	No

Integration Parameters

Sampling Rate (seconds).....	0.20
Peak Threshold (mV or uS/data pt interval).....	0.400
Starting Peak Width (seconds).....	10.0
Peak Area Reject.....	1000

Integration Timed Events

Time	Description
-----	-----

Calibration Parameters

External or Internal Calibration.....	External
Calibrate by Area or Height.....	Height
Replace Or Average Calibrations.....	Replace
Number Of Levels for Calibration.....	6
Calibration fit type.....	Quadratic
Response Factor for unknown peaks.....	0.0
Default Injection Volume.....	1.0
Default Dilution Factor.....	1.0
Area Reject for Reference Peaks.....	1000
Percent Retention Time Window for Reference Peaks.....	5.0

Component # 1 FLUORIDE Retention Time 0.98
 Reference Peak FLUORIDE Window Size 5.00%
 Least Squares Slope = 3.04147E-004
 Least Squares Intercept = 3.47570E-002
 Ka = -7.20865E-010

Level	Amount	Area	Height
1	9.90000E-002	1547	292
2	2.46800E-001	3814	699
3	4.91900E-001	8014	1443
4	9.72600E-001	18242	3026
5	1.90760E+000	37268	6331
6	3.67410E+000	82101	12310

Component # 2 CHLORIDE Retention Time 1.62
 Reference Peak FLUORIDE Window Size 7.00%
 Least Squares Slope = 5.28644E-004
 Least Squares Intercept = 2.16858E-002
 Ka = -4.60695E-009

Level	Amount	Area	Height
1	1.21800E-001	1174	211
2	3.03500E-001	3127	534
3	6.03960E-001	6384	1090
4	1.19610E+000	13048	2257
5	2.34610E+000	26949	4596
6	4.51840E+000	56310	9249

Component # 3 NITRITE Retention Time 2.00
 Reference Peak FLUORIDE Window Size 7.00%
 Least Squares Slope = 8.73191E-004
 Least Squares Intercept = 4.13109E-001
 Ka = 4.52304E-009

Level	Amount	Area	Height
1	1.06600E+000	6728	955
2	2.65700E+000	17398	2609
3	5.28700E+000	36144	5455
4	1.04706E+001	72479	10242
5	2.05385E+001	156757	21248
6	3.95550E+001	299406	37444

Component # 4 NITRATE Retention Time 3.72
 Reference Peak FLUORIDE Window Size 10.00%
 Least Squares Slope = 1.94357E-003
 Least Squares Intercept = -1.10952E-001
 Ka = 1.78504E-008

Level	Amount	Area	Height
1	9.99000E-001	5537	519
2	2.49000E+000	15382	1323
3	4.95550E+000	33066	2632
4	9.81380E+000	67697	4864
5	1.92492E+001	142777	9168
6	3.70730E+001	296013	16606

Component # 5 PHOSPHATE Retention Time 5.22
 Reference Peak FLUORIDE Window Size 7.00%
 Least Squares Slope = 4.25699E-003
 Least Squares Intercept = 1.29081E-001
 Ka = -2.86430E-008



Level	Amount	Area	Height
1	9.99200E-001	2588	195
2	2.49050E+000	11460	648
3	4.95640E+000	16063	1091
4	9.81580E+000	33316	2235
5	1.92531E+001	72002	4693
6	3.70804E+001	149467	9246



Component # 6 SULFATE Retention Time 7.10
 Reference Peak FLUORIDE Window Size 10.00%
 Least Squares Slope = 2.01226E-003
 Least Squares Intercept = 2.11736E-001
 Ka = -1.25909E-009



Level	Amount	Area	Height
1	9.99000E-001	8490	472
2	2.49000E+000	20937	1134
3	4.95550E+000	42330	2329
4	9.81380E+000	84876	4644
5	1.92492E+001	182078	9631
6	3.70730E+001	373505	18512

Component # 7 Oxalate Retention Time 9.77
 Reference Peak FLUORIDE Window Size 10.00%
 Least Squares Slope = 0.00000E+000
 Least Squares Intercept = 0.00000E+000
 Ka = 0.00000E+000

Level	Amount	Area	Height
1	0.00000E+000	0	0
2	0.00000E+000	0	0
3	0.00000E+000	0	0
4	0.00000E+000	0	0
5	0.00000E+000	0	0
6	0.00000E+000	98993	5848







***** AUTOMATIC CALIBRATION UPDATE *****

Method File: C:\DX\METHOD\Grout01.met

Calibration Level : 1

***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.77	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET. TIME	MEASURED RET. TIME	NEW RET. TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
1	FLUORIDE	1.00	1.00	1.00	2.917e+002	2.917e+002	2.917e+002
2	CHLORIDE	1.63	1.63	1.63	2.108e+002	2.108e+002	2.108e+002
3	NITRITE	2.00	2.00	2.00	9.545e+002	9.545e+002	9.545e+002
4	NITRATE	4.03	4.03	4.03	5.186e+002	5.186e+002	5.186e+002
5	PHOSPHATE	5.35	5.35	5.35	1.954e+002	1.954e+002	1.954e+002
6	SULFATE	7.13	7.13	7.13	4.723e+002	4.723e+002	4.723e+002

DATA REPROCESSED ON Wed Mar 21 13:29:28 1990

Sample Name: autocall1R Date: Fri Dec 01 12:44:23 1989
 Data File : c:\dx\data\89120102.d04
 Method : C:\DX\METHOD\Grout01.met
 Interface : 1 System : 1 Inject# : 4 Detector: CDM


***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2821

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	1.00	FLUORIDE	9.900e-002	1.547e+003	292	1	0 0.00%
2	1.63	CHLORIDE	1.218e-001	1.174e+003	211	1	0 0.00%
3	2.00	NITRITE	1.066e+000	6.728e+003	955	1	0 0.00%
4	4.03	NITRATE	9.990e-001	5.537e+003	519	1	0 0.00%
5	5.35	PHOSPHATE	9.992e-001	2.588e+003	195	1	0 0.00%
6	7.13	SULFATE	9.990e-001	8.490e+003	472	1	0 0.00%



***** AUTOMATIC CALIBRATION UPDATE *****

Method File: C:\DX\METHOD\Grout01.met

Calibration Level : 2

***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.61	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET. TIME	MEASURED RET. TIME	NEW RET. TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
----------	----------------	---------------	--------------------	---------------	------------	-----------------	------------

1	FLUORIDE	1.00	0.98	0.98	2.917e+002	6.990e+002	6.990e+002
2	CHLORIDE	1.63	1.62	1.62	2.108e+002	5.339e+002	5.339e+002
3	NITRITE	2.00	2.02	2.02	9.545e+002	2.609e+003	2.609e+003
4	NITRATE	4.03	4.02	4.02	5.186e+002	1.323e+003	1.323e+003
5	PHOSPHATE	5.35	5.40	5.40	1.954e+002	6.483e+002	6.483e+002
6	SULFATE	7.13	7.22	7.22	4.723e+002	1.134e+003	1.134e+003


DATA REPROCESSED ON Wed Mar 21 13:29:56 1990

Sample Name: autocal2R Date: Fri Dec 01 12:54:32 1989
 Data File : c:\dx\data\89120102.d05
 Method : C:\DX\METHOD\Grout01.met
 Interface : 1 System : 1 Inject# : 5 Detector: CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2821
 Area reject = 1000 One Data Point per 0.2 seconds
 Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	% DELTA BL PEAK RET TIME
1	0.98	FLUORIDE	2.468e-001	3.814e+003	699	1 0 0.00%
2	1.62	CHLORIDE	3.035e-001	3.127e+003	534	1 0 0.00%
3	2.02	NITRITE	2.657e+000	1.740e+004	2609	1 0 0.00%
4	4.02	NITRATE	2.490e+000	1.538e+004	1323	1 0 0.00%
5	5.40	PHOSPHATE	2.491e+000	1.146e+004	648	1 0 0.00%
6	7.22	SULFATE	2.490e+000	2.094e+004	1134	1 0 0.00%



***** AUTOMATIC CALIBRATION UPDATE *****
 Method File: C:\DX\METHOD\Grout01.met
 Calibration Level : 3
 ***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.77	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET.TIME	MEASURED RET.TIME	NEW RET.TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
1	FLUORIDE	0.98	0.98	0.98	2.917e+002	1.443e+003	1.443e+003
2	CHLORIDE	1.62	1.62	1.62	2.108e+002	1.090e+003	1.090e+003
3	NITRITE	2.02	2.02	2.02	9.545e+002	5.455e+003	5.455e+003
4	NITRATE	4.02	3.97	3.97	5.186e+002	2.632e+003	2.632e+003
5	PHOSPHATE	5.40	5.35	5.35	1.954e+002	1.091e+003	1.091e+003
6	SULFATE	7.22	7.22	7.22	4.723e+002	2.329e+003	2.329e+003


DATA REPROCESSED ON Wed Mar 21 13:30:26 1990

Sample Name: autocal3r Date: Fri Dec 01 13:04:42 1989
 Data File : c:\dx\data\89120102.d06
 Method : C:\DX\METHOD\Grout01.met
 Interface : 1 System : 1 Inject# : 6 Detector: CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2821
 Area reject = 1000 One Data Point per 0.2 seconds
 Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	% DELTA BL PEAK	RET TIME	
1	0.98	FLUORIDE	4.919e-001	8.014e+003	1443	1	0	0.00%
2	1.62	CHLORIDE	6.040e-001	6.384e+003	1090	1	0	0.00%
3	2.02	NITRITE	5.287e+000	3.614e+004	5455	1	0	0.00%
4	3.97	NITRATE	4.956e+000	3.307e+004	2632	1	0	0.00%
5	5.35	PHOSPHATE	4.956e+000	1.606e+004	1091	1	0	0.00%
6	7.22	SULFATE	4.955e+000	4.233e+004	2329	1	0	0.00%



***** AUTOMATIC CALIBRATION UPDATE *****
Method File: C:\DX\METHOD\Grout01.met
Calibration Level : 4
***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.77	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET.TIME	MEASURED RET.TIME	NEW RET.TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
1	FLUORIDE	0.98	0.98	0.98	2.917e+002	3.026e+003	3.026e+003
2	CHLORIDE	1.62	1.62	1.62	2.108e+002	2.257e+003	2.257e+003
3	NITRITE	2.02	2.00	2.00	9.545e+002	1.024e+004	1.024e+004
4	NITRATE	3.97	3.88	3.88	5.186e+002	4.864e+003	4.864e+003
5	PHOSPHATE	5.35	5.28	5.28	1.954e+002	2.235e+003	2.235e+003
6	SULFATE	7.22	7.13	7.13	4.723e+002	4.644e+003	4.644e+003

DATA REPROCESSED ON Wed Mar 21 13:31:01 1990

6
=====


Sample Name:	Autocal4r	Date:	Fri Dec 01 13:14:52 1989
Data File :	c:\dx\data\89120102.d07		
Method :	C:\DX\METHOD\Grout01.met		
Interface :	1	System :	1
		Inject#:	7
		Detector:	CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2820
Area reject = 1000 One Data Point per 0.2 seconds
Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	% DELTA BL PEAK RET TIME
1	0.98	FLUORIDE	9.726e-001	1.824e+004	3026	1 0 0.00%
2	1.62	CHLORIDE	1.196e+000	1.305e+004	2257	1 0 0.00%
3	2.00	NITRITE	1.047e+001	7.248e+004	10242	1 0 0.00%
4	3.88	NITRATE	9.814e+000	6.770e+004	4864	1 0 0.00%
5	5.28	PHOSPHATE	9.816e+000	3.332e+004	2235	1 0 0.00%
6	7.13	SULFATE	9.814e+000	8.488e+004	4644	1 0 0.00%

06

File: C:\DX\DATA\89120102.D07 Sample: AUTOCAL4R
9 1 1 2 0 5 2 0 7 0 7

***** AUTOMATIC CALIBRATION UPDATE *****
 Method File: C:\DX\METHOD\Grout01.met
 Calibration Level : 5
 ***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.77	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET. TIME	MEASURED RET. TIME	NEW RET. TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
1	FLUORIDE	0.98	0.98	0.98	2.917e+002	6.331e+003	6.331e+003
2	CHLORIDE	1.62	1.62	1.62	2.108e+002	4.596e+003	4.596e+003
3	NITRITE	2.00	2.00	2.00	9.545e+002	2.125e+004	2.125e+004
4	NITRATE	3.88	3.82	3.82	5.186e+002	9.168e+003	9.168e+003
5	PHOSPHATE	5.28	5.27	5.27	1.954e+002	4.693e+003	4.693e+003
6	SULFATE	7.13	7.15	7.15	4.723e+002	9.631e+003	9.631e+003


DATA REPROCESSED ON Wed Mar 21 13:31:36 1990

Sample Name: Autocal5r	Date: Fri Dec 01 13:25:03 1989
Data File : c:\dx\data\89120102.d08	
Method : C:\DX\METHOD\Grout01.met	
Interface : 1 System : 1	Inject# : 8 Detector: CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2820
 Area reject = 1000 One Data Point per 0.2 seconds
 Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.98	FLUORIDE	1.908e+000	3.727e+004	6331	1	0 0.00%
2	1.62	CHLORIDE	2.346e+000	2.695e+004	4596	2	0 0.00%
3	2.00	NITRITE	2.054e+001	1.568e+005	21248	2	0 0.00%
4	3.82	NITRATE	1.925e+001	1.428e+005	9168	1	0 0.00%
5	5.27	PHOSPHATE	1.925e+001	7.200e+004	4693	1	0 0.00%
6	7.15	SULFATE	1.925e+001	1.821e+005	9631	1	0 0.00%



***** AUTOMATIC CALIBRATION UPDATE *****

Method File: C:\DX\METHOD\Grout01.met

Calibration Level : 6

***** PEAKS NOT FOUND IN THIS RUN *****

Name	Adjusted Ret Time	Reference Peak
Oxalate	9.77	1

***** COMPONENTS FOUND IN THIS RUN *****

COMP NUM	COMPONENT NAME	OLD RET. TIME	MEASURED RET. TIME	NEW RET. TIME	OLD HEIGHT	MEASURED HEIGHT	NEW HEIGHT
----------	----------------	---------------	--------------------	---------------	------------	-----------------	------------

1	FLUORIDE	0.98	0.98	0.98	2.917e+002	1.231e+004	1.231e+004
2	CHLORIDE	1.62	1.62	1.62	2.108e+002	9.249e+003	9.249e+003
3	NITRITE	2.00	2.00	2.00	9.545e+002	3.744e+004	3.744e+004
4	NITRATE	3.82	3.72	3.72	5.186e+002	1.661e+004	1.661e+004
5	PHOSPHATE	5.27	5.22	5.22	1.954e+002	9.246e+003	9.246e+003
6	SULFATE	7.15	7.10	7.10	4.723e+002	1.851e+004	1.851e+004


DATA REPROCESSED ON Wed Mar 21 13:32:10 1990

Sample Name:	Autocal6r	Date:	Fri Dec 01 13:35:14 1989			
Data File :	c:\dx\data\89120102.d09					
Method :	C:\DX\METHOD\Grout01.met					
Interface :	1	System :	1 Inject#:	9	Detector:	CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2821
 Area reject = 1000 One Data Point per 0.2 seconds
 Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.98	FLUORIDE	3.674e+000	8.210e+004	12310	2	0 0.00%
2	1.62	CHLORIDE	4.518e+000	5.631e+004	9249	2	0 0.00%
3	2.00	NITRITE	3.956e+001	2.994e+005	37444	2	0 0.00%
4	3.72	NITRATE	3.707e+001	2.960e+005	16606	1	0 0.00%
5	5.22	PHOSPHATE	3.708e+001	1.495e+005	9246	1	0 0.00%
6	7.10	SULFATE	3.707e+001	3.735e+005	18512	1	0 0.00%



DIONEX SCHEDULE - C:\DX\SCHEDULE\89120401.SCH

Inj#	Sample Name	Method Name	Data File	Vol.	Dil.	Int.Std.
1	SETUP	...\\GROUT01	...\\891204011	1	0	
2	BLANK	...\\GROUT01	...\\891204011	1	0	
3	LMCS/6C11HC	...\\GROUT01	...\\891204011	101	0	
4	5054b	...\\GROUT01	...\\891204011	1	0	
5	5043	...\\GROUT01	...\\891204011	101	0	
6	5044d	...\\GROUT01	...\\891204011	101	0	
7	5045s	...\\GROUT01	...\\891204011	101	0	
8	LMCS/6C11HC	...\\GROUT01	...\\891204011	101	0	

Q
-
7
0
6
In
C
N
-
6

DATA REPROCESSED ON Mon Mar 26 14:33:50 1990

=====
Sample Name: BLANK Date: Mon Dec 04 14:55:38 1989
Data File : C:\DX\DATA\89120401.D02
Method : c:\windows\ai400\method\GROUT01.met
ACI Address: 1 System : 1 Inject#: 2 Detector: CDM
=====

***** EXTERNAL STANDARD REPORT *****


Stop time = 9.40 Minutes Number of Data Points = 2821

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL	PEAK	RET TIME
-------------	-------------	--------------	-------------------	------	-----	---------	----	------	----------

File: C:\DX\DATA\89120401.D02 Sample: BLANK



DATA REPROCESSED ON Mon Mar 26 14:56:13 1990

Sample Name: LMCS/6C11HC

Date: Mon Dec 04 15:05:40 1989

Data File : C:\DX\DATA\89120401.D03

Method : c:\windows\ai400\method\GROUT01.met

ACI Address: 1 System : 1 Inject#: 3 Detector: CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes

Number of Data Points = 2821

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF BL	% DELTA PEAK	RET TIME
1	1.00	FLUORIDE	6.298e+001	1.241e+004	1924	2	0	0.00%
2	1.35		0.000e+000	2.758e+003	277	2		
3	1.63	CHLORIDE	7.381e+001	8.452e+003	1378	2	0	0.00%
4	3.98	NITRATE	6.085e+002	3.925e+004	3068	1	0	7.17%
5	5.52	PHOSPHATE	5.873e+002	1.996e+004	1343	1	0	5.75%
6	7.57	SULFATE	5.759e+002	5.145e+004	2728	1	0	6.57%

File: C:\DX\DATA\89120401.D03 Sample: LMCS/6C11HC



DATA REPROCESSED ON Mon Mar 26 14:36:45 1990


=====
Sample Name: 5054B Date: Mon Dec 04 15:15:46 1989
Data File : C:\DX\DATA\89120401.D04
Method : c:\windows\ai400\method\GROUT01.met
ACI Address: 1 System : 1 Inject#: 4 Detector: CDM
=====

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2820
Area reject = 1000 One Data Point per 0.2 seconds
Amount Injected = 1 Dilution factor = 1

PEAK NUM.	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA
					HEIGHT	BL PEAK
					RET	TIME

File: C:\DX\DATA\89120401.D04 Sample: 5054B



DATA REPROCESSED ON Mon Mar 26 15:01:14 1990


=====
Sample Name: 5043 Date: Mon Dec 04 15:25:50 1989
Data File : C:\DX\DATA\89120401.D05
Method : c:\windows\ai400\method\GROUT01.met
ACI Address: 1 System : 1 Inject#: 5 Detector: CDM
=====

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes Number of Data Points = 2821
Area reject = 1000 One Data Point per 0.2 seconds
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	1.00	FLUORIDE	2.218e+002	3.814e+004	7153	1	0 0.00%
2	2.03		0.000e+000	1.004e+004	1419	1	
3	3.98	NITRATE	6.225e+002	4.016e+004	3135	1	0 7.17%
4	5.55	PHOSPHATE	4.773e+002	1.624e+004	1083	1	0 6.39%
5	7.67	SULFATE	5.367e+001	2.528e+003	151	1	0 7.98%

File: C:\DX\DATA\89120401.D05 Sample: 5043



DATA REPROCESSED ON Mon Mar 26 15:02:25 1990

=====
Sample Name: 5044D Date: Mon Dec 04 15:35:54 1989
Data File : C:\DX\DATA\89120401.D06
Method : c:\windows\ai400\method\GROUT01.met
ACI Address: 1 System : 1 Inject#: 6 Detector: CDM
=====

***** EXTERNAL STANDARD REPORT *****


Stop time = 9.40 Minutes Number of Data Points = 2821

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF BL	% DELTA PEAK	RET TIME
1	1.00	FLUORIDE	1.830e+002	3.166e+004	5863	1	0	0.00%
2	2.03		0.000e+000	8.677e+003	1287	1		
3	3.98	NITRATE	5.427e+002	3.459e+004	2750	1	0	7.17%
4	5.52	PHOSPHATE	3.551e+002	1.151e+004	795	1	0	5.75%
5	7.62	SULFATE	4.576e+001	1.737e+003	112	1	0	7.28%

File: C:\DX\DATA\89120401.D06 Sample: 5044D



DATA REPROCESSED ON Mon Mar 26 15:03:39 1990

=====
Sample Name: 5045S Date: Mon Dec 04 15:45:59 1989
Data File : C:\DX\DATA\89120401.D07
Method : c:\windows\ai400\method\GROUT01.met
ACI Address: 1 System : 1 Inject#: 7 Detector: CDM
=====

***** EXTERNAL STANDARD REPORT *****


Stop time = 9.40 Minutes Number of Data Points = 2820

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	1.00	FLUORIDE	2.226e+002	3.912e+004	7181	1	0 0.00%
2	1.63	CHLORIDE	6.858e+001	7.304e+003	1276	1	0 0.00%
3	2.03		0.000e+000	4.559e+004	6752	1	
4	3.90	NITRATE	1.027e+003	7.074e+004	5053	1	0 0.00%
5	5.47	PHOSPHATE	7.565e+002	2.582e+004	1746	1	0 0.00%
6	7.52	SULFATE	4.939e+002	4.362e+004	2322	1	0 5.87%

File: C:\DX\DATA\89120401.D07 Sample: 5045S



DATA REPROCESSED ON Mon Mar 26 15:05:05 1990

Sample Name: LMCS/6C11HC

Date: Mon Dec 04 15:56:04 1989

Data File : C:\DX\DATA\89120401.D08

Method : c:\windows\ai400\method\GROUT01.met

ACI Address: 1 System : 1 Inject#: 8 Detector: CDM

***** EXTERNAL STANDARD REPORT *****

Stop time = 9.40 Minutes

Number of Data Points = 2820

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF BL	% DELTA PEAK	RET TIME
1	0.98	FLUORIDE	5.439e+001	9.682e+003	1645	1	0	0.00%
2	1.63	CHLORIDE	7.652e+001	7.955e+003	1431	1	0	0.00%
3	3.93	NITRATE	5.927e+002	3.897e+004	2992	1	0	5.83%
4	5.42	PHOSPHATE	5.922e+002	2.019e+004	1355	1	0	0.00%
5	7.45	SULFATE	5.959e+002	5.281e+004	2826	1	0	0.00%

File: C:\DX\DATA\89120401.D08 Sample: LMCS/6C11HC



Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	WB39937
Procedure / Rev	LA-344-105/A-3
Technologist	80725 S. Cervantes
Date	12-4-89
Temperature	25 C
Starting Time	12-4-89/11:00
Ending Time	12-4-89/14:00
Chemist	R. E. Brandt

Total Organic Carbon

Water Digestion

Samples were not acidified prior to analysis.
Results reported as Total Organic Carbon
and Carbonate combined.

	Description	Lab. Id.
1	Reagent Blank	F5054
2	Initial Check Std	F5042
3	Sample 89-041	F5043
4	Duplicate 89-041	F5044
5	Spike 89-041	F5045
6	Final LMCS Check Std	F5046
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		


Interim

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Standard	70C11B/11 uL			2.2 mL
Spike	70C11B/50 ug	Sample/200 uL		2.4 mL

Rev.E
4/04/90

Approved by:

SST-102

Prepared by:  Signature	H. S. Rich Printed Name	Date: May 4, 1990
Verified by:  Signature	C. M. Seidel Printed Name	Date: May 4, 1990
Approved by: _____ Signature	Printed Name	Date: _____

This page intentionally left blank.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: BLK Date: 12-04-1989 Time: 10:32:41

Blank = N/A Sample Size = 200 Dilution Factor = 11
% Difference = 10 Min Readings = 7 Max Readings = 10

	Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00	
2	2.01	1.10	100.00	
3	3.01	1.50	26.67	
4	4.01	1.90	21.05	
5	5.01	2.20	13.64	
6	6.01	2.60	15.38	
7	7.01	2.90	10.34	
8	8.01	3.40	14.71	
9	9.01	3.80	10.53	
10	10.01	4.10	7.32	

BLANK VALUE = 4.1 / 10.00653 = .4097324 g/L/minute

Sample Run By: 80725 J.C.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: STD/5042 Date: 12-04-1989 Time: 11:08:15

Blank = .4097324 Sample Size = 200 Dilution Factor = 11
% Difference = 10 Min Readings = 7 Max Readings = 10

== Reading ===== Analysis Time ===== Coulometer ===== % Difference ==
1 1.01 0.00 0.00

2	2.01	49.00	100.00
3	3.01	53.60	6.58
4	4.01	54.80	2.19
5	5.01	55.50	1.26
6	6.01	56.10	1.07
7	7.01	56.60	0.88

$$(56.6 - 2.870828)(11)/(200) = 2.955104 \text{ g/L Carbon}$$

$$(56.6 - 2.870828)(11)/(200)(12) = .2462587 \text{ Molar Carbon}$$

Sample Run By: 80725_A.C.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: 5043

Date: 12-04-1989 Time: 11:57:24

Blank = .4097324
% Difference = 10

Sample Size = 200 Dilution Factor = 11
Min Readings = 7 Max Readings = 10

== Reading ==	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	1.30	100.00
3	3.01	1.80	27.78
4	4.01	2.20	18.18
5	5.01	2.70	18.52
6	6.01	3.20	15.63
7	7.01	3.70	13.51
8	8.01	4.10	9.76

A
c
7
0
6
L
C
V
—
I
0

(4.1 - 3.280035)(11) / (200) = 4.509807E-02 g/L Carbon

(4.1 - 3.280035)(11) / (200)(12) = 3.758172E-03 Molar Carbon

Sample Run By: 80725 J.C.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: 5044

Date: 12-04-1989

Time: 13:40:45

Blank = .4097324

% Difference = 10

Sample Size = 200 Dilution Factor = 11

Min Readings = 7 Max Readings = 10

== Reading ===== Analysis Time ===== Coulometer ===== % Difference ==
1 1.01 0.00 0.00

2 2.01 1.60 100.00

3 3.01 2.00 20.00

4 4.01 2.50 20.00

5 5.01 2.90 13.79

6 6.01 3.30 12.12

7 7.01 3.70 10.81

8 8.01 4.10 9.76

(4.1 - 3.28046)(11) / (200) = 4.507469E-02 g/L Carbon

(4.1 - 3.28046)(11) / (200)(12) = 3.756224E-03 Molar Carbon

Sample Run By: 80725 A.C.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: 5045

Date: 12-04-1989 Time: 13:52:28

Blank = .4097324
% Difference = 10

Sample Size = 200 Dilution Factor = 11
Min Readings = 7 Max Readings = 10

== Reading ===== Analysis Time ===== Coulometer ===== % Difference ==
1 1.01 0.00 0.00

2	2.01	47.50	100.00
3	3.01	49.40	3.85
4	4.01	50.30	1.79
5	5.01	50.70	0.79
6	6.01	51.20	0.98
7	7.01	51.70	0.97

6
7
0
0
6
5
7
1
2
1
7
8

$$(51.7 - 2.870453) * (11) / (200) = 2.685625 \text{ g/L Carbon}$$

$$(51.7 - 2.870453) * (11) / (200) * (12) = .2238021 \text{ Molar Carbon}$$

Sample Run By: 80725 J.C.

COULOMETER ANALYSIS REPORT
TICTOC Rev. 0

Sample: 5046

Date: 12-04-1989 Time: 14:03:49

Blank = .4097324
% Difference = 10

Sample Size = 200 Dilution Factor = 11
Min Readings = 7 Max Readings = 10

	Reading	Analysis Time	Coulometer	% Difference
1		1.01	0.00	0.00
2		2.01	53.30	100.00
3		3.01	56.20	5.16
4		4.01	57.20	1.75
5		5.01	57.80	1.04
6		6.01	58.20	0.69
7		7.01	58.70	0.85

7
0
7
0
5
6
2
1
6

$$(58.7 - 2.870528)(11)/(200) = 3.070621 \text{ g/L Carbon}$$

$$(58.7 - 2.870528)(11)/(200)(12) = .2558851 \text{ Molar Carbon}$$

Sample Run By: 80725 J.C.

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	WB39937
Procedure / Rev	LA-344-105/A-3
Technologist	E. Colvin
Date	06-29-90
Temperature	N/A
Starting Time	12:30
Ending Time	15:00
Chemist	R. Brandt

Total Organic Carbon from Water Digestion II.

Sample was acidified prior to analysis.

	Description	Lab. Id.
1	Reagent Blank	F5054
2	Initial LMCS Check Std	F5042
3	Sample 89-041	F5043
4	Duplicate 89-041	F5044
5	Spike 89-041	F5045
6	Final LMCS Check Std	F5046
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume of Standard
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std	70C11D /200uL			2.2 mL
Spike	80C11A /200uL	Sample/200 uL		.5 mL

Interim

4/04/90

Rev.E

SST-102

Prepared by: Shirley Cervantes S. A. Cervantes Date: 08/09/90
 Signature Printed Name

Verified by: C. M. Seidel C. M. Seidel Date: 08/09/90
 Signature Printed Name

Approved by: L. H. Taylor L. H. Taylor Date: 8/30/90
 Signature Printed Name

CULLOMETER ANALYSIS REPORT
FACTORY Rev. 0

Sample# F-50054 Date: October 1990 Time: 12:43:26

Blank = N/A Sample Size = 200 Dilution Factor = 1
% Difference = 10 Min Reading = / Max Reading = /

#	Reading	Analysis Time	Culometer	% Difference
1	1.01	0.00	0.00	0.00
2	2.01	1.00	1.00	100.00
3	3.01	1.60	1.60	57.10
4	4.01	2.10	2.10	23.81
5	5.01	2.60	2.60	19.23
6	6.01	3.00	3.00	15.33
7	7.01	3.50	3.50	14.29

BLANK VALUE = 3.03 / 2,006693 = .4995253 mg/minute.

Sample Run By: 800548...

POOR COPY RECEIVED

COURCUMETER ANALYSIS REPORT
TCLDC Rev. 0

Sample# E-SO42 Date: 06-29-1990 Time: 12:51:19

Blank = .4995253 Sample Size = 200 Dilution Factor = 11
 % Difference = 10 Min Reading = 7 Max Reading = 7

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	35.10	100.00
3	3.01	47.80	26.57
4	4.01	52.70	9.30
5	5.01	54.90	4.01
6	6.01	55.20	2.31
7	7.01	57.10	1.53

$$(57.1 - 3.500001) / (11) / (200) = 2.948 \text{ g/L Carbon}$$

$$(57.1 - 3.500001) / (11) / (200) (12) = .2456687 \text{ Molar Carbon}$$

Sample Run By: 60028

"POOR COPY RECEIVED"

COULOMETRIC ANALYSIS REPORT
TTCDOC Rev. 0

Sample ID: 1-10104-3

Date: 08-29-1990

Time: 14:59:55

Blank = .4995253
% Difference = 10

Sample Size = 200 Dilution Factor = 1.1
Min Readings = 7 Max Readings = 7

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	1.90	100.00
3	3.01	2.80	82.14
4	4.01	3.50	20.00
5	5.01	4.10	14.63
6	6.01	4.50	8.89
7	7.01	5.00	10.00

$$\text{C}_5 = (3.499604 \times 1.1) / (200) = 8.2352178E-03 \text{ g/l Carbon}$$

$$\text{C}_5 = (3.499604 \times 1.1) / (200 \times 1.1) = 6.876816E-04 \text{ g/lar Carbon}$$

Sample Run By: 60028

POOR COPY RECEIVED

Coulometer Analysis Report
TIC/TUC Rev. 0

Sample# = P-100041 Dilution Factor = 0.6-2.9-1.9%v/v Time# = 13s Observed

Blank = -1.4995253 Sample Size = 200 Dilution Factor = 1.1
 % Difference = 10 Min Reading = -2 Max Reading = -2

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	2.60	100.00
3	3.01	3.70	29.73
4	4.01	4.30	13.95
5	5.01	4.90	12.24
6	6.01	5.40	9.09
7	7.01	5.90	8.47

$$\text{C}_5\text{H}_9 = 3.500001 \times 1.1 \times (200) = .0132 \text{ g/L Carbon}$$

$$\text{C}_5\text{H}_9 = 3.500001 \times 1.1 \times (200) \times 120 = .0011 \text{ Molal Carbon}$$

Sample Run By: 8002B

cPOOR COPY RECEIVED

COULOMETER ANALYSIS REPORT
TICLOC Rev. 0

Sample #: R-50045 Date: 06-29-1990 Time: 1:31:16:48

Blank = .4995253 Sample Size = 200 Dilution Factor = 1.1
 % Difference = 10 Min Readings = 7 Max Readings = 7

#	Reading	Analysis Time	Coulometer	% Difference
1	1.01		0.00	0.00
2	2.01	80.10	100.00	
3	3.01	108.10	125.90	
4	4.01	118.70	135.93	
5	5.01	122.80	134.84	
6	6.01	124.70	134.52	
7	7.01	126.00	134.03	

$$\sqrt{126} = 3.500001 \times (1.1) / (200) = .67375 \text{ g/L Carbon}$$

$$(126 - 3.500001) \times (1.1) / (200) / (12) = 5.614583E-02 \text{ Molar Carbon}$$

Sample Run By: 60028

'POOR COPY RECEIVED'

COLORIMETER ANALYSIS REPORT
TICTOC Rev. 0

Sampler: R-5046 Dates: 06-29-1990 Liner: 1.3 ± 2% ± 0%

Batch #: 4995253 Sample Size = 200 Dilution Factor = 11
 % Difference = 10 Min Readings = 7 Max Readings = 7

	Reading	Analysis Time	Colorimeter	% Difference
1	1.01	0.00	0.00	0.00
2	2.01	36.40	100.00	
3	3.01	48.20	24.48	
4	4.01	52.70	8.54	
5	5.01	54.60	3.48	
6	6.01	55.80	2.15	
7	7.01	56.60	1.41	

$$C_{\text{MS,6}} = 3.500031 \times (11) / (200) = 2.920498 \text{ g/l Carbon}$$

$$(S_{\text{MS,6}} = 3.500031 \times (11) / (200) \times 12) = .2433749 \text{ Molar Carbon}$$

Sample Run By: 80028

"POOR COPY RECEIVED"

ACID DIGESTION TEST ANALYSIS

ICP Results Data Summary

Date Analyzed:	February 22, 1990	Acid Digested Standard	F5047
Procedure:	LA-505-151/A-0	Sample 89-041	F5048
Analyst:	D. M. Southwick	Duplicate of 89-041	F5049
Digestion	Acid Digestion	Spike of 89-041	F5050
Procedure:	LA-505-159/A-0		

	Instrument Starting Standard %	Acid Digest. LMCS Standard %	Wet Weight Sample ug/g	Wet Weight Sample Duplicate ug/g	Spike Recovery %	Closing LMCS Standard %
Aluminum	103.92%		34730	37858	NOT CALC.	109.27%
Barium	97.78%		22	16	94.30%	102.20%
Bismuth	100.00%	92.25%	20942	19470	NOT CALC.	96.06%
Boron	107.54%	101.60%	5 LT	-8 LT	53.64%	103.98%
Cadmium	103.12%	96.85%	3 LT	-3 LT	76.13%	99.23%
Calcium	104.97%	113.20%	155	152	7.28%	102.15%
Chromium	93.60%		446	429	53.42%	99.46%
Copper	104.29%	104.65%	27	15	75.98%	101.92%
Iron	96.80%		7927	8145	NOT CALC.	103.53%
Lanthanum	93.57%	90.75%	10 LT	-15 LT	55.81%	90.76%
Lead	93.02%	83.85%	131	-40 LT	87.36%	90.04%
Lithium	104.96%		4	-2 LT	101.89%	107.85%
Magnesium	103.54%	102.20%	103	116	99.67%	99.90%
Manganese	97.11%		1993	2048	58.27%	103.56%
Molybdenum	105.61%	106.45%	20	11	80.73%	114.90%
Nickel	96.54%		59	44	64.14%	103.22%
Phosphorous	108.41%	108.45%	9243	9534	NOT CALC.	101.74%
Potassium	99.29%	90.15%	198	-53 LT	43.29%	94.82%
Silver	105.57%	101.15%	-9 LT	-21 LT	NOT CALC.	103.69%
Sodium	99.57%	110.45%	65496	65413	NOT CALC.	97.18%
Strontium	104.12%	102.45%	243	245	94.41%	101.31%
Tantalum	100.59%		63	13 LT	-3.48%	106.69%
Tin	90.31%		31	17	47.32%	97.19%
Titanium	104.73%		11	2 LT	77.34%	110.47%
Zinc	102.20%	100.00%	57	50	94.70%	99.54%
Zirconium	102.37%		44	7 LT	-5.36%	108.13%

LT: Less Than

NC: Not Calibrated

NOT CALC: Not Calculated

Instrument Standards Outside Control Limits

Prepared by:

H. S. Rich

Date: May 2, 1990

Verified by:

C. M. Seidel

Date: May 2, 1990

Approved by:

L. H. Taylor

Date: May 9, 1990

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Acid Digestion

Instrument	N/A
Procedure / Rev	LA-505-159/A-0
Technologist	D.M. Southwick/69769
Date	12-6-89
Temperature	72 F
Starting Time	08:20
Ending Time	14:00
Chemist	S. A. Jones

Note: A reagent blank was not prepared.

*LMCS Check Standard was digested, but not used in the analysis.

	Description	Lab. Id.
1	LMCS Check Standard	F5047
2	Sample of 89-041	F5048
3	Duplicate of 89-041	F5049
4	Spike of 89-041	F5050
5	LMCS Check Standard	F5051*
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	of Standard
LMCS Check Std.	34C11CM/5mL			25 mL
Spike	103C15C/5mL	104C15D/5mL		25 mL

Interim

Rev E
4/04/90

SST-102

Prepared by: H.S. Rich H.S. Rich Date: May 04, 1990
Signature Printed Name

Verified by: C.M. Seidel C.M. Seidel Date: May 04, 1990
Signature Printed Name

Approved by: L.H. Taylor L.H. Taylor Date: May 9, 1990
Signature Printed Name

Analytical Batch

Lab Segment Serial No.: F5033

Customer ID.: 89-041

Instrument	WB39939
Procedure / Rev	LA-505-151/A-0
Technologist	D. M. Southwick
Date	Feb. 22, 1990
Temperature	70 F
Starting Time	16:00
Ending Time	21:30
Chemist	S. A. Jones

ICP analysis of sample 89-041.
No reagent blank was analyzed in this batch.

	Description	Lab. Id.
1	Initial LMCS Check Std.	na
2	Digested Standard 34C11C	F5047
3	Sample 89-041	F5048
4	Duplicate of 89-041	F5049
5	Spike of 89-041	F5050
6	Ending LMCS Check Std.	na
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	of Standard
LMCS Check Std.	81C11A/1.0 mL	82C11A/1.0 mL		11.0 mL
Digested LMCS Std.	34C11C/5.0 mL			50.0 mL
Spike of Sample	103C15C/5.mL	104C15D/5.0mL	Sample/.5838g	25.0 mL

Interim

4:04:90

Rev E

SST-102

Prepared by:	H. S. Rich Printed Name	Date: May 4, 1990
Verified by:	C. M. Seidel Printed Name	Date: May 4, 1990
Approved by:	L. H. Taylor Printed Name	Date: May 9, 1990

ICP Results Raw Data Summary

Date Analyzed:	February 22, 1990	Acid Digested Standard	F5047
Procedure:	LA-505-151/A-0	Reagent Blank	F5055
Analyst:	D. M. Southwick	Sample 89-041	F5048
Digestion	Acid Digestion	Duplicate of 89-041	F5049
Procedure:	LA-505-159/A-0	Spike of 89-041	F5050
		Acid Digested Standard	F5051

	Instrument Starting LMCS Standard %	Acid Digest. LMCS Standard %	Reagent BLANK ppm	Wet Weight Sample ug/g	Wet Weight Sample Duplicate ug/g	Spike Recovery %	LMCS ACID Digestion %	Closing LMCS Standard %
Aluminum	103.92%		NOT RUN	34730	37858	NOT CALC.		109.27%
Antimony				73	-2 LT		Not Run	
Arsenic				47	26			
Barium	97.78%			22	16	94.30%		102.20%
Beryllium				1	0 LT			
Bismuth	100.00%	92.25%		20942	19470	NOT CALC.		96.06%
Boron	107.54%	101.60%		5 LT	-8 LT	53.64%		103.98%
Cadmium	103.12%	96.85%		3 LT	-3 LT	76.13%		99.23%
Calcium	104.97%	113.20%		155	152	7.28%		102.15%
Cerium	86.74%			203	10 LT	1.66%		91.59%
Chromium	93.60%			446	429	53.42%		99.46%
Cobalt	78.65%			-76 LT	-96 LT	NOT CALC.		89.73%
Copper	104.29%	104.65%		27	15	75.98%		101.92%
Europium				3	-1 LT			
Iron	96.80%			7927	8145	NOT CALC.		103.53%
Lanthanum	93.57%	90.75%		10 LT	-15 LT	55.81%		90.76%
Lead	93.02%	83.85%		131	-40 LT	87.36%		90.04%
Lithium	104.96%			4	-2 LT	101.89%		107.85%
Magnesium	103.54%	102.20%		103	116	99.67%		99.90%
Manganese	97.11%			1993	2048	58.27%		103.56%
Mercury				18	8			
Molybdenum	105.61%	106.45%		20	11	80.73%		114.90%
Neodymium	88.30%			84 LT	-44 LT	68.73%		91.77%
Nickel	96.54%			59	44	64.14%		103.22%
Phosphorus	108.41%	108.45%		9243	9534	NOT CALC.		101.74%
Potassium	99.29%	90.15%		198	-53 LT	43.29%		94.82%
Samarium				119	-62 LT			
Selenium				448	337			
Silicon	78.22%	81.60%		707	961	42.33%		83.71%
Silver	105.57%	101.15%		-9 LT	-21 LT	NOT CALC.		103.69%
Sodium	99.57%	110.45%		65496	65413	NOT CALC.		97.18%
Strontium	104.12%	102.45%		243	245	94.41%		101.31%
Sulfur				921	707			
Tantalum	100.59%			63	13 LT	3.48%		106.69%
Thallium				266	16 LT			
Thorium				76	-62 LT			
Tin	90.31%			31	17	47.32%		97.19%
Titanium	104.73%			11	2 LT	77.34%		110.47%
Tungsten				100	58			
Uranium				2116	1025			
Vanadium				25	3 LT			
Zinc	102.20%	100.00%		57	50	94.70%		99.54%
Zirconium	102.37%			44	7 LT	-5.36%		108.13%

LT: Less Than

NC: Not Calibrated

NOT CALC: Not Calculated

Instrument Standards Outside Control Limits

ICP Results Raw Data

Date Analyzed: February 22, 1990
 Procedure: LA-505-151/A-0
 Analyst: D. M. Southwick
 Digestion: Acid Digestion
 Procedure: LA-505-159/A-0

Acid Digested Standard F5047
 Reagent Blank F5055
 Sample 89-041 F5048
 Duplicate of 89-041 F5049
 Spike of 89-041 F5050
 Acid Digested Standard F5051

					Digestion Weight	0.02173	Digestion Weight		
					Sample	0.5432 g	Volume		
	Starting LMCS Standard		LMCS Acid Digestion	Reagent Blank	Sample	25.00 mL	Sample		
	Instrument Standard	Recovery	Digestion Standard	Recovery	Dilution	Dilution	Duplicate Sample Dilution		
	ppm	%	ppm	%	Three ppm	Two ppm	One ppm		
Aluminum	103.92	103.92%			NOT RUN	754.61	765.73		
Antimony					0.00	1.58			
Arsenic					-0.70	1.03			
Barium	97.78	97.78%			0.02	0.48			
Beryllium					-0.03	0.01			
Bismuth	100.00	100.00%	9.23	92.25%	455.03	470.49			
Boron	107.54	107.54%	10.16	101.60%	-1.40	0.11 LT			
Cadmium	103.12	103.12%	9.69	96.85%	-1.05	0.06 LT			
Calcium	104.97	104.97%	11.32	113.20%	-6.00	3.37			
Cerium	86.74	86.74% #			-10.48	4.40			
Chromium	94.45	93.60%			9.69	11.99			
Gerbalt	78.65	78.65% #			-6.70	-1.64 LT			
Copper	104.29	104.29%	10.47	104.65%	-0.53	0.59			
Europium					-0.24	0.06			
Iron	96.80	96.80%			172.23	183.77			
Yttrium	93.57	93.57%	9.08	90.75%	-0.91	0.22 LT			
Danthanum					-3.74	2.85			
Lead	93.02	93.02%	8.39	83.85%	-0.38	0.09			
Lithium	104.96	104.96%			2.25	2.79			
Magnesium	103.54	103.54%	10.22	102.20%	43.31	45.06			
Manganese	97.11	97.11%			-0.47	0.38			
Mercury					-0.29	0.43			
Molybdenum	105.40	105.61%	10.65	106.45%	-1.47	1.83 LT			
Neodymium	88.30	88.30% #			-0.12	1.27			
Nickel	96.54	96.54%			200.83	209.40			
Phosphorous	108.41	108.41%	10.85	108.45%	-6.85	4.30			
Potassium	99.29	99.29%	9.02	90.15%	-13.63	2.58			
Samarium					-4.73	9.72			
Selenium					12.93	15.37			
Silicon	78.22	78.22% #	8.16	81.60%	-5.20	-0.20 LT			
Silver	105.57	105.57%	10.12	101.15%	1423.10	1403.50			
Sodium	99.57	99.57%	11.05	110.45%	5.27	5.49			
Strontium	104.12	104.12%	10.25	102.45%	20.02	18.73			
Sulfur					-3.63	1.38			
Tantalum	100.09	100.59%			-20.50	5.78			
Thallium					-11.30	1.65			
Thorium					-2.19	0.67			
Tin	90.31	90.31%			-0.58	0.24			
Titanium	104.83	104.73%			-1.64	2.17			
Tungsten					-56.10	45.98			
Uranium					-1.28	0.54			
Vanadium					0.50	1.24			
Zinc	102.20	102.20%	10.00	100.00%	-1.11	0.96			
Zirconium	102.17	102.37%							
Dilution Factor	11.00		10.00		1.00	1.00	101.00	11.00	1.00

LT: Less Than

NC: Not Calibrated

NOT CALC: Not Calculated

Instrument Standards Outside Control Limits

	0.02369	Digestion	0.02335								
	0.5922 g 25.00 mL	Weight Volume	0.5838 g 25.00 mL	Spike of Sample	Spike of Sample	Spike of Sample	Standard Recovery	Standard LMCS Acid	Acid Digestion Standard	Ending LMCS Standard	Standard Recovery
	Duplicate Dilution	Duplicate Dilution	Dilution	Three ppm	Dilution	Dilution	%	Digestion ppm	Recovery %	ppm	%
	Two ppm	One ppm			Two ppm	One ppm					
Aluminum	896.77	883.47			818.09	880.31	NOT CALC.			109.27	109.27%
Antimony	-8.00	-0.05 LT			1.00	11.06					
Arsenic	1.64	0.61			-3.82	1.07			NOT RUN		
Barium	0.25	0.39			9.45	10.89	94.30%			102.20	102.20%
Beryllium	-0.03	0.00 LT			-0.07	0.00 LT					
Bismuth	461.21	447.03			427.80	469.40	NOT CALC.			96.06	96.06%
Boron	-1.94	-0.20 LT			5.48	9.93	53.64%			103.98	103.98%
Cadmium	-0.92	-0.06 LT			7.68	10.02	76.13%			99.23	99.23%
Calcium	-5.76	3.60			4.35	14.70	7.28%			102.15	102.15%
Cerium	-2.96	0.23 LT			-28.65	0.17 LT	1.66%			91.59	91.59%
Chromium	10.16	11.78			15.76	22.03	53.42%			100.36	99.46%
Cobalt	-19.26	-2.28 LT			1.68	7.30	NOT CALC.			89.73	89.73%
Copper	0.16	0.35			8.23	11.00	75.98%			101.92	101.92%
Europium	-0.06	-0.02 LT			-0.58	0.00 LT					
Iron	192.94	184.98			178.96	195.55	NOT CALC.			103.53	103.53%
Lanthanum	-2.73	-0.36 LT			5.82	9.52	55.81%			90.76	90.76%
Lead	-31.13	-0.95 LT			-24.90	11.80	87.36%			90.04	90.04%
Lithium	-0.30	-0.06 LT			10.29	11.68	101.89%			107.85	107.85%
Magnesium	2.74	2.89			12.38	13.76	99.67%			99.90	99.90%
Manganese	48.50	46.23			52.37	56.34	58.27%			103.56	103.56%
Mercury	1.08	0.19			-1.55	0.17					
Molybdenum	0.01	0.26			8.54	11.23	80.73%			114.67	114.90%
Neodymium	2.24	-1.05 LT			-10.43	8.84	68.73%			91.77	91.77%
Nickel	1.42	1.04			7.78	11.15	64.14%			103.22	103.22%
Phosphorus	225.84	208.40			227.65	251.76	NOT CALC.			101.74	101.74%
Potassium	-12.08	-1.25 LT			-28.40	8.95	43.29%			94.82	94.82%
Samarium	-6.69	-1.48 LT			-31.96	-1.77 LT					
Selenium	3.63	7.98			-10.54	12.31					
Silicon	29.42	22.77			18.13	18.50	42.33%			83.71	83.71%
Silver	-4.39	-0.51 LT			3.61	10.28	NOT CALC.			103.69	103.69%
Sodium	1549.50	1534.90			1558.70	1642.10	NOT CALC.			97.18	97.18%
Strontium	5.81	5.67			15.11	16.36	94.41%			101.31	101.31%
Sulfur	16.76	19.02			16.18	21.26					
Tantalum	0.03	0.31 LT			-8.67	0.35 LT	3.48%				
Thallium	-6.68	0.37 LT			-58.85	-0.58 LT					
Thorium	-8.81	-1.48 LT			-21.01	-0.86 LT					
Tin	-0.53	0.40			5.45	9.93	47.32%			97.19	97.19%
Titanium	-0.13	0.05 LT			7.99	10.33	77.34%			110.58	110.47%
Tungsten	1.44	1.37			-4.78	1.56					
Uranium	-0.97	24.28			-166.80	24.81					
Vanadium	-0.39	0.06 LT			-3.73	0.06 LT					
Zinc	0.59	1.19			10.01	11.72	94.70%			99.54	99.54%
Zirconium	0.14	0.16 LT			-3.53	0.50	-5.36%			107.91	108.13%
Dilution Factor	101.00	11.00	1.00	101.00	11.00			10.00		11.00	

	Spike Standard	Spike Standard	LMCS Standards	LMCS Standard	ACID DIGESTION	ACID DIGEST.
	LMCS	ID	Values	IDs	LMCS STANDARD VALUES	LMCS IDs Book
	ppm added	5.0 mL	ppm	Book	ppm in Sample	#
Aluminum	10.00	5.0 mL	100.00			
Antimony		104C15D				
Arsenic						
Barium	10.00		100.00			
Beryllium						
Bismuth	10.00		100.00		100.00	
Boron	10.00		100.00		100.00	
Cadmium	10.00		100.00		100.00	
Calcium	10.00		100.00		100.00	
Cerium	10.00		100.00			
Chromium	10.00		100.90			
Cobalt #	10.00		100.00			
Copper	10.00		100.00		100.00	
Europium						
Iron	10.00		100.00			
Yttrium	10.00		100.00		100.00	
Lead	10.00		100.00		100.00	
Lithium	10.00		100.00			
Magnesium	10.00		100.00		100.00	
Manganese	10.00		100.00			
Mercury						
Molybdenum #	10.00		99.80		99.80	
Neodymium	10.00		100.00			
Nickel	10.00		100.00			
Phosphorus	10.00		100.00		100.00	
Potassium	10.00		100.00		100.00	
Samarium						
Selenium						
Silicon #	10.00		100.00		100.00	
Silver	10.00		100.00		100.00	
Sodium	10.00		100.00		100.00	
Strontium	10.00		100.00		100.00	
Sulfur						
Tantalum	9.95		99.50			
Thallium						
Thorium						
Tin	10.00		100.00			
Titanium #	10.00		100.10			
Tungsten						
Uranium						
Vanadium						
Zinc	10.00		100.00		100.00	
Zirconium	9.98		99.80			
Dilution Factor						

ICP Calibration Report

Procedure: LA-505-151 Revision: A-0
Instrument: WB39939
Technologist: D.M. Southwick
Date: Feb. 22, 1990
Time: 16:28

Calibration Standards for ICP Program "SST"

Element	Standard	Element	Standard
Aluminum	SST-3	Antimony	SST-4
Arsenic	SST-4	Barium	SST-2
Beryllium	SST-2	Bismuth	SST-4
Boron	SST-3	Cadmium	SST-2
Calcium	SST-2	Cerium	SST-5
Chromium	SST-2	Cobalt	SST-2
Copper	SST-2	Europium	SST-5
Iron	SST-2	Lanthanum	SST-5
Lead	SST-4	Lithium	SST-1
Magnesium	SST-2	Manganese	SST-2
Mercury	SST-3	Molybdenum	SST-3
Neodymium	SST-5	Nickel	SST-2
Phosphorous	SST-3	Potassium	SST-1
Samarium	SST-5	Selenium	SST-4
Silicon	SST-3	Silver	SST-2
Sodium	SST-1	Strontium	SST-2
Sulfur	SST-3	Tantalum	SST-3
Thallium	SST-4	Thorium	SST-4
Tin	SST-4	Titanium	SST-3
Tungsten	SST-3	Uranium	SST-4
Vanadium	SST-2	Zinc	SST-2
Zirconium	SST-3		

ICP Standard Formulations

SST-0:
Calibration blank, 1 M ultrex HNO3.

SST-1:

Stock solutions from AESAR/John Mathey Inc., Seabrook, NH 03874.
Individual element solutions as follows:

Li LiCO₃ 10,000 ppm in 5% HNO₃ Lot# 14394A

K KNO₃ 10,000 ppm in 5% HNO₃ Lot# 14379A

Na NaCO₃ 10,000 ppm in 5% HNO₃ Lot# 14400A

200 mL of standard made by combining 25 mL HCl/HNO₃ mixed acid, 1 mL each single element standards, and water.

SST-2:

Stock solutions from VHG labs, Inc., 180 Zachary Rd. #5, Manchester, NH 03103. Mixed element standards as follows:

SM-10 Li, Na, K, Rb, Cs, Be, Mg, Ca, Sr, & Ba 100 ppm
Lot# 0-119A

SM-20 V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ag, & Cd 100 ppm
Lot# 0-119B

50 mL of each mixed standard are added to a 250 mL volumetric flask and diluted to volume with 1 M HNO₃.

SST-3:

Stock solutions from AESAR/John Mathey Inc., Seabrook, NH 03874.
Individual element solutions as follows:

Al Al 10,000 ppm in 10% HCl Lot# 9-053A

B H₃BO₃ 10,000 ppm in 1% NH₄OH Lot# 9-335A

Hg Hg 10,000 ppm in 5% HNO₃ Lot# 8-656S

Mo Mo 10,000 ppm in 5% HCl Lot# 9-159T

P P 10,000 ppm in 5% HNO₃ Lot# 9-160A

Si Si 1000 ppm in KOH Lot# 086DM Spex Industries, Edison, NJ

S (NH₄)₂SO₄ in H₂O Lot# 9-231M

Ta TaCl₅ 10,000 ppm in 5% HCl/tr HF Lot# 9-335M

Ti Ti 10,000 ppm in 5% HF Lot# 9-079EE

W W 10,000 ppm in 5% HF/tr HNO₃ Lot# 8-685L

Zr ZrCl₂O 10,100 ppm in 5% HCl Lot# 9-078G

50 mL of each mixed standard are added to a 250 mL volumetric flask and diluted to volume with 1 M HNO₃.

SST-4:

Stock solution from VHG labs, Inc., 180 Zachary Rd. #5, Manchester, NH 03103. Mixed element standard as follows:

SM-50 Ga, In, Tl, Ge, Sn, Pb, As, Sb, Bi, Se, Te, Th, & U 100 ppm Lot# 0-119D

Solution is used directly for calibration.

SST-5:

Stock solution from VHG labs, Inc., 180 Zachary Rd. #5, Manchester, NH 03103. Mixed element standard as follows:

SM-60 Sc, Y, La, Ce, Pr, Nd, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, & Lu 100 ppm Lot# 7-165F

50 mL of SM-60 is added to a 250 mL volumetric flask and diluted to volume with 1 M HNO₃.

Sample name : SST0
Programme : SST 22-Feb-90 16:05:28

NAME	MV	INT	RSD
AL	1.34	0.74	
CB	0.33	0.35	
AS	0.72	0.93	
BA	2.57	0.82	
BE	0.51	0.60	
BT	2.48	1.49	
B-	3.07	1.30	
CD	1.52	1.03	
CA	2.20	0.98	
CR	3.51	0.69	
CR	1.01	0.31	
CO	0.23	0.43	
CU	2.00	0.77	
EU	2.71	0.89	
FE	1.14	0.82	
LA	0.32	0.18	
PB	0.25	0.00	
LT	2.77	0.63	
AG	0.48	0.46	
BN	0.56	0.31	
HG	2.52	0.94	
MO	1.13	0.37	
ND	3.92	1.63	
NI	2.25	0.53	
P	1.00	2.36	
K	2.46	0.44	
CM	3.45	0.76	
SE	1.16	1.04	
SI	2.19	0.54	
OG	9.52	0.77	
NA	3.43	0.63	
SR	2.68	0.70	
S	0.53	1.04	
TA	2.34	1.06	
TL	2.73	0.50	
TH	1.04	0.64	
CK	0.86	0.56	
TI	2.42	0.64	
U	0.89	0.67	
C	3.31	0.78	
VI	2.99	0.24	
ZN	1.66	0.47	
ZR	3.44	0.50	

Sample name : SST1
Programme : SST 22-Feb-90 16:09:21

NAME	MV	INT	RSD
LI	411.60	0.42	
K	10.82	0.51	
NA	47.15	0.43	

Programme : S9T3

NAME MV INT RSD

BA	227.52	0.94	CD	40.19	0.95	CA	310.01	0.93	CB	150.82	1.17	CD	150.82	0.93	CE	421.95	0.93	DE	421.95	0.93	EE	227.52	0.94	NAME	MV	INT	RSD
----	--------	------	----	-------	------	----	--------	------	----	--------	------	----	--------	------	----	--------	------	----	--------	------	----	--------	------	------	----	-----	-----

Programme : S9T3

NAME MV INT RSD

AL	17.25	2.95	AL	534.70	2.85	HG	383.52	9.48	HO	174.05	9.26	HI	52.47	8.97	IT	67.59	8.73	JI	30.20	8.83	JK	79.02	9.33	LM	344.49	9.56	LN	135.19	8.30	NAME	MV	INT	RSD
----	-------	------	----	--------	------	----	--------	------	----	--------	------	----	-------	------	----	-------	------	----	-------	------	----	-------	------	----	--------	------	----	--------	------	------	----	-----	-----

Programme : S9T4

NAME MV INT RSD

AS	7.06	2.54	AS	82.74	2.94	BS	98.03	2.55	PS	2.95	2.77	BS	28.95	2.58	ES	1.95	2.77	IS	1.95	2.77	NAME	MV	INT	RSD
----	------	------	----	-------	------	----	-------	------	----	------	------	----	-------	------	----	------	------	----	------	------	------	----	-----	-----

Programme : S9T5

NAME MV INT RSD

NAHIL	22-FEB-90	16:19:17	NAME	MV	INT	RSD																		
-------	-----------	----------	-------	-----------	----------	-------	-----------	----------	-------	-----------	----------	-------	-----------	----------	-------	-----------	----------	-------	-----------	----------	------	----	-----	-----

Sample name : SST5
Programme : SST 22-Feb-90 16:19:17

NAME	MV	INT	RSD
CE	11.93	1.44	
HJ	357.24	2.01	
LA	4.03	1.84	
ND	23.16	1.62	
SM	11.61	1.43	

Programme name : SST Channel name : AL Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	1.2768	18.115	-0.422429E+01	0.314307E+01		

Programme name : SST Channel name : SB1 Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	0.3091	7.4119	-0.483192E+01	0.148522E+02		

Programme name : SST Channel name : AS Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	0.6802	86.877	-0.672919E+00	0.121916E+01		

Programme name : SST Channel name : BA Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	2.4399	238.90	-0.228345E+00	0.889080E-01		

Programme name : SST Channel name : BE1 Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	0.4848	442.95	-0.242241E-01	0.424673E-01		

Programme name : SST Channel name : BI Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0		C1	C2		C3
CRV1	128.5620	92.115	-0.295602E+01	0.211707E+01		

CRV1 2.3538 92.491 -0.289610E+01 0.1116886E+01

Programme name : SST Channel name : B Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 2.9178 561.44 -0.288859E+00 0.940502E-01

Programme name : SST Channel name : CD Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 1.4434 158.36 -0.203531E+00 0.133961E+00

O

—

Programme name : SST Channel name : CA Polynomial type : CC

O

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 2.0925 325.51 -0.143119E+00 0.649749E-01

Q

Programme name : SST Channel name : CE Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 3.3358 12.522 -0.834575E+01 0.237680E+01

Programme name : SST Channel name : CR Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 0.9560 42.201 -0.513632E+00 0.510399E+00

Programme name : SST Channel name : CO Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 129.0 129.0 129.0

Process Data Output

CRV1 0.2185 1.9303 -0.286010E+01 0.124353E+02

Programme name : SST Channel name : CU Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 1.9041 81.555 -0.529777E+00 0.264316E+00

Programme name : SST Channel name : EU Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 2.5713 375.10 -0.152690E+00 0.564125E-01

Programme name : SST Channel name : FE Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 1.0798 68.763 -0.353265E+00 0.310791E+00

Programme name : SST Channel name : LA Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 0.3037 4.2263 -0.172544E+01 0.539763E+01

Programme name : SST Channel name : PB Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1 0.2366 3.1007 -0.920856E+01 0.369823E+02

Programme name : SST Channel name : LI Polynomial type : CC

Programme name : SST Channel name : B Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	2.9178	561.44	-0.288859E+00	0.940502E-01			

Programme name : SST Channel name : CD Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	1.4434	158.36	-0.203531E+00	0.133961E+00			

Programme name : SST Channel name : CA Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	2.0925	325.51	-0.143119E+00	0.649749E-01			

Programme name : SST Channel name : CE Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	3.3358	12.522	-0.834575E+01	0.237680E+01			

Programme name : SST Channel name : CR Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	0.9560	42.201	-0.513632E+00	0.510399E+00			

Programme name : SST Channel name : CO Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients	C0	C1	C2	C3
CRV1	0.13185	1.9303	-0.286019E+01	0.124252E-01			

Programme name : SST Channel name : CU Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefficients
C1 C2 C3
CRV1 1.9041 81.555 -0.529777E+00 0.264316E+00

Programme name : SST Channel name : EU Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefficients
C1 C2 C3
CRV1 2.5713 375.10 -0.152690E+00 0.564125E-01

CV

Programme name : SST Channel name : FE Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefricients
C1 C2 C3
CRV1 1.0798 68.763 -0.353265E+00 0.310791E+00

CF

Programme name : SST Channel name : LA Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefficients
C1 C2 C3
CRV1 0.3037 4.2263 -0.172544E+01 0.539763E+01

Programme name : SST Channel name : PB Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefficients
C1 C2 C3
CRV1 0.2366 3.1007 -0.920858E+01 0.369823E+02

Programme name : SST Channel name : LI Polynomial type : CC

Curve Min Int Max Int
CO Curve Coefficients
C1 C2 C3
CRV1 1.32...

Programme name : SST Channel name : MG Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	0.4126	316.82	-0.286306E-01	0.663789E-01	

Programme name : SST Channel name : MN Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	0.5330	172.97	-0.683423E-01	0.121823E+00	

3

4 Programme name : SST Channel name : HG Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	2.3959	402.70	-0.330974E+00	0.131235E+00	

5

6 Programme name : SST Channel name : MD Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	1.0757	183.38	-0.326289E+00	0.288157E+00	

Programme name : SST Channel name : ND Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	3.7246	24.319	-0.407554E+01	0.103950E+01	

Programme name : SST Channel name : NI Polynomial type : CC

Curve	Min Int	Max Int	CO	Curve Coefficients	C3
				C1	C2
CRV1	133	2.3242	0.000000E+00	0.000000E+00	

Programme name : SST Channel name : P Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	0.9459	55.093	-0.967154E+00	0.971364E+00	
------	--------	--------	---------------	--------------	--

Programme name : SST Channel name : K Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.3414	11.365	-0.147432E+02	0.598182E+01	
------	--------	--------	---------------	--------------	--

4

Programme name : SST Channel name : SM Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	3.2819	12.189	-0.847356E+01	0.245279E+01	
------	--------	--------	---------------	--------------	--

5

Programme name : SST Channel name : SE Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	1.1007	30.396	-0.416942E+01	0.359946E+01	
------	--------	--------	---------------	--------------	--

Programme name : SST Channel name : SI Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.0649	70.968	-0.167803E+01	0.764596E+00	
------	--------	--------	---------------	--------------	--

Programme name : SST Channel name : AG Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	9.0437	337.22	-0.605908E+00	0.324774E+00	
------	--------	--------	---------------	--------------	--

Programme name : SST Channel name : NA Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	3.3028	49.506	-0.398040E+01	0.114489E+01			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : SR Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.5498	439.59	-0.129049E+00	0.480605E-01			
------	--------	--------	---------------	--------------	--	--	--

L0

Programme name : SST Channel name : S Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	0.5516	31.710	-0.980227E+00	0.168811E+01			
------	--------	--------	---------------	--------------	--	--	--

L0

Programme name : SST Channel name : TA Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.2205	82.970	-0.152405E+01	0.652046E+00			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : TL2 Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.5916	25.632	-0.125813E+02	0.461191E+01			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : TH Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

135							
-----	--	--	--	--	--	--	--

Programme name : SST Channel name : SN Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

CRV1	0.8211	128.67	-0.710316E+00	0.821807E+00		
------	--------	--------	---------------	--------------	--	--

Programme name : SST Channel name : TI Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

CRV1	2.2946	361.71	-0.353046E+00	0.146169E+00		
------	--------	--------	---------------	--------------	--	--

6

Programme name : SST Channel name : W Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

CRV1	0.8474	41.128	-0.116518E+01	0.120626E+01		
------	--------	--------	---------------	--------------	--	--

7

Programme name : SST Channel name : U Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

CRV1	3.1483	10.790	-0.476013E+02	0.143637E+02		
------	--------	--------	---------------	--------------	--	--

Programme name : SST Channel name : V1 Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

CRV1	2.8383	73.530	-0.891291E+00	0.298324E+00		
------	--------	--------	---------------	--------------	--	--

Programme name : SST Channel name : ZN Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
C1 C2 C3

Programme name : SST Channel name : ZR Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients			
					C0	C1	C2	C3
CRV1	3.2680	141.94	-0.130561E+01	0.37953BE+00				

Sample name : SST3
Programme : SST 22-Feb-90 16:28:22

NAME	MV	INT	RSD
AL	15.96	0.30	
B	486.20	0.39	
HG	338.57	0.19	
MO	154.52	0.24	
P	47.31	0.46	
SI	62.10	0.33	
S	26.90	0.89	
TA	70.07	0.13	
TI	309.66	0.46	
W	34.70	0.82	
ZR	122.07	0.48	

Programme name : SST Channel name : AL Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients			
					C0	C1	C2	C3
CRV1	1.2768	16.753	-0.459917E+01	0.342200E+01				

Programme name : SST Channel name : SBl Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients			
					C0	C1	C2	C3
CRV1	0.3091	7.4113	-0.483192E+01	0.148522E+02				

Programme name : SST Channel name : AS Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients			
					C0	C1	C2	C3
CRV1	0.6002	86.877	-0.672919E+00	0.121916E+01				

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.4399	238.90	-0.228345E+00	0.889080E-01		

Programme name : SST Channel name : BE1 Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.4848	442.95	-0.242241E-01	0.474673E-01		

Programme name : SST Channel name : BI Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.3538	92.431	-0.289610E+01	0.116888E+01		

Programme name : SST Channel name : B Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.9178	510.51	-0.317860E+00	0.103492E+00		

Programme name : SST Channel name : CD Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.4434	158.36	-0.203531E+00	0.133961E+00		

Programme name : SST Channel name : CA Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.0925	325.51	-0.143118E+00	0.649749E-01		

Programme name : SST Channel name : CE Polynomial type : UC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	3.3358	12.522	-0.834575E+01	0.237680E+01		

Programme name : SST Channel name : CR Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.9560	42.201	-0.513632E+00	0.510399E+00		

Programme name : SST Channel name : CO Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.2185	1.9303	-0.286010E+01	0.124352E+02		

Programme name : SST Channel name : CU Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.9041	81.555	-0.529777E+00	0.264316E+00		

Programme name : SST Channel name : EU Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.5713	375.10	-0.152690E+00	0.564125E-01		

Programme name : SST Channel name : FE Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.0798	68.763	-0.353265E+00	0.310791E+00		

Programme name : SST Channel name : LA Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	0.3037	4.2263	-0.172544E+01	0.539763E+01			

Programme name : SST Channel name : PB Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	0.2366	3.1007	-0.920858E+01	0.369823E+02			

Programme name : SST Channel name : LI Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	2.6337	432.18	-0.339059E+00	0.122301E+00			

Programme name : SST Channel name : MG Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	0.4126	316.92	-0.288306E-01	0.663789E-01			

Programme name : SST Channel name : MN Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	0.5330	172.97	-0.683423E-01	0.121822E+00			

Programme name : SST Channel name : HG Polynomial type : CC

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	2140.59	355.50	-0.375247E+06	0.1403789E+06			

Programme name : SST Channel name : MU Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 1.0757 162.24 -0.369118E+00 0.325980E+00 C2 C3

Programme name : SST Channel name : ND Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 3.7246 24.319 -0.407554E+01 0.103950E+01 C1 C2 C3

Programme name : SST Channel name : NI Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 2.1343 92.746 -0.521979E+00 0.232335E+00 C1 C2 C3

Programme name : SST Channel name : P Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 0.9459 49.679 -0.107482E+01 0.107950E+01 C1 C2 C3

Programme name : SST Channel name : K Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 2.3414 11.365 -0.147432E+02 0.598182E+01 C1 C2 C3

Programme name : SST Channel name : SM Polynomial type : CC

Curve Min Int Max Int
C0 Curve Coefficients
CRV1 3.2819 12.189 -0.847356E+01 0.245279E+01 C1 C2 C3

Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	1.1007	30.396	-0.416942E+01	0.359846E+01			
Programme name : SST	Channel name : SI	Polynomial type : CC					
Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	2.0849	65.200	-0.183192E+01	0.834715E+00			
Programme name : SST	Channel name : AG	Polynomial type : CC					
Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	9.0437	337.72	-0.609998E+00	0.640776E-01			
Programme name : SST	Channel name : NA	Polynomial type : CC					
Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	3.3028	49.506	-0.398040E+01	0.114489E+01			
Programme name : SST	Channel name : SR	Polynomial type : CC					
Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	2.5498	439.59	-0.129048E+00	0.480805E-01			
Programme name : SST	Channel name : S	Polynomial type : CC					
Curve	Min Int	Max Int	C0	Curve Coefficients	C1	C2	C3
CRV1	0.5516	28.249	-0.110296E+01	0.189948E+01			

Programme name : SST Channel name : TA Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.2205	73.578	-0.172529E+01	0.728145E+00			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : TL2 Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.5916	25.632	-0.125813E+02	0.461191E+01			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : TH Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	0.9896	16.741	-0.699011E+01	0.671051E+01			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : SN Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	0.8211	128.67	-0.710316E+00	0.821807E+00			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : TI Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	2.2946	325.15	-0.393059E+00	0.162273E+00			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : W Polynomial type : CC

Curve Min Int Max Int Curve Coefficients
C0 C1 C2 C3

CRV1	0.8474	36.431	-0.131936E+01	0.147910E+01			
------	--------	--------	---------------	--------------	--	--	--

Programme name : SST Channel name : V0 Polynomial type : CC
 Curve Min Int Max Int Curve Coefficients
 CRV1 3.1483 10.790 -0.476013E+02 0.143637E+02 C0 C1 C2 C3

Programme name : SST Channel name : V1 Polynomial type : CC
 Curve Min Int Max Int Curve Coefficients
 CRV1 2.8383 73.530 -0.891291E+00 0.298324E+00 C0 C1 C2 C3

Programme name : SST Channel name : ZN Polynomial type : CC
 Curve Min Int Max Int Curve Coefficients
 CRV1 1.5723 351.23 -0.994446E-01 0.600874E-01 C0 C1 C2 C3

Programme name : SST Channel name : ZR Polynomial type : CC
 Curve Min Int Max Int Curve Coefficients
 CRV1 3.2680 128.17 -0.144987E+01 0.421475E+00 C0 C1 C2 C3

Sample name : SIC11A
Programme : SST 22-Feb-90 16:43:59

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.30	-0.153	-1.631	-10.58	
Sb	0.32	-0.084	-0.926	-44.41	
As	0.89	0.212	2.333	11.74	
Ba	2.57	0.001	0.006	131.47	
Be	0.50	-0.000	-0.005	-32.08	
Bi	10.26	9.091	100.00	2.37	
B	97.53	9.776	107.54	1.70	
Cd	71.50	9.375	103.12	2.90	
Ca	149.07	9.542	104.97	2.37	
Ce	3.49	-0.047	-0.514	-60.66	
Cr	1.12	0.059	0.646	22.68	
Co	0.21	(-0.203	(-2.234	-12.74	
Cu	37.87	9.481	104.29	2.27	
Eu	2.68	-0.001	-0.016	-35.29	
Fe	1.11	-0.008	-0.093	-79.51	
La	1.90	8.507	93.573	2.63	
Pb	0.48	8.457	93.023	3.07	
Li	2.71	-0.007	-0.079	-18.12	
Mg	142.24	9.413	103.54	2.76	
Mn	0.59	0.004	0.045	22.52	
Hg	3.24	0.107	1.173	12.72	
Mo	1.12	-0.003	-0.032	-138.92	
Nd	4.25	0.345	3.792	22.09	
Ni	2.23	-0.004	-0.043	-127.70	
P	10.13	9.855	108.41	2.80	
K	3.97	9.027	99.292	2.88	
Sm	3.37	-0.210	-2.311	-9.55	
Se	1.14	-0.073	-0.805	-36.25	
Si	2.25	0.043	0.468	38.74	
Ag	159.29	9.597	105.57	2.17	
Na	11.38	9.052	99.571	2.85	
Sr	199.55	9.465	104.12	2.56	
C	0.63	0.099	1.087	15.54	
Ta	2.30	-0.030	-0.330	-39.89	
Tl	2.68	-0.209	-2.300	-32.29	
Th	1.01	-0.183	-2.018	-5.59	
Sn	0.98	0.098	1.079	7.92	
Ti	2.36	-0.009	-0.097	-22.40	
W	0.97	0.120	1.318	7.41	
U	3.29	-0.393	-4.319	-47.18	
V	2.95	-0.011	-0.123	-33.92	
Zn	156.28	9.291	102.20	2.58	
Zr	3.39	-0.020	-0.216	-18.22	

Dilution factor : 11.0000

*Not OK
has
/*

Sample name : SIC11A
Programme : SST 22-Feb-90 16:48:03

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.35	0.032	0.351	368.42	
Sb	0.31	-0.168	-1.852	-55.10	
As	0.74	0.028	0.313	128.74	
Ba	2.54	-0.002	-0.024	-313.64	
Be	0.50	-0.001	-0.006	-52.92	
Bi	2.43	-0.056	-0.617	-151.79	
B	3.64	0.059	0.647	65.45	
Cd	1.44	-0.010	-0.110	-47.36	
Ca	0.85	(-0.088	(-0.968	-0.60	
Ce	3.49	-0.075	-0.828	-350.92	
Cr	1.02	0.009	0.103	232.09	
Co	0.21	(-0.191	(-2.097	-13.58	
Cu	2.02	0.004	0.045	408.54	
Eu	2.68	-0.002	-0.018	-282.49	
Fe	1.10	-0.010	-0.113	-81.39	
La	0.31	-0.036	-0.396	-68.74	
Pb	0.24	-0.419	-4.610	-5.09	
Li	2.75	-0.003	-0.035	-258.76	
Mg	0.53	0.006	0.071	11.81	
Mn	0.57	0.002	0.018	85.75	
Hg	2.86	0.050	0.553	40.73	
Ho	30.53	9.582	105.40	2.25	
Nd	3.96	0.039	0.427	379.65	
Ni	2.23	-0.004	-0.048	-461.18	
IP	0.98	-0.021	-0.234	-105.07	
K	2.45	-0.088	-0.965	-413.37	
CSm	3.41	-0.098	-1.079	-254.04	
Se	1.14	-0.053	-0.581	-226.03	
Si	10.71	7.111	78.217	2.00	
Ag	8.81	(-0.045	(-0.499	-45.28	
Na	3.47	-0.010	-0.109	-1102.2	
Sr	2.67	-0.001	-0.007	-521.97	
S	0.59	0.018	0.195	130.35	
Ta	2.31	-0.022	-0.241	-228.87	
Tl	2.65	-0.347	-3.822	-8.64	
Th	1.03	-0.049	-0.541	-338.82	
Sn	0.84	-0.019	-0.211	-55.71	
Ti	2.38	-0.005	-0.060	-213.23	
W	0.88	-0.016	-0.174	-200.49	
U	3.27	-0.608	-6.689	-242.14	
V	2.92	-0.021	-0.230	-15.23	
Zn	1.52	(-0.008	(-0.086	-37.76	
Zr	3.42	-0.007	-0.082	-405.39	

Dilution factor : 11.0000

Sample name : 82C11A
Programme : SST 22-Feb-90 16:52:38

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	4.10	9.447	103.92	1.14	
Sb	0.32	-0.045	-0.490	-19.25	
As	0.90	0.226	2.490	0.31	
Ba	102.55	8.889	97.778	1.07	
Be	0.55	0.002	0.021	8.73	
Bi	2.52	0.044	0.484	38.50	
B	3.13	0.006	0.069	76.55	
Cd	1.48	-0.005	-0.055	-44.76	
Ca	3.01	0.052	0.575	3.62	
Ce	6.83	7.885	86.740	1.61	
Cr	17.83	8.586	94.445	1.24	
Co	0.91	7.150	78.653	1.25	
Cu	1.98	-0.005	-0.060	-55.03	
Eu	3.06	0.020	0.222	3.66	
Fe	29.45	8.800	96.802	1.31	
La	0.32	0.018	0.198	30.00	
Pb	0.24	-0.357	-3.932	-15.80	
Li	80.79	9.542	104.96	0.65	
Mg	1.44	0.067	0.732	1.44	
Mn	73.02	8.828	97.105	1.25	
Hg	2.79	0.039	0.432	27.90	
Mo	1.19	0.020	0.218	21.04	
Nd	11.64	8.027	88.301	1.70	
Ni	40.02	8.776	96.540	1.20	
P	0.98	-0.019	-0.214	-128.66	
K	2.43	-0.217	-2.391	-29.69	
Sm	3.39	-0.161	-1.772	-21.16	
Se	2.16	3.612	39.728	0.40	
Si	2.65	0.379	4.166	2.87	
Ag	9.43	-0.006	-0.062	-67.00	
Na	3.46	-0.019	-0.214	-55.91	
Sr	2.70	0.001	0.009	51.28	
S	0.66	0.141	1.553	14.98	
Ta	14.66	9.099	100.09	1.56	
Tl	2.74	0.043	0.473	356.70	
Th	1.07	0.182	2.018	6.34	
Sn	10.86	8.210	90.314	1.96	
Ti	60.93	9.530	104.93	1.10	
W	0.89	0.000	0.005	2551.40	
U	3.63	4.606	50.665	5.66	
V	2.92	-0.020	-0.217	-15.33	
Zn	2.03	0.022	0.246	5.04	
Zr	25.48	9.288	102.17	1.01	

Dilution factor : 11.0000

Sample name : 34C11C
Sample code 1 : DIGEST
Programme : SST

F5047
22-Feb-90 16:56:57

NAME	MV	INT	CONCEN	RSD
Al	1.35	0.022	221.82	
Sb	0.32	-0.069	-49.49	
As	0.74	0.025	31.62	
Ba	2.59	0.002	177.62	
Be	0.50	-0.001	-19.61	
Bi	4.06	1.845	8.54	
B	22.71	2.032	8.09	
Cd	15.98	1.937	8.90	
Ca	37.05	2.264	8.58	
Ce	3.44	-0.166	-78.49	
Cr	1.04	0.020	12.33	
Co	0.23	-0.012	-199.99	
Cu	9.92	2.093	7.63	
Eu	2.64	-0.004	-70.20	
Fe	1.47	0.103	9.46	
La	0.66	1.815	8.58	
Pb	0.29	1.677	12.93	
Li	2.73	-0.005	-92.54	
Mg	31.22	2.044	8.44	
Mn	0.59	0.004	22.53	
Hg	2.49	-0.004	-151.82	
Mo	7.66	2.129	8.80	
Nd	3.92	0.000		
Ni	2.21	-0.008	-117.26	
P	3.01	2.169	9.21	
K	2.77	1.803	8.80	
Sm	3.37	-0.203	-70.91	
Se	1.13	-0.114	-79.45	
Sc	4.15	1.632	7.23	
Ag	41.08	2.023	7.63	
Na	5.41	2.209	6.33	
Sr	45.31	2.049	8.08	
S	0.67	0.172	6.91	
Ca	2.28	-0.041	-75.01	
Tl	2.62	-0.517	-7.94	
Th	1.03	-0.110	-84.47	
Sn	0.88	0.010	31.22	
Ti	3.24	0.134	7.77	
W	0.89	0.004	836.28	
U	3.24	-1.125	-71.96	
V	2.89	-0.030	-15.61	
Zn	34.93	2.000	7.85	
Zr	3.38	-0.026	-69.89	

Sample name : F5048
 Sample code 1 : 100-10
 Programme : SST 22-Feb-90 17:01:02

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	3.53	7.471	754.61	1.39	
Sb	0.33	-0.000	-0.000	*****	
As	0.71	-0.007	-0.698	-36.74	
Ba	2.57	0.000	0.021	173.20	
Be	0.50	-0.000	-0.034	-8.25	
Bi	6.33	4.505	455.03	1.84	
B	2.94	-0.014	-1.401	-36.96	
Cd	1.44	(-0.010	(-1.051	-15.36	
Ca	1.29	(-0.059	(-6.002	-1.29	
Ce	3.47	-0.104	-10.48	-4.77	
Cr	1.19	0.096	9.691	7.38	
Co	0.22	-0.066	-6.699	-39.03	
Cu	1.98	-0.005	-0.525	-11.74	
Eu	2.66	-0.002	-0.239	-4.96	
Fe	6.62	1.705	172.23	1.82	
La	0.32	-0.009	-0.909	-60.00	
Pb	0.25	-0.032	-3.735	-100.00	
Li	2.74	-0.004	-0.375	-25.75	
Mg	0.77	0.022	2.246	1.41	
Mn	4.08	0.429	43.306	0.70	
He	2.49	-0.005	-0.471	-100.54	
Mo	1.12	-0.003	-0.285	-58.07	
Nd	3.91	-0.015	-1.470	-536.30	
Ni	2.24	-0.001	-0.117	-406.12	
Si	2.84	1.988	200.83	1.79	
K	2.45	-0.068	-6.847	-22.21	
GSm	3.40	-0.135	-13.63	-2.10	
Se	1.15	-0.047	-4.725	-42.36	
Si	2.35	0.128	12.927	1.73	
Ag	8.72	(-0.051	(-5.199	-0.76	
Na	15.78	14.090	1423.1	1.79	
Sr	3.77	0.052	5.271	1.48	
S	0.69	0.198	20.016	5.34	
Ta	2.29	-0.036	-3.628	-43.89	
Tl	2.68	-0.203	-20.50	-55.30	
Th	1.03	-0.112	-11.30	-6.00	
Sn	0.84	-0.022	-2.186	-53.16	
Ti	2.38	-0.006	-0.581	-14.98	
W	0.88	-0.016	-1.643	-39.63	
U	3.28	-0.555	-56.10	-10.45	
V	2.95	-0.013	-1.276	-27.98	
Zn	1.74	0.005	0.500	25.02	
Zr	3.41	-0.011	-1.107	-7.69	

Dilution factor : 101.000

dilution factor : 11.0000

A1	21.69	069.612	0765.73	1.10
S15	0.34	0.144	1.579	98.68
A5	0.79	0.093	1.028	51.14
B3	3.06	0.044	0.481	27.58
B4	6.92	0.307	3.373	2.19
C3	3.68	0.400	4.401	122.19
C4	3.14	1.090	11.987	5.10
C5	0.22	(-0.149	(-1.641	-8.23
C6	2.21	0.054	0.599	47.86
C7	2.80	0.005	0.058	168.43
D1	54.89	16.706	183.77	1.94
D2	0.32	0.020	0.218	193.49
D3	0.84	0.008	0.099	209.99
D4	4.25	0.259	2.848	42.96
D5	2.76	0.025	0.384	47.41
D6	1.25	0.039	0.433	46.78
D7	4.08	0.167	1.833	137.53
D8	2.74	0.116	1.274	29.75
D9	18.63	19.027	209.40	1.50
E1	2.53	0.391	4.299	153.01
E2	1.40	0.235	2.581	198.07
E3	3.55	0.235	1.375	80.52
E4	1.48	1.703	18.726	0.78
E5	1.06	0.150	1.649	211.16
E6	2.84	0.526	5.783	141.36
E7	1.49	0.061	0.669	44.70
E8	0.94	0.061	1.649	211.16
E9	1.03	0.157	2.169	30.06
F1	3.61	4.180	45.978	62.46
F2	3.15	0.049	0.587	72.40
F3	3.53	0.113	1.238	4.10
F4	3.65	0.087	0.961	64.85

NAME	MU INT	CONCEN	OILCOPR	RSD
------	--------	--------	---------	-----

Programme : SSI

Sample Code 1 : I-10

Sample Name : F5048

Sample name : F5049
 Sample code 1 : 100-10
 Programme : SST 22-Feb-90 17:11:01

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	3.94	8.879	896.77	2.61	
Sb	0.32	-0.079	-2.000	-19.75	
As	0.73	0.016	1.642	56.79	
Ba	2.60	0.003	0.254	45.52	
Be	0.50	-0.000	-0.034	-21.82	
Bi	6.38	4.566	461.21	2.38	
R	2.89	(-0.019	(-1.937	-7.48	
Cd	1.45	-0.009	-0.920	-6.12	
Ca	1.32	(-0.057	(-5.760	-2.63	
Ce	3.50	-0.029	-2.961	-135.43	
Cr	1.20	0.101	10.155	19.69	
Co	0.21	(-0.191	(-19.26	-9.96	
Cu	2.01	0.002	0.160	201.60	
Eu	2.70	-0.001	-0.055	-145.94	
Fe	7.28	1.910	192.94	3.71	
La	0.31	-0.027	-2.726	-23.09	
Pb	0.24	-0.308	-31.13	-6.93	
Li	2.75	-0.003	-0.296	-52.10	
Mg	0.84	0.027	2.742	3.73	
Mn	4.50	0.480	48.500	2.34	
Hg	2.59	0.011	1.077	10.66	
Mo	1.13	0.000	0.011	3036.50	
Nd	3.94	0.022	2.240	245.90	
Ni	2.31	0.014	1.424	17.16	
P	3.07	2.236	225.84	4.94	
K	2.44	-0.120	-12.08	-28.43	
Sm	3.43	-0.066	-6.689	-57.65	
Se	1.17	0.036	3.634	101.16	
Si	2.54	0.291	29.423	4.91	
Ag	3.84	(-0.043	(-4.388	-5.02	
Na	16.88	15.342	1549.5	2.76	
Sr	3.88	0.058	5.810	2.93	
S	0.67	0.166	16.755	8.02	
Ta	2.34	0.000	0.025	2787.50	
Tl	2.71	-0.066	-6.677	-59.34	
Th	1.03	-0.087	-8.811	-39.47	
Sn	0.86	-0.005	-0.526	-98.61	
Ti	2.41	-0.001	-0.126	-126.68	
W	0.90	0.014	1.444	68.89	
U	3.31	-0.010	-0.969	-2471.7	
V	2.97	-0.004	-0.392	-169.64	
Zn	1.75	0.006	0.587	11.06	
Zr	3.44	0.001	0.142	467.07	

Dilution factor : 101.000

Sample name : F5049
 Sample code 1 : 1-10
 Programme : SST 22-Feb-90 17:15:40

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	24.81	180.316	1883.47	2.01	
Sb	0.33	-0.005	-0.054	-793.69	
As	0.76	0.056	0.612	20.23	
Ba	2.97	0.035	0.389	6.46	
Be	0.52	0.000	0.004	67.51	
Bi	37.24	40.639	447.03	2.26	
B	2.90	(-0.018	(-0.198	-43.33	
Cd	1.48	-0.005	-0.060	-50.65	
Ca	7.23	0.327	3.595	3.17	
Ce	3.52	0.021	0.227	375.14	
Cr	3.10	1.070	11.775	4.22	
Co	0.21	(-0.207	(-2.280	-9.17	
Cu	2.12	0.032	0.349	19.55	
Eu	2.67	-0.002	-0.023	-59.81	
Fe	55.24	16.816	184.98	3.09	
La	0.31	-0.032	-0.356	-41.94	
Pb	0.25	-0.086	-0.949	-24.74	
Li	2.73	-0.005	-0.058	-44.03	
Mg	4.39	0.262	2.887	3.08	
Mn	35.06	4.203	46.229	3.01	
Hg	2.64	0.017	0.192	27.70	
Mo	1.21	0.024	0.263	19.11	
Nd	3.83	-0.095	-1.049	-18.55	
Ni	2.65	0.095	1.041	11.51	
P	18.55	18.946	208.40	3.69	
K	2.45	-0.114	-1.250	-56.27	
Sm	3.40	-0.134	-1.475	-43.71	
Se	1.36	0.726	7.983	8.45	
Si	4.67	2.070	22.765	3.19	
Ag	8.80	(-0.046	(-0.505	-13.45	
Na	125.36	139.54	1534.9	1.95	
Sr	13.41	0.516	5.674	2.46	
Ta	1.49	1.729	19.031	3.95	
Tl	2.38	0.028	0.309	94.01	
Th	2.74	0.034	0.372	493.63	
Sn	1.02	-0.134	-1.476	-30.55	
Ti	0.91	0.036	0.401	40.37	
W	2.45	0.005	0.054	73.51	
U	0.98	0.125	1.372	25.61	
V	3.47	2.207	24.279	19.75	
Zn	3.01	0.005	0.060	203.80	
Zr	3.46	0.108	1.193	4.20	
	3.49	0.015	0.162	66.58	

Dilution factor : 11.0000

Sample name : HNO3
Programme : SST 22-Feb-90 17:19:54

NAME	MV	INT	CONCEN	RSD
Al	1.29	-0.194	-2.70	
Sb	0.31	-0.158	-5.41	
As	0.69	-0.036	-7.78	
Ba	2.50	-0.006	-8.41	
Be	0.50	-0.001	-19.24	
Bi	2.34	(-0.162	-7.31	
B	2.22	(-0.036	-2.74	
Cd	1.41	(-0.015	-6.97	
Ca	0.47	(-0.113	-0.12	
Ce	3.42	-0.211	-4.56	
Cr	0.92	(-0.045	-13.35	
Co	0.22	-0.112	-11.11	
Cu	1.94	-0.018	-12.21	
Eu	2.63	-0.005	-7.74	
Fe	1.09	-0.016	-36.28	
La	0.31	-0.038	-16.50	
Pb	0.24	(-0.493	-8.66	
Li	2.70	-0.008	-21.25	
Mg	0.34	(-0.006	-0.60	
Mn	0.55	-0.002	-18.41	
Hg	2.43	-0.013	-13.07	
Mo	1.09	-0.013	-15.86	
Nd	3.80	-0.126	-16.94	
Ni	2.17	-0.017	-30.38	
P	0.95	-0.049	-23.17	
K	2.42	-0.277	-7.58	
Sm	3.36	-0.230	-6.67	
Se	1.10	(-0.209	-5.54	
Si	2.09	-0.092	-5.55	
As	0.59	(-0.060	-3.07	
Na	3.36	-0.136	-5.26	
Sr	2.63	-0.003	-10.91	
S	0.56	-0.038	-49.56	
Ta	2.25	-0.068	-15.95	
Tl	2.62	-0.489	-11.59	
Th	1.02	-0.172	0.00	
Sr	0.82	(-0.039	-5.24	
Ti	2.33	-0.013	-6.04	
W	0.86	-0.048	-22.98	
U	3.22	-1.422	-2.07	
V	2.89	-0.081	-11.45	
Zn	1.49	(-0.010	-3.19	
Zr	3.37	-0.028	-9.08	

NAME	MV INT	CONCEN	DILCEN	RSD	Dilution Factor : 101.000
A1	3.71	8.100	819.09	1.65	
Sb	0.33	0.010	1.000	150.01	
As	0.69	-0.038	-3.817	-16.13	
Ba	3.62	0.094	9.453	1.38	
Be	6.10	4.236	427.80	0.99	
Ca	0.50	-0.001	-0.027	-16.50	
Ce	3.39	-0.284	-28.65	-6.54	
Cr	1.31	0.156	15.757	12.73	
Co	0.23	0.017	1.675	96.60	
Cu	2.31	0.091	8.331	3.50	
Eu	2.60	0.006	-0.583	-6.86	
Fb	6.84	1.772	178.96	2.70	
Fe	0.33	0.058	5.815	10.83	
Li	3.60	0.102	10.285	1.56	
Mg	2.28	0.123	12.381	2.20	
Mn	4.82	0.519	52.370	2.13	
Mo	1.39	0.015	-1.548	-12.90	
Nd	2.42	0.103	-10.43	-18.90	
P	3.09	2.254	227.65	4.32	
Pa	2.58	0.077	7.783	5.73	
Na	3.82	-0.103	-10.43	-18.90	
Ng	4.42	0.281	-28.40	-13.00	
Os	1.13	0.104	-10.54	-24.95	
Pa	3.33	0.316	-31.96	-10.29	
Pr	1.13	0.104	-10.54	-24.95	
Sm	3.33	0.316	-31.96	-10.29	
Ta	2.22	0.086	-8.673	-7.49	
Tb	1.01	-0.208	-21.01	-6.72	
Tm	1.01	-0.208	-21.01	-6.72	
U	3.20	-1.652	-166.8	-7.13	
V	0.86	-0.047	-4.790	-17.40	
W	2.90	0.079	7.988	1.60	
Zn	3.30	0.099	10.007	2.10	
Zr	3.96	-0.095	-3.593	-2.41	

Sample name : F5050
 Sample Code 1 : 100-10
 Programme : SST
 22-Feb-90 17:23:49

Sample name : F5050
 Sample code 1 : 1-10
 Programme : SST . 22-Feb-90 17:27:43

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	24.73	180.028	1880.31	1.87	
Sb	0.39	1.005	11.055	5.33	
As	0.80	0.098	1.073	3.31	
Ba	13.70	0.990	10.885	1.90	
Be	0.52	0.000	0.004	60.84	
Bi	38.99	42.673	469.40	1.84	
B	11.80	0.903	9.934	1.59	
Cd	8.32	0.911	10.020	2.24	
Ca	22.78	1.337	14.704	2.05	
Ce	3.52	0.015	0.166	114.22	
Cr	4.93	2.002	22.025	3.41	
Co	0.28	0.663	7.295	4.72	
Cu	5.79	1.000	10.998	1.95	
Eu	2.71	0.000	0.001	196.39	
Fe	58.34	17.777	195.55	1.95	
La	0.48	0.865	9.520	3.30	
Pb	0.29	1.072	11.797	6.90	
Li	11.45	1.061	11.675	2.00	
Mg	19.28	1.251	13.758	1.84	
Mn	42.60	5.121	56.336	1.96	
Hg	2.63	0.016	0.171	7.72	
Mo	4.26	1.021	11.233	1.57	
Nd	4.69	0.804	8.643	10.84	
Ni	6.61	1.014	11.151	1.55	
PP	22.20	22.888	251.76	3.36	
K	2.60	0.814	0.949	0.42	
Sm	3.39	-0.161	-1.772	-5.28	
Se	1.47	1.119	12.310	3.73	
Si	4.21	1.681	18.495	1.16	
Ag	24.11	0.935	10.283	1.99	
Na	133.37	149.28	1642.1	2.00	
Sr	33.61	1.487	16.358	1.81	
S	1.60	1.933	21.263	1.92	
Ta	2.38	0.031	0.346	16.24	
Tl	2.72	-0.052	-0.575	-89.26	
Th	1.03	-0.078	-0.861	-22.68	
Sn	1.96	0.902	9.926	2.51	
Ti	0.19	0.939	10.332	1.91	
W	0.99	0.142	1.562	10.05	
U	3.47	2.255	24.806	3.31	
V	3.01	0.005	0.057	5.77	
Zn	19.39	1.065	11.720	1.70	
Zr	3.55	0.045	0.498	3.77	

Dilution factor : 11.0000

Sample name : HN03
Programme : SST 22-Feb-90 17:31:35

NAME	MV	INT	CONCEN	RSD
Al	1.36	0.062	270.52	
Sb	0.32	-0.025	-150.99	
As	0.73	0.014	261.44	
Ba	2.63	0.005	169.88	
Be	0.51	0.000	414.31	
Bi	2.48	0.008	1629.55	
P	2.92	-0.016	-66.33	
Cd	1.49	-0.004	-169.22	
Ca	0.49	(-0.111	-0.93	
Ce	3.60	0.215	150.33	
Cr	0.95	(-0.026	-39.96	
Co	0.22	-0.120	-11.95	
Cu	2.05	0.011	178.60	
Eu	2.73	0.004	154.68	
Fe	1.13	-0.001	-469.98	
La	0.32	0.004	606.19	
Pb	0.25	-0.111	-66.67	
Li	2.82	0.006	191.52	
Mg	0.35	(-0.005	-8.82	
Mn	0.57	0.001	182.48	
Hg	2.58	0.009	143.04	
Mo	1.15	0.005	297.34	
Nd	3.96	0.039	459.61	
Ni	2.30	0.011	158.53	
P	1.00	0.000	10446.8	
K	2.51	0.263	169.66	
Sm	3.53	0.196	168.62	
Se	1.16	0.012	1213.98	
Si	2.22	0.018	357.81	
Ag	9.11	-0.026	-97.83	
Na	3.54	0.068	215.22	
Sr	2.73	0.002	177.97	
S	0.59	0.014	98.33	
Ta	2.39	0.038	173.59	
Tl	2.76	0.125	510.00	
Th	1.06	0.105	227.50	
Sn	0.87	0.005	495.63	
Ti	2.46	0.008	182.69	
W	0.91	0.028	188.37	
U	3.39	1.135	164.23	
V	3.09	0.030	130.45	
Zn	1.55	(-0.006	-32.57	
Zr	3.50	0.026	167.75	

Dilution factor : 11.0000

NAME	MV INT	CONCEN	DILCOR	RSI0
A1	1.37	0.079	0.966	8.70
A2	0.74	0.030	0.980	34.69
A3	0.52	0.032	0.989	34.69
A4	0.52	0.030	0.980	34.69
A5	0.52	0.031	0.980	34.69
A6	0.51	0.000	0.001	34.64
B1	0.51	0.003	0.030	37.07
B2	0.42	0.062	0.681	130.90
B3	0.85	0.023	0.252	14.48
C1	1.45	0.005	0.097	5.35
C2	0.89	0.085	0.940	1.69
C3	3.47	0.101	1.116	32.84
C4	0.96	0.023	0.351	79.16
C5	0.22	0.178	1.961	4.03
C6	3.01	0.001	0.010	217.01
C7	1.11	0.007	0.076	37.55
C8	2.67	0.002	0.021	34.09
C9	0.58	0.002	0.019	10.66
C10	0.55	0.008	0.086	7.87
C11	2.73	0.005	0.052	24.64
C12	0.24	0.284	3.119	7.53
C13	0.31	0.036	0.396	39.69
C14	1.11	0.007	0.076	37.55
C15	3.01	0.001	0.010	217.01
C16	0.96	0.023	0.351	79.16
C17	0.22	0.178	1.961	4.03
C18	1.11	0.001	0.010	217.01
C19	2.67	0.002	0.021	34.09
C20	0.58	0.002	0.019	10.66
C21	0.55	0.008	0.086	7.87
C22	2.73	0.005	0.052	24.64
C23	0.24	0.284	3.119	7.53
C24	0.31	0.036	0.396	39.69
C25	1.11	0.007	0.076	37.55
C26	3.01	0.001	0.010	217.01
C27	0.96	0.023	0.351	79.16
C28	0.22	0.178	1.961	4.03
C29	1.11	0.001	0.010	217.01
C30	2.67	0.002	0.021	34.09
C31	0.58	0.002	0.019	10.66
C32	0.55	0.008	0.086	7.87
C33	2.73	0.005	0.052	24.64
C34	0.24	0.284	3.119	7.53
C35	0.31	0.036	0.396	39.69
C36	1.11	0.007	0.076	37.55
C37	3.01	0.001	0.010	217.01
C38	0.96	0.023	0.351	79.16
C39	0.22	0.178	1.961	4.03
C40	1.11	0.001	0.010	217.01
C41	2.67	0.002	0.021	34.09
C42	0.58	0.002	0.019	10.66
C43	0.55	0.008	0.086	7.87
C44	2.73	0.005	0.052	24.64
C45	0.24	0.284	3.119	7.53
C46	0.31	0.036	0.396	39.69
C47	1.11	0.007	0.076	37.55
C48	3.01	0.001	0.010	217.01
C49	0.96	0.023	0.351	79.16
C50	0.22	0.178	1.961	4.03
C51	1.11	0.001	0.010	217.01
C52	2.67	0.002	0.021	34.09
C53	0.58	0.002	0.019	10.66
C54	0.55	0.008	0.086	7.87
C55	2.73	0.005	0.052	24.64
C56	0.24	0.284	3.119	7.53
C57	0.31	0.036	0.396	39.69
C58	1.11	0.007	0.076	37.55
C59	3.01	0.001	0.010	217.01
C60	0.96	0.023	0.351	79.16
C61	0.22	0.178	1.961	4.03
C62	1.11	0.001	0.010	217.01
C63	2.67	0.002	0.021	34.09
C64	0.58	0.002	0.019	10.66
C65	0.55	0.008	0.086	7.87
C66	2.73	0.005	0.052	24.64
C67	0.24	0.284	3.119	7.53
C68	0.31	0.036	0.396	39.69
C69	1.11	0.007	0.076	37.55
C70	3.01	0.001	0.010	217.01
C71	0.96	0.023	0.351	79.16
C72	0.22	0.178	1.961	4.03
C73	1.11	0.001	0.010	217.01
C74	2.67	0.002	0.021	34.09
C75	0.58	0.002	0.019	10.66
C76	0.55	0.008	0.086	7.87
C77	2.73	0.005	0.052	24.64
C78	0.24	0.284	3.119	7.53
C79	0.31	0.036	0.396	39.69
C80	1.11	0.007	0.076	37.55
C81	3.01	0.001	0.010	217.01
C82	0.96	0.023	0.351	79.16
C83	0.22	0.178	1.961	4.03
C84	1.11	0.001	0.010	217.01
C85	2.67	0.002	0.021	34.09
C86	0.58	0.002	0.019	10.66
C87	0.55	0.008	0.086	7.87
C88	2.73	0.005	0.052	24.64
C89	0.24	0.284	3.119	7.53
C90	0.31	0.036	0.396	39.69
C91	1.11	0.007	0.076	37.55
C92	3.01	0.001	0.010	217.01
C93	0.96	0.023	0.351	79.16
C94	0.22	0.178	1.961	4.03
C95	1.11	0.001	0.010	217.01
C96	2.67	0.002	0.021	34.09
C97	0.58	0.002	0.019	10.66
C98	0.55	0.008	0.086	7.87
C99	2.73	0.005	0.052	24.64
C100	0.24	0.284	3.119	7.53
C101	0.31	0.036	0.396	39.69
C102	1.11	0.007	0.076	37.55
C103	3.01	0.001	0.010	217.01
C104	0.96	0.023	0.351	79.16
C105	0.22	0.178	1.961	4.03
C106	1.11	0.001	0.010	217.01
C107	2.67	0.002	0.021	34.09
C108	0.58	0.002	0.019	10.66
C109	0.55	0.008	0.086	7.87
C110	2.73	0.005	0.052	24.64
C111	0.24	0.284	3.119	7.53
C112	0.31	0.036	0.396	39.69
C113	1.11	0.007	0.076	37.55
C114	3.01	0.001	0.010	217.01
C115	0.96	0.023	0.351	79.16
C116	0.22	0.178	1.961	4.03
C117	1.11	0.001	0.010	217.01
C118	2.67	0.002	0.021	34.09
C119	0.58	0.002	0.019	10.66
C120	0.55	0.008	0.086	7.87
C121	2.73	0.005	0.052	24.64
C122	0.24	0.284	3.119	7.53
C123	0.31	0.036	0.396	39.69
C124	1.11	0.007	0.076	37.55
C125	3.01	0.001	0.010	217.01
C126	0.96	0.023	0.351	79.16
C127	0.22	0.178	1.961	4.03
C128	1.11	0.001	0.010	217.01
C129	2.67	0.002	0.021	34.09
C130	0.58	0.002	0.019	10.66
C131	0.55	0.008	0.086	7.87
C132	2.73	0.005	0.052	24.64
C133	0.24	0.284	3.119	7.53
C134	0.31	0.036	0.396	39.69
C135	1.11	0.007	0.076	37.55
C136	3.01	0.001	0.010	217.01
C137	0.96	0.023	0.351	79.16
C138	0.22	0.178	1.961	4.03
C139	1.11	0.001	0.010	217.01
C140	2.67	0.002	0.021	34.09
C141	0.58	0.002	0.019	10.66
C142	0.55	0.008	0.086	7.87
C143	2.73	0.005	0.052	24.64
C144	0.24	0.284	3.119	7.53
C145	0.31	0.036	0.396	39.69
C146	1.11	0.007	0.076	37.55
C147	3.01	0.001	0.010	217.01
C148	0.96	0.023	0.351	79.16
C149	0.22	0.178	1.961	4.03
C150	1.11	0.001	0.010	217.01
C151	2.67	0.002	0.021	34.09
C152	0.58	0.002	0.019	10.66
C153	0.55	0.008	0.086	7.87
C154	2.73	0.005	0.052	24.64
C155	0.24	0.284	3.119	7.53
C156	0.31	0.036	0.396	39.69
C157	1.11	0.007	0.076	37.55
C158	3.01	0.001	0.010	217.01
C159	0.96	0.023	0.351	79.16
C160	0.22	0.178	1.961	4.03
C161	1.11	0.001	0.010	217.01
C162	2.67	0.002	0.021	34.09
C163	0.58	0.002	0.019	10.66
C164	0.55	0.008	0.086	7.87
C165	2.73	0.005	0.052	24.64
C166	0.24	0.284	3.119	7.53
C167	0.31	0.036	0.396	39.69
C168	1.11	0.007	0.076	37.55
C169	3.01	0.001	0.010	217.01
C170	0.96	0.023	0.351	79.16
C171	0.22	0.178	1.961	4.03
C172	1.11	0.001	0.010	217.01
C173	2.67	0.002	0.021	34.09
C174	0.58	0.002	0.019	10.66
C175	0.55	0.008	0.086	7.87
C176	2.73	0.005	0.052	24.64
C177	0.24	0.284	3.119	7.53
C178	0.31	0.036	0.396	39.69
C179	1.11	0.007	0.076	37.55
C180	3.01	0.001	0.010	217.01
C181	0.96	0.023	0.351	79.16
C182	0.22	0.178	1.961	4.03
C183	1.11	0.001	0.010	217.01
C184	2.67	0.002	0.021	34.09
C185	0.58	0.002	0.019	10.66
C186	0.55	0.008	0.086	7.87
C187	2.73	0.005	0.052	24.64
C188	0.24	0.284	3.119	7.53
C189	0.31	0.036	0.396	39.69
C190	1.11	0.007	0.076	37.55
C191	3.01	0.001	0.010	217.01
C192	0.96	0.023	0.351	79.16
C193	0.22	0.178	1.961	4.03
C194	1.11	0.001	0.010	217.01
C195	2.67	0.002	0.021	34.09
C196	0.58	0.002	0.019	10.66
C197	0.55	0.008	0.086	7.87
C198	2.73	0.005	0.052	24.64
C199	0.24	0.284	3.119	7.53
C200	0.31	0.036	0.396	39.69
C201	1.11	0.007	0.076	37.55
C202	3.01	0.001	0.010	217.01
C203	0.96	0.023	0.351	79.16
C204	0.22	0.178	1.961	4.03
C205	1.11	0.001	0.010	217.01
C206	2.67	0.002	0.021	34.09
C207	0.58	0.002	0.019	10.66
C208	0.55	0.008	0.086	7.87
C209	2.73	0.005	0.052	24.64
C210	0.24	0.284	3.119	7.53
C211	0.31	0.036	0.396	39.69
C212	1.11	0.007	0.076	37.55
C213</td				

Sample name : SICLIA
Programme : SST 22-Feb-90 17:40:41

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.29	-0.201	-2.208	-8.59	
Sb	0.32	-0.144	-1.579	-15.80	
As	0.87	0.193	2.119	2.76	
Ba	2.55	-0.002	-0.023	-33.47	
Be	0.50	-0.001	-0.008	-16.78	
Bi	9.95	8.733	96.060	0.92	
B	94.41	9.453	103.98	0.52	
Cd	68.86	9.021	99.234	1.43	
Ca	145.12	9.286	102.15	1.27	
Ce	3.45	-0.145	-1.595	-17.68	
Cr	0.94	(-0.032	(-0.032	-23.36	
Co	0.21	(-0.207	(-0.200	-9.17	
Cu	37.06	9.265	101.92	1.09	
Eu	2.65	-0.003	-0.036	-14.04	
Fe	1.09	-0.015	-0.163	-31.68	
La	1.85	8.251	90.763	1.13	
Pb	0.47	8.185	90.039	1.38	
Li	2.68	-0.011	-0.126	-12.33	
Mg	137.25	9.081	99.896	1.36	
Mn	0.58	0.002	0.026	13.02	
Hg	2.55	0.004	0.040	182.76	
Mo	1.11	-0.008	-0.090	-11.55	
Nd	3.99	0.074	0.812	35.22	
Ni	2.20	-0.010	-0.113	-31.09	
P	9.56	9.250	101.74	5.69	
K	3.91	8.620	94.818	1.18	
Sm	3.33	-0.301	-3.310	-7.78	
Se	1.12	-0.151	-1.662	-19.10	
Si	2.18	-0.014	-0.156	-88.30	
Ag	156.62	9.426	103.69	0.91	
Na	11.19	8.835	97.182	1.67	
Sr	194.23	9.210	101.31	1.22	
S	0.62	0.084	0.919	22.77	
Ta	2.28	-0.039	-0.430	-20.78	
Tl	2.63	-0.435	-4.786	-18.95	
Th	1.00	-0.257	-2.830	-9.16	
Sn	0.96	0.076	0.841	4.35	
Ti	2.34	-0.013	-0.143	-10.97	
W	0.95	0.092	1.009	7.39	
U	3.24	-1.029	-11.32	-15.68	
V	2.90	-0.027	-0.292	-6.77	
Zn	152.26	9.049	99.544	1.22	
Cr	3.36	-0.032	-0.352	-12.69	

Dilution factor : 11.0000

Sample name : 82C11A
Sample code 1 : 1-10
Programme : SST 22-Feb-90 17:44:52

NAME	MV	INT	CONCEN	DILCOR	RSB
Al	4.25	9.934	109.27	1.79	
Sb	0.32	-0.054	-0.599	-155.08	
As	0.92	0.252	2.767	18.90	
Ba	107.07	9.291	102.20	1.67	
Be	0.55	0.002	0.022	30.25	
Bi	2.55	0.080	0.083	194.95	
R	3.41	0.035	0.089	55.20	
Cd	1.50	-0.002	-0.027	-335.48	
Ca	3.23	0.067	0.733	3.77	
Ce	7.01	8.326	91.585	4.95	
Cr	18.88	9.124	100.36	3.03	
Co	0.89	8.158	89.733	6.92	
Cu	2.02	0.004	0.043	573.25	
Eu	3.12	0.023	0.256	30.59	
Fe	31.42	9.412	103.53	3.19	
La	0.33	0.029	0.317	117.09	
Pb	0.24	-0.303	-3.390	-36.66	
Li	82.94	9.804	107.85	3.75	
Mg	1.56	0.075	0.822	3.20	
Mn	77.84	9.415	103.56	3.62	
Hg	2.70	0.027	0.296	48.96	
Mo	1.21	0.025	0.275	47.94	
Nd	11.95	8.342	91.766	1.37	
Ni	42.64	9.384	103.22	3.65	
P	0.99	-0.001	-0.008	-6300.0	
K	2.44	-0.156	-1.711	-384.33	
Sm	3.43	-0.066	-0.728	-605.57	
Se	2.25	3.943	43.370	5.52	
Si	2.70	0.425	4.671	20.04	
Ag	9.61	0.006	0.067	500.33	
Na	3.50	0.023	0.248	769.75	
Sr	2.73	0.002	0.026	198.30	
S	0.66	0.153	1.679	28.14	
Ta	15.41	9.651	106.16	4.15	
Tl	2.87	0.670	7.373	102.78	
Th	1.08	0.248	2.721	108.66	
Sn	11.62	8.835	97.195	5.21	
Ti	64.19	10.053	110.58	1.46	
W	0.90	0.016	0.179	337.94	
U	3.70	5.554	61.093	41.07	
V	3.02	0.009	0.098	406.78	
Zn	2.11	0.027	0.299	15.97	
Zr	26.72	9.810	107.91	1.43	

dilution factor : 11.0000

Sample name : 81C11A
Programme : SST 22-Feb-90 18:59:39

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.29	-0.175	-1.920	-15.31	
Sb	0.31	-0.233	-2.500	-3.69	
As	0.70	-0.017	-0.183	-59.59	
Ba	2.42	(-0.013	(-0.144	-6.78	
Be	0.49	-0.001	-0.012	-10.94	
Bi	2.26	(-0.259	(-2.854	-16.35	
B	2.71	(-0.037	(-0.408	-14.59	
Cd	1.37	(-0.020	(-0.224	-7.39	
Ca	1.10	(-0.072	(-0.787	-0.82	
Ce	3.30	(-0.494	(-5.429	-7.83	
Cr	0.93	(-0.041	(-0.449	-56.25	
Co	0.22	(-0.178	(-1.901	-4.03	
Cu	1.92	-0.023	-0.256	-23.30	
Eu	2.54	(-0.009	(-0.102	-8.87	
Fe	1.06	(-0.025	(-0.276	-15.03	
La	0.31	-0.074	-0.811	0.00	
Pb	0.23	(-0.579	(-6.373	-9.75	
Li	2.62	(-0.018	(-0.203	-3.88	
Zg	0.67	0.016	0.171	2.02	
Mn	0.58	0.002	0.021	101.79	
Hg	2.46	-0.010	-0.109	-127.32	
No	28.94	9.065	99.719	0.93	
Od	3.69	(-0.236	(-2.592	-37.56	
Ni	2.10	(-0.033	(-0.368	-9.76	
P	0.90	(-0.099	(-1.088	-12.58	
K	2.35	-0.708	-7.786	-6.23	
Sm	3.25	(-0.504	(-5.549	-8.31	
Se	1.10	(-0.211	(-2.322	-18.98	
Si	10.27	6.740	74.144	0.48	
Hg	8.32	(-0.077	(-0.844	-5.93	
Na	3.29	(-0.210	(-2.313	-11.64	
Sr	2.57	-0.005	-0.060	-9.27	
C	0.55	-0.049	-0.543	-32.03	
Ta	2.20	(-0.103	(-1.134	-16.78	
Tl	2.57	(-0.723	(-7.948	-15.60	
Th	0.99	-0.336	-3.691	-7.57	
Sr.	0.80	(-0.054	(-0.591	-16.20	
Ti	2.27	(-0.024	(-0.261	-9.60	
W	0.84	(-0.077	(-0.846	-28.46	
U	3.10	(-3.021	(-33.23	-8.10	
V	2.83	(-0.047	(-0.522	-7.97	
Zn	1.43	(-0.013	(-0.147	-6.54	
Zr	3.31	-0.057	-0.624	-12.30	

Dilution factor : 11.0000

Sample name : 81C11A
Programme : SST 22-Feb-90 19:03:42

NAME	MV	INT	CONCEN	DILCOR	RSD
Al		1.28	-0.226	-2.484	-27.81
Sb		0.32	-0.099	-1.089	-31.22
As		0.88	0.194	2.132	10.39
Ba		2.54	-0.003	-0.029	-133.67
Be		0.49	-0.001	-0.011	-7.27
Bi		10.11	8.921	92.134	2.96
P		95.29	9.544	104.98	1.84
Cd		68.28	8.948	98.369	3.75
Ca		146.80	9.395	103.35	2.51
Ce		3.44	-0.166	-1.821	-79.65
Cr		1.01	0.004	0.041	618.98
Co		0.22	-0.083	-0.912	-31.22
Cu		37.73	9.442	103.87	2.27
Eu		2.64	-0.004	-0.042	-57.60
Fe		1.09	-0.016	-0.177	-35.21
La		1.87	8.343	91.772	2.79
Pb		0.47	8.173	89.904	3.70
Li		2.67	-0.012	-0.134	-32.67
Mg		139.01	9.132	100.45	2.78
Mn		0.59	0.004	0.044	3.03
Hg		2.46	-0.009	-0.104	-81.70
Mo		1.10	-0.012	-0.133	-45.25
Nd		4.03	0.116	1.277	117.43
Ni		2.17	-0.016	-0.199	-52.27
P		0.98	9.702	106.72	3.76
K		3.96	8.973	98.700	2.94
Sm		3.33	-0.312	-3.436	-40.43
Se		1.11	-0.193	-2.124	-6.45
Si		2.16	-0.027	-0.300	-122.56
Ag		159.92	9.638	106.01	2.53
Na		11.44	9.121	100.33	2.36
Sr		197.49	9.366	103.03	2.40
S		0.63	0.101	1.107	10.39
Ta		2.24	-0.068	-0.752	-28.68
Tl		2.58	(-0.664	(-7.305	-8.01
Tb		1.01	-0.226	-2.485	-38.84
Sn		0.94	0.066	0.726	22.28
Ti		2.31	-0.016	-0.181	-34.14
W		0.95	0.085	0.933	9.61
U		3.23	-1.207	-13.27	-61.58
V		2.84	(-0.045	(-0.491	-4.17
Zn		153.78	9.141	100.55	2.89
Zr		3.35	-0.039	-0.427	-40.22

Dilution factor : 11.0000

Sample name : 82C11A
Programme : SST 22-Feb-90 19:07:53

NAME	MV	INT	CONCEN	DILCOR	RSB
Al		4.09	9.395	103.24	4.09
Sb		0.32	-0.010	-0.109	-624.49
As		0.88	0.200	2.195	11.44
Ba	102.65	6.898		97.676	3.35
Be		0.54	0.002	0.012	28.64
Bi		2.44	-0.039	-0.429	-174.47
R		3.22	0.015	0.169	72.22
Cd		1.44	(-0.010)	(-0.112)	-35.34
Ca		3.00	0.052	0.567	11.97
Co		6.78	7.776	85.537	5.64
Cr	17.39	8.363		91.992	4.30
Co		0.83	7.523	82.756	5.37
Cu		1.95	-0.016	-0.123	-36.46
Eu		3.00	0.016	0.180	13.47
Fe	29.07	8.693		93.513	4.26
La		0.33	0.031	0.336	56.73
Pb		0.24	-0.160	-1.763	-26.65
Li	81.37	9.613		105.74	2.69
Hg		1.42	0.065	0.718	5.16
Mn	71.03	8.683		95.500	4.14
Hg		2.56	0.006	0.069	82.79
Mo		1.17	0.012	0.137	50.88
Nd	11.21	7.678		84.463	2.59
Ni	39.11	8.564		94.200	4.05
P		0.97	-0.029	-0.317	-77.94
K		2.28	-0.495	-5.440	-23.61
Sm		3.31	-0.343	-3.777	-25.58
Se		2.10	3.379	37.169	5.75
Si		2.59	0.331	3.636	12.60
		9.20	-0.021	-0.226	-45.08

20.45				
V		2.88	-0.034	-0.370
Zn		2.01	0.021	0.232
Zr		25.25	9.194	101.14

Dilution factor : 11.0000

APPENDIX A

ANALYTICAL ANALYSIS CARDS

卷之三

9 1 1 2 2 5 2 0 7 ? 4

Physical Properties

Serial No.	Sample Point	Date	Time Issued	Priority
F 5033.-5201	SEGMENT-4	11- 7-89	11:16	18
Determination	Method/Standard	Result Units	Charge Code	Reruns
VDA-SAMP	LI-000-200	NONE	WB75L	0
Sample Size		Customer ID		
?		84-041		
Remarks, Calculations, Results: THIRD SAMPLE Bottle ID # 030 GR=53 21.82 TARE 12.13 SAMPLE WT 4.39g.				
Analyst-1 SAS-11-8-84	Analyst-2 R2H 12-14-85	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date	Time Completed	Lab Unit Mgr	KP	

34-6800-061 (R-10-83)

Serial No.	Sample Point	Date	Time Issued	Priority
F 5033.-5000	SEGMENT-2	11- 7-89	11:15	18
Determination	Method/Standard	Result Units	Charge Code	Reruns
APPAR/DTIR	LI-000-200	NONE	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: a. JAR ID# 004 tip/100g Solids Pantrometer 60 b. JAR TARE WT. c. JAR TOTAL WT. d. C-BE e. EST. VOL./LENGTH 16" x 8" x 16" f. VISUAL REMARKS Bottom 4" very soft Brown - Runny Free running middle Firmer Brown clear liquid Top 5" firm Black dark exterior				
Analyst-1 SAS-11-8-84 16093	Analyst-2 R2H 12-14-85	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date	Time Completed	Lab Unit Mgr	KP	

34-6800-061 (R-10-83)

Serial No.	Sample Point	Date	Time Issued	Priority
F 5033.-5003	SEGMENT-2	11- 7-89	11:15	16
Determination	Method/Standard	Result Units	Charge Code	Reruns
HOMOGZT	LI-000-200	NONE	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: Homogenization Complete Completed				
Analyst-1 DJH 12-14-85	Analyst-2 R2H 12-14-85	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date	Time Completed	Lab Unit Mgr	KP	

34-6800-061 (R-10-83)


Serial No.	Sample Point	Date	Time Issued	Priority
F 5033.-5002	SEGMENT-2	11- 7-89	11:15	24
Determination	Method/Standard	Result Units	Charge Code	Reruns
PRT-SIZE	LI-000-200	NONE	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: PARTICLE SIZE DISTRIBUTION Bottle ID# 033 results: 1) particles dia. <150 μm 2) 5 sheets attached GR=53 18.63 TARE 12.85 SAMPLE WT 0.78g				
Analyst-1 SAS-11-8-84	Analyst-2 SBR 11-15-89	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date	Time Completed	Lab Unit Mgr	KP	

34-6800-061 (R-10-83)

Physical Properties

Liquids	Solids
SAR 60%	SAR 0%
239.50	527.40
222.73	222.87
16.77%	305.03%

F 5033,-5000



9 1 1 2 3 5 9 0 7 3 6

pH Analysis of Solid Samples

Serial No.	Sample Point		Date	Time issued	Priority
F 5035 -5015	SEGMENT-2		11- 7-89	11:15	15
Determination	Method/Standard	Result Units	Charge Code	Reruns	
pH	LA-212-103	NONE	WB75L	0	
Sample Size	Customer ID				
?	2.9389 / 2.935 ml		89-041		
Remarks, Calculations, Results:					
pH	12.71	Bottle #44	SQWAD		
SAMPLE TEMP	24.6				
	2.938				
		36.19			
		22.29			
		13.909			
	WHL-N-313-# Pg 3				
Analyst-1 JR Smith 65186	Analyst-2 CC269	Analyst-3	Analyst-4	Analyst-5 R.E. Bennett	
DL Herrington 66431	Hrs	Mins	Hrs	Hrs	
Date 11-22-89	Time Completed	Lab Unit Mgr	CJA	dl	
					54-6800-061 (R-10-83)

Serial No.	Sample Point		Date	Time issued	Priority
F 5035 -5015	SEGMENT-3		11- 7-89	11:15	15
Determination	Method/Standard	Result Units	Charge Code	Reruns	
pH	LA-212-103	% RECOVERY	WB75L	0	
Sample Size	Customer ID				
?	2.1029 / 2.102 ml		89-041		
Remarks, Calculations, Results:					
pH	12.92	Bottle #44	KP		
SAMPLE TEMP	24.6				
	2.1029				
		100.0%			
		10.00 / 10.00			
	Duplicate				
	2.1014				
Analyst-1 CC269	Analyst-2	Analyst-3	Analyst-4	Analyst-5 R.E. Bennett	
Hrs	Mins	Hrs	Hrs	Hrs	
Date 11-22-89	Time Completed	Lab Unit Mgr	CJA	dl	
					54-6800-061 (R-10-83)

9 1 1 2 0 5 0 0 7 ? 7

pH Analysis of Solid Sample

Serial No.	Sample Point	Date	Time Issued	Priority
F 5053 - 5315	SEGMENT-22	11-7-89	11:18	1B
Determination	Method/Standard	Result Units	Charge Code	Reruns
pH	LA-212-103	NONE	WB75L	0
Sample Size				Customer ID
?				
Remarks, Calculations, Results:				
LMCS CHECK SAMPLE "Q" pH FOUND <u>5.67</u> STD ID <u>H2O</u> SAMPLE TEMP <u>24.7</u> H ₂ O Blank <u>5.67 (OK)</u> CG				
Analyst - 1 <u>60269</u> Hrs <u>xx.7</u>	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5 <u>REB</u> Hrs
Date <u>11-27-89</u>	Time Completed	Lab Unit Mgr <u>CG</u>		Hrs <u>x</u>

SI-MOD-081 (R-10-83)

Percent Water Analysis

9 1 1 2 2 5 9 0 7 7 8

Serial No.	Sample Point		Date	Time Issued	Priority
F 5033.-5010	SEGMENT-2		11- 7-89	11:15	19
Determination	Method/Standard	Result Units	Charge Code	Recons	
% H ₂ O	LA-564-101	%	WB75L	0	
Sample Size	Customer ID				
?	89-041				
Remarks, Calculations, Results:					
350 mPa 24.5416 G 6.276 KT Bottle # 47 23.5805 T 38.89% KT 24.1682 W1 28.80% 24.1678 W2 22.45% WAC-N-313-4 Pg 3 6.30g					
Analyst-1 JR Smith 6526	Analyst-2 6B598	Analyst-3	Analyst-4	Analyst-5 RE Bennett	
Hrs R.D. Hartung 66931	Hrs R.D. Hartung	Hrs	Hrs	Hrs	
Date 11-22-89	Time Completed 11-28	Lab Unit Mgr Cga	OK		

\$4-6800-081 (R-10-83)

Serial No.	Sample Point		Date	Time Issued	Priority
F 5032.-3410	SEGMENT-1		11- 7-89	11:15	19
Determination	Method/Standard	Result Units	Charge Code	Recons	
% H ₂ O	LA-564-101	% RECOVERY	WB75L	0	
Sample Size	Customer ID				
?	1mL				
Remarks, Calculations, Results:					
LMCS CHECK SAMPLE LMCS ID 24.7941 G 99.10% 23.4130 T 23.9785 W1 23.9760 W2 59.05% / 59.610					
Analyst-1 6B598	Analyst-2	Analyst-3	Analyst-4	Analyst-5 RE Bennett	
Hrs R.D. Hartung	Hrs	Hrs	Hrs	Hrs	
Date 11-22-89	Time Completed 11-28	Lab Unit Mgr Cga	OK		

\$4-6800-081 (R-10-83)

Serial No.	Sample Point		Date	Time Issued	Priority
F 5036.-5510	SEGMENT-5		11- 7-89	11:15	19
Determination	Method/Standard	Result Units	Charge Code	Recons	
% H ₂ O	LA-564-101	% RECOVERY	WB75L	0	
Sample Size	Customer ID				
?	1mL				
Remarks, Calculations, Results:					
LMCS CHECK SAMPLE LMCS ID 25.3256 G 95.0% 23.9514 T 24.5430 W1 24.5343 W2 56.66 / 59.61					
Analyst-1 6B598	Analyst-2	Analyst-3	Analyst-4	Analyst-5 RE Bennett	
Hrs R.D. Hartung	Hrs	Hrs	Hrs	Hrs	
Date 11-28	Time Completed	Lab Unit Mgr Cga	OK		

\$4-6800-081 (R-10-83)

Serial No.	Sample Point		Date	Time Issued	Priority
F 5034.-5110	SEGMENT-3		11- 7-89	11:15	19
Determination	Method/Standard	Result Units	Charge Code	Recons	
% H ₂ O	LA-564-101	%	WB75L	0	
Sample Size	Customer ID				
?	89-041				
Remarks, Calculations, Results:					
DUPLICATE SAMPLE 24.5378 G Bottle # 47 23.7367 T 24.2248 W1 24.2247 W2 56.90% / 39.1%					
Analyst-1 6B598	Analyst-2	Analyst-3	Analyst-4	Analyst-5 RE Bennett	
Hrs R.D. Hartung	Hrs	Hrs	Hrs	Hrs	
Date 11-28	Time Completed	Lab Unit Mgr Cga	OK		

\$4-6800-081 (R-10-83)

9 1 1 2 0 5 9 0 7 1 3

Percent Water Analysis

Serial No. F 5053-5310	Sample Point SEGMENT-22		Date 11- 7-89	Time Issued 11:18	Priority 18
Determination % H ₂ O	Method/Standard LA-564-101	Result Units %	Charge Code WB75L	Returns 0	
Sample Size ?					Customer ID
Remarks, Calculations, Results: REAGENT BLANK 23.5575 G 23.5575 T - .0088 23.5487 W1 23.5488 W2					
Analyst-1 6.8598 R.O. Hale	Analyst-2	Analyst-3	Analyst-4	Analyst-5 XEBank	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 11-28	Time Completed	Lab Unit Mgr CG	DLM		

SI-8800-061 (R-10-83)

9 1 1 2 3 5 7 0 7 7 0

Fusion Dissolution

Serial No.	Sample Point	Date	Time Issued	Priority
F 5039.-6100	SEGMENT-8	11- 7-89	11:16	18
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results:				
DUPLICATE ANALYSIS GRAMS SAMPLE <u>6810</u> VOLUME ON COMPLETION <u>200 ml</u> 6810 200 = <u>3.40 g/ml</u> <u>3.40</u> -3 <u>200</u> → <u>0.6810 GRAMS</u>				
Analyst-1 Name	Analyst-2 Name	Analyst-3	11/27/89 8 12:19:48	
J.R.Smith	R.D.Hale	68598	11/27/89 8 12:19:48	
Hrs	Hrs	SL	Hrs	
11-22-89		R.D.Hale		
Date	Time Completed	Lab Unit Mgr		
C/S 54-800-081 (R-10-82)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 5038.-6000	SEGMENT-7	11- 7-89	11:16	18
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L-	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results:				
GRAMS SAMPLE 6059 VOLUME ON COMPLETION <u>200 ml</u> 6059 200 = <u>3.03 g/ml</u> <u>3.03</u> -3 <u>200</u> → <u>0.6059 GRAMS</u>				
Analyst-1 Name	Analyst-2 Name	Analyst-3	11/27/89 8 12:17:32	
J.R.Smith	R.D.Hale	68598	11/27/89 8 12:17:32	
Hrs	Hrs	SL	Hrs	
11-22-89		R.D.Hale		
Date	Time Completed	Lab Unit Mgr		
C/S 54-800-081 (R-10-82)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 5052.-6300	SEGMENT-21	11- 7-89	11:18	18
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L	WB75L	0
Sample Size		Customer ID		
?				
Remarks, Calculations, Results:				
REAGENT BLANK				
<i>Complete</i>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
68598				me
Hrs	Hrs	Hrs	Hrs	Hrs
R.D.Hale				
Date	Time Completed	Lab Unit Mgr		
11-22		C/S		
C/S 54-800-081 (R-10-82)				

9 1 1 2 3 5 3 0 7 7 1

Fusion Dissolution (second dissolution) Used for GEA Analysis

Serial No.	Sample Point	Date	Time issued	Priority
F 5039.-6100	SEGMENT-H	11- 7-89	11:16	26
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L	E21D1	0
Sample Size		Customer ID		
?				
Remarks, Calculations, Results:				
DUPLICATE ANALYSIS GRAMS SAMPLE <u>.4374</u> VOLUME ON COMPLETION <u>250 ml</u> SEQUENCE # : 138 MT 1: 41.7530 MT 2: 42.1344 NET WEIGHT: <u>1.73 x 10⁻³ g/ml</u> <u>1.73 g/l</u>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
<u>68590/kh</u>				
Hrs	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	<i>Zellie Paul Kathy Mandelke</i> SI-8000-001 (R-10-82)	
4-30-90				

Serial No.	Sample Point	Date	Time issued	Priority
F 5038.-6000	SEGMENT-B	11- 7-89	11:16	26
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L	E21D1	0
Sample Size		Customer ID		
?				
Remarks, Calculations, Results:				
GRAMS SAMPLE <u>.4796</u> VOLUME ON COMPLETION <u>250 ml</u> SEQUENCE # : 137 MT 1: 42.8147 MT 2: 43.2943 NET WEIGHT: <u>1.918 x 10⁻³ g/ml</u> <u>1.918 g/l</u>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
<u>68590/kh</u>				
Hrs	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	<i>Zellie Paul Kathy Mandelke</i> SI-8000-001 (R-10-82)	
4-30-90				

Serial No.	Sample Point	Date	Time issued	Priority
F 5052.-6300	SEGMENT-U	11- 7-89	11:18	26
Determination	Method/Standard	Result Units	Charge Code	Refuse
FUSION	LA-549-141	G/L	E21D1	0
Sample Size		Customer ID		
?				
REAGENT BLANK				
<i>complete</i>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
<u>68590/kh</u>				
Hrs	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	<i>Zellie Paul Kathy Mandelke</i> SI-8000-001 (R-10-82)	
4-30-90				

Total Alpha Analysis on the Fusion Dissolution

9 1 1 2 3 5 7 0 7 7 2

Serial No. F 5038.-6020	Sample Point SEGMENT-7		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination AT	Method/Standard LA-508-101	Result Units uCi/L	Charge Code WB75L	Returns 0	
Sample Size ?100-10 - 500					Customer ID
Remarks, Calculations, Results: 1.10 well					
Analyst-1 CC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 GW	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW			

SI-6000-061 (R-10-83)

Serial No. F 5037.-6520	Sample Point SEGMENT-6		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination AT	Method/Standard LA-508-101	Result Units % RECOVERY	Charge Code WB75L	Returns 0	
Sample Size ?10ml					Customer ID STO
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 83894					
Analyst-1 CC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 GW	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW			

SI-6000-061 (R-10-83)

9.511^{-3}
 $\sqrt{1.0003^{-2}}$
 K1
 95.10%

Serial No. F 5040.-6220	Sample Point SEGMENT-9		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination AT	Method/Standard LA-508-101	Result Units % RECOVERY	Charge Code WB75L	Returns 0	
Sample Size ?100-10 - 500					Customer ID
Remarks, Calculations, Results: SPIKE SAMPLE F5038 SPIKE ID 83894 SPIKE VOLUME 10ml Ratio of add to sample insuff.					
Analyst-1 CC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 GW	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW			

SI-6000-061 (R-10-83)

Serial No. F 5039.-6120	Sample Point SEGMENT-8		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination AT	Method/Standard LA-508-101	Result Units uCi/L	Charge Code WB75L	Returns 0	
Sample Size ?100-10 - 500					Customer ID GW
Remarks, Calculations, Results: DUPLICATE SAMPLE					
Analyst-1 CC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 GW	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW			

SI-6000-061 (R-10-83)

1.489 well

Total Alpha Analysis on the Fusion Dissolution

9 1 1 2 0 5 9 0 7 0 3

16/2

28 - 12
16

Alpha Calculation by EMB on 11-30-1989 at 03:05:21
Set #16 2-inch mount Alpha eff. : .22
Sample size : 1 mL Dilution : 202

29
16

Mount # 1

28

----- - 0.2 = 1.0753E+00 uCi/L alpha
10

Mount # 2

29

----- - 0.2 = 1.1167E+00 uCi/L alpha
10

16/3

454

Alpha Calculation by EMB on 11-30-1989 at 01:02:33
Set #16 2-inch mount Alpha eff. : .22
Sample size : 10 mL Dilution : 1

16 Mount # 1

479

454 ----- - 0.2 = 9.2547E-03 uCi/L alpha
10

Mount # 2

479

----- - 0.2 = 9.7666E-03 uCi/L alpha
10

F 5038.-6020

16/2

236 - 12
5

Alpha Calculation by EMB on 11-30-1989 at 03:50:51
Set #16 2-inch mount Alpha eff. : .22
Sample size : 1 mL Dilution : 202

244
5

Mount # 1

236 ----- - 0.2 = 1.9439E+01 uCi/L alpha
5

Mount # 2

244 ----- - 0.2 = 2.0101E+01 uCi/L alpha
5

16/2

15 - 12
5

Alpha Calculation by EMB on 11-30-1989 at 03:51:36
Set #16 2-inch mount Alpha eff. : .22
Sample size : 1 mL Dilution : 202

20
5

15 ----- - 0.2 = 1.4062E+00 uCi/L alpha
5

Mount # 2

20 ----- - 0.2 = 1.5717E+00 uCi/L alpha
5

F 5039.-6120

9 1 1 2 2 3 2 0 7 7 4

Total Alpha Analysis on the Fusion Dissolution

Serial No F 5052-6320	Sample Point SEGMENT-21		Date 11- 7-89	Time Issued 11:18	Priority 1B
Determination AT	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Reruns 0	
Sample Size ?					Customer ID
Remarks Calculations, Results: REAGENT BLANK					
$\sim 2.4 \times 10^{-4}$ mCi/l					
Analyst-1 CCZ6.5	Analyst-2	Analyst-3	Analyst-4	Analyst-5 GW	
CPA 444 = 5.0000	Hrs	Hrs	Hrs	Hrs	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW	KH		

54-6300-061 (R-1C-83)

Serial No F 5041-6320	Sample Point SEDIMENT 10		Date 11- 7-89	Time Issued 11:16	Priority 1B
Determination AT	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size ?					Customer ID
Remarks Calculations, Results: LMDS CHECK SAMPLE LMDS ID 83844					
86.4%					
Analyst-1 CCZ6.5	Analyst-2	Analyst-3	Analyst-4	Analyst-5 GW	
CPA 444 = 5.0000	Hrs	Hrs	Hrs	Hrs	
Date 11-29-89	Time Completed	Lab Unit Mgr CJW	KH		

54-6300-061 (R-1C-83)

Total Alpha Analysis on the Fusion Dissolution

9 1 1 2 7 5 9 0 7 7 5

16/2

$$\frac{8}{10} - .6$$

$$\frac{5}{10}$$

Alpha Calculation by DM on 11-29-1989 at 21:53:41
Det #16 2-inch mount Alpha eff. : .22
Sample size : 1 ea Dilution : 1

Mount # 1

$$\frac{8}{10} - 0.6 < 1.6736E-06 \mu\text{Ci}/\text{ea alpha}$$

Mount # 2

$$\frac{5}{10} - 0.6 < 1.1686E-06 \mu\text{Ci}/\text{ea alpha}$$

F 5052-6320

16/2

$$\frac{434}{10}$$

Alpha Calculation by EMR on 11-30-1989 at 02:52:09
Det #16 2-inch mount Alpha eff. : .22
Sample size : 10 mL Dilution : 1

Mount # 1

$$\frac{434}{10} - 0.2 = 8.8452E-03 \mu\text{Ci}/\text{L alpha}$$

$$\frac{414}{10}$$

Mount # 2

$$\frac{414}{10} - 0.2 = 8.4357E-03 \mu\text{Ci}/\text{L alpha}$$

F 5041-6520

Total Beta Analysis on the Fusion Dissolution

9 1 1 2 5 9 0 7 7 6

Serial No. F 5038.-6025	Sample Point SEGMENT-7	Date 11-7-89	Time Issued 11:16	Priority 19
Determination TB	Method/Standard LA-54B-101	Result Units uCi/L	Charge Code WB75L	Reruns 0
Sample Size ? 100-10 - 500	Customer ID 89-041			
Remarks, Calculations, Results: 3.283 ³ well				
Analyst-1 LC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 me
<i>Mayer</i> J. Frazee				
Date 11-7-89	Time Completed	Lab Unit Mgr Chay		

54-6000-061 (R-10-83)

Serial No. F 5037.-6525	Sample Point SEGMENT-6	Date 11-7-89	Time Issued 11:16	Priority 19
Determination TH	Method/Standard LA-54B-101	Result Units % RECOVERY	Charge Code WB75L	Reruns 0
Sample Size ? 100-10 - 500	Customer ID STD 89-041			
Remarks, Calculations, Results: LMDS CHECK SAMPLE LMDS ID 83842 K1				
Analyst-1 LC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 me
<i>Mayer</i> J. Frazee				
Date 11-7-89	Time Completed	Lab Unit Mgr Chay		

54-6000-061 (R-10-83)

1.346⁻¹ / 1.411 - 1 95.40%

Serial No. F 5040.-6225	Sample Point SEGMENT-9	Date 11-7-89	Time Issued 11:16	Priority 19
Determination TB	Method/Standard LA-54B-101	Result Units % RECOVERY	Charge Code WB75L	Reruns 0
Sample Size ? 100-10 - 500	Customer ID 89-041			
Remarks, Calculations, Results: SPIKE SAMPLE F 5038 SPIKE ID 83842 SPIKE VOLUME 1000				
Analyst-1 LC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 me
<i>Mayer</i> J. Frazee				
Date 11-7-89	Time Completed	Lab Unit Mgr Chay		

54-6000-061 (R-10-83)

3.0976³ ratio of std to sample concn?

Serial No. F 5039.-6125	Sample Point SEGMENT-8	Date 11-7-89	Time Issued 11:16	Priority 19
Determination TB	Method/Standard LA-54B-101	Result Units uCi/L	Charge Code WB75L	Reruns 0
Sample Size ? 100-10 - 500	Customer ID 89-041			
Remarks, Calculations, Results: DUPLICATE SAMPLE				
Analyst-1 LC269	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 me
<i>Mayer</i> J. Frazee				
Date 11-7-89	Time Completed	Lab Unit Mgr Chay		

54-6000-061 (R-10-83)

3.88³ well

Total Beta Analysis on the Fusion Dissolution

9 1 1 2 0 5 9 0 7 0 7

16/2

DF = 202

16/2

106821

10

Beta-Calculation by EMB on 11-30-1989 at 03:05:19
Bet #16 2-inch mount Beta eff.: .291
Sample size : 1 mL Dilution : 202

103512

10

Count # 1

106821
----- - 17.0 = 3.3348E+03 uCi/L beta

10

Count # 2

3.28³

103512

10

----- - 17.0 = 3.2313E+03 uCi/L beta

F 5038,-6025

61969

5

Beta-Calculation by EMB on 11-30-1989 at 03:51:34
Bet #16 2-inch mount Beta eff.: .291
Sample size : 1 mL Dilution : 202

62275

5

Count # 1

61969
----- - 17.0 = 3.8700E+03 uCi/L beta

5

Count # 2

62275
----- - 17.0 = 3.8892E+03 uCi/L beta

5

F 5039,-6125

16/2

497492

5

Beta-Calculation by EMB on 11-30-1989 at 03:50:49
Bet #16 2-inch mount Beta eff.: .291
Sample size : 1 mL Dilution : 202

Count # 1

497492

5

Count # 2

497492
----- - 17.0 = 3.0898E+03 uCi/L beta

5

F 5040 - 6225

16/2

8856

10

Beta-Calculation by EMB on 11-30-1989 at 01:02:31
Bet #16 2-inch mount Beta eff.: .291
Sample size : 10 mL Dilution : 1

Count # 1

8856
----- - 17.0 = 1.3445E-01 uCi/L beta

10

Count # 2

8879
----- - 17.0 = 1.3481E-01 uCi/L beta

10

F 5037,-6525

9 1 1 2 3 5 9 0 7 0 8

Total Beta Analysis on the Fusion Dissolution

5-16

Serial No. F 5052,-6325	Sample Point SEGMENT-21	Date 11- 7-89	Time issued 11:18	Priority 18
Determination TB	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Retuns 0
Sample Size ?	Customer ID 89-041			

Remarks, Calculations, Results:
REAGENT BLANK

~~25.58 -4~~
~~36.97~~ ~~8%~~
 will will

Analyst-1 4.C26G	Analyst-2	Analyst-3	Analyst-4	Analyst-5 're'
Hrs 4.00	Hrs	Hrs	Hrs	Hrs
Date 11-29-89	Time Completed	Lab Unit Mgr CJW		

SI-6500-061 (R-1C-83)

Serial No. F 5041,-6525	Sample Point SEGMENT-10	Date 11- 7-89	Priority 19	
Determination TB	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WR75L	
Sample Size ? 10ml	Customer ID 89-041			

Remarks, Calculations, Results:
LMCS CHECK SAMPLE
LMCS ID_83P44

Analyst-1 4.C26G	Analyst-2	Analyst-3	Analyst-4	Analyst-5 're'
Hrs 4.00	Hrs	Hrs	Hrs	Hrs
Date 11-29-89	Time Completed	Lab Unit Mgr CJW		

SI-6500-061 (R-1C-83)

Total Beta Analysis on the Fusion Dissolution

9 1 1 2 3 5 7 0 7 7 9

6/2

189
10 - 10

180
10

Beta Calculation by DM on 11-29-1989 at 21:53:39
Det #16 2-inch mount Beta eff. : .291
Sample size : 1 ea Dilution : 1

Mount # 1

189
----- - 18.0 < 6.2994E-06 uCi/ea beta
10

Mount # 2

180
----- - 18.0 < 4.8461E-06 uCi/ea beta
10

*Slo. 97 - S KF
new/0*

6/2

8760
10

- 8663 Calculation by EMB on 11-30-1989 at 02:52:07
Det #16 2-inch mount Beta eff. : .291
Sample size : 10 mL Dilution : 1

Mount # 1

8760
----- - 17.0 = 1.3297E-01 uCi/L beta
10

Mount # 2

8663
----- - 17.0 = 1.3456E-01 uCi/L beta
10

F 5052,-6325

9 1 1 2 2 5 7 0 3 7 0

Gamma Energy Analysis of the Fusion Dissolution

3617

Serial No. F 5038.-6030	Sample Point SEGMENT-7		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units UCl/L	Charge Code WB75L	Remarks, Calculations, Results: Second fusion $Cs^{137} 9.29 \times 10^1$ mCi/l RERUN Analyst-1 68598111 Hrs Date 5-2-90	
Sample Size ? 200 & 500		Customer ID			

3615

Serial No. F 5037.-6530	Sample Point SEGMENT-6		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units % RECOVERY	Charge Code WB75L	Remarks, Calculations, Results: LMCS CHECK, SAMPLE LMCS ID 2244 $Co^{60} 2.13 \times 10^1 / 2.0566 \times 10^1$ 103.6% RERUN $Cs^{137} 2.97 \times 10^1 / 2.889 \times 10^1$ 102.8% Analyst-1 68598111 Hrs Date 5-2-90	
Sample Size ? 500		Customer ID			

3619

Serial No. F 5040.-6230	Sample Point SEGMENT-9		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units % RECOVERY	Charge Code WB75L	Remarks, Calculations, Results: Second fusion SPIKE SAMPLE F5038 SPIKE ID 122044 SPIKE VOLUME 200 $Co^{60} (1.06 \times 10^2 - 0) \times 2 = 2.12 \times 10^1$ REP 106.0% $Cs^{137} (2.46 \times 10^2) - (9.29 \times 10^1) = 15.31 \times 2 = 30.62 \times 10^1$ Analyst-1 68598111 Hrs Date 5-2-90	
Sample Size ? 200 & 500 Spike		Customer ID			

3618

Serial No. F 5039.-6130	Sample Point SEGMENT-8		Date 11- 7-89	Time Issued 11:16	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units UCl/L	Charge Code WB75L	Remarks, Calculations, Results: DUPLICATE SAMPLE Second fusion $Cs^{137} 2.47 \times 10^1$ mCi/l RERUN Analyst-1 68598111 Hrs Date 5-2-90	
Sample Size ? 200 & 100		Customer ID			

9 1 1 2 3 5 9 0 3 0 1

Gamma Energetics Analysis of the Fusion Dissolution

3621																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Sample No. F 5041--6530</td> <td colspan="2">Sample Point SEGMENT-10</td> <td>Date 11-7-89</td> <td>Time Issued 11:16</td> <td>Priority 19</td> </tr> <tr> <td>Determination GEA</td> <td>Method/Standard LA-348-121</td> <td>Result Units % RECOVERY</td> <td>Charge Code WB7SL</td> <td colspan="2">Reruns 0</td> </tr> <tr> <td colspan="3">Sample Size ? 500L</td> <td colspan="3">Customer ID</td> </tr> <tr> <td colspan="6">Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID...122344</td> </tr> <tr> <td colspan="6" style="text-align: center;">Co^{60} $2.18 \times 10^3 / 2.0566 \times 10^3$ RERUN 106.0%</td> </tr> <tr> <td colspan="6" style="text-align: center;">Cs^{137} $2.95 \times 10^3 / 2.8888 \times 10^3$ 102.1%</td> </tr> <tr> <td>Analyst-1 HRS 6.8595/1011</td> <td>Analyst-2 HRS</td> <td>Analyst-3 HRS</td> <td>Analyst-4 HRS Kathy Mandible</td> <td>Analyst-5 HRS</td> <td></td> </tr> <tr> <td>Date 5-2-90</td> <td>Time Completed</td> <td>Lab Unit Mgr</td> <td colspan="3"><i>Kathy Mandible</i> SA-4000-001 (P-10-83)</td> </tr> </table>						Sample No. F 5041--6530	Sample Point SEGMENT-10		Date 11-7-89	Time Issued 11:16	Priority 19	Determination GEA	Method/Standard LA-348-121	Result Units % RECOVERY	Charge Code WB7SL	Reruns 0		Sample Size ? 500L			Customer ID			Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID...122344						Co^{60} $2.18 \times 10^3 / 2.0566 \times 10^3$ RERUN 106.0%						Cs^{137} $2.95 \times 10^3 / 2.8888 \times 10^3$ 102.1%						Analyst-1 HRS 6.8595/1011	Analyst-2 HRS	Analyst-3 HRS	Analyst-4 HRS Kathy Mandible	Analyst-5 HRS		Date 5-2-90	Time Completed	Lab Unit Mgr	<i>Kathy Mandible</i> SA-4000-001 (P-10-83)		
Sample No. F 5041--6530	Sample Point SEGMENT-10		Date 11-7-89	Time Issued 11:16	Priority 19																																																
Determination GEA	Method/Standard LA-348-121	Result Units % RECOVERY	Charge Code WB7SL	Reruns 0																																																	
Sample Size ? 500L			Customer ID																																																		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID...122344																																																					
Co^{60} $2.18 \times 10^3 / 2.0566 \times 10^3$ RERUN 106.0%																																																					
Cs^{137} $2.95 \times 10^3 / 2.8888 \times 10^3$ 102.1%																																																					
Analyst-1 HRS 6.8595/1011	Analyst-2 HRS	Analyst-3 HRS	Analyst-4 HRS Kathy Mandible	Analyst-5 HRS																																																	
Date 5-2-90	Time Completed	Lab Unit Mgr	<i>Kathy Mandible</i> SA-4000-001 (P-10-83)																																																		
3616																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Sample No. F 5052--6530</td> <td colspan="2">Sample Point SEGMENT-21</td> <td>Date 11-7-89</td> <td>Time Issued 11:18</td> <td>Priority 18</td> </tr> <tr> <td>Determination GEA</td> <td>Method/Standard LA-348-121</td> <td>Result Units uCi/L</td> <td>Charge Code WB7SL</td> <td colspan="2">Reruns 0</td> </tr> <tr> <td colspan="3">Sample Size ? 1ml</td> <td colspan="3">Customer ID</td> </tr> <tr> <td colspan="6">Remarks, Calculations, Results: REAGENT BLANK</td> </tr> <tr> <td colspan="6" style="text-align: center;">$\text{Cs}^{137} < 6.32 \times 10^{-2}$ mCi/L Second fusion RERUN</td> </tr> <tr> <td>Analyst-1 HRS 6.8595/1011</td> <td>Analyst-2 HRS</td> <td>Analyst-3 HRS</td> <td>Analyst-4 HRS</td> <td>Analyst-5 HRS</td> <td></td> </tr> <tr> <td>Date 5-2-90</td> <td>Time Completed</td> <td>Lab Unit Mgr</td> <td colspan="3"><i>Kathy Mandible</i> SA-4000-001 (P-10-83)</td> </tr> </table>						Sample No. F 5052--6530	Sample Point SEGMENT-21		Date 11-7-89	Time Issued 11:18	Priority 18	Determination GEA	Method/Standard LA-348-121	Result Units uCi/L	Charge Code WB7SL	Reruns 0		Sample Size ? 1ml			Customer ID			Remarks, Calculations, Results: REAGENT BLANK						$\text{Cs}^{137} < 6.32 \times 10^{-2}$ mCi/L Second fusion RERUN						Analyst-1 HRS 6.8595/1011	Analyst-2 HRS	Analyst-3 HRS	Analyst-4 HRS	Analyst-5 HRS		Date 5-2-90	Time Completed	Lab Unit Mgr	<i>Kathy Mandible</i> SA-4000-001 (P-10-83)								
Sample No. F 5052--6530	Sample Point SEGMENT-21		Date 11-7-89	Time Issued 11:18	Priority 18																																																
Determination GEA	Method/Standard LA-348-121	Result Units uCi/L	Charge Code WB7SL	Reruns 0																																																	
Sample Size ? 1ml			Customer ID																																																		
Remarks, Calculations, Results: REAGENT BLANK																																																					
$\text{Cs}^{137} < 6.32 \times 10^{-2}$ mCi/L Second fusion RERUN																																																					
Analyst-1 HRS 6.8595/1011	Analyst-2 HRS	Analyst-3 HRS	Analyst-4 HRS	Analyst-5 HRS																																																	
Date 5-2-90	Time Completed	Lab Unit Mgr	<i>Kathy Mandible</i> SA-4000-001 (P-10-83)																																																		

9 1 1 2 3 5 7 0 2 9 2

Uranium Analysis of the Fusion Dissolution

Serial No. F 5038-6040	Sample Point SEGMENT-7		Date 11-7-89	Time Issued 11:16	Priority 23
Determination U	Method/Standard LA-925-106	Result Units G/L	Charge Code WB75L	Return 0	Customer ID 89-041
Sample Size ? 100-10-100			Customer ID 89-041		
Remarks, Calculations, Results: $\text{Spk Vol: } 100\text{A}$ $\text{Spk ID: } 5.62^{-4}$ Sample: .04 Imp+Spk: .38 $\frac{(04)(5.62^{-4})(.1)(1010)}{.38-.04} = \cancel{10.08}^{3.33}$					
Analyst-1 6C269 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 11-28-89	Time Completed	Lab Unit Mgr CGV			
34-6800-061 (R-10-83)					

Serial No. F 5037-6040	Sample Point SEGMENT-6		Date 11-7-89	Time Issued 11:16	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Return 0	Customer ID 58B38
Sample Size ? 100-10-100			Customer ID 58B38		
Remarks, Calculations, Results: LMCS CHECK SAMPLE $\text{LMCS ID: } 58B38$ $\text{Spk Vol: } 100\text{A}$ $\text{Spk ID: } 5.62^{-4}$ Sample: .17 Imp+Spk: .46 $\frac{(.17)(5.62^{-4})(.1)(1010)}{.46-.17} = \cancel{3.33}^{2.99}$					
Analyst-1 6C269 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 11-28-89	Time Completed	Lab Unit Mgr CGV			
34-6800-061 (R-10-83)					

Serial No. F 5040-6240	Sample Point SEGMENT-9		Date 11-7-89	Time Issued 11:16	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Return 0	Customer ID 89-041
Sample Size ? 100-10-100			Customer ID 89-041		
Remarks, Calculations, Results: SPIKE SAMPLE $\text{SPIKE ID: } 58B38$ $\text{SPKE VOLUME } 100 \times 10^{-100}$ $\text{Spk Vol: } 100\text{A}$ $\text{Spk ID: } 5.62^{-4}$ Sample: .22 Imp+Spk: .52 $\frac{(.22)(5.62^{-4})(.1)(1010)}{.52-.22} = \cancel{4.16}^{3.49}$ $4.16^{-2} / 6.68^{-3} = \cancel{3.49}^{2.99} \cdot .52 - .22$					
Analyst-1 6C269 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 11-28-89	Time Completed	Lab Unit Mgr CGV			
34-6800-061 (R-10-83)					

Serial No. F 5039-6140	Sample Point SEGMENT-8		Date 11-7-89	Time Issued 11:16	Priority 23
Determination U	Method/Standard LA-925-106	Result Units G/L	Charge Code WB75L	Return 0	Customer ID 89-041
Sample Size ? 100-10-100			Customer ID 89-041		
Remarks, Calculations, Results: DUPLICATE SAMPLE $\text{Spk Vol: } 100\text{A}$ $\text{Spk ID: } 5.62^{-4}$ Sample: .06 Imp+Spk: .40 $\frac{(.06)(5.62^{-4})(.1)(1010)}{.40-.06} = \cancel{1.06}^{.81}$					
Analyst-1 6C269 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 11-28-89	Time Completed	Lab Unit Mgr CGV			
34-6800-061 (R-10-83)					

9 1 - 12 15 9 0 3 2 3

Uranium Analysis of the Fusion Dissolution

Serial No. F 5052--6340	Sample Point SEGMENT-21	Date 11- 7-89	Time Issued 11:18	Priority 18
Determination U	Method/Standard LA-925-106	Result Units G/L	Charge Code WB7EL	Recurve 0
Sample Size <i>? direct</i>				Customer ID
Remarks, Calculations, Results:				
REAGENT BLANK $\text{Spk Vol: } 100\lambda$ $\text{Spk TD: } 5.113-5$ $\frac{(0.02)(5.113-5)(1)(1)}{30} = 1.83 \text{ S/l}$ $\cancel{3.05} \rightarrow \cancel{81\lambda} \\ \cancel{81\lambda} : .02 \\ \cancel{B1\lambda + Spk} : .30$				
Analyst-1 <i>GC769</i> <i>11-28-89</i>	Analyst-2 <i>11-28-89</i>	Analyst-3 <i>Hrs</i>	Analyst-4 <i>Hrs</i>	Analyst-5 <i>Hrs</i>
Date 11-28-89	Time Completed	Lab Unit Mgr <i>Chase</i>		

Serial No.	Sample Point	Date	Time Issued	Priorit 23
F 5041-6540	SEGMENT-10	11- 7-89	11:15	
Determination	Method/Standard	Result Units	Charge Code	Requ 0
U	LA-925-106	% RECOVERY	WB75L	
Sample Size				Customer ID
? 100-10-100	58838		50	
Remarks, Calculations, Results	<p>LMCS CHECK SAMPLE</p> <p>LMCS ID 58838</p> <p>$\text{SPK Val: } 100 \lambda$</p> <p>$\text{SPK Val: } 5.62^{-4}$</p> <p>$\frac{(10)(5.62^{-4})(.1)(100)}{.46} = 30.2$</p> <p>$30.2 - .16 = 29.8$</p>			
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
LC265	Mrs	Mrs	Mrs	Mrs
11-28-89	Time Completed	Lab Unit Mgr	<i>John</i>	

Uranium Analysis of the Fusion Dissolution

F 5037.-6540

$$\cancel{3.34 \times 10^{-2}} / \cancel{3.99 \times 10^{-2}} = 108.37\%$$

F 5052.-6340

$$\cancel{(1.0)(9.0)} \cancel{(5.68^{-2})} / \cancel{(0.2)} = 3.59 \times 10^{-4}$$

F 5041.-6540

$$\cancel{3.39 \times 10^{-2}} / \cancel{3.99 \times 10^{-2}} = 98.67\%$$

F 5040.-6240

$$\cancel{3.39 \times 10^{-2}} / \cancel{3.99 \times 10^{-2}} = 113.37\%$$

$$\cancel{(9.9^{-4})(5.68^{-2})} / \cancel{(0.16)(9.9^{-4})(5.68^{-2})} = \cancel{2.95 \times 10^{-2}}$$

$$\cancel{(0.22)(9.9^{-4})} \cancel{(5.68^{-2})} - \cancel{(9.9^{-4})(6.55^{-3})} = \cancel{(9.9^{-4})}$$

2 1 1 2 0 5 9 0 3 7 5

Uranium Analysis of the Fusion Dissolution

$$\frac{(.06)(9.9^{-4})(5.68^{-2})}{(9.9^{-4})[(.40)(\frac{5.7}{5.6}) - (.06)]} = 9.82 \times 10^{-3} \text{ g/l}$$

CJL

$$\frac{(.04)(9.9^{-4})(5.68^{-2})}{(9.9^{-4})[(.38)(\frac{5.7}{5.6}) - (.04)]} = 6.55 \times 10^{-3} \text{ g/l}$$

CJL

F 5039.-6140

F 5038.-6040

9 1 1 2 7 5 7 0 9 7 6

Serial No.	Sample Point	Date	Time Issued	Priority
F 5054.-7100	SEGMENT-13	11-7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Returns
H ₂ O-DGST	LA-504-101	G/L	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: IMPLICATE ANALYSIS GRAMS SAMPLE 8971 VOLUME ON COMPLETION 100 ml DF 111.47 8.97 ³ g/ml				
Analyst-1 80725 J.C.	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Assigned Chemist J.P. 9171
Date 12-4-89	Time Completed	Lab Unit Mgr CJW		
34-000-061 (R-10-83)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 5047.-7000	SEGMENT-12	11-7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Returns
H ₂ O-DGST	LA-504-101	G/L	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: GRAMS SAMPLE 9587 VOLUME ON COMPLETION 100 ml Bottle # 080 350 ml 29.98 22.35 7.63 g				
Analyst-1 J.R. 5286	Analyst-2 D.K. 66731	Analyst-3 Hrs	Analyst-4 Hrs	Assigned Chemist J.P. 9171
Date 11-3-89	Time Completed	Lab Unit Mgr CJW		
34-000-061 (R-10-83)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 5054.-7300	SEGMENT-23	11-7-89	11:18	18
Determination	Method/Standard	Result Units	Charge Code	Returns
H ₂ O-DGST	LA-504-101	G/L	WR75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: REAGENT BLANK Complete				
Analyst-1 80725 J.C.	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Assigned Chemist J.P. 9171
Date 12-4-89	Time Completed	Lab Unit Mgr CJW		
34-000-061 (R-10-83)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 5045.-7200	SEGMENT-14	11-7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Returns
H ₂ O-DGST	LA-504-101	% RECOVERY	WB75L	0
Sample Size		Customer ID		
?		89-041		
Remarks, Calculations, Results: SPiked ANALYSIS GRAMS SAMPLE 8.73 VOLUME ON COMPLETION 100 ml VOLUME SPIKE SPIKE ID DF 114.53 8.73 ³ g/ml				
Analyst-1 80725 J.C.	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Assigned Chemist J.P. 9171
Date 12-4-89	Time Completed	Lab Unit Mgr CJW		
34-000-061 (R-10-83)				

F 5043.-7000

** JET_WT ** REV 0 **

SEQUENCE # : 6
WT 1: 100.0650
WT 2: 100.9360

NET WEIGHT:

--> 0.8731 GRAMS

12-01-89 @ 14:51:13

F 5044.-7100

** JET_WT ** REV 0 **

SEQUENCE # : 5
WT 1: 99.8121
WT 2: 99.7192

NET WEIGHT:

--> 0.8971 GRAMS

12-01-89 @ 14:43:06

F 5045.-7200

** REV 0 **

SEQUENCE # : 4
WT 1: 99.9723
WT 2: 99.9318

NET WEIGHT:

--> 0.9587 GRAMS

12-01-89 @ 14:53:14

2 0 8 0 6 8 0 7 1 1 6

Water Digestion II

9 1 1 2 3 5 7 0 3 0 8

Serial No. F 5044-7100	Sample Point SEGMENT-M	Date 11-7-89	Time Issued 11:17	Priority 26
Determination H2O-DGST	Method/Standard LA-504-101	Result Units G/L	Charge Code E21D1	Return 0
Sample Size DIRECT	Customer ID <i>2nd batch</i>			
Remarks, Calculations, Results: DUPLICATE ANALYSIS GRAMS SAMPLE <u>6013.9</u> VOLUME ON COMPLETION <u>50 mL</u> <u>1.203 E -2 g/ml</u> <u>1.203 E 1 g/l</u>				
Analyst-1 <u>80028</u> Hrs	Analyst-2 <i>Ed Cohn</i> Hrs	Analyst-3 <i>John Smith</i> Hrs	Analyst-4 <i>John Smith</i> Hrs	Analyst-5
Date <u>6-27-90</u>	Time Completed	Lab Unit/Mgr <i>Ed Cohn</i>	Lab Unit/Mgr <i>John Smith</i>	Lab Unit/Mgr <i>John Smith</i>

SI-8000-001 (R-10-82)

Serial No. F 5043-7000	Sample Point SEGMENT-L	Date 11-7-89	Time Issued 11:17	Priority 26
Determination H2O-DGST	Method/Standard LA-504-101	Result Units G/L	Charge Code E21D1	Return 0
Sample Size DIRECT	Customer ID <i>2nd batch</i>			
Remarks, Calculations, Results: GRAMS SAMPLE <u>.4802.9</u> VOLUME ON COMPLETION <u>50 mL</u> <u>9.604 E -3 g/ml</u> <u>9.604 g/l</u>				
Analyst-1 <u>80028</u> Hrs	Analyst-2 <i>Ed Cohn</i> Hrs	Analyst-3 <i>John Smith</i> Hrs	Analyst-4 <i>John Smith</i> Hrs	Analyst-5
Date <u>6-27-90</u>	Time Completed	Lab Unit/Mgr <i>Ed Cohn</i>	Lab Unit/Mgr <i>John Smith</i>	Lab Unit/Mgr <i>John Smith</i>

SI-8000-001 (R-10-82)

Serial No. F 5043-7200	Sample Point SEGMENT-N	Date 11-7-89	Time Issued 11:17	Priority 26
Determination H2O-DGST	Method/Standard LA-504-101	Result Units % RECOVERY	Charge Code E21D1	Return 0
Sample Size ?	Customer ID <i>2nd batch</i>			
Remarks, Calculations, Results: SPIKED ANALYSIS GRAMS SAMPLE <u>.5213g</u> VOLUME ON COMPLETION <u>50mL</u> VOLUME SPIKE <u>.50mL</u> SPIKE ID <u>1.043 E -2 g/ml</u> <u>1.043 E 1 g/l</u>				
Analyst-1 <u>80028</u> Hrs	Analyst-2 <i>Ed Cohn</i> Hrs	Analyst-3 <i>John Smith</i> Hrs	Analyst-4 <i>John Smith</i> Hrs	Analyst-5
Date <u>6-27-90</u>	Time Completed	Lab Unit/Mgr <i>Ed Cohn</i>	Lab Unit/Mgr <i>John Smith</i>	Lab Unit/Mgr <i>John Smith</i>

SI-8000-001 (R-10-82)

9 1 1 2 7 5 7 0 9 9 9

Ion Chromatographic Analysis of the Water Digestion

Serial No F 5045.-7271	Sample Point SEGMENT-14	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination F	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Returns 0
Sample Size ?	Customer ID 89-041			
100-10 ml				
Remarks, Calculations, Results: SPIKE SAMPLE R5043 F 49.79 ppm SPIKE ID 4009-A SPIKE VOLUME 50µL $\frac{(5.05)}{5.00} \left(222.6 \text{ ppm} \right) - (221.8) \left(\frac{1.73}{5.58} \right) \times 100 = 45.6\%$ $\frac{(.050)}{5.05} (49.79) (101)$ Spike L 20% Not calculated				
Analyst-1 6B107	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>Chambers</i>
14:30				
Date 12-5-89	Time Completed	Lab Unit Mgr <i>Tellie Paul Kishikawa</i>		
SI-6000-061 (R-10-83)				

Serial No F 5044.-7171	Sample Point SEGMENT-13	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination F	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Returns 0
Sample Size ?	Customer ID 89-041			
Remarks, Calculations, Results: DUPLICATE SAMPLE 1.83^2 ppm				
Analyst-1 6B107	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>Chambers</i>
14:30				
Date 12-5-89	Time Completed	Lab Unit Mgr <i>Chambers</i>		
SI-6000-061 (R-10-83)				

Serial No F 5042.-7571	Sample Point SEGMENT-11	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination F	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Returns 0
Sample Size 100-10	Customer ID 89-041			
Remarks, Calculations, Results: LMDS CHECK SAMPLE LMDS ID 6C11HC $62.98 / 62.0$ 101.6%				
Analyst-1 6B107	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>Chambers</i>
14:30				
Date 12/5/89	Time Completed	Lab Unit Mgr <i>Chambers</i>		
SI-6000-061 (R-10-83)				

Serial No F 5043.-7071	Sample Point SEGMENT-12	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination F	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Returns 0
Sample Size ?	Customer ID 89-041			
Remarks, Calculations, Results: 2.22^2 ppm				
Analyst-1 6B107	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>Chambers</i>
14:30				
Date 12/5/89	Time Completed	Lab Unit Mgr <i>Chambers</i>		
SI-6000-061 (R-10-83)				

Ion Chromatographic Analysis of the Water Digestion

9 1 1 2 7 5 9 0 3 1 1

Ion Chromatographic Analysis of the Water Digestion

Serial No. F 5045.-7272	Sample Point SEGMENT-14	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination CL	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 0
Sample Size ?	Customer ID 89-041			
<i>100-10 - A1 = 10ml</i>				
Remarks, Calculations, Results: SPIKE SAMPLE R5043 C1 = 59.78 ppm. SPIKE ID 40C9-A SPIKE VOLUME 5.04L				
$\frac{(1.01)(68.58) - (0.0)(91)}{(1.05)(59.78) / (1.01)} \times 100 = 115.9\%$				
Analyst-1 6B107 Hrs 14:30	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12-5-89	Time Completed	Lab Unit/Mgr <i>Kathy Mandelbo</i>	Signature/Changer <i>John H. Johnson</i>	

Serial No. F 5054.-7372	Sample Point SEGMENT-23	Date 11- 7-89	Time Issued 11:18	Priority 18
Determination CL	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 0
Sample Size ?	Customer ID 89-041			
<i>Direct</i>				
Remarks, Calculations, Results: REAGENT BLANK.				
<i>2.1 ppm</i>				
Analyst-1 6B107 Hrs 14:30	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12-5-89	Time Completed	Lab Unit/Mgr <i>Cja</i>	Signature/Changer <i>John H. Johnson</i>	

Serial No. F 5044.-7172	Sample Point SEGMENT-13	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination CL	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 0
Sample Size 100-10	Customer ID			
Remarks, Calculations, Results: DUPLICATE SAMPLE				
<i>1.01x10³ ppm</i>				
Analyst-1 6B107 Hrs 14:30	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12-5-89	Time Completed	Lab Unit/Mgr <i>Kathy Mandelbo</i>	Signature/Changer <i>John H. Johnson</i>	

Serial No. F 5046.-7572	Sample Point SEGMENT-15	Date 11- 7-89	Time Issued 11:18	Priority 19
Determination CL	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 0
Sample Size 100-10	Customer ID			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 6C 11 HC				
$\frac{7.65}{75} \times 100 = 102\%$				
Analyst-1 6B107 Hrs 14:30	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12-5-89	Time Completed	Lab Unit/Mgr <i>Cja</i>	Signature/Changer <i>John H. Johnson</i>	

9 1 1 2 2 5 3 0 0 1 2

Ion Chromatographic Analysis of the Water Digestion

Serial No.	Sample Point		Date	Time Issued	Priority
F 5043.-7073	SEGMENT-12		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
N03	LA-533-105	PPM	WB75L	0	
Sample Size			Customer ID		
<i>? 100-0</i>		<i>89-041</i>			
Remarks, Calculations, Results:					
<i>6.22² ppm</i>					
Analyst-1 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>6B107</i> <i>1430</i>				<i>6B107 7-17</i>	
Date 12/5/89	Time Completed	Lab Unit Mgr. <i>CJG</i>			
54-6000-001 (R-10-83)					

Serial No.	Sample Point		Date	Time Issued	Priority
F 5042.-7573	SEGMENT-11		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
N03	LA-533-105	% RECOVERY	WB75L	0	
Sample Size			Customer ID		
<i>100-10</i>		<i>89-041</i>			
Remarks, Calculations, Results:					
<i>LMCS CHECK SAMPLE</i> <i>LMCS ID 6C1LMC</i>					
Analyst-1 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>6B107</i> <i>1430</i>				<i>6B107 7-17</i>	
Date 12/5/89	Time Completed	Lab Unit Mgr. <i>CJG</i>		<i>bfp</i>	
54-6000-001 (R-10-83)					

Serial No.	Sample Point		Date	Time Issued	Priority
F 5045.-7273	SEGMENT-14		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
N03	LA-533-105	% RECOVERY	WB75L	0	
Sample Size			Customer ID		
<i>? 100uL-10 ml</i>		<i>89-041</i>			
Remarks, Calculations, Results:					
<i>R 5043 NO₃-500 ppm</i> <i>SPIKE SAMPLE R 5043 NO₃-500 ppm</i> <i>SPIKE ID 4DC9-A</i> <i>SPIKE VOLUME 5.00 ml</i>					
$\frac{(1.01)(10.2) - (6.22)(.91)}{(0.50)(500)} \times 100 = 94.25\%$					
Analyst-1 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>6B107</i> <i>1430</i>				<i>6B107 7-17</i>	
Date 12-5-89	Time Completed	Lab Unit Mgr. <i>CJG bfp</i>			
54-6000-001 (R-10-83)					

Serial No.	Sample Point		Date	Time Issued	Priority
F 5044.-7173	SEGMENT-13		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
N03	LA-533-105	PPM	WB75L	0	
Sample Size			Customer ID		
<i>? 100-0</i>		<i>89-041</i>			
Remarks, Calculations, Results:					
<i>DUPLICATE SAMPLE</i>					
Analyst-1 Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>6B107</i> <i>1430</i>				<i>6B107 7-17</i>	
Date 12/5/89	Time Completed	Lab Unit Mgr. <i>CJG bfp</i>			
54-6000-001 (R-10-83)					

Ion Chromatographic Analysis of the Water Digestion

REAGENT BLANK									
Sample No.	Sample Point	Date	Time issued	Procedure	Charge Code	Charger ID	Sample Size	Method/Standard	Determination
F 5054-A-7572	SEGMENT-23	11-7-89	11:18	NDS	LA-525-105	MR75L	0	NDS	F 5046-A-7572 SEGMENT-15
F 5043-A-7074 SEGMENT-12									
Sample No.	Sample Point	Date	Time issued	Procedure	Charge Code	Charger ID	Sample Size	Determination	
E 5042-A-7574	SEGMENT-11	11-7-89	11:17	DBS	LA-525-105	WB75L	0	PQ-A	F 5042-A-7574 SEGMENT-11
F 5043-A-7074 SEGMENT-12									
Sample No.	Sample Point	Date	Time issued	Procedure	Charge Code	Charger ID	Sample Size	Determination	
E 5042-A-7574	SEGMENT-11	11-7-89	11:17	DBS	LA-525-105	WB75L	0	PQ-A	F 5042-A-7574 SEGMENT-11
4.77 ppm									
1430	12/5/89	hrs	hrs	hrs	hrs	hrs	hrs	hrs	hrs
168107	12/5/89	hrs	hrs	hrs	hrs	hrs	hrs	hrs	hrs
1430	12/5/89	hrs	hrs	hrs	hrs	hrs	hrs	hrs	hrs
168107	12/5/89	hrs	hrs	hrs	hrs	hrs	hrs	hrs	hrs

3 | 2 | 3 | 9 | 1 | 6 | 0 | 3 | 1 | 2 | 1 | 3 | 0 | 9 | 1 | 6

9 1 1 2 3 5 3 0 3 1 4

Ion Chromatographic Analysis of the Water Digestion

Serial No.	Sample Point	Date	Time Issued	Priority
F-5034.-7374	SEGMENT-23	11-7-89	11:18	18
Determination	Method/Standard	Result Units	Charge Code	Runno
P014	LA-503-105	FPM	WB75L	0
Sample Size				Customer ID
? Direct			89-041	
Remarks, Calculations, Results:				
REAGENT BLANK				
<i><1 ppm</i>				
Analyst-1 68107	Analyst-2	Analyst-3	Analyst-4 <i>79-3411</i>	Analyst-5 <i>79-3411</i>
Hrs 1430	Hrs	Hrs	Hrs	Hrs
Date 12/5/89	Time Completed	LAD UNH MGT <i>CGP</i>	<i>YR</i>	<i>off</i>

Serial No F 5045.-7274	Sample Point SEGMENT-14	Date 11- 7-89	Time Issued 11:17	Priority 19
Determination PO4	Method: Standard LA-333-105	Result Units % RECOVERY	Charge Code WB75L	Recurse 0
Sample Size ?				Customer ID 89-041
100mL = 10mL		PO4 - 501.69 ppm		
Remarks, Calculations, Results: SPike Sample SPike ID 4009-a SPike Volume .504L				
$\frac{(1.01) 765.5 - (477)(.91)}{(.050)(501.69)} \times 100 = 67.6\%$ $\frac{5.05}{(101)}$				
Analyst-1 10307 Hrs 14:30	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>[Signature]</i> Hrs
Date 12-5-89	Time Completed	Lab Unit Mgr <i>[Signature]</i>		

Serial No.	Sample Point	Date	Time Issued	Priority
F 5044-7174	SEGMENT-13	11-7-89	11:17	19
Determination	Method-Standard	Result Units	Charge Code	Review
PO4	LA-533-105	PPM	WE75L	0
Sample Size				Customer ID
? 100-10				89-041
Remarks, Calculations, Results:				
DUPLICATE SAMPLE				
3.55^2 ppm				
Analyst - 1 68107 HRS	Analyst - 2 HRS	Analyst - 3 HRS	Analyst - 4 HRS	Review Signature Date 440
1430	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>DP</i>	
Date 12/5/89				SI-8000-061 (R-10-0)

Serial No	Sample Point	Date	Time issued	Priority
F 5046,-7574	SEGMENT-15	11- 7-89	11:18	19
Determination	Method/Standard	Result/Units	Charge Code	Runns
P04	LA-533-105	% RECOVERY	WB75L	0
Sample Size			Customer ID	
100-10			89-041	
Remarks, Calculations, Results:				
LMCS CHECK SAMPLE LMCS ID <u>60G11H2</u> 2 + ^{10%} 2 HEP <u>592.2</u> 95.4				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5 Check
<u>60G107</u>	Hrs	Hrs	Hrs	<u>60G11711</u> Hrs
<u>1420</u>				
Date	Time Completed	Lab Unit Mgr		
<u>12/5/89</u>		<u>Cope</u>	<u>KP</u>	

9 1 1 2 2 5 9 0 8 1 5

Ion Chromatographic Analysis of the Water Digestion

Sample No F 5042-7575	Sample Point SEGMENT-11	Date 11-7-89	Time Issued 11:17	Priority 19
Determination S04	Method/Standard LA-535-105	Result Units % RECOVERY	Charge Code WB75L	Ref ID C
Sample Size 100-10				Customer ID 80-041
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <u>6C11AC</u>				
<p style="text-align: center;">575.9 / 621 92.1%</p>				
Analyst-1 <u>68107</u> Hrs <u>1430</u>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 Hrs <u>232121</u>
Date <u>12/5/89</u>	Time Completed	Lab Unit Mgr <u>CJG</u>		QD

Serial No.	Sample Point	Date	Time Issued	Priority
F 5054-7375	SEGMENT-23	11-7-89	11:18	16
Determination	Method/Standard	Result Units	Charge Code	Run No.
SU4	LA-533-105	PPM	WB75L	61
Sample Size				Customer ID
? Direct				89-041
Remarks, Calculations, Results:				
REAGENT BLANK				
<i><1 ppm</i>				
Analyst-1 <i>6B107</i>	Analyst - 2 <i>1430</i>	Analyst - 3 <i>1430</i>	Analyst - 4 <i>1430</i>	Analyst - 5 <i>1430</i>
Hrs <i>1430</i>	Hrs <i>1430</i>	Hrs <i>1430</i>	Hrs <i>1430</i>	Hrs <i>1430</i>
Date <i>12/5/89</i>	Time Completed	Labs Unit Mgr <i>CJA</i>	Q/C <i>OK</i>	

Serial No.	Sample Point	Date	Time Issued	Priority
F 5043-7075	SEGMENT-12	11-7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return
SO4	LA-533-105	PPM	WB75L	0
Sample Size				Customer ID
100-R				89-041
Remarks, Calculations, Results				
5.37 ppm				
Analyst -1	Analyst -2	Analyst -3	Analyst -4	Review by
46107				9421
Hrs	Hrs	Hrs	Hrs	Hrs
1420				
Date	Time Computed		Lab Unit Mgr	<i>Cga.</i>
12/5/87				<i>ff</i>

9 1 1 2 7 5 9 0 8 1 6

Ion Chromatographic Analysis of the Water Digestion

Serial No. F 5045.-7275	Sample Point SEGMENT-14	Date 11- 7-89	Time issued 11:17	Priority 19
Determination SO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB7SL	Reruns 0
Sample Size ?	Customer ID 100 uL - 10 mL 89-041			
Remarks, Calculations, Results: SPIKE SAMPLE SO ₄ - 501.03 SPIKE ID 40CG-19 SPIKE VOLUME 50uL $\frac{(1.0)(403.9) - (53.7)(.7)}{(1.050)(501.03)} \times 100 = 89.8\%$ $\frac{1.050}{5.05} (101)$				
Analyst-1 6B107 Hrs 14:30	Analyst-2	Analyst-3	Analyst-4	Analyst-5
Date 12-5-89	Time Completed	Lab Unit No. C-10000-001 (A-10-00)	<i>[Signature]</i>	

Serial No. F 5046.-7375	Sample Point SEGMENT-15	Date 11- 7-89	Time issued 11:18	Priority 19
Determination SO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB7SL	Reruns 0
Sample Size 100-10	Customer ID 89-041			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 6C11HC				
Analyst-1 6B107 Hrs 14:30	Analyst-2	Analyst-3	Analyst-4	Analyst-5
Date 12-5-89	Time Completed	Lab Unit No. C-10000-001 (A-10-00)	<i>[Signature]</i>	H1

5959 96.07

Total Organic Carbon Analysis of the Water Digestion I - Not Acidified

Serial No.	Sample Point		Date	Time issued	Priority
F 5042.-7526	SEGMENT-11		11- 7-89	11:16	19
Determination	Method/Standard	Result Units	Charge Code	Return	
TOC	LA-344-105	% RECOVERY	WB75L	0	
Sample Size				Customer ID	
? 200-2-200				89-041	
Remarks, Calculations, Results:					
LMCS CHECK SAMPLE LMCS ID 70C11B					
2.96 g/L $\cancel{3.03}$ 48.80					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
80725 L.C.				R/Element	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr			
12-4-89		Cyg			
64-8000-061 (R-10-82)					

Serial No.	Sample Point		Date	Time issued	Priority
F 5044.-7126	SEGMENT-13		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
TOC	LA-344-105	G/L	WB75L	0	
Sample Size				Customer ID	
? 200-2-200				89-041	
Remarks, Calculations, Results:					
DUPLICATE SAMPLE					
4.50 g/L 4.51 g/L					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
80725 L.C.				R/Element	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr			
12-4-89		Cyg			
64-8000-061 (R-10-82)					

Serial No.	Sample Point		Date	Time issued	Priority
F 5043.-7026	SEGMENT-12		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
TOC	LA-344-105	G/L	WB75L	0	
Sample Size				Customer ID	
? 200-2-200				89-041	
Remarks, Calculations, Results:					
4.50 g/L 4.51 g/L					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
80725 L.C.				R/Element	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr			
12-4-89		Cyg			
64-8000-061 (R-10-82)					

Serial No.	Sample Point		Date	Time issued	Priority
F 5045.-7226	SEGMENT-14		11- 7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Return	
TOC	LA-344-105	% RECOVERY	WB75L	0	
Sample Size				Customer ID	
? 200-2-200 + 200 spike				89-041	
Remarks, Calculations, Results:					
$(54.7 \text{ ppm} - 4.1 \text{ ppm}) / (4.1 \text{ ppm} \cdot 4.18 \text{ mL}) = 95.2\%$ $\text{SPIKE SAMPLE } 95.2\% \text{ K2HPO}_4$ $\text{SPIKE ID } 70C11B$ $\text{SPIKE VOLUME } 200 \text{ mL}$					
2.69 g/L $2.69 \text{ g/L} \cdot 95.2\% = 86.69 \text{ g/L}$ $86.69 \text{ g/L} \cdot 4.18 \text{ mL} = 360 \text{ mL}$ $360 \text{ mL} / 2.69 \text{ g/L} = 134 \text{ g/L}$					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
80725 L.C.				R/Element	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr			
12-4-89		Cyg			
64-8000-061 (R-10-82)					

9 1 1 2 7 5 9 0 8 1 0

Total Organic Carbon Analysis of the Water Digestion I - Not Acidified

Serial No.	Sample Point		Date	Time issued	Priority
F 5046.-7526	SEGMENT-15		11-7-89	11:17	19
Determination	Method/Standard	Result Units	Charge Code	Refine	
TOC	LA-344-105	% RECOVERY	WB75L	0	
Sample Size				Customer ID	
? 89-041					
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 70C11B					
$\frac{3.07 \text{ g/l}}{3.00} = 102\%$					
Analyst-1 80725 J.C. Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 R/Reanalyzed Hrs	
Date 12-4-89	Time Completed	Lab Unit No.	<i>Chas</i>	<i>100</i>	
S4-9800-081 (R-10-82)					

Serial No.	Sample Point		Date	Time issued	Priority
F 5054.-7326	SEGMENT-23		11-7-89	11:18	18
Determination	Method/Standard	Result Units	Charge Code	Refine	
TOC	LA-344-105	G/L	WB75L	0	
Sample Size				Customer ID	
? 200-2-200				89-041	
Remarks, Calculations, Results: REAGENT BLANK					
$\frac{2.25 \text{ g/l}}{2.00} = 112\%$					
Analyst-1 80725 J.C. Hrs	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 R/Reanalyzed Hrs	
Date 12-4-89	Time Completed	Lab Unit No.	<i>Chas</i>	<i>100</i>	
S4-9800-081 (R-10-82)					

Total Organic Carbon Analysis of the Water Digestion II- Acidified

9 1 1 2 3 5 7 0 9 1 9

Serial No. F 5043.-7026	Sample Point SEGMENT-12		Date 11- 7-89	Time issued 11:17	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units G/L	Charge Code WB75L	Refers 1	
Sample Size ? 1ml + 100ul 5m H ₂ SO ₄ - 200ul			Customer ID 89041		
Remarks, Calculations, Results: Second leach 8.25E -3 g/l 8.25E -3 g/l					
Analyst-1 80028 Hrs Ed Loh	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	Analyst-6 Hrs
Date 6-29-90	Time Completed	Lab Unit Mgr Zell M. Rohr Kelly Homelbeck 64-6800-081 (R-10-83)			

Serial No. F 5042.-7526	Sample Point SEGMENT-11		Date 11- 7-89	Time issued 11:16	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Refers 1	
Sample Size ? 200ul - 2ml - 200ul			Customer ID 89041		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 20611D					
Analyst-1 80028 Hrs Ed Loh	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	Analyst-6 Hrs
Date 6-29-90	Time Completed	Lab Unit Mgr Zell M. Rohr Kelly Homelbeck 64-6800-081 (R-10-83)			

Serial No. F 5045.-7226	Sample Point SEGMENT-14		Date 11- 7-89	Time issued 11:17	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Refers 1	
Sample Size ? 200ul + 100ul 5m H ₂ SO ₄ - 200ul			Customer ID 89041		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 20611A SPIKE VOLUME 200ul Second leach $(126 - 3.5) - (50.35) \times 100 = 100.9\%$					
Analyst-1 80028 Hrs Ed Loh	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	Analyst-6 Hrs
Date 6-29-90	Time Completed	Lab Unit Mgr Zell M. Rohr Kelly Homelbeck 64-6800-081 (R-10-83)			

Serial No. F 5044.-7126	Sample Point SEGMENT-13		Date 11- 7-89	Time issued 11:17	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units G/L	Charge Code WB75L	Refers 1	
Sample Size ? 1ml + 100ul 5m H ₂ SO ₄ - 200ul			Customer ID 89041		
Remarks, Calculations, Results: DUPLICATE SAMPLE Second leach $1.32 \times 10^{-2} g/l$					
Analyst-1 80028 Hrs Ed Loh	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	Analyst-6 Hrs
Date 6-29-90	Time Completed	Lab Unit Mgr Zell M. Rohr Kelly Homelbeck 64-6800-081 (R-10-83)			

9 1 1 2 2 5 2 0 3 2 0

Total Organic Carbon Analysis of the Water Digestion II - Acidified

Serial No. F 5054.-7326	Sample Point SEGMENT-23	Date 11- 7-89	Time issued 11:18	Priority 18
Determination TOC	Method/Standard LA-344-105	Result Units G/L	Charge Code WB75L	Run no. 1
Sample Size ? 200ul		Customer ID 89041		
Remarks, Calculations, Results: REAGENT BLANK				
<i>Second leach</i>				
3.5ug 4995ug/min				
Analyst-1 80028 <i>Ed Cohn</i>	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 R Bennett
Date 6-29-90	Time Completed	Lab Unit/Mgr. 7000-001 (R-10-83)	<i>Kathy Hendrickson</i>	

Serial No. F 5046.-7526	Sample Point SEGMENT-15	Date 11- 7-89	Time issued 11:17	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Run no. 1
Sample Size ? 200ul - 2ml - 200ul		Customer ID 89041		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 700110				
<i>2nd leach</i>				
2.920 gpl / 3.000sl 97.30%				
Analyst-1 80028 <i>Ed Cohn</i>	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 R Bennett
Date 6-29-90	Time Completed	Lab Unit/Mgr. 7000-001 (R-10-83)	<i>Kathy Hendrickson</i>	

2 1 1 2 0 5 9 0 9 ?

Acid Digestion

Serial No F-5047-8500	Sample Point SEGMENT-P	Date 11-7-89	Time Issued 11:18	Priority 26
Determination ACD-DGST	Method/Standard LA-505-159	Result Units G/L	Charge Code E21D1	Reruns 0
Sample Size 10 ml	Customer ID 89-041			
Remarks, Calculations, Results: GRAMS SAMPLE .5432 VOLUME ON COMPLETION 20ml 25 1.09 g/ml				
Analyst-1 JASmit 65286	Analyst-2 D6 Hartung 66481	Analyst-3 69769/EMR	Analyst-4 KJYML	Analyst-5 KJYML
Hrs 11/6/89	Hrs 11/6/89	Hrs 11/6/89	Hrs 11/6/89	Hrs 11/6/89
Date 12/6/89	Time Completed 11:18	Lab Unit Mgr Corynne	Signature	
54-6800-061 (R-10-83)				

Serial No F-5050-8500	Sample Point SEGMENT-19	Date 11-7-89	Time Issued 11:18	Priority 26
Determination ACD-DGST	Method/Standard LA-505-159	Result Units G/L	Charge Code W.E.751	Reruns 0
Sample Size 10 ml	Customer ID 89-041			
Remarks, Calculations, Results: GRAMS SAMPLE .5838 VOLUME ON COMPLETION 20ml 103C15C + 104C15D 1.17 g/ml				
Analyst-1 69769/EMR	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12/6/89	Time Completed 11:18	Lab Unit Mgr Corynne	Signature	
54-6800-061 (R-10-83)				

Serial No F-5047-8500	Sample Point SEGMENT-P	Date 11-7-89	Time Issued 11:18	Priority 26
Determination ACD-DGST	Method/Standard LA-505-159	Result Units G/L	Charge Code E21D1	Reruns 0
Sample Size 5 ml	Customer ID 89-041			
Remarks, Calculations, Results: GRAMS SAMPLE .5922 VOLUME ON COMPLETION 20ml 25 1.18 g/ml				
Analyst-1 69769/EMR	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs
Date 12/6/89	Time Completed 11:18	Lab Unit Mgr Corynne	Signature	
54-6800-061 (R-10-83)				

9 1 1 2 3 5 7 0 8 ? 2

TCP Analysis

Serial No.	Sample Point		Date	Time Issued	Priority
F 5049.-8050	SEGMENT-17		11- 7-89	11:18	23
Determination	Method/Standard	Result Units	Charge Code	Remarks	
ICP	LA-505-151	PPM	WB75L	0	
Sample Size			Customer ID		
? 100-10	1-10		89-041		
Remarks, Calculations, Results:					
<i>Completed</i>					
Analyst-1 69769	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>R. Southchurch</i>					
Date 2/22/90	Time Completed	Lab Unit Mgr <i>QD</i>	Comments 84-000-081 (R-10-83)		

Serial No.	Sample Point		Date	Time Issued	Priority
F 5047.-8150	SEGMENT-1c		11- 7-89	11:18	23
Determination	Method/Standard	Result Units	Charge Code	Remarks	
ICP	LA-505-151	% RECOVERY	WB75L	0	
Sample Size			Customer ID		
? 1-10			89-041		
Remarks, Calculations, Results:					
<i>LMCS CHECK SAMPLE Check Standards</i> <i>LMCS ID E-111 (1,18 (2) E-111A 28</i> <i>34C11C = digested LMCS Check</i>					
Analyst-1 69769	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>R. Southchurch</i>					
Date 2/22/90	Time Completed	Lab Unit Mgr <i>QD</i>	Comments 84-000-081 (R-10-83)		

Serial No.	Sample Point		Date	Time Issued	Priority
F 5050.-8250	SEGMENT-19		11- 7-89	11:18	23
Determination	Method/Standard	Result Units	Charge Code	Remarks	
ICP	LA-505-151	PPM	WB75L	0	
Sample Size			Customer ID		
? 100-10	1-10		89-041		
Remarks, Calculations, Results:					
<i>103C15C = 104C15D</i> <i>5.0 ml each</i> <i>SPIKE</i> <i>Completion 50.0ml</i> <i>1.17⁻² g/ml Sample</i> <i>89-041</i>					
Analyst-1 69769	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>R. Southchurch</i>					
Date 2/22/90	Time Completed	Lab Unit Mgr <i>QD</i>	Comments 84-000-081 (R-10-83)		

Serial No.	Sample Point		Date	Time Issued	Priority
F 5049.-8150	SEGMENT-18		11- 7-89	11:18	23
Determination	Method/Standard	Result Units	Charge Code	Remarks	
ICP	LA-505-151	PPM	WB75L	0	
Sample Size			Customer ID		
? 100-10	1-10		89-041		
Remarks, Calculations, Results:					
<i>DUPLICATE SAMPLE</i> <i>Completed</i>					
Analyst-1 69769	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
<i>R. Southchurch</i>					
Date 2/22/90	Time Completed	Lab Unit Mgr <i>QD</i>	Comments 84-000-081 (R-10-83)		

9 1 1 2 7 5 9 0 3 ? 3

Acid Digestion

Serial No.	Sample Point		Date	Time Issued	Priority
F 5051 -8500	SEGMENT-T		11- 7-89	11:18	26
Determination	Method/Standard	Result Units	Charge Code	Reruns	
ACD-DGST	LA-505-159	G/L	E21D1	0	
Sample Size				Customer ID	
? .5ml				ICP 2	
Remarks, Calculations, Results: GRAMS SAMPLE _____ VOLUME ON _____ COMPLETION <u>sample</u> 35 C11CJ					
Analyst-1 69769 /SMG	Analyst-2	Analyst-3	Analyst-4	Analyst-5	Analyst-6
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 10/16/89	Time Completed	Lab Unit Mgr <i>John C. Smith</i>	S-0000-001 (R-10-63)		