



Westinghouse  
Hanford Company

P.O. Box 1970 Richland, WA 99352

012193 857  
4 of 10

## 222-S/RCRA ANALYTICAL LABORATORIES

PROJECT: SINGLE-SHELL TANK WASTE  
CHARACTERIZATION

TANK: 241-U-110

CORE: 6

SEGMENT: 4

CUSTOMER ID. NUMBER:  
89-045



REPORT REVISION: 1

DATE PRINTED: AUGUST 24, 1990

## Table of Contents

Signature Form.....	2
Introduction.....	3
Sampling and Custody Data.....	6
Sample Data Summary.....	9
Physical Test Results.....	11
Extrusion of Segment.....	12
Appearance of Solid.....	12
Homogenization.....	12
Single Shell Tank Segment - Subsamples.....	13
VOA Sampling.....	13
Photograph of Tank 241-U-11 Core 6 Segment 4.....	14
Particle Size Analysis I .....	15
Particle Size Analysis II.....	20
Undigested Sample Analysis.....	26
Untreated Sample Results Summary.....	27
pH Analysis.....	28
Percent Water.....	29
KOH Fusion Analysis.....	30
Fusion Analysis Results Summary .. .	31
Fusion Analysis of Laboratory Solution. ....	32
Fusion Dissolution.....	35
Total Alpha and Beta Analysis I .....	36
Total Alpha and Beta Analysis II.....	37
GEA Analysis.....	44
Uranium Analysis I .....	84
Uranium Analysis II.....	85
Water Digestion Test Analysis.....	86
Water Digestion Results Summary.....	87
Water Digestion Analysis of Laboratory Solutions.....	88
Water Digestion .....	89
Ion Chromatographic Analysis I .....	90
Ion Chromatographic Analysis II.....	91
Total Organic Carbon Analysis..(Acidified).....	123

Acid Digestion Test Analysis.....	130
ICP Data Summary.....	131
Acid Digestion Analysis.....	132
Inductively Coupled Plasma Spectrographic Analysis.....	133
ICP Raw Data Summary.....	134

Appendix A

Analytical Analysis Cards

9  
8  
7  
6  
5  
4  
3  
2  
1

I have reviewed this report and certify that the package is in compliance with Quality Assurance Project Plan - WHC-SD-CP-QAPP-002. I found it to be a true and accurate accounting both technically and for completeness of the laboratory analyses performed on this sample.

Shirley Cervantes  
Shirley A. Cervantes  
Data Coordinator

Date August 24, 1990

Cary M Seidel  
Cary M. Seidel  
Unit Manager

Date Aug 24, 1990

L.H. Taylor  
Larry H. Taylor  
Laboratory Q.A. Officer

Date Aug 30, 1990

# INTRODUCTION

## INTRODUCTION

Westinghouse Hanford Company 222-S/RCRA Analytical Laboratories are supporting the characterization efforts of the single shell tanks. The characterization of tank 241-U-110 was performed under Phase 1A and 1B of the Waste Characterization Plan for the Hanford Site Single-Shelled Tanks (WHC-EP-0210).

Tank 241-U-110 has a 500,000 gallon capacity, construction was completed in 1944. The tank received first cycle waste, REDOX high-level waste, coating waste, and laboratory waste until 1975. Between July 7, 1975 and February 2, 1976, P-10 pumps were installed, and 41,700 gallons of liquid waste were pumped from the tank. Tank 241-U-110 still contains an estimated 195,000 gallons of waste.

Analytical Laboratories performs all analytical analysis to the specifications of the Quality Assurance Project Plan, WHC-SD-CP-QAPP-002. In accordance with WHC-SD-CP-QAPP-002 the following laboratory policies are being followed.

Spikes are performed on either the undissolved sample, or the sample after dissolution as directed by the chemist. If the spike addition is found to be less than 20% of an analyte concentration, the spike recovery is not reported due to errors introduced by the precision of the sample analysis. The concentration of spike additions will be re-evaluated before the start of phase 1C. Two spiking routines are being used during phase 1A and 1B. For the following analyses, Ion Chromatography, Inductively Coupled Plasma, Mercury Hydride, Total Organic Carbon, and Carbonate analyses the solid sample is spiked independently from the sample digestion. Any non-homogeneity of the sample could adversely affect the spike recoveries. For the radio-isotopic analysis and other analyses not specified above the spikes were preformed by spiking an aliquot of sample after digestion.

The laboratory does not report sample results from batch analyses that are questionable. The results from questionable batches are discarded and the analysis is repeated. Sample cards (laboratory travelers) for the repeated analysis are reissued for analysis after they have been stamped "rerun". Laboratory travelers are issued using a computerized routine according to a "sample point". This sample point label (segment-n) on the Laboratory travelers and on the GEA analysis reports has no relationship to the sampling activities or the sample identification. All results in this data package relate only to the sample identified as segment 4 from core 6 taken from tank 241-U-110.

The organic analysis of this sample will be performed by Pacific Northwest Laboratories (PNL). Due to instrument and procedure problems, PNL has been unable to separate organics from the normal paraffin hydrocarbon present in the samples. The results from the organic analysis will be provided when available.

Samples analyzed for Total Organic Carbon between November 1, 1989 and February 22, 1990 were not acidified. The results from these analyses include total organic carbon, carbonate, and dissolved carbon dioxide from the air. The validity of these analysis are subject to interpretation. The total organic carbon procedure was corrected and these analyses will be repeated wherever possible.

All sample results reported here by weight are reported as the "wet weight" of the sample. Some samples did noticeably lose moisture during the process of aliquoting and weighing the sample for digestion. In order to minimize errors due to loss of moisture, the percent moisture was determined at the earliest opportunity. Attempts to dry the sample before analysis resulted in approximately a ten fold increase in radiation levels. In order to reduce and control radiation exposure to laboratory personnel the samples were not dried before aliquouting and digestion. This may result in some laboratory results being biased high.

This report is formatted into sections corresponding to the type of dissolutions performed prior to analysis. A brief summary of analytical results is reported, followed by calibration data and an analysis batch report. Any notable observations regarding an analysis are noted on the batch report for that analysis. Copies of laboratory travelers can be found in Appendix A.

9 1 1 2 2 3 0 1 0 7 4

## SAMPLING AND CUSTODY DATA

Lab Segment Serial #  
F0101

9 1 1 2 2 5 0 1 0 7 5

CHAIN-OF-CUSTODY RECORD FOR CORE SAMPLING

(1) Shipment Number S-024-89 (2) Sample Number 89-045 (3) Supervisor D.C. Hartley  
(4) Tank 104 (5) Riser 17 (6) Segment #4 (7) Cask Serial Number 1003C

Radiation Survey Data:	(8) FIELD	(20) LABORATORY	(9) Shipment Description:
Over Top Dose Rate	<u>1.5 Mr/hr</u>	<u>1.5 MR/hr</u>	A. Work Package Number <u>2W-89-0095RW</u>
Side Dose Rate	<u>2 mR/hr</u>	<u>2 mR/hr</u>	B. Cask Seal Number <u>For Future Use</u>
Bottom Dose Rate	<u>1 mR/hr</u>	<u>1 mR/hr</u>	C. Sampler Number Used <u>#49</u>
Smearable Contamination	<u>1 Det.</u> (alpha)	<u>1 Det</u> (alpha)	D. Date and Time Sampler Unseated <u>11-14-89, 1730</u>
	<u>&lt; DET.</u> (beta-gamma)	<u>&lt; Det</u> (beta-gamma)	E. Expected Liquid Content <u>20%</u>
	RPT <u>R. Sain</u> (Signature)	RPT <u>D. Arnold</u> (Signature)	F. Expected Solid Content <u>80%</u>
			G. Dose Rate Through Drill String <u>70 mR/hr.</u>
			H. Expected Sample Length <u>19"</u>

(10) INFORMATION (Include statement of laboratory tests to be performed.\*)

Core #006, WHE-EP-0210, Waste characterization Plan for the  
Hanford Site Single Shell Tank

\*Reference laboratory work request, if available.

Comments:

(11) POINT OF ORIGIN <u>241-4</u> <u>116</u>	(12) SENDER NAME <u>D.C. Hartley</u> SENDER SIGNATURE <u>DCHartley</u>	(13) DATE AND TIME RELEASED <u>11-15-89</u> <u>1025</u>	(14) DESTINATION <u>222 S</u> <u>LABS.</u> <u>200 West</u>	(16) RECIPIENT NAME <u>C. M. Seidel</u> RECIPIENT SIGNATURE <u>Craig M Seidel</u>	(17) DATE AND TIME RECEIVED <u>1045</u> <u>11-15-89</u>
(15) Seal Intact Upon Release? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	(18) Seal Intact Upon Receipt? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	(19) Seal Data Consistent with this Record? Shipment No. <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No      Sample No. <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No			

**Single Shell Tank Waste Characterization  
Summary of Core Sample**

**Phase  
I-A**

Tank ID:	241-U-110
Riser ID:	17
Core ID:	006

Date Sampling Initiated:	11/10/90
Date Sampling Completed:	11/14/90

Segment 1	Lab Serial No.	F0029	Segment 8	Lab Serial No.
	Customer ID. No.	89-042		Customer ID. No.
	Last Segment?	No		Last Segment?
Segment 2	Lab Serial No.	F0053	Segment 9	Lab Serial No.
	Customer ID. No.	89-043		Customer ID. No.
	Last Segment?	No		Last Segment?
Segment 3	Lab Serial No.	F0077	Segment 10	Lab Serial No.
	Customer ID. No.	89-044		Customer ID. No..
	Last Segment?	No		Last Segment?
Segment 4	Lab Serial No.	F0101	Segment 11	Lab Serial No.
	Customer ID. No.	89-045		Customer ID. No.
	Last Segment?	Yes		Last Segment?
Segment 5	Lab Serial No.		Segment 12	Lab Serial No.
	Customer ID. No.			Customer ID. No.
	Last Segment?			Last Segment?
Segment 6	Lab Serial No.		Segment 13	Lab Serial No.
	Customer ID. No.			Customer ID. No.
	Last Segment?			Last Segment?
Segment 7	Lab Serial No.		Segment 14	Lab Serial No.
	Customer ID. No.			Customer ID. No.
	Last Segment?			Last Segment?

REMARKS: CUSTOMER ID# 89-042  
WAS RECEIVED EMPTY.

Interim

Prepared by: Shirley Cervantes S. A. Cervantes Date: 07/07/90  
Signature Printed Name

Verified by: Cary M Seidel C. M. Seidel Date: 07/07/90  
Signature Printed Name

Approved by: J.H. Taylor L.H. Taylor Date: 08-30-90  
Signature Printed Name

# SAMPLE DATA SUMMARY

9 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 0 | 7 | 7

## SUMMARY DATA REPORT

Tank	241-U-110			Acid Digestion		
Core	6			Wet Weight Sample	Wet Weight Sample	Duplicate
Segment	4					
Customer ID	89-045					
<b>Untreated Sample</b>						
	Sample	Duplicate		Aluminum	45520 ug/g	52568 ug/g
pH	11.96	12.22		Antimony	LT	LT
Percent Water	38.50 %	37.30 %		Arsenic	LT	LT
<b>Fusion Analysis</b>						
Total Alpha	<0.517 uci/g	<0.873	uci/g	Barium	LT	LT
Total Beta	84.6 uci/g	86.0	uci/g	Beryllium	LT	LT
GEA				Bismuth	21590 ug/g	26620 ug/g
Cs-137	18.7 uci/g	23.4	uci/g	Boron	LT	LT
Uranium	1390 ug/g	1310	ug/g	Cadmium	LT	LT
	1080 ug/g			Calcium	295 ug/g	1113 ug/g
				Cerium	LT	LT
				Chromium	119 ug/g	188 ug/g
				Copper	LT	LT
				Iron	14837 ug/g	16885 ug/g
				Lanthanum	LT	LT
				Lead	LT	602 ug/g
				Magnesium	479 ug/g	449 ug/g
				Manganese	156 ug/g	269 ug/g
				Mercury	LT	LT
				Molybdenum	LT	LT
				Nickel	LT	LT
<b>Water Digestion</b>						
Fluoride	17200 ug/g	18600	ug/g	Phosphorous	53626 ug/g	47284 ug/g
Chloride	<953 ug/g	<979	ug/g	Potassium	LT	LT
Nitrate	25500 ug/g	26600	ug/g	Samarium	LT	LT
Phosphate	146000 ug/g	159000	ug/g	Selenium	LT	LT
Sulfate	4060 ug/g	<9790	ug/g	Silver	LT	LT
Total Organic Carbon *	726 ug/g	693	ug/g	Sodium	187906 ug/g	168975 ug/g
				Strontium	120 ug/g	106 ug/g
				Sulfur	113 ug/g	278 ug/g
				Tantalum	LT	LT
				Thallium	LT	LT
				Thorium	LT	LT
				Tin	LT	LT
				Titanium	LT	LT
				Uranium	LT	LT
				Vanadium	LT	LT
				Zinc	142 ug/g	5380 ug/g
				Zirconium	LT	97 ug/g

LT      Less Than Detection Limit

revised August 24, 1990

6  
4  
C  
I  
C  
S  
C  
2  
—  
1  
9

## PHYSICAL TEST RESULTS

**Single Shell Tank  
Extrusion of Segment -- Physical Tests**

**Phase  
I-A**

Lab Segment Serial No.: F-0101

Customer ID: 89-045

Analyst: RICHARD L. WEISS

Date Extruded: 11-16-89

**Drainable Liquid**

Liquid Submitted for Segment Analysis? -- No

Gross 25mL	Tare	Net
Serial	Date/Time _____ / _____	Estimated
Specific	Calculated	

Appearance of Liquid:

**Dimensions of Segment**

Complete Segment Obtained? NO	Length: 6.0 in	Calculated Volume: 4.71 cubic in
-------------------------------	----------------	----------------------------------

Remarks

Appearance of Solid:

SAMPLE DARK BROWN GRADING UPWARD TO MEDIUM BROWN. GRANULAR LOOKING THROUGHOUT SAMPLE. MOIST AND STICKY APPEARANCE AT THE BOTTOM; GRADING UPWARDS TO CRUMBLY AND NON-COHESIVE AT THE TOP.

Penetrometer

10

lbs/sq in

Remarks:

**Homogenization**

Procedure: T038A-00712 Revision: F	Quantity of Material 138.63 grams
Date Homogenized: 12-29-89	Time Homogenized: 5 Minutes
Operator: JOHN R. SMITH (65286)	

**Laboratory Notebook Reference**

WHC-N-313-4

Notebook No.

13

Page No.

Prepared by: <i>Herlene S. Rich</i> Signature	HERLENE S. RICH Printed Name	Date: 6-01-90
Verified by: <i>Cary M. Seidel</i> Signature	C.M. SEIDEL Printed Name	Date: 6-01-90
Approved by: <i>L.H. Taylor</i> Signature	L.H. Taylor Printed Name	Date: 08-30-90

Interim

3/27/90

Rev.C

SST-3

12

**Single Shell Tank  
Segment -- Subsamples**

**Phase  
1A**

Customer ID: 89-045

Lab Segment Serial No. F0101

**Volatile Organic Analysis**

VOA Sample

Laboratory Serial Number: 89-045-36

Date Sampled: 11/15/89

**Particle Size Distribution Analysis**

Particle Size Sample

Laboratory Serial Number: F0101

Date Sampled: 11/16/89

**Homogenized Solids**

**Undigested Solids Analysis**

Laboratory Serial Number for Sample: F0101

Date Sampled: 11/16/89

Laboratory Serial Number of Duplicate Sample: F0102

**Fusion Analysis of Solids**

Laboratory Serial Number for Sample: F0107

Date Sampled: 11/16/89

Laboratory Serial Number of Duplicate Sample: F0108

Laboratory Serial Number of Spiked Sample:

**Acid Digestion Analysis of Solids**

Laboratory Serial Number for Sample: F0116

Date Sampled: 11/16/89

Laboratory Serial Number of Duplicate Sample: F0117

Laboratory Serial Number of Spiked Sample: F0118

**Water Digestion Analysis of Solids**

Laboratory Serial Number for Sample: F0111

Date Sampled: 11/16/89

Laboratory Serial Number of Duplicate Sample: F0112

Laboratory Serial Number of Spiked Sample: F0113

**Laboratory Notebook Reference**

WHC-N-313-4  
Notebook No.

13  
Page No.

04/24/90  
Rev.A  
SST-17


Prepared by: Shirley Cervantes  
Signature S. A. Cervantes Date: 06-06-90

Verified by: Cary M. Seidel  
Signature C. M. Seidel Date: 06-06-90

Approved by: L.H. Taylor  
Signature L.H. Taylor Date: 08-30-90

9 1 1 2 3 5 2 1 0 3 2

14



TANK 241-II-110 CORE 6 SEGMENT 4

## Brinkmann

## Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010  
STATISTICS

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.001

DATE	:	04/12/1989	ACQ. RANGE	:	0.5-150	COUNTS	:	379228
TIME	:	08:14	ACQ. MODE	:	SAMPLE	S.N.F.	:	0.72
CONFIG.	:	1 (0.7 S1)	ACQ. TIME	:	1873 SEC	S.D.U.	:	4694
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	:	5	CONCENTR.	:	6.5E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	:	95.00%(V)	SOLIDS	:	2.9E-03 %

	MEAN Diameter	S.D.
Number, Length	1.17 $\mu$ m	0.89 $\mu$ m
Number, Area	1.47 $\mu$ m	0.94 $\mu$ m
Number, Volume	2.04 $\mu$ m	1.25 $\mu$ m
Length, Area	1.85 $\mu$ m	1.95 $\mu$ m
Length, Volume	2.69 $\mu$ m	2.13 $\mu$ m
Area, Volume	3.91 $\mu$ m	5.23 $\mu$ m
Volume, Moment	10.90 $\mu$ m	14.11 $\mu$ m

	MEDIAN Diameter	MODE	CONFIDENCE
Number	0.90 $\mu$ m	0.75 $\mu$ m	100.00%
Area	2.65 $\mu$ m	4.75 $\mu$ m	100.00%
Volume	5.47 $\mu$ m	4.75 $\mu$ m	99.92%

Sample red-brown, mushy, with hard lumps  
 Dispersed well in  $H_2O$ , no agglomeration  
 Dispersed particles < 150  $\mu$ m

## Brinkmann

## Particle Size Analyzer

## PROCESS CHEMISTRY LABS PARTICLE ANALYSIS

VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK


FILE NAME : F0101.001

DATE	:	04/12/1989	ACQ. RANGE	:	0.5-150	COUNTS	:	379228
TIME	:	08:14	ACQ. MODE	:	SAMPLE	S.N.F.	:	0.72
CONFIG.	:	1 (0.7 s1)	ACQ. TIME	:	1873 SEC	S.D.U.	:	4694
CELL TYPE	:	MAGNETIC (3)	SAMPLE SIZE	:	5	CONCENTR.	:	6.5E+06 #/ml
SAMPLE TYPE	:	REGULAR	REQ. CONF.	:	95.00%(v)	SOLIDS	:	2.9E-03 %

## PROBABILITY VOLUME DENSITY GRAPH

Name: SST,B000029,F0101,H2O.SBK

2.9E-05 cc/ml(100.0%)

Mode at 5.00  $\mu\text{m}$ ▼ << SCALE RANGE ( $\mu\text{m}$ ): ADJUSTED >>Median : 5.47  $\mu\text{m}$ Mean( $\text{vv}$ ): 18.90  $\mu\text{m}$ S.D.( $\text{vv}$ ): 14.11  $\mu\text{m}$ Conf( $\text{vv}$ ): 99.92 %

## Brinkmann

## Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.001

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-150	COUNTS	: 379228
TIME	: 08:14	ACQ. MODE	: SAMPLE	S.N.F.	: 0.72
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 1873 SEC	S.D.U.	: 4694
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 5	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(v)	SOLIDS	: 2.9E-03 %


PROBABILITY VOLUME DISTRIBUTION GRAPH

Name: SST,B000029,F0101,H2O.SBK

2.9E-05 cc/ml(100.0%)

Mean(nv): 2.04 $\mu$ mMedian : 5.47 $\mu$ mS.D.(nv): 1.25 $\mu$ mMean(vm): 10.90 $\mu$ mS.D.(vm): 14.11 $\mu$ mS.D.(vm): 14.11 $\mu$ m

Conf(vm): 99.92 %

L9 << SCALE RANGE ( $\mu$ m): ADJUSTED >>

## Brinkmann

## Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.001

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-150	COUNTS	: 379228
TIME	: 08:14	ACQ. MODE	: SAMPLE	S.N.F.	: 0.72
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 1873 SEC	S.D.U.	: 4694
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 5	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 2.9E-03 %


## PROBABILITY VOLUME DENSITY GRAPH

Name: SST,B000029,F0101,H2O.SBK

2.9E-05 cc/ml(100.0%)

Mode at 5.00  $\mu\text{m}$ Mean(nv): 2.04  $\mu\text{m}$ S.D.(nv): 1.25  $\mu\text{m}$ Median : 5.47  $\mu\text{m}$ Mean(vm): 10.98  $\mu\text{m}$ S.D.(vm): 14.11  $\mu\text{m}$ 

Conf(vm): 99.92 %

<< SCALE RANGE ( $\mu\text{m}$ ): ADJUSTED >>

Brinkmann

Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.001

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-150	COUNTS	: 379220
TIME	: 08:14	ACQ. MODE	: SAMPLE	S.N.F.	: 0.72
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 1873 SEC	S.D.U.	: 4694
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 5	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(v)	SOLIDS	: 2.9E-03 %

PROBABILITY VOLUME DISTRIBUTION GRAPH

Name: SST,B000029,F0101,H2O.SBK

2.9E-05 cc/ml(100.0%)

Mean(nv): 2.04 $\mu$ m

Median : 5.47 $\mu$ m

S.D.(nv): 1.25 $\mu$ m


Mean(vm): 10.90 $\mu$ m

S.D.(vm): 14.11 $\mu$ m

S.D.(vm): 14.11 $\mu$ m

Conf(vm): 99.92 %

<< SCALE RANGE ( $\mu$ m): ADJUSTED >>



PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010  
STATISTICS

SAMPLE NAME : SST.B000029,F0101,H2O.SBK

FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.:	6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %

	MEAN Diameter	S.D.
--	---------------	------

Number, Length	: 0.96 $\mu$ m	0.75 $\mu$ m
Number, Area	: 1.22 $\mu$ m	0.80 $\mu$ m
Number, Volume	: 1.69 $\mu$ m	1.05 $\mu$ m
Length, Area	: 1.55 $\mu$ m	1.63 $\mu$ m
Length, Volume	: 2.25 $\mu$ m	1.78 $\mu$ m
Area, Volume	: 3.27 $\mu$ m	3.41 $\mu$ m
Volume, Moment	: 6.83 $\mu$ m	5.53 $\mu$ m

	MEDIAN Diameter	MODE	CONFIDENCE
--	-----------------	------	------------

Number	: 0.75 $\mu$ m	0.55 $\mu$ m	100.00%
Area	: 2.13 $\mu$ m	4.49 $\mu$ m	98.89%
Volume	: 4.90 $\mu$ m	18.08 $\mu$ m	99.25%

## Brinkmann

## Particle Size Analyzer

## PROCESS CHEMISTRY LABS PARTICLE ANALYSIS

VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H20.SBK

FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %

## PROBABILITY NUMBER DENSITY GRAPH


Name: SST,B000029,F0101,H20.SBK

Median : 0.75 $\mu$ m

6.5E+06 #/ml(100.0%)

Mean(nI): 0.95 $\mu$ mMode at 0.70  $\mu$ mS.D.(nI): 0.71 $\mu$ m<< SCALE RANGE ( $\mu$ m): 0 - 10 >>

Conf(nI):100.00 %



## Brinkmann

## Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %


## PROBABILITY NUMBER DISTRIBUTION GRAPH

Name: SST,B000029,F0101,H2O.SBK

6.5E+06 #/ml(100.0%)

Median : 0.75 $\mu$ mMean(nl): 0.95 $\mu$ mS.D.(nl): 0.71 $\mu$ m

Conf(nl):100.00 %

<< SCALE RANGE ( $\mu$ m): 0 - 10 >>

## Brinkmann

## Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010  
STATISTICS

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %

MEAN Diameter	S.D.
---------------	------

Brinkmann  
Particle Size Analyzer

PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
VIA BRINKMANN 2010  
STATISTICS

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.001

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-150	COUNTS	: 379228
TIME	: 08:14	ACQ. MODE	: SAMPLE	S.N.F.	: 0.72
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 1873 SEC	S.D.U.	: 4694
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 5	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 2.9E-03 %

MEAN Diameter	S.D.
---------------	------

Number, Length	:	1.17 $\mu$ m	0.89 $\mu$ m
Number, Area	:	1.47 $\mu$ m	0.94 $\mu$ m
Number, Volume	:	2.04 $\mu$ m	1.25 $\mu$ m
Length, Area	:	1.85 $\mu$ m	1.95 $\mu$ m
Length, Volume	:	2.69 $\mu$ m	2.13 $\mu$ m
Area, Volume	:	3.91 $\mu$ m	5.23 $\mu$ m
Volume, Moment	:	10.90 $\mu$ m	14.11 $\mu$ m

MEDIAN Diameter	MODE	CONFIDENCE		
Number	:	0.90 $\mu$ m	0.75 $\mu$ m	100.00%
Area	:	2.65 $\mu$ m	4.75 $\mu$ m	100.00%
Volume	:	5.47 $\mu$ m	4.75 $\mu$ m	99.92%

23 Sample red-brown, mushy with last line lumps

## Brinkmann

## Particle Size Analyzer

## PROCESS CHEMISTRY LABS PARTICLE ANALYSIS

VIA BRINKMANN 2010

SAMPLE NAME : SST,B000029,F0101,H2O.SBK

FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %


## PROBABILITY NUMBER DENSITY GRAPH

Name: SST,B000029,F0101,H2O.SBK

6.5E+06 #/ml(100.0%)

Mode at 0.70  $\mu\text{m}$ C << SCALE RANGE ( $\mu\text{m}$ ): 0 - 10 >>Median : 0.75  $\mu\text{m}$ Mean(n1): 0.95  $\mu\text{m}$ S.D.(n1): 0.71  $\mu\text{m}$ 

Conf(n1): 100.00 %



**Brinkmann**  
**Particle Size Analyzer**

**PROCESS CHEMISTRY LABS PARTICLE ANALYSIS  
 VIA BRINKMANN 2010**

SAMPLE NAME : SST,B000029,F0101,H2O.SBK  
 FILE NAME : F0101.002

DATE	: 04/12/1989	ACQ. RANGE	: 0.5-60	COUNTS	: 96549
TIME	: 08:49	ACQ. MODE	: SAMPLE	S.N.F.	: 0.84
CONFIG.	: 1 (0.7 S1)	ACQ. TIME	: 494 SEC	S.D.U.	: 5496
CELL TYPE	: MAGNETIC (3)	SAMPLE SIZE	: 4	CONCENTR.	: 6.5E+06 #/ml
SAMPLE TYPE	: REGULAR	REQ. CONF.	: 95.00%(V)	SOLIDS	: 1.6E-03 %

**PROBABILITY NUMBER DISTRIBUTION GRAPH**

Name: SST,B000029,F0101,H2O.SBK

6.5E+06 #/ml(100.0%)


Median : 0.75 $\mu$ m

Mean(n1): 0.95 $\mu$ m

S.D.(n1): 0.71 $\mu$ m

Conf(n1): 100.00 %

« SCALE RANGE ( $\mu$ m): 0 - 10 »



UNDIGESTED SAMPLE ANALYSIS

9 1 1 2 0 0 1 0 0 4

## Single Shell Tank Project

## Untreated Sample Results

Tank: 241-U-110

Core: 6

Segment: 4

Customer ID: 89-045

	Check Standard	Blank	Sample	Sample Duplicate	Check Standard
Laboratory ID:	F0100	F0121	F0101	F0102	F0292
pH	101.00%	6.83	11.96	12.22	100.90%
Laboratory ID:	F0100	F0309	F0101	F0102	F0292
%Water	96.63%	N/A	38.50%	37.30%	96.80%

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	AL10653
Procedure / Rev	LA-212-103 / A-3
Technologist	M. Franz
Date	01-02-90
Temperature	23.6
Starting Time	15:30
Ending Time	20:00
Chemist	R. E. Brandt

pH Analysis of the Solid Sample.

	Description	Lab. Id.
1	Initial LMCS Check Std	F0100
2	Reagent Blank	F0121
3	Sample 89-045	F0101
4	Duplicate 89-045	F0102
5	Sample 89-047	F0125
6	Duplicate 89-047	F0126
7	Sample 89-048	F0149
8	Duplicate 89-048	F0150
9	Sample 89-050	F0289
10	Duplicate 89-050	F0290
11	Final LMCS Check Std	F0292

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim Standard Type	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std	72C11-A/5.0mL			5.0 mL

Rev.E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 6-01-90
SST-102	Verified by: <u>Cary M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 6-01-90
	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Percent Water in Sample

Instrument	N/A
Procedure / Rev	LA-564-101 / D-1
Technologist	R. D. Hale
Date	01-03-90
Temperature	120 C
Starting Time	11:00; 01-02-90
Ending Time	11:00; 01-03-90
Chemist	R. E. Brandt

	Description	Lab. Id.
1	Initial LMCS Check Std	F0100
2	Reagent Blank	F0309
3	Sample 89-047	F0125
4	Duplicate 89-047	F0126
5	Sample 89-048	F0149
6	Duplicate 89-048	F0150
7	Sample 89-050	F0289
8	Duplicate 89-050	F0290
9	LMCS Check Std	F0292
10	Sample 89-045	F0101
11	Duplicate 89-045	F0102

	Description	Lab. Id.
12	Ending LMCS Check Std	F0292
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim 4/04/90	Primary Book	Second Book	Third Book	Final Volume	
	Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	of Standard
	LMCS Check Std	11C11AG/1.0g			1 gram

Rev.E SST-102	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 6-01-90
	Verified by: <u>Cathy M Seidel</u> Signature	C. M. Seidel Printed Name	Date: 6-01-90
	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

## KOH FUSION ANALYSIS

## Single Shell Tank Project

## Fusion Analysis

## Results On Laboratory Digestion

Tank: 241-U-110  
 Core: 6  
 Segment: 4  
 Customer ID: 89-045

	Check Standard F0105	Blank F0308	Sample F0106	Sample Duplicate F0107	Spike of Sample F0296	Check Standard F0297
Laboratory ID:						
Fusion Digestion			2.40	g/L	2.29	g/L
Total Alpha	111.90%	<1.00E-04 uci/L	*	uci/L	<2.00	uci/L
Total Beta	98.80%	<2.58E-04 uci/L	*	uci/L	1.97E+02	uci/L
Laboratory ID:	F0105	F0120	F0106	F0107	F0108	F0109
Total Alpha	113.70%	<1.79E-04 uci/L	<1.24	uci/L	**	uci/L
Total Beta	105.20%	6.65E-04 uci/L	2.03E+02	uci/L	**	uci/L
Laboratory ID:	F0081	F0192	F0106	F0107	F0084	F0181
GEA Cs-137	99.90%	<4.77E-02 uci/L	4.48E+01	uci/L	5.36E+01	uci/L
Laboratory ID:	F0105	F0120	F0106	F0107	F0108	F0297
Uranium (First Run)	98.70%	<1.04E+04 ug/L	3.34E-03	g/L	3.01E-03	g/L
Uranium (Second Run)	90.30%	<1.34E+02 ug/L	2.59E-03	g/L	**	95.00%
						92.30%

\* Invalid Results Due To High Deviation.

\*\* F0107 Not Rerun Due To Insufficient Quantity of Sample.

\*\*\* Ratio Of Spike To Sample Insufficient To Calculate Spike Recovery.

## Single Shell Tank Project

**Fusion Analysis**  
**Laboratory Results of Solids**  
**Units are Sample Wet Weight**

Tank:	241-U-110						
Core:	6						
Segment:	4						
Customer ID:	89-045						
	Check Standard	Blank	Sample		Sample Duplicate	Spike of Sample	Check Standard
Laboratory ID:	F0105	F0308	F0106		F0107	F0296	F0297
Fusion Digestion			2.40	g/L	2.29	g/L	
Total Alpha	111.90%	<1.00E-04 uci/L	*	uci/g	<8.73E-01 uci/g	98.00%	100.30%
Total Beta	98.80%	<2.58E-04 uci/L	*	uci/g	8.60E+01 uci/g	***	96.50%
Laboratory ID:	F0105	F0120	F0106		F0107	F0108	F0109
Total Alpha	113.70%	<1.79E-04 uci/L	<5.17E-01 uci/g	uci/g	**	uci/g	***
Total Beta	105.20%	6.65E-04 uci/L	8.46E+01 uci/g	uci/g	**	uci/g	***
Laboratory ID:	F0081	F0192	F0106		F0107	F0084	F0181
GEA Cs-137	99.90%	<4.77E-02 uci/L	1.87E+01 uci/g	uci/g	2.34E+01 uci/g	103.90%	96.00%
Laboratory ID:	F0105	F0120	F0106		F0107	F0108	F0297
Uranium (First Run)	98.70%	<1.04E+04 ug/L	1.39E+03 ug/g	ug/g	1.31E+03 ug/g	***	108.30%
Uranium (Second Run)	90.30%	<1.34E+02 ug/L	1.08E+03 ug/g	ug/g	**	95.00%	92.30%

\* Invalid Due To high Deviation.

\*\* F0107 Was Not Rerun Due To Insufficient Quantity of Sample.

\*\*\* Ratio Of Spike To Sample Insufficient to Calculate Spike Recovery.

This page intentionally left blank.

This page intentionally left blank.

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

## Fusion Dissolution

Instrument	N/A
Procedure / Rev	LA-549-141/A-0
Technologist	R. D. Hale
Date	01/03/90
Temperature	23 C
Starting Time	10:00
Ending Time	12:00
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Reagent Blank	F0168
2	Sample 89-045	F0106
3	Duplicate 89-045	F0107
4	Sample 89-047	F0130
5	Duplicate 89-047	F0131
6	Sample 89-048	F0154
7	Duplicate 89-048	F0155
8	Sample 89-050	F0294
9	Duplicate 89-050	F0295
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim Standard Type	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
N/A				

Rev E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 06/01/90
SST-102	Verified by: <u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 06/01/90
	Approved by: <u>L. H. Taylor</u> Signature	L. H. Taylor Printed Name	Date: 08-30-90

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	WA93415
Procedure / Rev	LA-508-101/C-2
Technologist	J. A. Hopkins
Date	01-05-90
Temperature	70 F
Starting Time	09:30
Ending Time	14:00
Chemist	S. A. Catlow

Total Alpha

Total Beta

Fusion Dissolution

Detector #18

	Description	Lab. Id.
1	Initial LMCS Check Std	F0105
2	Reagent Blank	F0308
3	Sample 89-045	F0106
4	Duplicate 89-045	F0107
5	Sample 89-047	F0130
6	Duplicate 89-047	F0131
7	Sample 89-048	F0154
8	Duplicate 89-048	F0155
9	Sample 89-050	F0294
10	Duplicate 89-050	F0295
11	Spike 89-050	F0296

	Description	Lab. Id.
12	Final LMCS Check Std	F0297
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim Rev E 4/04/90	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
Standard Type				
LMCS Check Std	83B44/10mL			10mL
Spike	83B44/10mL	F0296/0.100mL		10.1mL

Prepared by: S. A. Cervantes S. A. Cervantes Date: 06-01-90  
Signature Printed Name

Verified by: C. M. Seidel C. M. Seidel Date: 06-01-90  
Signature Printed Name

Approved by: L. H. Taylor L. H. Taylor Date: 08-30-90  
Signature Printed Name

## Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	N/A
Procedure / Rev	LA-508-101
Technologist	M. Franz
Date	01/25/90
Temperature	N/A
Starting Time	08:00
Ending Time	09:00
Chemist	S. A. Catlow

Total Alpha

Total Beta

Fusion Dissolution

Detector #18

F0107 was not rerun because of insufficient amount of sample.

	Description	Lab. Id.
1	Initial LMCS Check Std	F0105
2	Reagent Blank	F0120
3	Sample 89-045	F0106
4	Spike 89-045	F0108
5	Final LMCS Check Std	F0109
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim 4/04/90	Standard Type	Primary Book	Second Book	Third Book	Final Volume
		No. & Aliquot	No. & Aliquot	No. & Aliquot	of Standard
	LMCS Check Std	102B44/10 mL			10 mL
	Spike	102B44/10 mL	F0108/200 uL		10.2 mL

Rev.E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 07/07/90
	Verified by: <u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 07/07/90
SST-102	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

**Single Shell Tank  
Calibration Record**

**Phase  
I-A**

Analyte: **Am<sup>241</sup>**

Procedure **L0-508-002**

Revision: **A-0**

Instrument: **Detector #18**

Property Number: **WA93415**

Technologist: **R. A. Jones**

Payroll Number: **65801**

Date: **06/28/89**

Calibration Standard ID: **36B40A3; 36B40B3, 36B40C3**

Analyte Concentration: **61800, 110700, 161400 cpm**

Type of Calibration: **Efficiency**

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3			
4	<b>SEE ATTACHED SHEETS</b>		
5			
6			
7			
8			
9			
10			

Interim

Comments:

Rev. (Draft) **1/18/89**

Prepared by:




Signature

H. S. Rich

Printed Name

Date: **06-01-90**

Verified by:




Signature

C. M. Seidel

Printed Name

Date: **6-27-90**

Approved by:



Signature

L.H. Taylor

Printed Name

Date: **08-30-90**

## CALIBRATION SHEET FOR ALPHA/BETA SYSTEMS: USING PROCEDURE LQ-508-002

DETECTOR No. 18  
 RADIONUCLIDE: Am-241  
 HALF LIFE: 154497  
 COUNT TIME: 5  
 CPM BKG: 0.2  
 TIME ZERO DATE (HD): 15897  
 DATE COUNTED (HD): 16347  
 CALIBRATED BY: RA JONES HD 0 = 09/25/44

STANDARD ID	SIZE	DATE	TIME	COUNTS @ 0 DEG.	COUNTS @ 90 DEG.	COUNTS @ 180 DEG.	COUNTS @ 270 DEG.
36B40A8	1						
36B40B7	1						
36B40C7	1						
36B40A3	2	06/28/80	1542	67207	66768	67025	66645
36B40B3	2	06/28/80	1547	115573	116337	116289	116143
36B40C3	2	06/28/80	1552	162269	162819	162370	161593
36B40A6	5	06/28/80	1558	61627	62404	61970	61272
36B40B6	5	06/28/80	1603	118582	119217	118566	119430
36B40C5	5	06/28/80	1608	164322	165699	166216	166176

STANDARD ID	SIZE	STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
36B40A8	1"	60570	-0	1.00	-0	-0.0000
36B40B7	1"	109900	-0	1.00	-0	-0.0000
36B40C7	1"	159700	-0	1.00	-0	-0.0000
AVERAGE, 1" =		-0.0000 +/- @95%	0.0000	-97.62 %	ON	06/28/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
36B40A3	2"	61800	13382	1.00	13409	0.2170
36B40B3	2"	110700	23217	1.00	23264	0.2102
36B40C3	2"	161400	32452	1.00	32518	0.2015
AVERAGE, 2" =		0.2095 +/- @95%	0.0152	7.27 %	ON	06/28/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
36B40A6	5"	59470	12363	1.00	12388	0.2083
36B40B6	5"	109800	23790	1.00	23838	0.2171
36B40C5	5"	160100	33120	1.00	33187	0.2073
AVERAGE, 5" =		0.2109 +/- @95%	0.0106	5.01 %	ON	06/28/89
NEW EFFS FOR DET		18 Am-241	1" =	-0.0000	2" =	0.2095
				5" =		0.2109

9 1 1 2 2 3 3 1 1 1 9

Single Shell Tank Calibration Record		Phase I-A

Analyte:  $\text{Co}^{60}$

Procedure LQ-508-002

Revision: A-0

Instrument: Detector #18

Property Number: WA93415

Technologist: R. A. Jones

Payroll Number: 65801

Date: 06/28/89

Calibration Standard ID: 32B40A4, 32B40B3, 32B40C4

Analyte Concentration: 70480, 135100, 202400

Type of Calibration: Efficiency

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3			
4	SEE ATTACHED SHEETS		
5			
6			
7			
8			
9			
10			

Interim

Comments:

Rev. (Draft)

17889

Prepared by:




H. S. Rich

Printed Name

Date: 06-01-90

Verified by:




C. M. Seidel

Printed Name

Date: 6-27-90

Approved by:



L.H. Taylor

Date: 08-30-90

## CALIBRATION SHEET FOR ALPHA/BETA SYSTEMS: USING PROCEDURE LQ-508-002

DETECTOR No. 18  
 RADIONUCLIDE: Co-60  
 HALF LIFE: 1925  
 COUNT TIME: 5  
 CPM BKG: 5  
 CPM 1" BKG:  
 CALIBRATED BY: RA JONES HD 0 = 09/25/44

STANDARD ID	SIZE	DATE	TIME	COUNTS @ 0 DEG.	COUNTS @ 90 DEG.	COUNTS @ 180 DEG.	COUNTS @ 270 DEG.
100B40A2	1						
100B40B1	1						
100B40C1	1						
32B40A4	2	06/28/89	1510	95552	95030	96367	94943
32B40B3	2	06/28/89	1515	179993	179923	180564	179845
32B40C4	2	06/28/89	1521	266251	266109	266791	262848
32B40A5	5	06/28/89	1526	80056	79664	81559	79720
32B40B6	5	06/28/89	1531	159760	162820	161429	163674
32B40C5	5	06/28/89	1536	234482	235955	237348	236432

STANDARD ID	SIZE	STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
100B40A2	1"	67290	0	0.00	0	0.0000
100B40B1	1"	137800	0	0.00	0	0.0000
100B40C1	1"	199700	0	0.00	0	0.0000
AVERAGE, 1" =		0.0000 +/- @95%	0.0000	ERR % ..	ON	06/28/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
32B40A4	2"	70480	19090	1.18	22561	0.3201
32B40B3	2"	135100	36011	1.18	42560	0.3150
32B40C4	2"	202400	53095	1.18	62750	0.3100
AVERAGE, 2" =		0.3151 +/- @95%	0.0099	3.13 %	ON	06/28/89
STANDARD ID		STD VALUE	AVE CPM	DECAY CORR	DECAY CORR CPM	EFFICIENCY
32B40A5	5"	70160	16045	1.18	18963	0.2703
32B40B6	5"	135700	32379	1.18	38267	0.2820
32B40C5	5"	201900	47206	1.18	55790	0.2763
AVERAGE, 5" =		0.2762 +/- @95%	0.0115	4.16 %	ON	06/28/89
NEW EFFS FOR DET		18 Co-60	1" =	0.0000	2" =	0.3151
				5" =		0.2762

## Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	WA
Procedure / Rev	LA-548-121/C-3
Technologist	D.M. Southwick
Date	01/09/90
Temperature	72 F
Starting Time	12:30
Ending Time	14:00
Chemist	S. A. Catlow

GEA Analysis  
Fusion Dissolution

Samples were prepared in batch, but counted randomly.

Detectors 1, 2, and 4.

	Description	Lab. Id.
1	Initial LMCS Check Std.	F0081
2	Blank	F0192
3	Sample 89-044	F0082
4	Duplicate 89-044	F0083
5	Spike 89-044	F0084
6	Sample 89-045	F0106
7	Duplicate 89-045	F0107
8	Sample 89-049	F0178
9	Duplicate 89-049	F0179
10	Final LMCS Check Std.	F0181
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
LMCS Check Std.	89B44/500 uL			22 mL
Spike	89B44/100 uL	F0084/100 uL		22 mL

Prepared by:	<u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 07/07/90
Verified by:	<u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 07/07/90
Approved by:	<u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

Interim

Rev E 4/04/90

SSST-102  
44

**Single Shell Tank  
Calibration Record**

**Phase  
I-A**

**Analyte:** Isotope, Mixed Gamma

**Procedure:** LQ-508-003

**Revision:** A-0

**Instrument:** GEA Detector #1

**Property Number:** 401934

**Technologist:** JL Anderson

**Payroll Number:** 61413

**Date:** 3/2/89

**Calibration Standard ID:** 56B40 D1

**Analyte Concentration:** N/A

**Type of Calibration:** Gamma Energy Analysis (Efficiency)

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3	SEE ATTACHED		
4			
5			
6			
7			
8			
9			
10			

**Comments:**

Rev. (Draft) 1/18/89  
Prepared by: J.S. Rich  
Signature

H. S. Rich  
Printed Name

Date: 06-01-90

Verified by: C. M. Seidel  
Signature

C. M. Seidel  
Printed Name

Date: 6-27-90

Approved by: L.H. Taylor  
Signature

L.H. Taylor  
Printed Name

Date: 08-30-90

DETECTOR: 1  
 GEOMETRY CODE: 42  
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 2  
 CALIBRATION DATE: 14-Feb-89  
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL  
 STANDARD ID: 56840 D1

ENERGY (KEV)                    EFFICIENCY (COUNTS/GAMMA)

59.536	5.721347E-03
88.032	1.512568E-02
122.0614	2.041958E-02
165.853	1.856472E-02
279.1967	
391.668	1.042777E-02
513.99	7.856059E-03
661.65	6.838966E-03
898.021	5.300244E-03
1173.237	4.218416E-03
1332.501	3.785537E-03
1836.129	2.931033E-03

EQUATION 0-165 KEV

$$\begin{aligned} \text{LOG(EFF)} = & -5.343694\text{E+01} \\ & + 2.034704\text{E+01} * \text{LOG(ENERGY)} \\ & + -2.088264\text{E+00} * \text{LOG(ENERGY)}^2 \end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned} \text{LOG(EFF)} = & 8.372735\text{E+00} \\ & + -7.762489\text{E+00} * \text{LOG(ENERGY)} \\ & + 2.017698\text{E+00} * \text{LOG(ENERGY)}^2 \\ & + -2.447560\text{E-01} * \text{LOG(ENERGY)}^3 \\ & + 1.067720\text{E-02} * \text{LOG(ENERGY)}^4 \end{aligned}$$

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

DETECTOR: 1  
 GEOMETRY CODE: 43  
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 3  
 CALIBRATION DATE: 16-Feb-89  
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL  
 STANDARD ID: 56840 D1

ENERGY (KEV)                    EFFICIENCY (COUNTS/GAMMA)

59.536	1.397695E-03
88.032	3.641448E-03
122.0614	5.035820E-03
165.853	4.620516E-03
279.1967	
391.668	2.619018E-03
513.99	1.890740E-03
661.65	1.782478E-02
898.021	1.392563E-03
1173.237	1.117189E-03
1332.501	1.007670E-03
1836.129	7.782502E-04

EQUATION 0-165 KEV

$$\text{LOG(EFF)} = -5.354869\text{E+01}$$

+ 1.975356E+01 \*LOG(ENERGY)  
+ -2.020858E+00 \*LOG(ENERGY)^2

EQUATION 165-1836 KEV

LOG(EFF) = 4.001880E+01  
+ -2.857555E+01 \*LOG(ENERGY)  
+ 6.748440E+00 \*LOG(ENERGY)^2  
+ 7.173093E-01 \*LOG(ENERGY)^3  
+ 2.821780E-02 \*LOG(ENERGY)^4

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

LQ  
-----  
-  
-  
-  
-  
C  
C  
C  
2  
-  
9  
1  
6

**Single Shell Tank  
Calibration Record**

**Phase  
I-A**

**Analyte:** Mixed Isotope Standards

**Procedure** L0-508-003

**Revision:** A-0

**Instrument:** GEA Detector #2

**Property Number:** 401934

**Technologist:** JL Anderson

**Payroll Number:** 61413

**Date:** 9-1-88

**Calibration Standard ID:** 56B40 D1

**Analyte Concentration:** N/A


**Type of Calibration:** Gamma Energy Analysis (Efficiency)

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3	SEE ATTACHED		
4			
5			
6			
7			
8			
9			
10			

Interim


**Comments:**

Rev. (Draft) 1/18/89

**Prepared by:**   
Signature


H. S. Rich  
Printed Name

Date: 06-01-90

**Verified by:**   
Signature

C. M. Seidel  
Printed Name

Date: 6-27-90

**Approved by:**   
Signature

L.H. Taylor  
Printed Name

Date: 08-30-90

DETECTOR: 2  
 GEOMETRY CODE: 42  
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 2  
 CALIBRATION DATE: 21-Oct-88  
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL  
 STANDARD ID: 56B40 D1

ENERGY (KEV) EFFICIENCY (COUNTS/GAMMA)

59.536	3.417000E-03
88.032	1.090000E-02
122.0614	1.408000E-02
165.853	1.516000E-02
279.1967	9.929000E-03
391.668	7.578000E-03
513.99	5.875000E-03
661.65	4.927000E-03
898.021	3.727000E-03
1173.237	3.085000E-03
1332.501	2.683000E-03
1836.129	2.102000E-03

EQUATION 0-122 KEV

$$\begin{aligned} \text{LOG(EFF)} = & -6.654070\text{E+01} \\ & + 2.583780\text{E+01} * \text{LOG(ENERGY)} \\ & + -2.677550\text{E+00} * \text{LOG(ENERGY)}^2 \end{aligned}$$

EQUATION 122-1836 KEV

$$\begin{aligned} \text{LOG(EFF)} = & -1.050740\text{E+02} \\ & + 6.428950\text{E+01} * \text{LOG(ENERGY)} \\ & + -1.503170\text{E+01} * \text{LOG(ENERGY)}^2 \\ & + 1.533670\text{E+00} * \text{LOG(ENERGY)}^3 \\ & + -5.838530\text{E-02} * \text{LOG(ENERGY)}^4 \end{aligned}$$

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

DETECTOR: 2  
 GEOMETRY CODE: 43  
 GEOMETRY DECSRIPTION: 22 ML LIQUID, POS 3  
 CALIBRATION DATE: 28-Sep-88  
 ANALYST(S): J. L. ANDERSON/M. R. DOWELL  
 STANDARD ID: 56B40 D1

ENERGY (KEV) EFFICIENCY (COUNTS/GAMMA)

59.536	1.476000E-03
88.032	4.721000E-03
122.0614	6.589000E-03
165.853	6.613000E-03
279.1967	4.692000E-03
391.668	3.542000E-03
513.99	2.810000E-03
661.65	2.327000E-03
898.021	1.790000E-03
1173.237	1.437000E-03
1332.501	1.277000E-03
1836.129	9.824000E-04

EQUATION 0-165 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -5.826830\text{E+01} \\ & + 2.165450\text{E+01} * \text{LOG(ENERGY)} \\ & + -2.198930\text{E+00} * \text{LOG(ENERGY)}^2\end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -2.233890\text{E+01} \\ & + 1.174520\text{E+01} * \text{LOG(ENERGY)} \\ & + -2.739550\text{E+00} * \text{LOG(ENERGY)}^2 \\ & + 2.655450\text{E-01} * \text{LOG(ENERGY)}^3 \\ & + -9.668420\text{E-03} * \text{LOG(ENERGY)}^4\end{aligned}$$

# Single Shell Tank Calibration Record

Phase  
I-A

Analyte: Mixed Isotope Standards

Procedure	L0-508-003	Revision:	A-0
Instrument:	GEA Detector #4	Property Number:	401934
Technologist:	J. L. Anderson	Payroll Number:	61913
Date:	2-07-89		

Calibration Standard ID: 56B40 D1

Analyte Concentration: N/A

Type of Calibration: Gamma Energy Analysis (Efficiency)

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3	SEE ATTACHED		
4			
5			
6			
7			
8			
9			
10			

Interim

Comments:

1/18/89

Prepared by: H. S. Rich  
Signature

H. S. Rich

Printed Name

Date: 06-01-90

Rev. (Draft)

Verified by: C. M. Seidel  
Signature

C. M. Seidel

Printed Name

Date: 6-24-90

SST-103

Approved by: L.H. Taylor  
Signature

L.H. Taylor

Printed Name

Date: 08-30-90

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

DETECTOR: 4  
GEOMETRY CODE: 41  
GEOMETRY DESCRIPTION: 22 ML LIQUID, POS 1  
CALIBRATION DATE: 1-Sep-89  
ANALYST(S): J. L. ANDERSON/M. R. DOWELL  
STANDARD ID: 56B40 D1

ENERGY (KEV) EFFICIENCY (COUNTS/GAMMA)

59.536	2.682446E-02
88.032	8.210956E-02
122.0614	1.118411E-01
165.853	1.066653E-01
279.1967	
391.668	5.704220E-02
513.99	
661.65	3.685958E-02
898.021	2.541629E-02
1173.237	2.161710E-02
1332.501	1.973393E-02
1836.129	1.484468E-02

EQUATION 0-165 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -5.844056\text{E+01} \\ & + 2.310700\text{E+01} * \text{LOG(ENERGY)} \\ & + 2.371355\text{E+00} * \text{LOG(ENERGY)}^2\end{aligned}$$

EQUATION 165-1836 KEV

$$\begin{aligned}\text{LOG(EFF)} = & -1.718967\text{E+01} \\ & + 8.164155\text{E+00} * \text{LOG(ENERGY)} \\ & + -1.384196\text{E+00} * \text{LOG(ENERGY)}^2 \\ & + 7.025985\text{E-02} * \text{LOG(ENERGY)}^3\end{aligned}$$

GEA CALIBRATION RECORD

PROCEDURE LQ-508-003

\*  
\*  
\*       G A M M A   S P E C T R U M   A N A L Y S I S                   \*  
\*  
\* \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM

30-MAY-90 13:47:25

A N A L Y S I S   P A R A M E T E R S

MCA UNIT NUMBER: 2   /   ADC UNIT NUMBER: 1.0  
DETECTOR NUMBER: 1   /   GEOMETRY NUMBER: 42  
SPECTRUM SIZE: 4096 CHANNELS  
ORDER OF SMOOTHING FUNCTION: 5  
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK  
PEAK CONFIDENCE FACTOR: 85.0%  
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV  
-ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY

- ENVIRONMENTAL BACKGROUND SUBTRACTED  
  LLD CALCULATION PERFORMED  
- MEASURED ENERGY DIFFERENCES LISTED  
  MULTIPLET ANALYSIS PERFORMED

- ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD1998  
  ANALYZED BY:           AJ

- SAMPLE DESCRIPTION: E81-942.1L F81  
  GEOMETRY DESCRIPTION: 5/30/90  
- SAMPLE SIZE: 1.0000E-03 LI           / CONVERSION FACTOR: 5.0000E-01  
  STANDARD SIZE: 1.0000E+00 EA  
- ANALYSIS LIBRARY FILE: ANL000

- COLLECT STARTED ON 9-JAN-90 AT 13:41:59

- COLLECT LIVE TIME:   3000. SECONDS  
  REAL TIME:          3006. SECONDS  
  DEAD TIME:          0.20 %

DECAYED TO        0. DAYS,  0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 23-NOV-89  
EFFICIENCY CALIBRATION PERFORMED 2-MAR-89

SAMPLE: E81

DATA COLLECTED ON 9-JAN-90 AT 13:41:59

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

E-0081

## RADIONUCLIDE ANALYSIS REPORT

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI			ENERGY COMPARISON (KEV)	
	MEASURED	ERROR	DECAY CORRECTED	ERROR	EXPECT
AC-228	LLD<7.90E-01		LLD<7.90E-01		911.07
AG-108M	LLD<2.25E-01		LLD<2.25E-01		433.94
AG-110M	LLD<1.05E+00		LLD<1.05E+00		657.76
AM-241	LLD<9.30E-01		LLD<9.30E-01		59.54
AM-243	LLD<2.44E-01		LLD<2.44E-01		74.67
AR-41	LLD<1.44E-01		LLD<1.44E-01		1293.64
AU-198	LLD<2.00E-01		LLD<2.00E-01		411.80
BA-133	LLD<3.00E-01		LLD<3.00E-01		356.02
BA-139	LLD<5.85E-01		LLD<5.85E-01		165.85
BA-140	LLD<7.76E-01		LLD<7.76E-01		537.27
BA-141	LLD<6.03E-01		LLD<6.03E-01		190.23
BE-7	LLD<2.07E+00		LLD<2.07E+00		477.59
BI-207	LLD<1.98E-01		LLD<1.98E-01		569.70
BT-212	LLD<2.58E+00		LLD<2.58E+00		727.27
BI-214	LLD<8.68E-01		LLD<8.68E-01		609.32
CD-109	LLD<3.47E+00		LLD<3.47E+00		88.03
CE-139	LLD<1.32E-01		LLD<1.32E-01		165.85
CE-141	LLD<2.02E-01		LLD<2.02E-01		145.44
CEPR144	LLD<1.75E+00		LLD<1.75E+00		133.51
CO-56	LLD<1.93E-01		LLD<1.93E-01		846.76
CO-57	LLD<1.08E-01		LLD<1.08E-01		122.06
CO-58	LLD<1.86E-01		LLD<1.86E-01		810.75
CO-60	2.34E+01	+7.09E-01	2.34E+01	+7.09E-01	1332.50 -0.53
CR-51	LLD<1.51E+00		LLD<1.51E+00		1173.24 -0.62
CS-134	2.11E+01	+6.82E-01	2.11E+01	+6.82E-01	795.84 -0.60
					604.70 -0.52
CS-136	LLD<1.85E-01		LLD<1.85E-01		818.51
CS-137	3.81E+01	+8.42E-01	3.81E+01	+8.42E-01	661.65 -0.57
CS-138	LLD<1.83E-01		LLD<1.83E-01		1435.86
EU-152	LLD<2.61E-01		LLD<2.61E-01		1408.01
EU-154	LLD<3.59E-01		LLD<3.59E-01		1274.45
EU-155	LLD<4.22E-01		LLD<4.22E-01		105.31
FE-59	LLD<4.31E-01		LLD<4.31E-01		1099.25
HF-181	LLD<2.29E-01		LLD<2.29E-01		482.20
HG-203	LLD<1.81E-01		LLD<1.81E-01		279.20
I-131	LLD<2.32E-01		LLD<2.32E-01		364.48
I-132	LLD<2.14E-01		LLD<2.14E-01		667.69
I-133	LLD<2.29E-01		LLD<2.29E-01		529.69
I-134	LLD<2.70E-01		LLD<2.70E-01		847.03
I-135	LLD<4.50E-01		LLD<4.50E-01		1260.41
K-40	LLD<1.88E+00		LLD<1.88E+00		1460.75
KR-85	LLD<4.52E+01		LLD<4.52E+01		513.99
KR-85M	LLD<1.30E-01		LLD<1.30E-01		151.17
KR-87	LLD<4.97E-01		LLD<4.97E-01		402.58
KR-89	LLD<7.42E+00		LLD<7.42E+00		220.90
LA-140	LLD<9.05E-02		LLD<9.05E-02		1596.20

LA-142	LLD<4.45E-01	LLD<4.45E-01	641.83
MN-54	LLD<1.81E-01	LLD<1.81E-01	834.83
MN-56	LLD<2.17E-01	LLD<2.17E-01	846.76
NA-22	LLD<1.10E-01	LLD<1.10E-01	1274.55
NA-24	LLD<2.33E-01	LLD<2.33E-01	1368.60
NB-94	LLD<1.59E-01	LLD<1.59E-01	702.63
NB-95	LLD<1.57E-01	LLD<1.57E-01	765.78
NB-97	LLD<1.28E+00	LLD<1.28E+00	657.92
NP-238	LLD<8.29E-01	LLD<8.29E-01	984.45
NP-239	LLD<1.04E+00	LLD<1.04E+00	277.60
PA-233	LLD<4.70E-01	LLD<4.70E-01	311.98
PA-234M	LLD<3.78E+01	LLD<3.78E+01	1001.03
PB-210	LLD<5.44E+00	LLD<5.44E+00	465.03
PB-212	LLD<3.30E-01	LLD<3.30E-01	239.00
PB-214	LLD<5.01E-01	LLD<5.01E-01	351.92
PO-210	LLD<2.06E+04	LLD<2.06E+04	804.00
PO-214	LLD<8.05E+03	LLD<8.05E+03	799.70
PO-216	LLD<1.67E+04	LLD<1.67E+04	804.90
PU-239	LLD<1.45E+03	LLD<1.45E+03	129.30
PU-241	LLD<5.32E+04	LLD<5.32E+04	148.57
RA-224	LLD<3.57E+00	LLD<3.57E+00	240.99
RA-226	LLD<3.27E+00	LLD<3.27E+00	186.10
RB-88	LLD<6.33E-01	LLD<6.33E-01	1836.00
RB-89	LLD<1.04E+00	LLD<1.04E+00	1031.88
RN-220	LLD<1.69E+02	LLD<1.69E+02	549.73
RU-103	LLD<2.11E-01	LLD<2.11E-01	497.08
RURH106	LLD<3.62E+00	LLD<3.62E+00	621.80
SB-124	LLD<2.52E-01	LLD<2.52E-01	602.72
SB-125	LLD<1.70E+00	LLD<1.70E+00	176.33
SC-46	LLD<2.13E-01	LLD<2.13E-01	1120.45
SE-75	LLD<2.38E-01	LLD<2.38E-01	264.66
SN-113	LLD<2.89E-01	LLD<2.89E-01	391.67
SR-85	LLD<1.98E-01	LLD<1.98E-01	513.99
SR-91	LLD<3.54E-01	LLD<3.54E-01	555.60
SR-92	LLD<9.07E-02	LLD<9.07E-02	1383.94
TA-182	LLD<6.07E-01	LLD<6.07E-01	1121.30
TC-99M	LLD<1.16E-01	LLD<1.16E-01	140.51
TE-123M	LLD<1.25E-01	LLD<1.25E-01	159.00
TE-125M	LLD<3.28E+01	LLD<3.28E+01	109.27
TE-132	LLD<1.53E-01	LLD<1.53E-01	228.16
TH-228	LLD<1.09E+01	LLD<1.09E+01	84.37
TL-208	LLD<2.47E-01	LLD<2.47E-01	583.14
U-235	LLD<2.17E-01	LLD<2.17E-01	185.71
U-237	LLD<6.23E-01	LLD<6.23E-01	208.00
W-187	LLD<5.94E-01	LLD<5.94E-01	685.74
XE-131M	LLD<5.77E+00	LLD<5.77E+00	163.98
XE-133	LLD<3.81E-01	LLD<3.81E-01	81.00
XE-133M	LLD<1.35E+00	LLD<1.35E+00	233.21
XE-135	LLD<1.60E-01	LLD<1.60E-01	249.79
XE-138	LLD<1.21E+00	LLD<1.21E+00	258.41
Y-88	LLD<6.00E-02	LLD<6.00E-02	1836.06
Y-91	LLD<4.73E+01	LLD<4.73E+01	1204.90
Y-91M	LLD<2.68E-01	LLD<2.68E-01	555.60
ZN-65	LLD<5.03E-01	LLD<5.03E-01	1115.55
ZR-95	LLD<3.46E-01	LLD<3.46E-01	756.73
ZR-97	LLD<1.84E-01	LLD<1.84E-01	743.33

TOTAL      8.25E+01 +-1.29E+00      8.25E+01 +-1.29E+00

STANDARD DEVIATION = 0.04

## P E A K , A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1C	1126.46	562.70	1.41	645.	764.	10.1	CS-134, EU-152
2C	1138.69	568.82	1.41	640.	1301.	9.0	CS-134, BI-207
3	1209.44	604.18	1.47	583.	8369.	2.3	SB-124, CS-134
4	1323.29	661.08	1.54	419.	12295.	1.8	CS-137
4B		661.82			35.	46.4	
5C	1591.64	795.25	1.56	321.	5855.	3.0	CS-134
6C	1603.76	801.31	1.56	274.	553.	9.3	CS-134
7	2346.14	1172.62	1.68	255.	5341.	2.8	CO-60
8	2664.59	1331.97	1.88	55.	4914.	2.8	CO-60
9	2729.47	1364.44	2.40	13.	144.	18.2	CS-134
10	2800.60	1400.04	1.92	13.	59.	32.1	I-132,BI-214
11	2921.54	1460.58	1.84	11.	172.	16.2	K-40
c 11B		1461.77			182.	11.2	

— ERROR QUOTATION AT 1.96 SIGMA  
 — PEAK CONFIDENCE LEVEL AT 85.0%

C - MULTIPLET ANALYSIS CONVERGED NORMALLY

✓ B - ENVIRONMENTAL BACKGROUND PEAK

✓ BACKGROUND SUBTRACTION PERFORMED USING FILE BK0011

BACKGROUND DESCRIPTION: BK0011

BACKGROUND COLLECT STARTED ON 10-JAN-85 AT 12:00:00

BACKGROUND LIVE TIME: 6000. SECONDS

EBAR = \*\*\*\*\* MEV/DISINTEGRATION  
MAXIMUM PERMISSABLE ACTIVITY = 1.45E-09 UC/LI  
TOTAL MEASURED ACTIVITY = 8.25E+01 (+-1.29E+00) UC/LI  
% TECH. SPEC. = \*\*\*\*\* (+-\*\*\*\*\*)

ERROR QUOTATION AT 1.96 SIGMA  
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
1126.46	562.70	764.	10.1	3.27E+01
1138.69	568.82	1301.	9.0	5.62E+01
1603.76	801.31	553.	9.3	3.17E+01
2729.47	1364.44	144.	18.2	1.29E+01
2800.60	1400.04	59.	32.1	5.40E+00

PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
2921.54	1460.58	172.	16.2	1.63E+01

## \* \* \* \* \* GAMMA SPECTRUM ANALYSIS \* \* \* \* \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM

22-MAY-90 12:38:39

## A N A L Y S I S   P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 1.0  
DETECTOR NUMBER: 1 / GEOMETRY NUMBER: 42  
SPECTRUM SIZE: 4096 CHANNELS  
ORDER OF SMOOTHING FUNCTION: 5  
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK  
PEAK CONFIDENCE FACTOR: 85.0%  
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV  
ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY

ENVIRONMENTAL BACKGROUND SUBTRACTED  
LLD CALCULATION PERFORMED  
MEASURED ENERGY DIFFERENCES LISTED  
MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD1002  
ANALYZED BY: AJ

SAMPLE DESCRIPTION: F192  
GEOMETRY DESCRIPTION:  
SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E+00  
STANDARD SIZE: 1.0000E+00 EA  
ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 9-JAN-90 AT 19:58:58

COLLECT LIVE TIME: 3000. SECONDS  
REAL TIME: 3002. SECONDS  
DEAD TIME: 0.07 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 23-NOV-89  
EFFICIENCY CALIBRATION PERFORMED 2-MAR-89

222-S COUNTING ROOM

22-MAY-90 12:38:39

F-0192

P E A K   A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	1322.80	660.84	1.47	33.	33.	62.0	CS-137
1B		661.82			35.	46.4	
2	2921.42	1460.52	2.04	3.	146.	16.7	K-40
2B		1461.77			182.	11.2	

ERROR QUOTATION AT 1.96 SIGMA  
PEAK CONFIDENCE LEVEL AT 85.0%

B - ENVIRONMENTAL BACKGROUND PEAK

BACKGROUND SUBTRACTION PERFORMED USING FILE BK0011

BACKGROUND DESCRIPTION: BK0011

BACKGROUND COLLECT STARTED ON 10-JAN-85 AT 12:00:00

BACKGROUND LIVE TIME: 6000. SECONDS

222-S COUNTING ROOM

22-MAY-90 12:38:39

SAMPLE: F192

DATA COLLECTED ON 9-JAN-90 AT 19:58:58

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

## R A D I O N U C L I D E A N A L Y S I S R E P O R T

NUCLIDE	ACTIVITY CONCENTRATION IN $\mu\text{Ci}/\text{LI}$			ENERGY COMPARISON (KEV)		
	MEASURED	ERROR	DECAY CORRECTED	ERROR	EXPECT	DIFF
AC-228	LLD<1.03E-01		LLD<1.03E-01		911.07	
AG-108M	LLD<2.60E-02		LLD<2.60E-02		433.94	
AG-110M	LLD<3.82E-02		LLD<3.82E-02		657.76	
AM-241	LLD<1.54E-01		LLD<1.54E-01		59.54	
AM-243	LLD<4.29E-02		LLD<4.29E-02		74.67	
AR-41	LLD<3.48E-02		LLD<3.48E-02		1293.64	
AU-198	LLD<2.19E-02		LLD<2.19E-02		411.80	
BA-133	LLD<4.59E-02		LLD<4.59E-02		356.02	
BA-139	LLD<9.13E-02		LLD<9.13E-02		165.85	
BA-140	LLD<1.04E-01		LLD<1.04E-01		537.27	
BA-141	LLD<8.40E-02		LLD<8.40E-02		190.23	
BE-7	LLD<2.45E-01		LLD<2.45E-01		477.59	
BI-207	LLD<2.76E-02		LLD<2.76E-02		569.70	
BI-212	LLD<4.40E-01		LLD<4.40E-01		727.27	
BI-214	LLD<8.03E-02		LLD<8.03E-02		609.32	
CD-109	LLD<5.91E-01		LLD<5.91E-01		88.03	
CE-139	LLD<2.07E-02		LLD<2.07E-02		165.85	
CE-141	LLD<3.22E-02		LLD<3.22E-02		145.44	
CEPR144	LLD<2.84E-01		LLD<2.84E-01		133.51	
CO-56	LLD<2.34E-02		LLD<2.34E-02		846.76	
CO-57	LLD<1.83E-02		LLD<1.83E-02		122.06	
CO-58	LLD<2.81E-02		LLD<2.81E-02		810.75	
CO-60	LLD<1.31E-02		LLD<1.31E-02		1332.50	
CR-51	LLD<2.12E-01		LLD<2.12E-01		320.09	
CS-134	LLD<3.07E-02		LLD<3.07E-02		795.84	
CS-136	LLD<2.53E-02		LLD<2.53E-02		818.51	
CS-137	LLD<4.77E-02		LLD<4.77E-02		661.65	
CS-138	LLD<6.89E-02		LLD<6.89E-02		1435.86	
EU-152	LLD<1.30E-01		LLD<1.30E-01		1408.01	
EU-154	LLD<1.13E-01		LLD<1.13E-01		1274.45	
EU-155	LLD<6.59E-02		LLD<6.59E-02		105.31	
FE-59	LLD<5.71E-02		LLD<5.71E-02		1099.25	
HF-181	LLD<2.39E-02		LLD<2.39E-02		482.20	
HG-203	LLD<2.46E-02		LLD<2.46E-02		279.20	
I-131	LLD<2.90E-02		LLD<2.90E-02		364.48	
I-132	LLD<3.48E-02		LLD<3.48E-02		667.69	
I-133	LLD<2.85E-02		LLD<2.85E-02		529.69	
I-134	LLD<3.71E-02		LLD<3.71E-02		847.03	
I-135	LLD<1.29E-01		LLD<1.29E-01		1260.41	
K-40	LLD<9.20E-01		LLD<9.20E-01		1460.75	
R-85	LLD<8.59E+00		LLD<8.59E+00		513.99	
R-85M	LLD<2.16E-02		LLD<2.16E-02		151.17	
R-87	LLD<6.41E-02		LLD<6.41E-02		402.58	
R-89	LLD<1.09E+00		LLD<1.09E+00		220.90	
A-140	LLD<3.16E-02		LLD<3.16E-02		1596.20	
A-142	LLD<6.40E-02		LLD<6.40E-02		641.83	
N-54	LLD<2.94E-02		LLD<2.94E-02		834.83	

LA-142	LLD<4.45E-01	LLD<4.45E-01	641.83
MN-54	LLD<1.81E-01	LLD<1.81E-01	834.83
MN-56	LLD<2.17E-01	LLD<2.17E-01	846.76
NA-22	LLD<1.10E-01	LLD<1.10E-01	1274.55
NA-24	LLD<2.33E-01	LLD<2.33E-01	1368.60
NB-94	LLD<1.59E-01	LLD<1.59E-01	702.63
NB-95	LLD<1.57E-01	LLD<1.57E-01	765.78
NB-97	LLD<1.28E+00	LLD<1.28E+00	657.92
NP-238	LLD<8.29E-01	LLD<8.29E-01	984.45
NP-239	LLD<1.04E+00	LLD<1.04E+00	277.60
PA-233	LLD<4.70E-01	LLD<4.70E-01	311.98
PA-234M	LLD<3.78E+01	LLD<3.78E+01	1001.03
PB-210	LLD<5.44E+00	LLD<5.44E+00	465.03
PB-212	LLD<3.30E-01	LLD<3.30E-01	239.00
PB-214	LLD<5.01E-01	LLD<5.01E-01	351.92
PO-210	LLD<2.06E+04	LLD<2.06E+04	804.00
PO-214	LLD<8.05E+03	LLD<8.05E+03	799.70
PO-216	LLD<1.67E+04	LLD<1.67E+04	804.90
PU-239	LLD<1.45E+03	LLD<1.45E+03	129.30
PU-241	LLD<5.32E+04	LLD<5.32E+04	148.57
RA-224	LLD<3.57E+00	LLD<3.57E+00	240.99
RA-226	LLD<3.27E+00	LLD<3.27E+00	186.10
RB-88	LLD<6.33E-01	LLD<6.33E-01	1836.00
RB-89	LLD<1.04E+00	LLD<1.04E+00	1031.88
RN-220	LLD<1.69E+02	LLD<1.69E+02	549.73
RU-103	LLD<2.11E-01	LLD<2.11E-01	497.08
RURH106	LLD<3.62E+00	LLD<3.62E+00	621.80
SB-124	LLD<2.52E-01	LLD<2.52E-01	602.72
SB-125	LLD<1.70E+00	LLD<1.70E+00	176.33
SC-46	LLD<2.13E-01	LLD<2.13E-01	1120.45
SE-75	LLD<2.38E-01	LLD<2.38E-01	264.66
SN-113	LLD<2.89E-01	LLD<2.89E-01	391.67
SR-85	LLD<1.98E-01	LLD<1.98E-01	513.99
SR-91	LLD<3.54E-01	LLD<3.54E-01	555.60
SR-92	LLD<9.07E-02	LLD<9.07E-02	1383.94
TA-182	LLD<6.07E-01	LLD<6.07E-01	1121.30
TC-99M	LLD<1.16E-01	LLD<1.16E-01	140.51
TE-123M	LLD<1.25E-01	LLD<1.25E-01	159.00
TE-125M	LLD<3.28E+01	LLD<3.28E+01	109.27
TE-132	LLD<1.53E-01	LLD<1.53E-01	228.16
TH-228	LLD<1.09E+01	LLD<1.09E+01	84.37
TL-208	LLD<2.47E-01	LLD<2.47E-01	583.14
U-235	LLD<2.17E-01	LLD<2.17E-01	185.71
U-237	LLD<6.23E-01	LLD<6.23E-01	208.00
W-187	LLD<5.94E-01	LLD<5.94E-01	685.74
XE-131M	LLD<5.77E+00	LLD<5.77E+00	163.98
XE-133	LLD<3.81E-01	LLD<3.81E-01	81.00
XE-133M	LLD<1.35E+00	LLD<1.35E+00	233.21
XE-135	LLD<1.60E-01	LLD<1.60E-01	249.79
XE-138	LLD<1.21E+00	LLD<1.21E+00	258.41
Y-88	LLD<6.00E-02	LLD<6.00E-02	1836.06
Y-91	LLD<4.73E+01	LLD<4.73E+01	1204.90
Y-91M	LLD<2.68E-01	LLD<2.68E-01	555.60
ZN-65	LLD<5.03E-01	LLD<5.03E-01	1115.55
ZR-95	LLD<3.46E-01	LLD<3.46E-01	756.73
ZR-97	LLD<1.84E-01	LLD<1.84E-01	743.33
TOTAL	8.25E+01 +-1.29E+00	8.25E+01 +-1.29E+00	

STANDARD DEVIATION = 0.04

MN-56	LLD<2.64E-02	LLD<2.64E-02	846.76	F-0192
NA-22	LLD<3.52E-02	LLD<3.52E-02	1274.55	
NA-24	LLD<2.71E-02	LLD<2.71E-02	1368.60	
NB-94	LLD<2.93E-02	LLD<2.93E-02	702.63	
NB-95	LLD<2.36E-02	LLD<2.36E-02	765.78	
NB-97	LLD<4.63E-02	LLD<4.63E-02	657.92	
NP-238	LLD<1.09E-01	LLD<1.09E-01	984.45	
NP-239	LLD<1.43E-01	LLD<1.43E-01	277.60	
PA-233	LLD<6.07E-02	LLD<6.07E-02	311.98	
PA-234M	LLD<2.96E+00	LLD<2.96E+00	1001.03	
PB-210	LLD<6.08E-01	LLD<6.08E-01	465.03	
PB-212	LLD<4.85E-02	LLD<4.85E-02	239.00	
PB-214	LLD<6.66E-02	LLD<6.66E-02	351.92	
PO-210	LLD<2.43E+03	LLD<2.43E+03	804.00	
PO-214	LLD<2.78E+02	LLD<2.78E+02	799.70	
PO-216	LLD<1.27E+03	LLD<1.27E+03	804.90	
PU-239	LLD<2.51E+02	LLD<2.51E+02	129.30	
PU-241	LLD<9.09E+03	LLD<9.09E+03	148.57	
RA-224	LLD<5.50E-01	LLD<5.50E-01	240.99	
RA-226	LLD<4.77E-01	LLD<4.77E-01	186.10	
RB-88	LLD<1.98E-01	LLD<1.98E-01	1836.00	
RB-89	LLD<1.44E-01	LLD<1.44E-01	1031.88	
RN-220	LLD<2.35E+01	LLD<2.35E+01	549.73	
RU-103	LLD<2.63E-02	LLD<2.63E-02	497.08	
RURH106	LLD<5.53E-01	LLD<5.53E-01	621.80	
SB-124	LLD<2.50E-02	LLD<2.50E-02	602.72	
SB-125	LLD<2.71E-01	LLD<2.71E-01	176.33	
SC-46	LLD<3.04E-02	LLD<3.04E-02	1120.45	
SE-75	LLD<3.41E-02	LLD<3.41E-02	264.66	
SN-113	LLD<3.28E-02	LLD<3.28E-02	391.67	
SR-85	LLD<3.77E-02	LLD<3.77E-02	513.99	
SR-91	LLD<4.97E-02	LLD<4.97E-02	555.60	
SR-92	LLD<2.99E-02	LLD<2.99E-02	1383.94	
TA-182	LLD<9.06E-02	LLD<9.06E-02	1121.30	
TC-99M	LLD<1.94E-02	LLD<1.94E-02	140.51	
TE-123M	LLD<2.01E-02	LLD<2.01E-02	159.00	
TE-125M	LLD<5.19E+00	LLD<5.19E+00	109.27	
TE-132	LLD<2.28E-02	LLD<2.28E-02	228.16	
TH-228	LLD<1.90E+00	LLD<1.90E+00	84.37	
TE-208	LLD<2.89E-02	LLD<2.89E-02	583.14	
U-235	LLD<3.28E-02	LLD<3.28E-02	185.71	
U-237	LLD<8.54E-02	LLD<8.54E-02	208.00	
W-187	LLD<7.53E-02	LLD<7.53E-02	685.74	
XE-131M	LLD<8.52E-01	LLD<8.52E-01	163.98	
XE-133	LLD<6.62E-02	LLD<6.62E-02	81.00	
XE-133M	LLD<1.87E-01	LLD<1.87E-01	233.21	
XE-135	LLD<2.24E-02	LLD<2.24E-02	249.79	
XE-138	LLD<1.65E-01	LLD<1.65E-01	258.41	
Y-88	LLD<1.88E-02	LLD<1.88E-02	1836.06	
Y-91	LLD<1.28E+01	LLD<1.28E+01	1204.90	
Y-91M	LLD<3.76E-02	LLD<3.76E-02	555.60	
ZN-65	LLD<1.12E-01	LLD<1.12E-01	1115.55	
ZR-95	LLD<4.95E-02	LLD<4.95E-02	756.73	
ZR-97	LLD<2.33E-02	LLD<2.33E-02	743.33	
<hr/>				
TOTAL	0.00E-01 +-0.00E-01	0.00E-01 +-0.00E-01		

ERROR QUOTATION AT 1.96 SIGMA  
LLD CONFIDENCE LEVEL AT 85.0%

ALL DETECTED PEAKS WERE USED IN THE ANALYSIS

PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
1322.80	660.84	33.	62.0	1.59E+00
2921.42	1460.52	146.	16.7	1.38E+01

## \* \* \* \* \* GAMMA SPECTRUM ANALYSIS \* \* \* \* \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM

22-MAY-90 12:57:40

## A N A L Y S I S P A R A M E T E R S

MCA UNIT NUMBER: 2 / ADC UNIT NUMBER: 1.0  
DETECTOR NUMBER: 1 / GEOMETRY NUMBER: 42  
SPECTRUM SIZE: 4096 CHANNELS  
ORDER OF SMOOTHING FUNCTION: 5  
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK  
PEAK CONFIDENCE FACTOR: 85.0%  
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV  
ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY

ENVIRONMENTAL BACKGROUND SUBTRACTED

LLD CALCULATION PERFORMED

MEASURED ENERGY DIFFERENCES LISTED

MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD1000

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F84

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 1.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 9-JAN-90 AT 17:38:48

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3003. SECONDS

DEAD TIME: 0.10 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 23-NOV-89

EFFICIENCY CALIBRATION PERFORMED 2-MAR-89

222-S COUNTING ROOM

22-MAY-90 12:57:40

## P E A K   A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	93.93	46.86	1.14	351.	105.	53.0	
2C	1126.47	562.71	1.54	183.	186.	27.2	CS-134, EU-152
3C	1138.56	568.75	1.54	210.	285.	23.1	CS-134, BI-207
4	1209.44	604.18	1.49	221.	1797.	5.2	SB-124, CS-134
5	1323.32	661.10	1.47	149.	6884.	2.4	CS-137
5B		661.82			35.	46.4	
6C	1591.62	795.24	1.46	92.	1210.	7.1	CS-134
7C	1603.60	801.23	1.46	87.	134.	21.5	CS-134
8	2346.20	1172.65	1.87	56.	1128.	6.2	CO-60
9	2664.57	1331.96	1.86	13.	1060.	6.1	CO-60
10	2921.50	1460.56	1.63	5.	162.	16.0	K-40
10B		1461.77			182.	11.2	

—  
ERROR QUOTATION AT 1.96 SIGMA  
PEAK CONFIDENCE LEVEL AT 85.0%

C - MULTIPLET ANALYSIS CONVERGED NORMALLY  
B - ENVIRONMENTAL BACKGROUND PEAK

—  
BACKGROUND SUBTRACTION PERFORMED USING FILE BK0011  
BACKGROUND DESCRIPTION: BK0011  
BACKGROUND COLLECT STARTED ON 10-JAN-85 AT 12:00:00  
BACKGROUND LIVE TIME: 6000. SECONDS

222-S COUNTING ROOM

22-MAY-90 12:57:40

SAMPLE: F84

DATA COLLECTED ON 9-JAN-90 AT 17:38:48

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

## RADIONUCLIDE ANALYSIS REPORT

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI			ENERGY COMPARISON		
	MEASURED	ERROR	DECAY CORRECTED	ERROR	(KEV) EXPECT	DIFF
AC-228	LLD<2.06E+00		LLD<2.06E+00		911.07	
AG-108M	LLD<7.20E-01		LLD<7.20E-01		433.94	
AG-110M	LLD<3.92E+00		LLD<3.92E+00		657.76	
AM-241	LLD<3.84E+00		LLD<3.84E+00		59.54	
AM-243	LLD<1.05E+00		LLD<1.05E+00		74.67	
AR-41	LLD<4.93E-01		LLD<4.93E-01		1293.64	
AU-198	LLD<6.21E-01		LLD<6.21E-01		411.80	
BA-133	LLD<9.63E-01		LLD<9.63E-01		356.02	
BA-139	LLD<2.18E+00		LLD<2.18E+00		165.85	
BA-140	LLD<2.14E+00		LLD<2.14E+00		537.27	
BA-141	LLD<2.15E+00		LLD<2.15E+00		190.23	
BE-7	LLD<6.79E+00		LLD<6.79E+00		477.59	
BI-207	LLD<5.86E-01		LLD<5.86E-01		569.70	
BI-212	LLD<7.22E+00		LLD<7.22E+00		727.27	
Bi-214	LLD<2.25E+00		LLD<2.25E+00		609.32	
CD-109	LLD<1.50E+01		LLD<1.50E+01		88.03	
CE-139	LLD<4.94E-01		LLD<4.94E-01		165.85	
CE-141	LLD<7.99E-01		LLD<7.99E-01		145.44	
CEPR144	LLD<6.74E+00		LLD<6.74E+00		133.51	
CO-56	LLD<5.21E-01		LLD<5.21E-01		846.76	
CO-57	LLD<4.35E-01		LLD<4.35E-01		122.06	
CO-58	LLD<4.57E-01		LLD<4.57E-01		810.75	
CO-60	2.52E+01	+1.56E+00	2.52E+01	+1.56E+00	1332.50	-0.54
					1173.24	-0.59
CR-51	LLD<5.03E+00		LLD<5.03E+00		320.09	
CS-134	2.18E+01	+1.57E+00	2.18E+01	+1.57E+00	795.84	-0.61
					604.70	-0.52
CS-136	LLD<4.96E-01		LLD<4.96E-01		818.51	
CS-137	1.06E+02	+2.90E+00	1.06E+02	+2.90E+00	661.65	-0.55
CS-138	LLD<6.89E-01		LLD<6.89E-01		1435.86	
EU-152	LLD<1.18E+00		LLD<1.18E+00		1408.01	
EU-154	LLD<9.67E-01		LLD<9.67E-01		1274.45	
EU-155	LLD<1.78E+00		LLD<1.78E+00		105.31	
FE-59	LLD<1.15E+00		LLD<1.15E+00		1099.25	
HF-181	LLD<7.35E-01		LLD<7.35E-01		482.20	
HG-203	LLD<6.16E-01		LLD<6.16E-01		279.20	
I-131	LLD<7.08E-01		LLD<7.08E-01		364.48	
I-132	LLD<6.05E-01		LLD<6.05E-01		667.69	
I-133	LLD<6.64E-01		LLD<6.64E-01		529.69	
I-134	LLD<7.39E-01		LLD<7.39E-01		847.03	
I-135	LLD<1.39E+00		LLD<1.39E+00		1260.41	
K-40	LLD<9.13E+00		LLD<9.13E+00		1460.75	
KR-85	LLD<1.52E+02		LLD<1.52E+02		513.99	
KR-85M	LLD<4.99E-01		LLD<4.99E-01		151.17	
KR-87	LLD<1.49E+00		LLD<1.49E+00		402.58	
KR-89	LLD<2.49E+01		LLD<2.49E+01		220.90	
LA-140	LLD<3.68E-01		LLD<3.68E-01		1596.20	

LA-142	LLD<1.38E+00	LLD<1.38E+00	641.83
MN-54	LLD<5.70E-01	LLD<5.70E-01	834.83
MN-56	LLD<5.87E-01	LLD<5.87E-01	846.76
NA-22	LLD<3.91E-01	LLD<3.91E-01	1274.55
NA-24	LLD<6.13E-01	LLD<6.13E-01	1368.60
NB-94	LLD<4.96E-01	LLD<4.96E-01	702.63
NB-95	LLD<5.02E-01	LLD<5.02E-01	765.78
NB-97	LLD<4.75E+00	LLD<4.75E+00	657.92
NP-238	LLD<2.24E+00	LLD<2.24E+00	984.45
NP-239	LLD<3.38E+00	LLD<3.38E+00	277.60
PA-233	LLD<1.48E+00	LLD<1.48E+00	311.98
PA-234M	LLD<9.84E+01	LLD<9.84E+01	1001.03
PB-210	LLD<1.76E+01	LLD<1.76E+01	465.03
PB-212	LLD<1.13E+00	LLD<1.13E+00	239.00
PB-214	LLD<1.59E+00	LLD<1.59E+00	351.92
PO-210	LLD<5.80E+04	LLD<5.80E+04	804.00
PO-214	LLD<1.89E+04	LLD<1.89E+04	799.70
PO-216	LLD<4.45E+04	LLD<4.45E+04	804.90
PU-239	LLD<5.92E+03	LLD<5.92E+03	129.30
PU-241	LLD<2.02E+05	LLD<2.02E+05	148.57
RA-224	LLD<1.21E+01	LLD<1.21E+01	240.99
RA-226	LLD<1.13E+01	LLD<1.13E+01	186.10
RB-88	LLD<3.16E+00	LLD<3.16E+00	1836.00
RB-89	LLD<2.42E+00	LLD<2.42E+00	1031.88
RN-220	LLD<4.89E+02	LLD<4.89E+02	549.73
RU-103	LLD<6.34E-01	LLD<6.34E-01	497.08
RURH106	LLD<1.04E+01	LLD<1.04E+01	621.80
SB-124	LLD<6.53E-01	LLD<6.53E-01	602.72
SB-125	LLD<5.89E+00	LLD<5.89E+00	176.33
SC-46	LLD<5.37E-01	LLD<5.37E-01	1120.45
SE-75	LLD<7.77E-01	LLD<7.77E-01	264.66
SN-113	LLD<8.83E-01	LLD<8.83E-01	391.67
SR-85	LLD<6.69E-01	LLD<6.69E-01	513.99
SR-91	LLD<9.44E-01	LLD<9.44E-01	555.60
SR-92	LLD<4.94E-01	LLD<4.94E-01	1383.94
TA-182	LLD<1.73E+00	LLD<1.73E+00	1121.30
TC-99M	LLD<4.47E-01	LLD<4.47E-01	140.51
TE-123M	LLD<4.48E-01	LLD<4.48E-01	159.00
TE-125M	LLD<1.36E+02	LLD<1.36E+02	109.27
TE-132	LLD<5.24E-01	LLD<5.24E-01	228.16
TH-228	LLD<4.75E+01	LLD<4.75E+01	84.37
TL-208	LLD<6.98E-01	LLD<6.98E-01	583.14
U-235	LLD<7.44E-01	LLD<7.44E-01	185.71
U-237	LLD<2.17E+00	LLD<2.17E+00	208.00
W-187	LLD<1.66E+00	LLD<1.66E+00	685.74
XE-131M	LLD<2.00E+01	LLD<2.00E+01	163.98
XE-133	LLD<1.72E+00	LLD<1.72E+00	81.00
XE-133M	LLD<4.69E+00	LLD<4.69E+00	233.21
XE-135	LLD<5.41E-01	LLD<5.41E-01	249.79
XE-138	LLD<4.03E+00	LLD<4.03E+00	258.41
Y-88	LLD<3.00E-01	LLD<3.00E-01	1836.06
Y-91	LLD<1.52E+02	LLD<1.52E+02	1204.90
Y-91M	LLD<7.14E-01	LLD<7.14E-01	555.60
ZN-65	LLD<1.31E+00	LLD<1.31E+00	1115.55
ZR-95	LLD<8.81E-01	LLD<8.81E-01	756.73
ZR-97	LLD<4.69E-01	LLD<4.69E-01	743.33
TOTAL	1.53E+02 +-3.65E+00	1.53E+02 +-3.65E+00	

STANDARD DEVIATION = 0.04

EBAR = \*\*\*\*\* MEV/DISINTEGRATION  
 MAXIMUM PERMISSABLE ACTIVITY = 2.33E-09 UC/LI  
 TOTAL MEASURED ACTIVITY = 1.53E+02 (+-3.65E+00) UC/LI  
 % TECH. SPEC. = \*\*\*\*\* (+-\*\*\*\*\*)

ERROR QUOTATION AT 1.96 SIGMA  
 LLD CONFIDENCE LEVEL AT 85.0%

#### PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
93.93	46.86	105.	53.0	1.51E+01
1126.47	562.71	186.	27.2	7.94E+00
1138.56	568.75	285.	23.1	1.23E+01
1603.60	801.23	134.	21.5	7.67E+00

P

#### PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
2921.50	1460.56	162.	16.0	1.54E+01

## \* \* \* \* \* GAMMA SPECTRUM ANALYSIS \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM

22-MAY-90 13:13:25

## A N A L Y S I S   P A R A M E T E R S

MCA UNIT NUMBER: 1 / ADC UNIT NUMBER: 2.0  
DETECTOR NUMBER: 2 / GEOMETRY NUMBER: 42  
SPECTRUM SIZE: 4096 CHANNELS  
ORDER OF SMOOTHING FUNCTION: 5  
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK  
PEAK CONFIDENCE FACTOR: 85.0%  
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV  
~~ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY~~

ENVIRONMENTAL BACKGROUND SUBTRACTED

LLD CALCULATION PERFORMED

MEASURED ENERGY DIFFERENCES LISTED

MULTIPLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD2742

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F106

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 5.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 9-JAN-90 AT 17:40:21

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3003. SECONDS

DEAD TIME: 0.10 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 17-MAR-89

EFFICIENCY CALIBRATION PERFORMED 21-OCT-88

222-S COUNTING ROOM

22-MAY-90 13:13:25

SAMPLE: F106

DATA COLLECTED ON 9-JAN-90 AT 17:40:21

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

## R A D I O N U C L I D E A N A L Y S I S R E P O R T

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI			ENERGY COMPARISON (KEV)	
	MEASURED	ERROR	CORRECTED	ERROR	EXPECT
AC-228	LLD<3.51E-01		LLD<3.51E-01		911.07
AG-108M	LLD<2.03E-01		LLD<2.03E-01		433.94
AG-110M	LLD<1.45E+00		LLD<1.45E+00		657.76
AM-241	LLD<9.56E-01		LLD<9.56E-01		59.54
AM-243	LLD<2.46E-01		LLD<2.46E-01		74.67
AR-41	LLD<1.26E-01		LLD<1.26E-01		1293.64
AU-198	LLD<1.75E-01		LLD<1.75E-01		411.80
BA-133	LLD<2.31E-01		LLD<2.31E-01		356.02
BA-139	LLD<4.90E-01		LLD<4.90E-01		165.85
BA-140	LLD<5.28E-01		LLD<5.28E-01		537.27
BA-141	LLD<5.18E-01		LLD<5.18E-01		190.23
BE-7	LLD<1.79E+00		LLD<1.79E+00		477.59
BI-207	LLD<1.36E-01		LLD<1.36E-01		569.70
BI-212	LLD<1.38E+00		LLD<1.38E+00		727.27
BI-214	LLD<3.05E-01		LLD<3.05E-01		609.32
GD-109	LLD<3.12E+00		LLD<3.12E+00		88.03
CE-139	LLD<1.11E-01		LLD<1.11E-01		165.85
CE-141	LLD<1.83E-01		LLD<1.83E-01		145.44
CEPR144	LLD<1.44E+00		LLD<1.44E+00		133.51
CO-56	LLD<8.20E-02		LLD<8.20E-02		846.76
CO-57	LLD<9.31E-02		LLD<9.31E-02		122.06
CO-58	LLD<7.72E-02		LLD<7.72E-02		810.75
CO-60	LLD<8.02E-02		LLD<8.02E-02		1332.50
CR-51	LLD<1.23E+00		LLD<1.23E+00		320.09
CS-134	LLD<9.59E-02		LLD<9.59E-02		795.84
CS-136	LLD<6.92E-02		LLD<6.92E-02		818.51
CS-137	4.48E+01	+ -1.03E+00	4.48E+01	+ -1.03E+00	661.65
CS-138	LLD<2.37E-01		LLD<2.37E-01		1435.86
EU-152	LLD<6.05E-01		LLD<6.05E-01		1408.01
EU-154	LLD<2.52E-01		LLD<2.52E-01		1274.45
EU-155	LLD<4.35E-01		LLD<4.35E-01		105.31
FE-59	LLD<1.84E-01		LLD<1.84E-01		1099.25
HF-181	LLD<2.17E-01		LLD<2.17E-01		482.20
HG-203	LLD<1.37E-01		LLD<1.37E-01		279.20
I-131	LLD<1.72E-01		LLD<1.72E-01		364.48
I-132	LLD<7.62E-01		LLD<7.62E-01		667.69
I-133	LLD<1.48E-01		LLD<1.48E-01		529.69
I-134	LLD<1.14E-01		LLD<1.14E-01		847.03
I-135	LLD<3.96E-01		LLD<3.96E-01		1260.41
K-40	LLD<2.04E+00		LLD<2.04E+00		1460.75
KR-85	LLD<3.89E+01		LLD<3.89E+01		513.99
KR-85M	LLD<1.11E-01		LLD<1.11E-01		151.17
KR-87	LLD<3.88E-01		LLD<3.88E-01		402.58
KR-89	LLD<5.86E+00		LLD<5.86E+00		220.90
LA-140	LLD<8.48E-02		LLD<8.48E-02		1596.20
LA-142	LLD<3.34E-01		LLD<3.34E-01		641.83
MN-54	LLD<8.33E-02		LLD<8.33E-02		834.83

F-0106

MN-56	LLD<9.25E-02	LLD<9.25E-02	846.76
NA-22	LLD<8.97E-02	LLD<8.97E-02	1274.55
NA-24	LLD<6.61E-02	LLD<6.61E-02	1368.60
NB-94	LLD<8.40E-02	LLD<8.40E-02	702.63
NB-95	LLD<8.14E-02	LLD<8.14E-02	765.78
NB-97	LLD<1.65E+00	LLD<1.65E+00	657.92
NP-238	LLD<2.64E-01	LLD<2.64E-01	984.45
NP-239	LLD<8.51E-01	LLD<8.51E-01	277.60
PA-233	LLD<3.13E-01	LLD<3.13E-01	311.98
PA-234M	LLD<1.55E+01	LLD<1.55E+01	1001.03
PB-210	LLD<4.84E+00	LLD<4.84E+00	465.03
PB-212	LLD<2.61E-01	LLD<2.61E-01	239.00
PB-214	LLD<3.59E-01	LLD<3.59E-01	351.92
PO-210	LLD<8.20E+03	LLD<8.20E+03	804.00
PO-214	LLD<8.48E+02	LLD<8.48E+02	799.70
PO-216	LLD<5.56E+03	LLD<5.56E+03	804.90
PU-239	LLD<1.40E+03	LLD<1.40E+03	129.30
PU-241	LLD<4.31E+04	LLD<4.31E+04	148.57
RA-224	LLD<2.68E+00	LLD<2.68E+00	240.99
RA-226	LLD<2.60E+00	LLD<2.60E+00	186.10
RB-88	LLD<8.48E-01	LLD<8.48E-01	1836.00
RB-89	LLD<5.02E-01	LLD<5.02E-01	1031.88
RN-220	LLD<1.28E+02	LLD<1.28E+02	549.73
RU-103	LLD<1.67E-01	LLD<1.67E-01	497.08
RURH106	LLD<2.67E+00	LLD<2.67E+00	621.80
SB-124	LLD<1.29E-01	LLD<1.29E-01	602.72
SB-125	LLD<1.32E+00	LLD<1.32E+00	176.33
SC-46	LLD<1.00E-01	LLD<1.00E-01	1120.45
SE-75	LLD<1.99E-01	LLD<1.99E-01	264.66
SN-113	LLD<2.29E-01	LLD<2.29E-01	391.67
SR-85	LLD<1.71E-01	LLD<1.71E-01	513.99
SR-91	LLD<2.50E-01	LLD<2.50E-01	555.60
SR-92	LLD<1.47E-01	LLD<1.47E-01	1383.94
TA-182	LLD<3.17E-01	LLD<3.17E-01	1121.30
TC-99M	LLD<9.77E-02	LLD<9.77E-02	140.51
TE-123M	LLD<1.06E-01	LLD<1.06E-01	159.00
TE-125M	LLD<3.06E+01	LLD<3.06E+01	109.27
TE-132	LLD<1.24E-01	LLD<1.24E-01	228.16
TH-228	LLD<1.05E+01	LLD<1.05E+01	84.37
TC-208	LLD<1.63E-01	LLD<1.63E-01	583.14
U-235	LLD<1.71E-01	LLD<1.71E-01	185.71
U-237	LLD<4.87E-01	LLD<4.87E-01	208.00
W-187	LLD<2.83E-01	LLD<2.83E-01	685.74
XE-131M	LLD<4.75E+00	LLD<4.75E+00	163.98
XE-133	LLD<3.47E-01	LLD<3.47E-01	81.00
XE-133M	LLD<1.05E+00	LLD<1.05E+00	233.21
XE-135	LLD<1.24E-01	LLD<1.24E-01	249.79
XE-138	LLD<9.11E-01	LLD<9.11E-01	258.41
Y-88	LLD<8.05E-02	LLD<8.05E-02	1836.06
Y-91	LLD<2.71E+01	LLD<2.71E+01	1204.90
Y-91M	LLD<1.89E-01	LLD<1.89E-01	555.60
ZN-65	LLD<2.33E-01	LLD<2.33E-01	1115.55
ZR-95	LLD<1.44E-01	LLD<1.44E-01	756.73
ZR-97	LLD<9.09E-02	LLD<9.09E-02	743.33
<hr/>			
TOTAL	4.48E+01 +-1.03E+00	4.48E+01 +-1.03E+00	

EBAR = \*\*\*\*\* MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 1.16E-08 UC/LI

TOTAL MEASURED ACTIVITY = 4.48E+01 (+-1.03E+00) UC/LI

% TECH. SPEC. = \*\*\*\*\* (+-\*\*\*\*)

F-0106

ERROR QUOTATION AT 1.96 SIGMA  
LLD CONFIDENCE LEVEL AT 85.0%

ALL DETECTED PEAKS WERE USED IN THE ANALYSIS

PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
2921.98	1460.53	130.	19.1	1.70E+01

222-S COUNTING ROOM

22-MAY-90 13:13:25

F-0106

P E A K   A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	1324.33	661.78	1.63	84.	10396.	1.9	CS-137
1B		661.85			36.	13.9	
2	2921.98	1460.53	2.51	10.	130.	19.1	K-40
2B		1460.85			156.	3.8	

ERROR QUOTATION AT 1.96 SIGMA  
PEAK CONFIDENCE LEVEL AT 85.0%

B - ENVIRONMENTAL BACKGROUND PEAK

BACKGROUND SUBTRACTION PERFORMED USING FILE BK0012

BACKGROUND DESCRIPTION: BKG

BACKGROUND COLLECT STARTED ON 30-AUG-88 AT 16:46:00

BACKGROUND LIVE TIME: 60000. SECONDS

## \* \* \* \* \* GAMMA SPECTRUM ANALYSIS \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

22-MAY-90 12:39:33

## ANALYSIS PARAMETERS

MCA UNIT NUMBER: 1 / ADC UNIT NUMBER: 1.0

DETECTOR NUMBER: 4 / GEOMETRY NUMBER: 41

SPECTRUM SIZE: 4096 CHANNELS

ORDER OF SMOOTHING FUNCTION: 5

NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK

PEAK CONFIDENCE FACTOR: 85.0%

IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV

ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY

~~ENVIRONMENTAL BACKGROUND SUBTRACTED~~~~LLD CALCULATION PERFORMED~~~~MEASURED ENERGY DIFFERENCES LISTED~~~~MULTIPLET ANALYSIS PERFORMED~~

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD4880

ANALYZED BY: AJ

SAMPLE DESCRIPTION: F107

GEOMETRY DESCRIPTION:

SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 5.0000E-01

STANDARD SIZE: 1.0000E+00 EA

ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 9-JAN-90 AT 17:37:32

COLLECT LIVE TIME: 3000. SECONDS

REAL TIME: 3015. SECONDS

DEAD TIME: 0.50 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 26-DEC-89

EFFICIENCY CALIBRATION PERFORMED 1-SEP-89

## P E A K   A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	53.60	26.98	1.09	1497.	794.	15.2	
1B		27.06			123.	34.3	
2	1219.03	609.41	2.29	715.	271.	31.4	BI-214,
2B		609.25			197.	22.8	RU-103
3	1323.70	661.74	1.63	641.	93246.	0.6	CS-137
3B		661.35			379.	12.7	
4	2241.14	1120.54	1.09	110.	42.	87.2	BI-214, SC-46,
							TA-182
5	2921.77	1461.06	2.33	38.	801.	7.4	K-40
5B		1460.80			854.	7.1	

ERROR QUOTATION AT 1.96 SIGMA  
 PEAK CONFIDENCE LEVEL AT 85.0%

B<sup>-</sup> - ENVIRONMENTAL BACKGROUND PEAK

BACKGROUND SUBTRACTION PERFORMED USING FILE BK0014  
 BACKGROUND DESCRIPTION: BKG  
 BACKGROUND COLLECT STARTED ON 8-SEP-89 AT 12:00:00  
 BACKGROUND LIVE TIME: 3000. SECONDS

222-S COUNTING ROOM WESTINGHOUSE HANFORD

22-MAY-90 12:39:33

SAMPLE: F107

DATA COLLECTED ON 9-JAN-90 AT 17:37:32

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

## R A D I O N U C L I D E A N A L Y S I S R E P O R T

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI			ENERGY COMPARISON (KEV)		
	MEASURED	ERROR	CORRECTED	ERROR	EXPECT	DIFF
AC-228	LLD<1.45E-01		LLD<1.45E-01		911.07	
AG-108M	LLD<6.67E-02		LLD<6.67E-02		433.94	
AG-110M	LLD<4.79E-01		LLD<4.79E-01		657.76	
AM-241	LLD<2.98E-01		LLD<2.98E-01		59.54	
AM-243	LLD<7.18E-02		LLD<7.18E-02		74.67	
AR-41	LLD<3.20E-02		LLD<3.20E-02		1293.64	
AU-198	LLD<6.29E-02		LLD<6.29E-02		411.80	
BA-133	LLD<7.58E-02		LLD<7.58E-02		356.02	
BA-139	LLD<1.61E-01		LLD<1.61E-01		165.85	
BA-140	LLD<1.91E-01		LLD<1.91E-01		537.27	
BA-141	LLD<1.60E-01		LLD<1.60E-01		190.23	
BE-7	LLD<6.59E-01		LLD<6.59E-01		477.59	
BI-207	LLD<3.93E-02		LLD<3.93E-02		569.70	
BI-212	LLD<3.28E-01		LLD<3.28E-01		727.27	
BI-214	LLD<1.05E-01		LLD<1.05E-01		609.32	
CD-109	LLD<1.02E+00		LLD<1.02E+00		88.03	
CE-139	LLD<3.63E-02		LLD<3.63E-02		165.85	
CE-141	LLD<5.39E-02		LLD<5.39E-02		145.44	
CEPR144	LLD<4.62E-01		LLD<4.62E-01		133.51	
CO-56	LLD<2.45E-02		LLD<2.45E-02		846.76	
CO-57	LLD<3.01E-02		LLD<3.01E-02		122.06	
CO-58	LLD<2.35E-02		LLD<2.35E-02		810.75	
CO-60	LLD<2.80E-02		LLD<2.80E-02		1332.50	
CR-51	LLD<4.26E-01		LLD<4.26E-01		320.09	
CS-134	LLD<2.69E-02		LLD<2.69E-02		795.84	
CS-136	LLD<2.19E-02		LLD<2.19E-02		818.51	
CS-137	5.36E+01	+-5.35E-01	5.36E+01	+-5.35E-01	661.65	0.09
CS-138	LLD<6.14E-02		LLD<6.14E-02		1435.86	
EU-152	LLD<1.36E-01		LLD<1.36E-01		1408.01	
EU-154	LLD<8.28E-02		LLD<8.28E-02		1274.45	
EU-155	LLD<1.27E-01		LLD<1.27E-01		105.31	
FE-59	LLD<5.04E-02		LLD<5.04E-02		1099.25	
HF-181	LLD<7.86E-02		LLD<7.86E-02		482.20	
HG-203	LLD<4.81E-02		LLD<4.81E-02		279.20	
I-131	LLD<5.87E-02		LLD<5.87E-02		364.48	
I-132	LLD<7.62E-02		LLD<7.62E-02		667.69	
I-133	LLD<5.56E-02		LLD<5.56E-02		529.69	
I-134	LLD<3.21E-02		LLD<3.21E-02		847.03	
I-135	LLD<1.09E-01		LLD<1.09E-01		1260.41	
K-40	LLD<8.78E-01		LLD<8.78E-01		1460.75	
KR-85	LLD<1.23E+01		LLD<1.23E+01		513.99	
KR-85M	LLD<3.75E-02		LLD<3.75E-02		151.17	
KR-87	LLD<1.33E-01		LLD<1.33E-01		402.58	
KR-89	LLD<1.89E+00		LLD<1.89E+00		220.90	
LA-140	LLD<2.46E-02		LLD<2.46E-02		1596.20	
LA-142	LLD<1.07E-01		LLD<1.07E-01		641.83	
MN-54	LLD<2.44E-02		LLD<2.44E-02		834.83	

MN-56	LLD<2.77E-02	LLD<2.77E-02	846.76
NA-22	LLD<2.85E-02	LLD<2.85E-02	1274.55
NA-24	LLD<2.66E-02	LLD<2.66E-02	1368.60
NB-94	LLD<2.15E-02	LLD<2.15E-02	702.63
NB-95	LLD<2.37E-02	LLD<2.37E-02	765.78
NB-97	LLD<6.84E-01	LLD<6.84E-01	657.92
NP-238	LLD<9.11E-02	LLD<9.11E-02	984.45
NP-239	LLD<2.79E-01	LLD<2.79E-01	277.60
PA-233	LLD<1.16E-01	LLD<1.16E-01	311.98
PA-234M	LLD<4.06E+00	LLD<4.06E+00	1001.03
PB-210	LLD<1.73E+00	LLD<1.73E+00	465.03
PB-212	LLD<9.83E-02	LLD<9.83E-02	239.00
PB-214	LLD<1.29E-01	LLD<1.29E-01	351.92
PO-210	LLD<1.87E+03	LLD<1.87E+03	804.00
PO-214	LLD<2.25E+02	LLD<2.25E+02	799.70
PO-216	LLD<1.13E+03	LLD<1.13E+03	804.90
PU-239	LLD<3.99E+02	LLD<3.99E+02	129.30
PU-241	LLD<1.41E+04	LLD<1.41E+04	148.57
RA-224	LLD<9.56E-01	LLD<9.56E-01	240.99
RA-226	LLD<1.01E+00	LLD<1.01E+00	186.10
RB-88	LLD<1.77E-01	LLD<1.77E-01	1836.00
RB-89	LLD<1.05E-01	LLD<1.05E-01	1031.88
RN-220	LLD<3.96E+01	LLD<3.96E+01	549.73
RU-103	LLD<6.35E-02	LLD<6.35E-02	497.08
RURH106	LLD<8.36E-01	LLD<8.36E-01	621.80
SB-124	LLD<4.06E-02	LLD<4.06E-02	602.72
SB-125	LLD<4.21E-01	LLD<4.21E-01	176.33
SC-46	LLD<3.03E-02	LLD<3.03E-02	1120.45
SE-75	LLD<6.55E-02	LLD<6.55E-02	264.66
SN-113	LLD<8.01E-02	LLD<8.01E-02	391.67
SR-85	LLD<5.42E-02	LLD<5.42E-02	513.99
SR-91	LLD<7.82E-02	LLD<7.82E-02	555.60
SR-92	LLD<3.56E-02	LLD<3.56E-02	1383.94
TA-182	LLD<7.60E-02	LLD<7.60E-02	1121.30
TC-99M	LLD<2.99E-02	LLD<2.99E-02	140.51
TE-123M	LLD<3.34E-02	LLD<3.34E-02	159.00
TE-125M	LLD<9.72E+00	LLD<9.72E+00	109.27
TE-132	LLD<4.28E-02	LLD<4.28E-02	228.16
TH-228	LLD<3.17E+00	LLD<3.17E+00	84.37
TL-208	LLD<5.68E-02	LLD<5.68E-02	583.14
U-235	LLD<5.60E-02	LLD<5.60E-02	185.71
U-237	LLD<1.75E-01	LLD<1.75E-01	208.00
W-187	LLD<7.05E-02	LLD<7.05E-02	685.74
XE-131M	LLD<1.52E+00	LLD<1.52E+00	163.98
XE-133	LLD<1.15E-01	LLD<1.15E-01	81.00
XE-133M	LLD<3.42E-01	LLD<3.42E-01	233.21
XE-135	LLD<3.98E-02	LLD<3.98E-02	249.79
XE-138	LLD<3.31E-01	LLD<3.31E-01	258.41
Y-88	LLD<1.67E-02	LLD<1.67E-02	1836.06
Y-91	LLD<9.71E+00	LLD<9.71E+00	1204.90
Y-91M	LLD<5.92E-02	LLD<5.92E-02	555.60
ZN-65	LLD<7.59E-02	LLD<7.59E-02	1115.55
ZR-95	LLD<3.77E-02	LLD<3.77E-02	756.73
ZR-97	LLD<2.02E-02	LLD<2.02E-02	743.33
<hr/>			
TOTAL	5.36E+01 +-5.35E-01	5.36E+01 +-5.35E-01	

EBAR = \*\*\*\*\* MEV/DISINTEGRATION

MAXIMUM PERMISSABLE ACTIVITY = 1.16E-08 UC/LI

TOTAL MEASURED ACTIVITY = 5.36E+01 (+-5.35E-01) UC/LI

% TECH. SPEC. = \*\*\*\*\* (+-\*\*\*\*)

F-0107

ERROR QUOTATION AT 1.96 SIGMA  
LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
53.60	26.98	671.	19.0	6.99E+02
2241.14	1120.54	42.	87.2	6.14E-01

PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
1219.03	609.41	271.	31.4	2.28E+00
2921.77	1461.06	801.	7.4	1.48E+01

## \* \* \* \* \* GAMMA SPECTRUM ANALYSIS \* \* \* \* \*

CANBERRA SPECTRAN-F V2.06 SOFTWARE

222-S COUNTING ROOM WESTINGHOUSE HANFORD

22-MAY-90 12:58:18

## A N A L Y S I S   P A R A M E T E R S

MCA UNIT NUMBER: 1 / ADC UNIT NUMBER: 1.0  
DETECTOR NUMBER: 4 / GEOMETRY NUMBER: 41  
SPECTRUM SIZE: 4096 CHANNELS  
ORDER OF SMOOTHING FUNCTION: 5  
NUMBER OF BACKGROUND CHANNELS: 4 ON EACH SIDE OF PEAK  
PEAK CONFIDENCE FACTOR: 85.0%  
IDENTIFICATION ENERGY WINDOW: +- 1.50 KEV  
ERROR QUOTATION: 1.96 SIGMA UNCERTAINTY

ENVIRONMENTAL BACKGROUND SUBTRACTED  
LLD CALCULATION PERFORMED  
MEASURED ENERGY DIFFERENCES LISTED  
MULTIPLLET ANALYSIS PERFORMED

ANALYSIS OF SPECTRUM SAVED IN DISK FILE: SD4881  
ANALYZED BY: AJ

SAMPLE DESCRIPTION: F181  
GEOMETRY DESCRIPTION:  
SAMPLE SIZE: 1.0000E-03 LI / CONVERSION FACTOR: 5.0000E-01  
STANDARD SIZE: 1.0000E+00 EA  
ANALYSIS LIBRARY FILE: ANL000

COLLECT STARTED ON 9-JAN-90 AT 21:47:12

COLLECT LIVE TIME: 3000. SECONDS  
REAL TIME: 3031. SECONDS  
DEAD TIME: 1.02 %

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT

ENERGY CALIBRATION PERFORMED 26-DEC-89  
EFFICIENCY CALIBRATION PERFORMED 1-SEP-89

222-S COUNTING ROOM WESTINGHOUSE HANFORD

22-MAY-90 12:58:18

## P E A K   A N A L Y S I S

PK	CENTROID CHANNEL	ENERGY KEV	FWHM KEV	BACKGND COUNTS	NET AREA COUNTS	ERROR %	NUCLIDES
1	53.52	26.94	1.09	2647.	1374.	11.6	
1B		27.06			123.	34.3	
2	951.30	475.58	1.50	4649.	925.	22.0	CS-134
3C	1126.85	563.33	1.57	3190.	3491.	6.5	CS-134, EU-152
4C	1138.98	569.39	1.57	3111.	6628.	4.7	CS-134, BI-207
5	1209.75	604.77	1.55	3004.	41038.	1.0	CS-134
6	1323.62	661.70	1.66	2076.	63912.	0.8	CS-137
6B		661.35			379.	12.7	
7C	1591.87	795.82	1.72	1666.	29757.	1.5	CS-134
8C	1604.14	801.96	1.72	1599.	2941.	7.8	CS-134
9	2077.47	1038.67	1.87	1554.	402.	32.7	CS-134
10C	2335.39	1167.68	2.06	968.	634.	19.9	CS-134
11C	2346.28	1173.13	2.06	950.	27131.	1.4	CO-60
12	2664.80	1332.48	2.26	274.	24465.	1.3	CO-60
13	2729.85	1365.03	2.22	113.	815.	8.1	CS-134
14	2801.21	1400.74	2.60	92.	374.	13.2	BI-214
T5	2921.41	1460.88	2.44	62.	818.	7.6	K-40
15B		1460.80			854.	7.1	

ERROR QUOTATION AT 1.96 SIGMA  
 PEAK CONFIDENCE LEVEL AT 85.0%

MULTIPLET ANALYSIS CONVERGED NORMALLY  
 ENVIRONMENTAL BACKGROUND PEAK

BACKGROUND SUBTRACTION PERFORMED USING FILE BK0014  
 BACKGROUND DESCRIPTION: BKG  
 BACKGROUND COLLECT STARTED ON 8-SEP-89 AT 12:00:00  
 BACKGROUND LIVE TIME: 3000. SECONDS

222-S COUNTING ROOM WESTINGHOUSE HANFORD

22-MAY-90 12:58:18

SAMPLE: F181

DATA COLLECTED ON 9-JAN-90 AT 21:47:12

DECAYED TO 0. DAYS, 0.0000 HOURS BEFORE THE START OF COLLECT.

## RADIONUCLIDE ANALYSIS REPORT

NUCLIDE	ACTIVITY CONCENTRATION IN uCi/LI			ENERGY COMPARISON (KEV)		
	MEASURED	ERROR	DECAY CORRECTED	ERROR	EXPECT	DIFF
AC-228	LLD<3.84E-01		LLD<3.84E-01		911.07	
AG-108M	LLD<8.23E-02		LLD<8.23E-02		433.94	
AG-110M	LLD<4.07E-01		LLD<4.07E-01		657.76	
AM-241	LLD<3.90E-01		LLD<3.90E-01		59.54	
AM-243	LLD<9.32E-02		LLD<9.32E-02		74.67	
AR-41	LLD<5.96E-02		LLD<5.96E-02		1293.64	
AU-198	LLD<8.38E-02		LLD<8.38E-02		411.80	
BA-133	LLD<1.06E-01		LLD<1.06E-01		356.02	
BA-139	LLD<2.11E-01		LLD<2.11E-01		165.85	
BA-140	LLD<3.10E-01		LLD<3.10E-01		537.27	
BA-141	LLD<2.06E-01		LLD<2.06E-01		190.23	
BE-7	LLD<8.20E-01		LLD<8.20E-01		477.59	
BI-207	LLD<7.94E-02		LLD<7.94E-02		569.70	
BT-212	LLD<1.10E+00		LLD<1.10E+00		727.27	
BI-214	LLD<5.84E-01		LLD<5.84E-01		609.32	
CD-109	LLD<1.27E+00		LLD<1.27E+00		88.03	
CE-139	LLD<4.76E-02		LLD<4.76E-02		165.85	
CE-141	LLD<7.06E-02		LLD<7.06E-02		145.44	
CEPR144	LLD<5.97E-01		LLD<5.97E-01		133.51	
CO-56	LLD<8.92E-02		LLD<8.92E-02		846.76	
CO-57	LLD<3.81E-02		LLD<3.81E-02		122.06	
CO-58	LLD<8.29E-02		LLD<8.29E-02		810.75	
CO-60	2.26E+01	+ -3.28E-01	2.26E+01	+ -3.28E-01	1332.50	-0.02
—					1173.24	-0.11
CR-51	LLD<5.75E-01		LLD<5.75E-01		320.09	
CS-134	2.03E+01	+ -3.39E-01	2.03E+01	+ -3.39E-01	795.84	-0.03
					604.70	0.07
CS-136	LLD<7.91E-02		LLD<7.91E-02		818.51	
CS-137	3.67E+01	+ -4.07E-01	3.67E+01	+ -4.07E-01	661.65	0.05
CS-138	LLD<8.67E-02		LLD<8.67E-02		1435.86	
EU-152	LLD<2.07E-01		LLD<2.07E-01		1408.01	
EU-154	LLD<1.65E-01		LLD<1.65E-01		1274.45	
EU-155	LLD<1.63E-01		LLD<1.63E-01		105.31	
FE-59	LLD<1.96E-01		LLD<1.96E-01		1099.25	
HF-181	LLD<9.90E-02		LLD<9.90E-02		482.20	
HG-203	LLD<6.64E-02		LLD<6.64E-02		279.20	
I-131	LLD<8.20E-02		LLD<8.20E-02		364.48	
I-132	LLD<9.17E-02		LLD<9.17E-02		667.69	
I-133	LLD<8.80E-02		LLD<8.80E-02		529.69	
I-134	LLD<1.18E-01		LLD<1.18E-01		847.03	
I-135	LLD<2.10E-01		LLD<2.10E-01		1260.41	
K-40	LLD<8.92E-01		LLD<8.92E-01		1460.75	
KR-85	LLD<1.69E+01		LLD<1.69E+01		513.99	
KR-85M	LLD<4.91E-02		LLD<4.91E-02		151.17	
KR-87	LLD<1.81E-01		LLD<1.81E-01		402.58	
KR-89	LLD<2.55E+00		LLD<2.55E+00		220.90	
LA-140	LLD<3.53E-02		LLD<3.53E-02		1596.20	

LA-142	LLD<1.83E-01	LLD<1.83E-01	641.83
MN-54	LLD<8.51E-02	LLD<8.51E-02	834.83
MN-56	LLD<1.01E-01	LLD<1.01E-01	846.76
NA-22	LLD<5.34E-02	LLD<5.34E-02	1274.55
NA-24	LLD<7.43E-02	LLD<7.43E-02	1368.60
NB-94	LLD<6.84E-02	LLD<6.84E-02	702.63
NB-95	LLD<8.00E-02	LLD<8.00E-02	765.78
NB-97	LLD<5.74E-01	LLD<5.74E-01	657.92
NP-238	LLD<3.59E-01	LLD<3.59E-01	984.45
NP-239	LLD<3.79E-01	LLD<3.79E-01	277.60
PA-233	LLD<1.60E-01	LLD<1.60E-01	311.98
PA-234M	LLD<1.81E+01	LLD<1.81E+01	1001.03
PB-210	LLD<1.97E+00	LLD<1.97E+00	465.03
PB-212	LLD<1.27E-01	LLD<1.27E-01	239.00
PB-214	LLD<1.76E-01	LLD<1.76E-01	351.92
PO-210	LLD<7.22E+03	LLD<7.22E+03	804.00
PO-214	LLD<3.73E+03	LLD<3.73E+03	799.70
PO-216	LLD<6.39E+03	LLD<6.39E+03	804.90
PU-239	LLD<5.10E+02	LLD<5.10E+02	129.30
PU-241	LLD<1.82E+04	LLD<1.82E+04	148.57
RA-224	LLD<1.29E+00	LLD<1.29E+00	240.99
RA-226	LLD<1.29E+00	LLD<1.29E+00	186.10
RB-88	LLD<3.41E-01	LLD<3.41E-01	1836.00
RB-89	LLD<4.32E-01	LLD<4.32E-01	1031.88
RN-220	LLD<6.96E+01	LLD<6.96E+01	549.73
RU-103	LLD<8.38E-02	LLD<8.38E-02	497.08
RURH106	LLD<1.46E+00	LLD<1.46E+00	621.80
SB-124	LLD<1.92E-01	LLD<1.92E-01	602.72
SB-125	LLD<5.59E-01	LLD<5.59E-01	176.33
SG-46	LLD<1.12E-01	LLD<1.12E-01	1120.45
SE-75	LLD<8.99E-02	LLD<8.99E-02	264.66
SN-113	LLD<1.10E-01	LLD<1.10E-01	391.67
SR-85	LLD<7.40E-02	LLD<7.40E-02	513.99
SR-91	LLD<1.40E-01	LLD<1.40E-01	555.60
SR-92	LLD<4.77E-02	LLD<4.77E-02	1383.94
TA-182	LLD<3.00E-01	LLD<3.00E-01	1121.30
TC-99M	LLD<3.90E-02	LLD<3.90E-02	140.51
TE-123M	LLD<4.36E-02	LLD<4.36E-02	159.00
TE-125M	LLD<1.21E+01	LLD<1.21E+01	109.27
TE-132	LLD<5.78E-02	LLD<5.78E-02	228.16
TH-228	LLD<3.95E+00	LLD<3.95E+00	84.37
TL-208	LLD<1.00E-01	LLD<1.00E-01	583.14
U-235	LLD<7.18E-02	LLD<7.18E-02	185.71
U-237	LLD<2.32E-01	LLD<2.32E-01	208.00
W-187	LLD<2.39E-01	LLD<2.39E-01	685.74
XE-131M	LLD<1.98E+00	LLD<1.98E+00	163.98
XE-133	LLD<1.45E-01	LLD<1.45E-01	81.00
XE-133M	LLD<4.62E-01	LLD<4.62E-01	233.21
XE-135	LLD<5.41E-02	LLD<5.41E-02	249.79
XE-138	LLD<4.54E-01	LLD<4.54E-01	258.41
Y-88	LLD<3.21E-02	LLD<3.21E-02	1836.06
Y-91	LLD<2.36E+01	LLD<2.36E+01	1204.90
Y-91M	LLD<1.06E-01	LLD<1.06E-01	555.60
ZN-65	LLD<2.27E-01	LLD<2.27E-01	1115.55
ZR-95	LLD<1.37E-01	LLD<1.37E-01	756.73
ZR-97	LLD<7.73E-02	LLD<7.73E-02	743.33

TOTAL      7.96E+01 + -6.23E-01      7.96E+01 + -6.23E-01

STANDARD DEVIATION = 0.07

EBAR = \*\*\*\*\* MEV/DISINTEGRATION  
 MAXIMUM PERMISSABLE ACTIVITY = 1.45E-09 UC/LI  
 TOTAL MEASURED ACTIVITY = 7.96E+01 (+-6.23E-01) UC/LI  
 % TECH. SPEC. = \*\*\*\*\* (+-\*\*\*\*\*)

ERROR QUOTATION AT 1.96 SIGMA  
 LLD CONFIDENCE LEVEL AT 85.0%

PEAKS NOT USED IN ANALYSIS

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
53.52	26.94	1251.	13.2	1.32E+03
951.30	475.58	925.	22.0	6.21E+00
1126.85	563.33	3491.	6.5	2.73E+01
1138.98	569.39	6628.	4.7	5.24E+01
1604.14	801.96	2941.	7.8	3.19E+01
2077.47	1038.67	402.	32.7	5.52E+00
2335.39	1167.68	634.	19.9	9.66E+00
2729.85	1365.03	815.	8.1	1.42E+01
2801.21	1400.74	374.	13.2	6.66E+00

PEAKS ELIMINATED BY BACKGROUND SUBTRACTION

CENTROID CHANNEL	ENERGY KEV	NET AREA COUNTS	ERROR %	GAMMAS/SEC
2921.41	1460.88	818.	7.6	1.51E+01

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

## Uranium Analysis

Instrument	WA77344
Procedure / Rev	LA-925-106 / A-2
Technologist	M. Franz
Date	01-05-90
Temperature	N/A
Starting Time	08:00
Ending Time	15:00
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Initial LMCS Check Std	F0105
2	Reagent Blank	F0120
3	Sample 89-045	F0106
4	Duplicate 89-045	F0107
5	Spike 89-045	F0108
6	Sample 89-047	F0130
7	Duplicate 89-047	F0131
8	Sample 89-048	F0154
9	Duplicate 89-048	F0155
10	Sample 89-050	F0294
11	Duplicate 89-050	F0295

	Description	Lab. Id.
12	Final LMCS Check Std	F0297
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std	58B38 / 1 <u>L</u>			5.7 mL
Spike	58B83 /100 <u>uL</u>	F0108 / 1 <u>uL</u>		5.80 mL

Rev.E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 06-01-90
	Verified by: <u>Cary M Seidel</u> Signature	C.M. Seidel Printed Name	Date: 06-01-90
SST-102	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

## Uranium Analysis

Instrument	WA77344
Procedure / Rev	LA-925-106 / A-2
Technologist	S. Lai
Date	07-12-90
Temperature	N/A
Starting Time	Not Reported.
Ending Time	Not Reported.
Chemist	S. A. Catlow

	Description	Lab. Id.
1	Initial LMCS Check Std	F0105
2	Reagent Blank	F0120
3	Sample 89-045	F0106
4	Duplicate 89-045	F0107
5	Spike 89-045	F0108
6	Final LMCS Check Std	F0297
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book	Second Book	Third Book	Final Volume of Standard
	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std	58B38 / 1 uL			5.7 mL
Spike	54B38 / 100 uL	F0106 / 1 uL		5.8 mL

Rev.E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 8-09-90
	Verified by: <u>Cary M Seidel</u> Signature	C. M. Seidel Printed Name	Date: 8-09-90
SST-102	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

WATER DIGESTION TEST ANALYSIS

9 1 1 2 9 6 9 1 1 7 5

L8

**Single Shell Tank Project**

**Water Digestion**  
**Laboratory Results of Solids**  
**Units are Sample Wet Weight**

Tank 241-U-110

Core 6

Segment 4

Customer ID: 89-045

Laboratory Segment Serial No: F0101

Laboratory ID:	Check Standard F0110	Blank F0122	Sample F0111	Sample Duplicate F0112	Spike of Sample F0113	Check Standard F0066
Water Digestion			10.6 g/L	10.32 g/L	10.05 g/L	
Ion Chromatograph						
Fluoride	93.10%	<0.1 ppm	1.72E+04 ug/g	1.86E+04 ug/g	96.30%	93.10%
Chloride	97.90%	0.15 ppm	<9.53E+02 ug/g	<9.79E+02 ug/g	107.20%	96.30%
Nitrate	98.30%	<1.0 ppm	2.55E+04 ug/g	2.66E+04 ug/g	104.40%	105.20%
Sulfate	92.90%	<1.0 ppm	4.06E+03 ug/g	<9.79E+03 ug/g	101.30%	98.90%
Laboratory ID:	F0110	F0122	F0111	F0112	F0113	F0114
Phosphate	94.30%	<1.0 ppm	1.46E+05 ug/g	1.59E+05 ug/g	103.20%	101.00%
Laboratory ID:	F0110	F0122	F0111	F0112	F0113	F0114
Total Organic Carbon	97.40%	6.4 ug	7.26E+02 ug/g	6.93E+02 ug/g	94.50%	96.10%

9 | 1 2 3 5 0 | 1 5 6

8

## Single Shell Tank Project

Water Digestion  
Sample Results on Laboratory Digestion

Tank 241-U-110

Core 6

Segment 4

Customer ID: 89-045

Laboratory Segment Serial No. F0101

Laboratory	ID:	Check Standard F0110	Blank F0122	Sample F0111	Sample duplicate F0112	Spike of Sample F0113	Check Standard F0066
Water Digestion				10.60 g/L	10.32 g/L	10.05 g/L	

## Ion Chromatograph

Fluoride	93.10%	<0.1 ppm	1.82E+02 ppm	1.92E+02 ppm	96.30%	93.10%
Chloride	97.90%	0.15 ppm	<10.1 ppm	<10.1 ppm	107.20%	96.30%
Nitrate	98.30%	<1.0 ppm	2.70E+02 ppm	2.74E+02 ppm	104.40%	105.20%
Sulfate	92.90%	<1.0 ppm	4.31E+01 ppm	1.01E+02 ppm	101.30%	98.90%
Laboratory ID:	F0110	F0122	F0111	F0112	F0113	F0114
Phosphate	94.30%	<1.0 ppm	1.55E+03 ppm	1.64E+03 ppm	103.20%	101.00%
Laboratory ID:	F0110	F0122	F0111	F0112	F0113	F0114
Total Organic Carbon	97.40%	6.4 ug	7.70E-03 g/L	7.15E-03 g/L	94.50%	96.10%

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	N/A
Procedure / Rev	LA-504-101/A-2
Technologist	N.E. Wright
Date	01-09-90
Temperature	25 C
Starting Time	12:00; 01-08-90
Ending Time	09:00; 01-09-90
Chemist	H. S. Rich

## Water Digestion

Note: Sample is not spiked prior to digestion.

This procedure provides a sample to be spiked later with the appropriate elements.

	Description	Lab. Id.
1	Reagent Blank	F0122
2	Sample 89-045	F0111
3	Duplicate 89-045	F0112
4	Spike 89-045	F0113
5		
6		
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
N/A				
Spike (See Note)				

Prepared by:	<u>S. A. Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 06-01-90
Verified by:	<u>C. M. Seidel</u> Signature	C. M. Seidel Printed Name	Date: 06-01-90
Approved by:	<u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	AL10653
Procedure / Rev	LA-533-105/A-3
Technologist	6B107/N.E. Wright
Date	02/16/90
Temperature	24 C
Starting Time	11:30
Ending Time	14:00
Chemist	H. S. Rich

Ion Chromatograph Analysis

Water Digestion

\*Chromatogram Only

	Description	Lab. Id.
1	Eluent Blank	*
2	Initial LMCS Check Std.	F0110
3	Reagent Blank 89-045	F0122
4	Sample 89-045	F0111
5	Duplicate 89-045	F0112
6	Spike 89-045	F0113
7	Sample 89-043	F0063
8	Duplicate 89-043	F0064
9	Final LMCS Check Std.	F0066
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim		Primary Book	Second Book	Third Book	Final Volume of Standard
	Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
	LMCS Check Std.	6C11H1/100uL			10.10mL
	Spike	35C9-67/300uL	F0113/3.02mg		5.3mL

Rev. E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 06/06/90
SST-102	Verified by: <u>Cary M. Seidel</u> Signature	C.M. Seidel Printed Name	Date: 06/06/90
90	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

## Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	AL10653
Procedure / Rev	LA-533-105/A-3
Technologist	6B107/N.E. Wright
Date	02/22/90
Temperature	24 C
Starting Time	12:00
Ending Time	14:40
Chemist	H. S. Rich

Ion Chromatograph Analysis

Water Digestion

Re-Ran For Phosphate Only

\*Chromatogram Only

	Description	Lab. Id.
1	Eluent Blank	*
2	Initial LMCS Check Std.	F0110
3	Reagent Blank 89-045	F0122
4	Sample 89-045	F0111
5	Duplicate 89-045	F0112
6	Spike 89-045	F0113
7	Sample 89-043	F0063
8	Duplicate 89-043	F0064
9	Final LMCS Check Std.	F0114
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume of Standard
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std.	6C11HI/100ul			10.10mL
Spike	35C9-67/300ul	F0113/3.02mg		5.3mL

Interim

4/04/90

Rev.E

SST-102

Prepared by:	<u>S. A. Cervantes</u> Signature	S. A. Cervantes Printed Name	Date: 06/06/90
Verified by:	<u>C.M. Seidel</u> Signature	C.M. Seidel Printed Name	Date: 06/06/90
Approved by:	<u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

Single Shell Tank  
Calibration Record

Phase  
I-A

Analyte: Ion Chromatograph

Procedure LA-533-105

Revision: A-3

Instrument: Dionex

Property Number: WB24721

Technologist: Nora Wright

Payroll Number: 6B107

Date: 2-15-90; 2-22-90

Calibration Standard ID: 35C9-67 issued 02-06-90

Analyte Concentration: F=47; Cl=60.0; NO<sub>3</sub>=481.0; PO<sub>4</sub>=492.0; SO<sub>4</sub>=482.0 in ppm

Type of Calibration: Quadratic least squares

	Dilution	Concentration	Instrument Reading Units =
1			
2			
3			
4			
5	See Attached Calibration Sheets		
6			
7			
8			
9			
10			

Initials

Comments:

Rev. (Draft) 1/15/89

Prepared by: Shirley Cervantes S. A. Cervantes Date: 07-17-90  
Signature Printed Name

Rev. (Draft) 1/15/89

Verified by: Cary M Seidel C. M. Seidel Date: 07-17-90  
Signature Printed Name

ST-103

Approved by: L.H. Taylor L.H. Taylor Date: 08-30-90

## DIONEX METHOD PARAMETERS - SST.MET

## Detector Parameters

Number of Detectors.....	1
Detector 1 Type.....	CDM-1

## Report Options

Run Time (minutes).....	11.50
Detector 1 real time plot scale.....	20.00
Print Report.....	Yes
Print Replot.....	Yes
AutoScale Replot to Highest Peak.....	Yes
Print Retention Times on Chromatogram.....	Yes
List Peaks Not Found in this run.....	No
Report Unknowns found in run.....	No
Record Raw Data.....	Yes
Raw Data File Name: c:\WINDOWS\AI400\DATA\90021401.D08	
Record Result Data.....	No

## Integration Parameters

Sampling Rate (seconds).....	0.20
Peak Threshold (mV or uS/data pt interval).....	0.400
Starting Peak Width (seconds).....	10.0
Peak Area Reject.....	1000

## Integration Timed Events

Time	Description
------	-------------

## Calibration Parameters

External or Internal Calibration.....	External
Calibrate by Area or Height.....	Height
Replace Or Average Calibrations.....	Replace
Number Of Levels for Calibration.....	6
Calibration fit type.....	Quadratic
Response Factor for unknown peaks.....	0.0
Default Injection Volume.....	1.0
Default Dilution Factor.....	1.0
Area Reject for Reference Peaks.....	1000
Percent Retention Time Window for Reference Peaks.....	5.0

Component # 1      FLUORIDE      Retention Time    0.97  
 Reference Peak      FLUORIDE      Window Size       2.00%  
 Least Squares Slope = 3.60283E-004  
 Least Squares Intercept = 1.37365E-002  
 Ka =                          -4.29870E-009

Level	Amount	Area	Height
1	9.38100E-002	1336	246
2	2.33800E-001	3478	616
3	4.65300E-001	7271	1226
4	9.21500E-001	15622	2625
5	1.80760E+000	37180	5315
6	3.48130E+000	74116	11092

Component # 2      CHLORIDE      Retention Time    1.53  
 Reference Peak      FLUORIDE      Window Size       2.50%  
 Least Squares Slope = 6.25472E-004  
 Least Squares Intercept = -1.74138E-002  
 Ka =                          -1.51561E-008

Level	Amount	Area	Height
1	1.19800E-001	1267	200
2	2.98500E-001	2903	489
3	5.94100E-001	6430	1017
4	1.17650E+000	13713	2088
5	2.30760E+000	28370	4066
6	4.44430E+000	59291	9183

Component # 3      NITRITE      Retention Time    1.90  
 Reference Peak      FLUORIDE      Window Size       5.00%  
 Least Squares Slope = 8.73497E-004  
 Least Squares Intercept = 5.05282E-001  
 Ka =                          4.08720E-009

Level	Amount	Area	Height
1	1.28940E+000	8214	1202
2	3.21850E+000	20050	2969
3	6.39600E+000	43200	6258
4	1.26667E+001	91140	13081
5	2.48451E+001	181100	25092
6	4.78505E+001	353600	44776

Component # 4      NITRATE      Retention Time    3.28  
 Reference Peak      FLUORIDE      Window Size       7.00%  
 Least Squares Slope = 1.95501E-003  
 Least Squares Intercept = -6.94769E-003  
 Ka =                          4.87686E-009

Level	Amount	Area	Height
1	9.60100E-001	5253	468
2	2.39300E+000	15077	1234
3	4.76240E+000	30767	2460
4	9.43140E+000	63964	4753
5	1.84992E+001	132213	9250
6	3.56286E+001	274103	17468

Component # 5 PHOSPHATE Retention Time 5.10  
 Reference Peak FLUORIDE Window Size 7.00%  
 Least Squares Slope = 5.16373E-003  
 Least Squares Intercept = 1.15944E-001  
 Ka = -7.81115E-008

Level	Amount	Area	Height
1	9.82000E-001	2854	180
2	2.44770E+000	7216	437
3	4.87130E+000	16848	948
4	9.64710E+000	32539	1886
5	1.89223E+001	69406	3876
6	3.64434E+001	142872	8003

Component # 6 SULFATE Retention Time 6.43  
 Reference Peak FLUORIDE Window Size 5.00%  
 Least Squares Slope = 2.38873E-003  
 Least Squares Intercept = 2.12534E-001  
 Ka = -1.74248E-008


Level	Amount	Area	Height
1	9.62100E-001	6077	338
2	2.39800E+000	17483	922
3	4.77230E+000	36861	1890
4	9.45100E+000	79004	4005
5	1.85377E+001	163581	8158
6	3.57027E+001	343157	16953

IC Control File: C:\WINDOWS\AI400\METHOD\GROUT01.TE


Step	Time	Description
Init		CDM AutoOffset Off
Init		CDM Recorder Mark OFF
Init		CDM Temp. Comp. = 1.7 / Deg C
Init		CDM Recorder Range = 1.000 uS
Init		CDM Cell ON
Init		CMA Heater = 25 Deg. C
Init		Valve A ON
Init		Valve B ON
Init		Inject Valve OFF
Init		CIM Relay 1 OFF
Init		CIM Relay 2 OFF
Init		CIM AC 1 OFF
Init		CIM AC 2 OFF
Init		GPM Start
Init		GPM Hold Gradient Clock
Init		GPM Reset ON
1	0.0	CDM AutoOffset ON
1	0.0	GPM Reset OFF
2	0.1	Inject Valve ON
2	0.1	GPM Run Gradient Clock
3	3.0	Inject Valve OFF
4	3.5	CIM Relay 1 ON
5	4.0	CIM Relay 1 OFF

GpmFile: C:\WINDOWS\AI400\METHOD\GROUT01.GPM  
Lo Pressure Limit = 200  
Hi Pressure Limit = 2000  
Eluant 1 - DI WATER  
Eluant 2 - BICARBONATE  
Eluant 3 - CARBONATE  
Eluant 4 -


Time	Flow	%1	%2	%3	%4	Comment
0.0	2.0	84	8	8	0	
15.8	2.0	84	8	8	0	




Component: FLUORIDE  
Fit Type: Quadratic  
Conc = ( -4.298701e-009 \* Resp\* )  
( 3.602831e-004 \* Resp )+ 0.0137  
Standardization: Ext  
Calibration: Height




Component: CHLORIDE  
Fit Type: Quadratic  
Conc = ( -1.515615e-008 \* Resp\* )  
( 6.254716e-004 \* Resp )+ -0.0174  
Standardization: Ext  
Calibration: Height




Component: NITRITE  
Fit Type: Quadratic  
Conc = ( 4.087201e-009 \* Resp\* )+  
( 8.734971e-004 \* Resp )+ 0.5053  
Standardization: Ext  
Calibration: Height



Component: NITRATE  
Fit Type: Quadratic  
Conc = ( 4.876863e-009 \* Resp<sup>2</sup> ) +  
( 1.955009e-003 \* Resp ) + -0.0069  
Standardization: Ext  
Calibration: Height



Component: PHOSPHATE  
Fit Type: Quadratic  
Conc = ( -7.811152e-008 \* Resp<sup>2</sup> ) +  
( 5.163732e-003 \* Resp ) + 0.1159  
Standardization: Ext  
Calibration: Height



Component: SULFATE  
Fit Type: Quadratic  
Conc = ( -1.742481e-008 \* Resp<sup>2</sup> ) +  
( 2.388730e-003 \* Resp ) + 0.2125  
Standardization: Ext  
Calibration: Height

DATA REPROCESSED ON Wed Jun 06 11:31:37 1990


Sample Name: AUTOCALLR Date: Wed Feb 14 14:51:13 1990  
Data File : A:\90021401.D03  
Method : c:\windows\ai400\method\SST.met  
ACI Address: 1 System : 1 Inject#: 3 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3450  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	BL	REF PEAK	% DELTA RET TIME
1	0.98	FLUORIDE	9.381e-002	1.336e+003	246	1	0	0.00%
2	1.55	CHLORIDE	1.198e-001	1.267e+003	200	1	0	0.00%
3	3.52	NITRATE	9.601e-001	5.253e+003	468	1	0	0.00%
4	5.40	PHOSPHATE	9.820e-001	2.854e+003	180	1	0	0.00%
5	6.60	SULFATE	9.621e-001	6.077e+003	338	1	0	0.00%

File: A:\90021401.D03 Sample: AUTOCAL1R



DATA REPROCESSED ON Wed Jun 06 11:34:04 1990

Sample Name: AUTOCAL2R

Date: Wed Feb 14 15:03:32 1990

Data File : A:\90021401.D04

Method : c:\windows\ai400\method\SST.met

ACI Address: 1 System : 1 Inject #: 4 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes

Number of Data Points = 3451

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	2.338e-001	3.478e+003	616	1	0
2	1.53	CHLORIDE	2.985e-001	2.903e+003	489	1	0
3	3.48	NITRATE	2.393e+000	1.508e+004	1234	1	0
4	5.30	PHOSPHATE	2.448e+000	7.216e+003	437	1	0
5	6.55	SULFATE	2.398e+000	1.748e+004	922	1	0

File: A:\90021401.D04 Sample: AUTOCAL2R



DATA REPROCESSED ON Wed Jun 06 11:36:14 1990

Sample Name: AUTOCAL3R

Date: Wed Feb 14 15:15:50 1990

Data File : A:\90021401.D05

Method : c:\windows\ai400\method\SST.met

ACI Address: 1 System : 1 Inject#: 5 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes

Number of Data Points = 3451

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	4.653e-001	7.271e+003	1226	1	0 0.00%
2	1.53	CHLORIDE	5.941e-001	6.430e+003	1017	1	0 0.00%
3	3.45	NITRATE	4.762e+000	3.077e+004	2460	1	0 0.00%
4	5.32	PHOSPHATE	4.871e+000	1.685e+004	948	1	0 0.00%
5	6.57	SULFATE	4.772e+000	3.686e+004	1890	1	0 0.00%

File: A:\90021401.D05 Sample: AUTOCAL3R



DATA REPROCESSED ON Wed Jun 06 11:39:05 1990

=====  
Sample Name: AUTOCAL4R Date: Wed Feb 14 15:28:08 1990  
Data File : A:\90021401.D06  
Method : c:\windows\ai400\method\SST.met  
ACI Address: 1 System : 1 Inject#: 6 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3451

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	9.215e-001	1.562e+004	2625	1	0 0.00%
2	1.53	CHLORIDE	1.177e+000	1.371e+004	2088	1	0 0.00%
3	3.40	NITRATE	9.431e+000	6.396e+004	4753	1	0 0.00%
4	5.33	PHOSPHATE	9.647e+000	3.254e+004	1886	1	0 0.00%
5	6.57	SULFATE	9.451e+000	7.900e+004	4005	1	0 0.00%

File: A:\90021401.D06 Sample: AUTOCAL4R



DATA REPROCESSED ON Wed Jun 06 11:41:31 1990


Sample Name: AUTOCAL5R Date: Wed Feb 14 15:40:27 1990  
Data File : A:\90021401.D07  
Method : c:\windows\ai400\method\SST.met  
ACI Address: 1 System : 1 Inject#: 7 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	1.808e+000	3.718e+004	5315	2	0 0.00%
2	1.53	CHLORIDE	2.308e+000	2.837e+004	4066	2	0 0.00%
3	3.37	NITRATE	1.850e+001	1.322e+005	9250	1	0 0.00%
4	5.37	PHOSPHATE	1.892e+001	6.941e+004	3876	1	0 0.00%
5	6.58	SULFATE	1.854e+001	1.636e+005	8158	1	0 0.00%

File: A:\90021401.D07 Sample: AUTOCAL5R



DATA REPROCESSED ON Wed Jun 06 11:44:00 1990

Sample Name: AUTOCAL6R Date: Wed Feb 14 15:52:46 1990  
Data File : A:\90021401.D08  
Method : c:\windows\ai400\method\SST.met  
ACI Address: 1 System : 1 Inject#: 8 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3451

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA REF TIME
1	0.97	FLUORIDE	3.481e+000	7.412e+004	11092	2	0 0.00%
2	1.53	CHLORIDE	4.444e+000	5.927e+004	9183	2	0 0.00%
3	3.28	NITRATE	3.563e+001	2.741e+005	17468	1	0 0.00%
4	5.15	PHOSPHATE	3.644e+001	1.429e+005	8003	1	0 0.00%
5	6.43	SULFATE	3.570e+001	3.432e+005	16953	1	0 0.00%

File: A:\90021401.D08 Sample: AUTOCAL6R



## DIONEX SCHEDULE - A:\90021500.SCH

Inj #	Sample Name	Method Name	Data File	Vol.	Dil.	Int.Std.
1	SETUP	...\\sst	...\\900215001	1	0	
2	BLANK	...\\sst	...\\900215001	1	0	
3	LMCS/6C11HI	...\\sst	...\\900215001	101	0	
4	6294	...\\sst	...\\900215001	1111	0	
5	LMCS/6C11HI	...\\sst	...\\900215001	101	0	
6	LMCS/73C11J	...\\sst	...\\900215001	101	0	
7	122B	...\\sst	...\\900215001	1	0	
8	111	...\\sst	...\\900215001	101	0	
9	112D	...\\sst	...\\900215001	101	0	
10	113S	...\\sst	...\\900215001	101	0	
11	63	...\\sst	...\\900215001	101	0	
12	64D	...\\sst	...\\900215001	101	0	
13	LMCS/6C11HI	...\\sst	...\\900215001	101	0	
14	LMCS/73C11J	...\\sst	...\\900215001	101	0	

9 1 1 2 2 2 1 1 7 3

DATA REPROCESSED ON Wed Jun 06 15:29:58 1990


Sample Name: BLANK Date: Thu Feb 15 11:42:40 1990  
Data File : A:\90021500.D02  
Method : c:\windows\ai400\method\sst.met  
ACI Address: 1 System : 1 Inject #: 2 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA
					HEIGHT	BL PEAK
					RET TIME	
1	2.53		0.000e+000	6.737e+004	1200	1

File: A:\90021500.D02 Sample: BLANK



DATA REPROCESSED ON Wed Jun 06 15:32:41 1990


=====  
Sample Name: LMCS/6C11HI Date: Thu Feb 15 12:19:28 1990  
Data File : A:\90021500.D05  
Method : c:\windows\ai400\method\sst.met  
ACI Address: 1 System : 1 Inject#: 5 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3450  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	6.704e+001	1.122e+004	1845	1	0 0.00%
2	1.53	CHLORIDE	8.524e+001	9.110e+003	1426	1	0 0.00%
3	3.40	NITRATE	7.100e+002	4.452e+004	3567	1	0 0.00%
4	5.62	PHOSPHATE	6.809e+002	2.401e+004	1309	1	0 0.00%
5	6.80	SULFATE	6.960e+002	5.899e+004	2855	1	0 5.70%

File: A:\90021500.D05 Sample: LMCS/6C11HI



DATA REPROCESSED ON Wed Jun 06 15:35:12 1990


Sample Name: LMCS/73C11J	Date: Thu Feb 15 12:31:44 1990
Data File : A:\90021500.D06	
Method : c:\windows\ai400\method\sst.met	
ACI Address: 1	System : 1      Inject#: 6      Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes                  Number of Data Points = 3451  
Area reject = 1000                  One Data Point per 0.2 seconds  
Amount Injected = 1                  Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF BL	% DELTA PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	6.607e+001	1.145e+004	1817	1	0	0.00%
2	1.53	CHLORIDE	8.177e+001	8.913e+003	1367	2	0	0.00%
3	1.87	NITRITE	6.406e+002	5.139e+004	6485	2	0	0.00%
4	3.40	NITRATE	7.045e+002	4.629e+004	3540	1	0	0.00%
5	5.62	PHOSPHATE	6.738e+002	2.391e+004	1295	1	0	0.00%
6	6.78	SULFATE	6.847e+002	5.741e+004	2806	1	0	5.44%

File: A:\90021500.D06 Sample: LMCS/73C11J



DATA REPROCESSED ON Wed Jun 06 15:39:12 1990


=====  
Sample Name: 122B Date: Thu Feb 15 12:44:01 1990  
Data File : A:\90021500.D07  
Method : c:\windows\ai400\method\sst.met  
ACI Address: 1 System : 1 Inject#: 7 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3450  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
-------------	-------------	--------------	-------------------	------	---------------	------------	---------------------

File: A:\90021500.D07 Sample: 122B



DATA REPROCESSED ON Wed Jun 06 15:41:55 1990


Sample Name: 111	Date: Thu Feb 15 12:56:16 1990
Data File : A:\90021500.D08	
Method : c:\windows\ai400\method\sst.met	
ACI Address: 1 System : 1 Inject#: 8 Detector: CDM	

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3450  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF	% DELTA	BL PEAK	RET TIME
1	0.97	FLUORIDE	1.819e+002	3.635e+004	5295	1	0	0.00%	
2	1.87	NITRITE	9.151e+001	3.375e+003	458	1	0	0.00%	
3	3.43	NITRATE	2.703e+002	1.737e+004	1368	1	0	0.00%	
4	5.57	PHOSPHATE	1.411e+003	5.179e+004	2803	1	0	0.00%	
5	6.80	SULFATE	4.313e+001	1.605e+003	90	1	0	5.70%	

File: A:\90021500.D08 Sample: 111



DATA REPROCESSED ON Wed Jun 06 15:44:44 1990


=====  
Sample Name: 112D Date: Thu Feb 15 13:08:33 1990  
Data File : A:\90021500.D09  
Method : c:\windows\ai400\method\sst.met  
ACI Address: 1 System : 1 Inject#: 9 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3450  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF	% DELTA	BL PEAK	RET TIME
1	0.97	FLUORIDE	1.920e+002	3.594e+004	5614	1	0	0.00%	
2	1.87	NITRITE	9.335e+001	3.123e+003	479	1	0	0.00%	
3	3.40	NITRATE	2.740e+002	1.730e+004	1387	1	0	0.00%	
4	5.47	PHOSPHATE	1.539e+003	5.677e+004	3071	1	0	0.00%	

File: A:\90021500.D09 Sample: 112D



DATA REPROCESSED ON Fri Jun 08 15:46:32 1990


=====  
Sample Name: 113S Date: Thu Feb 15 13:20:49 1990  
Data File : A:\90021500.D10  
Method : c:\windows\ai400\method\sst.met  
ACI Address: 1 System : 1 Inject#: 10 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL	PEAK	RET TIME
1	0.97	FLUORIDE	4.069e+002	1.000e+005	13231	2	0	0.00%	
2	1.53	CHLORIDE	3.470e+002	4.432e+004	6565	3	0	0.00%	
3	1.87	NITRITE	9.066e+001	3.394e+003	448	4	0	0.00%	
4	3.27	NITRATE	2.951e+003	2.218e+005	14428	1	0	0.00%	
5	5.47	PHOSPHATE	3.877e+003	1.676e+005	8506	2	0	0.00%	
6	6.72	SULFATE	2.672e+003	2.529e+005	12046	2	0	0.00%	

File: A:\90021500.D10 Sample: 113S



DATA REPROCESSED ON Wed Jun 06 15:52:24 1990

Sample Name: LMCS/6C11HI

Date: Thu Feb 15 13:57:37 1990

Data File : A:\90021500.D13

Method : c:\windows\ai400\method\sst.met

ACI Address: 1 System : 1 Inject#: 13 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes

Number of Data Points = 3450

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL	PEAK	RET TIME
1	0.97	FLUORIDE	6.706e+001	1.166e+004	1845	1	0	0.00%	
2	1.52	CHLORIDE	8.379e+001	8.963e+003	1402	1	0	0.00%	
3	3.32	NITRATE	7.593e+002	4.930e+004	3813	1	0	0.00%	
4	5.12	PHOSPHATE	7.649e+002	2.497e+004	1477	1	0	-6.97%	
5	6.47	SULFATE	7.410e+002	5.990e+004	3050	1	0	0.00%	

File: A:\90021500.D13 Sample: LMCS/6C11HI



DATA REPROCESSED ON Wed Jun 06 15:54:23 1990


Sample Name: LMCS/73C11J	Date: Thu Feb 15 14:09:53 1990
Data File : A:\90021500.D14	
Method : c:\windows\ai400\method\sst.met	
ACI Address: 1 System : 1 Inject#: 14	Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL	PEAK	RET TIME
1	0.97	FLUORIDE	6.845e+001	1.126e+004	1885	1	0	0.00%	
2	1.52	CHLORIDE	8.427e+001	8.820e+003	1410	2	0	0.00%	
3	1.85	NITRITE	6.674e+002	5.062e+004	6772	2	0	0.00%	
4	3.33	NITRATE	7.171e+002	4.605e+004	3603	1	0	-0.00%	
5	5.17	PHOSPHATE	7.487e+002	2.728e+004	1445	1	0	-6.06%	
6	6.48	SULFATE	7.071e+002	5.762e+004	2903	1	0	0.00%	

File: A:\90021500.D14 Sample: LMCS/73C11J



## DIONEX SCHEDULE - A:\90022200.SCH

Inj #	Sample Name	Method Name	Data File	Vol.	Dil.	Int.Std.
1	SETUP	...\\sst	...\\900222001	1	0	
2	BLANK	...\\sst	...\\900222001	1	0	
3	LMCS/6C11HI	...\\sst	...\\900222001	101	0	
4	LMCS/73C11J	...\\sst	...\\900222001	101	0	
5	122B	...\\sst	...\\900222001	1	0	
6	111	...\\sst	...\\900222001	101	0	
7	112D	...\\sst	...\\900222001	101	0	
8	113S	...\\sst	...\\900222001	101	0	
9	63	...\\sst	...\\900222001	101	0	
10	64D	...\\sst	...\\900222001	101	0	
11	LMCS/6C11HI	...\\sst	...\\900222001	101	0	
12	LMCS/73C11J	...\\sst	...\\900222001	101	0	

DATA REPROCESSED ON Thu Feb 22 14:46:32 1990

=====

| Sample Name: LMCS/6C11HI Date: Thu Feb 22 12:00:58 1990|  
| Data File : C:\WINDOWS\AI400\DATA\90-24263.D03 |  
| Method : C:\WINDOWS\AI400\METHOD\SST.MET |  
| CIM Address: 1 System : 1 Cycle#: 3 Detector: CDM |  
=====


\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL	PEAK	RET TIME
1	0.97	FLUORIDE	103.7%	7.464e+001	1.175e+004	2064	1	0	-0.34%
2	1.52	CHLORIDE	100.1%	8.707e+001	9.453e+003	1457	1	0	-0.87%
3	3.35	NITRATE	104.3%	7.532e+002	4.905e+004	3782	1	0	2.13%
4	5.03	PHOSPHATE	102.1%	7.374e+002	2.374e+004	1422	1	0	-1.31%
5	6.30	SULFATE	101.9%	7.630e+002	6.028e+004	3146	1	0	-2.07%



DATA REPROCESSED ON Wed Jun 06 15:24:51 1990

Sample Name: LMCS/73C11J

Date: Thu Feb 22 12:13:15 1990

Data File : A:\90-24263.D04

Method : c:\windows\ai400\method\sst.met

ACI Address: 1 System : 1 Inject#: 4 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes

Number of Data Points = 3451

Area reject = 1000


One Data Point per 0.2 seconds

Amount Injected = 1

Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL PEAK	RET TIME
1	0.97	FLUORIDE	6.938e+001	1.095e+004	1912	1	0	0.00%
2	1.52	CHLORIDE	8.285e+001	9.123e+003	1386	2	0	0.00%
3	1.85	NITRITE	6.540e+002	5.051e+004	6629	2	0	0.00%
4	3.37	NITRATE	7.114e+002	4.580e+004	3575	1	0	0.00%
5	5.12	PHOSPHATE	6.922e+002	2.241e+004	1331	1	0	-6.97%
6	6.37	SULFATE	7.285e+002	5.774e+004	2996	1	0	0.00%

File: A:\90-24263.D04 Sample: LMCS/73C11J



DATA REPROCESSED ON Wed Jun 06 15:01:10 1990

=====

Sample Name: 122B	Date: Thu Feb 22 12:25:34 1990
Data File : A:\90022200.D05	
Method : c:\windows\ai400\method\sst.met	
ACI Address: 1	System : 1 Inject#: 5 Detector: CDM

=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3450

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 1

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	
					HEIGHT	BL PEAK	RET TIME

File: A:\90022200.D05 Sample: 122B



DATA REPROCESSED ON Thu Jun 07 10:43:16 1990


Sample Name: 111	Date: Thu Feb 22 12:37:50 1990
Data File : A:\90022200.D06	
Method : C:\WINDOWS\AI400\METHOD\SST.MET	
ACI Address: 1 System : 1 Inject#: 6	Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	HEIGHT	REF	% DELTA	RET TIME
1	0.97	FLUORIDE	2.135e+002	3.996e+004	6302	1	0	0.00%
2	1.85	NITRITE	9.251e+001	3.218e+003	469	1	0	0.00%
3	3.40	NITRATE	2.813e+002	1.800e+004	1423	1	0	0.00%
4	5.08	PHOSPHATE	1.550e+003	5.343e+004	3095	1	0	0.00%

File: A:\90022200.D06 Sample: 111



DATA REPROCESSED ON Thu Jun 07 10:41:35 1990

=====  
Sample Name: 112D Date: Thu Feb 22 12:50:06 1990  
Data File : A:\90022200.D07  
Method : C:\WINDOWS\AI400\METHOD\SST.MET  
ACI Address: 1 System : 1 Inject#: 7 Detector: CDM  
=====

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3450

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	RET TIME
1	0.97	FLUORIDE	2.148e+002	3.860e+004	6345	1	0 0.00%
2	1.83	NITRITE	9.352e+001	3.519e+003	480	1	0 0.00%
3	3.37	NITRATE	2.784e+002	1.712e+004	1408	1	0 0.00%
4	4.98	PHOSPHATE	1.639e+003	5.547e+004	3282	1	0 0.00%

File: A:\90022200.D07 Sample: 112D



DATA REPROCESSED ON Thu Jun 07 10:39:35 1990

Sample Name: 113S Date: Thu Feb 22 13:02:23 1990  
Data File : A:\90022200.D08  
Method : C:\WINDOWS\AI400\METHOD\SST.MET  
ACI Address: 1 System : 1 Inject #: 8 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3451

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF HEIGHT	BL PEAK	% DELTA RET TIME
1	0.97	FLUORIDE	4.640e+002	1.078e+005	15623	2	0 0.00%
2	1.52	CHLORIDE	3.794e+002	4.754e+004	7338	3	0 0.00%
3	1.85	NITRITE	9.137e+001	3.247e+003	456	4	0 0.00%
4	3.20	NITRATE	3.645e+003	2.721e+005	17684	1	0 0.00%
5	4.98	PHOSPHATE	4.739e+003	1.959e+005	10839	1	0 0.00%
6	6.28	SULFATE	3.501e+003	3.236e+005	16378	1	0 0.00%

File: A:\90022200.D08 Sample: 113S



DATA REPROCESSED ON Thu Jun 07 10:37:20 1990


Sample Name: LMCS/6C11HI Date: Thu Feb 22 13:39:15 1990  
Data File : A:\90022200.D11  
Method : C:\WINDOWS\AI400\METHOD\SST.MET  
ACI Address: 1 System : 1 Inject#: 11 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*

Stop time = 11.50 Minutes Number of Data Points = 3451  
Area reject = 1000 One Data Point per 0.2 seconds  
Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	RET TIME
1	0.97	FLUORIDE	7.285e+001	1.157e+004	2012	1	0 0.00%
2	1.52	CHLORIDE	8.801e+001	9.378e+003	1474	1	0 0.00%
3	3.33	NITRATE	7.474e+002	4.753e+004	3753	1	0 0.00%
4	5.10	PHOSPHATE	7.289e+002	2.432e+004	1405	1	0 0.00%
5	6.35	SULFATE	7.537e+002	5.982e+004	3105	1	0 0.00%

File: A:\90022200.D11 Sample: LMCS/6C11HI



DATA REPROCESSED ON Thu Jun 07 10:35:30 1990

Sample Name: LMCS/73C11J	Date: Thu Feb 22 13:51:32 1990
Data File : A:\90022200.D12	
Method : C:\WINDOWS\AI400\METHOD\SST.MET	
ACI Address: 1	System : 1 Inject#: 12 Detector: CDM

\*\*\*\*\* EXTERNAL STANDARD REPORT \*\*\*\*\*


Stop time = 11.50 Minutes Number of Data Points = 3450

Area reject = 1000 One Data Point per 0.2 seconds

Amount Injected = 1 Dilution factor = 101

PEAK NUM	RET TIME	PEAK NAME	CONC. in ug/ml	AREA	REF	% DELTA	BL PEAK	RET TIME
1	0.97	FLUORIDE	7.164e+001	1.124e+004	1977	1	0	0.00%
2	1.52	CHLORIDE	8.536e+001	8.720e+003	1428	2	0	0.00%
3	1.83	NITRITE	6.043e+002	4.953e+004	6097	2	0	0.00%
4	3.33	NITRATE	7.644e+002	5.078e+004	3838	1	0	0.00%
5	5.07	PHOSPHATE	7.552e+002	2.494e+004	1458	1	0	0.00%
6	6.35	SULFATE	7.506e+002	5.970e+004	3092	1	0	0.00%

File: A:\90022200.D12 Sample: LMCS/73C11J



# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Total Organic Carbon From Water Digestion

Instrument	WA399937
Procedure / Rev	LA-344-105/A-3
Technologist	E. Colvin
Date	05/30/90
Temperature	N/A
Starting Time	12:30
Ending Time	15:00
Chemist	R. E. Brandt

	Description	Lab. Id.
1	Reagent Blank	F0122
2	Initial LMCS Check Std	F0110
3	Sample 89-045	F0111
4	Duplicate 89-045	F0112
5	Spike 89-045	F0113
6	Final LMCS Check Std	F0114
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

	Primary Book	Second Book	Third Book	Final Volume of Standard
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
LMCS Check Std	80C11B/200uL			2.2 mL
Spike	80C11B/200uL	F0113/2.01mg		.5 mL

Prepared by: S. A. Cervantes S. A. Cervantes Date: 07/13/90  
 Signature Printed Name

Verified by: C. M. Seidel C. M. Seidel Date: 07/13/90  
 Signature Printed Name

Approved by: L. H. Taylor L. H. Taylor Date: 08-30-90  
 Signature Printed Name

Interim

Rev. E 4/04/90

SST-102

COULOMETER ANALYSIS REPORT  
TICTOC Rev. 0

Sample: F-122

Date: 05-30-1990

Time: 09:42:54

Blank = N/A  
% Difference = 10

Sample Size = 200 Dilution Factor = 1  
Min Readings = 7 Max Readings = 7

== Reading == Analysis Time == Coulometer == % Difference ==  
1 1.01 0.00 0.00

2	2.01	2.00	100.00
3	3.01	2.90	31.03
4	4.01	3.80	23.68
5	5.01	4.70	19.15
6	6.01	5.40	12.96
7	7.01	6.40	15.63

BLANK VALUE = 6.4 / 7.005615 = .9135529 ua/minute

Sample Run By: 8002B

COULOMETER ANALYSIS REPORT  
TICTOC Rev. 0

Sampler: F-110      Date: 05-30-1990      Timer: 10:03:59

Blank = .9135529      Sample Size = 200      Dilution Factor = 11  
% Difference = 10      Min Readings = 7      Max Readings = 7

== Reading == Analysis Time == Coulometer == % Difference ==  
1                1.01                0.00                0.00

2	2.01	46.70	100.00
3	3.01	52.60	11.22
4	4.01	55.40	5.05
5	5.01	57.20	3.15
6	6.01	58.40	2.05
7	7.01	59.50	1.85

$$(59.5 - 6.400056) / (11) / (200) = 2.920497 \text{ avl. Carbon}$$

$$(59.5 - 6.400056) / (11) / (200) (12) = .2433747 \text{ Molar Carbon}$$

Sample Run By: S0029

COULOMETER ANALYSIS REPORT  
EICTOC Rev. 0

Sample #: F-111

Date: 05-30-1990 Time: 10:30:26

Blank = .9133329  
% Difference = 10

Sample Size = 200 Dilution Factor = 1.1  
Min Readings = 7 Max Readings = 7

Readings ===== Analysis Time ===== Coulometer ===== % Difference =====  
J 1.01 0.00 0.00

Z	2.01	3.00	100.00
3	3.01	4.20	28.57
4	4.01	5.00	16.00
5	5.01	6.10	18.03
6	6.01	6.90	11.59
7	7.01	7.80	11.54

$$C_{ATM} = 6.4 \times (1.1) / (200) = 7.700001E-03 \text{ g/L Carbon}$$

$$C_{MOL} = 6.4 \times (1.1) / (200) (12) = 6.416667E-04 \text{ Molar Carbon}$$

Sample Run By: 80028

POOR COPY RECEIVED

COULOMETER ANALYSIS REPORT  
TIC TOC Rev. 0

Sample# F-112

Date: 05-30-1990 Time: 10:38:32

Blank = .9135529  
% Difference = 10

Sample Size = 200 Dilution Factor = 1.1  
Min Readings = 7 Max Readings = 2

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	32.01	3.10	100.00
3	3.01	4.30	27.91
4	4.01	5.10	15.69
5	5.01	6.10	16.39
6	6.01	6.90	11.59
7	7.01	7.70	10.39

$$C_{\text{vol}} = 6.400948 \times (1.1) \times (200) = 7.144785E-03 \text{ g/L Carbon}$$

$$C_{\text{mol}} = 6.400948 \times (1.1) \times (200) (12) = 5.953968E-04 \text{ Molar Carbon}$$

Sample Run# 80028

POOR COPY RECEIVED

COULOMETER ANALYSIS REPORT  
EIC10C Rev. 0

Sample #: E-313      Date: 05-30-1990      Time: 10:46:30

Blank = .9135529      Sample Size = 200      Dilution Factor = 1.1  
 % Difference = 10      Min Readings = ?      Max Readings = ?

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	89.80	100.00
3	3.01	106.40	15.60
4	4.01	113.70	6.42
5	5.01	117.50	3.23
6	6.01	119.60	1.76
7	7.01	121.10	1.29

**POOR COPY RECEIVED**

$$(121.1 - 6.400948) \times 1.1 / (200) = .6308448 \text{ a.t. Carbon}$$

$$(121.1 - 6.400948) \times 1.1 / (200) / (12) = .0525704 \text{ Molar Carbon}$$

Sample Run By: 80028

COULOMETRIC ANALYSIS REPORT  
TCTOC Rev. 0

Sample: F-114      Date: 05-30-1990      Time: 10:54:54

Blank = .9135529      Sample Size = 200      Dilution Factor = 11  
% Difference = 10      Min Readings = /      Max Readings = /

Reading	Analysis Time	Coulometer	% Difference
1	1.01	0.00	0.00
2	2.01	39.40	-
3	3.01	49.40	-20.00
4	4.01	53.70	8.01
5	5.01	56.20	4.05
6	6.01	57.80	2.77
7	7.01	58.80	1.70

( 58.8 - 6.400948 ) ( 11 ) / ( 200 ) = 2.881948 a/l Carbon

( 58.8 - 6.400948 ) ( 11 ) / ( 200 ) ( 12 ) = .2401623 Molar Carbon

Sample Run By: 80028

## ACID DIGESTION TEST ANALYSIS

## ICP Results

## ICP DATA SUMMARY

Date Analyzed:	March 12, 1990	Acid Digested Standard	F0115
Procedure:	LA-505-151/A-0	Reagent Blank	F0123
Analyst:	D. M. Southwick	Sample 89-045	F0116
Digestion	Acid Digestion	Duplicate 89-045	F0117
Procedure:	LA-505-159/A-0	Spike of F0116	F0118
		Acid Digested Standard	F0119

	Instrument Starting Standard	Acid Digest. Standard	Reagent BLANK	Wet Weight Sample	Wet Weight Sample Duplicate	Spike Recovery	LMCS ACID Digestion	Closing LMCS Standard
	%	%	ppm	ug/g	ug/g	%	%	%
Aluminum	100.47%		-0.71 LT	45520	52568	NOT CALC.	95.63%	106.46%
Antimony	97.68%		-0.59 LT	-751 LT	10 LT			101.29%
Arsenic	104.97%		-0.14 LT	-172 LT	3 LT			103.90%
Barium	103.89%		-0.04 LT	-58 LT	-6 LT	NOT CALC.	87.43%	93.37%
Beryllium	100.50%		0.00 LT	-6 LT	-2 LT			105.22%
Bismuth	102.33%	83.52%	-0.66 LT	21590	26620	NOT CALC.		106.21%
Boron	99.53%	90.89%	-0.03 LT	-65 LT	14 LT	NOT CALC.		93.06%
Cadmium	96.13%	93.88%	-0.02 LT	-35 LT	-8 LT	NOT CALC.		98.54%
Calcium	105.75%	115.16%	0.19	295	1113	111.89%		95.68%
Cerium	93.76%		-1.74 LT	-2483 LT	-288 LT	NOT CALC.	83.71%	97.80%
Chromium	100.34%		-0.18 LT	119	188	98.46%	86.40%	97.42%
Copper	106.21%	97.85%	-0.08 LT	-112 LT	2 LT	NOT CALC.		98.57%
Iron	101.13%		0.02	14837	16885	NOT CALC.	95.97%	97.66%
Lanthanum	99.52%	95.72%	-0.13 LT	-193 LT	-40 LT	NOT CALC.		108.92%
Lead	102.54%	91.72%	-0.35 LT	-83 LT	602	NOT CALC.		100.30%
Magnesium	100.98%	96.86%	0.02	479	449	87.17%		98.04%
Manganese	101.63%		0.00 LT	156	269	107.06%	94.31%	97.98%
Mercury	104.69%		-0.05 LT	-47 LT	9 LT			99.42%
Molybdenum	101.07%	95.80%	-0.04 LT	-59 LT	3 LT	NOT CALC.		101.72%
Nickel	99.60%		-0.07 LT	-104 LT	-5 LT	NOT CALC.	95.42%	99.59%
Phosphorus	105.42%	97.37%	-0.06 LT	53626	47284	NOT CALC.		99.90%
Potassium	102.78%		-2.53 LT	-3986 LT	-1211 LT	NOT CALC.		101.04%
Samarium	102.38%		-1.95 LT	-3036 LT	-675 LT			104.47%
Selenium	104.47%		-0.62 LT	-609 LT	99 LT			103.29%
Silver	100.85%		-0.11 LT	-186 LT	-48 LT	NOT CALC.		98.23%
Sodium	108.57%	114.60%	-0.50 LT	187906	168975	NOT CALC.		101.63%
Strontium	103.85%	94.59%	-0.02 LT	120	106	95.52%		94.36%
Sulfur	102.58%		-0.05 LT	113	278			103.02%
Tantalum	95.45%		-0.28 LT	-417 LT	-73 LT	NOT CALC.	84.22%	93.36%
Thallium	101.40%		-2.02 LT	-3223 LT	-906 LT			104.50%
Thorium	97.10%		-1.11 LT	-1702 LT	-403 LT			104.62%
Tin	100.56%		-0.03 LT	-112 LT	-11 LT	NOT CALC.	101.08%	104.92%
Titanium	96.10%		0.03	-70 LT	8 LT	NOT CALC.	89.66%	101.92%
Uranium	98.51%		-10.90 LT	-15293 LT	-2485 LT			99.44%
Vanadium	102.61%		-0.09 LT	-128 LT	-46 LT			105.49%
Zinc	98.20%	94.36%	0.00	142	5380	98.35%		98.06%
Zirconium	96.38%		-0.19 LT	-276 LT	97	NOT CALC.	86.94%	101.68%

LT: Less Than

NC: Not Calibrated

NOT CALC: Not Calculated

# Instrument Standards Outside Control Limits

# Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

## Acid Digestion

Instrument	N/A
Procedure / Rev	LA-505-158/A-1
Technologist	D. Southwick
Date	01/31/90
Temperature	72 C
Starting Time	08:00
Ending Time	14:45
Chemist	S.A. Jones

	Description	Lab. Id.
1	Reagent Blank	F0075
2	Sample 89-043	F0068
3	Duplicate 89-043	F0069
4	Sample 89-045	F0116
5	Duplicate 89-045	F0117
6	Spike 89-045	F0118
7		
8		
9		
10		
11		

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Interim	Primary Book	Second Book	Third Book	Final Volume of Standard
Standard Type	No. & Aliquot	No. & Aliquot	No. & Aliquot	
Spike	104C15B/5mL	103C15D/5mL		50 mL

Rev. E 4/04/90	Prepared by: <u>Shirley Cervantes</u> Signature	S.A. Cervantes Printed Name	Date: 07/13/90
SST-102	Verified by: <u>Cathy M Seidel</u> Signature	C.M. Seidel Printed Name	Date: 07/13/90
132	Approved by: <u>L.H. Taylor</u> Signature	L.H. Taylor Printed Name	Date: 08-30-90

## Analytical Batch

Lab Segment Serial No.: F0101

Customer ID.: 89-045

Instrument	Not reported.
Procedure / Rev	LA-505-151
Technologist	D. M. Southwick
Date	03-12-90
Temperature	70 F
Starting Time	12:30
Ending Time	15:11
Chemist	S. A. Jones

ICP Analysis of Sample 89-045.

Acid Digestion.

No inter-element corrections have been made on this data.

	Description	Lab. Id.
1	Initial LMCS Check Std	N/A
2	Reagent Blank	F0123
3	Digested LMCS Std	F0115
4	Sample 89-045	F0116
5	Duplicate 89-045	F0117
6	Spike 89-045	F0118
7	Sample 89-077	F0600
8	Duplicate 89-077	F0601
9	Sample 89-049	F0188
10	Duplicate 89-049	F0189
11	Final LMCS Check Std	N/A

	Description	Lab. Id.
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

Standard Type	Primary Book No. & Aliquot	Second Book No. & Aliquot	Third Book No. & Aliquot	Final Volume of Standard
Spike 89-045	103C15C/5.0mL	104CV15D/5.0mL		50.0 mL
LMCS Check Std	78C11C/1.0 mL	82B38A/1.0mL	77C11C/1.0 mL	11.0 mL
Acid Digested Std	81C11A/5.0mL	82C11A/5.0mL		50.0 mL

Interim

4/04/90

Prepared by: Shirley Ceravntes S. A. Ceravntes Date: 08-05-90  
Signature Printed Name

Verified by: Craig M Seidel C. M. Seidel Date: 08-05-90  
Signature Printed Name

SST-102

Approved by: L.H. Taylor L.H. Taylor Date: 08-30-90  
Signature Printed Name

## ICP Results      RAW DATA SUMMARY

Date Analyzed:	March 12, 1990	Acid Digested Standard	F0115
Procedure:	LA-505-151/A-0	Reagent Blank	F0123
Analyst:	D. M. Southwick	Sample 89-045	F0116
Digestion	Acid Digestion	Duplicate 89-045	F0117
Procedure:	LA-505-159/A-0	Spike of F0116	F0118
		Acid Digested Standard	F0119

	Instrument Starting LMCS Standard %	Acid Digest. LMCS Standard %	Reagent BLANK ppm	Wet Weight Sample ug/g	Wet Weight Sample ug/g	Spike Recovery %	LMCS ACID Digestion %	Closing LMCS Standard %
Aluminum	100.47%		-0.71 LT	45520	52568	NOT CALC.	95.63%	106.46%
Antimony	97.68%		-0.59 LT	-751 LT	10 LT			101.29%
Arsenic	104.97%		-0.14 LT	-172 LT	3 LT			103.90%
Barium	103.89%		-0.04 LT	-58 LT	-6 LT	NOT CALC.	87.43%	93.37%
Beryllium	100.50%		0.00 LT	-6 LT	-2 LT			105.22%
Bismuth	102.33%	83.52%	-0.66 LT	21590	26620	NOT CALC.		106.21%
Boron	99.53%	90.89%	-0.03 LT	-65 LT	14 LT	NOT CALC.		93.06%
Cadmium	96.13%	93.88%	-0.02 LT	-35 LT	-8 LT	NOT CALC.		98.54%
Calcium	105.75%	115.16%	0.19	295	1113	111.89%		95.68%
Cerium	93.76%		-1.74 LT	-2483 LT	-288 LT	NOT CALC.	83.71%	97.80%
Chromium	100.34%		-0.18 LT	119	188	98.46%	86.40%	97.42%
Cobalt	70.93%		-0.19 LT	-265 LT	-117 LT	NOT CALC.	72.36%	76.26%
Copper	106.21%	97.85%	-0.08 LT	-112 LT	2 LT	NOT CALC.		98.57%
Europium	99.63%		-0.03 LT	-49 LT	-18 LT			116.33%
Iron	101.13%		0.02	14837	16885	NOT CALC.	95.97%	97.66%
Lanthanum	99.52%	95.72%	-0.13 LT	-193 LT	-40 LT	NOT CALC.		108.92%
Lead	102.54%	91.72%	-0.35 LT	-83 LT	602	NOT CALC.		100.30%
Lithium	116.33%		-0.06 LT	-114 LT	-27 LT	NOT CALC.	105.88%	104.53%
Magnesium	100.98%	96.86%	0.02	479	449	87.17%		98.04%
Manganese	101.63%		0.00 LT	156	269	107.06%	94.31%	97.98%
Mercury	104.69%		-0.05 LT	-47 LT	9 LT			99.42%
Molybdenum	101.07%	95.80%	-0.04 LT	-59 LT	3 LT	NOT CALC.		101.72%
Neodymium	86.98%		-2.79 LT	-4417 LT	-3004 LT			85.68%
Nickel	99.60%		-0.07 LT	-104 LT	-5 LT	NOT CALC.	95.42%	99.59%
Phosphorus	105.42%	97.37%	-0.06 LT	53626	47284	NOT CALC.		99.90%
Potassium	102.78%		-2.53 LT	-3986 LT	-1211 LT	NOT CALC.		101.04%
Samarium	102.38%		-1.95 LT	-3036 LT	-675 LT			104.47%
Selenium	104.47%		-0.62 LT	-609 LT	99 LT			103.29%
Silicon	82.09%	86.92%	2.70	1369	3584	33.15%		84.14%
Silver	100.85%		-0.11 LT	-186 LT	-48 LT	NOT CALC.		98.23%
Sodium	108.57%	114.60%	-0.50 LT	187906	168975	NOT CALC.		101.63%
Strontium	103.85%	94.59%	-0.02 LT	120	106	95.52%		94.36%
Sulfur	102.58%		-0.05 LT	113	278			103.02%
Tantalum	95.45%		-0.28 LT	-417 LT	-73 LT	NOT CALC.	84.22%	93.36%
Thallium	101.40%		-2.02 LT	-3223 LT	-906 LT			104.50%
Thorium	97.10%		-1.11 LT	-1702 LT	-403 LT			104.62%
Tin	100.56%		-0.03 LT	-112 LT	-11 LT	NOT CALC.	101.08%	104.92%
Titanium	96.10%		0.03	-70 LT	8 LT	NOT CALC.	89.66%	101.92%
Tungsten	85.32%		-0.16 LT	-199 LT	30 LT			84.78%
Uranium	98.51%		-10.90 LT	-15293 LT	-2485 LT			99.44%
Vanadium	102.61%		-0.09 LT	-128 LT	-46 LT			105.49%
Zinc	98.20%	94.36%	0.00	142	5380	98.35%		98.06%
Zirconium	96.38%		-0.19 LT	-276 LT	97	NOT CALC.	86.94%	101.68%

LT: Less Than

NC: Not Calibrated

NOT CALC: Not Calculated

# Instrument Standards Outside Control Limits

## ICP Results

Date Analyzed: March 12, 1990  
 Procedure: LA-505-151/A-0  
 Analyst: D. M. Southwick  
 Digestion Procedure: Acid Digestion  
 LA-505-159/A-0

Acid Digested Standard  
 Reagent Blank  
 Sample 89-045  
 Duplicate 89-045  
 Spike of F0116  
 Acid Digested Standard

F0115  
 F0123  
 F0116  
 F0117  
 F0118  
 F0119

					Digestion Weight	0.01301	Digestion Weight
Starting LMCS Standard	Instrument Recovery Standard	LMCS Acid Digestion Standard	Acid Digestion Standard Recovery	Reagent Blank	Volume	Sample	Volume
ppm	%	ppm	%	ppm	ppm	ppm	ppm
Aluminum	50.24	100.47%		-0.71 LT	592.03	632.45	
Antimony	9.77	97.68%		-0.59 LT	-43.16	-9.77 LT	
Arsenic	52.48	104.97%		-0.14 LT	-9.29	-2.23 LT	
Barium	10.39	103.89%		-0.04 LT	-2.59	-0.75 LT	
Beryllium	10.05	100.50%		0.00 LT	-0.27	-0.08 LT	
Bismuth	51.27	102.33%	8.35	83.52%	-0.66 LT	280.80	315.97
Boron	9.95	99.53%	9.09	90.89%	-0.03 LT	-2.66	-0.84 LT
Gadolinium	9.61	96.13%	9.39	93.88%	-0.02 LT	-1.82	-0.46 LT
Calcium	10.58	105.75%	11.52	115.16%	0.19	3.84	2.90
Cerium	9.38	93.76%		-1.74 LT	-118.80	-32.30 LT	
Chromium	10.03	100.34%		-0.18 LT	-12.25	1.55	
Cobalt	7.09	70.93% #		-0.19 LT	-11.68	-3.45 LT	
Copper	10.62	106.21%	9.79	97.85%	-0.08 LT	-6.01	-1.45 LT
Europium	9.96	99.63%		-0.03 LT	-2.22	-0.64 LT	
Iron	10.11	101.13%		0.02	192.97	192.07	
Lanthanum	49.86	99.52%	9.57	95.72%	-0.13 LT	-8.41	-2.51 LT
Lead	51.37	102.54%	9.17	91.72%	-0.35 LT	-8.45	-1.08 LT
Lithium	11.63	116.33% #		-0.06 LT	-4.95	-1.48 LT	
Magnesium	10.10	100.98%	9.69	96.86%	0.02	6.23	2.68
Manganese	10.16	101.63%		0.00 LT	2.03	2.35	
Mercury	26.17	104.69%		-0.05 LT	-2.75	-0.61 LT	
Molybdenum	50.54	101.07%	9.56	95.80%	-0.04 LT	-3.13	-0.77 LT
Neodymium	8.70	86.98% #		-2.79 LT	-227.00	-57.45 LT	
Nickel	9.96	99.60%		-0.07 LT	-5.21	-1.35 LT	
Phosphorus	52.71	105.42%	9.74	97.37%	-0.06 LT	697.46	671.68
Potassium	25.70	102.78%		-2.53 LT	-182.60	-51.84 LT	
Samarium	10.24	102.38%		-1.95 LT	-135.00	-39.48 LT	
Selenium	52.24	104.47%		-0.62 LT	-41.70	-7.92 LT	
Silicon	41.04	82.09% #	8.69	86.92%	2.70	8.72	17.80
Silver	10.09	100.85%	9.42	-0.11 LT	-8.54	-2.42 LT	
Sodium	27.14	108.57%	11.46	114.60%	-0.50 LT	2443.90	2510.10
Strontium	10.39	103.85%	9.46	94.59%	-0.02 LT	0.65	1.56
Sulfur	51.29	102.58%		-0.05 LT	-1.78	1.48	
Tantalum	47.73	95.45%		-0.28 LT	-19.86	-5.43 LT	
Thallium	50.70	101.40%		-2.02 LT	-139.40	-41.92 LT	
Thorium	48.65	97.10%		-1.11 LT	-76.34	-22.13 LT	
Tin	50.28	100.56%		-0.03 LT	-6.21	-1.46 LT	
Titanium	48.05	96.10%		0.03	-3.59	-0.92 LT	
Tungsten	21.33	85.32% #		-0.16 LT	-11.93	-2.58 LT	
Uranium	49.35	98.51%		-10.90 LT	-730.60	-198.90 LT	
Vanadium	10.26	102.61%		-0.09 LT	-6.38	-1.67 LT	
Zinc	9.82	98.20%	9.44	94.36%	0.00	1.85	2.20
Zirconium	48.19	96.38%		-0.19 LT	-12.77	-3.59 LT	
Dilution Factor	1.00		10.00		1.00	101.00	21.00
							1.00

	0.01011		0.01007					
	0.5054 g	50.00 mL	Weight	0.5033 g	50.00 mL			
	Duplicate Sample Dilution	Duplicate Sample Dilution	Spike of Sample Dilution	Spike of Sample Dilution	Spike of Sample Dilution	Spike Recovery	Standard LMCS Acid	Acid Digestion Standard Ending LMCS Standard
	Two ppm	One ppm	Three ppm	Two ppm	One ppm	%	Digestion ppm	Recovery %
Aluminum	531.36	531.78		493.71	529.23	NOT CALC.	9.56	95.63%
Antimony	-21.58	0.10 LT		-51.80	2.49 LT			10.13
Arsenic	-4.38	0.03 LT		-13.46	-0.84 LT			51.95
Barium	-1.40	-0.07 LT		6.40	9.48	NOT CALC.	8.74	87.43%
Beryllium	-0.13	-0.02 LT		-0.49	-0.06 LT			9.34
Bismuth	269.08	282.99		234.64	287.99	NOT CALC.		10.52
Boron	-0.95	0.14 LT		9.85	10.37	NOT CALC.		53.21
Cadmium	-1.00	-0.08 LT		7.00	9.41	NOT CALC.		9.31
Calcium	11.25	2.57		14.16	13.14	111.89%		9.85
Cerium	-56.48	-2.91 LT		-195.10	-21.42 LT	NOT CALC.	8.37	83.71%
Chromium	-10.87	1.90		-5.57	11.04	98.46%	8.72	86.40%
Cobalt	-8.69	-1.19 LT		-13.31	4.52	NOT CALC.	7.24	72.36%
Copper	-2.80	0.02 LT		1.69	9.71	NOT CALC.		7.63
Europium	-1.11	-0.18 LT		-3.48	-0.42 LT			9.57
Iron	170.67	164.51		170.47	166.29	NOT CALC.	9.60	95.97%
Lanthanum	-4.46	-0.40 LT		-2.74	8.63	NOT CALC.		9.78
Lead	-5.85	6.08		-24.05	9.46	NOT CALC.		54.57
Lithium	-2.47	-0.27 LT		5.38	11.08	NOT CALC.	10.59	105.88%
Magnesium	4.54	2.63		13.54	12.73	87.17%		10.45
Manganese	2.72	2.56		12.28	12.10	107.06%	9.43	94.31%
Mercury	-0.78	0.10 LT		-1.14	0.69	NOT CALC.		9.80
Molybdenum	-1.52	0.03 LT		5.17	9.47	NOT CALC.		9.80
Neodymium	-165.20	-30.36 LT		-293.60	-36.64 LT			50.25
Nickel	-2.38	-0.05 LT		1.55	9.00	NOT CALC.	9.54	95.42%
Phosphorus	477.95	467.13		512.80	507.61	NOT CALC.		8.57
Potassium	-100.70	-12.24 LT		-272.20	-26.98 LT	NOT CALC.		49.95
Samarium	-65.86	-6.83 LT		-217.00	-27.39 LT			25.26
Selenium	-20.38	1.00 LT		-60.15	-1.63 LT			101.04%
Silicon	19.16	36.23		-12.79	17.09	33.15%		104.47%
Silver	-4.39	-0.48 LT		-12.17	-0.51 LT	NOT CALC.		51.64
Sodium	1708.00	1630.80		1925.20	1840.90	NOT CALC.		103.29%
Strontium	1.07	1.52		10.76	11.47	95.52%		52.07
Sulfur	1.19	2.81		-7.97	2.04			84.14%
Tantalum	-10.26	-0.74 LT		-30.01	-3.59 LT	NOT CALC.	8.38	84.22%
Thallium	-60.70	-9.15 LT		-212.60	-25.67 LT			46.68
Thorium	-37.29	-4.07 LT		-121.20	-14.94 LT			93.36%
Tin	-3.01	-0.11 LT		-0.87	8.29	NOT CALC.	10.11	52.25
Titanium	-1.64	0.08 LT		2.86	7.57	NOT CALC.	8.98	52.41
Tungsten	-5.34	0.30 LT		-18.45	-1.51 LT			104.50%
Uranium	-356.40	-25.12 LT		-1175.00	-135.20 LT			52.41
Vanadium	-2.10	-0.46 LT		-9.22	-0.95 LT	NOT CALC.		104.62%
Zinc	54.38	2.32		11.26	11.96	98.35%		52.46
Zirconium	-5.20	0.98		-20.23	-2.38 LT	NOT CALC.	8.68	104.92%
Dilution Factor	101.00	21.00	1.00	101.00	21.00		10.00	105.55
								105.49%
								9.81
								98.06%
								50.84
								101.68%

	Spike Standard LMCS	Spike Standard ID Book	LMCS Standards Values	LMCS Standard IDs Book	ACID DIGESTION LMCS STANDARD VALUES	ACID DIGEST. LMCS	Calc. Sample	Calc. Duplicate	Calc. Spike
	ppm added	#	ppm	#	ppm in Sample	#			
Aluminum	10.00	5.0 mL	50.00		100.00		592	531	494
Antimony		104C15D	10.00				-10	0	2
Arsenic			50.00				-2	0	-1
Barium	10.00		10.00		100.00		-1	0	6
Beryllium			10.00				0	0	0
Bismuth	10.00		50.10		100.00		281	269	235
Boron	10.00		10.00		100.00		-1	0	10
Cadmium	10.00		10.00		100.00		0	0	7
Calcium	10.00		10.00		100.00		4	11	14
Cerium	10.00		10.00		100.00		-32	-3	-21
Chromium	10.00		10.00		100.90		2	2	11
Cobalt #	10.00		10.00		100.00		-3	-1	5
Copper	10.00		10.00		100.00		-1	0	10
Europium #			10.00				-1	0	0
Iron	10.00		10.00		100.00		193	171	170
Lanthanum	10.00		50.10		100.00		-3	0	9
Lead	10.00		50.10		100.00		-1	6	9
Lithium	10.00		10.00		100.00		-1	0	5
Magnesium	10.00		10.00		100.00		6	5	14
Manganese	10.00		10.00		100.00		2	3	12
Mercury			25.00				-1	0	1
Molybdenum	10.00		50.00		99.80		-1	0	5
Neodymium #			10.00		100.00		-57	-30	-37
Nickel	10.00		10.00		100.00		-1	0	9
Phosphorus	10.00		50.00		100.00		697	478	513
Potassium	10.00		25.00		100.00		-52	-12	-27
Samarium			10.00		100.00		-39	-7	-27
Selenium			50.00				-8	1	-2
Silicon #	10.00		50.00		100.00		18	36	17
Silver	10.00		10.00				-2	0	-1
Sodium	10.00		25.00		100.00		2444	1708	1925
Strontium	10.00		10.00		100.00		2	1	11
Sulfur			50.00				1	3	2
Tantalum	9.95		50.00		99.50		-5	-1	-4
Thallium			50.00				-42	-9	-26
Thorium			50.10				-22	-4	-15
Tin	10.00		50.00		100.00		-1	0	8
Titanium	10.00		50.00		100.10		-1	0	8
Tungsten #			25.00				-3	0	-2
Uranium			50.10				-199	-25	-135
Vanadium			10.00				-2	0	-1
Zinc	10.00		10.00		100.00		2	54	11
Zirconium	9.98		50.00		99.80		-4	1	-2
Dilution Factor					10.00				

## ICP Calibration Report

7  
C  
C  
—  
C  
—  
6

Procedure: LA-505-151 Revision: A-0  
Instrument: WB39939  
Technologist: D. M. Southwick  
Date: March 12, 1990  
Time: 08:36

### Calibration Standards for ICP Program "SST"

Element	Standard	Element	Standard
Aluminum	SST-3	Antimony	SST-4
Arsenic	SST-4	Barium	SST-2
Beryllium	SST-2	Bismuth	SST-4
Boron	SST-3	Cadmium	SST-2
Calcium	SST-2	Cerium	SST-5
Chromium	SST-2	Cobalt	SST-2
Copper	SST-2	Europium	SST-5
Iron	SST-4	Lanthanum	SST-5
Lead	SST-4	Lithanum	SST-1
Magnesium	SST-2	Manganese	SST-2
Mercury	SST-2	Molybdenum	SST-3
Neodymium	SST-5	Nickel	SST-2
Phosphorous	SST-3	Potassium	SST-1
Samarium	SST-5	Selenium	SST-4
Silicon	SST-3	Silver	SST-2
Sodium	SST-1	Strontium	SST-2
Sulfur	SST-3	Tantalum	SST-3
Thallium	SST-4	Thorium	SST-4
Tin	SST-4	Titanium	SST-3
Tungsten	SST-3	Uranium	SST-4
Vanadium	SST-2	Zinc	SST-2
Zirconium	SST-3		

### ICP Standard Formulations

SST-0:  
Calibration blank, 1 M ultrx HNO<sub>3</sub>.

12-Mar-90 08:36:20

Condition	Value	Min	/	Max
VACUUM	= 14.52	7.000	/	50.00
SPTEMP	= 38.75	37.00	/	39.00
MAINS	= 234.8	220.0	/	247.0
-1000V	= -1004	-1010	/	-990
CTEMP	= 23.25	19.00	/	35.00
+5V	= 5.150	4.750	/	5.250
+12V	= 12.14	11.70	/	12.30
-12V	= -12.2	-12.3	/	-11.7
+24V	= 23.12	22.50	/	26.50
-100V	= -100	-101	/	-99.0
+5VSQ	= 5.146	4.750	/	5.250
+15VSQ	= 15.13	14.70	/	15.30
-15VSQ	= -15.2	-15.3	/	-14.7

Position Calibration in Progress

SLIT	PM	ALPHA	BETA	ALPHA	BETA	ALPHA	BETA
POS'N		SLIT	SLIT	LAMBDA1	LAMBDA1	LAMBDA2	LAMBDA2

Previous data :

INSTR 0.00000 586.625 1.00093 -0.2905 1.00009 -0.0682 0.00000 0.00000

Current data :

INSTR 0.00000 586.663 1.00092 -0.2222 1.00009 -0.0677 0.00000 0.00000

START THE PLASMA NOW, PLEASE. 12-Mar-90 08:44:12




Figure 1 ICP Calibration March 12, 1990




Figure 2ICP Calibration March 12, 1990

```
Profile position = -1.03
Option :
1. Accept new position
2. Re-run sample
3. Enter new profile position via keyboard
<CR> Exit

[AAAAAA] Input not numeric - Please retype

Option = 1
SAMI positioning on profile
>

Profile position = -1.03
Option :
1. Accept new position
2. Re-run sample
3. Enter new profile position via keyboard
<CR> Exit

[AAAAAA] Input not numeric - Please retype

Option = 1
SAMI positioning on profile
>
```

Figure 3 ICP Analysis March 12, 1990

Sample name : SST0  
Programme : SST 12-Mar-90 10:59:13

NAME	MV	INT	RSD
AL	1.74	2.00	
SB	0.40	1.48	
AS	1.00	2.36	
BA	3.12	1.99	
BE	0.68	1.30	
BI	3.72	1.54	
B	3.92	1.28	
CD	2.14	1.44	
CA	0.51	0.56	
CE	4.33	2.05	
CR	1.62	35.05	
CO	0.23	0.91	
CU	2.49	1.79	
EU	3.36	2.14	
FE	1.48	0.98	
LA	0.35	0.91	
PB	0.26	1.09	
LI	3.24	1.96	
HG	0.40	1.03	
MN	0.73	0.97	
HG	3.46	1.54	
MO	1.67	2.20	
ND	4.85	10.38	
NI	3.28	2.43	
P	1.35	1.22	
K	2.81	1.34	
SM	4.15	2.08	
SE	1.64	1.71	
SI	2.79	1.97	
AG	11.90	2.09	
NA	4.19	2.27	
SR	3.11	1.63	
S	0.71	1.85	
IA	3.20	1.76	
TL	3.68	1.52	
TH	1.29	1.78	
SN	1.24	2.94	
TI	3.09	1.74	
W	1.31	2.24	
U	4.11	2.07	
V1	3.81	2.05	
ZN	2.51	1.34	
ZR	4.04	1.53	

Sample name : SST1  
Programme : SST 12-Mar-90 11:03:33

NAME	MV	INT	RSD
LI	477.08	0.73	
K	13.99	0.55	

Figure 4 ICP Calibration March 12, 1990

Sample name : SST1  
Programme : SST 12-Mar-90 11:03:33

NAME	MV	INT	RSD
LI	477.08	0.73	
K	13.99	0.55	
NA	60.88	0.70	

Sample name : SST2  
Programme : SST 12-Mar-90 11:05:42

NAME	MV	INT	RSD
BA	348.09	1.28	
BE	554.46	0.91	
CD	333.58	0.45	
CA	509.69	1.11	
CR	70.72	1.60	
CO	2.71	1.59	
CU	112.27	1.45	
EE	158.44	0.85	
HG	493.87	0.62	
MN	325.52	0.50	
NI	178.09	0.47	
AG	496.65	0.87	
SR	624.33	1.23	
V1	216.30	0.54	
ZN	611.79	0.37	

Sample name : SST3  
Programme : SST 12-Mar-90 11:10:28

NAME	MV	INT	RSD
AL	24.63	0.18	
B	771.85	0.51	
HG	671.49	0.18	
MO	317.71	0.17	
P	90.17	0.98	
SI	91.60	0.23	
S	48.86	0.81	
TA	131.98	0.22	
TI	567.83	0.24	
W	72.55	0.34	
ZR	185.89	0.21	

Sample name : SST4  
Programme : SST 12-Mar-90 11:15:15

NAME	MV	INT	RSD
SB	7.42	0.79	
AS	133.32	0.54	
BI	113.47	0.43	
PB	5.44	0.99	
SE	52.01	0.11	
TL	52.32	0.22	

Figure 5 ICP Analysis March 12, 1990

NA 60.88 0.70

Sample name : SST2  
Programme : SST

12-Mar-90 11:05:42

NAME	MV	INT	RSD
BA	348.09	1.28	
BE	554.46	0.91	
CD	333.58	0.45	
CA	509.69	1.11	
CR	70.72	1.60	
CO	2.71	1.59	
CU	112.27	1.45	
FE	158.44	0.85	
HG	493.87	0.62	
NN	325.52	0.50	
NI	178.09	0.47	
AG	496.65	0.87	
SR	624.33	1.23	
VI	216.30	0.54	
ZN	611.79	0.37	

Sample name : SST3  
Programme : SST

12-Mar-90 11:10:28

NAME	MV	INT	RSD
AL	24.63	0.18	
B	771.85	0.51	
HG	671.49	0.18	
MO	317.71	0.17	
P	90.17	0.98	
SI	91.60	0.23	
S	48.86	0.81	
TA	131.98	0.22	
TI	567.03	0.24	
W	72.55	0.34	
ZR	185.89	0.21	

Sample name : SST4  
Programme : SST

12-Mar-90 11:15:15

NAME	MV	INT	RSD
SB	7.42	0.79	
AS	133.32	0.54	
BI	113.47	0.43	
PB	3.44	0.99	
SE	52.01	0.11	
TL	52.32	0.22	
TH	26.03	0.54	
SN	238.24	0.76	
U	13.00	0.28	

Sample name : SST5  
Programme : SST

12-Mar-90 11:18:24

Figure 6 ICP Calibration March 12, 1990

NAME	MV INT	RSD
CE	15.49	0.31
EU	513.43	0.37
LA	7.00	0.09
ND	15.22	0.68
SM	13.93	0.05

Programme name : SST      Channel name : AL      Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.6521	25.862	-0.379832E+01	0.216420E+01		

Programme name : SST      Channel name : SB1      Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.3768	7.7872	-0.565079E+01	0.142457E+02		

Programme name : SST      Channel name : AB      Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.9510	139.99	-0.756486E+00	0.755730E+00		

Programme name : SST      Channel name : BA      Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.9640	365.49	-0.180888E+00	0.579768E-01		

Programme name : SST      Channel name : BE1      Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.6489	582.18	-0.246672E-01	0.361159E-01		

Figure 7 ICP Calibration March 12, 1990

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	3.5346	119.15	-0.339009E+01	0.911151E+00		
 Programme name : SST Channel name : B Polynomial type : CC						
Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	3.7243	810.44	-0.255254E+00	0.651102E-01		
 Programme name : SST Channel name : CD Polynomial type : CC						
Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.0286	350.25	-0.128852E+00	0.603427E-01		
 Programme name : SST Channel name : CA Polynomial type : CC						
Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	0.4861	535.17	-0.200977E-01	0.392789E-01		
 Programme name : SST Channel name : CE Polynomial type : CC						
Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	4.1122	16.261	-0.775909E+01	0.179249E+01		
 Programme name : SST Channel name : CR Polynomial type : CC						
Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.5381	74.252	-0.468614E+00	0.289447E+00		

Figure 8 ICP Calibration March 12, 1990

Programme name : SST Channel name : CU Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	0.2163	2.8413	-0.183726E+01	0.806994E+01		

Programme name : SST Channel name : CU Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	2.3608	117.88	-0.452699E+00	0.182173E+00		

Programme name : SST Channel name : EU Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	3.2116	539.10	-0.132562E+00	0.392118E-01		

Programme name : SST Channel name : FE Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	1.4082	166.36	-0.188884E+00	0.127423E+00		

Programme name : SST Channel name : LA Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	0.3338	7.3493	-0.105696E+01	0.300842E+01		

Programme name : SST Channel name : PB Polynomial type : CC

Curve	Min Int	Max Int		Curve Coefficients		
	C0		C1	C2		C3
CRV1	0.2503	5.7144	-0.509140E+01	0.193100E+02		

Figure 9 ICP Calibration March 12, 1990

Programme name : SST	Channel name : LI	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	3.0786 500.93 -0.341959E+00 0.103521E+00		
Programme name : SST	Channel name : HG	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	0.3841 518.56 -0.163876E-01 0.403298E-01		
Programme name : SST	Channel name : MN	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	0.6919 341.80 -0.448494E-01 0.615781E-01		
Programme name : SST	Channel name : HG	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	3.2899 705.07 -0.259196E+00 0.748472E-01		
Programme name : SST	Channel name : MO	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	1.3881 333.59 -0.264476E+00 0.158211E+00		
Programme name : SST	Channel name : ND	Polynomial type : CC	
Curve	Min Int Max Int	Curve Coefficients	
	C0 C1 C2 C3		
CRV1	4.6046 15.979 -0.934691E+01 0.192839E+01		

Figure 10ICP Calibration March 12, 1990

Programme name : SST Channel name : NI Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	3.1122	186.99	-0.374799E+00	0.114408E+00		

Programme name : SST Channel name : P Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.2796	94.683	-0.758210E+00	0.562888E+00		

Programme name : SST Channel name : K Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.6701	14.691	-0.125689E+02	0.447187E+01		

Programme name : SST Channel name : SM Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	3.9384	14.627	-0.847380E+01	0.204401E+01		

Programme name : SST Channel name : SE Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	1.5342	54.614	-0.324747E+01	0.198501E+01		

Programme name : SST Channel name : SI Polynomial type : CC

Curve	Min Int	Max Int	C0	C1	C2	C3
CRV1	2.6521	96.182	-0.157171E+01	0.563000E+00		

Figure 11 ICP Calibration March 12, 1990

Programme name : SST Channel name : AG Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	11.302	521.48	-0.490852E+00	0.412585E-01				

Programme name : SST Channel name : NA Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	3.9786	63.921	-0.369384E+01	0.882006E+00				

Programme name : SST Channel name : SR Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	2.9551	655.55	-0.100146E+00	0.321945E-01				

Programme name : SST Channel name : S Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	0.6793	51.305	-0.742512E+00	0.103848E+01				

Programme name : SST Channel name : TA Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	3.0416	138.58	-0.124308E+01	0.388260E+00				

Programme name : SST Channel name : IL2 Polynomial type : CC

Curve	Min	Int	Max	Int	Coefficients			
					C0	C1	C2	C3
CRV1	3.4969	54.936	-0.756795E+01	0.205595E+01				

Figure 12 ICP Calibration March 12, 1990

Programme name : SST Channel name : TH Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	1.2223	27.328	-0.520068E+01	0.404198E+01		

Programme name : SST Channel name : SN Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	1.1745	250.15	-0.521655E+00	0.421937E+00		

Programme name : SST Channel name : TI Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	2.9323	596.23	-0.273279E+00	0.885352E-01		

Programme name : SST Channel name : W Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	1.2458	76.178	-0.920380E+00	0.701866E+00		

Programme name : SST Channel name : U Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	3.9020	13.646	-0.462069E+02	0.112499E+02		

Programme name : SST Channel name : V1 Polynomial type : CC

Curve	Min	Int	Max	Int	Curve Coefficients	
	C0	C1	C2	C3		
CRV1	3.6208	227.12	-0.358727E+00	0.941210E-01		

Figure 13 ICP Calibration March 12, 1990

Programme name : SST Channel name : ZN Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients		
	C0	C1	C2	C3	
CRV1	2.3986	642.38	-0.825355E-01	0.328260E-01	

Programme name : SST Channel name : ZR Polynomial type : CC

Curve	Min Int	Max Int	Curve Coefficients		
	C0	C1	C2	C3	
CRV1	3.8367	193.19	-0.111041E+01	0.274945E+00	

Figure 14 ICP Calibration March 12, 1990

Sample name : BIC11A  
 Sample code 1 : ICP-1B  
 Sample code 2 : MULTI  
 Sample code 3 : 1-10  
 Programme : SSI                    12-Mar-90 12:39:20

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.57	(-0.369	(-4.060	-19.95	
Sb	0.38	-0.280	-3.082	-18.33	
As	1.27	0.203	2.231	7.97	
Se	2.88	(-0.014	(-0.151	-29.52	
Be	0.63	(-0.002	(-0.023	-23.73	
Bi	13.79	9.170	100.87	0.41	
B	137.40	8.691	95.603	0.97	
Cd	150.36	8.944	98.386	0.51	
Ca	231.19	9.061	99.669	1.66	
Ce	3.96	(-0.655	(-7.203	-25.86	
Cr	1.17	(-0.131	(-1.441	-3.38	
Co	0.22	-0.089	-0.976	-10.50	
Cu	52.86	9.176	100.94	1.14	
Eu	3.06	(-0.013	(-0.139	-23.58	
Fe	1.46	-0.002	-0.026	-179.79	
La	3.40	9.186	101.04	1.29	
Pb	0.74	9.166	100.82	1.20	
Li	2.92	(-0.034	(-0.371	-21.54	
Hg	221.00	8.941	98.346	0.74	
Mn	0.73	0.000	0.001	793.68	
Hg	3.35	-0.009	-0.094	-46.99	
Mo	1.56	(-0.017	(-0.190	-25.19	
Nd	4.24	(-1.176	(-12.93	-10.01	
Ni	3.06	(-0.025	(-0.270	-42.23	
P	18.23	9.500	104.50	0.67	
K	4.67	8.325	91.576	1.61	
Sm	3.72	(-0.878	(-9.661	-22.03	
Se	1.50	(-0.261	(-2.868	-21.01	
Si	2.58	(-0.121	(-1.329	-26.18	
Ag	231.60	9.065	99.711	0.73	
Na	14.56	9.149	100.64	0.47	
Sr	280.11	8.918	98.097	1.54	
S	0.83	0.122	1.344	12.90	
Ts	2.90	(-0.119	(-1.304	-17.38	
Tl	3.32	(-0.746	(-8.209	-6.20	
Th	1.17	(-0.476	(-5.232	-21.45	
Sn	1.44	0.086	0.948	12.42	
Ti	2.81	(-0.025	(-0.270	-21.65	
W	1.40	0.062	0.682	13.22	
U	3.75	(-4.031	(-44.34	-25.01	
V	3.45	(-0.034	(-0.378	-23.89	
Zn	274.68	8.934	98.274	0.51	
Zr	3.73	(-0.084	(-0.928	-21.73	

Dilution factor : 11.0000

*enter*

Figure 1 ICP Analysis March 12, 1990

Sample name : S1C11A  
 Sample code 1 : ICP-1B  
 Sample code 2 : MO-SI  
 Sample code 3 : 1-10  
 Programme : SST                    12-Mar-90 12:44:01

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.80	0.125	1.370	61.55	
Sb	0.39	-0.043	-0.470	-258.91	
As	1.05	0.036	0.391	52.68	
Ba	3.10	-0.001	-0.013	-311.10	
Bm	0.68	-0.000	-0.002	-225.08	
Bi	3.70	-0.019	-0.214	-297.06	
B	4.24	0.021	0.228	18.41	
Cd	2.09	-0.003	-0.028	-98.50	
Ca	0.72	0.008	0.091	1.92	
Cr	4.31	-0.034	-0.375	-501.37	
Cr	1.17	(-0.129	(-1.415	-8.68	
Co	0.23	-0.016	-0.178	-76.38	
Cu	2.55	0.012	0.130	54.37	
Bu	3.36	-0.001	-0.009	-362.90	
Fs	1.54	0.007	0.077	89.19	
La	0.35	-0.004	-0.044	-270.43	
Pb	0.27	0.071	0.779	36.77	
Li	3.23	-0.001	-0.014	-520.10	
Mg	0.44	0.001	0.016	15.18	
Mn	0.74	0.001	0.010	70.23	
Hg	3.50	0.003	0.032	186.54	
Mo	58.51	8.992	98.916	1.59	
Nd	4.26	(-1.129	(-12.42	-11.11	
Hi	3.28	0.000	0.005	1589.85	
P	1.38	0.016	0.173	52.61	
K	2.75	-0.258	-2.837	-87.53	
Sm	4.13	-0.036	-0.397	-486.91	
Se	1.63	-0.003	-0.036	-1317.7	
Si	15.46	7.135	76.480	1.16	
Ag	11.92	0.001	0.011	556.41	
Na	4.20	0.010	0.107	771.49	
Er	3.11	-0.000	-0.001	-1726.2	
S	0.76	0.050	0.548	20.83	
Ta	3.17	-0.012	-0.131	-216.62	
Tl	3.71	0.064	0.709	288.37	
Th	1.29	0.027	0.296	334.50	
Sn	1.30	0.025	0.277	24.69	
Ti	3.07	-0.001	-0.015	-377.56	
W	1.31	0.000	0.005	5338.04	
U	4.08	-0.266	-2.929	-363.34	
V	3.87	0.006	0.061	25.24	
Zn	2.57	0.002	0.021	41.98	
Zr	4.02	-0.005	-0.059	-309.41	

Dilution factor : 11.0000

Sample name : S2C11A  
 Sample code 1 : ICP-2B  
 Sample code 2 : 1-10  
 Programme : SST                    12-Mar-90 12:48:22

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	5.85	8.974	98.716	0.60	
Sb	0.38	-0.223	-2.455	-35.54	

enter

Figure 2 ICP Analysis March 12, 1990

Sample name : 82C11A  
 Sample code 1 : ICP-2B  
 Sample code 2 : 1-10  
 Programme : SSX                    12-Mar-90 12:48:22

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	5.85	8.974	98.716	0.60	
Br	0.38	-0.223	-2.455	-35.54	
As	1.27	0.202	2.222	7.33	
Se	150.72	9.021	99.235	1.81	
Be	0.70	0.001	0.008	64.38	
Bi	3.32	(-0.186	(-2.041	-39.81	
B	4.10	0.012	0.128	60.83	
Cd	2.05	-0.005	-0.056	-23.31	
Ca	4.76	0.167	1.837	1.76	
Ce	9.25	8.830	97.128	1.23	
Cr	32.77	9.016	99.175	0.90	
Co	1.15	7.446	81.905	3.43	
Cu	2.31	(-0.033	(-0.360	-32.94	
Eu	3.66	0.011	0.121	24.36	
Fm	74.02	9.243	101.67	0.67	
La	0.36	0.015	0.165	70.24	
Pb	0.26	0.019	0.212	288.67	
Li	100.19	10.230	112.53	2.26	
Mg	2.90	0.101	1.111	0.64	
Mn	150.49	9.222	101.44	0.91	
Hg	3.50	0.010	0.111	57.16	
Mo	1.72	0.008	0.085	92.13	
Nd	8.46	6.970	76.675	1.58	
Ni	83.98	9.233	101.56	0.74	
P	1.31	-0.022	-0.241	-119.75	
K	2.50	(-1.297	(-15.36	-18.41	
Sm	3.71	(-0.892	(-9.811	-23.53	
Se	3.21	3.130	34.427	1.13	
Si	3.43	0.358	3.941	8.72	
Ag	11.85	-0.002	-0.020	-731.86	
Na	3.84	(-0.305	(-2.350	-31.86	
Br	2.93	(-0.006	(-0.063	-34.55	
S	0.82	0.109	1.196	19.52	
Ta	25.66	8.718	95.902	1.36	
Tl	3.44	(-0.491	(-5.405	-62.95	
Th	1.24	-0.191	-2.105	-60.66	
Sn	24.00	9.604	105.65	1.27	
Ti	101.14	8.681	95.490	1.34	
W	1.24	(-0.049	(-0.540	-49.12	
U	4.16	0.540	5.940	176.85	
V	3.37	(-0.042	(-0.460	-19.74	
Zn	3.50	0.032	0.355	3.27	
Zr	36.09	8.812	96.934	1.31	

Dilution factor : 11.0000

Figure 3 ICP Analysis March 12, 1990

Sample name : 78C11A  
Sample code 1 : SST-1  
Sample code 2 : DIRECT  
Programme : SST                    12-Mar-90 12:57:18

NAME	MV	INI	CONCEN	RSD
A1	1.53	(-0.461	-15.36	
Sb	1.08	9.768	1.22	
As	1.01	0.004	320.00	
Ba	182.31	10.389	1.85	
Be	0.61	(-0.003	-27.05	
Bi	3.18	(-0.489	-16.79	
B	156.79	9.953	1.27	
Cd	161.44	9.613	0.65	
Ca	269.73	10.575	1.53	
Cr	9.56	9.376	1.24	
Cr	36.29	10.034	0.50	
Co	1.11	7.093	3.41	
Cu	60.79	10.621	1.48	
Eu	3.50	0.005	48.13	
Fe	80.84	10.113	0.36	
La	0.34	-0.028	-21.43	
Pb	0.25	-0.212	-13.89	
Li	113.48	11.633	1.88	
Mg	249.56	10.098	0.66	
Mn	163.76	10.163	0.67	
Hg	3.26	(-0.015	-22.00	
Mo	1.54	(-0.020	-15.90	
Nd	9.36	8.698	4.10	
Ni	90.33	9.960	0.54	
P	1.29	-0.034	-31.51	
K	8.56	25.695	0.73	
Sm	3.35	(-1.637	-11.90	
Sm	3.40	3.504	0.83	
Bi	2.39	(-0.225	-13.44	
Ag	9.66	(-0.092	-11.49	
Na	34.96	27.143	1.35	
Sr	323.68	10.385	1.51	
B	0.90	0.189	7.27	
Ia	2.69	(-0.200	-16.78	
II	3.07	(-1.263	-20.12	
Th	1.13	(-0.632	-16.20	
Sn	120.40	50.281	0.99	
Ti	2.39	(-0.044	-11.48	
W	1.36	0.031	62.67	
U	3.72	(-4.395	-23.18	
V	3.17	(-0.061	-7.41	
Zn	301.67	9.820	0.61	
Zr	3.54	(-0.138	-10.76	

XAT  
2/21

Figure 4 ICP Analysis March 12, 1990

Sample name : 81B38A *should be 82B38A*  
 Sample code 1 : SST-2  
 Sample code 2 : DIRECT  
 Programme : SST                    12-Mar-90 13:02:00

KAT 36/21

NAME	MV	INT	CUNCEN	RSD
A1	3.53	3.920	1.95	
Bb	0.42	0.313	25.85	
An	2.87	1.410	1.71	
Ba	3.23	0.006	53.81	
Be	0.67	-0.000	-100.00	
Bi	59.99	51.265	1.04	
B	4.69	0.050	7.98	
Cd	2.15	0.001	137.94	
Ca	1.17	0.036	1.52	
Cr	4.43	0.182	78.46	
Cr	1.45	(-0.050	-14.86	
Co	0.22	-0.062	-13.04	
Cu	3.71	0.223	2.78	
Iu	269.99	10.454	1.22	
Fe	1.76	0.036	7.21	
La	16.93	49.061	1.28	
Pb	2.92	51.371	1.18	
Li	3.10	-0.015	-42.15	
Hg	0.60	0.008	2.89	
Mn	0.87	0.009	7.74	
Hg	3.91	0.033	15.63	
No	1.73	0.009	51.17	
Nd	4.50	(-0.471	-47.24	
Mi	3.35	0.009	73.37	
P	1.66	0.174	8.83	
K	2.65	(-0.733	-24.70	
Gm	9.15	10.238	1.26	
Gm	1.71	0.137	48.95	
Si	3.65	0.483	5.62	
Ag	256.34	10.085	0.85	
Na	4.09	-0.084	-87.43	
Sr	3.20	0.003	44.45	
B	0.84	0.130	7.87	
Ia	3.49	0.110	23.00	
II	4.33	5.456	3.14	
Ih	13.32	48.649	1.00	
Sn	1.41	0.071	5.03	
II	3.50	0.037	12.36	
W	1.29	-0.014	-76.33	
U	6.49	49.353	1.39	
V	6.06	0.212	2.50	
Zn	2.70	0.006	17.17	
Zr	4.31	0.075	17.68	

Figure 5 ICP Analysis March 12, 1990

Sample name :: 77C11C  
 Sample code 1 : SST-3  
 Sample code 2 : DIRECT  
 Programme : SST

12-Mar-90 13:07:56

NAME	MV	INT	CONCEN	RSD
Al	24.74	50.235	1.41✓	
Sb	0.46	0.878	8.93	
As	70.45	52.483	0.91✓	
Ba	3.03	-0.005	-111.90	
Be	278.97	10.050	1.06✓	
Bi	4.29	0.514	32.36	
B	4.36	0.029	29.22	
Cd	2.16	0.001	292.09	
Ca	0.81	0.012	1.22	
Ce	3.86	(-0.836	-29.08	
Cr	1.20	(-0.120	-10.47	
Co	0.23	-0.019	-89.21	
Cu	2.55	0.012	112.49	
Bu	3.01	(-0.014	-33.68	
Fe	1.72	0.031	17.18	
La	0.33	(-0.058	-31.03	
Pb	0.37	0.129	52.68	
Li	2.95	(-0.030	-40.27	
Hg	0.50	0.004	8.08	
Mn	0.95	0.014	6.59	
Hg	353.14	26.172	1.21	
Mo	321.10	50.537✓	0.86	
Nd	3.93	(-1.772	-12.11	
Ni	7.11	0.439	2.42	
P	94.99	)52.712✓	0.33	
K	2.50	(-1.389	-24.44	
Sm	3.73	(-0.856	-32.24	
Se	27.95	52.236✓	0.69	
Si	75.69	41.044✓	1.25	
Ag	18.15	0.258	8.05	
Na	4.02	-0.147	-81.63	
Sr	2.93	(-0.006	-45.29	
S	50.10	51.288✓	1.35	
Ta	126.12	47.723✓	0.63	
Tl	28.34	50.702✓	1.03	
Th	1.34	0.209	74.97	
Sn	1.67	0.184	16.00	
Ti	545.83	48.051✓	1.29	
W	31.70	21.330✓	1.59	
U	4.78	7.541	22.62	
V	112.83	10.261✓	2.34	
Zn	3.45	0.031	9.14	
Zr	179.31	48.190✓	1.20	

Figure 6 ICP Analysis March 12, 1990

Sample name : F123  
Sample code 1 : REAGEN  
Sample code 2 : DIRECT  
Programme : SST

12-Mar-90 13:14:26

NAME	MV	INT	CONCEN	RSD
Al	1.41	(-0.712	-19.22	
Sb	0.36	(-0.594	-11.00	
As	0.81	(-0.142	-27.54	
Ba	2.44	(-0.040	-14.87	
Be	0.59	(-0.004	-24.37	
Bi	3.00	(-0.659	-17.79	
B	3.53	(-0.025	-12.59	
Cd	1.77	(-0.022	-7.66	
Ca	5.23	0.186	86.27	
Ce	3.36	(-1.736	-20.65	
Cr	1.01	(-0.176	-12.77	
Co	0.20	(-0.191	-7.32	
Cu	2.06	(-0.078	-29.11	
Eu	2.55	(-0.032	-14.81	
Fe	1.61	0.017	199.07	
La	0.31	(-0.132	-14.61	
Pb	0.25	(-0.348	-12.83	
Li	2.63	(-0.063	-19.62	
Mg	0.84	0.018	90.26	
Mn	0.69	(-0.003	-84.08	
Hg	2.85	(-0.046	-22.13	
Mo	1.41	(-0.042	-21.87	
Nd	3.40	(-2.791	-10.24	
Ni	2.70	(-0.066	-9.84	
P	1.24	(-0.063	-159.67	
K	2.24	(-2.531	-18.14	
Sm	3.19	(-1.952	-16.81	
Se	1.32	(-0.621	-18.96	
Si	7.58	2.697	18.29	
Ag	9.30	(-0.107	-19.31	
Ns	3.63	(-0.496	-46.42	
Br	2.59	(-0.017	-18.79	
B	0.67	(-0.045	-110.12	
Tz	2.47	(-0.284	-9.47	
Tl	2.70	(-2.024	-13.62	
Ih	1.01	(-1.114	-12.13	
Sn	1.16	(-0.033	-224.37	
Ti	3.37	0.025	256.18	
W	1.09	(-0.155	-25.36	
U	3.14	(-10.90	-14.96	
V	2.90	(-0.085	-15.27	
Zn	2.63	0.004	281.98	
Zr	3.36	(-0.186	-14.61	

Figure 7 ICP Analysis March 12, 1990

Sample name : 34C11C  
 Sample code 1 : ICPI  
 Sample code 2 : 5-25  
 Programme : SST

12-Mar-90 13:20:48

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	1.70	-0.095	-0.477	-92.75	
Sb	0.38	-0.247	-1.235	-20.26	
As	1.01	0.004	0.019	541.47	
Ba	2.86	(-0.015	(-0.075	-32.37	
Bm	0.62	(-0.002	(-0.011	-11.04	
Bi	5.53	1.670	8.352	6.91	
B	31.84	1.818	9.089	1.39	
Cd	33.25	1.878	9.388	1.50	
Ca	59.15	2.303	11.516	0.19	
Ce	3.85	(-0.858	(-4.290	-25.25	
Cr	1.10	(-0.152	(-0.758	-4.58	
Co	0.22	-0.086	-0.430	-9.38	
Cu	13.23	1.957	9.785	1.09	
Bu	2.97	(-0.016	(-0.080	-23.01	
Pm	1.92	0.055	0.277	5.45	
La	0.99	1.914	9.572	1.35	
Pb	0.36	1.834	9.172	5.99	
Li	2.97	(-0.028	(-0.142	-28.39	
Hg	48.20	1.937	9.686	1.24	
Mn	0.86	0.008	0.040	13.41	
Hg	3.26	(-0.015	(-0.077	-20.97	
Mo	13.76	1.912	9.561	1.98	
Nd	3.92	(-1.795	(-8.973	-13.75	
W	3.06	(-0.024	(-0.121	-42.52	
P	4.81	1.947	9.737	1.94	
K	2.96	0.660	3.302	45.97	
Sm	3.67	(-0.975	(-4.875	-23.97	
Ba	1.49	(-0.295	(-1.476	-24.96	
Si	5.88	1.738	8.692	3.54	
Ag	57.58	1.885	9.424	1.56	
Na	6.79	2.292	11.460	5.31	
Sr	61.87	1.892	9.459	0.38	
S	0.83	0.117	0.583	12.66	
Ta	2.85	(-0.136	(-0.678	-21.37	
Il	3.25	(-0.879	(-4.396	-12.76	
Th	1.16	(-0.524	(-2.621	-22.72	
Sn	1.30	0.027	0.136	15.57	
Ti	3.45	0.032	0.161	17.79	
W	1.25	-0.044	-0.320	-41.65	
U	3.64	(-3.272	(-26.36	-23.93	
V	3.43	(-0.036	(-0.181	-18.65	
Zn	60.00	1.887	9.436	1.85	
Zr	3.70	(-0.092	(-0.461	-21.75	

Dilution factor : 5.00000

Figure 8 ICP Analysis March 12, 1990

Sample name : 35C11C  
 Sample code 1 : ICP2  
 Sample code 2 : 5-25  
 Programme : 88T      12-Mar-90 13:25:09

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	2.61	1.913	9.563	2.24	
Sb	0.39	-0.037	-0.285	-125.83	
As	1.06	0.047	0.236	13.64	
Ba	33.28	1.749	8.743	0.53	
Be	0.68	0.000	0.000	800.03	
Bi	3.69	-0.023	-0.117	-174.06	
B	4.56	0.048	0.240	7.30	
Cd	2.14	0.000	0.001	497.81	
Ca	11.04	0.414	2.068	5.16	
Cr	5.26	1.674	8.371	5.41	
Cr	7.64	1.744	8.718	0.58	
Co	0.41	1.447	7.236	1.48	
Cu	2.51	0.005	0.026	80.80	
Eu	3.39	0.000	0.002	462.42	
Fe	16.55	1.919	9.597	0.77	
La	0.35	-0.002	-0.010	-86.61	
Pb	0.28	0.225	1.126	21.57	
Li	23.31	2.118	10.588	1.33	
Mg	2.19	0.072	0.362	1.33	
Mn	31.36	1.886	9.431	0.59	
Hg	3.45	-0.001	-0.006	-100.66	
Mo	1.70	0.005	0.024	80.58	
Nd	5.11	0.514	2.568	8.33	
Ni	19.96	1.908	9.542	1.37	
P	1.49	0.078	0.388	27.18	
K	2.70	-0.507	-2.534	-25.61	
Sm	4.02	-0.253	-1.264	-41.79	
Se	1.94	0.603	3.027	7.01	
Si	3.56	0.433	2.163	4.62	
Ag	11.80	-0.004	-0.021	-142.27	
Na	4.80	0.543	2.717	7.63	
Br	3.16	0.001	0.007	74.46	
S	0.85	0.135	0.677	14.88	
Ta	7.52	1.676	8.380	2.05	
Tl	3.55	-0.273	-1.364	-42.68	
Ih	1.27	-0.084	-0.418	-63.46	
Sn	6.03	2.022	10.108	1.47	
Ti	23.36	1.795	8.975	0.66	
W	1.29	-0.013	-0.066	-46.19	
U	4.07	-0.446	-2.231	-121.44	
V	3.65	-0.015	-0.075	-15.94	
Zn	3.56	0.034	0.171	2.11	
Zr	10.35	1.735	8.677	1.15	

Dilution factor : 5.00000

Figure 9 ICP Analysis March 12, 1990

Sample name : E116  
 Sample code 1 : SAMPLE  
 Sample code 2 : 100-10  
 Programme : SST , 12-Mar-90 13:29:28

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	4.42	5.862	592.03	0.51	
Sb	0.37	(-0.427	(-43.16	-9.62	
As	0.86	(-0.092	(-9.287	-8.27	
Ba	2.68	(-0.026	(-2.588	-4.82	
Bm	0.61	(-0.003	(-0.271	-23.20	
Bi	6.77	2.780	280.80	1.31	
B	3.52	(-0.026	(-2.663	-12.24	
Cd	1.84	(-0.018	(-1.816	-4.70	
Ca	1.48	0.038	3.836	0.60	
Ce	3.67	(-1.176	(-118.8	-3.72	
Cr	1.20	(-0.121	(-12.25	-0.41	
Co	0.21	(-0.116	(-11.68	-4.03	
Cu	2.16	(-0.060	(-6.010	-4.30	
Eu	2.82	(-0.022	(-2.222	-2.43	
Fe	16.48	1.911	192.97	0.28	
La	0.32	(-0.083	(-8.407	-7.52	
Pb	0.26	-0.084	-8.451	-13.32	
Li	2.78	(-0.049	(-4.949	-4.31	
Mg	1.93	0.062	6.233	0.69	
Mn	1.05	0.020	2.032	3.84	
Hg	3.10	(-0.027	(-2.747	-4.46	
Mo	1.48	(-0.031	(-3.132	-6.78	
Nd	3.68	(-2.248	(-227.0	-2.14	
Ni	2.83	(-0.052	(-5.211	-6.54	
P	13.62	6.906	697.46	0.04	
K	2.41	(-1.808	(-182.6	-4.10	
Sm	3.49	(-1.337	(-135.0	-4.12	
Se	1.43	(-0.413	(-41.70	-1.44	
Si	2.94	0.086	8.719	5.09	
Ag	9.85	(-0.085	(-8.537	-3.24	
Na	31.62	24.197	2443.9	0.53	
Br	3.31	0.006	0.654	5.46	
S	0.70	-0.018	-1.783	-31.13	
Ta	2.70	(-0.197	(-19.86	-1.53	
Tl	3.01	(-1.380	(-139.4	-6.26	
Th	1.10	(-0.756	(-76.34	-4.85	
Sn	1.09	(-0.061	(-6.208	-5.20	
Ti	2.69	(-0.036	(-3.592	-3.73	
W	1.14	(-0.118	(-11.93	-7.20	
U	3.46	(-7.234	(-730.6	-3.41	
V	3.14	(-0.063	(-6.376	-6.20	
Zn	3.07	0.018	1.847	5.30	
Zr	3.58	(-0.126	(-12.77	-4.24	

Dilution factor : 101.000

Figure 10 ICP Analysis March 12, 1990

Sample name : F116  
 Sample code 1 : SAMPLE  
 Sample code 2 : 500-10  
 Programme : SSI                    12-Mar-90 13:33:45

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	15.53	30.116	632.45	0.17	
Sb	0.36	(-0.465	(-9.773	-3.06	
As	0.86	(-0.106	(-2.232	-12.80	
Ba	2.50	(-0.036	(-0.732	-9.87	
Be	0.58	(-0.004	(-0.080	-4.11	
Bi	20.23	15.046	315.97	0.60	
B	3.31	(-0.040	(-0.840	-9.90	
Cd	1.77	(-0.022	(-0.461	-4.94	
Ca	4.03	0.138	2.898	0.41	
Cr	3.47	(-1.538	(-32.30	-9.97	
Cr	1.87	0.074	1.548	6.92	
Co	0.21	(-0.164	(-3.446	-2.84	
Cu	2.11	(-0.069	(-1.452	-11.21	
Hu	2.60	(-0.031	(-0.641	-8.97	
Fe	73.26	9.146	192.07	0.62	
La	0.31	(-0.119	(-2.506	-9.34	
Pb	0.26	-0.051	-1.081	-99.22	
Li	2.57	(-0.071	(-1.483	-9.03	
Hg	3.56	0.128	2.684	0.84	
Mn	2.55	0.112	2.351	0.72	
Hg	3.07	(-0.029	(-0.612	-10.35	
Ho	1.44	(-0.037	(-0.770	-10.48	
Nd	3.43	(-2.736	(-57.45	-6.88	
Ni	2.72	(-0.064	(-1.348	-9.28	
P	58.17	31.985	671.68	1.30	
K	2.26	(-2.468	(-51.84	-9.04	
Sm	3.23	(-1.880	(-39.48	-9.08	
Se	1.45	(-0.377	(-7.920	-19.64	
Si	4.30	0.848	17.801	3.21	
Ag	9.10	(-0.115	(-2.421	-8.99	
Na	139.71	)119.53	)2510.1	0.42	
Sr	5.41	0.074	1.553	1.71	
S	0.78	0.070	1.476	6.66	
Ta	2.54	(-0.258	(-5.425	-11.41	
Tl	2.71	(-1.996	(-41.92	-6.42	
Th	1.03	(-1.054	(-22.13	-8.36	
Gn	1.07	(-0.069	(-1.459	-11.94	
Ti	2.59	(-0.044	(-0.915	-9.68	
W	1.14	(-0.123	(-2.584	-10.88	
U	3.27	(-9.472	(-198.9	-9.50	
V	2.97	(-0.080	(-1.670	-1.93	
Zn	3.70	0.105	2.199	0.92	
Zr	3.42	(-0.171	(-3.587	-9.47	

Dilution factor : 21.0000

Figure 11 ICP Analysis March 12, 1990

Sample name : 1MHN03  
Sample code 1 : DIRECT  
Programme : SST                    12-Mar-90 13:37:50

NAME	MV	INT	CONCEN	RSD
Al	1.37	(-0.810	-8.81	
Sb	0.35	(-0.627	-7.98	
As	0.80	(-0.149	-7.86	
Ba	2.38	(-0.043	-9.50	
Be	0.56	(-0.004	-8.77	
Bi	2.93	(-0.717	-9.94	
B	3.14	(-0.051	-10.90	
Cd	1.69	(-0.027	-8.62	
Ca	0.49	(-0.001	-24.72	
Ce	3.31	(-1.831	-10.17	
Cr	0.91	(-0.205	-2.36	
Co	0.21	(-0.167	-9.68	
Cu	1.96	(-0.096	-10.92	
Hg	2.55	(-0.033	-9.90	
Fe	1.15	(-0.042	-9.60	
La	0.31	(-0.136	-11.46	
Pb	0.24	(-0.399	-12.80	
Li	2.53	(-0.075	-10.70	
Mg	0.36	(-0.002	-16.60	
Mn	0.60	(-0.008	-10.44	
Hg	2.97	(-0.037	-12.40	
Mo	1.36	(-0.050	-9.31	
Nd	3.36	(-2.871	-6.05	
Ni	2.57	(-0.080	-12.41	
P	1.10	(-0.138	-12.61	
K	2.22	(-2.655	-9.93	
Sm	3.17	(-2.004	-10.46	
Se	1.30	(-0.664	-7.80	
Si	2.20	(-0.334	-10.22	
Ag	8.88	(-0.125	-9.61	
Na	3.22	(-0.850	-9.90	
Sr	2.53	(-0.019	-10.03	
S	0.61	(-0.107	-9.85	
Ta	2.49	(-0.278	-11.59	
Tl	2.82	(-1.780	-17.15	
Th	1.01	(-1.120	-10.45	
Sn	1.01	(-0.095	-12.30	
Ti	2.45	(-0.057	-9.91	
W	1.05	(-0.186	-7.65	
U	3.12	(-11.07	-10.05	
V	3.00	(-0.077	-13.74	
Zn	2.08	(-0.014	-7.80	
Zr	3.35	(-0.190	-10.71	

Figure 12 ICP Analysis March 12, 1990

Sample name : F117  
 Sample code 1 : SAM  
 Sample code 2 : 100-10  
 Sample code 3 : 089045  
 Programme : SST                    12-Mar-90 13:42:45

NAME	HV	INT	CONCEN	DILCOR	RSD
Al	4.15	5.261	531.36	0.25	
Bb	0.38	-0.214	-21.58	-7.70	
Br	0.94	(-0.043	(-4.376	-15.63	
Ca	2.88	(-0.014	(-1.398	-11.40	
Cr	0.65	(-0.001	(-0.130	-38.88	
Di	6.64	2.664	269.08	1.80	
Si	3.78	-0.009	-0.954	-15.43	
Cd	1.97	(-0.010	(-1.004	-10.56	
Ca	3.35	0.111	11.250	0.51	
Ce	4.02	(-0.559	(-56.48	-14.39	
Cr	1.25	(-0.108	(-10.87	-5.27	
Co	0.22	-0.086	-8.694	-9.38	
Cu	2.33	(-0.028	(-2.797	-12.60	
Eu	3.10	(-0.011	(-1.109	-11.66	
Fe	14.74	1.690	170.67	0.51	
La	0.34	-0.044	-4.456	-10.41	
Pb	0.26	-0.058	-5.851	-50.92	
Li	3.01	(-0.024	(-2.465	-10.38	
Mg	1.51	0.045	4.540	0.21	
Mn	1.17	0.027	2.722	1.87	
Hg	3.36	-0.008	-0.779	-6.06	
Mo	1.58	(-0.015	(-1.523	-18.33	
Nd	4.00	(-1.635	(-165.2	-10.32	
Ni	3.07	(-0.024	(-2.384	-20.77	
P	9.75	4.732	477.95	1.46	
K	2.59	(-0.997	(-100.7	-7.53	
Sr	3.83	(-0.652	(-65.86	-12.81	
Se	1.53	(-0.202	(-20.38	-21.47	
Si	3.13	0.190	19.163	4.76	
Ag	10.84	(-0.044	(-4.394	-10.93	
Na	23.36	16.911	1708.0	0.24	
Br	3.44	0.011	1.069	6.94	
S	0.73	0.012	1.189	33.41	
Ta	2.94	(-0.102	(-10.26	-14.34	
Tl	3.39	(-0.601	(-60.70	-30.49	
Th	1.20	(-0.369	(-37.29	-13.60	
Sn	1.17	(-0.030	(-3.012	-28.60	
Ti	2.90	(-0.016	(-1.642	-13.30	
W	1.24	(-0.053	(-5.340	-25.43	
U	3.79	(-3.529	(-356.4	-12.62	
V	3.59	(-0.021	(-2.095	-25.31	
Zn	18.92	0.538	54.380	0.06	
Zr	3.85	-0.052	-5.202	-13.82	

Dilution factor : 101.000

Figure 13 ICP Analysis March 12, 1990

Sample name : F117  
 Sample code 1 : 8AM  
 Sample code 2 : 500-10  
 Sample code 3 : 089045  
 Programme : 88T      12-Mar-90 13:46:53

NAME	MV	INT	CONCEN	DILCOR	RSD
A1	13.33	25.323	531.78	1.94	
Sb	0.40	0.005	0.100	1039.39	
As	1.00	0.001	0.026	330.45	
Ba	3.07	-0.003	-0.065	-106.92	
Be	0.66	-0.001	-0.017	-56.70	
Bi	18.51	13.476	282.99	1.94	
B	4.02	0.007	0.141	94.50	
Cd	2.07	-0.004	-0.077	-43.12	
Ca	3.63	0.122	2.568	2.55	
Ce	4.25	-0.139	-2.911	-105.25	
Cr	1.93	0.090	1.896	5.78	
Co	0.22	-0.056	-1.186	-21.82	
Cu	2.49	0.001	0.024	595.09	
Eu	3.24	-0.006	-0.117	-41.07	
Fe	62.96	7.834	164.51	2.12	
La	0.35	-0.019	-0.400	-41.78	
Pb	0.28	0.290	6.083	3.85	
Li	3.12	-0.013	-0.270	-46.63	
Mg	3.49	0.125	2.629	1.87	
Mn	2.71	0.122	2.558	1.61	
Hg	3.52	0.005	0.095	92.44	
Ho	1.68	0.002	0.032	391.51	
Nd	4.10	(-1.446	(-30.36	-15.11	
Hi	3.26	-0.002	-0.047	-261.86	
P	40.87	22.244	467.13	1.08	
K	2.68	-0.583	-12.24	-29.36	
Sm	3.99	-0.325	-6.825	-47.67	
Se	1.66	0.048	1.000	64.95	
Si	5.86	1.725	36.230	8.31	
Ag	11.34	-0.023	-0.483	-40.09	
Na	92.23	)77.655	)1630.8	2.18	
Br	5.36	0.072	1.518	0.59	
S	0.84	0.134	2.806	1.19	
Ta	3.11	-0.035	-0.737	-59.74	
Tl	3.47	(-0.436	(-9.153	-74.44	
Th	1.24	-0.194	-4.074	-41.16	
Sn	1.22	-0.005	-0.112	-207.73	
Ti	3.13	0.004	0.078	120.29	
W	1.33	0.014	0.300	83.65	
U	4.00	-1.196	-25.12	-72.38	
V	3.58	(-0.022	(-0.461	-55.93	
Zn	5.88	0.110	2.318	1.40	
Zr	4.21	0.047	0.983	28.59	

Dilution factor : 21.0000

Figure 14 ICP Analysis March 12, 1990

Sample name : 1MHN03  
Sample code 1 : DIRECT  
Programme : SST                    12-Mar-90 13:50:44

NAME	MV	INT	CONCEN	RSD
Al	1.46	(-0.609	-7.88	
Sb	0.36	(-0.470	-6.31	
As	0.87	(-0.096	-12.79	
Ba	2.56	(-0.032	-8.84	
Br	0.58	(-0.004	-8.03	
Bi	3.12	(-0.551	-1.90	
B	3.38	(-0.035	-19.96	
Cd	1.81	(-0.020	-7.34	
Ca	0.51	-0.000	-57.28	
Ce	3.57	(-1.363	-8.65	
Cr	0.95	(-0.192	-1.08	
Co	0.21	(-0.148	-11.35	
Cu	2.09	(-0.072	-9.04	
Hg	2.76	(-0.024	-8.40	
Fe	1.24	(-0.031	-11.56	
La	0.32	(-0.091	-6.59	
Pb	0.25	-0.187	-17.92	
Li	2.70	(-0.057	-9.29	
Hg	0.38	-0.001	-17.92	
Hn	0.63	(-0.006	-6.43	
Hg	3.11	(-0.026	-16.05	
Ho	1.43	(-0.038	-11.00	
Nd	3.60	(-2.405	-6.67	
Hi	2.75	(-0.060	-8.75	
P	1.14	(-0.114	-20.71	
K	2.36	(-2.035	-6.69	
Sm	3.41	(-1.496	-9.05	
Sm	1.39	(-0.486	-7.05	
Si	2.34	(-0.252	-9.19	
Ag	9.62	(-0.094	-9.24	
Na	3.47	(-0.629	-8.11	
Br	2.67	(-0.014	-8.96	
S	0.65	(-0.067	-14.69	
Ta	2.65	(-0.215	-7.99	
Tl	2.82	(-1.775	-13.14	
Th	1.08	(-0.843	-8.42	
Sn	1.07	(-0.070	-0.92	
Ti	2.61	(-0.043	-7.36	
W	1.12	(-0.136	-10.33	
U	3.37	(-8.284	-8.56	
V	3.04	(-0.073	-11.84	
Zn	2.17	(-0.011	-5.94	
Zr	3.52	(-0.143	-8.37	

Figure 15 ICP Analysis March 12, 1990

Sample name : F118  
Sample code 1 : SAM  
Sample code 2 : 100-10  
Sample code 3 : 089045  
Programme : SST

12-Mar-90 13:55:27

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	3.98	4.888	493.71	1.59	
Bb	0.36	(-0.513	(-51.80	-8.93	
As	0.82	(-0.133	(-13.46	-8.79	
Ba	4.21	0.063	6.400	3.87	
Be	0.55	(-0.005	(-0.485	-6.68	
Bi	6.27	2.323	234.64	2.86	
B	5.42	0.098	9.853	41.01	
Cd	3.28	0.069	6.993	4.14	
Ca	4.08	0.140	14.158	0.58	
Cr	3.25	(-1.932	(-195.1	-6.34	
Co	1.43	(-0.055	(-5.574	-9.84	
Co	0.21	(-0.132	(-13.31	-7.07	
Cu	2.58	0.017	1.693	30.93	
Eu	2.50	(-0.034	(-3.481	-5.64	
Fe	14.73	1.688	170.47	0.64	
La	0.34	-0.027	-2.735	-44.91	
Pb	0.25	-0.238	-24.05	-20.41	
Li	3.73	0.053	5.379	8.96	
Mg	3.71	0.134	13.541	1.01	
Mn	2.70	0.122	12.279	1.46	
Hg	3.31	-0.011	-1.144	-210.83	
Mo	2.00	0.051	5.172	5.44	
Nd	3.34	(-2.907	(-293.6	-4.07	
Ni	3.41	0.015	1.545	37.21	
P	10.37	5.077	512.80	2.84	
K	2.21	(-2.693	(-272.2	-5.08	
Sm	3.09	(-2.149	(-217.0	-5.17	
Ge	1.34	(-0.596	(-60.15	-6.57	
Si	2.57	(-0.127	(-12.79	-15.66	
Ag	8.98	(-0.121	(-12.17	-5.35	
Na	25.80	19.062	1925.2	0.70	
Sr	6.42	0.106	10.755	1.33	
S	0.64	(-0.079	(-7.971	-6.96	
Ta	2.44	(-0.297	(-30.01	-6.33	
Tl	2.56	(-2.105	(-212.6	-11.68	
Th	0.99	(-1.200	(-121.2	-5.44	
Sn	1.22	-0.009	-0.867	-17.73	
Ti	3.41	0.028	2.864	12.39	
W	1.05	(-0.183	(-18.45	-5.20	
U	3.07	(-11.63	(-1175	-5.60	
V	2.84	(-0.091	(-9.218	-2.48	
Zn	5.91	0.112	11.264	2.35	
Zr	3.31	(-0.200	(-20.23	-5.16	

Dilution factor : 101,000

Figure 16 ICP Analysis March 12, 1990

Sample name : F118  
 Sample code 1 : SAM  
 Sample code 2 : 500-10  
 Sample code 3 : 089045  
 Programme : SST                            12-Mar-90 13:59:44

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	13.28	25.201	529.23	1.31	
Br	0.41	0.119	2.493	63.50	
As	0.95	(-0.040	(-0.841	-57.10	
Ba	10.90	0.451	9.475	1.20	
Be	0.61	(-0.003	(-0.058	-20.90	
Bi	18.77	13.714	287.93	0.32	
B	11.50	0.494	10.370	4.31	
Cd	9.56	0.448	9.407	2.32	
Ca	16.44	0.626	13.139	2.12	
Cr	3.76	(-1.020	(-21.42	-25.20	
Co	3.44	0.526	11.044	1.26	
Co	0.25	0.215	4.519	28.64	
Cu	5.02	0.462	9.709	0.79	
Eu	2.87	(-0.020	(-0.419	-22.63	
Fe	63.63	7.919	166.29	0.54	
La	0.49	0.411	8.634	2.64	
Pb	0.29	0.451	9.462	15.45	
Li	8.24	0.528	11.080	1.44	
Mg	15.37	0.606	12.734	0.59	
Mn	10.09	0.576	12.100	0.24	
Hg	3.90	0.033	0.685	16.67	
No	4.32	0.451	9.472	1.38	
Nd	3.94	(-1.743	(-36.64	-18.85	
Ni	7.02	0.429	9.003	3.14	
P	44.29	24.172	507.61	0.90	
K	2.52	(-1.285	(-26.98	-26.51	
Sn	3.51	(-1.304	(-27.39	-22.16	
Se	1.60	-0.077	-1.626	-105.69	
Si	4.24	0.814	17.092	5.11	
Ag	11.31	-0.024	-0.511	-69.94	
Na	103.58	187.663	1840.9	2.11	
Sr	20.07	0.546	11.468	1.71	
S	0.81	0.097	2.035	15.61	
Ta	2.76	(-0.171	(-3.585	-22.26	
Tl	3.09	(-1.223	(-25.67	-22.26	
Th	1.11	(-0.711	(-14.94	-23.59	
Sn	2.17	0.395	8.291	7.09	
Ti	7.16	0.361	7.573	0.84	
W	1.21	(-0.072	(-1.508	-29.66	
U	3.54	(-6.439	(-135.2	-23.62	
V	3.33	(-0.045	(-0.947	-16.64	
Zn	19.86	0.570	11.960	0.64	
Zr	3.63	(-0.113	(-2.381	-24.15	

Bilution factor : 21.0000

Figure 17 ICP Analysis March 12, 1990

Sample name : 1MHN03  
Sample code 1 : DIRECT  
Programme : SST                    12-Mar-90 14:03:36

NAME	MV	INT	CONCEN	RSD
Al	1.63	(-0.244	-20.86	
Sb	0.38	-0.233	-24.74	
As	0.54	(-0.048	-26.06	
Ba	2.91	(-0.012	-21.67	
Be	0.64	(-0.002	-5.00	
Bi	3.51	(-0.194	-3.73	
Ca	4.10	0.012	16.57	
Cd	1.98	(-0.009	-10.60	
Co	0.52	0.000	26.02	
Cs	4.06	(-0.480	-23.94	
Cr	1.06	(-0.163	-2.53	
Co	0.22	-0.059	-20.83	
Cu	2.34	(-0.026	-21.04	
Dy	3.15	(-0.009	-23.74	
Fe	1.38	(-0.013	-23.55	
La	0.34	-0.035	-13.09	
Pb	0.26	-0.084	-81.04	
Li	3.05	(-0.020	-21.42	
Mg	0.41	0.000	48.50	
Mn	0.69	-0.002	-27.52	
Hg	3.88	0.031	4.77	
Mo	1.58	(-0.015	-24.79	
Nd	4.02	(-1.601	-7.56	
Ni	3.09	(-0.021	-12.86	
P	1.26	(-0.047	-28.49	
K	2.62	(-0.854	-21.05	
Sm	3.88	(-0.538	-24.36	
Sc	1.53	(-0.210	-28.43	
Si	2.63	(-0.093	-19.02	
Ag	11.00	(-0.037	-17.00	
Na	3.94	(-0.220	-22.27	
Sr	2.96	-0.005	-25.09	
S	0.69	-0.029	-22.08	
Ta	2.98	(-0.086	-21.10	
Tl	3.35	(-0.678	-18.96	
Th	1.21	(-0.317	-21.26	
Sn	1.16	(-0.033	-25.17	
Fr	2.90	(-0.016	-18.31	
W	1.23	(-0.055	-29.64	
U	3.84	(-3.022	-21.43	
V	3.51	(-0.028	-36.20	
Zn	2.37	(-0.005	-16.80	
Zr	3.85	-0.053	-20.01	

Figure 18 ICP Analysis March 12, 1990

Sample name : FG00  
 Sample code 1 : SAM  
 Sample code 2 : 100-10  
 Programme : SST                            12-Mar-90 14:09:12

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	5.02	7.164	723.58	0.93	
Sb	0.38	-0.266	-26.86	-32.59	
As	0.92	(-0.060	(-6.106	-47.15	
Ba	2.81	(-0.018	(-1.805	-37.86	
Be	0.62	(-0.002	(-0.224	-30.34	
Bi	3.37	(-0.324	(-32.70	-34.25	
B	3.87	-0.004	-0.357	-307.74	
Cd	1.93	(-0.012	(-1.243	-27.44	
Ca	0.86	0.014	1.382	1.04	
Ce	3.91	(-0.752	(-75.98	-38.64	
Cr	1.03	(-0.169	(-17.11	-5.31	
Co	0.22	-0.081	-8.151	-35.12	
Cu	2.26	(-0.041	(-4.152	-37.36	
Eu	3.03	(-0.014	(-1.395	-35.19	
Fe	1.62	0.017	1.725	36.46	
La	0.34	-0.049	-4.963	-42.42	
Pb	0.26	-0.064	-6.501	-86.60	
Li	2.94	(-0.032	(-3.201	-37.05	
Mg	0.83	0.017	1.752	1.29	
Mn	0.69	(-0.002	(-0.236	-41.80	
Hg	3.71	0.019	1.885	40.43	
Mo	1.55	(-0.019	(-1.875	-35.98	
Nd	3.89	(-1.848	(-186.7	-16.40	
Ni	2.96	(-0.036	(-3.598	-43.14	
P	1.25	(-0.057	(-5.723	-57.72	
K	2.53	(-1.240	(-125.3	-29.99	
Sm	3.74	(-0.839	(-84.71	-38.54	
Se	1.51	(-0.252	(-25.46	-42.57	
Si	2.70	-0.053	-5.402	-97.26	
Ag	10.58	(-0.054	(-5.470	-35.61	
Na	3.82	(-0.325	(-32.87	-43.23	
Br	2.88	(-0.007	(-0.755	-40.81	
S	0.68	(-0.039	(-3.951	-49.41	
Ta	2.89	(-0.123	(-12.39	-41.18	
Tl	3.15	(-1.088	(-109.8	-35.76	
Th	1.17	(-0.474	(-47.90	-38.82	
Sn	1.13	(-0.044	(-4.416	-22.96	
Ti	2.83	(-0.023	(-2.277	-37.44	
W	1.21	(-0.069	(-6.971	-44.61	
U	3.69	(-4.650	(-469.6	-36.96	
V	3.34	(-0.044	(-4.436	-32.08	
Zn	2.68	0.005	0.538	42.22	
Zr	3.75	(-0.080	(-8.118	-37.49	

Dilution factor : 101.000

Figure 19 ICP Analysis March 12, 1990

Sample name : F600  
Sample code 1 : SAM  
Sample code 2 : 500-10  
Programme : SST

12-Mar-90 14:13:20

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	19.04	37.794	793.67	0.66	
Bb	0.36	(-0.480	(-10.07	-10.71	
As	0.87	(-0.097	(-2.042	-11.24	
Ba	2.46	(-0.038	(-0.801	-6.91	
Bm	0.57	(-0.004	(-0.083	-3.28	
Bi	3.01	(-0.544	(-13.53	-5.74	
B	3.48	(-0.029	(-0.604	-14.29	
Cd	1.75	(-0.023	(-0.484	-3.71	
Ca	1.20	0.027	0.571	1.06	
Ce	3.42	(-1.638	(-34.39	-7.04	
Cr	0.92	(-0.202	(-4.233	-1.82	
Co	0.21	(-0.148	(-3.107	-3.15	
Cu	2.02	(-0.085	(-1.789	-6.83	
Eu	2.64	(-0.029	(-0.612	-6.82	
Fe	2.54	0.134	2.821	1.86	
La	0.31	(-0.113	(-2.380	-5.33	
Pb	0.26	-0.129	-2.703	-25.98	
Li	2.61	(-0.067	(-1.403	-7.44	
Hg	0.60	0.008	0.165	3.02	
Mn	0.72	-0.001	-0.013	-81.45	
Hg	3.42	-0.004	-0.074	-129.03	
No	1.46	(-0.034	(-0.707	-11.20	
Nd	3.44	(-2.711	(-56.94	-4.69	
Ni	2.65	(-0.072	(-1.503	-3.97	
P	1.17	(-0.099	(-2.073	-15.57	
K	2.28	(-2.376	(-49.90	-6.09	
Sm	3.27	(-1.789	(-37.56	-7.20	
Se	1.42	(-0.423	(-8.893	-8.82	
Bi	2.68	-0.060	-1.265	-34.77	
Ag	9.20	(-0.111	(-2.335	-7.36	
Na	3.43	(-0.667	(-14.02	-7.39	
Sr	2.60	(-0.016	(-0.344	-7.40	
S	0.64	(-0.077	(-1.614	-1.35	
Ta	2.54	(-0.258	(-5.417	-9.70	
Tl	2.74	(-1.934	(-40.61	-3.01	
Th	1.04	(-1.005	(-21.11	-8.45	
Sn	1.08	(-0.067	(-1.415	-16.09	
Ti	2.60	(-0.043	(-0.900	-8.47	
W	1.14	(-0.122	(-2.565	-10.29	
U	3.23	(-9.859	(-207.0	-7.25	
V	2.94	(-0.082	(-1.732	-3.07	
Zn	2.67	0.005	0.109	14.77	
Zr	3.42	(-0.169	(-3.557	-6.73	

Dilution factor : 21.0000

Figure 20 ICP Analysis March 12, 1990

Sample name : F601  
 Sample code 1 : SAM  
 Sample code 2 : 500-10  
 Programme : SST                    12-Mar-90 14:18:34

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	15.82	30.750	645.75	1.84	
Sb	0.35	(-0.627	(-13.16	-7.98	
As	0.84	(-0.121	(-2.550	-17.39	
Bm	2.37	(-0.044	(-0.918	-12.17	
Bm	0.56	(-0.005	(-0.096	-7.95	
Bi	2.92	(-0.730	(-15.34	-13.86	
B	3.39	(-0.035	(-0.726	-25.33	
Cd	1.71	(-0.026	(-0.344	-10.69	
Ca	1.33	0.032	0.674	2.06	
Co	3.28	(-1.880	(-39.49	-12.15	
Cr	0.90	(-0.209	(-4.387	-2.99	
Co	0.21	(-0.156	(-3.276	-11.95	
Cu	1.94	(-0.099	(-2.072	-12.48	
Bu	2.54	(-0.033	(-0.696	-12.32	
Te	3.96	0.315	6.619	3.55	
La	0.31	(-0.139	(-2.927	-9.41	
Pb	0.25	(-0.270	(-5.677	-22.96	
Li	2.51	(-0.077	(-1.615	-11.93	
Hg	0.58	0.007	0.150	5.12	
Mn	0.69	(-0.002	(-0.050	-43.07	
Hg	3.34	-0.009	-0.195	-86.40	
Mo	1.41	(-0.041	(-0.870	-15.89	
Nd	3.31	(-2.959	(-62.15	-6.60	
Ni	2.60	(-0.077	(-1.627	-13.81	
P	1.15	(-0.112	(-2.344	-21.25	
K	2.20	(-2.713	(-56.97	-11.36	
Sm	3.14	(-2.056	(-43.17	-12.11	
Sr	1.37	(-0.523	(-10.99	-13.26	
Si	2.54	(-0.142	(-2.979	-30.13	
Ag	8.82	(-0.127	(-2.665	-12.03	
Na	3.30	(-0.783	(-16.44	-12.95	
Sr	2.53	(-0.019	(-0.394	-12.67	
S	0.62	(-0.094	(-1.970	-20.30	
Ta	2.46	(-0.287	(-6.034	-10.95	
Il	2.73	(-1.950	(-40.96	-14.55	
Th	1.00	(-1.147	(-24.08	-11.41	
Sn	1.02	(-0.090	(-1.884	-12.24	
Ti	2.51	(-0.051	(-1.073	-12.91	
W	1.10	(-0.149	(-3.135	-16.04	
U	3.10	(-11.33	(-237.9	-11.84	
V	2.95	(-0.081	(-1.694	-17.71	
Zn	2.28	(-0.008	(-0.164	-14.59	
Zr	3.33	(-0.194	(-4.076	-12.37	

Dilution factor : 21.0000

Figure 21 ICP Analysis March 12, 1990

Sample name : F188  
 Sample code 1 : SAM  
 Sample code 2 : 100-10  
 Sample code 3 : 089049  
 Programme : SST                    12-Mar-90 14:23:10

NAME	MV	INI	CONCEN	DILCOR	RSD
Al	5.06	7.256	732.85	0.34	
Sb	0.37	(-0.328	(-33.09	-17.57	
As	0.92	(-0.058	(-5.852	-19.10	
Ba	2.82	(-0.017	(-1.729	-20.82	
Be	0.62	(-0.002	(-0.226	-21.15	
Bi	8.32	4.186	422.83	1.87	
B	3.87	-0.003	-0.322	-110.06	
Cd	1.93	(-0.013	(-1.282	-20.51	
Ca	1.71	0.047	4.738	0.81	
Ce	3.92	(-0.736	(-74.29	-20.59	
Cr	1.23	(-0.113	(-11.38	-7.38	
Co	0.22	-0.065	-6.521	-19.09	
Cu	2.31	(-0.032	(-3.202	-25.00	
Bu	3.02	(-0.014	(-1.430	-19.25	
Fe	22.02	2.616	264.25	0.62	
La	0.33	(-0.061	(-6.178	-22.54	
Pb	0.26	-0.109	-11.05	-46.69	
Li	2.94	(-0.032	(-3.254	-21.00	
Hg	1.41	0.041	4.119	0.57	
Mn	1.20	0.029	2.923	2.90	
Hg	3.69	0.017	1.701	38.25	
No	1.54	(-0.021	(-2.077	-21.13	
Nd	3.86	(-1.913	(-193.2	-10.08	
Ni	3.02	(-0.029	(-2.970	-21.96	
P	11.32	5.612	566.81	0.98	
K	2.54	(-1.233	(-124.5	-13.19	
Sm	3.73	(-0.860	(-86.84	-20.10	
Se	1.52	(-0.232	(-23.46	-22.66	
Si	3.13	0.188	19.030	13.75	
Ag	10.54	(-0.056	(-5.658	-19.24	
Na	28.39	21.348	2156.1	0.68	
Sr	3.89	0.025	2.520	5.14	
S	0.69	-0.023	-2.307	-24.05	
Ta	2.87	(-0.128	(-12.90	-23.30	
Tl	3.24	(-0.898	(-90.74	-19.40	
Th	1.16	(-0.497	(-50.21	-18.72	
Sn	1.13	(-0.044	(-4.489	-17.58	
Ti	2.83	(-0.022	(-2.253	-22.43	
W	1.21	(-0.068	(-6.900	-13.36	
U	3.70	(-4.616	(-466.2	-19.51	
V	3.43	(-0.036	(-3.600	-11.59	
Zn	3.05	0.018	1.780	9.90	
Zr	3.74	(-0.082	(-8.322	-19.86	

Dilution factor : 101.000

Figure 22 ICP Analysis March 12, 1990

Sample name : F189  
 Sample code 1 : SAM  
 Sample code 2 : 500-10  
 Sample code 3 : 089049  
 Programme : SST

12-Mar-90 14:27:39

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	16.61	32.489	682.27	2.29	
Sb	0.39	-0.066	-1.396	-119.31	
As	1.00	0.000	0.000	8784300	
Ba	3.08	-0.003	-0.054	-169.54	
Be	0.66	-0.001	-0.017	-79.64	
Bi	24.98	19.372	406.81	1.27	
B	4.02	0.007	0.143	87.87	
Cd	2.07	-0.004	-0.084	-63.15	
Ca	4.92	0.173	3.636	2.91	
Cr	4.22	-0.193	-4.053	-95.73	
Co	2.01	0.112	2.358	3.81	
Co	0.22	-0.046	-0.960	-46.69	
Cu	2.62	0.025	0.516	37.89	
Hu	3.19	(-0.007	(-0.154	-43.19	
Fe	90.40	11.331	237.95	1.98	
La	0.34	-0.021	-0.442	-54.09	
Pb	0.28	0.296	6.218	11.30	
Li	3.08	(-0.017	(-0.363	-47.72	
Mg	4.26	0.156	3.279	2.02	
Mn	2.97	0.138	2.901	1.03	
Hg	3.88	0.031	0.657	17.99	
Mo	1.67	-0.000	-0.002	-5626.7	
Nd	4.04	(-1.548	(-32.50	-14.09	
Ni	3.24	-0.004	-0.077	-276.68	
P	45.52	24.863	522.13	2.18	
K	2.64	(-0.747	(-15.68	-40.28	
Sm	3.92	(-0.454	(-9.543	-47.66	
Se	1.66	0.047	0.987	143.56	
Si	5.01	1.249	26.219	0.90	
Ag	11.15	(-0.031	(-0.644	-41.15	
Na	108.45	91.955	1931.1	2.61	
Sr	7.34	0.136	2.861	1.55	
S	0.80	0.084	1.766	17.15	
Ta	3.06	-0.053	-1.122	-55.36	
Tl	3.42	(-0.535	(-11.24	-34.86	
Th	1.22	(-0.267	(-5.602	-41.56	
Sn	1.23	-0.001	-0.027	-485.74	
Ti	3.04	-0.004	-0.084	-142.04	
W	1.38	0.050	1.056	54.38	
U	3.96	-1.710	-35.91	-65.14	
V	3.60	(-0.020	(-0.420	-39.56	
Zn	5.69	0.104	2.190	0.75	
Zr	3.88	-0.044	-0.918	-43.39	

Dilution factor : 21.0000

Figure 23 ICP Analysis March 12, 1990

Sample name : IMHN03  
Sample code 1 : DIRECT  
Programme : SST 12-Mar-90 14:31:38

NAME	MV	INT	CONCEN	RSD
Al	1.44	(-0.654	-14.42	
Sb	0.36	(-0.494	-10.40	
As	0.85	(-0.117	-13.27	
Ba	2.51	(-0.036	-13.91	
Ba	0.58	(-0.004	-9.46	
Bi	3.07	(-0.597	-14.92	
B	3.47	(-0.030	-23.54	
Cd	1.77	(-0.022	-11.47	
Ca	0.51	(-0.000	-118.88	
Cr	3.50	(-1.492	-13.92	
Cr	0.93	(-0.199	-4.53	
Co	0.21	(-0.121	-7.70	
Cu	2.05	(-0.080	-13.69	
Bu	2.70	(-0.027	-14.08	
Fe	1.21	(-0.035	-15.86	
La	0.32	(-0.099	-14.95	
Pb	0.25	(-0.270	-10.91	
Li	2.66	(-0.052	-14.31	
Mg	0.38	(-0.001	-21.82	
Mn	0.62	(-0.006	-10.49	
Hg	3.43	(-0.002	-266.04	
Ho	1.42	(-0.040	-13.87	
Nd	3.48	(-2.633	-12.17	
Ni	2.70	(-0.066	-12.84	
P	1.14	(-0.117	-13.94	
K	2.31	(-2.220	-11.56	
Sm	3.34	(-1.642	-13.84	
Se	1.37	(-0.524	-11.97	
Si	2.30	(-0.279	-13.09	
Ag	9.39	(-0.103	-13.08	
Na	3.40	(-0.697	-15.15	
Sr	2.63	(-0.016	-14.21	
S	0.63	(-0.090	-5.92	
Ta	2.60	(-0.232	-13.16	
Tl	2.81	(-1.787	-16.24	
Th	1.06	(-0.923	-13.97	
Sn	1.03	(-0.089	-15.78	
Ti	2.56	(-0.046	-13.43	
W	1.08	(-0.161	-12.63	
U	3.30	(-9.049	-13.83	
V	3.01	(-0.075	-8.17	
Zn	2.13	(-0.013	-10.55	
Zr	3.47	(-0.156	-13.92	

Figure 24 ICP Analysis March 12, 1990

Sample name : F189  
 Sample code 1 : SAM  
 Sample code 2 : 100-10  
 Sample code 3 : 089049  
 Programme : SST                            12-Mar-90 14:35:30

NAME	MV	INT	CONCEN	DILCOR	RSD
Al	4.24	5.461	551.58	1.30	
Bb	0.35	(-0.627	(-63.31	-2.62	
As	0.81	(-0.145	(-14.63	-2.63	
Ba	2.36	(-0.044	(-4.433	-1.21	
Be	0.55	(-0.005	(-0.490	-3.01	
Bi	7.22	3.191	322.25	0.90	
B	3.29	(-0.041	(-4.154	-4.70	
Cd	1.69	(-0.027	(-2.734	-2.31	
Ca	2.19	0.066	6.664	0.70	
Cr	3.27	(-1.904	(-192.3	-0.89	
Cr	1.09	(-0.154	(-15.57	-2.02	
Co	0.21	(-0.156	(-15.76	-5.97	
Cu	2.36	(-0.023	(-2.349	-4.59	
Eu	2.50	(-0.034	(-3.469	-1.43	
Fe	17.75	2.073	209.38	0.72	
La	0.31	(-0.137	(-13.88	-1.26	
Pb	0.25	(-0.341	(-34.46	-11.32	
Li	2.49	(-0.080	(-8.043	-1.82	
Mg	1.31	0.037	3.710	0.96	
Mn	1.60	0.054	5.413	0.78	
Hg	3.33	-0.010	-1.031	-34.59	
Mo	1.35	(-0.050	(-5.081	-4.04	
Nd	3.28	(-3.015	(-304.6	-3.49	
Ni	2.57	(-0.081	(-8.139	-1.94	
P	11.33	5.620	567.63	0.93	
K	2.18	(-2.808	(-283.6	-1.04	
Sm	3.11	(-2.121	(-214.2	-1.06	
Se	1.31	(-0.645	(-65.16	-3.93	
Si	2.77	-0.010	-1.005	-59.64	
A <sub>9</sub>	8.71	(-0.131	(-13.28	-1.60	
Na	27.77	20.804	2101.2	0.68	
Sr	3.38	0.009	0.864	4.36	
S	0.62	(-0.097	(-9.789	-2.23	
Ta	2.43	(-0.298	(-30.12	-1.91	
Tl	2.55	(-2.329	(-235.3	-3.59	
Th	0.99	(-1.186	(-119.8	-1.10	
Sn	1.00	(-0.102	(-10.26	-5.73	
Ti	2.44	(-0.057	(-5.756	-1.33	
W	1.06	(-0.173	(-17.51	-8.80	
U	3.09	(-11.50	(-1161	-1.18	
V	2.86	(-0.090	(-9.059	-2.38	
Zn	2.75	0.008	0.768	7.07	
Zr	3.30	(-0.202	(-20.45	-1.26	

Dilution factor : 101.000

Figure 25 ICP Analysis March 12, 1990

Sample name : F189  
 Sample code 1 : SAM  
 Sample code 2 : 500-10  
 Sample code 3 : 089049  
 Programme : SST                    12-Mar-90 14:39:48

NAME	MV	INT	CONCEN	DILCOR	RSN
Al	14.61	28.115	590.42	1.67	
Sb	0.36	(-0.465	(-9.773	-20.07	
As	0.87	(-0.095	(-2.005	-24.71	
Ba	2.67	(-0.026	(-0.545	-25.57	
Bm	0.60	(-0.003	(-0.064	-21.52	
Bi	23.24	17.789	373.58	0.50	
B	3.59	(-0.022	(-0.453	-34.56	
Cd	1.82	(-0.019	(-0.403	-13.60	
Ca	7.42	0.271	5.696	2.49	
Ce	3.62	(-1.271	(-26.70	-22.81	
Cr	1.88	0.076	1.591	17.61	
Co	0.21	(-0.126	(-2.655	-12.77	
Cu	4.05	0.284	5.971	2.24	
Eu	2.73	(-0.025	(-0.532	-19.64	
Fe	78.18	9.773	205.24	0.64	
La	0.32	(-0.103	(-2.211	-17.22	
Pb	0.26	-0.097	-2.028	-50.33	
Li	2.68	(-0.060	(-1.251	-20.75	
Mg	4.00	0.146	3.062	0.57	
Mn	5.30	0.281	5.908	0.50	
Hg	3.59	0.010	0.200	94.18	
Mo	1.48	(-0.030	(-0.627	-24.75	
Nd	3.56	(-2.475	(-51.98	-9.93	
Ni	2.96	(-0.037	(-0.769	-35.06	
P	49.70	27.220	571.62	0.31	
K	2.33	(-2.142	(-44.98	-20.15	
Sm	3.38	(-1.575	(-33.07	-20.37	
Se	1.47	(-0.321	(-6.739	-18.83	
Si	4.90	1.188	24.947	3.75	
Ag	9.55	(-0.097	(-2.036	-19.02	
Na	114.60	)97.385	)2045.1	2.57	
Br	6.65	0.114	2.396	0.35	
S	0.73	0.017	0.363	143.16	
Ta	2.67	(-0.207	(-4.331	-23.15	
Fl	2.99	(-1.420	(-29.82	-10.18	
In	1.07	(-0.684	(-18.56	-19.47	
Sn	1.12	(-0.051	(-1.072	-32.67	
Tl	2.73	(-0.032	(-0.671	-27.09	
W	1.24	(-0.050	(-1.046	-44.68	
U	3.42	(-7.770	(-163.2	-22.26	
V	3.25	(-0.053	(-1.118	-20.06	
Zn	5.50	0.098	2.056	1.25	
Zr	3.50	(-0.149	(-3.126	-19.49	

Dilution factor : 21.0000

Figure 26 ICP Analysis March 12, 1990

Sample name : 1MHN03  
Sample code 1 : DIRECI  
Programme : SST      12-Mar-90 14:43:37

NAME	MV	INT	CONCEN	RSD
Al	1.58	(-0.347	-14.21	
Sb	0.37	(-0.313	-6.94	
As	0.92	(-0.064	-10.91	
Ba	2.81	(-0.018	-13.45	
Be	0.62	(-0.002	-14.01	
Bi	3.40	(-0.294	-9.80	
B	3.85	-0.005	-81.17	
Cd	1.92	(-0.013	-11.99	
Ca	0.52	0.000	43.27	
Ce	3.93	(-0.719	-14.70	
Cr	1.02	(-0.174	-2.31	
Co	0.22	-0.075	-22.30	
Cu	2.27	(-0.039	-13.51	
Eu	3.05	(-0.013	-13.02	
Fe	1.34	(-0.018	-11.27	
La	0.33	(-0.055	-23.78	
Pb	0.26	-0.161	-18.33	
Li	2.96	(-0.030	-15.87	
Mg	0.40	-0.000	-150.31	
Mn	0.67	(-0.004	-10.43	
Hg	3.81	0.026	5.46	
Mo	1.54	(-0.021	-14.05	
Nd	3.88	(-1.858	-4.17	
Hi	3.00	(-0.031	-11.55	
P	1.23	(-0.067	-16.31	
K	2.55	(-1.170	-13.22	
Sm	3.76	(-0.788	-14.29	
Sm	1.50	(-0.278	-12.64	
Si	2.54	(-0.142	-12.44	
Ag	10.63	(-0.052	-14.32	
Na	3.81	(-0.332	-12.71	
Sr	2.88	(-0.007	-13.70	
S	0.68	-0.036	-15.89	
Ta	2.90	(-0.119	-11.59	
Tl	3.24	(-0.914	-15.76	
Th	1.17	(-0.454	-13.51	
Sn	1.13	(-0.047	-11.79	
Ti	2.82	(-0.023	-12.16	
W	1.19	(-0.084	-1.94	
U	3.71	(-4.477	-13.25	
V	3.45	(-0.034	-24.03	
Zn	2.30	(-0.007	-7.62	
Zr	3.76	(-0.077	-11.61	

Figure 27 ICP Analysis March 12, 1990

Sample name : 1MHN03  
Sample code 1 : DIRECT  
Programme : SST . 12-Mar-90 14:43:37

NAME	MV	INT	CONCEN	RSD
Al	1.58	(-0.347	-14.21	
Sb	0.37	(-0.313	-6.94	
As	0.92	(-0.054	-10.91	
Ba	2.81	(-0.018	-13.45	
Br	0.62	(-0.002	-14.01	
Bi	3.40	(-0.294	-9.80	
B	3.85	-0.005	-81.17	
Cd	1.92	(-0.013	-11.99	
Ca	0.52	0.000	43.27	
Ce	3.93	(-0.719	-14.70	
Cr	1.02	(-0.174	-2.31	
Co	0.22	-0.075	-22.30	
Cu	2.27	(-0.039	-13.51	
Eu	3.05	(-0.013	-13.02	
Fe	1.34	(-0.018	-11.27	
La	0.33	(-0.055	-23.78	
Pb	0.26	-0.161	-18.33	
Li	2.96	(-0.030	-15.87	
Mg	0.40	-0.000	-150.31	
Mn	0.67	(-0.004	-10.43	
Hg	3.81	0.026	5.46	
Ho	1.54	(-0.021	-14.05	
Nd	3.88	(-1.858	-4.17	
Ni	3.00	(-0.031	-11.55	
P	1.23	(-0.067	-16.31	
K	2.55	(-1.170	-13.22	
Sr	3.76	(-0.788	-14.29	
Se	1.50	(-0.278	-12.64	
Si	2.54	(-0.142	-12.44	
Ag	10.63	(-0.052	-14.32	
Na	3.81	(-0.332	-12.71	
Sr	2.88	(-0.007	-13.70	
S	0.68	-0.036	-15.89	
Ta	2.90	(-0.119	-11.59	
Tl	3.24	(-0.914	-15.76	
Th	1.17	(-0.454	-13.51	
Sn	1.13	(-0.047	-11.79	
Ti	2.82	(-0.023	-12.16	
W	1.19	(-0.084	-1.94	
U	3.71	(-4.477	-13.25	
V	3.43	(-0.034	-24.03	
Zn	2.30	(-0.007	-7.62	
Zr	3.76	(-0.077	-11.61	

Figure 27 ICP Analysis March 12, 1990

Sample name : 78C11C  
Sample code 1 : SST1  
Sample code 2 : DIRECT  
Programme : SST

16:00  
11

12-Mar-90 T114743

NAME	MV	INT	CONCEN	RSD
Al	1.76	0.046	139.56	
Sb	1.11	10.129	✓ 1.63	
As	1.11	0.083	6.42	
Ba	164.17	9.337	✓ 2.69	
Be	0.68	-0.000	-346.41	
Bi	3.70	-0.015	-483.10	
B	146.85	9.306	✓ 2.70	
Cd	165.44	9.854	✓ 1.21	
Ca	244.10	9.568	✓ 2.68	
Cr	9.78	9.780	✓ 0.44	
Cr	35.27	9.742	✓ 2.00	
Co	1.17	7.626	✓ 4.35	
Cu	56.59	9.837	✓ 2.47	
Eu	3.96	0.023	7.83	
Fe	78.13	9.766	✓ 1.80	
La	0.37	0.060	14.43	
Pb	0.27	0.122	27.35	
Li	102.30	10.453	✓ 2.73	
Mg	242.31	9.804	✓ 1.72	
Mn	159.84	9.798	✓ 1.94	
Hg	3.77	0.023	31.33	
Ho	1.74	0.011	39.70	
Nd	9.29	8.568	✓ 2.70	
Ni	90.33	9.939	✓ 1.43	
P	1.45	0.057	18.18	
K	8.46	25.260	✓ 1.29	
Sm	4.00	-0.289	-66.99	
Se	3.52	3.730	0.98	
Si	2.78	-0.008	-347.93	
Ag	11.50	-0.016	-66.14	
Nm	32.99	25.407	✓ 1.98	
Sr	296.21	9.436	✓ 2.61	
S	0.95	0.246	2.72	
Ia	3.15	-0.019	-150.23	
Il	3.68	-0.006	-3906.7	
Th	1.31	0.098	103.72	
Sn	125.57	52.459	✓ 1.12	
Ti	3.00	-0.007	-62.07	
W	1.55	0.169	14.24	
U	4.33	2.501	36.13	
V	3.62	(-0.018	-58.25	
Zn	301.23	9.806	✓ 1.54	
Zr	3.99	-0.014	-114.13	

Figure 28 ICP Analysis March 12, 1990

*82930A*

Sample name : B1B38C  
 Sample code 1 : SST2  
 Sample code 2 : DIRECT  
 Programme : SST                    12-Mar-90 14:52:12

NAME	MV	INT	CONCEN	RSD
Al	3.43	3.691	1.73	
Sb	0.40	0.043	183.59	
As	2.74	1.317	2.72	
Ba	2.06	(-0.015	-26.08	
Be	0.61	(-0.003	-9.04	
Bi	62.12	53.211✓	0.44	
B	4.44	0.034	24.42	
Cd	1.91	(-0.014	-15.03	
Ca	1.26	0.029	2.55	
Ce	3.87	(-0.817	-24.22	
Cr	1.31	(-0.090	-5.86	
Co	0.22	-0.073	-50.10	
Cu	3.46	0.178	5.50	
Eu	300.04	11.633✓	1.33	
Fe	1.69	0.026	13.22	
La	18.49	)54.571✓	1.03	
Pb	2.87	50.251✓	1.24	
Li	2.70	(-0.057	-13.49	
Mg	0.63	0.009	2.73	
Mn	0.81	0.003	2.84	
Hg	3.79	0.025	29.59	
Mo	1.53	(-0.022	-30.50	
Nd	4.07	(-1.500	-13.46	
Ni	2.96	(-0.037	-27.78	
P	1.48	0.075	21.93	
K	2.34	(-2.118	-11.26	
Tm	9.26	10.447✓	0.23	
Se	1.54	(-0.189	-45.12	
Si	3.39	0.337	7.66	
Ag	249.98	9.823✓	5.98	
Na	3.56	(-0.553	-16.43	
Sr	2.91	(-0.007	-26.54	
S	0.77	0.060	20.10	
Ta	3.08	-0.046	-74.08	
Tl	5.88	4.524✓	2.49	
Ih	14.25	52.414✓	0.89	
Sn	1.29	0.022	58.65	
Ti	3.16	0.007	84.41	
W	1.13	(-0.126	-24.82	
U	8.54	49.818✓	0.18	
V	5.67	0.175	2.11	
Zn	2.51	-0.000	-314.77	
Zr	9.95	-0.023	-88.48	

Figure 29 ICP Analysis March 12, 1990

Sample name : 77C11C  
Sample code 1 : SST3  
Sample code 2 : DIRECT  
Programme : SST

12-Mar-90 14:56:16

*Al*  
*14:56*

NAME	MV	INT	CONCEN	RSD
Al	26.11	153.231✓	1.20	
Si	0.44	0.584	18.31	
As	69.75	51.952✓	1.17	
Ba	2.67	(-0.026	-9.14	
Be	292.02	10.532✓	0.36	
Bi	3.85	0.118	40.39	
B	4.15	0.015	26.36	
Cd	1.93	(-0.012	-14.57	
Ca	0.81	0.012	1.84	
Cr	3.30	(-1.846	-5.58	
Co	1.07	(-0.160	-4.14	
Eu	0.23	-0.016	-57.73	
Cu	2.26	(-0.041	-15.07	
Eu	2.57	(-0.032	-6.17	
Fe	1.57	0.012	8.82	
La	0.31	(-0.128	-7.16	
Pb	0.26	-0.148	0.00	
Li	2.52	(-0.076	-7.00	
Mg	0.47	0.003	3.63	
Mn	0.88	0.009	5.08	
Hg	335.55	24.856✓	1.13	
Mo	323.15	50.861✓	1.20	
Nd	3.40	(-2.790	-3.25	
Ni	6.65	0.386	2.68	
P	90.09	49.952✓	1.09	
K	2.19	(-2.774	-5.69	
Sm	3.19	(-1.949	-5.99	
Se	27.65	51.643✓	1.18	
Si	77.51	42.068✓	1.77	
Ag	34.57	0.936	22.99	
Na	3.51	(-0.595	-6.87	
Sr	2.62	(-0.016	-6.88	
S	50.32	51.509✓	1.88	
Ts	123.44	46.682✓	1.35	
Tl	29.09	52.249✓	0.19	
Th	1.20	(-0.368	-20.39	
Sn	1.52	0.120	5.59	
Ti	578.70	50.962✓	1.36	
W	31.51	21.196✓	2.87	
U	4.41	3.349	31.44	
V	115.89	10.549✓	1.00	
Zn	3.17	0.022	5.73	
Zr	188.94	50.839	1.45	

Figure 30 ICP Analysis March 12, 1990

M  
T  
C  
I  
C  
S  
C  
E  
N  
—  
9

## APPENDIX A ANALYTICAL ANALYSIS CARDS

9 1 1 2 2 6 7 1 2 7 4

## Physical Analysis Sample 89-045

Serial No. F 101.-5000	Sample Point <b>SEGMENT-2</b>	Date 11-15-89	Time Issued 10:55	Priority 18
Determination APP/OTR	Method/Standard LI-000-200	Result Units NONE	Charge Code WB75L	Reruns 0
Sample Size ?	Customer ID 104C N 3133			Customer ID 89-045
Remarks, Calculations, Results: A. JAR ID#051 B. JAR TARE Wt. 221.45 C. JAR TOTAL WT. 356.40 D. C-B= 134.94 E. EST. VOL./LENGTH in inches of solids F. VISUAL REMARKS: Colors tan beige. Contamination on inside on cask and on outside of Vial body. SAMPLE Dark Brown grading to Medium Brown in color. Granular looking throughout. moist and sticky at bottom grading up to crumbly at top.				
Analyst-1 RLW	Analyst-2 A297	Analyst-3	Analyst-4	Analyst-5
KJT OTS	61300/16297			
Hrs	Hrs	Hrs	Hrs	Hrs
11-16-89				
Date	Time Completed	Lab Unit Mgr <i>CJM</i> SJS		

SA-8900-045 (R-10-82)

Serial No. F 101.-5003	Sample Point <b>SEGMENT-2</b>	Date 11-15-89	Time Issued 10:55	Priority 18
Determination HOMO97T	Method/Standard LI-000-200	Result Units NONE	Charge Code WB75L	Reruns 0
Sample Size ?	Customer ID 89-045			
Remarks, Calculations, Results: <i>Homogenization fact# 135</i> <i>Complete 12-29-89</i> 104C N 3134 pg 13				
Analyst-1 RLW	Analyst-2 A297	Analyst-3	Analyst-4	Analyst-5
KJT OTS	61300/16297			
Hrs	Hrs	Hrs	Hrs	Hrs
12-29-89				
Date	Time Completed	Lab Unit Mgr <i>BJM</i> S		

SA-8900-045 (R-10-82)

Serial No. F 101.-5001	Sample Point <b>SEGMENT-2</b>	Date 11-15-89	Time Issued 10:55	Priority 18
Determination VOA SAMPL	Method/Standard LI-000-200	Result Units NONE	Charge Code WB75L	Reruns 0
Sample Size ?	Customer ID 2.5%			Customer ID 89-045
Remarks, Calculations, Results: DUPLICATE SAMPLE Total wt. 25.0 g Tare wt. 22.49 g Sample wt. 2.5 g <i>Sent to VNL</i> 104C N 313 3				
Analyst-1 RLW	Analyst-2	Analyst-3	Analyst-4	Analyst-5
KJT OTS				
Hrs	Hrs	Hrs	Hrs	Hrs
11-16-89				
Date	Time Completed	Lab Unit Mgr <i>JRW</i> S		

SA-8900-045 (R-10-82)

9 1 1 2 9 6 7 1 2 5 5

## pH Analysis of Solid Sample

Serial No. F 100.-5515	Sample Point SEGMENT-1		Date 11-15-89	Time Issued 10:55	Priority 19
Determination pH	Method/Standard LA-212-103	Result Units % RECOVERY	Charge Code WB75L	Remarks, Calculations, Results: LMCS CHECK SAMPLE OH FOUND <u>10.10</u> STD ID <u>12011-A</u> SAMPLE TEMP <u>23.6</u>	
Sample Size ?			Customer ID <u>089045</u>		
<i>10.10 / 1000 101.0%</i>					
Analyst - 1 <u>GC269</u>	Analyst - 2 <u>Mary Tracy</u>	Analyst - 3	Analyst - 4	Analyst - 5 <u>REB</u>	
Hrs	Hrs	Hrs	Hrs		
Date 1-2-90	Time Completed <u>CJA</u>	Lab Unit Mgr <u>bff</u>			
54-6800-061 (R-10-83)					

Serial No. F 102.-5115	Sample Point SEGMENT-3		Date 11-15-89	Time Issued 10:55	Priority 19
Determination pH	Method/Standard LA-212-103	Result Units NONE	Charge Code WB75L	Remarks, Calculations, Results: pH <u>12.02</u> SAMPLE TEMP <u>23.7</u>	
Sample Size ?			Customer ID <u>089045</u>		
<i>2.581 g / 2.581 ml</i>					
Analyst - 1 <u>GC269</u>	Analyst - 2 <u>Mary Tracy</u>	Analyst - 3	Analyst - 4	Analyst - 5 <u>REB</u>	
Hrs	Hrs	Hrs	Hrs		
Date 1-2-90	Time Completed <u>CJA</u>	Lab Unit Mgr <u>bff</u>			
54-6800-061 (R-10-83)					

Serial No. F 101.-5015	Sample Point SEGMENT-2		Date 11-15-89	Time Issued 10:55	Priority 19
Determination pH	Method/Standard LA-212-103	Result Units NONE	Charge Code WB75L	Remarks, Calculations, Results: pH <u>11.96</u> SAMPLE TEMP <u>23.7</u>	
Sample Size ?			Customer ID <u>089045</u>		
<i>2.5000 g / 2.5000 ml #216 8.13 grams</i>					
Analyst - 1 <u>GC269</u>	Analyst - 2 <u>Mary Tracy</u>	Analyst - 3	Analyst - 4	Analyst - 5 <u>REB</u>	
Hrs	Hrs	Hrs	Hrs		
Date 1-2-90	Time Completed <u>CJA</u>	Lab Unit Mgr <u>bff</u>			
54-6800-061 (R-10-83)					

Serial No. F 292.-5515	Sample Point SEGMENT-E		Date 11-21-89	Time Issued 8:30	Priority 19
Determination pH	Method/Standard LA-212-103	Result Units % RECOVERY	Charge Code WB75L	Remarks, Calculations, Results: LMCS CHECK SAMPLE OH FOUND <u>10.09</u> STD ID <u>12011-A</u> SAMPLE TEMP <u>23.8</u>	
Sample Size ?			Customer ID <u>089050</u>		
<i>10.09 / 10.0 (100.9%)</i>					
Analyst - 1 <u>GC269</u>	Analyst - 2 <u>Mary Tracy</u>	Analyst - 3	Analyst - 4	Analyst - 5 <u>REB</u>	
Hrs	Hrs	Hrs	Hrs		
Date 1-2-90	Time Completed <u>CJA</u>	Lab Unit Mgr <u>bff</u>			
54-6800-061 (R-10-83)					

9 1 1 2 0 6 0 1 2 7 6

## pH Analysis of Solid Sample

Serial No. • F 121-5315	Sample Point SEGMENT-22	Date 11-15-89	Time Issued 10:58	Priority 16
Determination pH	Method/Standard LA-212-103	Result Units NONE	Charge Code WB75L	Reurns O
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results: LMC'S CHECK SAMPLE pH FOUND _____ STD AD SAMPLE TEMP _____ <i>Reagent Blank = 6.83</i>				
Analyst - 1 <i>GC269</i>	Analyst - 2 <i>Mary</i>	Analyst - 3 <i>Mary</i>	Analyst - 4 <i>b1</i>	Analyst - 5 <i>RE</i>
Date 1-2-90	Time Completed <i>8pm</i>	Lab Unit/Mgr <i>b1</i>		

84-0800-061 (P-10-62)

Percent Water Analysis

9 1 1 2 2 5 7 1 2 7 7

<i>Really F940 See Comp #8</i>				
Serial No. F 309.-5310	Sample Point SEGMENT-V	Date 11-21-89	Time issued 8:32	Priority 18
Determination % H <sub>2</sub> O	Method/Standard LA-564-101	Result Units <i>(X)</i>	Charge Code WB75L	Runno 0
Sample Size ?	Customer ID <b>089050</b>			
Remarks, Calculations, Results:  REAGENT BLANK 0 22.0121 G 21.4901 22.0121 T 21.4901 22.0054 W1 21.4836 22.0051 W2 21.4836  Analyst-1 Analyst-2 Analyst-3 Analyst-4 Analyst-5 <i>63598/RH</i> <i>  </i> <i>  </i> <i>  </i> <i>  </i> Hrs      Hrs      Hrs      Hrs      Hrs  Date 1-3-90 Time Completed <i>Craig Siedle</i> Lab Unit Mgr				
54-8600-061 (R-10-83)				

<i>Really F940 See Comp #8</i>				
Serial No. F 101.-5010	Sample Point SEGMENT-2	Date 11-15-89	Time issued 10:55	Priority 19
Determination % H <sub>2</sub> O	Method/Standard LA-564-101	Result Units %	Charge Code WB75L	Runno 0
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results:  #225 5.86 grams G 22.7747 T 21.7273 W1 22.3713 W2 22.3715 10HC N 3134  Analyst-1 Analyst-2 Analyst-3 Analyst-4 Analyst-5 <i>63598/RH</i> <i>  </i> <i>  </i> <i>  </i> <i>  </i> Hrs      Hrs      Hrs      Hrs      Hrs  Date 1-3-90 Time Completed <i>Craig</i> Lab Unit Mgr				
54-8600-061 (R-10-83)				

<i>Really F939 See Comp</i>				
Serial No. F 100.-5510	Sample Point SEGMENT-1	Date 11-15-89	Time issued 10:54	Priority 19
Determination % H <sub>2</sub> O	Method/Standard LA-564-101	Result Units % RECOVERY	Charge Code WB75L	Runno 0
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results:  LMCS CHECK SAMPLE LMCS ID <i>11C11A6</i> G 23.1834 G 23.2273 57.90 96.63% T 21.8182 T 21.8485 <i>today</i> W1 22.3970 W2 22.4333 57.60 59.61 W2 22.3860 W2 22.4227  Analyst-1 Analyst-2 Analyst-3 Analyst-4 Analyst-5 <i>63598/RH</i> <i>  </i> <i>  </i> <i>  </i> <i>  </i> Hrs      Hrs      Hrs      Hrs      Hrs  Date 1-3-90 Time Completed <i>Craig</i> Lab Unit Mgr				
54-8600-061 (R-10-83)				

<i>Really F939 See Comp</i>				
Serial No. F 292.-5510	Sample Point SEGMENT-E	Date 11-21-89	Time issued 8:29	Priority 19
Determination % H <sub>2</sub> O	Method/Standard LA-564-101	Result Units % RECOVERY	Charge Code WB75L	Runno 0
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results:  LMCS CHECK SAMPLE LMCS ID <i>11C11A6</i> G 22.9125 G 23.2933 58.296 57.20% T 21.5294 T 21.9143 W1 22.1071 W2 22.5044 57.70% W2 22.1011 W2 22.4941 59.61  Analyst-1 Analyst-2 Analyst-3 Analyst-4 Analyst-5 <i>63598/RH</i> <i>  </i> <i>  </i> <i>  </i> <i>  </i> Hrs      Hrs      Hrs      Hrs      Hrs  Date 1-3-90 Time Completed <i>Craig</i> Lab Unit Mgr				
54-8600-061 (R-10-83)				

9 1 1 2 0 6 0 1 2 5 3

## Percent Water Analysis

Serial No. F 102-5110	Sample Point SEGMENT-3	Date 11-15-89	Time Issued 10:55	Priority 19
Determination % H <sub>2</sub> O	Method/Standard LA-564-101	Result Units %	Charge Code WB75L	Reruns 0
Sample Size ?				Customer ID <b>089045</b>
Remarks, Calculations, Results: <b>DUPLICATE SAMPLE</b>				
G 22.8091	T. 21.6142	W1 22.3629	W2 22.3623	37.38
Analyst - 1 <i>69598/21</i>	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 <i>R E Lando</i> Hrs
Date 1-3-90	Time Completed	Lab Unit/Mgr <i>Cope</i>	<i>11</i>	

84-9900-081 (R-10-82)

9 1 1 2 3 5 9 1 2 7 9

## Fusion Dissolution

Serial No. <b>F 106.-6000</b>	Sample Point <b>SEGMENT-7</b>	Date <b>11-15-89</b>	Time Issued <b>10:55</b>	Priority <b>18</b>
Determination <b>FUSION</b>	Method/Standard <b>LA-549-141</b>	Result Units <b>G/L</b>	Charge Code <b>WB75L</b>	Returns <b>0</b>
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results: GRAMS SAMPLE <u>.5995</u> VOLUME ON COMPLETION <u>250 ml</u>  <u>2.40 -3 g/ml</u> <del><u>2.40 -3 g/ml</u></del> <u>2.40 -3 g/l</u>				
Analyst-1 <i>63598/90</i>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>cpe</i>
Hrs	Hrs	Hrs	Hrs	Hrs
Date <b>1-3-90</b>	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>DL</i>	
54-6000-061 (R-10-82)				

Serial No. <b>F 168.-6300</b>	Sample Point <b>SEGMENT-U</b>	Date <b>11-17-89</b>	Time Issued <b>10:18</b>	Priority <b>18</b>
Determination <b>FUSION</b>	Method/Standard <b>LA-549-141</b>	Result Units <b>G/L</b>	Charge Code <b>WB75L</b>	Returns <b>0</b>
Sample Size ?	Customer ID <b>089048</b>			
Remarks, Calculations, Results: REAGENT BLANK  <i>Complete</i>				
Analyst-1 <i>63598/90</i>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>cpe</i>
Hrs	Hrs	Hrs	Hrs	Hrs
Date <b>1-3-90</b>	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>DL</i>	
54-6000-061 (R-10-82)				

Serial No. <b>F 107.-6100</b>	Sample Point <b>SEGMENT-8</b>	Date <b>11-15-89</b>	Time Issued <b>10:55</b>	Priority <b>18</b>
Determination <b>FUSION</b>	Method/Standard <b>LA-549-141</b>	Result Units <b>G/L</b>	Charge Code <b>WB75L</b>	Returns <b>0</b>
Sample Size ?	Customer ID <b>089045</b>			
Remarks, Calculations, Results: DUPLICATE ANALYSIS #2 GRAMS SAMPLE <u>.5731</u> VOLUME ON COMPLETION <u>250 ml</u>  <u>2.29 -3 g/ml</u> <del><u>2.29 -3 g/ml</u></del> <u>2.29 -3 g/l</u> <i>NH4 NO3/T</i>				
Analyst-1 <i>63598/90</i>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>cpe</i>
Hrs	Hrs	Hrs	Hrs	Hrs
Date <b>1-3-90</b>	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>DL</i>	
54-6000-061 (R-10-82)				

Total Alpha Analysis on the Fusion Dissolution

9 1 1 2 7 6 7 1 3 : 0

Serial No. F 107.-6120	Sample Point SEGMENT-8		Date 11-15-89	Time Issued 10:55	Priority 19
Determination AT	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Runno. 0	
Sample Size ? 100-10-100			Customer ID 089045		
Remarks, Calculations, Results: DUPLICATE SAMPLE  <i>&lt;200<sup>-2</sup></i> <i>J. Hopkins</i>					
Analyst-1 6A543	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Cja</i>	Q.M.	G.M.	
54-8800-061 (R-10-63)					

Serial No. F 105.-6520	Sample Point SEGMENT-6		Date 11-15-89	Time Issued 10:55	Priority 19
Determination AT	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WB75L	Runno. 0	
Sample Size ? 10ml 508			Customer ID 089045		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID _____					
Analyst-1 6A543	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Cja</i>	Q.M.	G.M.	
54-8800-061 (R-10-63)					

Serial No. F 106.-6020	Sample Point SEGMENT-7		Date 11-15-89	Time Issued 10:55	Priority 19
Determination AT	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Runno. 0	
Sample Size ? 100-10-100			Customer ID 089045		
Remarks, Calculations, Results:					
Analyst-1 6A543	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 1-5-90	Time Completed	Lab Unit Mgr			
54-8800-061 (R-10-63)					

Serial No. F 296.-6220	Sample Point SEGMENT-1		Date 11-21-89	Time Issued 8:30	Priority 19
Determination AT	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WB75L	Runno. 0	
Sample Size 508			Customer ID 089050		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 83044 SPIKE VOLUME 10ml 97.1% 9.71% / 1.000 98.0% 1.000 / 1.000 9.901 <sup>-1</sup> - 1.935 <sup>-2</sup> = 9.708 <sup>-1</sup> / 100 See back of Card <i>Cja</i> <i>0.000</i> →					
Analyst-1 6A543	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Cja</i>	Q.M.	G.M.	
54-8800-061 (R-10-63)					

Total Alpha Analysis on the Fusion Dissolution

$$\frac{18}{2} \frac{[(9.90)^1](.0001) - (1.935)^2(.0001)(100)}{(1.010)(1.001)^2} = \underline{\underline{92.8870}} = 97.0\% \text{ alpha}$$

482  
10 - .3  
Alpha Calculation by AJ on 01-05-1990 at 18:21:53  
Net #18 2-inch mount Alpha eff. : .2095  
Sample size : .1 mL Dilution : 1

Mount # 1

$$\frac{482}{10} - 0.3 = 1.0299E+00 \mu\text{Ci/L alpha}$$

Mount # 2

$$\frac{445}{10} - 0.3 = 4.5035E-01 \mu\text{Ci/L alpha}$$

F 296.-6220 AT

18 2" 1-6-90 VR

502  
10 - .3  
Alpha Calculation by VR on 01-06-1990 at 00:41:44  
Net #18 2-inch mount Alpha eff. : .2095  
Sample size : 10 mL Dilution : 1

Mount # 1

$$\frac{500}{10} - 0.3 = 1.0696E-02 \mu\text{Ci/L alpha}$$

Mount # 2

$$\frac{547}{10} - 0.3 = 1.3697E-02 \mu\text{Ci/L alpha}$$

F 105.-6520

18/2

17  
10 - 3

Alpha Calculation by AJ on 01-05-1990 at 18:17:07  
Net #18 2-inch mount Alpha eff. : .2095  
Sample size : .1 mL Dilution : 101

Mount # 1

$$\frac{17}{10} - 0.3 = 3.0403E+00 \mu\text{Ci/L alpha}$$

Mount # 2

$$\frac{9}{10} - 0.3 = 1.3030E+00 \mu\text{Ci/L alpha}$$

F 106.-6020

18 2" 1-6-90 VR

7  
10 - .3  
2 - 10

Alpha Calculation by VR on 01-06-1990 at 00:47:16  
Net #18 2-inch mount Alpha eff. : .2095  
Sample size : 1 mL Dilution : 1010

Mount # 1

$$\frac{7}{10} - 0.3 = 2.0017E+00 \mu\text{Ci/L alpha}$$

Mount # 2

$$\frac{7}{10} - 0.3 = 2.0017E+00 \mu\text{Ci/L alpha}$$

F 107.-6120

9 1 1 2 0 6 7 1 2 1 2

## Total Alpha Analysis on the Fusion Dissolution

Serial No.	Sample Point	Date	Time Issued	Priority
F 308.-6320	SEGMENT-U	11-21-89	8:32	18
Determination	Method/Standard	Result Units	Charge Code	Permit
AT	LA-548-101	uCi/L	WB75L	0
Sample Size				Customer ID
? 16ml				089050
Remarks, Calculations, Results: REAGENT BLANK  4.100 <sup>-4</sup> mili				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
6A543 Hrs D. Hopkins 1-5-90	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	<i>Cyan</i>	<i>GM</i>
54-6800-081 (R-10-83)				

Serial No.	Sample Point	Date	Time Issued	Priority
F 297.-6520	SEGMENT-J	11-21-89	8:30	19
Determination	Method/Standard	Result Units	Charge Code	Permit
AT	LA-548-101	% RECOVERY	WB75L	0
Sample Size				Customer ID
? 10ml				089050
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID _____  100.3% 1.003 <sup>-2</sup> / 1.0001 <sup>-2</sup>				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
6A543 Hrs D. Hopkins 1-5-90	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	<i>Cyan</i>	<i>GM</i>
54-6800-081 (R-10-83)				

Total Alpha Analysis on the Fusion Dissolution

8 2" 1-6-90 VR

473

466  
10

Alpha Calculation by VR on 01-06-1990 at 00:46:14  
Det #18 2-inch mount Alpha eff.: .2095  
Sample size: 10 mL Dilution: 1

Mount # 1

473  
-----  
10

0.3 = 1.0104E-02 uCi/L alpha

Mount # 2

466  
-----  
10

0.3 = 9.9551E-03 uCi/L alpha

297.-6520 AT

18 2" 1-6-90 VR

4  
10  
1  
10

Alpha Calculation by VR on 01-06-1990 at 02:38:47  
Det #18 2-inch mount Alpha eff.: .2095  
Sample size: 1 mL Dilution: 1

Mount # 1

4  
-----  
10

0.4 < 1.0034E-06 uCi/mL alpha

100

Mount # 2

1  
-----  
10

0.4 < 1.0020E-06 uCi/mL alpha

100

F 308.-6320 AT

Total Alpha Analysis on the Fusion Dissolution

9 1 1 2 7 6 7 1 2 3 4

Serial No.	Sample Point	Date	Time Issued	Priority
F 108.-6220	SEGMENT-I	11-15-89	10:56	26
Determination	Method/Standard	Result Units	Charge Code	Reruns
AT	LA-508-101	% RECOVERY	F215C	0
Sample Size	Customer ID			
? 100-10-200				
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 102B44 SPIKE VOLUME 10ml				
To follow				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
GC267				
Mandy Tracy	Hrs	Hrs	Hrs	Hrs
1-25-90	Time Completed	Lab Unit Mgr.	DM S	

54-000-081 (R-10-83)

Serial No.	Sample Point	Date	Time Issued	Priority
F 109.-6520	SEGMENT-J	11-15-89	10:56	26
Determination	Method/Standard	Result Units	Charge Code	Reruns
AT	LA-508-101	% RECOVERY	F215C	0
Sample Size	Customer ID			
? 10 ml				
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 102B44				
108.8%				
<del>1.09<sup>-2</sup></del> <del>1.0022<sup>-2</sup></del>				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
GC267				
Mandy Tracy	Hrs	Hrs	Hrs	Hrs
1-25-90	Time Completed	Lab Unit Mgr.	DM S	

54-000-081 (R-10-83)

Serial No.	Sample Point	Date	Time Issued	Priority
F 105.-6520	SEGMENT-6	11-15-89	10:55	19
Determination	Method/Standard	Result Units	Charge Code	Reruns
AT	LA-548-101	% RECOVERY	WB75L	1
Sample Size	Customer ID			
? 10 ml				
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 102B44				
113.7%				
<del>1.14<sup>-2</sup></del> <del>1.0022<sup>-2</sup></del>				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
GC267				
Mandy Tracy	Hrs	Hrs	Hrs	Hrs
1-25-90	Time Completed	Lab Unit Mgr.	DM S	

54-000-081 (R-10-83)

Serial No.	Sample Point	Date	Time Issued	Priority
F 106.-6020	SEGMENT-7	11-15-89	10:55	19
Determination	Method/Standard	Result Units	Charge Code	Reruns
AT	LA-548-101	uCi/L	WB75L	0
Sample Size	Customer ID			
? 100-10-200	089045			
Remarks, Calculations, Results:				
RERUN 1.024 will				
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
GC267				
Mandy Tracy	Hrs	Hrs	Hrs	Hrs
1-25-90	Time Completed	Lab Unit Mgr.	DM S	

54-000-081 (R-10-83)

Total Alpha Analysis on the Fusion Dissolution

<p><b>18 2" 1-25-80 VR</b></p> <p><del>10</del> <del>10</del> <del>10</del> <del>10</del> - .5</p> <p>Alpha Calculation by VR on 01-25-1990 at 18:40:25 Det #18 2-inch mount Alpha eff. : .2095 Sample size : 1 mL Dilution : 505</p> <p>Mount # 1</p> <p><del>10</del> <del>10</del> - 0.5 = 1.2368E+00 uCi/L alpha</p> <p>Mount # 2</p> <p><del>10</del> <del>10</del> - 0.5 = 1.2368E+00 uCi/L alpha</p>	<p><b>18 2" 1-25-90 VR</b></p> <p><del>540</del> <del>10</del> <del>527</del> <del>10</del> - .5</p> <p>Alpha Calculation by VR on 01-25-1990 at 19:14:10 Det #18 2-inch mount Alpha eff. : .2095 Sample size : 10 mL Dilution : 1</p> <p>Mount # 1</p> <p><del>540</del> <del>10</del> - 0.5 = 1.1503E-02 uCi/L alpha</p> <p>Mount # 2</p> <p><del>527</del> <del>10</del> - 0.5 = 1.1224E-02 uCi/L alpha</p>
<p>F 106.-6020</p>	<p>F 105.-6520</p>
<p><b>18 2" 1-25-90 VR</b></p> <p><del>538</del> <del>10</del> <del>429</del> <del>10</del> - .5</p> <p>Alpha Calculation by VR on 01-25-1990 at 18:43:00 Det #18 2-inch mount Alpha eff. : .2095 Sample size : 10 mL Dilution : 1</p> <p>Mount # 1</p> <p><del>538</del> <del>10</del> - 0.5 = 1.1460E-02 uCi/L alpha</p> <p>Mount # 2</p> <p><del>489</del> <del>10</del> - 0.5 = 1.0407E-02 uCi/L alpha</p>	<p><b>18 2" 1-25-90 VR</b></p> <p><del>506</del> <del>10</del> <del>551</del> <del>10</del> - .5</p> <p>Alpha Calculation by VR on 01-25-1990 at 18:41:18 Det #18 2-inch mount Alpha eff. : .2095 Sample size : 1 mL Dilution : 505</p> <p>Mount # 1</p> <p><del>506</del> <del>10</del> - 0.5 = 5.4399E+01 uCi/L alpha</p> <p>Mount # 2</p> <p><del>551</del> <del>10</del> - 0.5 = 5.7285E+01 uCi/L alpha</p>
<p>F 109.-6520</p>	<p>F 108.-6220</p>

9 1 1 2 2 6 9 1 2 5 6

## Total Alpha Analysis on the Fusion Dissolution

Serial No.	Sample Point		Date	Time Issued	Priority
F 120.-6320	SEGMENT-II		11-15-89	10:58	26
Determination	Method/Standard	Result Units	Charge Code	Returns	
AT	LA-508-101	uCi/L	F215C	0	
Sample Size	<i>? 10ml</i>				Customer ID
Remarks, Calculations, Results: REAGENT BLANK					
Analyst - 1 <i>GC269</i>	Analyst - 2 <i>Mary</i>	Analyst - 3 <i>Mary</i>	Analyst - 4 <i>John</i>	Analyst - 5 <i>John</i>	
Date 1-25-90	Time Completed	Lab Unit No. <i>CPL</i>		Comments <i>None</i>	

54-0000-061 (R-10-62)

9 1 1 2 2 5 2 1 2 5 7

Total Alpha Analysis on the Fusion Dissolution

Route # 1  
 Alpha Factor 7  
 Beta 18  
 Sample size 7  
 - - -  
 - - - 5  
 - - - - -  
 10 10 10  
 Route # 2  
 10  
 5

10	Sample size : 10 ml	Dilution : 1
10	0.5 ( 2.4591E-04 MCL alpha	
2	Count # 2	
10	0.5 ( 1.1219E-04 MCL alpha	
10	Sample size : 10 ml	Dilution : 1

Total Beta Analysis on the Fusion Dissolution

9 1 1 2 0 6 7 1 2 3 3

7-3 F948 Seg Comp 8

Serial No. F 308.-6325	Sample Point SEGMENT-U	Date 11-21-89	Time Issued 8:32	Priority 1B
Determination TB	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Retuns 0
Sample Size ? 10ml	Customer ID 089050			
Remarks, Calculations, Results: REAGENT BLANK				
<i>&lt;2.58 H</i> <i>well</i>				
Analyst -1 6A543	Analyst -2	Analyst -3	Analyst -4	Analyst -5
Hrs	Hrs	Hrs	Hrs	Hrs
<i>J. Hopkins</i>				
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Carter</i>	g/m <sup>2</sup>	
64-000-061 (R-10-82)				

7-8 F949 Seg Comp 8

Serial No. F 296.-6225	Sample Point SEGMENT-I	Date 11-21-89	Time Issued 8:30	Priority 19
Determination TB	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WB75L	Retuns 0
Sample Size 508	Customer ID 089050			
Sample Size ? 1001	Customer ID 089050			
Remarks, Calculations, Results: SPIKE SAMPLE F294 SPIKE ID 83844 SPIKE VOLUME 10ml				
<i>2.154' - 5.493 = 1.60'</i> <i>100 - 1.60' / 1.390' = 100</i> <i>See back of card</i>				
Analyst -1 6A543	Analyst -2	Analyst -3	Analyst -4	Analyst -5
Hrs	Hrs	Hrs	Hrs	Hrs
<i>J. Hopkins</i>				
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Carter</i>	g/m <sup>2</sup>	
64-000-061 (R-10-82)				

9 9 10

Serial No. F 105.-6525	Sample Point SEGMENT-6	Date 11-15-89	Time Issued 10:55	Priority 19
Determination TB	Method/Standard LA-548-101	Result Units % RECOVERY	Charge Code WB75L	Retuns 0
Sample Size ? 10 ml	Customer ID 089045 83844 End C			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 83844				
<i>1.376' / 1.3923' = 98.8%</i>				
Analyst -1 6A543	Analyst -2	Analyst -3	Analyst -4	Analyst -5
Hrs	Hrs	Hrs	Hrs	Hrs
<i>J. Hopkins</i>				
Date 1-5-90	Time Completed	Lab Unit Mgr <i>Carter</i>	g/m <sup>2</sup>	
64-000-061 (R-10-82)				

11 2

Serial No. F 106.-6025	Sample Point SEGMENT-7	Date 11-15-89	Time Issued 10:55	Priority 19
Determination TB	Method/Standard LA-548-101	Result Units uCi/L	Charge Code WB75L	Retuns 0
Sample Size ? 100-10-100	Customer ID 089045			
Remarks, Calculations, Results: <i>Bad checker</i>				
Analyst -1 6A543	Analyst -2	Analyst -3	Analyst -4	Analyst -5
Hrs	Hrs	Hrs	Hrs	Hrs
<i>J. Hopkins</i>				
Date 1-5-90	Time Completed	Lab Unit Mgr		
64-000-061 (R-10-82)				

9 1 1 2 0 6 0 1 2 1 9

## Total Beta Analysis on the Fusion Dissolution

1x 2" 1-6-90 VR

$$\frac{56}{10} - 6$$

$$\frac{50}{10} -$$

Beta Calculation by VR on 01-06-1990 at 02:38:45  
 Det #18 2-inch count Beta eff. : .3151  
 Sample size : 1 mL Dilution : 1

Mount # 1

$$\frac{56}{10} - 6.0 < 2.5801E-06 \mu\text{Ci}/\text{mL beta}$$

Mount # 2

$$\frac{50}{10} - 6.0 < 2.5801E-06 \mu\text{Ci}/\text{mL beta}$$

F 308.-6325 TB

16/2

Spots too low to calculate

$$\frac{15124}{10} - 6$$

Beta Calculation by AJ on 01-05-1990 at 18:21:51  
 Det #18 2-inch mount Beta eff. : .3151  
 Sample size : .1 mL Dilution : 1

$$\frac{14892}{10}$$

Mount # 1

$$\frac{15124}{10} - 6.0 < 2.1535E+01 \mu\text{Ci}/\text{L beta}$$

Mount # 2

$$\frac{14892}{10} - 6.0 < 2.1205E+01 \mu\text{Ci}/\text{L beta}$$

F 296.-6225 TB

16 2" 1-6-90 VR

$$\frac{9487}{10} -$$

$$\frac{9879}{10} -$$

Beta Calculation by VR on 01-06-1990 at 00:43:41  
 Det #18 2-inch count Beta eff. : .3151  
 Sample size : 10 mL Dilution : 1

Mount # 1

$$\frac{9489}{10} - 6.0 = 1.3479E-01 \mu\text{Ci}/\text{L beta}$$

Mount # 2

$$\frac{9879}{10} - 6.0 = 1.4057E-01 \mu\text{Ci}/\text{L beta}$$

PF 101

,100

18/2

$$\frac{3660}{10} - 6$$

Beta Calculation by AJ on 01-05-1990 at 18:17:05  
 Det #18 2-inch mount Beta eff. : .3151  
 Sample size : .1 mL Dilution : 101

$$\frac{1595}{10}$$

Mount # 1

$$\frac{3660}{10} - 6.0 = 5.1978E+02 \mu\text{Ci}/\text{L beta}$$

Mount # 2

$$\frac{1595}{10} - 6.0 = 2.2163E+02 \mu\text{Ci}/\text{L beta}$$

F 106.-6025

9 1 1 2 3 6 7 1 2 7 0

## Total Beta Analysis on the Fusion Dissolution

3 X 4				
Serial No. F 107-6125	Sample Point SEGMENT-8	Date 11-15-89	Time Issued 10:56	Priority 19
Determination TB	Method/Standard LA-5100-101 <i>50B</i>	Result Units uCi/L	Charge Code WB75L	Refers 0
Sample Size ? 100 -10-100				Customer ID 089045
Remarks, Calculations, Results: DUPLICATE SAMPLE				
<i>1.97<sup>2</sup> mif</i>				
Analyst - 1 <i>6A543</i> <i>J. Hopkins</i>	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 Hrs
Date 1-5-90	Time Completed	Lab Unit Mon <i>Copier</i>	Comments <i>YAM's</i>	

SA-0800-083 (R-10-82)

9 1 1 2 2 5 7 1 2 7 1

Total Beta Analysis on the Fusion Dissolution

1.01 ± 3

18 2" 1-6-90 uR

$$\begin{array}{r} 1432 \\ - 10 \\ \hline 1421 \\ - 10 \\ \hline \end{array}$$

Beta Calculation by VR on 01-04-1990 at 00:47:14  
Det #18 2-inch mount Beta eff. : .3151  
Sample size : 1 mL Dilution : 1010

Mount # 1

$$\begin{array}{r} 1432 \\ - 10 \\ \hline \end{array} = 6.0 = 1.9810E+02 \mu\text{Ci/L beta}$$

Mount # 2

$$\begin{array}{r} 1421 \\ - 10 \\ \hline \end{array} = 6.0 = 1.9651E+02 \mu\text{Ci/L beta}$$

F 107.-6125

9 1 1 2 2 6 7 1 2 7 2

## Total Beta Analysis on the Fusion Dissolution

152 F950 Soc Com S

Serial No. F 297.-6525	Sample Point SEGMENT-J	Date 11-21-89	Time Issued 8:30	Priority 19
Determination TB	Method/Standard LA-54B-101	Result Units % RECOVERY	Charge Code WB75L	Refuge 0
Sample Size ? 10 ml	Customer ID 83844 089050			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID _____				
Analyst - 1 68543	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
Hrs	Hrs	Hrs	Hrs	Hrs
Q. Hopkins	Date 1-5-90	Time Completed	Lab Unit Mgr C. Taylor	AM 8

1.344<sup>-1</sup>  
1.39231  
96.5%

34-9800-061 (R-10-63)

Total Beta Analysis on the Fusion Dissolution

18 2" 1-6-90 VR

9410 - 6  
10  
9517 -  
10

Beta Calculation by VR on 01-06-1990 at 00:46:12  
Bet #18 2-inch mount Beta eff. : .3151  
Sample size : 10 mL Dilution : 1

Mount # 1

9410 ----- 6.0 = 1.3366E-01 uCi/L beta  
10

Mount # 2

9517 ----- 6.0 = 1.3519E-01 uCi/L beta  
10

F 297.-6525 TB

2 7 1 1 2 3 4 5 6

Total Beta Analysis on the Fusion Dissolution

9 1 1 2 2 5 7 1 2 7 4

9-10

Serial No.	Sample Point	Date	Time Issued	Priority
F 106.-6025	SEGMENT-7	11-15-89	10:55	19
Determination	Method/Standard	Result Units	Charge Code	Reruns
TB	LA-508-101	uCi/L	WB75L	0
Sample Size	Customer ID			
? 100-10-200	089045			
Remarks, Calculations, Results:				
<p><i>2.03 well</i></p> <p><b>RERUN</b></p>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
GC269 Crauny Frank	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	RMS	
1-25-90		<i>CRA</i>	RMS	

54-8800-081 (R-10-82)

3-4

Serial No.	Sample Point	Date	Time Issued	Priority
F 105.-6525	SEGMENT-6	11-15-89	10:55	19
Determination	Method/Standard	Result Units	Charge Code	Reruns
TB	LA-548-101	% RECOVERY	WB75L	1
Sample Size	Customer ID			
? 10ml				
Remarks, Calculations, Results:				
<p>LMCS CHECK SAMPLE</p> <p>LMCS ID 102.B44</p>				
<p><i>1.457<sup>-1</sup></i></p> <p><i>1.3851<sup>-1</sup></i></p> <p><i>105.20</i></p>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
GC269 Crauny Frank	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	RMS	
1-25-90		<i>CRA</i>	RMS	

54-8800-081 (R-10-82)

7-3

Serial No.	Sample Point	Date	Time Issued	Priority
F 108.-6225	SEGMENT-I	11-15-89	10:56	26
Determination	Method/Standard	Result Units	Charge Code	Reruns
TB	LA-508-101	% RECOVERY	F215C	0
Sample Size	Customer ID			
? 100-10-200				
Remarks, Calculations, Results:				
<p>SPIKE SAMPLE</p> <p>SPIKE ID 102.B44</p> <p>SPIKE VOLUME 10ml</p>				
<p><i>To Low</i></p>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
GC269 Crauny Frank	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	RMS	
1-25-90		<i>CRA</i>	RMS	

54-8800-081 (R-10-83)

5-6

Serial No.	Sample Point	Date	Time Issued	Priority
F 120.-6325	SEGMENT-U	11-15-89	10:58	26
Determination	Method/Standard	Result Units	Charge Code	Reruns
TB	LA-508-101	uCi/L	F215C	0
Sample Size	Customer ID			
? 10ml				
Remarks, Calculations, Results:				
<p>REAGENT BLANK</p>				
<p><i>6.65<sup>-4</sup></i></p>				
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5
GC269 Crauny Frank	Hrs	Hrs	Hrs	Hrs
Date	Time Completed	Lab Unit Mgr	RMS	
1-25-90		<i>CRA</i>	RMS	

54-8800-081 (R-10-83)

Total Beta Analysis on the Fusion Dissolution

9   1   2   2   5   ?   2   7   5								
18 2" 1-25-90 VR	D.F. 505	18 2" 1-25-90 VR						
<u>2922</u> - 13		<u>10362</u> - 13						
<u>2958</u> 10		<u>10276</u> 10						
Beta Calculation by VR on 01-25-1990 at 18:40:22 Det #18 2-inch mount Beta eff. : .3151 Sample size : 1 uL Dilution : 505		Beta Calculation by VR on 01-25-1990 at 19:13:07 Det #18 2-inch mount Beta eff. : .3151 Sample size : 10 uL Dilution : 1						
Mount # 1		Mount # 1						
2922 10	- 13.0 = 2.0156E+02 uCi/L beta	10362 10	- 13.0 = 1.4627E-01 uCi/L beta					
Mount # 2		Mount # 2						
2958 10	- 13.0 = 2.0416E+02 uCi/L beta	10276 10	- 13.0 = 1.4504E-01 uCi/L beta					
	F 106.-6025							F 105.-6525
18 2" 1-25-90 VR	D.F. 505	18 2" 1-25-90 VR						
<u>12652</u> - 13		<u>187</u> - 13						
<u>13054</u> 10		<u>139</u> 10						
Beta Calculation by VR on 01-25-1990 at 18:41:16 Det #18 2-inch mount Beta eff. : .3151 Sample size : 1 uL Dilution : 505		Beta Calculation by VR on 01-25-1990 at 18:42:12 Det #18 2-inch mount Beta eff. : .3151 Sample size : 10 uL Dilution : 1						
Mount # 1		Mount # 1						
12652 10	- 13.0 = 9.0399E+02 uCi/L beta	187 10	- 13.0 = 8.1484E-04 uCi/L beta					
Mount # 2		Mount # 2						
13054 10	- 13.0 = 9.3301E+02 uCi/L beta	139 10	- 13.0 = 5.1552E-04 uCi/L beta					
	F 108.-6225							F 120.-6325

9 1 1 2 7 6 0 1 2 7 6

## Total Beta Analysis on the Fusion Dissolution

3-4				
Serial No. F 109.-6525	Sample Point SEGMENT-J	Date 11-15-89	Time Issued 10:56	Priority 26
Determination TB	Method/Standard LA-508-101	Result Units % RECOVERY	Charge Code F215C	Reruns 0
Sample Size ? 10 ml				Customer ID
Remarks Calculations Results: LMCS CHECK SAMPLE LMCS ID 102B44				
Analyst -1 6C269 Mason Greene	Analyst -2	Analyst -3	Analyst -4	Analyst -5
Date 1-25-90	Hrs	Hrs	Hrs	Hrs
	Time Completed	Lab Unit Mgr <i>EPA</i>	Pm 5	

100.4%

$1.391^{-1}$  /  $1.3851^{-1}$

54-0800-081 (R-10-85)

9 1 1 2 3 6 0 1 2 7 7

Total Beta Analysis on the Fusion Dissolution

18 2" 1-25-80 LR

$$\frac{9824}{10} - 13 \\ \frac{9902}{10}$$

Beta Calculation by VR on 01-25-1980 at 10:42:57  
Det #18 2-inch mount Beta eff. : .3151  
Sample size : 10 ml Dilution : 1

Mount # 1

$$\frac{9824}{10} - 13.0 = 1.3858E-01 \mu\text{Ci/L beta}$$

Mount # 2

$$\frac{9902}{10} - 13.0 = 1.3979E-01 \mu\text{Ci/L beta}$$

F 109.-6525

Gamma Energy Analysis of the Fusion Dissolution

9 1 1 2 2 6 0 1 2 7 3

Serial No. F 81.-6530	Sample Point SEGMENT-6		Date 11-15-89	Time Issued 10:22	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size ? 500 <sup>1</sup>					Customer ID Sld 89844
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 89844					
$\text{Cs}^{137}$ 2.34' / 2.2255' 105.1% $\text{Cs}^{137}$ 3.81' / 3.813' 99.9%					
Analyst-1 PMS/69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5 ghe	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 1-9-90	Time Completed	Lab Unit Mgr Cja	KI		
SI-4800-081 (R-10-83)					

**2742**

Serial No. F 106.-6030	Sample Point SEGMENT-7		Date 11-15-89	Time Issued 10:55	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units uCi/L	Charge Code WB75L	Reruns 0	
Sample Size ? 500 <sup>1</sup>					Customer ID 089045
Remarks, Calculations, Results:					
$\text{Cs}^{137}$ 4.48' ueifl or 1.87' ueifg					
Analyst-1 PMS/69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5 ghe	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 1-9-90	Time Completed	Lab Unit Mgr Cja	KI		
SI-4800-081 (R-10-83)					

Serial No. F 107.-6130	Sample Point SEGMENT-8		Date 11-15-89	Time Issued 10:56	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units uCi/L	Charge Code WB75L	Reruns 0	
Sample Size ? 500 <sup>1</sup>					Customer ID 089045
Remarks, Calculations, Results: DUPLICATE SAMPLE					
$\text{Cs}^{137}$ 5.36' ueifl or 2.34' ueifg					
Analyst-1 69769 PMS	Analyst-2	Analyst-3	Analyst-4	Analyst-5 ghe	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 1-9-90	Time Completed	Lab Unit Mgr Cja	KI		
SI-4800-081 (R-10-83)					

**1002**

Serial No. F 192.-6330	Sample Point SEGMENT-U		Date 11-17-89	Time Issued 10:34	Priority 18
Determination GEA	Method/Standard LA-548-121	Result Units uCi/L	Charge Code WB75L	Reruns 0	
Sample Size ? uml					Customer ID 089049
Remarks, Calculations, Results: REAGENT BLANK					
$\text{Cs}^{137}$ <4.77 <sup>-2</sup> ueifl					
Analyst-1 PMS/69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5 ghe	
Hrs	Hrs	Hrs	Hrs	Hrs	
Date 1-9-90	Time Completed	Lab Unit Mgr Cja	KI		
SI-4800-081 (R-10-83)					

Gamma Energy Analysis of the Fusion Dissolution

1000

Serial No. F 84-6230		Sample Point SEGMENT-9		Date 11-15-89	Time Issued 10:23	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units % RECOVERY	Charge Code WB75L	Re runs O		
Sample Size ? <del>500</del> 1001			Customer ID CS9044			
Remarks, Calculations, Results: SPIKE SAMPLE #139 F82 SPIKE ID 89844 SPIKE VOLUME 1001						
$1.06^2 - 6.64 = 3.96 / 3.813$ <span style="float: right;">104.07051/500 103.9</span>						
Analyst - 1 DMS/69765	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 Hrs		
Date 1-9-90	Time Completed	Lab Unit Mgr <i>CJA</i>	<i>TJW</i>	<i>DK</i>		

84-0000-061 (R-10-63)

#881

Serial No. F 181-6330		Sample Point SEGMENT-J		Date 11-17-89	Time Issued 10:32	Priority 19
Determination GEA	Method/Standard LA-548-121	Result Units % RECOVERY	Charge Code WB75L	Re runs O		
Sample Size ? 5001			Customer ID SLI 89844			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID _____						
$\text{Co}^{60} \quad 2.26 \text{ weight} / 2.2255 \quad 101.6\%$ $\text{Cs}^{137} \quad 3.66 \text{ weight} / 3.813 \quad 96.0\%$						
Analyst - 1 DMS/69765	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 Hrs		
Date 1-9-90	Time Completed	Lab Unit Mgr <i>CJA</i>	<i>DK</i>			

84-0000-061 (R-10-63)

9 1 1 2 2 6 7 1 2 0

Uranium Analysis of the Fusion Dissolution

Serial No.	Sample Point		Date	Time issued	Priority
F 120.-6340	SEGMENT-21		11-15-89	10:58	1B
Determination	Method/Standard	Result Units	Charge Code	Runno	
U	LA-925-106	G/L	WB75L	0	
Sample Size				Customer ID	
? 100.0 - 10ml - 100ml				089045	
Remarks, Calculations, Results: REAGENT BLANK Spike 54.838 Value 5.62E-4 Volume .100ml					
$\begin{aligned} & (.982)(.06)(5.62E-4)(1)(1010) \quad .06 \\ & .38 - [(.982)(.06)] \quad 2.104E-3 \end{aligned}$					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
GC269					
MMoney	Hrs	Hrs	f. thg	Hrs	
Date	Time Completed	Lab Unit Mgr	Signature		
1-5-90		MMoney	MMoney		
SI-6000-061 (R-10-83)					

Serial No.	Sample Point		Date	Time issued	Priority
F 105.-6540	SEGMENT-6		11-15-89	10:55	23
Determination	Method/Standard	Result Units	Charge Code	Runno	
U	LA-925-106	% RECOVERY	WB75L	0	
Sample Size				Customer ID	
? 100-10-100				089045	
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 58838					
$\begin{aligned} & 97.0\% \quad \text{Sample: 16} \\ & 98.2\% \quad \text{Spike + Spk: 44} \\ & \text{Spk Vol: 100} \\ & \text{Spk ID: 5.62E-4} \end{aligned}$					
$\begin{aligned} & (.16)(5.62E-4)(1)(1010) = 2.916E-3 \\ & (\text{over for recalculation}) \quad .30 \quad 2.916E-3 \end{aligned}$					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
GC269					
MMoney	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr	Signature		
1-5-90		MMoney	MMoney		
SI-6000-061 (R-10-83)					

Serial No.	Sample Point		Date	Time issued	Priority
F 107.-6140	SEGMENT-8		11-15-89	10:56	23
Determination	Method/Standard	Result Units	Charge Code	Runno	
U	LA-925-106	G/L	WB75L	0	
Sample Size				Customer ID	
? 100-10-100				089045	
Remarks, Calculations, Results: DUPLICATE SAMPLE					
$\begin{aligned} & \text{Spike Vol: 100} \quad \text{Sample: .18} \\ & \text{Spike + Spk: .48} \\ & \text{Spk ID: 5.113E-5} \end{aligned}$					
$\begin{aligned} & (.18)(5.113E-5)(1)(1010) = 3.01E-3 \end{aligned}$					
$\begin{aligned} & (\text{over for recalculation}) \quad 3.01E-3 \end{aligned}$					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
GC269					
MMoney	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr	Signature		
1-5-90		MMoney	MMoney		
SI-6000-061 (R-10-83)					

Serial No.	Sample Point		Date	Time issued	Priority
F 106.-6040	SEGMENT-7		11-15-89	10:55	23
Determination	Method/Standard	Result Units	Charge Code	Runno	
U	LA-925-106	G/L	WB75L	0	
Sample Size				Customer ID	
? 100-10-100				089045	
Remarks, Calculations, Results:					
$\begin{aligned} & \text{Spike Vol: 100} \quad \text{Sample: .20} \\ & \text{Spike + Spk: .50} \\ & \text{Spk ID: 5.113E-5} \end{aligned}$					
$\begin{aligned} & (.20)(5.113E-5)(1)(1010) = 3.4E-3 \end{aligned}$					
$\begin{aligned} & (\text{over for recalculation}) \quad 3.4E-3 \end{aligned}$					
Analyst-1	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
GC269					
MMoney	Hrs	Hrs	Hrs	Hrs	
Date	Time Completed	Lab Unit Mgr	Signature		
1-5-90		MMoney	MMoney		
SI-6000-061 (R-10-83)					

Uranium Analysis of the Fusion Dissolution

F 106.-6040

$\Delta^{18}\text{E} = 3.34 \times 10^{-3}$

$$\frac{\Delta^{18}\text{E}}{\Delta^{18}\text{O}_{\text{std}}} = \frac{[0.02 - \left(\frac{16}{18}\right)0.5] \cdot 66000}{(0.00099)(0.00516)(0.00001)}$$

F 105.-6540

$$\frac{\Delta^{18}\text{E}}{\Delta^{18}\text{O}_{\text{std}}} = \frac{0.02 - 0.016}{0.00099} = 58.79\%$$

$$\frac{\Delta^{18}\text{E}}{\Delta^{18}\text{O}_{\text{std}}} = \frac{[0.02 - \left(\frac{16}{18}\right)0.5] \cdot 66000}{(0.00099)(0.00516)(0.00001)} = 5.95 \times 10^{-2}$$

F 107.-6140

$\Delta^{18}\text{E} = 3.01 \times 10^{-3}$

$$\frac{\Delta^{18}\text{E}}{\Delta^{18}\text{O}_{\text{std}}} = \frac{[0.02 - \left(\frac{16}{18}\right)0.5] \cdot 66000}{(0.00099)(0.00516)(0.00001)} = 5.95 \times 10^{-2}$$

9 1 1 2 0 3 7 1 2 ? 2

## Uranium Analysis of the Fusion Dissolution

Serial No. F 297.-6540	Sample Point SEGMENT-J	Date 11-21-89	Time Issued 8:30	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Re runs 0
Sample Size ? 100-10-100	Customer ID 089050			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 58838  Spike Vol: 100x Spike ID: 5.62 <sup>-4</sup> (.1)(5.62 <sup>-4</sup> )(.1)(100) = 3.33 <sup>-3</sup> (over for recalculation) 46-17  Analyst -1 GC269 C. M. Gandy Date 1-5-90				
Analyst -2	Analyst -3	Analyst -4	Analyst -5	Analyst -6
Hrs	Hrs	Hrs	Hrs	
Date 1-5-90	Time Completed CJP	Lab Unit/Mgr GRW		

54-9800-081 (R-10-80)

Serial No. F 108.-6240	Sample Point SEGMENT-9	Date 11-15-89	Time Issued 10:56	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Re runs 0
Sample Size ? 100-10-100	Customer ID 089015			
Remarks, Calculations, Results: Spike Sample F108 Spike ID 58838 Spike Volume 100-10-100 ((20)(518-4))(100) Sample Vol: .20 Spike Vol: 100A (20)(518-4)(100) Sample Vol: .46 Spike ID: 5.62 <sup>-4</sup> .46-.20 = 3.33 <sup>-3</sup> (over for recalculation) 5.3 <sup>-4</sup> Spike To 10 Calculated  Analyst -1 GC269 C. M. Gandy Date 1-5-90				
Analyst -2	Analyst -3	Analyst -4	Analyst -5	Analyst -6
Hrs	Hrs	Hrs	Hrs	Hrs
Date 1-5-90	Time Completed CJP	Lab Unit/Mgr GRW	100	

54-9800-081 (R-10-80)

Uranium Analysis of the Fusion Dissolution

F 108.-6240

F 297.6540

$$\frac{2.99 \times 10^{-2}}{5.11 \times 10^{-4}} = 1.71 \%$$

$$\frac{(0.00059)(0.46(\frac{5.7}{5.8}) - 0.20)}{(0.00059)(0.46(\frac{5.7}{5.8}) - 0.20) + (0.000516)(0.20)} = \frac{0.00059}{0.00059 + 0.000516} = 0.001032$$

$$\frac{3.99 \times 10^{-2}}{3.24 \times 10^{-2}} = \frac{0.00059[(0.46(\frac{5.7}{5.8}) - 0.17)]}{(0.00059)(0.46(\frac{5.7}{5.8}) - 0.17) + (0.000516)(0.20)} = \frac{0.00059}{0.00059 + 0.000516} = 0.001032$$

9 | 1 | 2 | 0 | 3 | 7 | 1 | 2 | 7 | 4

Uranium Analysis of the Fusion Dissolution

Serial No. F 106.-6040		Sample Point SEGMENT-7		Date 11-15-89	Time Issued 10:55	Priority 23
Determination U	Method/Standard LA-925-106	Result Units G/L	Charge Code WB75L	Rerun		
Sample Size $100 - 10 - 500$ $100 - 10\text{ml} - .100$				Customer ID 089045		
Remarks, Calculations, Results:						
$\text{spike ID/vol: } 54838 \quad 2.59E-3 \text{ g/l RERUN}$ $5.62^{-4} \quad 0.050$						
$\text{spike vol: } 100\text{ml} \quad \frac{(6.0)}{6.1} \times 0.05 \times (5.62E-4) \times 10^2 \times (5.02) =$ $- 265 - \left( \frac{5.6}{5.7} \times 0.05 \right)$						
Analyst-1 Sue Lai	Analyst-2 Hrs 60916	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	<i>S. M. Scott</i>	
Date 7-12-90	Time Completed	Lab Unit Mgr.	<i>Tall Bob Dymitri Sankar</i> 54-0000-081 (R-10-83)			

Serial No. F 105.-6540		Sample Point SEGMENT-6		Date 11-15-89	Time Issued 10:55	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Rerun		
Sample Size .100 - 10ml - .100				Customer ID 089045		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 58838						
$\text{spike ID/vol: } 54838 \quad \frac{5.6}{5.7} \times 0.10 \times (5.62E-4) \times 10^2 / 100 = 2.7E-2$ $5.62^{-4} \quad 0.335 - \left( \frac{5.6}{5.7} \times 0.10 \right) = 0.335$						
$\text{spike vol: } 100\text{ml} \quad 90.3\% \quad 299\%$						
Analyst-1 Sue Lai	Analyst-2 Hrs 60916	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	<i>S. M. Scott</i>	
Date 7-12-90	Time Completed	Lab Unit Mgr.	<i>Tall Bob Dymitri Sankar</i> 54-0000-081 (R-10-83)			

Serial No. F 108.-6240		Sample Point SEGMENT-9		Date 11-15-89	Time Issued 10:56	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Rerun		
Sample Size .100 - 10ml - .100				Customer ID 089045		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 54838 F106 Sample SPIKE VOLUME 100ml RERUN $5.6 \times 10^{-5} \times (5.62E-4) \times 10^2$ $\frac{5.6}{5.7} \times 10^{-5} \times 10^2 = 0.350 = 0.125 = 95.0\%$						
Analyst-1 Sue Lai	Analyst-2 Hrs 60916	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	<i>S. M. Scott</i>	
Date 7-12-90	Time Completed	Lab Unit Mgr.	<i>Tall Bob Dymitri Sankar</i> 54-0000-081 (R-10-83)			

Serial No. F 297.-6540		Sample Point SEGMENT-J		Date 11-21-89	Time Issued 8:30	Priority 23
Determination U	Method/Standard LA-925-106	Result Units % RECOVERY	Charge Code WB75L	Rerun		
Sample Size .100ml - 10ml - .100ml				Customer ID 089050		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 58838 $3.10E-2 - 2.7E-3 =$ $\text{spike ID/vol: } 54838 \quad 2.84E-2 / 2.99E-2 \times 10 = 92.9\%$ $5.62^{-4} \quad \frac{5.6}{5.7} \times 0.335 = 0.1 \times 10^2 = 100\text{ml}$						
Analyst-1 Sue Lai	Analyst-2 Hrs 60916	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	<i>S. M. Scott</i>	
Date 7-12-90	Time Completed	Lab Unit Mgr.	<i>Tall Bob Dymitri Sankar</i> 54-0000-081 (R-10-83)			

9 1 - 2 2 3 2 - 2 2 5

$$\frac{\left(\frac{5.6}{5.7}\right) \times .110 \times (5.62e-4) \times (.1) \times (1010)}{.530 - \left[\left(\frac{5.6}{5.7}\right) \times .110\right]} = 2.76e-2 / 2.99e-2 = 90.9\%$$

93.9%

F297.-6540

Serial No <b>F 120.-6340</b>	Sample Point <b>SEGMENT-21</b>	Date <b>11-15-89</b>	Time Issued <b>10:38</b>	Priority <b>18</b>
Determination <b>U</b>	Method/Standard <b>LA-720-106</b>	Result Units <b>G/L</b>	Caliber Code <b>WB75E</b>	Remarks
Sample Size 1,400-16mls	50mL (S)			Customer ID <b>089045</b>
Preparation Calculations Results: <b>REAGENT BLANK</b>		<b>F 168 Blank</b> <b>RERUN</b> 1.346 0.000 <del>Pass</del>		
SPIKE: 54838 5.624.		$\frac{100}{( \frac{5.55}{5.62} ) - 1}$ $\frac{0.185}{( \frac{5.55}{5.62} ) - 1}$ $\frac{20}{( \frac{5.55}{5.62} ) - 1}$ $\frac{812}{( \frac{5.55}{5.62} ) - 1}$		
Analyst-1 <b>Sue Yon</b>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <b>Meter</b>
Hrs <b>6C716</b>	Hrs	Hrs	Hrs	Hrs
Date <b>7-12-90</b>	Time Completed	Lab Unit Mgr <b>Zell P. D. Jones</b>	Signature	

2

9 1 1 2 3 5 2 1 2 3 6

Water Digestion

Serial No. F 111.-7000	Sample Point SEGMENT-12	Date 11-15-89	Time issued 10:56	Priority 19
Determination H2O-DGST	Method/Standard LA-504-101	Result Units G/L	Charge Code WB75L	Runno 0
Sample Size ?		Customer ID 089045		
Remarks, Calculations, Results: GRAMS SAMPLE <u>5307</u> #233 VOLUME ON COMPLETION <u>50 ml</u> <u>1.06<sup>-2</sup> g/ml</u> <u>1.06<sup>-5</sup> g/l</u> <del>1.06<sup>-2</sup> g/ml</del> <del>1.06<sup>-5</sup> g/l</del> <u>10.6</u> <u>WTC N 313 4</u>		2.92 grams		
Analyst - 1 #1098	Analyst - 2 6B07	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 (Chemist) CJF 9471
20 min	0 min	<i>Tad M. Bal</i>	5/20/90	
Date 1/8/90	Time Completed 0900	Lab Unit Mgr <i>CJF</i>		

Serial No. F 112.-7100	Sample Point <b>SEGMENT-13</b>	Date 11-15-89	Time Issued 10:56	Priority 19
Determination H2O-DGST	Method/Standard LA-504-101	Result Units G/L	Charge Code WB75L	Runno 0
Sample Size ?		Customer ID <b>089045</b>		
Remarks, Calculations, Results: <b>DUPLICATE ANALYSIS</b> GRAMS SAMPLE <u>.5161</u> VOLUME ON COMPLETION <u>50ml</u>  <u>1.032 - 2 g/ml</u> <u>1.032 - 5 g/l</u> <sup>WB75L</sup> <u>5/15/90</u> <u>10.32</u> <u>1.032 - 5 g/l</u> <sup>WB75L</sup> <u>5/15/90</u> <u>1.032 g/l</u>				
Analyst - 1 81098	Analyst - 2 6B07	Analyst - 3	Analyst - 4 5/15/90	Analyst - 5 Completed
Hrs 20 min	Hrs 10 min.	Hrs Tall. Prol	Hrs 5/15/90	Hrs Completed
Date 1/8/90	Time Completed 0900	Lab Unit Mgr CJ		

Water Digestion

9 1 1 2 0 6 0 1 2 7 7

WT 1: 130.7603

SEQUENCE #: 100

WT 1: 129.5140  
WT 2: 129.3450

NET WEIGHT:

--> 0.5687 GRAMS

01/02/90 8 11:53:34

WT 2: 131.2660

NET WEIGHT:  
--> 0.5824 GRAMS

01/02/90 8 11:53:32

SPIKE...

SAMPLE

SEQUENCE #: 99

WT 1: 129.1260  
WT 2: 129.6420

NET WEIGHT:

--> 0.5161 GRAMS

01/02/90 8 11:43:03

DUPPLICATE

9 1 1 2 7 6 7 1 2 7 8

## Ion Chromatographic Analysis of the Water Digestion - Floride

Serial No. F 122.-7371	Sample Point SEGMENT-23		Date 11-15-89	Time Issued 10:58	Priority 18
Determination F	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? 700-10 Direct			Customer ID 089045		
Remarks, Calculations, Results: REAGENT BLANK  ~.1 ppm					
Analyst - 1 6B107/new	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5/checked	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr CJW	Customer ID 5215		

54-000-061 (R-10-83)

Serial No. F 110.-7571	Sample Point SEGMENT-11		Date 11-15-89	Time Issued 10:56	Priority 19
Determination F	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size 100-10			Customer ID 6CHIE		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 6C11P1					
Analyst - 1 6B107/new	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5/checked	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr CJW	Customer ID 5215		

54-000-061 (R-10-83)

$\frac{67}{72} \times 100 = 93.1\%$

67.04 5215

Serial No. F 112.-7171	Sample Point SEGMENT-13		Date 11-15-89	Time Issued 10:57	Priority 19
Determination F	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results: DUPLICATE SAMPLE  1.92 <sup>2</sup> ppm					
Analyst - 1 6B107/new	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5/checked	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr CJW	Customer ID 5215		

54-000-061 (R-10-83)

Serial No. F 111.-7071	Sample Point SEGMENT-12		Date 11-15-89	Time Issued 10:56	Priority 19
Determination F	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results:  1.82 <sup>2</sup>					
Analyst - 1 6B107/new	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5/checked	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr CJW	Customer ID 5215		

54-000-061 (R-10-83)

9 1 1 2 2 6 0 1 2 2 9

## Ion Chromatographic Analysis of the Water Digestion - Florida

Serial No. F 66.-7571	Sample Point SEGMENT-15		Date 11-15-89	Time Issued 10:21	Priority 19
Determination F	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Runnoe 0	
Sample Size 100-10			Customer ID		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <u>6C11A</u>					
$100 \left( \frac{67.06}{72} \right) = 93.1\%$ <p style="text-align: center;">gms</p>					
Analyst-1 <i>6B107/ew</i>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>SPK/2111</i>	Analyst-6
Hrs .5	Hrs	Hrs	Hrs	Hrs	Hrs
Date 2/15/90	Time Completed	Lab Unit Mgr <i>CJW</i>			
54-000-081 (R-10-83)					

Serial No. F 113.-7271	Sample Point SEGMENT-14		Date 11-15-89	Time Issued 10:57	Priority 19
Determination F	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Runnoe 0	
Sample Size ? 100 uL - 10 ml			Customer ID 89-045		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 35C9-67 SPIKE VOLUME .300/5mL					
$\left( \frac{5.3 \text{ mL}}{5.0 \text{ mL}} \right) \left( \frac{406.9 \text{ ppm}}{47 \text{ ppm}} \right) - \left( \frac{181.9 \text{ ppm}}{10.6 \text{ ppm}} \right) \left( \frac{10.05 \text{ mL}}{10.6 \text{ mL}} \right) \times 100 = 96.3\%$ <p style="text-align: center;">(.30 mL)(47 ppm) / (101) / 5.3 mL</p>					
Analyst-1 <i>6B107</i>	Analyst-2	Analyst-3	Analyst-4	Analyst-5 <i>SPK/2111</i>	Analyst-6
Hrs .5	Hrs	Hrs	Hrs	Hrs	Hrs
Date 2/15/90	Time Completed	Lab Unit Mgr <i>SPK/2111</i>			
54-000-081 (R-10-83)					

9 1 1 2 3 5 7 1 2 3 0

Ion Chromatographic Analysis of the Water Digestion - Chloride

Serial No. F 122.-7372 Sample Point SEGMENT-23 Date 11-15-89 Time Issued 10:58 Priority 18 Determination CL Method/Standard LA-533-105 Result Units PPM Charge Code WB75L Returns 0 Sample Size ? Direct Customer ID US89045					
Remarks, Calculations, Results: REAGENT BLANK  <i>L.1 ppm</i> <i>.15 ppm</i>					
Analyst-1 6B107/NEW	Analyst-2	Analyst-3	Analyst-4	Analyst-5 Chemist CLP/254171	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr. CGN	SEY		
54-0800-061 (R-10-82)					

Serial No. F 110.-7572 Sample Point SEGMENT-11 Date 11-15-89 Time Issued 10:56 Priority 19 Determination CL Method/Standard LA-533-105 Result Units % RECOVERY Charge Code WB75L Returns 0 Sample Size 100-10 Customer ID US89045					
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 6C11RT					
Analyst-1 6B107/NEW	Analyst-2	Analyst-3	Analyst-4	Analyst-5 Chemist CLP/254171	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr. CLP/254171	SEY		
54-0800-061 (R-10-82)					

Serial No. F 112.-7172 Sample Point SEGMENT-13 Date 11-15-89 Time Issued 10:57 Priority 19 Determination CL Method/Standard LA-533-105 Result Units PPM Charge Code WB75L Returns 0 Sample Size ? 100-10 Customer ID US89045					
Remarks, Calculations, Results: DUPLICATE SAMPLE  <i>&lt;10.1 ppm</i>					
Analyst-1 6B107/NEW	Analyst-2	Analyst-3	Analyst-4	Analyst-5 Chemist CLP/254171	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr. CGN	SEY		
54-0800-061 (R-10-82)					

Serial No. F 111.-7072 Sample Point SEGMENT-12 Date 11-15-89 Time Issued 10:56 Priority 19 Determination CL Method/Standard LA-533-105 Result Units PPM Charge Code WB75L Returns 0 Sample Size ? 100-10 Customer ID US89045					
Remarks, Calculations, Results:  <i>&lt;10.1 ppm</i>					
Analyst-1 6B107/NEW	Analyst-2	Analyst-3	Analyst-4	Analyst-5 Chemist CLP/254171	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr. CGN	SEY		
54-0800-061 (R-10-82)					

9 1 1 2 2 3 0 1 2 2 1

## Ion Chromatographic Analysis of the Water Digestion - chloride

Serial No. F 66-7572	Sample Point SEGMENT-15		Date 11-15-89	Time issued 10:21	Priority 19
Determination CL	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size 100-10			Customer ID		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <u>GC1112</u>					
$\left(\frac{53.8}{87}\right) 100 = 61.3\%$ <p style="text-align: center;">AM 5</p>					
Analyst - 1 <u>6B1C7/WB7</u>	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5 <u>CLP/24171</u>	Comments
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr <u>CLP</u>	<u>SH</u>		

54-0000-001 (R-10-83)

Serial No. F 113-7272	Sample Point SEGMENT-14		Date 11-15-89	Time issued 10:57	Priority 19
Determination CL	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size ? 100 μL - 10 mL			Customer ID 89-045		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID <u>3509-67</u> SPIKE VOLUME <u>.300/5mL</u>					
$\frac{(1.06)(347) - (0)(.948)}{(.300)(60)} \times 100 = 107.2\%$ <p style="text-align: center;">5.3</p>					
Analyst - 1 <u>6B1C7</u>	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5 <u>CLP</u>	Comments
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/15/90	Time Completed	Lab Unit Mgr <u>CLP</u>	<u>Kathy Hambley</u>		

54-0000-001 (R-10-83)

9 1 1 2 0 5 0 1 2 0 2

## Ion Chromatographic Analysis of the Water Digestion - Nitrate

Serial No. F 122.-7373	Sample Point SEGMENT-23		Date 11-15-89	Time Issued 10:58	Priority 18
Determination NO3	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Runno 0	
Sample Size ? Direct			Customer ID 089045		
Remarks, Calculations, Results: REAGENT BLANK					
$\leq 1 \text{ ppm}$					
Analyst-1 6B107/wk4	Analyst-2	Analyst-3	Analyst-4 <i>6B107/wk4</i>	Analyst-5 <i>6B107/wk4</i>	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/18/80	Time Completed	Lab Unit Mgr <i>CJW</i>	SLX		
SI-0000-001 (R-10-03)					

Serial No. F 110.-7573	Sample Point SEGMENT-11		Date 11-15-89	Time Issued 10:56	Priority 19
Determination NO3	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Runno 0	
Sample Size 100-10			Customer ID 089045		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 6C17A1					
$(\frac{710}{722})_{\text{YOC}} = 98.3\%$					
Analyst-1 6B107/wk4	Analyst-2	Analyst-3	Analyst-4 <i>6B107/wk4</i>	Analyst-5 <i>6B107/wk4</i>	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/18/80	Time Completed	Lab Unit Mgr <i>CJW</i>	SLX		
SI-0000-001 (R-10-03)					

Serial No. F 112.-7173	Sample Point SEGMENT-13		Date 11-15-89	Time Issued 10:57	Priority 19
Determination NO3	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Runno 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results: DUPLICATE SAMPLE					
$2.74^{\text{a}}$ ppm					
Analyst-1 6B107/wk4	Analyst-2	Analyst-3	Analyst-4 <i>6B107/wk4</i>	Analyst-5 <i>6B107/wk4</i>	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/18/80	Time Completed	Lab Unit Mgr <i>CJW</i>	SLX		
SI-0000-001 (R-10-03)					

Serial No. F 111.-7073	Sample Point SEGMENT-12		Date 11-15-89	Time Issued 10:56	Priority 19
Determination NO3	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Runno 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results:					
$2.70^{\text{a}}$ ppm					
Analyst-1 6B107/wk4	Analyst-2	Analyst-3	Analyst-4 <i>6B107/wk4</i>	Analyst-5 <i>6B107/wk4</i>	
Hrs .5	Hrs	Hrs	Hrs	Hrs	
Date 2/18/80	Time Completed	Lab Unit Mgr <i>CJW</i>	SLX		
SI-0000-001 (R-10-03)					

9 1 1 2 7 6 7 1 2 7 3

## Ion Chromatographic Analysis of the Water Digestion - Nitrate

Serial No. F 66-7573	Sample Point SEGMENT-15	Date 11-15-89	Time Issued 10:21	Priority 19
Determination NO <sub>3</sub>	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Renew 0
Sample Size 100-10		Customer ID		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID EC11H1				
$\left(\frac{759.3}{722}\right)_{10} = 105.3\%$ <p style="text-align: center;">gms</p>				
Analyst - 1 603107/NWJ	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2/15/90	Time Completed	Lab Unit Mgr CJW	SLV	Analyst - 5 Signature
54-9000-001 (R-10-32)				

Serial No. F 113-7273	Sample Point SEGMENT-14	Date 11-15-89	Time Issued 10:57	Priority 19
Determination NO <sub>3</sub>	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Renew 0
Sample Size ? 100 $\mu$ L - 10 ml		Customer ID 89-C45		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 7509-67 SPIKE VOLUME, 300/5mL				
$\frac{(1.06)(2951) - (270.3)(.948)}{(1.300)(481)} \times 100 = 104.4\%$ $\frac{5.3}{(101)}$				
Analyst - 1 603107	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2-15-90	Time Completed	Lab Unit Mgr Kathy Remender	Analyst - 5 Signature	
54-9000-001 (R-10-32)				

9 1 1 2 2 3 0 1 2 2 4

## Ion Chromatographic Analysis of the Water Digestion - Phosphate

Serial No F 122.-7374	Sample Point SEGMENT-23	Date 11-15-89	Time Issued 10:58	Priority 18
Determination PO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 1
Sample Size ? Direct	Customer ID			
Remarks, Calculations, Results: REAGENT BLANK				
RERUN				
<i>1 ppm</i>				
Analyst - 1 <i>108107/100</i>	Analyst - 2	Analyst - 3	Analyst - 4	<i>Analyst - 5 Chemist</i>
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2/22/90	Time Completed	Lab Unit/Mgr <i>Cga</i>		
54-6800-061 (R-10-53)				

Serial No. F 110.-7574	Sample Point SEGMENT-11	Date 11-15-89	Time Issued 10:56	Priority 19
Determination PO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 1
Sample Size 100-10	Customer ID			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <i>1001142</i>				
RERUN				
<i>(680.9)/100 = 94.3 %</i>				
<i>737.1 / 722.0 = 102.1 % 5/14/90</i>				
Analyst - 1 <i>108107/100</i>	Analyst - 2	Analyst - 3	Analyst - 4	<i>Analyst - 5 Chemist</i>
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2/22/90	Time Completed	Lab Unit/Mgr <i>Cga</i>	4MS	
54-6800-061 (R-10-53)				

Serial No. F 112.-7174	Sample Point SEGMENT-13	Date 11-15-89	Time Issued 10:57	Priority 19
Determination PO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 1
Sample Size ? 100-10	Customer ID			
Remarks, Calculations, Results: DUPLICATE SAMPLE				
RERUN				
<i>1.64<sup>3</sup> ppm</i>				
Analyst - 1 <i>108107/100</i>	Analyst - 2	Analyst - 3	Analyst - 4	<i>Analyst - 5 Chemist</i>
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2/22/90	Time Completed	Lab Unit/Mgr <i>Cga</i>	4MS	
54-6800-061 (R-10-53)				

Serial No. F 111.-7074	Sample Point SEGMENT-12	Date 11-15-89	Time Issued 10:56	Priority 19
Determination PO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 1
Sample Size ? 100-10	Customer ID			
Remarks, Calculations, Results:				
RERUN				
<i>1.55<sup>3</sup> ppm</i>				
Analyst - 1 <i>108107/100</i>	Analyst - 2	Analyst - 3	Analyst - 4	<i>Analyst - 5 Chemist</i>
Hrs .5	Hrs	Hrs	Hrs	Hrs
Date 2/22/90	Time Completed	Lab Unit/Mgr <i>Cga</i>	4MS	
54-6800-061 (R-10-53)				

9 1 1 2 0 5 9 1 2 2 5

## Ion Chromatographic Analysis of the Water Digestion - Phosphate

Sample No. F 114.-7574	Sample Point SEGMENT-0	Date 11-15-89	Time Issued 10:57	Priority 26
Determination PO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code E21D1	Returns 0
Sample Size ?	Customer ID			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID GC11A1				
101.0 %				
Analyst -1 46107 Hrs Date 2-22-90	Analyst -2 Hrs	Analyst -3 Hrs	Analyst -4 Hrs	Analyst -5 Chemist Signature Cover)
(S-3000-061 (R-10-42))				

Sample No. F 113.-7274	Sample Point SEGMENT-14	Date 11-15-89	Time Issued 10:57	Priority 19
Determination PO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Returns 1
Sample Size ? 100-10	Customer ID			
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID 33C9-69 SPIKE VOLUME .300/5ml				
<del>4.738<sup>3</sup>-1.55<sup>3</sup>= 3188 / 2888 = 103.2%</del>				
Analyst -1 63107/260 Hrs .5	Analyst -2 Hrs	Analyst -3 Hrs	Analyst -4 Hrs	Analyst -5 Chemist Signature Cover)
Date 2/22/90	Time Completed Lab Unit Mgr Cogn	Lab Unit Mgr Cogn	Time Completed Pm 2	(S-3000-061 (R-10-42))

9 1 1 2 0 5 7 1 2 0 6

Ion Chromatographic Analysis of the Water Digestion - Phosphate

F 113.-7274

$$\frac{\left(\frac{5.3}{5}\right)(4.738) - (1553)\left(\frac{10.03}{10.6}\right)}{(6.300)(6.01)(1.01)} \times 100 = 103.2\%$$

9 1 1 2 0 6 0 1 2 7 7

## Ion Chromatographic Analysis of the Water Digestion - Sulfate

Serial No. F 122.-7375	Sample Point SEGMENT-23	Date 11-15-89	Time Issued 10:58	Priority 18
Determination SO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 0
Sample Size ? Direct		Customer ID <b>089045</b>		
Remarks, Calculations, Results: REAGENT BLANK				
<i>&lt;1 ppm</i>				
Analyst-1 <i>6B107/100</i> Hrs .5	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>6B107/100</i> Hrs
Date <i>2/15/90</i>	Time Completed	Lab Unit Mgr <i>CJW</i>		<i>ESR</i>
54-6800-081 (R-10-82)				

Serial No. F 110.-7575	Sample Point SEGMENT-11	Date 11-15-89	Time Issued 10:56	Priority 19
Determination SO4	Method/Standard LA-533-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 0
Sample Size 100-10		Customer ID <b>089045</b>		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <i>PC1111</i>				
$\frac{100}{(696/749)} = 92.9\% \text{ BMST}$				
Analyst-1 <i>6B107/100</i> Hrs .5	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>6B107/100</i> Hrs
Date <i>2/15/90</i>	Time Completed	Lab Unit Mgr <i>CJW</i>		<i>ESR</i>
54-6800-081 (R-10-82)				

Serial No. F 112.-7175	Sample Point SEGMENT-13	Date 11-15-89	Time Issued 10:57	Priority 19
Determination SO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 0
Sample Size ? 100-10		Customer ID <b>089045</b>		
Remarks, Calculations, Results: DUPLICATE SAMPLE				
<i>&lt;10.1 ppm</i>				
Analyst-1 <i>6B107/100</i> Hrs .5	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>6B107/100</i> Hrs
Date <i>2/15/90</i>	Time Completed	Lab Unit Mgr <i>CJW</i>		<i>ESR</i>
54-6800-081 (R-10-82)				

Serial No. F 111.-7075	Sample Point SEGMENT-12	Date 11-15-89	Time Issued 10:56	Priority 19
Determination SO4	Method/Standard LA-533-105	Result Units PPM	Charge Code WB75L	Retruns 0
Sample Size ? 100-10		Customer ID <b>089045</b>		
Remarks, Calculations, Results:				
<i>&lt;1.31 ppm</i>				
Analyst-1 <i>6B107/100</i> Hrs .5	Analyst-2 Hrs	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 <i>6B107/100</i> Hrs
Date <i>2/15/90</i>	Time Completed	Lab Unit Mgr <i>CJW</i>		<i>ESR</i>
54-6800-081 (R-10-82)				

9 1 1 2 7 6 7 1 2 7 3

Ion Chromatographic Analysis of the Water Digestion - Sulfate

Serial No.	Sample Point		Date	Time Issued	Priority
F 113.-7275	SEGMENT-14		11-15-89	10:57	19
Determination	Method/Standard	Result Units			
SO4	LA-533-105	% RECOVERY	WB75L	Refuge	0
Sample Size			Customer ID		
? 100 uL - 10 mL			089045		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID <u>35C9-67</u> SPIKE VOLUME <u>.300/5 mL</u> $\frac{(1.06)(2.672) - (4.31)(1.948)}{(1.300)(4.82)} \times 100 = 101.3\%$ <u>5.3</u>					
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5	Analyst - 6
6B107				<i>CHP</i>	
Hrs	Hrs	Hrs	Hrs	Hrs	
.5					
Date	Time Completed	Lab Unit/Mer			
2/15/90		<i>Tell. Park</i>	<i>Kathy Womble</i>		
54-5800-081 (R-10-63)					

Serial No.	Sample Point		Date	Time Issued	Priority
F 66.-7575	SEGMENT-15		11-15-89	10:21	19
Determination	Method/Standard	Result Units			
SO4	LA-533-105	% RECOVERY	WB75L	Refuge	0
Sample Size			Customer ID		
100-10			Customer ID		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <u>6C11H1</u>					
$\left(\frac{7.41}{7.49}\right) 100 = 98.9\% \text{ PMS}$					
Analyst - 1	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5	Analyst - 6
6B107/ <i>new</i>				<i>CHP</i>	<i>CHP</i>
Hrs	Hrs	Hrs	Hrs	Hrs	
.5					
Date	Time Completed	Lab Unit/Mer			
2/15/90		<i>CJN</i>	<i>KZL</i>		
54-5800-081 (R-10-63)					

9 1 1 2 0 6 7 1 2 0 9

## Total Organic Carbon Analysis of the Water Digestion

Serial No. F 113.-7226	Sample Point SEGMENT-14		Date 11-15-89	Time Issued 10:57	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 1	
Sample Size ? 200ul +100ul .5m H <sub>2</sub> SO <sub>4</sub> - 200ul			Customer ID 89045		
Remarks, Calculations, Results: SPIKE SAMPLE SPIKE ID <u>TOC11</u> SPIKE VOLUME <u>200 ul</u>					
$\frac{(21.1-4)_{\text{lab}} - (7.8-6.4)_{\text{lab}}}{119.9 \text{ ug}} \times 100 = 94.5\%$					
Analyst-1 <u>80028</u>	Analyst-2 <u>Ed Cohn</u>	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 5-30-90	Time Completed	Lab Unit Mgr <u>Tell M. P. D. C. H. L. B. H. B. H. B.</u> S4000-001 (R-10-82)			

Serial No. F 114.-7526	Sample Point SEGMENT-15		Date 11-15-89	Time Issued 10:57	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Retruns 1	
Sample Size ? 200ul - 2nd - 200ul			Customer ID 89045		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID TOC11D					
$2.882 \text{ g/l}$ 96.1% 3.00 g/l					
Analyst-1 <u>80028</u>	Analyst-2 <u>Ed Cohn</u>	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 5-30-90	Time Completed	Lab Unit Mgr <u>Tell M. P. D. C. H. L. B. H. B. H. B.</u> S4000-001 (R-10-82)			

Serial No. F 111.-7026	Sample Point SEGMENT-12		Date 11-15-89	Time Issued 10:56	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units G/L	Charge Code WB75L	Retruns 1	
Sample Size ? 1ml +100ul .5m H <sub>2</sub> SO <sub>4</sub> - 200ul			Customer ID 89045		
Remarks, Calculations, Results:					
$7.7 \times 10^{-3} \text{ g/l}$					
Analyst-1 <u>80028</u>	Analyst-2 <u>Ed Cohn</u>	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 5-30-90	Time Completed	Lab Unit Mgr <u>Tell M. P. D. C. H. L. B. H. B. H. B.</u> S4000-001 (R-10-82)			

Serial No. F 112.-7126	Sample Point SEGMENT-13		Date 11-15-89	Time Issued 10:56	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units G/L	Charge Code WB75L	Retruns 1	
Sample Size ? 1ml +100ul .5m H <sub>2</sub> SO <sub>4</sub> - 200ul			Customer ID 89045		
Remarks, Calculations, Results: DUPLICATE SAMPLE					
$7.15 \times 10^{-3} \text{ g/l}$					
Analyst-1 <u>80028</u>	Analyst-2 <u>Ed Cohn</u>	Analyst-3 Hrs	Analyst-4 Hrs	Analyst-5 Hrs	
Date 5-30-90	Time Completed	Lab Unit Mgr <u>Tell M. P. D. C. H. L. B. H. B. H. B.</u> S4000-001 (R-10-82)			

9 1 1 2 0 6 7 1 3 0 0

## Total Organic Carbon Analysis of the Water Digestion

Serial No. F 110.-7526	Sample Point SEGMENT-11	Date 11-15-89	Time Needed 10:56	Priority 19
Determination TOC	Method/Standard LA-344-105	Result Units % RECOVERY	Charge Code WB75L	Run No. 1
Sample Size ? 200ml - 2ml - 200ul	Customer ID 89045			
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID <u>70611D</u>				
<p style="text-align: center;">97.4%</p> <p style="text-align: center;"><del>2.9205</del> <del>3.00</del></p>				
Analyst - 1 <u>80028</u> <i>Ed Cohen</i>	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 <u>78104</u> <i>Ed Cohen</i>
Date 5-30-90	Time Completed <i>Tell M. Paul H. H. H. H. H.</i>	Lab Unit Mgr <i>Tell M. Paul H. H. H. H. H.</i>		

Serial No. F 122.-7326	Sample Point SEGMENT-23	Date 11-15-89	Time Issued 10:58	Priority 18
Determination TOC	Method/Standard LA-344-105	Result Units <u>TOC ug</u>	Charge Code WB75L	Run No. 1
Sample Size ? 200ul	Customer ID 89045			
Remarks, Calculations, Results: REAGENT BLANK				
<p style="text-align: center;">.9120 ug/min</p> <p style="text-align: center;">6.4 ug</p>				
Analyst - 1 <u>80028</u> <i>Ed Cohen</i>	Analyst - 2 Hrs	Analyst - 3 Hrs	Analyst - 4 Hrs	Analyst - 5 <u>78104</u> <i>Ed Cohen</i>
Date 5-30-90	Time Completed <i>Tell M. Paul H. H. H. H. H.</i>	Lab Unit Mgr <i>Tell M. Paul H. H. H. H. H.</i>		

Acid Digestion

9   1   2   0   6   9   3   0																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Serial No. F 117.-8100</td> <td>Sample Point SEGMENT-18</td> <td>Date 11-15-89</td> <td>Time Issued 10:57</td> <td>Priority 23</td> <td></td> </tr> <tr> <td>Determination ACD-DGST</td> <td>Method/Standard LA-505-159</td> <td>Result Units g/ml sj</td> <td>Charge Code WB75L</td> <td>Reruns 0</td> <td></td> </tr> <tr> <td>Sample Size ?</td> <td colspan="4">Customer ID <b>089045</b></td> <td></td> </tr> <tr> <td colspan="6">           Remarks, Calculations, Results:  <b>DUPLICATE ANALYSIS</b>  <b>GRAMS SAMPLE</b>  <b>VOLUME ON COMPLETION SJml</b>  <math>1.01^{-2}</math> g/ml         </td> </tr> <tr> <td>Analyst - 1 69769</td> <td>Analyst - 2</td> <td>Analyst - 3</td> <td>Analyst - 4</td> <td>Analyst - 5</td> <td></td> </tr> <tr> <td>K. Saubelick</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td></td> </tr> <tr> <td>Date 1/3/90</td> <td>Time Completed</td> <td>Lab Unit Mgr <i>CJW</i></td> <td><i>OK</i></td> <td><i>SJONES</i></td> <td></td> </tr> <tr> <td colspan="6" style="text-align: right;">SJ-8908-081 (R-10-83)</td> </tr> </table>						Serial No. F 117.-8100	Sample Point SEGMENT-18	Date 11-15-89	Time Issued 10:57	Priority 23		Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0		Sample Size ?	Customer ID <b>089045</b>					Remarks, Calculations, Results: <b>DUPLICATE ANALYSIS</b> <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.01^{-2}$ g/ml						Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5		K. Saubelick	Hrs	Hrs	Hrs	Hrs		Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>		SJ-8908-081 (R-10-83)					
Serial No. F 117.-8100	Sample Point SEGMENT-18	Date 11-15-89	Time Issued 10:57	Priority 23																																																	
Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0																																																	
Sample Size ?	Customer ID <b>089045</b>																																																				
Remarks, Calculations, Results: <b>DUPLICATE ANALYSIS</b> <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.01^{-2}$ g/ml																																																					
Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5																																																	
K. Saubelick	Hrs	Hrs	Hrs	Hrs																																																	
Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>																																																	
SJ-8908-081 (R-10-83)																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Serial No. F 75.-8300</td> <td>Sample Point SEGMENT-24</td> <td>Date 11-15-89</td> <td>Time Issued 10:22</td> <td>Priority 18</td> <td></td> </tr> <tr> <td>Determination ACD-DGST</td> <td>Method/Standard LA-505-159</td> <td>Result Units g/ml sj</td> <td>Charge Code WB75L</td> <td>Reruns 0</td> <td></td> </tr> <tr> <td>Sample Size ?</td> <td colspan="4">Customer ID <b>089045</b></td> <td></td> </tr> <tr> <td colspan="6">           Remarks, Calculations, Results:  <b>REAGENT BLANK</b>  <b>VOLUME ON COMPLETION SJml</b>  <math>Stab-34C11Co</math>  <math>35C11CK</math> </td> </tr> <tr> <td>Analyst - 1 69769</td> <td>Analyst - 2</td> <td>Analyst - 3</td> <td>Analyst - 4</td> <td>Analyst - 5</td> <td></td> </tr> <tr> <td>K. Saubelick</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td></td> </tr> <tr> <td>Date 1/3/90</td> <td>Time Completed</td> <td>Lab Unit Mgr <i>CJW</i></td> <td><i>OK</i></td> <td><i>SJONES</i></td> <td></td> </tr> <tr> <td colspan="6" style="text-align: right;">SJ-8908-081 (R-10-83)</td> </tr> </table>						Serial No. F 75.-8300	Sample Point SEGMENT-24	Date 11-15-89	Time Issued 10:22	Priority 18		Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0		Sample Size ?	Customer ID <b>089045</b>					Remarks, Calculations, Results: <b>REAGENT BLANK</b> <b>VOLUME ON COMPLETION SJml</b> $Stab-34C11Co$ $35C11CK$						Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5		K. Saubelick	Hrs	Hrs	Hrs	Hrs		Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>		SJ-8908-081 (R-10-83)					
Serial No. F 75.-8300	Sample Point SEGMENT-24	Date 11-15-89	Time Issued 10:22	Priority 18																																																	
Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0																																																	
Sample Size ?	Customer ID <b>089045</b>																																																				
Remarks, Calculations, Results: <b>REAGENT BLANK</b> <b>VOLUME ON COMPLETION SJml</b> $Stab-34C11Co$ $35C11CK$																																																					
Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5																																																	
K. Saubelick	Hrs	Hrs	Hrs	Hrs																																																	
Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>																																																	
SJ-8908-081 (R-10-83)																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Serial No. F 118.-8200</td> <td>Sample Point SEGMENT-19</td> <td>Date 11-15-89</td> <td>Time Issued 10:57</td> <td>Priority 23</td> <td></td> </tr> <tr> <td>Determination ACD-DGST</td> <td>Method/Standard LA-505-158</td> <td>Result Units % RECOVERY</td> <td>Charge Code WB75L</td> <td>Reruns 0</td> <td></td> </tr> <tr> <td>Sample Size ?</td> <td colspan="4">Customer ID <b>089045</b></td> <td></td> </tr> <tr> <td colspan="6">           Remarks, Calculations, Results:  <b>GRAMS SAMPLE</b>  <b>VOLUME ON COMPLETION SJml</b>  <math>1.01^{-2}</math> g/ml         </td> </tr> <tr> <td>Analyst - 1 69769</td> <td>Analyst - 2</td> <td>Analyst - 3</td> <td>Analyst - 4</td> <td>Analyst - 5</td> <td></td> </tr> <tr> <td>K. Saubelick</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td></td> </tr> <tr> <td>Date 1/3/90</td> <td>Time Completed</td> <td>Lab Unit Mgr <i>CJW</i></td> <td><i>OK</i></td> <td><i>SJONES</i></td> <td></td> </tr> <tr> <td colspan="6" style="text-align: right;">SJ-8908-081 (R-10-83)</td> </tr> </table>						Serial No. F 118.-8200	Sample Point SEGMENT-19	Date 11-15-89	Time Issued 10:57	Priority 23		Determination ACD-DGST	Method/Standard LA-505-158	Result Units % RECOVERY	Charge Code WB75L	Reruns 0		Sample Size ?	Customer ID <b>089045</b>					Remarks, Calculations, Results: <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.01^{-2}$ g/ml						Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5		K. Saubelick	Hrs	Hrs	Hrs	Hrs		Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>		SJ-8908-081 (R-10-83)					
Serial No. F 118.-8200	Sample Point SEGMENT-19	Date 11-15-89	Time Issued 10:57	Priority 23																																																	
Determination ACD-DGST	Method/Standard LA-505-158	Result Units % RECOVERY	Charge Code WB75L	Reruns 0																																																	
Sample Size ?	Customer ID <b>089045</b>																																																				
Remarks, Calculations, Results: <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.01^{-2}$ g/ml																																																					
Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5																																																	
K. Saubelick	Hrs	Hrs	Hrs	Hrs																																																	
Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>																																																	
SJ-8908-081 (R-10-83)																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td>Serial No. F 116.-8000</td> <td>Sample Point SEGMENT-17</td> <td>Date 11-15-89</td> <td>Time Issued 10:57</td> <td>Priority 23</td> <td></td> </tr> <tr> <td>Determination ACD-DGST</td> <td>Method/Standard LA-505-159</td> <td>Result Units g/ml sj</td> <td>Charge Code WB75L</td> <td>Reruns 0</td> <td></td> </tr> <tr> <td>Sample Size ?</td> <td colspan="4">Customer ID <b>089045</b></td> <td></td> </tr> <tr> <td colspan="6">           Remarks, Calculations, Results:  <b>GRAMS SAMPLE</b>  <b>VOLUME ON COMPLETION SJml</b>  <math>1.30^{-2}</math> g/ml  <math>10HC N : 1.27</math> </td> </tr> <tr> <td>Analyst - 1 69769</td> <td>Analyst - 2</td> <td>Analyst - 3</td> <td>Analyst - 4</td> <td>Analyst - 5</td> <td></td> </tr> <tr> <td>K. Saubelick</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td>Hrs</td> <td></td> </tr> <tr> <td>Date 1/3/90</td> <td>Time Completed</td> <td>Lab Unit Mgr <i>CJW</i></td> <td><i>OK</i></td> <td><i>SJONES</i></td> <td></td> </tr> <tr> <td colspan="6" style="text-align: right;">SJ-8908-081 (R-10-83)</td> </tr> </table>						Serial No. F 116.-8000	Sample Point SEGMENT-17	Date 11-15-89	Time Issued 10:57	Priority 23		Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0		Sample Size ?	Customer ID <b>089045</b>					Remarks, Calculations, Results: <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.30^{-2}$ g/ml $10HC N : 1.27$						Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5		K. Saubelick	Hrs	Hrs	Hrs	Hrs		Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>		SJ-8908-081 (R-10-83)					
Serial No. F 116.-8000	Sample Point SEGMENT-17	Date 11-15-89	Time Issued 10:57	Priority 23																																																	
Determination ACD-DGST	Method/Standard LA-505-159	Result Units g/ml sj	Charge Code WB75L	Reruns 0																																																	
Sample Size ?	Customer ID <b>089045</b>																																																				
Remarks, Calculations, Results: <b>GRAMS SAMPLE</b> <b>VOLUME ON COMPLETION SJml</b> $1.30^{-2}$ g/ml $10HC N : 1.27$																																																					
Analyst - 1 69769	Analyst - 2	Analyst - 3	Analyst - 4	Analyst - 5																																																	
K. Saubelick	Hrs	Hrs	Hrs	Hrs																																																	
Date 1/3/90	Time Completed	Lab Unit Mgr <i>CJW</i>	<i>OK</i>	<i>SJONES</i>																																																	
SJ-8908-081 (R-10-83)																																																					

9 1 1 2 3 9 1 3 0 2

## ICP Analysis

Serial No. F 123.-8350	Sample Point SEGMENT-24		Date 11-15-89	Time issued 10:58	Priority 13
Determination ICP	Method/Standard LA-505-151	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? Direct			Customer ID 089045		
Remarks, Calculations, Results: REAGENT BLANK  <i>Complete</i>					
Analyst-1 69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
Hrs R. Southwick	Hrs	Hrs	Hrs	Hrs	
Date 3/12/90	Time Completed	Lab Unit Mgr OP			
54-8800-081 (R-10-83)					

Serial No. F 115.-8550	Sample Point SEGMENT-16		Date 11-15-89	Time issued 10:57	Priority 23
Determination ICP	Method/Standard LA-505-151	Result Units % RECOVERY	Charge Code WB75L	Reruns 0	
Sample Size ? Direct			Customer ID 089045		
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID SIC11A Mag. 5 ICP18 78C1C SST1 SIC11A Multi 77C1C SST3 SIC11A - ICP 25 81B38 SST2 <i>Complete</i>					
Analyst-1 69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
Hrs R. Southwick	Hrs	Hrs	Hrs	Hrs	
Date 3/12/90	Time Completed	Lab Unit Mgr OP			
54-8800-081 (R-10-83)					

Serial No. F 116.-8050	Sample Point SEGMENT-17		Date 11-15-89	Time issued 10:57	Priority 23
Determination ICP	Method/Standard LA-505-151	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results:  <i>Complete</i>					
Analyst-1 69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
Hrs R. Southwick	Hrs	Hrs	Hrs	Hrs	
Date 3/12/90	Time Completed	Lab Unit Mgr OP			
54-8800-081 (R-10-83)					

Serial No. F 117.-8150	Sample Point SEGMENT-18		Date 11-15-89	Time issued 10:57	Priority 23
Determination ICP	Method/Standard LA-505-151	Result Units PPM	Charge Code WB75L	Reruns 0	
Sample Size ? 100-10			Customer ID 089045		
Remarks, Calculations, Results: DUPLICATE SAMPLE  <i>Complete</i>					
Analyst-1 69769	Analyst-2	Analyst-3	Analyst-4	Analyst-5	
Hrs R. Southwick	Hrs	Hrs	Hrs	Hrs	
Date 3/12/90	Time Completed	Lab Unit Mgr OP			
54-8800-081 (R-10-83)					

9 1 1 2 0 6 9 1 3 0 3

## ICP Analysis

Serial No F 119.-8550      Sample Point SEGMENT-20      Date 11-15-89      Time Issued 10:58      Priority 23 Determination ICP      Method/Standard LA-505-151      Result Units % RECOVERY      Charge Code WB75L      Returns 0 Sample Size ? Direct      Custom ID 089045				
Remarks, Calculations, Results: LMCS CHECK SAMPLE LMCS ID 7801C SST 1 7741C SST 3 81838A SST 2 <i>Complete</i>				
Analyst - 1 69769 <i>D. Southwick</i>	Analyst - 2 Hrs <i>00</i>	Analyst - 3 Hrs <i>00</i>	Analyst - 4 Hrs <i>00</i>	Analyst - 5 Hrs <i>00</i>
Date 3/12/90	Time Completed Lab Unit Mgr			
SA-8400-061 (R-10-83)				

Serial No F 118.-8250      Sample Point SEGMENT-19      Date 11-15-89      Time Issued 10:57      Priority 23 Determination ICP      Method/Standard LA-505-151      Result Units PPM      Charge Code WB75L      Returns 0 Sample Size ? 100-10      Custom ID 089045				
Remarks, Calculations, Results: DUPLICATE SAMPLE <i>Spike sample</i> <i>Complete</i>				
Analyst - 1 69769 <i>D. Southwick</i>	Analyst - 2 Hrs <i>00</i>	Analyst - 3 Hrs <i>00</i>	Analyst - 4 Hrs <i>00</i>	Analyst - 5 Hrs <i>00</i>
Date 3/12/90	Time Completed Lab Unit Mgr			
SA-8400-061 (R-10-83)				