

Lockheed Environmental Systems & Technologies Co.
Lockheed Analytical Services
975 Kelly Johnson Drive Las Vegas, Nevada 89119-3705
Telephone 702-361-0220 800-582-7605 Facsimile 702-361-8146

LK7008-LAS
0045077

LOCKHEED MARTIN

May 23, 1996

Ms. Joan Kessner
Bechtel Hanford, Inc.
3350 George Washington Way
MS B1-35
Richland, WA 99352



RE: Log-in No.:
Quotation No.:
SAF:
Document File No.:
BHI Document Control No.:
SDG No.:

L7008
Q400000-B
B96-036
0510596
363
LK7008



The attached data report contains the analytical results of samples that were submitted to Lockheed Analytical Services on 10 May 1996. The temperature of the cooler upon receipt was 2°C. Sample containers received agree with the chain-of-custody documentation. Sample containers were received intact. Samples were received in time to meet the analytical holding time requirements.

The case narratives included in the following attachments provide a detailed description of all events that occurred during sample preparation, analysis, and data review specific to the samples and analytical methods requested.

A list of data qualifiers, chain-of-custody forms, sample receiving checklist, and log-in report are also enclosed representing the samples received within this group.

If you have any questions concerning the analysis or the data please call Kathleen M. Hall at (509) 375-4741.

"I certify that this data package is in compliance with the SOW, both technically and for completeness, for other than the conditions detailed above. Release of the data contained in this hard copy data package has been authorized by the Laboratory Manager or a designee, as verified by the following signature."

Sincerely,

Kathleen M. Hall
Client Services Representative

cc: Client Services
Document Control

0003

CASE NARRATIVE RADIOCHEMICAL ANALYSES

The routine calibration and quality control (QC) analyses performed for this batch include as applicable: instrument calibration, initial and continuing calibration verification, quench monitoring standards, instrument background analysis, method blanks, yield tracer, laboratory control samples, matrix spike samples, and duplicate samples.

NOTE: Chemical recoveries and minimum detectable activities (MDAs) can be found on the preparation sheets and calculation sheets of the attached raw data.

Holding Time Requirements

All holding time requirements were met.

Analytical Method Total Strontium

The total strontium analysis was performed using standard operating procedure, LAL-92-SOP-0196. The samples were analyzed in workgroup 36899. The instrument calibration verification met criteria. The method blank was within QC criteria. The laboratory control sample recovery was within QC criteria. The duplicate (DUP) recoveries were within QC criteria. Sample BOHDY4 (L7008-2) and the DUP (36899DUP1) had chemical recoveries that were slightly out of criteria. Since all other QC criteria were within limits and all MDAs were below the reporting detection limit, data quality is not believed to be adversely affected. The client sample was flagged with a "Y" qualifier. No re-analyses were performed.

Yvonne M. Jacoby
Prepared By

May 22, 1996
Date

Lockheed Analytical Services
DATA QUALIFIERS FOR RADIOCHEMICAL ANALYSES
[Revised 04/05/96]

For Use on the Analytical Data Reporting Forms	
B	Any constituent that was detected in the associated method blank at a concentration was greater than the reporting detection limit (RDL).
C	The minimum detectable activity exceeded the RDL due to the residue weight limitations forcing a volume reduction.
D	Constituent detected in the diluted sample.
E	Constituent concentration exceeded the calibration or attenuation curve range.
F	<i>For Alpha Spectrometry Only</i> -- Full width half max exceeded the acceptance limits.
H	Sample analysis performed outside of method-specified maximum holding time requirement.
Y	Chemical yield exceeded acceptance limits.
For Use on the QC Data Reporting Forms	
*	QC data (i.e., percent recovery data for laboratory control standard and matrix spike; and RPD for replicate analyses) exceeded acceptance limits.
a¹	The spike recovery and/or RPD for matrix spike and duplicates cannot be evaluated due to insufficient spiking level compared to the elevated sample analyte concentration.
b¹	The RPD cannot be computed because the sample and/or duplicate concentration was below the MDA.

¹ Used as foot note designations on the QC summary form.

LOCKHEED ANALYTICAL SERVICES
 LOGIN CHAIN OF CUSTODY REPORT (ln01)
 May 10 1996, 05:07 pm

Login Number: L7008
 Account: 596 Bechtel Hanford, Inc. * Richland, WA
 Project: BECHTEL-HANFORD Bechtel Hanford Project

Laboratory Sample Number	Client Sample Number	Collect Date	Receive Date	Due PR Date
L7008-1 TEMP 2 Location: 157 Water 1 S SCREENING	BOHDY4	08-MAY-96	10-MAY-96	25-MAY-96
		Hold:04-NOV-96		
L7008-2 TEMP 2 Location: 157 Water 1 S SR-90 LAL-0196	BOHDY4	08-MAY-96	10-MAY-96	25-MAY-96
		Hold:04-NOV-96		
L7008-3 TEMP 2 Location: 157	BOHDY4	08-MAY-96	10-MAY-96	25-MAY-96
L7008-4 TEMP 2 Location: 157	BOHDY4	08-MAY-96	10-MAY-96	25-MAY-96
L7008-5 TEMP 2 Location: 157	BOHDY4	08-MAY-96	10-MAY-96	25-MAY-96
L7008-6 Location: Water 1 S EDD - DISK DEL. Water 1 S RAD RPT TYPE 2	REPORT TYPE	10-MAY-96	10-MAY-96	25-MAY-96

Page 1

Signature: Paul J. [Signature]

Date: 5-10-96

0006

050596

Bechtel Hanford, Inc.		CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST						L7008		Page <u>1</u> of <u>1</u>	
Collector <i>A. Rizzo, M. Mehlhorn</i>				Company Contact J.V. Borghese				Telephone (509) 373-4790			
Project Designation 100-NR-2 Monthly Performance Monitoring - May				Sampling Location 100 N				SAF No. B96-036			
Ice Chest No. SML-569				Field Logbook No. EL-1288				Method of Shipment Hand Delivered			
Shipped To Lockheed				Offsite Property No. W96-0-0640-44				Bill of Lading/Air Bill No. 2904656526			
Possible Sample Hazards/Remarks				Preservation		HNO ₃		Cool 4°C			
				Type of Container		P/G		P/G			
				No. of Container(s)		4		1			
Special Handling and/or Storage Maintain samples between 2°C and 6°C.				Volume		1L		20mL			
SAMPLE ANALYSIS				Sr-90		Activity Scan					
Sample No.	Matrix*	Date Sampled	Time Sampled								
BOHDY4	W	5-8-96	1025				X	X			
CHAIN OF POSSESSION				SPECIAL INSTRUCTIONS				Matrix*			
Relinquished By		Date/Time		Received By		Date/Time		<div style="font-size: 0.8em;"> S = Soil SE = Sediment SO = Solid SL = Sludge W = Water O = Oil A = Air DS = Drum Solids DL = Drum Liquids T = Tissue WI = Wipe L = Liquid V = Vegetation X = Other </div>			
<i>Monty Mehlhorn</i>		5-8-96		<i>Eric</i>		1300					
<i>Monty Mehlhorn</i>		1300		<i>Eric</i>		5-8-96					
<i>Eric</i>		0900		<i>Eric</i>		5-9-96					
Relinquished By		Date/Time		Received By		Date/Time					
Relinquished By		Date/Time		Received By		Date/Time					
Relinquished By		Date/Time		Received By		Date/Time					
LABORATORY SECTION		Received By		Title		Date/Time					
<i>Eric</i>		<i>Eric</i>		Sample Custodian		5-10-96/1000					
FINAL SAMPLE DISPOSITION		Disposal Method		Disposed By		Date/Time					

LOCKHEED MARTIN

Sample Login Login Review Checklist

Lot Number L2008

The login review should be conducted by that person logging in the samples as well as a peer. Please use this checklist to ensure that such reviews occur in a uniform basis. Please sign and date below to verify that a login review has occurred. This checklist should be affixed to each login package prior to distribution.

For effective login review, at a minimum, five reports from the login process are required. These are the COC (or equivalent), the login COC report, the sample summary report, the sample receiving checklist, and the login quotation. Before beginning review, ensure that these five components are available. Jobs with single component samples, the sample summary report may be omitted.

SAMPLE SUMMARY REPORT

	<u>YES</u>	<u>NO</u>	<u>N/A</u>	<u>Comment</u>
1. Are all sample ID's correct?	<u>X</u>	—	—	_____
2. Are all samples present?	<u>X</u>	—	—	_____
3. Are all matrices indicated correctly?	<u>X</u>	—	—	_____
4. Are all analyses on the COC logged in for the appropriate samples?	<u>X</u>	—	—	_____
5. Are all analyses logged in for the correct container?	<u>X</u>	—	—	_____
6. Are samples logged in according to LAS batching procedures?	<u>X</u>	—	—	_____

LOGIN CHAIN OF CUSTODY

	<u>YES</u>	<u>NO</u>	<u>N/A</u>	<u>Comment</u>
1. Are the collect, receive, and due dates correct for every sample?	<u>X</u>	—	—	_____
2. Have all appropriate comments been indicated in the comment section?	<u>X</u>	—	—	_____

SAMPLE RECEIVING CHECKLIST

	<u>YES</u>	<u>NO</u>	<u>N/A</u>	<u>Comment</u>
1. Are all discrepancies between the COC and the login noted (if applicable)?	—	—	<u>X</u>	_____

Paul D. Davis
primary review signature

5-10-96
date

John M. Miller
secondary review signature

5-10-96
date

0008

0570596

Lockheed Analytical Services
Sample Receiving Checklist

Page 01

Client Name: Bentley/Hanford

Job No. L7008

Cooler ID: mt

COOLER CONDITION UPON RECEIPT

Temperature of cooler upon receipt: 20

temperature of temp. blank upon receipt: -

	Yes	No	* Comments/Discrepancies
custody seals intact	X		
chain of custody present	X		
blue ice (or equiv.) present/frozen	X		
rad survey completed	X		

SAMPLE CONDITION UPON RECEIPT

	Yes	No	* Comments/Discrepancies
all bottles labeled	X		
samples intact	X		
proper container used for sample type	X		
sample volume sufficient for analysis	X		
proper pres. indicated on the COC	X		
VOA's contain headspace			<u>N/A</u>
are samples bi-phasic (if so, indicate sample ID'S):			<u>N/A</u>

MISCELLANEOUS ITEMS

	Yes	No	* Comments/Discrepancies
samples with short holding times		X	
samples to subcontract			<u>N/A</u>

ADDITIONAL COMMENTS/DISCREPANCIES

Completed by / date: Paul Davis 5-10-96

Sent to the client (date/initials):

** Client's signature upon receipt:

Notes: * = contact the appropriate CSR of any discrepancies immediately upon receipt

** = please review this information and return via facsimile to the appropriate CSR (702) 361-8146

05105
609

Lockheed Analytical Laboratory
SAMPLE SUMMARY REPORT (su02)
Bechtel Hanford, Inc. * Richland, WA

Client Sample Number	LAL Sample Number	SDG Number	Matrix	Method
BOHDY4 -	L7008-1 L7008-2		Water Water	SCREENING SR-90 LAL-0196
REPORT TYPE -	L7008-6 L7008-6		Water Water	EDD - DISK DEL. RAD RPT TYPE 2

0010

051054

LOCKHEED ANALYTICAL SERVICES

RADIOCHEMISTRY DATA REPORT

Account Name: Bechtel Hanford, Inc. * Richland, WA

Project Name: BECHTEL-HANFORD

Project Desc: Bechtel Hanford Project

Client Sample ID: B0H DY4

Date Collected: 08-MAY-96

Matrix: Water

Login Number: L7008

Date Received: 10-MAY-96

Constituent	Method	Batch	Activity	Error	MDA	Qualifier	Units	Analyzed	Lab ID
Sr-89,90	LAL-0196	36899	2.7	1.0	1.5	Y	pCi/L	20-MAY-96	L7008-2

LOCKHEED ANALYTICAL LABORATORY

SAMPLE PREPARATION LOG FOR STRONTIUM ANALYSIS TOTAL RADIOSTRONTIUM - LAL-91-SOP-0196

Date Prep Started : 5-16-96
WorkGroup : SR-90 LAL-0196 36899

Matrix : Water
Prep Due Date : 05/20/96

CUSTOMER ID	PARENT LAL ID	NO	QC	CHILD LAL ID	ALIUQUOT VOLUME (g or L)	SR CARRIER (mL)	YTTRIUM SEP DATE	YTTRIUM SEP TIME	PLANCHET TARE WT (grams)	PLANCHET GROSS WT (grams)	* * *	RESIDUE WEIGHT (grams)	COMMENTS						
L7008-2	36899DUP1	1	DUP1	36899-01	500ml	0.5	5-20-96	4:52pm	6.91969	6.92574									
Lab Ctrl Sample	36899LCS1	2	LCS1	36899-02	↓	↓	↓	↓	6.53564	6.54192									
Method Blank	36899MBB	3	MBB1	36899-03	↓	↓	↓	↓	6.85662	6.86281									
BOHDY4	L7008-2	4	SMP1	36899-04	↓	↓	↓	↓	6.85878	6.86441									
		5																	
		6																	
		7																	
		8																	
		9																	
		10																	
		11																	
		12																	
		13																	
		14																	
		15																	
		16																	
		17																	
		18																	
		19																	
		20																	
		21																	
		22																	
		23																	
		24																	
Conc&Vol of Carrier					0.01225g / 0.5ml ; 0.5ml					Act & Vol of LCS					26.78µCi/ml ; 1.0ml				
Carrier Exp Date					9-30-96					LCS Ref Date					4-1-94				
Carrier ID#					94-658-20					LCS ID#					94-677-44-1				
										Prep Anlst		HW							
										Start Date		5-16-96							
										Count Anlst									

Balance Number : 40020021 (✓)
()

Pipette Number : 134488 (✓)
139746 (✓)

Carrier and LCS added by: HW 5-16-96
Witnessed by: _____

Comments :

Analyst : HW 5-20-96

Checked by : HW

Cnt Rm Custody\Date : HW 5/20/96

V96024

0017

LOCKHEED ANALYTICAL LABORATORY

SAMPLE PREPARATION LOG FOR STRONTIUM ANALYSIS

TOTAL RADIOSTRONTIUM - LAL-91-SOP-0196

Date Prep Started : 5/16/96

Matrix : Water

WorkGroup : SR-90 LAL-0196 36899

Prep Due Date : 05/20/96

CUSTOMER ID	PARENT LAL ID	NO	QC	CHILD LAL ID	ALIQOT VOLUME (L)	SR CARRIER (mL)	YTTRIUM SEP DATE	YTTRIUM SEP TIME	PLANCHET TARE WT (grams)	PLANCHET GROSS WT (grams)	* * *	RESIDUE WEIGHT (grams)	COMMENTS
L7008-2	36899DUP1	1	DUP1	36899-01	0.5	0.5	5/20/96	4:52	6.91969	6.92574		0.00605	
Lab Ctrl Sample	36899LCS1	2	LCS1	36899-02	0.5	0.5	5/20/96	4:52	6.53564	6.54192		0.00628	
Method Blank	36899MBB	3	MBB1	36899-03	0.5	0.5	5/20/96	4:52	6.85662	6.86281		0.00619	
BOHDY4	L7008-2	4	SMP1	36899-04	0.5	0.5	5/20/96	4:52	6.85878	6.86441		0.00563	
		5											
		6											
		7											
		8											
		9											
		10											
		11											
		12											
		13											
		14											
		15											
		16											
		17											
		18											
		19											
		20											
		21											
		22											
		23											
		24											
Conc&Vol of Carrier	24.50 mg/mL; 0.5 mL				Act & Vol of LCS		26.78 pCi/mL; 1.0 mL				Prep Anlst		AW
Carrier Exp Date	30-Sep-96				LCS Ref Date		01-Apr-94				Start Date		5/16/96
Carrier ID#	94-658-20				LCS ID#		94-677-44-1				Count Anlst		LV

Balance Number : 40020021 ()

Pipette Number : 134488 ()

Carrier and LCS added by: AW

0018

139746 ()

Witnessed by: MY

Comments :

Analyst : SR FOR AW
DATA ENTRY

Checked by : 

Cnt Rm Custody\Date : _____

V96024

U.S. Environmental Protection Agency
Environmental Monitoring Systems Laboratory-Las Vegas
Nuclear Radiation Assessment Division

Calibration Certificate

Description

Principal radionuclide Strontium-90 Half life 28.6 years
Nominal activity 27 nano curies
Nominal volume 5 ml in ampoule/bottle number 94003-1

Measurement Activity of principal radionuclide

Activity per gram of this solution

5.40 nano curies of Strontium-90
at 0400 hours PST on April 1, 1994

Activity of daughter radionuclide

The principal activity was accompanied at the quoted time by

5.40 nano curies Per gram
of the daughter nuclide Yttrium-90

Total mass of this solution

Approximately 5.0 grams

Method of measurement

The activity of the primary solution was measured
by liquid scintillation counting.

The activity of the dilution was measured by
liquid scintillation counting.

Useful Life

This radionuclide has decayed through 0.0 half lives since it was obtained by EMSL-LV

We recommend that this solution should not be used after

August 1994

This dilution was prepared for the 1994 ASTM
Collaborative Study of a test method for the
determination of Sr-90 in water.

0021

CST 3

Purity

The manufacturer states that activities other than that of the principal nuclide and of its daughter nuclides, if any, were estimated/known to be:

(1) <input type="text"/>	less than equal to <input type="text"/> %	of the principal activity
(2) <input type="text"/>	less than equal to <input type="text"/> %	of the principal activity
(3) <input type="text"/>	less than equal to <input type="text"/> %	of the principal activity

The activity of impurity (1) is not (2) is not (3) is not included in the quoted figures of the principal activity.

Random Errors

The precision of this standard was such that the certified value of the radioactive concentration of the principal activity had a standard error (sm) not greater than $\pm 0.1\%$ (The 99.7% confidence limits are given by $t(sm)$ where t is obtained from the student t factor for the degree of freedom ($n-1$)).

The maximum uncertainty due to the assessable systematic errors (dilution, counting, and known uncertainty of the standard) is obtained by the separate arithmetic summation of the positive and negative systematic error ($+\delta - \delta'$). These have been estimated not to exceed

$+3.8\%$ or -3.8%

the overall uncertainty (often called accuracy) is an estimate of the possible divergence of the quoted result from the true value. It is a combination of random error $[t(sm)]$ at the 99.7% confidence limits and the worst case estimate of the systematic errors ($+\delta - \delta'$). The overall uncertainty is therefore calculated on the basis of $+[t(sm) + \delta] - [t(sm) + \delta']$ and is $+4.0\%$ or -4.0% of the quoted radioactive concentration.

Decay Schemes

This standardization is based on the following assumptions of the principle nuclide, its daughter nuclides and impurities (no allowance for error in these assumptions or the assumption of quoted half-life have been included in the statement of accuracy above).

Strontium-90 decays 100 percent by beta emission to yttrium-90. Yttrium-90 also decays 100 percent by beta emission.

Chemical Composition of Solution

Carrier content per gram of solution:
30 micrograms strontium

Other components:
0.1 M HCl

Preservative:

Remarks

Date Certificate Prepared

April 26, 1994

Approval Signature

0022
Paul B. Hahn

Sr-90
INITIAL STANDARD DILUTION RECORD

Standard Information:

Isotope:

Sr-90

Vendor:

EPA

Activity of Standard Received:

 2.7×10^4 uCi

Vendor I.D. #

94003-1

Weight of Standard Received (g):

50 g

LAL I.D. #:

AC5281

Standard Activity (pCi/g):

 5.4×10^3 pCi/g

NIST Traceable ?

yes

Half-life in Years or Days:

28.6 yrs

Certificate #:

94003-1

Reference Date:

4-1-1994

Receiver's Name:

K. Free

Date Received:

5-3-94

Primary Dilution

Balance Verification?:

yes

Diluent Used:

0.1 M HCl

a: Decay Corrected Standard Activity (pCi/g):

 5.4×10^3 pCi/g

b: Weight of the Source Transferred (g):

4.9670 g

c: Total diluted weight (g):

49.91 g

d: Total Diluted Volume (mL)

50 mL

e: Activity of Dilution by Weight (pCi/g) (a * b / c):

537.4 pCi/g

f: Calculated Density of Solution (g/mL) (c / d):

0.9982 g/mL

g: Activity of Dilution by Volume (pCi/mL) (e * f):

536.44 pCi/mL

h. Dilution Logbook I.D. #:

~~93-474~~ = 81 = 1 ⁹³⁻⁴⁷⁴⁻⁸²⁻¹ _{6/4/95}

Prepared By:

Dyane Wong

Preparation Date:

6-15-94

Reviewed By:

Joe Hutchinson

Review Date:

6/30/94

Purity/Cross Check Performed By:

Check Date:

Signed

Date

Signed

Date

0023

SECONDARY/WORKING LEVEL
STANDARD DILUTION RECORD

Dilution Source Information	
Isotope:	<u>Sr-90</u>
Ref. <u>4-1-94</u>	
Parent Barcode Number	<u>AC5281</u>
Vendor or Certificate I.D. # of Parent Standard:	<u>EPA 94003 - 1</u>
Diluted Source Logbook I.D. #:	<u>93-474 - 82 - 1</u>
Balance Verification?:	<u>Yes</u>
Diluent Used:	<u>0.1 M HCl</u>

Dilution	
*Diluent:	<u>0.1 M HCl</u>
*Density of diluent (g/ml):	<u>N/A</u>
a: Parent Specific Activity:	<u>536.44 pCi/mls</u>
b: Amount of Source Transferred:	<u>5.0018</u> g
c: Total amount of Dilution:	<u>100.20</u> g
d: Total Volume of Dilution:	<u>N/A</u>
e: Activity of Dilution (a * b / c):	<u>N/A</u>
f: Activity of Dilution (a * b / d):	<u>26.78 pCi/ml</u>
Dilution Logbook I.D. #:	<u>94-677-44 - 1</u>
Prepared By: <u>Dynes Wong</u>	Preparation Date: <u>3-2-95</u>
Reviewed By: <u>Joe H. Hines</u>	Review Date: <u>3/3/95</u>

*If the diluent remains unchanged from the diluent used for the dilution source, then a weight dilution of a volume unit source can be performed without a density conversion. If the diluent changes, a weighted proportion density conversion is necessary.

0024

Signed

Date

Signed

Date

Strontium Carrier Standardization

Strontium Carrier (10 mg/mL):

Use commercially available 10,000 μg Sr/mL ICP Standard or equivalent. Alternately, Dissolve 24.16 g of $\text{Sr}(\text{NO}_3)_2$ in water and dilute to 1 L in a volumetric flask with water.

Perform calibration check on a 0.5 mL pipet and then carefully pipet 3 - 0.5 mL portions of the strontium carrier solution into separate cleaned dried and tared planchets. Dry the planchet under a drying lamp. Cool the planchets in a desiccator and weigh.

	Calib # 1	Calib # 2	Calib # 3
Carrier plus planchet wt. (mg)	6865.84	6562.93	6540.56
Tare wt. of planchet (mg)	6853.57	68 6550.76	6528.26 6528.26
Net wt. of carrier added (mg)	12.27	12.17	12.30

AVERAGE $\text{Sr}(\text{NO}_3)_2 \pm \text{STD DEV.} = 12.25 \pm 0.0681$ $\left(0.01225\text{g} \pm 0.0001\text{g} \right)$
Per 0.5 mL

Expected mg of $\text{Sr}(\text{NO}_3)_2 = \text{cert. value} (=10\text{mg of Sr/mL}) * 0.5 \text{ mL} * 2.41$

Within 3% of expected (12.08 mg/0.5 mL) value (yes/no) yes

Initial and Date: W.L. 10-5-95

94-643-20

Walter L. 10-11-95

Continued on Page

Read and Understood By

0025

W.L.
Signed

10-4-95
Date

Dynus Wong
Signed

1-16-96
Date



CERTIFICATE OF ANALYSIS

Catalog Number: PLSR2-3X Lot No. G3-153SR
Element and Matrix: Sr/HNO₃/H₂O
Starting Material: Strontium Carbonate SrCO₃
Starting Material Lot Number: 05941D

Trace Metallic Impurities in the actual solution via ICP of the concentrate.

Element	PPM	Element	PPM	Element	PPM
Ag	<.05	Cu	<.05	P	<.05
Al	<.05	Fe	<.05	Pb	.06
As	<.05	Ga	<.05	Sb	<.10
B	<.05	In	<.05	Si	<.05
Ba	<.05	K	<.05	Sn	<.10
Be	<.05	Li	<.05	Ti	<.05
Bi	<.05	Mg	<.05	Tl	<.05
Ca	<.05	Mn	<.05	V	<.05
Cd	<.05	Mo	<.10	Zn	<.05
Co	<.05	Na	<.05	Zr	<.10
Cr	<.05	Ni	<.05		

Traceability Documentation For Solution Standard:

1. Classical Wet Assay: 10,050 ppm.

Titrimetry: EDTA titration using Methyl Thymol Blue as indicator. EDTA standardized against NIST Pb(NO₃)₂ SRM 928.

2. Instrumentation Analysis By Inductively Coupled Plasma Spectrometer[ICP]: 10,009 ppm via NIST SRM 3153a.

3. Balances are calibrated with NIST weight sets N.J. #92589 and #92550, according to NIST circular 547 3.4.3.

SPEX plasma solution standards are guaranteed stable and accurate to $\pm 0.5\%$ of labeled concentration for one year from date of shipment. This value is the sum of cumulative errors associated with analytical determinations, pipetting and diluting to final volume. For these solutions we use high purity acids, 18 megohm double deionized water and triple rinsed bottles. All glassware used is class A.

Signed by:

N. Kochendakota

Chemical Production Manager,

Date:

SEP 95



Printed on Recycled Paper