

DEC 07 '95

## ENGINEERING DATA TRANSMITTAL

Page 1 of 1

1. EDT NO 613475

2. To: (Receiving Organization) Distribution		3. From: (Originating Organization) Characterization Plans, Coordination and Reports		4. Related EDT No.: N/A							
5. Proj./Prog./Dept./Div.: Tank 241-BX-112/Waste Management/CPCR/Technical Basis Characterization		6. Cog. Engr.: John M. Conner		7. Purchase Order No.: N/A							
8. Originator Remarks: This document is being released into the supporting document system for retrievability purposes.				9. Equip./Component No.: N/A							
				10. System/Bldg./Facility: N/A							
11. Receiver Remarks: For release.				12. Major Assm. Dwg. No.: N/A							
				13. Permit/Permit Application No.: N/A							
				14. Required Response Date: 12/01/95							
15. DATA TRANSMITTED											
(A) Item No.	(B) Document/Drawing No.	(C) Sheet No.	(D) Rev. No.	(E) Title or Description of Data Transmitted	(F) Approval Designator						
1	WHC-SD-WM-DP-157	N/A	0	45-Day Safety Screening Results for Tank 241-BX-112, Auger Samples 95-AUG- 047 and 95-AUG-048	Q						
16. KEY											
Approval Designator (F)		Reason for Transmittal (G)		Disposition (H) & (I)							
E, S, Q, D or N/A (see WHC-CM-3-5, Sec.12.7)		1. Approval 2. Release 3. Information 4. Review 5. Post-Review 6. Dist. (Receipt Acknow. Required)		1. Approved 2. Approved w/comment 3. Disapproved w/comment 4. Reviewed no/comment 5. Reviewed w/comment 6. Receipt acknowledged							
(G)	(H)	17. SIGNATURE/DISTRIBUTION (See Approval Designator for required signatures)				(G)	(H)				
Reason	Disp.	(J) Name	(K) Signature	(L) Date	(M) MSIN	(J) Name	(K) Signature	(L) Date	(M) MSIN	Reason	Disp.
2	1	Cog.Eng. J.M. Conner	<i>J.M. Conner</i>	12-5-95							
2	1	Cog. Mgr. J.G. Kristofzski	<i>J.G. Kristofzski</i>	12-7-95							
2	1	QA T.L. Tung	<i>T.L. Tung</i>	12/6/95							
		Safety									
		Env.									
2	1	Program Support A.D.Rice	<i>A.D. Rice</i>	12-6-95							
18. A.E. Young <i>A.E. Young</i> Signature of EDT Originator		19. N/A Authorized Representative Date for Receiving Organization		20. <i>J.G. Kristofzski</i> J.G. Kristofzski Cognizant Manager		21. DOE APPROVAL (if required) Ctrl. No. <input type="checkbox"/> Approved <input type="checkbox"/> Approved w/comments <input type="checkbox"/> Disapproved w/comments					

## 45-Day Safety Screening Results for Tank 241-BX-112, Auger Samples 95-AUG-047 and 95-AUG-048

John M. Conner

Westinghouse Hanford Company, Richland, WA 99352  
U.S. Department of Energy Contract DE-AC06-87RL10930

EDT/ECN: EDT-613475 UC: 2070  
Org Code: 75310 Charge Code: MD378  
B&R Code: EW 3120074 Total Pages: 38


Key Words: 45-Day, Safety Screening Results, Safety, Screening, Results, Tank 241-BX-112, Tank BX-112, BX-112, BX Farm, Auger Samples, 95-AUG-047, 95-AUG-048

Abstract: N/A


---

TRADEMARK DISCLAIMER. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof or its contractors or subcontractors.

Printed in the United States of America. To obtain copies of this document, contact: WHC/BCS Document Control Services, P.O. Box 1970, Mailstop H6-08, Richland WA 99352, Phone (509) 372-2420; Fax (509) 376-4989.

  
Release Approval

  
Date



**Approved for Public Release**



Westinghouse  
Hanford Company

---

P.O. Box 1970 Richland, WA 99352

WHC-SD-WM-DP-157, REV. 0

**ANALYTICAL SERVICES**

**45-DAY SAFETY SCREENING RESULTS FOR TANK 241-BX-112,  
AUGER SAMPLES 95-AUG-047 AND 95-AUG 048**

**DATE PRINTED:**

**NOVEMBER 30, 1995**

WHC-SD-WM-DP-157, REV. 0

TABLE OF CONTENTS

Narrative . . . . .	1
Sample Data Summary . . . . .	6
Inorganic Analyses . . . . .	9
Differential Scanning Calorimetry (DSC)	
DSC Worklist # 3729 . . . . .	10
DSC Worklist # 3735 . . . . .	14
DSC Worklist # 4000 . . . . .	20
Thermogravimetric Analysis (TGA)	
TGA Worklist # 3730 . . . . .	21
TGA Worklist # 3734 . . . . .	25
TGA Worklist # 4003 . . . . .	31

This document consists of pages 1 through 35.

**WHC-SD-WM-DP-157, REV. 0**

**NARRATIVE**

**45-DAY SAFETY SCREENING REPORT FOR TANK 241-BX-112,  
AUGER SAMPLES 95-AUG-047 AND 95-AUG-048**

ANALYTICAL SUMMARY

Two auger samples were taken from tank 241-BX-112 (BX-112). The samples were received at the 222-S Laboratories and underwent safety screening analyses, consisting of differential scanning calorimetry (DSC), thermogravimetric analysis (TGA), bulk density, and determination of total alpha activity. As appropriate, the results were compared to the safety screening limits at a confidence level of 95%. All analytical results were within the action limits stated in the SAP.

Based on the results of vapor monitoring prior to sampling, the BX-112 vapor space is far below the lower explosive limit (LEL). None of the data indicate that the tank is "unsafe" when compared to the criteria (energetics, criticality, and flammability) in the *Safety Screening Data Quality Objective* (Dukelow, et al., 1995). However, the tank cannot be declared "safe," as two full length profiles were not obtained by auger sampling. Core sampling will be necessary to fully satisfy the DQO.

Any additional analyses on the BX-112 auger samples will be included in a revision to this report.

SCOPE

This document serves as the 45-day report deliverable for the tank BX-112 auger samples collected on November 16 and 17, 1995 (samples 95-AUG-047 and 95-AUG-048). The 222-S Laboratories received, extruded, and analyzed each sample in accordance with the SAP [1]. Included in this report are the primary safety screening results obtained from the analyses, and copies of all DSC and TGA raw data scans as requested in the SAP. The results of tank dome space flammability screening are also included.

No additional testing to support safety screening analyses for these auger samples is required. Any additional analyses conducted by the 222-S Laboratories on these auger samples will be included in a revision to this report.

TANK DOME SPACE FLAMMABILITY SCREENING

Prior to auger sampling, the vapor space of tank BX-112 was screened for flammability issues. The results of combustible gas monitoring inside riser 3 of tank BX-112 are presented in Table 1. This measurement is conducted in the field and recorded in the work package (work package for BX-112 auger sampling is #ES-95-00217). The results indicated that the tank vapor space was at 0% of the LEL, far below the action limit of 10% stated in the DQO [2].

Table 1. Characteristics of BX-112 Vapor Space as Determined by Combustible Gas Monitoring.

Measurement	Result
Lower Explosive Limit (LEL)	0%
Oxygen (O <sub>2</sub> )	20.9%
Total Organic Carbon (TOC)	11.3 ppm
Ammonia (NH <sub>3</sub> )	125 ppm

SAMPLE RECEIPT, EXTRUSION, AND SUBSAMPLING95-AUG-047

Auger sample 95-AUG-047 was collected from riser 3 of tank BX-112 on November 16, 1995, and extruded on November 21, 1995. This was a 50 cm (20 inch) auger sampler, with 39 flutes (this was the first time that this auger was used. It has twice as many flutes as the typical auger). The sample appeared fairly homogeneous. The waste was a very wet, light-brown sludge, which tended to drip from the auger onto the extrusion tray. Flutes 1-8 at the top of the auger were bare. The material on flutes 9-16 appeared slightly less wet than the lower portion of the sample. The material on flutes 9-16, totaling 45.8 grams, was segregated as the upper half-segment solids. Flutes 17-39 contained lumpy sludge (lumps disappeared upon subsampling). A total of 165.2 grams was subsampled as the lower half-segment solids. The half segment subsamples were homogenized and subsampled for further laboratory analyses, bulk density determination, and archiving.

95-AUG-048

Auger sample 95-AUG-048 was collected from riser 2 of Tank BX-112 on November 17, 1995, and extruded on November 21. This was a 50 cm (20 inch) auger sampler, with 19 flutes. The sample appeared to be a homogeneous, medium-brown, very wet sludge. Flutes 1-6 were bare. The sample was recovered on flutes 7-19, and mostly dripped off of the auger onto the extrusion tray. Due to the apparent length of the sample recovered [30 cm (12 inches)], the moderate-to-low recovery (81.3 g), and the apparent homogeneity, the sample was not subsampled into half segments, but homogenized and subsampled at the whole-segment level. Portions were then subsampled for bulk density determination and further laboratory analyses and archiving.

## ANALYTICAL RESULTS

### BULK DENSITY

Three subsamples were submitted for bulk density determination by centrifugation in a tared, graduated, vial per procedure LA-160-103, Rev. A-7. The results ranged from 1.31 to 1.35 g/cm<sup>3</sup>. These results are presented in the summary tables. In order to conserve sample, duplicate analyses were not conducted.

### THERMOGRAVIMETRIC ANALYSIS (TGA)

Three samples were submitted for moisture content determination by TGA per procedures LA-560-112, Rev. B-2, or LA-514-114, Rev. C-1 (a different procedure is used for each instrument). The samples were analyzed in duplicate. The results are presented in the summary tables, and the raw data scans are attached. All results were between 55.59 and 65.50 percent moisture. The relative percent difference (RPD) between sample and duplicate results for sample S95T003746 was 11.4%, which slightly exceeded the criterion of less than 10% given in the SAP. Inspection of the raw data (attached) indicates that the sample and duplicate scans are similar in shape, except that the weight loss for the sample result appears to begin at approximately 100 °C, instead of at ambient temperature. The chemist attributed this to static charge holding the sample tray to the side of the furnace [3]. Once the static charge was overcome, the weight loss scan appears very similar in shape to the scan of the duplicate sample (although the endpoint differs by almost 7 weight percent). The sample was rerun in duplicate. The RPD for the rerun was 1.71%. Both the original results and the results of the rerun are included in the summary tables. The results of the rerun are noted by a "1" next to the sample number.

### DIFFERENTIAL SCANNING CALORIMETRY (DSC)

Three samples were submitted for determination of energetics by DSC per procedure LA-514-113, Rev. C-1 or procedure LA-514-114, Rev. C-1. The samples were analyzed in duplicate. The results are presented in the summary tables, and the raw data scans are attached. None of the samples exhibited exotherms. Since none of the samples exhibited any exotherms, the statistical calculation of an upper 95% confidence level for each sample is unnecessary.

### ALPHA TOTAL

Three solids samples were submitted for total alpha analysis per procedure LA-508-101, Rev. D-2. The samples were fused per procedure LA-549-141, Rev. D-0 prior to analysis. Two fusions were prepared per sample (for duplicate results). Each fused dilution was analyzed twice; the results were averaged and reported as one value. The highest result returned was 0.219  $\mu\text{Ci/g}$ , more than two orders of magnitude below the action limit of 41  $\mu\text{Ci/g}$ . The upper 95% confidence level for each sample has been calculated and is presented in Table 2. All of the adjusted results are far below the action limit of 41



$\mu\text{Ci/g}$  stated in the SAP. The RPD for sample S95T003751 was 22.3%. However, no rerun was deemed necessary, as the 95% confidence limit upper value for this sample was 0.336  $\mu\text{Ci/g}$ , far below the action limit of 41  $\mu\text{Ci/g}$ .

The alpha results reported by the lab are calculated assuming a density of 1.5  $\text{g/cm}^3$  for solid samples. As the bulk density results recorded for these samples are all below 1.5  $\text{g/cm}^3$ , the alpha results reported remain conservative per the calculation described in the SAP [1].

One of the two standards run with these samples exhibited a recovery slightly outside the range specified in the SAP (110.2%). Since the result was so close to being within range, and the sample results were far below the limit, a rerun was deemed unnecessary. This result was well within the method control limits of 72.3-125.9%. All quality control results are presented in the summary tables.

Table 2. Comparison of Total Alpha Results at a Confidence Level of 95%.

Sample Description/ Sample Number	Sample Result	Duplicate Result	Mean	Var(Mean)	Upper 95% Confidence Limit
AUG-047 UH S95T003747	0.187	0.178	0.182	2.025E-05	0.211
AUG-047 LH S95T003751	0.175	0.219	0.197	4.84E-04	0.336
AUG-048 WS S95T003755	0.183	0.170	0.176	4.23E-05	0.218
AUG-047, AUG-048 combined	-	-	0.183 (weighted mean)	5.12E-05	0.228

Notes: var(mean) - variance of the mean; UH - upper half;  
LH - lower half; WS - whole segment;  
weighted mean - average for each auger given equal weight

## REFERENCES

- [1] J. M. Conner, *Tank 241-BX-112 Auger Sampling and Analysis Plan*, WHC-SD-WM-TSAP-051, Rev. 0A, Westinghouse Hanford Company, Richland, Washington, November 15, 1995.
- [2] G. T. Dukelow, et al., *Tank Safety Screening Data Quality Objective*, WHC-SD-WM-SP-004, Rev. 2, Westinghouse Hanford Company, Richland, Washington, August 31, 1995.
- [3] Personal Communication with B. D. Valenzuela, December 1, 1995.

WHC-SD-WM-DP-157, REV. 0

#### SAMPLE DATA SUMMARY

# INTERIM

45-Day Report for Auger Samples 95-AUG-047, 95-AUG-048  
BX-112

CORE NUMBER: 95-AUG-047, 95-AUG-048  
SEGMENT #: 95-AUG-047

SEGMENT PORTION: U Upper Half of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T003745			Bulk Density of Sample	g/mL	None	None	n/a	n/a	1.350	n/a	n/a	n/a	n/a	5.00e-01	n/a
S95T003746	1		% Water by TGA on Perkin Elmer	%	None	None	101.0	n/a	60.72	61.77	61.25	1.71	n/a	n/a	n/a
S95T003746			% Water by TGA on Perkin Elmer	%	None	None	100.5	n/a	55.59	62.30	58.95	11.4	n/a	n/a	n/a
S95T003746			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+01	480.0	99.97	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003746			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+01	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003747	F		Alpha of Digested Solid	uCi/g	-1.0e+01	41.01	92.97	<4.33e-03	1.87e-01	1.78e-01	1.82e-01	4.93	n/a	5.47e-03	8.30E+00

L Lower Half of Segment: L Lower Half of Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T003749			Bulk Density of Sample	g/mL	None	None	n/a	n/a	1.310	n/a	n/a	n/a	n/a	5.00e-01	n/a
S95T003750			% Water by TGA on Perkin Elmer	%	None	None	100.5	n/a	63.37	63.52	63.45	0.24	n/a	n/a	n/a
S95T003750			DSC Exotherm on Perkin Elmer	Joules/g	-1.0e+01	480.0	99.97	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003750			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+01	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003751	F		Alpha of Digested Solid	uCi/g	-1.0e+01	41.01	110.2	<3.12e-03	1.75e-01	2.19e-01	1.97e-01	22.3	98.47	3.80e-03	8.51E+00

=> Limit violated  
=> Selected Limit

WHC-SD-WM-DP-157, REV. 0

# INTERIM

45-Day Report for Auger Samples 95-AUG-047, 95-AUG-048  
BX-112

CORE NUMBER: 95-AUG-047, 95-AUG-048  
SEGMENT #: 95-AUG-048

SEGMENT PORTION: W Whole Segment

Sample#	R	A#	Analyte	Unit	Action Limits		Standard %	Blank	Result	Duplicate	Average	RPD %	Spk Rec %	Det Limit	Count Err%
					Lower	Upper									
S95T003753			Bulk Density of Sample	g/mL	None	None	n/a	n/a	1.310	n/a	n/a	n/a	n/a	5.00e-01	n/a
S95T003754			% Water by TGA using Mettler	%	None	None	100.5	n/a	65.42	65.50	65.46	0.12	n/a	n/a	n/a
S95T003754			DSC Exotherm Dry Calculated	Joules/g Dry	-1.0e+01	480.0	n/a	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003754			DSC Exotherm using Mettler	Joules/g	-1.0e+01	480.0	92.44	n/a	0.00e+00	0.00e+00	0.00e+00	0.00	n/a	n/a	n/a
S95T003755	F		Alpha of Digested Solid	uCi/g	-1.0e+01	41.01	92.97	<4.24e-03	1.83e-01	1.70e-01	1.76e-01	7.37	n/a	4.98e-03	7.73E+00

=> Limit violated

=> Selected Limit

α

WHC-SD-WM-DP-157, REV. 0

**WHC-SD-WM-DP-157, REV. 0**

**INORGANIC ANALYSES**

# LABCORE Data Entry Template for Worklist#

3729

Analyst: SME Instrument: DSC01 DSC01 Book # 12N149

Method: LA-514-113 Rev/Mod C-1

Worklist Comment: BX-112 DSC PLease run under N2. PRIORITY.

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-01	SOLID	<u>28.45</u>	<u>26.3</u>	<u>N/A</u>	Joules/g
95000202	BX-112	2 SAMPLE	S95T003754	0	DSC-01	SOLID	<u>N/A</u>	<u>0</u>		Joules/g
95000202	BX-112	3 DUP	S95T003754	0	DSC-01	SOLID	<u>0</u>	<u>0</u>	<u>N/A</u>	Joules/g

Final page for worklist #

3729

[Signature] [Signature] 11/27/95  
Analyst Signature Date

[Signature] 11/27/95  
Analyst Signature Date

Verified by Blandina Valenzuela  
11/27/95

## Data Entry Comments:

Sample has another endotherm beside the endotherm at 249.5 mJ/g  
nature of 3rd g in sample 9 Aug SME

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT  
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 11 TO 13.

DSC STD 12N14A

4.000 mg

Rate: 10.0 °C/min

File: 00015.001

Ident: 0.0

DSC METTLER 21-Nov-95

222-S Laboratory

ΔX

10.00

Integration

Delta H 211 mJ

20.8 J/g

Peak 157.7 °C

-15.5 mW

140.

140.

150.


180.

*Signature* 11.21.95

File: 00082.001 DSC METTLEH 24-Nov-95  
 Ident: 0.0 222-S Laboratory

000T0002754 N2  
 20.920 mg

Rate: 10.0 °C/min




Integration  
 Delta H 112.00 mJ  
 3.8 J/g  
 Peak 249.4°C  
 -1.5 mW





File: 00024.001 DSC METTLER 21-Nov-95  
 Ident: 0.0 222-0 Laboratory

Rate: 10.0 °C/min



# LABCORE Data Entry Template for Worklist#

3735

Analyst: SME Instrument: DSC01 DSC03 Book # 12044

Method: LA-514-114 Rev/Mod C-1

Worklist Comment: PE BX112 DSC

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			DSC-03	SOLID	<u>2845</u>	<u>2844</u>	<u>N/A</u>	Joules/g
95000202	BX-112	2 SAMPLE	S95T003746	0	DSC-03	SOLID	<u>N/A</u>	<u>0</u>		Joules/g
95000202	BX-112	3 DUP	S95T003746	0	DSC-03	SOLID	<u>0</u>	<u>0</u>	<u>N/A</u>	Joules/g
95000202	BX-112	4 SAMPLE	S95T003750	0	DSC-03	SOLID	<u>N/A</u>	<u>0</u>		Joules/g
95000202	BX-112	5 DUP	S95T003750	0	DSC-03	SOLID	<u>0</u>	<u>0</u>	<u>N/A</u>	Joules/g

Final page for worklist #

3735

[Signature] 11/21/95  
Analyst Signature Date

[Signature]  
Analyst Signature Date

Verified by Blandina  
Valenzuela  
11/27/95

S95T003746 produced an endotherm at 119.6°C with a delta H of 1464.2 J/g

Data Entry Comments:

11/27/95  
S95T003750 produced an endotherm at 118.9 113.8°C  
with a delta H of 1333.3 J/g

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,  
R = Replicate Number, A = Aliquot Code.


Curve 1: DSC

File info: IND112101 Tue Nov 21 07:55:04 1995

Sample Weight: 6.860 mg

12N14A Indium at 10C/min

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT  
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 15 TO 19.



N2, EXOTHERM DOWN

TEMP 1: 150.0 °C TIME 1: 0.0 min RATE 1: 10.0 °C/min  
TEMP 2: 170.0 °C

JD SPELLMAN  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 17:32:38 1995


Curve 1: DSC

File info: SAM112103 Tue Nov 21 18:26:33 1995

Sample Weight: 32.490 mg

S95T003746

91  
WHC-SD-WM-DP-157, REV. 0



exotherm down, N2 purge gas

TEMP: 35.0 °C    TIME: 0.0 min    RATE: 10.0 °/min  
TEMP: 500.0 °C

Temperature (°C)

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 20:14:51 1995


WHC-SD-WM-DP-157, Rev. 0

Curve 1: DSC

File info: SAM112104 Tue Nov 21 21:05:32 1995

Sample Weight: 11.160 mg

S95TC03746 DUP



WHC-SD-WM-DP-157, REV.0

exotherm down, N2 purge gas

TEMP: 25.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON

PERKIN-ELMER

7 Series Thermal Analysis System


Tue Nov 21 21:09:18 1995

Curve 1: DSC

File info: SAM112105 Tue Nov 21 22:19:21 1995

Sample Weight: 24.290 mg

S95T003750



exotherm down, N2 purge gas

TEMP: 35.0 °C TIME: 0.0 min RATE: 10.0 °/min

Temperature (°C)

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 22:24:53 1995


Curve 1: DSC

File info: SAM112106 Tue Nov 21 23:21:17 1995

Sample Weight: 41.360 mg

S95T003750 DUP

19



WHC-SD-WM-DP-157, REV. 0

exotherm down, N2 purge gas

TEMP1: 50.0 °C TEMP2: 500.0 °C TIME1: 0.0 min RATE1: 10.0 °C/min

Temperature (°C)

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Wed Nov 22 00:25:42 1995

# LABCORE Data Entry Template for Worklist#

4000

Analyst: BDV Instrument: DSC01 Book #       

Method: LA-514-113 Rev/Mod       

Worklist Comment: Dry DSCs for BX-112. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
95000202	BX-112	1 SAMPLE	S95T003746	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000202	BX-112	2 DUP	S95T003746	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000202	BX-112	3 SAMPLE	S95T003750	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000202	BX-112	4 DUP	S95T003750	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry
95000202	BX-112	5 SAMPLE	S95T003754	0	DSC-02	SOLID	N/A	Ø		Joules/g Dry
95000202	BX-112	6 DUP	S95T003754	0	DSC-02	SOLID	Ø	Ø	N/A	Joules/g Dry

Final page for worklist #

4000

Data entered + verified by  
Blandina Valenzuela  
Analyst Signature

12/1/95  
Date

Analyst Signature

Date

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number, R = Replicate Number, A = Aliquot Code.



# LABCORE Data Entry Template for Worklist#

3730

Analyst: SME Instrument: TGA01 TGA01 Book # 65N8A

Method: LA-560-112 Rev/Mod B-1

Worklist Comment: BX-112 TGA. Please run under N2. PRIORITY.

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	<u>59.74</u>	<u>60.05</u>	<u>N/A</u>	%
95000202	BX-112	2 SAMPLE	S95T003754	0	TGA-01	SOLID	<u>N/A</u>	<u>65.42</u>		%
95000202	BX-112	3 DUP	S95T003754	0	TGA-01	SOLID	<u>65.42</u>	<u>65.50</u>	<u>N/A</u>	%

Final page for worklist #

3730

[Signature] 11/20/95  
Analyst Signature Date

[Signature] 11/22/95  
Analyst Signature Date

Verified by Blandina Valenzuela  
11/27/95

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number,  
R = Replicate Number, A = Aliquot Code.

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT  
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 22 TO 24 .

TGA STD GENEA

PC.764 mg

Rate: 10.0 °C/min

File: 99914.001 TG METTLER 21-Nov-95

Ident: 9.0 222-B Laboratory

Step Analysis

Height: 36.49 mg

-60.05 %

Resid: 24.27 mg

39.95 %

Onset: 138.3 °C

WHC-SD-WM-DP-157, REV. 0

22

100.00

100.

100.

150.

Signature 11-21-95 °C

Q000000754 ME

24.484 mg

Rate: 10.0 °C/min

File: 00023.001

TG

METTLEH

21-Nov-95

Ident: 0.0

222-3 Laboratory

Step Analysis


Height=14.54 mg

65.42 %

Height= 7.69 mg

64.58 %

Drop= 95.0 °C





# LABCORE Data Entry Template for Worklist#

3734

Analyst: SmE Instrument: TGA01 TGA03 Book # 65NSA

Method: LA-514-114 Rev/Mod C-1

Worklist Comment: PE BX-112 tga

GROUP	PROJECT	S TYPE	SAMPLE#	R A	-----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.04</u>	<u>N/A</u>	%
95000202	BX-112	2 SAMPLE	S95T003746	0	TGA-03	SOLID	<u>N/A</u>	<u>55.57</u>		%
95000202	BX-112	3 DUP	S95T003746	0	TGA-03	SOLID	<u>55.57</u>	<u>63.30</u>	<u>N/A</u>	%
95000202	BX-112	4 SAMPLE	S95T003750	0	TGA-03	SOLID	<u>N/A</u>	<u>63.37</u>		%
95000202	BX-112	5 DUP	S95T003750	0	TGA-03	SOLID	<u>63.37</u>	<u>63.52</u>	<u>N/A</u>	%

Final page for worklist #

3734

For signature for Sign Section 11/21/95  
Analyst Signature Date

11/27/95  
Analyst Signature Date

Verified by Blandina Valenzuela  
11/27/95

S95T003746 the flat portion of the thermogram is due to static  
Data Entry Comments: electricity between the sample pan and the furnace.  
the sample will be run again.


Curve 1: TGA

File info: TER112101 Tue Nov 21 07:47:10 1995

Sample Weight: 24.173 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT  
COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 26 TO 30.



N2 100/MIN

TEMP: 35.0 °C TIME: 0.0 min RATE: 10.0 °C/min

JD SPELLMAN  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 17:39:44 1995


WHC-SD-WM-DP-157, REV. 0

Curve 1: TGA

File info: SAM112103 Tue Nov 21 18:30:51 1995

Sample Weight: 32.917 mg

S95T003746



WHC-SD-WM-DP-157, REV.0

100/MIN N2

TEMP: 50.0 °C TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)


SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 19:04:17 1995

Curve 1: TGA

File info: SAM112104 Tue Nov 21 21:08:37 1995

Sample Weight: 17.745 mg

S95T003746 DUP



TIME: 95.0 S TIME: 0.0 min RATE: 10.0 C/min

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 21:14:57 1995

WHC-SD-WM-DP-157 REV 1




Curve 1: TGA

File info: SAM112105 Tue Nov 21 22:15:25 1995

Sample Weight: 16.036 mg

S95T003750



WHC-SD-WM-DP-157, REV. 0

10C/MIN N2  
TEMP: 25.0 °C  
TIME: 0.0 min RATE: 10.0 °C/min

Temperature (°C)

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Nov 21 22:31:24 1995


Curve 1: TGA

File info: SAM112106 Tue Nov 21 23:26:24 1995

Sample Weight: 28.341 mg

S95T003750 DuP

*SMF*  
*11.21.95*



WHC-SD-WM-DP-157, REV. 0

100/MIN N2

TEMP: 35.0 °C TIR#1: 0.0 min RATE: 10.0 °/min

Temperature (°C)

SM FULTON  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Wed Nov 22 00:32:42 1995

**LABCORE Data Entry Template for Worklist#**

**4003**

**Analyst:** JDS **Instrument:** TGA0 3 **Book #** 65N8-A

**Method:** LA-514-114 Rev/Mod C-1

**Worklist Comment:** Please run BX-112 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	TEST	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-03	SOLID	<u>59.74</u>	<u>60.32</u>	<u>N/A</u>	%
95000202	BX-112	2 SAMPLE	S95T003746	1	TGA-03	SOLID	<u>N/A</u>	<u>60.72</u>		%
95000202	BX-112	3 DUP	S95T003746	1	TGA-03	SOLID	<u>60.72</u>	<u>61.77</u>	<u>N/A</u>	%

**Final page for worklist # 4003**

See attached for signatures  
**Analyst Signature**                      **Date** 12/5/95

                     12-6-95  
**Analyst Signature**                      **Date**                     

Verified by Blundina  
Valenzuela  
12/6/95

**Data Entry Comments:** Sample produced a second weight loss step of  
3.976%

# LABCORE Data Entry Template for Worklist#

4003

Analyst: Jds Instrument: TGA0 12/5/95 BBN Book # 65N8-A

Method: LA-560-112 Rev/Mod \_\_\_\_\_

Worklist Comment: Please run BX-112 TGAs under N2. bdv

GROUP	PROJECT	S TYPE	SAMPLE#	R A	-----TEST-----	MATRIX	ACTUAL	FOUND	DL	UNIT
		1 STD			TGA-01	SOLID	_____	_____	N/A	%
95000202	BX-112	2 SAMPLE	S95T003746	0	TGA-01	SOLID	N/A	_____	_____	%
95000202	BX-112	3 DUP	S95T003746	0	TGA-01	SOLID	_____	_____	N/A	%

Final page for worklist #

4003

Jah Spk 12-5-95  
Analyst Signature Date

\_\_\_\_\_  
Analyst Signature Date

Other instrument  
was used.  
12/5/95  
BBN

Data Entry Comments:

Units shown for QC (SPK & STD) may not reflect the actual units. DL = Detection Limit, S = Worklist Slot Number.  
R = Replicate Number, A = Aliquot Code.

Curve 1: TGA


File info: TER120401 Mon Dec 4 18:36:35 1995

Sample Weight: 24.774 mg

65N8-A Terliq

SIGNATURE BELOW REPRESENTS CHEMICAL TECHNOLOGIST/CHEMIST THAT

COMPLETED/VERIFIED THE CALIBRATION/ANALYSIS ON PAGES 33 TO 35.



WHC-SD-WM-DP-157, REV. 0

N2 10C/MIN

TEMP1: 30.0 C TIME1: 0.0 min RATE1: 10.0 C/min  
TEMP2: 300.0 C

Temperature (°C)


JD SPELLMAN  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Dec 5 04:36:24 1995

Curve 1: TGA

File info: SAM120401 Mon Dec 4 20:32:45 1995

Sample Weight: 17.990 mg

S95T003746



WHC-SD-WM-DP-157, REV. 0

10C/MIN N2

TEMP: 25.0 C TIME: 0.0 min RATE: 10.0 C/min

Temperature (°C)


JD SPELLMAN  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Mon Dec 4 20:36:09 1995

Curve 1: TGA

File info: SAM120502 Tue Dec 5 03:42:41 1995

Sample Weight: 21.542 mg

S95T003746 DUP



10C/MIN N2  
TEMP: 35.0 C  
TEMP: 500.0 C  
TIME: 0.0 min RATE: 10.0 C/min

Temperature (°C)

JD SPELLMAN  
PERKIN-ELMER  
7 Series Thermal Analysis System  
Tue Dec 5 04:19:48 1995

WHC-SD-WM-DP-157, REV.0

# DISTRIBUTION SHEET

To Distribution	From Characterization Plans, Coordination and Reports	Page 1 of 1	
		Date: 12/01/95	
Project Title/Work Order WHC-SD-WM-DP-157, Rev. 0, "45-Day Safety Screening Results for Tank 241-BX-112, Auger Samples 95-AUG-047 and 95-AUG-048"		EDT NO.: EDT-613475	
		ECN NO.: N/A	
Name	MSIN	Text With all Attach	EDT/ECN ONLY

Pacific Northwest Laboratory

J. R. Gormsen	K7-28		X
S. J. Harris	K7-22	X	
K. L. Silvers	P7-27		X

U.S. Department of Energy, RL

C. A. Babel	S7-54	X	
-------------	-------	---	--

Westinghouse Hanford Company

J. N. Appel	G3-21		X
H. Babad	S7-30	X	
J. M. Conner	R2-12	X	
G. D. Forehand	S7-31		X
V. W. Hall	T6-03	X	
D. C. Hetzer	S6-31		X
L. Jensen	R2-12	X	
N. W. Kirch	R2-11	X	
M. J. Kupfer	H5-49	X	
J. E. Meacham	S7-15	X	
P. M. Morant	H4-25	X	
K. L. Powell	T6-04		X
L. W. Shelton	H5-49	X	
B. C. Simpson	R2-12		X
J. A. Voogd	H5-03		X
Central Files	A3-88	2	
EDMC	H6-08	X	
LTIC	T6-03		X
TCRC	R2-12	X	
TFIC (Tank Farm Information Center)	R1-20		X

Washington State Department of Ecology

Single-Shell Tank Unit Manager			
A. B. Stone	B5-18	X	

U. S. Department of Energy

Jim Poppiti			X
12800 Middlebrook Rd.			
Trevion II, EM-36			
Germantown, MD 20874			

Los Alamos Technical Associates

A. T. DiCenso		X	
309 Bradley Blvd.			
Richland, WA 99352			

*TTTAA 1*