

**Recra LabNet Philadelphia
Analytical Report
REVISION**

Client : TNU-HANFORD B99-078
RFW# : 9909L101
SDG# : H0536
SAF# : B99-078

W.O. # : 10985-001-001-9999-00
Date Received: 09-03-99


INORGANIC CASE NARRATIVE

RECEIVED
AUG 7 2000

This package was revised to add Appendix A.


EDMC

1. This narrative covers the analyses of 7 soil samples. A method detection limit (MDL) determination was performed for Hydrazine. The information pertaining to the MDL study is provided within the data package.
2. The samples were prepared and analyzed in accordance with the methods checked on the attached glossary.
3. Sample holding times as required by the method and/or contract were met.
4. The method blank for Hydrazine was within method criteria.
5. The Laboratory Control Samples (LCS) for Hydrazine were within the laboratory control limits. The duplicate LCS was within the 20% Relative Percent Difference (RPD) control limit.
6. The replicate analysis for Percent Solids was within the 20% RPD control limit.
7. Results for the MDL study were reported on a wet weight basis.


J. Michael Taylor
Vice President
Philadelphia Analytical Laboratory

njp&pef\09-101

The results presented in this report relate only to the analytical testing and conditions of the samples received. All pages of this report are integral parts of the analytical data. Therefore, this report should only be reproduced in its entirety. 21 pages.



WET CHEMISTRY
METHODS GLOSSARY FOR SOIL/SOLIDS SAMPLE ANALYSIS

	<u>ASTM</u>	<u>SW846</u>	<u>OTHER</u>
% Ash	___ D2216-80		
% Moisture	___ D2216-80		___ ILMO4.0 (e)
% Solids	✓ D2216-80		___ ILMO4.0 (e)
% Volatile Solids	___ D2216-80		
ASTM Extraction in Water	___ D3987-81/85		
BTU	___ D240-87		
CEC		___ 9081	___ c
Chromium VI		___ 3060A/7196A	
Corrosivity ___ by coupon ___ by pH		___ 1110(mod) ___ 9045C	
Cyanide, Total		___ 9010B	___ ILMO4.0 (e)
Cyanide, Reactive		___ Section 7.3	
Halides, Extractable Organic		___ 9020B	___ EPA 600/4/84-008
Halides, Total		___ 9020B	___ EPA 600/4/84-008
EP Toxicity		___ 1310A	
Flash Point		___ 1010	
Ignitability		___ 1010	
Oil & Grease		___ 9071A	
Carbon, Total Organic		___ 9060	___ Lloyd Kahn (mod)
Oxygen Bomb Prep for Anions	___ D240-87(mod)	___ 5050	
Petroleum Hydrocarbons, Total Recoverable		___ 9071	___ EPA 418.1
pH, Soil		___ 9045C	
Sulfide, Reactive		___ Section 7.3	
Sulfide		___ 9030B(mod)	
Specific Gravity	___ D1429-76C/	___ D5057-90	
Sulfur, Total		___ 9056	
Synthetic Preparation Leach		___ 1312	
Paint Filter		___ 9095A	

Other: *dehydration*

Method: *Aspirator Method: F33615-84-D-4400/0016*

Other:

Method

Recra LabNet Philadelphia
METHOD REFERENCES AND DATA QUALIFIERS

DATA QUALIFIERS

- U = Indicates that the parameter was not detected at or above the reported limit. The associated numerical value is the sample detection limit.
- * = Indicates that the original sample result is greater than 4x the spike amount added.

ABBREVIATIONS

- MB = Method or Preparation Blank.
MS = Matrix Spike.
MSD = Matrix Spike Duplicate.
REP = Sample Replicate
LC = Laboratory Control Sample.
NC = Not calculated.

A suffix of -R, -S, or -T following these codes indicate a replicate, spike or sample duplicate analysis respectively.

ANALYTICAL WET CHEMISTRY METHODS

1. ASTM Standard Methods.
2. USEPA Methods for Chemical Analysis of Water and Wastes (USEPA 600/4-79-020).
3. Test Methods for Evaluating Solid Waste (USEPA SW-846).
 - a. Standard Methods for the Examination of Water and Waste, 16 ed, (1983).
 - b. Standard Methods for the Examination of Water and Waste, 17 ed, (1989)/18ed (1992).
 - c. Method of Soil Analysis, Part 1, Physical and Mineralogical Methods, 2nd ed, (1986).
 - d. Method of Soil Analysis, Part 2, Chemical and Microbiological Properties, Am. Soc. Agron., Madison, WI (1965).
 - e. USEPA Contract Laboratory Program, Statement of Work for Inorganic Analysis.
 - f. Code of Federal Regulations.

Recra LabNet - Lionville

INORGANICS DATA SUMMARY REPORT 12/28/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
-----	-----	-----	-----	-----	-----	-----
-001	BOW9P9	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.0	u MG/KG	1.0	1.0
-002	BOW9R0	% Solids	88.6	%	0.01	1.0
		Hydrazine	3.0	MG/KG	1.0	1.0
-003	BOW9R1	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.9	MG/KG	1.0	1.0
-004	BOW9R2	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.9	MG/KG	1.0	1.0
-005	BOW9R3	% Solids	88.6	%	0.01	1.0
		Hydrazine	2.8	MG/KG	1.0	1.0
-006	BOW9R4	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.8	MG/KG	1.0	1.0
-007	BOW9R5	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.2	MG/KG	1.0	1.0

Recra LabNet - Lionville

INORGANICS METHOD BLANK DATA SUMMARY PAGE 12/28/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10945-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
*****	*****	*****	*****	*****	*****	*****
BLANK10	99LHZ002-MB1	Hydrazine	1.0	u MG/KG	1.0	1.0

Recra LabNet - Lionville

INORGANICS ACCURACY REPORT 12/28/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	INITIAL RESULT	SPIKED AMOUNT	%RECOV	DILUTION FACTOR(SPK)
BLANK10	99LHZ002-MB1	Hydrazine	5.0	1.0 u	5.0	99.1	1.0
		Hydrazine MSD	5.0	1.0 u	5.0	100.5	1.0

Recra LabNet - Lionville

INORGANICS DUPLICATE SPIKE REPORT 12/28/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKE#1 %RECOV	SPIKE#2 %RECOV	%DIFF
BLANK10	99LHZ002-MB1	Hydrazine	99.1	100.5	1.4

Recra LabNet - Lionville

INORGANICS PRECISION REPORT 12/28/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	INITIAL RESULT	REPLICATE	RPD	DILUTION FACTOR (REP)
-----	-----	-----	-----	-----	-----	-----
-007REP	BOW9RS	% Solids	88.6	89.5	1.0	1.0

Recra LabNet - Lionville Laboratory
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNU-HANFORD B99-078

DATE RECEIVED: 09/03/99

RFW LOT # :9909L101

CLIENT ID /ANALYSIS	RFW #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
BOW9P9						
% SOLIDS	001	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	001	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R0						
% SOLIDS	002	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	002	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R1						
% SOLIDS	003	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	003	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R2						
% SOLIDS	004	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	004	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R3						
% SOLIDS	005	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	005	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R4						
% SOLIDS	006	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	006	S	99LHZ002	09/01/99	10/13/99	10/13/99
BOW9R5						
% SOLIDS	007	S	99L&S124	09/01/99	09/17/99	09/20/99
% SOLIDS	007 REP	S	99L&S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	007	S	99LHZ002	09/01/99	10/13/99	10/13/99

LAB QC:

HYDRAZINE	MB1	S	99LHZ002	N/A	10/13/99	10/13/99
-----------	-----	---	----------	-----	----------	----------

Recra LabNet - Lionville Laboratory
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNU-HANFORD B99-078

DATE RECEIVED: 09/03/99

RFW LOT # :9909L101

CLIENT ID /ANALYSIS	RFW #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
HYDRAZINE	MB1 BS	S	99LHZ002	N/A	10/13/99	10/13/99
HYDRAZINE	MB1 BSD	S	99LHZ002	N/A	10/13/99	10/13/99



99096 101

[illegible]

Appendix A

H0536 Hydrazine MDL Study: Raw Data and Analytical Method Information

19 October 1999

Joan Kessner
Bechtel Hanford, Inc.
2355 Stevens Drive
Building 1162
Richland, WA 99352

Reference: Hydrazine MDL Study
SAF# B99-078, SDG# H0536

Dear Ms. Kessner:


Attached is a summary of results and raw data for the referenced project. Also included is the final hydrazine method used for the MDL study and notes on the minor method modifications needed to overcome matrix problems observed in preliminary tests with the site soils.

It appears that the modified U.S. Air Force method will measure hydrazine concentrations in your soils down to approximately 2 ppm (calculated MDL = 2.3 mg/kg).

Please call me with questions and/or comments.

Sincerely,

Recra Environmental, Inc.


for Carter P. Nulton, Ph.D.
Vice President

RESULTS SUMMARY

1999 Wet Chemistry MDLs

PARAMETER	Methods	Date	Analyt	spike ppm	1	2	3	4	5	6	7	units	Ave	SD	MDL WATER (mg/L)	SOIL FACTOR	MDL SOIL (mg/Kg)
Hydrazine (color) soil ***	Air force method	10/13/99	RDB	2.50	(0.00)	2.95	1.92	1.88	2.81	1.81	1.17	mg/Kg	2.09	0.672	NA	1.0	2.3
*** For Hydrazine, sample #1 appears to have not been spiked, and has been excluded from the statistics as an outlier. The MDL is calculated from the other 6 values, using the student t-factor for 6 measurements.																	

RAW DATA

WESTON Analytical Inorganics Section

Date of Prep: 10/13/99
 Date of Analysis: 10/13/99
 Worksheet: HY2N1013
 Computer #: DAWIN
 Directory: INORGANICS
 Run Batch: 99LH2002
 Method: AIRFORCE
 Analyst: ROB
 Instrument: LAMDA B
 CALIB DATA
 SLOPE: 0.44
 INTERCEPT: 0.63
 CORRELATION COEFF.: 0.9977
 DETECTION LIMIT
 0.10 UNITS

RFV	SAMPLE ID	TEST	INST.	ABS	INITIAL	DILUTION	C/D	RESULT	FINAL	DETECTION	REPLICATE	ORIG	REP	SPIKE	SPIKE	LEVEL	AMOUNT	REDOV	LCS	ANT	LCS	?	REDOV	INITIAL	SAMPLE	FINAL	VOL.	SOLIDS	TIME	PREP
99LH2002-H81	HY2NS		0.244	0.496	1	10.0000		4.956	1.000	MS/MS	0.000	0.000		0.5	5.000	99.1								99LH2002	2	20	100.00			
99LH2002-H81	HY2N		0.017	-0.022	1	10.0000		-0.219	0	1.000	MS/MS	0.000												99LH2002	2	20	100.00			
MOL-1	HY2NS		0.017	-0.022	1	11.2455		-0.246	0	1.125	MS/MS	0.000		0.25	2.811	-0.7								99LH2002	2.000	20	88.57			
MOL-2	HY2NS		0.156	0.295	1	16.9937		3.243		1.079	MS/MS	0.000		0.25	2.748	118.0								99LH2002	2.054	20	88.57			
MOL-3	HY2NS		0.111	0.192	1	10.9776		2.113		1.078	MS/MS	0.000		0.25	2.744	77.0								99LH2002	2.057	20	88.57			
MOL-4	HY2NS		0.107	0.180	1	10.7837		2.026		1.078	MS/MS	0.000		0.25	2.696	75.2								99LH2002	2.094	20	88.57			
MOL-5	HY2NS		0.150	0.281	1	11.0000		3.117		1.100	MS/MS	0.000		0.25	2.770	112.5								99LH2002	2.038	20	88.57			
MOL-6	HY2NS		0.106	0.181	1	10.8719		1.968		1.067	MS/MS	0.000		0.25	2.710	72.4								99LH2002	2.077	20	88.57			
MOL-7	HY2NS		0.078	0.117	1	11.0529		1.296		1.105	MS/MS	0.000		0.25	2.763	46.9								99LH2002	2.043	20	88.57			
99LH2002-H81	HY2NT		0.247	0.502	1	10.0000		5.025		1.000	MS/MS	0.000	1.4	0.5	5.000	100.5								99LH2002	2	20	100.00			
CC-H8	HY2N		0.028	0.063	1	10.0000		0.032	0	1.000	MS/MS	0.000												99LH2002	2	20	100.00			
CC-LC BLK	HY2N		0.005	-0.049	1	9.0063		-0.407	0	0.907	MS/MS	0.000												99LH2002	2.023	20	100.00			

ANALYSIS OF HYDRAZINE IN SOILS

1. Reagents (All chemical are ACS reagent grade or better)

- 1.1 Concentrated sulfuric acid
- 1.2 Sulfuric acid, 0.1 N
- 1.3 Hydrazine sulfate
- 1.4 Acetic acid, glacial
- 1.5 Water, distilled or deionized
- 1.6 Hydrazine reagent, p-dimethylaminobenzaldehyde (Hydraver II may be purchased from Hach Chemical Co., P.O. Box 907, Ames, IA 50010; Catalog No. 1790 or prepare a 2.5% solution of p-dimethylaminobenzaldehyde in methanol.

2. Procedure

- 2.1 On soils, weigh out 2.0 g and add 20 mL 0.1 N H_2SO_4 . Tumble for 30 min. For water samples start here. Filter samples through 0.45 μm Millipore® and take 5 mL for analysis. Place in 25 mL volumetric flask.
- 2.2 Add 1 mL of Hydraver II reagent or hydrazine reagent.
- 2.3 Swirl the contents of the flasks intermittently for 8 min.
- 2.4 Adjust flask volume to 25 mL with the glacial acetic acid reagent.
- 2.5 Place stoppers in flask and invert bottle 5-6 times.
- 2.6 Allow to sit for 4 min.
- 2.7 **Filter samples through 0.45 μm Millipore® using 458 nm on the spectrophotometer.**
- 2.8 Read ABS against blank using 458 nm on the spectrophotometer.
- 2.9 Calculations, as on UDMH in soils.

3. Calibrations and Standards

- 3.1 Weigh out 0.4060 g of hydrazine sulfate ($\text{N}_2\text{H}_4 \cdot \text{H}_2\text{SO}_4$). Dissolve in 500

mL of 0.1 N H_2SO_4 in a 1 L volumetric flask. Bring to volume with 0.1 N H_2SO_4 . The resulting solution is 100 mg/mL in hydrazine. Prepare calibration curves using appropriate μL pipettes and the following procedure:

- 3.1.1 Make 10 ppm hydrazine sulfate daily (10 mL - 100 mg/mL to 100 mL with 0.1 N H_2SO_4).
- 3.1.2 Pipet 10 mL of 0.1 N sulfuric acid into each of six volumetric flasks.
- 3.1.3 Carefully pipet 0.05, 0.10, 0.20, 0.40 and 0.50 mL (.1, .2, .5, .8, and 1.0 ppm respectively) of the standard hydrazine solution into the flasks. Process one flask as a blank.
- 3.1.4 Add 1.0 mL of Hydraver II or hydrazine reagent to each flask.
- 3.1.5 Set time for 8 min. and swirl each flask intermittently.
- 3.1.6 After the 8-min. reaction period, bring each flask to 25 mL total volume with glacial acetic acid.
- 3.1.7 Place stoppers in flasks and invert 5-6 times.
- 3.1.8 Set timer for 4 min. to allow bubbles to disappear (tap flasks lightly).
- 3.1.9 Adjust spectrophotometer (458 nm) to 100%T using the blank solution (prepared in the same manner as the unknown samples).
- 3.1.10 Read sample or absorbance within 3 hours of Hydraver II addition.
- 3.1.11 Construct a calibration curve by plotting absorbance against total μg of hydrazine in solution.

4. References

See UDMH in Soils Analysis. USAFSAM Report TR-82-29, Field Sampling and Analysis of Hydrazine and UDMH Vapors in Air: The Firebrick Method, USAF School of Aerospace Medicine, Brooks AFB TX 78235-5501. Method modified for soil analysis by Tom Thomas, USAFOEHL, Brooks AFB TX 78235-5501.

NOTES ON MATRIX EFFECTS AND METHOD MODIFICATIONS

Preliminary tests of the method as written and performed by Recra on past projects indicated two matrix specific problems. The first being that the leachate became very turbid after the addition of glacial acetic acid (Step 2.4 in the method). We found that the turbidity could be removed by filtering through a 4.5 μm Millipore® filter (Step 2.7).

A second issue, which became apparent after the turbidity problem had been solved, was that color formation was not observed with a 0.2 ppm spike (a spike level which can be easily observed in clean water and soil samples). Apparently some co-extracted matrix constituents were interfering with color formation. After trying several dilutions and looking at different leaching solutions, we settled on a 10-fold dilution (2 g leached with 20 mL) for the MDL study.

FAX COVER SHEET**CONFIDENTIALITY NOTICE**

This facsimile transmission is intended only for the use of the individual or entity to which it is addressed, and may contain confidential information belonging to the sender. If you are not the intended recipient, you are hereby notified that any disclosure, copying, distribution, or the taking of any action in reliance on the contents of this information is strictly prohibited. If you have received this transmission in error, please immediately notify us by telephone to arrange for the return of these documents.

**RECRA
LabNet***a division of Recra Environmental, Inc.*

208 Welsh Pool Road
Lionville, PA 19341-1333
Phone: (610) 280-3000
Fax: (610) 280-3041

DATE:

6/6/00

TIME:**NO. PAGES:**

14

(including cover sheet)

**TO:**

Joan Kessner

OF:**FAX #:****FROM:**

Orletta

COPY TO:**SPECIAL INSTRUCTIONS:**

- | | |
|---------------------------------------|---|
| <input type="checkbox"/> Confidential | <input type="checkbox"/> Please Reply / Call |
| <input type="checkbox"/> Urgent | <input type="checkbox"/> For Your Information |
| <input type="checkbox"/> Other: | |

Send Original? ☐ No

- | |
|--|
| <input type="checkbox"/> Regular Mail / Fed Ex / UPS |
| <input type="checkbox"/> Overnight Express |

MESSAGE:

#0536 - originally mailed 1/21/2000
please let me know if this
is sufficient to be mailed to
you today



Chemical and Environmental Measurement Information

21 January 2000

Mr. Kevin Johnson
Thermo NUtech
2030 Wright Avenue
Richmond, California 94804

**Subject: TNU-Hanford Contract N501118
Analytical Data Package**

Dear Mr. Johnson:

Enclosed are the hard copy analytical reports (original plus 1 copy) for the referenced project listed above and the batch number/fraction indicated (marked X) in the ensuing table:

Recra Batch #	9909L101
SDG #	H0536
SAF #	B99-078
Date Received	09-03-99
# Samples	7
Matrix	Soil
Volatiles	
Semivolatiles	
Pest/PCB	
DRO	
GC Scan	
Metals	
Inorganics	X

The electronic data deliverable (EDD) will be sent separately by way of e-mail to Rich Weiss at Bechtel/Hanford. If you have any questions, please don't hesitate to contact me at (610) 280-3012.

Very truly yours,
RECRA LabNet Philadelphia



Orlette S. Johnson
Project Manager

Enclosure

c:\eroup\pm\orlette\tnu-hanford\data\bc ltr.doc

**RECRA
LabNet***a division of Recra Environmental, Inc.**Virtual Laboratories Everywhere***Recra LabNet Philadelphia
Analytical Report****Client : TNU-HANFORD B99-078**
RFW# : 9909L101
SDG# : H0536
SAF# : B99-078**W.O. # : 10985-001-001-9999-00**
Date Received: 09-03-99**INORGANIC CASE NARRATIVE**

1. This narrative covers the analyses of 7 soil samples. A method detection limit (MDL) determination was performed for Hydrazine. The information pertaining to the MDL study is provided within the data package.
2. The samples were prepared and analyzed in accordance with the methods checked on the attached glossary.
3. Sample holding times as required by the method and/or contract were met.
4. The method blank for Hydrazine was within method criteria.
5. The Laboratory Control Samples (LCS) for Hydrazine were within the laboratory control limits. The duplicate LCS was within the 20% Relative Percent Difference (RPD) control limit.
6. The replicate analysis for Percent Solids was within the 20% RPD control limit.
7. Results for the MDL study were reported on a wet weight basis.



J. Michael Taylor
Vice President
Philadelphia Analytical Laboratory1-20-00
Date

nlp/09-101

The results presented in this report relate only to the analytical testing and conditions of the samples at receipt and during storage. All pages of this report are integral parts of the analytical data. Therefore, this report should only be reproduced in its entirety of 12 pages.

Recra LabNet Philadelphia

WET CHEMISTRY

METHODS GLOSSARY FOR SOIL/SOLIDS SAMPLE ANALYSIS

	<u>ASTM</u>	<u>SW846</u>	<u>OTHER</u>
% Ash	___ D2216-80		
% Moisture	___ D2216-80		___ ILMO4.0 (e)
% Solids	✓ D2216-80		___ ILMO4.0 (e)
% Volatile Solids	___ D2216-80		
ASTM Extraction in Water	___ D3987-81/85		
BTU	___ D240-87		
CEC		___ 9081	___ c
Chromium VI		___ 3060A/7196A	
Corrosivity ___ by coupon ___ by pH		___ 1110(mod) ___ 9045C	
Cyanide, Total		___ 9010B	___ ILMO4.0 (e)
Cyanide, Reactive		___ Section 7.3	
Halides, Extractable Organic		___ 9020B	___ EPA 600/4/84-008
Halides, Total		___ 9020B	___ EPA 600/4/84-008
EP Toxicity		___ 1310A	
Flash Point		___ 1010	
Ignitability		___ 1010	
Oil & Grease		___ 9071A	
Carbon, Total Organic		___ 9060	___ Lloyd Kahn (mod)
Oxygne Bomb Prep for Anions	___ D240-87(mod)	___ 5050	
Petroleum Hydrocarbons, Total Recoverable		___ 9071	___ EPA 418.1
pH, Soil		___ 9045C	
Sulfide, Reactive		___ Section 7.3	
Sulfide		___ 9030B(mod)	
Specific Gravity	___ D1429-76C/	___ D5057-90	
Sulfur, Total		___ 9056	
Synthetic Prpearation Leach		___ 1312	
Paint Filter		___ 9095A	
Other: <i>dehydration</i>	Method: <i>As per Method</i>	Method: <i>F33615-84-D-4400/0010</i>	
Other:	Method		

Recra LabNet Philadelphia
METHOD REFERENCES AND DATA QUALIFIERS

DATA QUALIFIERS

- U = Indicates that the parameter was not detected at or above the reported limit. The associated numerical value is the sample detection limit.
- * = Indicates that the original sample result is greater than 4x the spike amount added.

ABBREVIATIONS

- MB = Method or Preparation Blank.
MS = Matrix Spike.
MSD = Matrix Spike Duplicate.
REP = Sample Replicate
LC = Laboratory Control Sample.
NC = Not calculated.

A suffix of -R, -S, or -T following these codes indicate a replicate, spike or sample duplicate analysis respectively.

ANALYTICAL WET CHEMISTRY METHODS

1. ASTM Standard Methods.
2. USEPA Methods for Chemical Analysis of Water and Wastes (USEPA 600/4-79-020).
3. Test Methods for Evaluating Solid Waste (USEPA SW-846).
 - a. Standard Methods for the Examination of Water and Waste, 16 ed, (1983).
 - b. Standard Methods for the Examination of Water and Waste, 17 ed, (1989)/18ed (1992).
 - c. Method of Soil Analysis, Part 1, Physical and Mineralogical Methods, 2nd ed, (1986).
 - d. Method of Soil Analysis, Part 2, Chemical and Microbiological Properties, Am. Soc. Agron., Madison, WI (1965).
 - e. USEPA Contract Laboratory Program, Statement of Work for Inorganic Analysis.
 - f. Code of Federal Regulations.

L-WI-034/D-6/99

1999 Wet Chemistry MDLs

PARAMETER	Method	Date	Analyst	spike ppm	1	2	3	4	5	6	7	units	Ave	SD	MDL WATER (mg/L)	SOIL FACTOR	MDL SOIL (mg/Kg)
zinc (collo) soil ***	Air force method	10/13/99	RDB	2.50	(0.00)	2.95	1.92	1.88	2.81	1.81	1.17	mg/Kg	2.09	0.672	NA	1.0	2.3
Hydrex, sample #1 appears to have not been spiked, and has been excluded from the analysis by an outlier. The MDL is calculated from the other 6 values, using the standard t-factor for 6 measurements.																	

Recra LabNet - Lionville

INORGANICS DATA SUMMARY REPORT 12/28/99

CLIENT: TNU-HANFORD 899-076

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
-001	B0W9P9	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.0	u MG/KG	1.0	1.0
-002	B0W9R0	% Solids	88.6	%	0.01	1.0
		Hydrazine	3.0	MG/KG	1.0	1.0
-003	B0W9R1	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.9	MG/KG	1.0	1.0
-004	B0W9R2	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.9	MG/KG	1.0	1.0
-005	B0W9R3	% Solids	88.6	%	0.01	1.0
		Hydrazine	2.8	MG/KG	1.0	1.0
-006	B0W9R4	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.8	MG/KG	1.0	1.0
-007	B0W9R5	% Solids	88.6	%	0.01	1.0
		Hydrazine	1.2	MG/KG	1.0	1.0

Jun-06-00 13:11

From-RECRA ENVIRONMENTAL INC

6102803041

T-076 P.08/14 F-561

Recra LabNet - Livertville

INORGANICS METHOD BLANK DATA SUMMARY PAGE 12/28/99

CLIENT: TNU-HANFORD W99-078

WORK ORDER: 10985-001-001-9999-00

RECRA LOT #: 99091101

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
BLANK10	99LHZ002-NB1	Hydrazine	1.0 u	MG/KG	1.0	1.0

Recra LabNet - Lionville

INORGANICS ACCURACY REPORT 12/28/99

CLIENT: TNU-HANFORD 899-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	INITIAL RESULT	SPIKED AMOUNT	%RECOV	DILUTION FACTOR (SPK)
BLANK10	99LHZ002-MB1	Hydrazine	5.0	1.0 u	5.0	99.1	1.0
		Hydrazine MSD	5.0	1.0 u	5.0	100.5	1.0

Recra LabNet - Lionville

INORGANICS DUPLICATE SPIKE REPORT 12/28/99

CLIENT: TNU-BANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKE#1 SPIKE#2		%DIFF
			%RECOV	%RECOV	
BLANK10	991H2002-MB1	Hydrazine	99.1	100.5	1.4

Recra LabNet - Lionville

INORGANICS PRECISION REPORT 12/22/99

CLIENT: TNU-HANFORD B99-078

RECRA LOT #: 9909L101

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	INITIAL RESULT	REPLICATE	RPD	DILUTION FACTOR (REP)
-007RES	BOWERS	* Solids	88.6	89.6	1.0	1.0

Recra LabNet - Lionville Laboratory
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNU-HANFORD B99-078

DATE RECEIVED: 09/03/99

RFW LOT # :9909L101

CLIENT ID /ANALYSIS	RFW #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
B0W9P9						
* SOLIDS	001	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	001	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R0						
* SOLIDS	002	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	002	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R1						
* SOLIDS	003	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	003	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R2						
* SOLIDS	004	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	004	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R3						
* SOLIDS	005	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	005	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R4						
* SOLIDS	006	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	006	S	99LHZ002	09/01/99	10/13/99	10/13/99
B0W9R5						
* SOLIDS	007	S	99L*S124	09/01/99	09/17/99	09/20/99
* SOLIDS	007 REP	S	99L*S124	09/01/99	09/17/99	09/20/99
HYDRAZINE	007	S	99LHZ002	09/01/99	10/13/99	10/13/99
LAB QC:						
HYDRAZINE	MB1	S	99LHZ002	N/A	10/13/99	10/13/99

Recra LabNet - Lionville Laboratory
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNU-HANFORD B99-078

DATE RECEIVED: 09/03/99

RFW LOT # :9909L101

CLIENT ID /ANALYSIS	RFW #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
HYDRAZINE	MB1 BS	S	99LHZ002	N/A	10/13/99	10/13/99
HYDRAZINE	MB1 BSD	S	99LHZ002	N/A	10/13/99	10/13/99

⑧ review

RECRA
LabNet

010

RECRA LabNet Use Only
99091101

Custody Transfer Record/Lab Work Request

Client <u>TNU-Hanford</u> <u>1599-078</u>				Refrigerator # <u>6</u>			
Est. Final Proj. Sampling Date				#Type Container			
Project # <u>10485-001-001-9999-00</u>				Liquid			
Project Contact/Phone #				Solid			
RECRA Project Manager <u>AS</u>				Volume			
QC <u>APCC</u> Del <u>AS</u> TAT <u>30 days</u>				Liquid			
Date Rec'd <u>9/15/99</u> Date Due <u>10/15/99</u>				Solid			
Account # <u>9/15/99</u> <u>10/13/99</u>				Preservatives			
ANALYSES REQUESTED →				ORGANIC			
				VOA BNA Pest/PCB Herb			
				INORG			
				Metal CN			

MATRIX CODES: S - Soil SE - Sediment SO - Solid SL - Sludge W - Water O - Oil A - Air DS - Drum Solids DL - Drum Liquids L - EP/TCLP Leachate WI - Wipe X - Other F - Fish	Lab ID	Client ID/Description	Matrix QC Chosen (✓)		Matrix	Date Collected	Time Collected	RECRA LabNet Use Only	
			MS	MSD					
	001	Baw 9P9			S	9/1/99	0736		✓
	2	R0							
	3	R1							
	4	R2							
	5	R3							
	6	R4							
	7	R5							

FIELD PERSONNEL: COMPLETE ONLY SHADED AREAS				DATE/REVISIONS:				RECRA LabNet Use Only			
Special Instructions: <u>Saf# 1599-078</u>				1. <u>Run matrix QC</u>				Samples were:			
				2. <u>9/13/99 Date rec'd - 9/13/99</u>				1) Shipped <u> </u> or			
				3. <u>Hydrazine MDL Study</u>				Hand Delivered <u> </u>			
				4. <u> </u>				Airfill # <u> </u>			
				5. <u>Samples relogged from</u>				2) Ambient or Chilled			
				6. <u>batch 990911007-99PM229</u>				3) Received in Good Condition Y or N			
								4) Labels Indicate Properly Preserved Y or N			
								5) Received Within Holding Times Y or N			

Relinquished by	Received by	Date	Time	Relinquished by	Received by	Date	Time
<u>Relog</u>	<u>Janson</u>	<u>9/15/99</u>	<u>-</u>				

Discrepancies Between Samples Labels and COC Record? Y or N	NOTES:

COC Tape was:
1) Present on Outer Package Y or N
2) Unbroken on Outer Package Y or N
3) Present on Sample Y or N
4) Unbroken on Sample Y or N
COC Record Present Upon Sample Rec't Y or N