

0054376

H/031

Recra LabNet - Lionville Laboratory
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNUHANFORD B00-068 H1031

DATE RECEIVED: 10/21/00

RFW LOT # :0010L008

CLIENT ID /ANALYSIS	RFW #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
---------------------	-------	-----	--------	------------	-----------	----------

BOYW00

TCLP	001	S	00LTO123	08/28/00	10/25/00	10/26/00
CHROMIUM, TCLP LEACH	002	W	99L1655	10/26/00	10/27/00	10/29/00
CHROMIUM, TCLP LEACH	002 REP	W	99L1655	10/26/00	10/27/00	10/29/00
CHROMIUM, TCLP LEACH	002 MS	W	99L1655	10/26/00	10/27/00	10/29/00

LAB QC:

CHROMIUM LABORATORY	LC1 BS	W	99L1655	N/A	10/27/00	10/29/00
CHROMIUM, TCLP LEACH	MB1	W	99L1655	N/A	10/27/00	10/29/00
CHROMIUM, TCLP LEACH	MB2	W	99L1655	N/A	10/27/00	10/29/00

RECEIVED
JAN 22 2001

EDMC

001



**Recra LabNet Philadelphia
Analytical Report**

Client: TNU-HANFORD B00-068
RFW#: 0010L008
SDG/SAF#: H1031/B00-068

W.O.#: 10985-001-001-9999-00
Date Received: 10-21-00


METALS CASE NARRATIVE

1. This narrative covers the analyses of 1 TCLP leachate sample.
2. The sample was prepared and analyzed in accordance with methods checked on the attached glossary. This is a re-log of 0009L519-001.
3. All analyses were performed within the required holding times.
4. The cooler temperature has been recorded on the original Chain of Custody.
5. All Initial and Continuing Calibration Verifications (ICV/CCVs) were within the 90-110% control limits (80-120% for Mercury).
6. All Initial and Continuing Calibration Blanks (ICB/CCBs) were within control limits (less than the PQL).
7. All preparation/method blanks (MB) were within method criteria {less than the Practical Quantitation Limit (3X the IDL), MB value less than 5% of the RCRA limit, or samples greater than 20X MB value}. Refer to the Inorganics Method Blank Data Summary.
8. All ICP Interference Check Standards were within control limits.
9. The laboratory control sample (LCS) was within the 80-120% control limits. Refer to form 7.
10. The duplicate analysis was outside the 20% Relative Percent Difference (RPD) control limits. Refer to the Inorganics Precision Report.
11. The TCLP extract from sample B0YW00 was selected for the matrix spike (MS) for this analytical batch. The MS recovery was greater than 50% as per method criteria.

The results presented in this report relate only to the analytical testing and conditions of the samples at receipt and during storage. All pages of this report are integral parts of the analytical data. Therefore, this report should only be reproduced in its entirety of 11 pages.

002

12. For the purposes of this report, the data has been reported to the Instrument Detection Limit (IDL). Values between the IDL and the Practical Quantitation Limit (PQL) are acquired in a region of less-certain quantification.
13. I certify that this sample data package is in compliance with SOW requirements, both technically and for completeness, other than the conditions detailed above. Release of the data contained in this hard-copy data package has been authorized by the Laboratory Manager or a designee, as verified by the following signature.



J. Michael Taylor
VP, Laboratory General Manager
Lionville Laboratory
gmb/m10-008

12-6-00
Date



METALS METHOD GLOSSARY

The following methods are used as reference for the digestion and analysis of samples contained within this

Recra Lot#:

00102008

Leaching Procedure: ☐ 1310 ☒ 1311 ☐ 1312 ☐ Other: _____

CLP Metals ☐ Digestion and ☐ Analysis Methods: ☐ ILM03.0 ☐ ILM04.0

Metals Digestion Methods: ☐ 3005A ☒ 3010A ☐ 3015 ☐ 3020A ☐ 3050B ☐ 3051 ☐ 200.7 ☐ SS17
☐ Other: _____

Metals Analysis Methods

	SW846	EPA	STD MTD	EPA OSWR	USATHAMA
Aluminum	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Antimony	<input type="checkbox"/> 6010B <input type="checkbox"/> 7041 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 204.2			<input type="checkbox"/> 99
Arsenic	<input type="checkbox"/> 6010B <input type="checkbox"/> 7060A ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 206.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Barium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Beryllium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Bismuth	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Boron	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Cadmium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7131A ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 213.2			<input type="checkbox"/> 99
Calcium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Chromium	<input checked="" type="checkbox"/> 6010B <input type="checkbox"/> 7191 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 218.2			<input type="checkbox"/> SS17
Cobalt	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Copper	<input type="checkbox"/> 6010B <input type="checkbox"/> 7211 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 220.2			<input type="checkbox"/> 99
Iron	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Lead	<input type="checkbox"/> 6010B <input type="checkbox"/> 7421 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 239.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Lithium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7430 ⁴	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Magnesium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Manganese	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Mercury	<input type="checkbox"/> 7470A ³ <input type="checkbox"/> 7471A ³	<input type="checkbox"/> 245.1 ² <input type="checkbox"/> 245.5 ²			<input type="checkbox"/> 99
Molybdenum	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Nickel	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Potassium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7610 ⁴	<input type="checkbox"/> 200.7 <input type="checkbox"/> 258.1 ⁴			<input type="checkbox"/> 99
Rare Earths	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Selenium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7740 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 270.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Silicon	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Silica	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Silver	<input type="checkbox"/> 6010B <input type="checkbox"/> 7761 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 272.2			<input type="checkbox"/> 99
Sodium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7770 ⁴	<input type="checkbox"/> 200.7 <input type="checkbox"/> 273.1 ⁴			<input type="checkbox"/> 99
Strontium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Thallium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7841 ^s	<input type="checkbox"/> 200.7 <input type="checkbox"/> 279.2 <input type="checkbox"/> 200.9			<input type="checkbox"/> 99
Tin	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Titanium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Uranium	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Vanadium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Zinc	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Zirconium	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99

Other: _____

Method: _____

METHOD REFERENCES AND DATA QUALIFIERS

DATA QUALIFIERS

U = Indicates that the parameter was not detected at or above the reported limit. The associated numerical value is the sample detection limit.

* = Indicates that the original sample result is greater than 4x the spike amount added.

ABBREVIATIONS

MB = Method or Preparation Blank.

MS = Matrix Spike.

MSD = Matrix Spike Duplicate.

REP = Sample Replicate

LCS = Laboratory Control Sample.

NC = Not calculated.

ANALYTICAL METAL METHODS

1. Not included in the method element list.
2. Modified Hg: Hg1 and Hg2 require less total volume of digestate due to the autosampler analysis. Sample volumes and reagents for mercury determinations in water and soil have been proportionately scaled down to adapt to this semi-automated technique. The sample volume used for water analysis is 33 mL. For soils, 0.1 grams of sample is taken to a final volume of 50 mL (including all reagents).
3. Modified Hg: Hg1 and Hg2 require less total volume of digestate due to the autosampler analysis. Sample volumes and reagents for mercury determinations in water and soil have been proportionately scaled down to adapt to this semi-automated technique. The sample volume used for water analysis is 33 mL. For soils, three 0.1 gram of sample is taken to a final volume of 50 mL (including all reagents).
4. Flame AA.
5. Graphite Furnace AA.

RFW 21-21L-033/N-10/96

Recra LabNet - Lionville

INORGANICS DATA SUMMARY REPORT 12/06/00

CLIENT: TNUHANFORD B00-068 H1031
WORK ORDER: 10985-001-001-9999-00

RECRA LOT #: 0010L008

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
-002	BOYW00	Chromium, TCLP Leachate	5.0	UG/L	3.5	1.0

Recra LabNet - Lionville

INORGANICS METHOD BLANK DATA SUMMARY PAGE 12/06/00

CLIENT: TNUHANFORD B00-068 H1031
WORK ORDER: 10985-001-001-9999-00

RECRA LOT #: 0010L008

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
*****	*****	*****	*****	*****	*****	*****
BLANK1	99L1655-MB1	Chromium, TCLP Leachate	3.5	u UG/L	3.5	1.0
BLANK2	99L1655-MB2	Chromium, TCLP Leachate	3.9	UG/L	3.5	1.0

Recra LabNet - Lionville

INORGANICS ACCURACY REPORT 12/06/00

CLIENT: TNUHANFORD B00-068 H1031

RECRA LOT #: 0010L008

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	INITIAL RESULT	SPIKED AMOUNT	%RECOV	DILUTION FACTOR (SPK)
-002	BOYW00	Chromium, TCLP Leachat	4390	5.0	5000	87.8	1.0

Recra LabNet - Lionville

INORGANICS PRECISION REPORT 12/06/00

CLIENT: TNUHANFORD B00-068 H1031

RECRA LOT #: 0010L008

WORK ORDER: 10985-001-001-9999-00

SAMPLE	SITE ID	ANALYTE	INITIAL RESULT	REPLICATE	RPD	DILUTION FACTOR (REP)
*****	*****	*****	*****	*****	*****	*****
-002REP	BOYW00	Chromium, TCLP Leachate	5.0	8.1	47.3	1.0

Recra LabNet - Lionville

INORGANICS LABORATORY CONTROL STANDARDS REPORT 12/06/00

CLIENT: TNUHANFORD B00-068 H1031
WORK ORDER: 10985-001-001-9999-00

RECRA LOT #: 0010L008

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	SPIKED AMOUNT	UNITS	%RECOV
-----	-----	-----	-----	-----	-----	-----
LCS1	99L1655-LC1	Chromium, LCS	497	500	UG/L	99.4

0010L008

Custody Transfer Record/Lab Work Request Page 1 of 1

Page 1 of 1

FIELD PERSONNEL: COMPLETE ONLY SHADED AREAS[illegible]