

H1227

0055514

Eberline Services  
W.O. No. R1-01-147-7620

Bechtel Hanford Inc.  
SDG H1227

Case Narrative

Page 1 of 1

## 1.0 GENERAL

Bechtel Hanford Inc. (BHI) Sample Delivery Group H1227 was composed of one water sample designated under SAF No. C01-015 with a Project Designation of: 200 UP1 IAM GW Monitoring, January 2001.

The sample was received as stated on the Chain-of-Custody document. Any discrepancies are noted on the Eberline Services Sample Receipt Checklist.

## 2.0 ANALYSIS NOTES

### 2.1 Technetium-99 Analyses

No problems were encountered during the course of the analyses.

#### Case Narrative Certification Statement

"I certify that this data package is in compliance with the SOW, both technically and for completeness, for other than the conditions detailed above. Release of the data obtained in this hard copy data package has been authorized by the Laboratory Manager or a designee, as verified by the following signature."

Melissa C. Mannion  
Melissa C. Mannion  
Program Manager

3/9/01  
Date

RECEIVED  
OCT 23 2001  
EDMC

TMA/RICHMOND  
SAMPLE DELIVERY GROUP H1227

SDG 7620  
Contact Melissa C. Mannion

SAMPLE SUMMARY

Client Hanford  
Contract TRC-SBB-207925  
Case no SDG H1227

CLIENT SAMPLE ID	LOCATION	MATRIX	LEVEL	LAB SAMPLE ID	SAF NO	CHAIN OF CUSTODY	COLLECTED
B113Y2	Hanford Site	WATER		R101147-01	C01-015	C01-015-6	01/23/01 12:43
Method Blank		WATER		R101147-03	C01-015		
Lab Control Sample		WATER		R101147-02	C01-015		
Duplicate (R101147-01)	Hanford Site	WATER		R101147-04	C01-015		01/23/01 12:43

**TMA/RICHMOND**  
**SAMPLE DELIVERY GROUP H1227**

SDG 7620  
 Contact Melissa C. Mannion

**QC SUMMARY**

Client Hanford  
 Contract TRC-SBB-207925  
 Case no SDG H1227

QC BATCH	CHAIN OF CUSTODY	CLIENT SAMPLE ID	MATRIX	% SOLIDS	SAMPLE AMOUNT	BASIS AMOUNT	DAYS SINCE RECEIVED	LAB COLL	LAB SAMPLE ID	DEPARTMENT SAMPLE ID
7620	C01-015-6	B113Y2	WATER				01/25/01	2	R101147-01	7620-001
		Method Blank	WATER						R101147-03	7620-003
		Lab Control Sample	WATER						R101147-02	7620-002
		Duplicate (R101147-01)	WATER				01/25/01	2	R101147-04	7620-004

QC SUMMARY

Page 1

SUMMARY DATA SECTION

Page 4

Lab id TMANC  
 Protocol Hanford  
 Version Ver 1.0  
 Form DVD-QS  
 Version 3.06  
 Report date 03/09/01

TMA/RICHMOND  
SAMPLE DELIVERY GROUP H1227

SDG 7620  
Contact Melissa C. Mannion

PREP BATCH SUMMARY

Client Hanford  
Contract TRC-SBB-207925  
Case no SDG H1227

TEST	MATRIX	METHOD	PREPARATION BATCH	ERROR 2 $\sigma$ %	CLIENT MORE	PLANCHETS ANALYZED RE BLANK	LCS	DUP/ORIG	MS/ORIG	QUALI- FIERS
Beta Counting										
TC	WATER	Technetium 99 in Water	6963-108	10.0	1	1	1	1/1		

Duplicates and Matrix Spikes are those with original (Client) sample in this Sample Delivery Group.  
Blank and LCS planchets are those in the same preparation batch as some Client, Duplicate or Spike sample.

**TMA/RICHMOND**  
**SAMPLE DELIVERY GROUP H1227**

**WORK SUMMARY**

SDG 7620  
 Contact Melissa C. Mannion

Client Hanford  
 Contract TRC-SBB-207925  
 Case no SDG H1227

CLIENT SAMPLE ID	LAB SAMPLE ID									
LOCATION	MATRIX	COLLECTED		SUF-						
CUSTODY	SAF No	RECEIVED	PLANCHET	TEST	FIX	ANALYZED	REVIEWED	BY	METHOD	
B113Y2		R101147-01	7620-001	TC		02/08/01	03/09/01	MCM	Technetium 99 in Water	
Hanford Site	WATER	01/23/01								
C01-015-6	C01-015	01/25/01								
Method Blank		R101147-03	7620-003	TC		02/06/01	03/09/01	MCM	Technetium 99 in Water	
	WATER									
	C01-015									
Lab Control Sample		R101147-02	7620-002	TC		02/05/01	03/09/01	MCM	Technetium 99 in Water	
	WATER									
	C01-015									
Duplicate (R101147-01)		R101147-04	7620-004	TC		02/06/01	03/09/01	MCM	Technetium 99 in Water	
Hanford Site	WATER	01/23/01								
	C01-015	01/25/01								

COUNTS OF TESTS BY SAMPLE TYPE										
TEST	SAF No	METHOD	REFERENCE	CLIENT	MORE	RE	BLANK	LCS	DUP	SPIKE TOTAL
TC	C01-015	Technetium 99 in Water	TC99_TR_SEP_LSC	1			1	1	1	4
TOTALS				1			1	1	1	4

WORK SUMMARY

Page 1

SUMMARY DATA SECTION

Page 6

Lab id TMANC  
 Protocol Hanford  
 Version Ver 1.0  
 Form DVD-CWS  
 Version 3.06  
 Report date 03/09/01

**R101147-03**

### Method Blank

## METHOD BLANK

SDG <u>7620</u>	Client/Case no <u>Hanford</u>	SDG <u>H1227</u>
Contact <u>Melissa C. Mannion</u>	Contract <u>TRC-SBB-207925</u>	
Lab sample id <u>R101147-03</u>	Client sample id <u>Method Blank</u>	
Dept sample id <u>7620-003</u>	Material/Matrix <u>WATER</u>	
	SAF No <u>C01-015</u>	

ANALYTE	CAS NO	RESULT pCi/L	2σ ERR (COUNT)	MDA pCi/L	RDL pCi/L	QUALI- FIERS	TEST
Technetium 99	14133-76-7	-1.26	4.8	14	15	U	TC

200 UP1 IAM GW Mon., Jan. 2001

QC-BLANK 37425

Lab id TMANC  
Protocol Hanford  
Version Ver 1.0  
Form DVD-DS  
Version 3.06  
Report date 03/09/01

TMA/RICHMOND  
SAMPLE DELIVERY GROUP H1227

R101147-02

Lab Control Sample

LAB CONTROL SAMPLE

SDG <u>7620</u>	Client/Case no <u>Hanford</u>	SDG <u>H1227</u>
Contact <u>Melissa C. Mannion</u>	Case no <u>TRC-SBB-207925</u>	
Lab sample id <u>R101147-02</u>	Client sample id <u>Lab Control Sample</u>	
Dept sample id <u>7620-002</u>	Material/Matrix <u>WATER</u>	
	SAF No <u>C01-015</u>	

ANALYTE	RESULT pCi/L	2 $\sigma$ ERR (COUNT)	MDA pCi/L	RDL pCi/L	QUALI- FIERS	TEST	ADDED pCi/L	2 $\sigma$ ERR pCi/L	REC %	3 $\sigma$ LMTS (TOTAL)	PROTOCOL LIMITS
Technetium 99	2440	56	13	15		TC	2360	94	103	83-117	80-120

200 UP1 IAM GW Mon., Jan. 2001

QC-LCS 37424
--------------

Lab id <u>TMANC</u>
Protocol <u>Hanford</u>
Version <u>Ver 1.0</u>
Form <u>DVD-LCS</u>
Version <u>3.06</u>
Report date <u>03/09/01</u>

**TMA/RICHMOND**  
**SAMPLE DELIVERY GROUP H1227**

R101147-04

B113Y2

**DUPLICATE**

SDG <u>7620</u>		Client/Case no <u>Hanford</u> SDG <u>H1227</u>	
Contact <u>Melissa C. Mannion</u>		Case no <u>TRC-SBB-207925</u>	
<b>DUPLICATE</b>		<b>ORIGINAL</b>	
Lab sample id <u>R101147-04</u>	Lab sample id <u>R101147-01</u>	Client sample id <u>B113Y2</u>	
Dept sample id <u>7620-004</u>	Dept sample id <u>7620-001</u>	Location/Matrix <u>Hanford Site</u> <u>WATER</u>	
	Received <u>01/25/01</u>	Collected <u>01/23/01 12:43</u>	
		Custody/SAF No <u>C01-015-6</u> <u>C01-015</u>	

ANALYTE	DUPLICATE pCi/L	2σ ERR (COUNT)	MDA pCi/L	RDL pCi/L	QUALI- FIERS	TEST	ORIGINAL pCi/L	2σ ERR (COUNT)	MDA pCi/L	QUALI- FIERS	RPD %	3σ TOT	PROT LIMIT
Technetium 99	138	18	11	15		TC	138	10	12		0	31	

200 UP1 IAM GW Mon., Jan. 2001

QC-DUP#1 37426



TMA / RICHMOND  
SAMPLE DELIVERY GROUP H1227

R101147-01

B113Y2

DATA SHEET

SDG <u>7620</u>	Client/Case no <u>Hanford</u>	SDG <u>H1227</u>
Contact <u>Melissa C. Mannion</u>	Contract <u>TRC-SBB-207925</u>	
Lab sample id <u>R101147-01</u>	Client sample id <u>B113Y2</u>	
Dept sample id <u>7620-001</u>	Location/Matrix <u>Hanford Site</u>	<u>WATER</u>
Received <u>01/25/01</u>	Collected <u>01/23/01 12:43</u>	
	Custody/SAF No <u>C01-015-6</u>	<u>C01-015</u>

ANALYTE	CAS NO	RESULT pCi/L	2σ ERR (COUNT)	MDA pCi/L	RDL pCi/L	QUALI- FIERS	TEST
Technetium 99	14133-76-7	138	10	12	15		TC

200 UP1 IAM GW Mon., Jan. 2001

DATA SHEETS

Page 1

SUMMARY DATA SECTION

Page 10

Lab id <u>TMANC</u>
Protocol <u>Hanford</u>
Version <u>Ver 1.0</u>
Form <u>DVD-DS</u>
Version <u>3.06</u>
Report date <u>03/09/01</u>

**TMA/RICHMOND**  
SAMPLE DELIVERY GROUP H1227

Test TC Matrix WATER  
SDG 7620  
Contact Melissa C. Mannion

**METHOD SUMMARY**  
TECHNETIUM 99 IN WATER  
BETA COUNTING

Client Hanford  
Contract TRC-SBB-207925  
Contract SDG H1227

**RESULTS**

CLIENT SAMPLE ID	LAB SAMPLE ID	RAW TEST	SUF- FIX	Technetium PLANCHET	99
Preparation batch 6963-108					
B113Y2	R101147-01			7620-001	138
BLK (QC ID=37425)	R101147-03			7620-003	U
LCS (QC ID=37424)	R101147-02			7620-002	ok
Duplicate (R101147-01)	R101147-04			7620-004	ok

Nominal values and limits from method RDLs (pCi/L) 15  
200 UP1 IAM GW Mon., Jan. 2001

**METHOD PERFORMANCE**

CLIENT SAMPLE ID	LAB SAMPLE ID	RAW TEST	SUF- FIX	MDA pCi/L	ALIQ L	PREP FAC	DILU- TION	YIELD %	EFF %	COUNT min	FWHM keV	DRIFT KeV	DAYS HELD	ANAL- PREPARED	YZED	DETECTOR
Preparation batch 6963-108 2σ prep error 10.0 % Reference Lab Notebook 6963 pg. 108																
B113Y2	R101147-01			12	0.0500			83		50			16	02/01/01	02/08	GRB-222
BLK (QC ID=37425)	R101147-03			14	0.0500			71		50				02/01/01	02/06	GRB-219
LCS (QC ID=37424)	R101147-02			13	0.0500			78		50				02/01/01	02/05	GRB-230
Duplicate (R101147-01) (QC ID=37426)	R101147-04			11	0.0500			87		50			14	02/01/01	02/06	GRB-220
Nominal values and limits from method				15	0.0500			20-105		50			180			

PROCEDURES REFERENCE TC99\_TR\_SEP\_LSC  
CP-021 Preparation of Tc-99m Tracer, rev 0  
CP-002 Q.C. Preparation, rev 2  
CP-003 Tracing, rev 2  
CP-541 Technetium-99 Purification (Water) by Extraction  
Chromatography, rev 0  
CP-008 Heavy Element Electroplating, rev 3

AVERAGES ± 2 SD MDA 12 ± 2.6  
FOR 4 SAMPLES YIELD 80 ± 14

**METHOD SUMMARIES**

Page 1

**SUMMARY DATA SECTION**

Page 11

Lab id TMANC  
Protocol Hanford  
Version Ver 1.0  
Form DVD-CMS  
Version 3.06  
Report date 03/09/01

PNNL		CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST						C.O.C. # C01-015-6  Page 1 of 1	
Collector D.R. BREWINGTON			Contact/Requester JH KESSNER			Telephone No. MSIN FAX (509) 375-4688			
SAF No. C01-015			Sampling Origin HANFORD SITE H1227 (7620)			Purchase Order/Charge Code SAL-79			
Project Title 200 UPLI AM GW MONITORING JANUARY 2001			Logbook No. Wm-SAWS-H41			Ice Chest No. SAC-25 Temp.			
Shipped To (Lab) TMA/RECRA			Method of Shipment GOVT VEHICLE			Bill of Lading/Air Bill No. 4235-7954-1838			
Protocol CERCLA			Data Turnaround 45 Days			Offsite Property No.			
POSSIBLE SAMPLE HAZARDS/REMARKS .. ..						SPECIAL INSTRUCTIONS Hold Time Fax TMA log-in to JH Kessner (372-9487) & DL Stewart (372-1704) Submit deliverable & invoices to JH Kessner, BHI			
						Total Activity Exemption: Yes [X] No [ ]			
Sample No.	Lab ID	*	Date	Time	No/Type Container	Sample Analysis		Preservative	
B113Y2 /		W	1/23/01	1243	1x20-mL P /	Activity Scan		None	
B113Y2 /		W	L	L	1x1000-mL G/P /	Technetium-99 /		HCl to pH <2	
Relinquished By D.R. BREWINGTON	Sign <i>[Signature]</i>	Date/Time 1/24/01	Received By FED EX	Print E. Leguro	Sign <i>[Signature]</i>	Date/Time 1-24-01	Matrix *		
Relinquished By FED EXPRESS		Date/Time 1-25-01	Received By 			Date/Time 1-25-01	S = Soil SE = Sediment SL = Sludge W = Water O = Oil A = Air	DS = Drum Solid DL = Drum Liqui T = Tissue WI = Wipe L = Liquid V = Vegetation X = Other	
Relinquished By		Date/Time	Received By			Date/Time			
Relinquished By		Date/Time	Received By			Date/Time			
FINAL SAMPLE DISPOSITION		Disposal Method (e.g., Return to customer, per lab procedure, used in process)				Disposed By		Date/Time	

ICE CHEST RECEIPT LOG

Use one form per shipment. Refer to Thermometer Correction Log for correction factor.

Customer: (WMTS) DECHTEL HANFORD (PNNL) Date: 1-25-01 10:AM

Ice chest # or description	SML-79					
Thermometer: time in	1:00 PM					
Thermometer: time out						
Thermometer reading						
Thermometer number	2132					
Correction factor	NONE					
Actual temperature*	—					
Custody seals on ice chest intact?	Yes					
Custody seals dated?	Yes					
Custody seals signed?	Yes					
Custody seals on samples?	Yes					
Ice chest scanned for activity?	Yes					

\* Temperature is in degrees centigrade.

Technician: \_\_\_\_\_

Comments: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

SAMPLE RECEIPT CHECKLIST

SAMPLE RECEIPT			
Client:	(WMTS) BECHTEL HANFORD (PNNL)	Date/Time received	1-25-01 10:00 AM
CoC No.	COI-015-6		
Container I.D. No.	SML-79	Requested TAT (Days)	45 P.O. Received Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>
INSPECTION			
1.	Custody seals on shipping container intact?	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/> N/A <input type="checkbox"/>	
2.	Custody seals on shipping container dated & signed?	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/> N/A <input type="checkbox"/>	
3.	Custody seals on sample containers intact?	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/> N/A <input type="checkbox"/>	
4.	Custody seals on sample containers dated & signed?	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/> N/A <input type="checkbox"/>	
5.	Cooler Temperature: _____	Packing material is:	Wet <input type="checkbox"/> Dry <input checked="" type="checkbox"/>
6.	Number of samples in shipping container: 1		
7.	Number of containers per sample: (2 EACH) (Or see CoC _____)		
8.	Paperwork agrees with samples?	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	
9.	Samples have: Tape <input type="checkbox"/> Hazard labels <input type="checkbox"/> Rad labels <input type="checkbox"/> Appropriate sample labels <input checked="" type="checkbox"/>		
10.	Samples are: In good condition <input checked="" type="checkbox"/> Leaking <input type="checkbox"/> Broken Container <input type="checkbox"/> Missing <input type="checkbox"/>		
11.	Describe any anomalies: _____		
13.	Was P.M. notified of any anomalies? Yes <input type="checkbox"/> No <input type="checkbox"/> Date _____		
14.	Received by E. Segura Date: 1-25-01 Time: 10:00 AM		

Customer Sample No.	cpm	mr/hr	Customer Sample No.	Cpm	mr/hr

Ion Chamber Ser. No. \_\_\_\_\_ Calibration date \_\_\_\_\_

Survey Meter Ser. No. \_\_\_\_\_ Calibration date \_\_\_\_\_