

W03760

0057843

Analytical Data Package Prepared For

Bechtel Hanford

Radiochemical Analysis By

STL Richland*2800 G.W. Way, Richland, Wa 99352, (509) 375-3131*

Assigned Laboratory Code: STLRL

Data Package Contains 30 Pages

Report No.: 19552

SDG No.	Order No.	Client Sample ID (List Order)	Lot-Sa No.	Work Order	Report DB ID	Batch No.
W03760	B02-032	B14JB2	J2D230207-1	E0AM91AC	9E0AM910	2113394
		B14JB2	J2D230207-1	E0AM91AA	9E0AM910	2113395

RECEIVED
AUG 1 2 2002
EDMC



CERTIFICATE OF ANALYSIS

Bechtel Hanford, Inc.
3350 George Washington Way
Richland, WA 99352

May 13, 2002

Attention: Joan Kessner

SAF Number	:	B02-032
Date SDG Closed	:	April 23, 2002
Number of Samples	:	One
Sample Type	:	Air Filter (Solid)
SDG Number	:	W03760
Data Deliverable	:	21-Day / Summary

I. Introduction

On April 23, 2002, one air filter (solid) sample was received at STL Richland (STLR) for radiochemical analysis. Upon receipt, the sample was assigned the following laboratory ID number to correspond with the Bechtel Hanford, Inc. (BHI) specific ID:

<u>STLR ID#</u>	<u>BHI ID#</u>	<u>MATRIX</u>	<u>DATE OF RECEIPT</u>
E0AM9	B14JB2	SOLID	4/23/02

II. Analytical Results/Methodology

The analytical results for this report are presented by laboratory sample ID. Each set of data includes sample identification information, analytical results and the appropriate associated statistical errors.

The requested analyses were:

- Gamma Spectroscopy**
Gamma Scan by method RICH-RC-5017
- Gas Proportional Counting**
Total Strontium by method RICH-RC-5006

III. Quality Control

The analytical results for each analysis performed under SDG W03760 include a minimum of one Laboratory Control Sample (LCS), one method (reagent) blank, and one duplicate sample analysis. Any exceptions have been noted in the "Comments" section.

Bechtel Hanford, Inc.
May 13, 2002
Page 2

QC and sample results are reported in the same units.

IV. Comments

Gamma Spectroscopy

Gamma Scan by method RICH-RC-5017:

There was insufficient sample material to meet the CRDLs. Except as noted, the LCS, batch blank, samples and sample duplicate (B14JB2) results are within contractual requirements.


Gas Proportional Counting

Total Strontium by method RICH-RC-5006:

There was insufficient sample material to process a batch sample duplicate. Except as noted, the LCS, batch blank, and sample results are within contractual requirements.

I certify that this Certificate of Analysis is in compliance with the SOW, both technically and for completeness, for other than the conditions detailed above. Release of the data contained in this hard copy data package has been authorized by the Laboratory Manager, or a designee as verified by the following signature.

Reviewed and approved:



Barbara M. Gillespie
Project Manager

0003

Drinking Water Method Cross References

DRINKING WATER ASTM METHOD CROSS REFERENCES		
Referenced Method	Isotope(s)	STL Richland's SOP number
EPA 901.1	Cs-134, I-131	RICH-RC-5017
EPA 900.0	Alpha & Beta	RICH-RC-5014
EPA 903.1	Ra-226	RICH-RC-5005
EPA 904.0	Ra-228	RICH-RC-5005
EPA 905.0	Sr89/90	RICH-RC-5006
ASTM D2460	Total Radium	RICH-RC-5027
Standard Method 7500-U-C & ASTM D5174	Uranium	RICH-RC-5058
EPA 906.0	Tritium	RICH-RC-5007
NOTE:		
The Gross Alpha LCS is prepared with Am-241 (unless otherwise specified in the case narrative)		
The Gross Beta LCS is prepared with Sr/Y-90 (unless otherwise specified in the case narrative)		

Uncertainty Estimation

STL Richland has adopted the internationally accepted approach to estimating uncertainties described in "NIST Technical Note 1297, 1994 Edition". The approach, "Law of Propagation of Errors", involves the identification of all variables in an analytical method which are used to derive a result. These variables are related to the analytical result (R) by some functional relationship, $R = \text{constants} * f(x, y, z, \dots)$. The components (x, y, z) are evaluated to determine their contribution to the overall method uncertainty. The individual component uncertainties (u_i) are then combined using a statistical model that provides the most probable overall uncertainty value. All component uncertainties are categorized as type A, evaluated by statistical methods, or type B, evaluated by other means. Uncertainties not included in the components, such as sample homogeneity, are combined with the component uncertainty as the square root of the sum-of-the-squares of the individual uncertainties. The uncertainty associated with the derived result is the combined uncertainty (u_c) multiplied by the coverage factor (1, 2, or 3).

When three or more sample replicates are used to derive the analytical result, the type A uncertainty is the standard deviation of the mean value (S/\sqrt{n}), where S is the standard deviation of the derived results. The type B uncertainties are all other random or non-random components that are not included in the standard deviation.

The derivation of the general "Law of Propagation of Errors" equations and specific example are available on request.

Report Definitions

Action Lev	An agreed upon activity level used to trigger some action when the final result is greater than or equal to the Action Level. Often the Action Level is related to the Decision Limit.
Batch	The QC preparation batch number that relates laboratory samples to QC samples that were prepared and analyzed together.
Bias	Defined by the equation $(\text{Result}/\text{Expected}) - 1$ as defined by ANSI N13.30.
COC No	Chain of Custody Number assigned by the Client or STL Richland.
Count Error (#s)	Poisson counting statistics of the gross sample count and background. The uncertainty is absolute and in the same units as the result. For Liquid Scintillation Counting (LSC) the batch blank count is the background.
Total Uncert (#s) u_c - Combined Uncertainty.	All known uncertainties associated with the preparation and analysis of the sample are propagated to give a measure of the uncertainty associated with the result, u_c the combined uncertainty. The uncertainty is absolute and in the same units as the result.
(#s), Coverage Factor	The coverage factor defines the width of the confidence interval, 1, 2 or 3 standard deviations.
CRDL (RL)	Contractual Required Detection Limit as defined in the Client's Statement Of Work or STL Richland "default" nominal detection limit. Often referred to the reporting level (RL)
Lc	Decision Level based on instrument background or blank, adjusted by the Efficiency, Chemical Yield, and Volume associated with the sample. The Type I error probability is approximately 5%. $L_c = (1.645 * \text{Sqrt}(2 * (\text{BkgrndCnt}/\text{BkgrndCntMin}) / \text{SCntMin})) * (\text{ConvFct}/(\text{Eff} * \text{Yld} * \text{Abn} * \text{Vol}) * \text{IngrFct})$. For LSC methods the batch blank is used as a measure of the background variability. Lc cannot be calculated when the background count is zero.
Lot-Sample No	The number assigned by the LIMS software to track samples received on the same day for a given client. The sample number is a sequential number assigned to each sample in the Lot.
MDC MDA	Detection Level based on instrument background or blank, adjusted by the Efficiency, Chemical Yield, and Volume with a Type I and II error probability of approximately 5%. $MDC = (4.65 * \text{Sqrt}((\text{BkgrndCnt}/\text{BkgrndCntMin}) / \text{SCntMin}) + 2.71 / \text{SCntMin}) * (\text{ConvFct}/(\text{Eff} * \text{Yld} * \text{Abn} * \text{Vol}) * \text{IngrFct})$. For LSC methods the batch blank is used as a measure of the background variability.
Primary Detector	The instrument identifier associated with the analysis of the sample aliquot.
Ratio U-234/U-238	The U-234 result divided by the U-238 result. The U-234/U-238 ratio for natural uranium in NIST SRM 4321C is 1.038.
Rst/MDC	Ratio of the Result to the MDC. A value greater than 1 may indicate activity above background at a high level of confidence. Caution should be used when applying this factor and it should be used in concert with the qualifiers associated with the result.
Rst/TotUcert	Ratio of the Result to the Total Uncertainty. If the uncertainty has a coverage factor of 2 a value greater than 1 may indicate activity above background at approximately the 95% level of confidence assuming a two-sided confidence interval. Caution should be used when applying this factor and it should be used in concert with the qualifiers associated with the result.
Report DB No	Sample Identifier used by the report system. The number is based upon the first five digits of the Work Order Number.
RER	The equation Replicate Error Ratio = $(S - D) / [\text{sqrt}(\text{TPUs}^2 + \text{TPUd}^2)]$ as defined by ICPT BOA where S is the original sample result, D is the result of the duplicate, TPUs is the total uncertainty of the original sample and TPUd is the total uncertainty of the duplicate sample.
SDG	Sample Delivery Group Number assigned by the Client or assigned by STL Richland upon sample receipt.
Sum Rpt Alpha Spec Rst(s)	The sum of the reported alpha spec results for tests derived from the same sample excluding duplicate result where the results are in the same units.
Work Order	The LIMS software assign test specific identifier.
Yield	The recovery of the tracer added to the sample such as Pu-242 used to trace a Pu-239/40 method.

Sample Results Summary

Date: 14-May-02

STL Richland STLRL

Ordered by Client Sample ID, Batch No.

Report No. : 19552

SDG No: W03760

Client ID	Work Order Number	Parameter	Result +- Uncertainty (2s)	Qual	Units	Yield	MDC MDA	RER
B14JB2	E0AM91AC	STRONTIUM	-7.46E-02 +- 2.5E-01	U	pCi/sample	86.90%	5.67E-01	
B14JB2	E0AM91AA	CO-60	1.42E+00 +- 2.3E+00	U	pCi/sample		4.58E+00	
		CS-137	3.67E-01 +- 1.9E+00	U	pCi/sample		3.45E+00	
		EU-152	1.12E+00 +- 4.7E+00	U	pCi/sample		8.46E+00	
		EU-154	1.81E+00 +- 5.3E+00	U	pCi/sample		1.06E+01	
		EU-155	1.74E+00 +- 4.4E+00	U	pCi/sample		8.12E+00	
B14JB2 DUP	E0AM91AE	CO-60	1.44E+00 +- 1.8E+00	U	pCi/sample		3.87E+00	
		CS-137	6.79E-02 +- 1.7E+00	U	pCi/sample		3.11E+00	
		EU-152	-9.37E-01 +- 3.9E+00	U	pCi/sample		6.74E+00	
		EU-154	3.82E-01 +- 5.4E+00	U	pCi/sample		1.03E+01	
		EU-155	6.43E-01 +- 2.7E+00	U	pCi/sample		4.71E+00	

Number of Results: 11

QC Results Summary
STL Richland STLRL
 Ordered by QC Type, Batch No.

Date: 14-May-02

Report No. : 19552

SDG No.: W03760

QC Type	Work Order Number	Parameter	Result +- Uncertainty (2s)	Qual	Units	Yield	Recovery	Bias	MDC MDA
BLANK QC	E0AW41AA	STRONTIUM	-2.44E-01 +- 2.3E-01	U	pCi/sample	88.60%			5.40E-01
BLANK QC	E0AW61AA	CO-60	-6.18E-01 +- 2.0E+00	U	pCi/sample				3.56E+00
		CS-137	5.62E-01 +- 1.6E+00	U	pCi/sample				3.10E+00
		EU-152	2.90E-01 +- 4.2E+00	U	pCi/sample				7.41E+00
		EU-154	5.44E-03 +- 6.1E+00	U	pCi/sample				1.13E+01
		EU-155	-5.81E-01 +- 2.8E+00	U	pCi/sample				4.82E+00
LCS	E0AW41AC	STRONTIUM	6.02E+00 +- 1.8E+00		pCi/sample	83.90%	87.33%	-0.1	6.03E-01
LCS	E0AW61AC	CO-60	7.56E+01 +- 1.2E+01		pCi/sample		97.63%	0.0	4.95E+00
		CS-137	4.65E+01 +- 9.5E+00		pCi/sample		92.79%	-0.1	4.64E+00
		EU-152	1.61E+02 +- 2.3E+01		pCi/sample		104.79%	0.0	1.01E+01

Number of Results: 10

FORM I

SAMPLE RESULTS

Date: 14-May-02

Lab Name: STL Richland

SDG: W03760

Collection Date: 4/23/2002 10:30:00 AM

Lot-Sample No.: J2D230207-1

Report No. : 19552

Received Date: 4/23/2002 10:45:00 AM

Client Sample ID: B14JB2

COC No. :

Matrix: SOLID

Ordered by Client Sample ID, Batch No.

Parameter	Result	Qual	Count Error (2 s)	Total Uncert(2 s)	MDC MDA, Action Lev	Rpt Unit, Lc	Yield CRDL(RL)	Rst/MDC, Rst/TotUncert	Analysis, Prep Date	Total Sa Size	Aliquot Size	Analy Method, Primary Detector
Batch: 2113394	Work Order: E0AM91AC			Report DB ID: 9E0AM910								
STRONTIUM	-7.46E-02	U	2.5E-01	2.5E-01	5.67E-01	pCi/sample 2.70E-01	86.90%	-0.13 -0.6	5/10/02 07:06 p	1.0 Sample	1.0 Sample	SRISO_SEP_PRECIP GPC28A
Batch: 2113395	Work Order: E0AM91AA			Report DB ID: 9E0AM910								
CO-60	1.42E+00	U	2.3E+00	2.3E+00	4.58E+00	pCi/sample		0.31 (1.3)	4/30/02 05:40 a		1.0 SA	GAMMA_GS GER3\$1
CS-137	3.67E-01	U	1.9E+00	1.9E+00	3.45E+00	pCi/sample		0.11 0.39	4/30/02 05:40 a		1.0 SA	GAMMA_GS GER3\$1
EU-152	1.12E+00	U	4.7E+00	4.7E+00	8.46E+00	pCi/sample		0.13 0.47	4/30/02 05:40 a		1.0 SA	GAMMA_GS GER3\$1
EU-154	1.81E+00	U	5.3E+00	5.3E+00	1.06E+01	pCi/sample		0.17 0.69	4/30/02 05:40 a		1.0 SA	GAMMA_GS GER3\$1
EU-155	1.74E+00	U	4.4E+00	4.4E+00	8.12E+00	pCi/sample		0.21 0.78	4/30/02 05:40 a		1.0 SA	GAMMA_GS GER3\$1

Number of Results: 6

Comments:

0008

FORM II

Date: 14-May-02

DUPLICATE RESULTS

Lab Name: STL Richland

SDG: W03760

Collection Date: 4/23/2002 10:30:00 AM

Lot-Sample No.: J2D230207-1

Report No.: 19552

Received Date: 4/23/2002 10:45:00 AM

Client Sample ID: B14JB2 DUP

COC No.:

Matrix: SOLID

Parameter	Result, Orig Rst	Qual	Count Error (2 s)	Total Uncert(2 s)	MDC MDA, Action Lev	Rpt Unit, CRDL	Yield	Rst/MDC, Rst/TotUncert	Analysis, Prep Date	Total Sa Size	Aliquot Size	Analy Method, Primary Detector
Batch: 2113395	Work Order: E0AM91AE			Report DB ID: E0AM91ER		Orig Sa DB ID: 9E0AM910						
CO-60	1.44E+00	U	1.8E+00	1.8E+00	3.87E+00	pCi/sample		0.37	5/1/02 05:43 a		1.0	GAMMA_GS
	1.42E+00	RER	0.0					(1.6)			SA	GER8\$1
CS-137	6.79E-02	U	1.7E+00	1.7E+00	3.11E+00	pCi/sample		0.02	5/1/02 05:43 a		1.0	GAMMA_GS
	3.67E-01	RER	0.2					0.08			SA	GER8\$1
EU-152	-9.37E-01	U	3.9E+00	3.9E+00	6.74E+00	pCi/sample		-0.14	5/1/02 05:43 a		1.0	GAMMA_GS
	1.12E+00	RER	0.7					-0.48			SA	GER8\$1
EU-154	3.82E-01	U	5.4E+00	5.4E+00	1.03E+01	pCi/sample		0.04	5/1/02 05:43 a		1.0	GAMMA_GS
	1.81E+00	RER	0.4					0.14			SA	GER8\$1
EU-155	6.43E-01	U	2.7E+00	2.7E+00	4.71E+00	pCi/sample		0.14	5/1/02 05:43 a		1.0	GAMMA_GS
	1.74E+00	RER	0.4					0.48			SA	GER8\$1

Number of Results: 5

Comments:

6009

FORM II BLANK RESULTS

Date: 14-May-02

Lab Name: STL Richland

SDG: W03760

Lot-Sample No.: J2D230000-394

Report No. : 19552

Matrix: SOLID

Parameter	Result	Qual	Count Error (2 s)	Total Uncert(2 s)	MDC MD A,	Rpt Unit, CRDL	Yield	Rst/MDC, Rst/TotUcert	Analysis, Prep Date	Total Sa Size	Allquot Size	Analy Method, Primary Detector
Batch: 2113394	Work Order: E0AW41AA			Report DB ID: E0AW41AB								
STRONTIUM	-2.44E-01	U	2.2E-01	2.3E-01	5.40E-01	pCi/sample	88.60%	-0.45	5/10/02 07:06 p	1.0	1.0	SRISO_SEP_PRECIP
					2.56E-01			-(2.1)		Sample	Sample	GPC28B

Number of Results: 1

Comments:

0010

FORM II BLANK RESULTS

Date: 14-May-02

Lab Name: STL Richland

SDG: W03760

Lot-Sample No.: J2D230000-395

Report No. : 19552

Matrix: SOLID

Parameter	Result	Qual	Count Error (2 s)	Total Uncert(2 s)	MDC MD A,	Rpt Unit, CRDL	Yield	Rst/MDC, Rst/TotUncert	Analysis, Prep Date	Total Sa Size	Aliquot Size	Analy Method, Primary Detector
Batch: 2113395	Work Order: E0AW61AA				Report DB ID: E0AW61AB							
CO-60	-6.18E-01	U	2.0E+00	2.0E+00	3.56E+00	pCi/sample		-0.17	4/30/02 05:40 a		1.0	GAMMA_GS
								-0.62			SA	GER8\$1
CS-137	5.62E-01	U	1.6E+00	1.6E+00	3.10E+00	pCi/sample		0.18	4/30/02 05:40 a		1.0	GAMMA_GS
								0.69			SA	GER8\$1
EU-152	2.90E-01	U	4.2E+00	4.2E+00	7.41E+00	pCi/sample		0.04	4/30/02 05:40 a		1.0	GAMMA_GS
								0.14			SA	GER8\$1
EU-154	5.44E-03	U	6.1E+00	6.1E+00	1.13E+01	pCi/sample		0.	4/30/02 05:40 a		1.0	GAMMA_GS
								0.			SA	GER8\$1
EU-155	-5.81E-01	U	2.8E+00	2.8E+00	4.82E+00	pCi/sample		-0.12	4/30/02 05:40 a		1.0	GAMMA_GS
								-0.41			SA	GER8\$1

Number of Results: 5

Comments:

0011

FORM II
LCS RESULTS

Date: 14-May-02

Lab Name: STL Richland

SDG: W03760

Lot-Sample No.: J2D230000-394

Report No. : 19552

Matrix: SOLID

Parameter	Result	Count Qual Error (2 s)	Total Uncert(2 s)	MDC MD	Report Unit	Yield	Expected	Expected Uncert	Recovery, Bias	Analysis, Prep Date	Aliquot Size	Analy Method, Primary Detector
Batch: 2113394	Work Order: E0AW41AC	Report DB ID: E0AW41CS										
STRONTIUM	6.02E+00	5.7E-01	1.8E+00	6.03E-01	pCi/sample	83.90%	6.89E+00	8.5E-02	87.33%	5/10/02 07:06 p	1.0	SRISO_SEP_PRECIP
Rec Limits:									-0.1		Sample	GPC28C

Number of Results: 1

Comments:

0012

FORM II
LCS RESULTS

Date: 14-May-02

Lab Name: STL Richland

SDG: W03760

Lot-Sample No.: J2D230000-395

Report No. : 19552

Matrix: SOLID

Parameter	Result	Count Qual Error (2 s)	Total Uncert(2 s)	MDC MD	Report Unit	Yield	Expected	Expected Uncert	Recovery, Bias	Analysis, Prep Date	Aliquot Size	Analy Method, Primary Detector
Batch: 2113395	Work Order: E0AW61AC				Report DB ID: E0AW61CS							
CO-60	7.56E+01	1.2E+01	1.2E+01	4.95E+00	pCi/sample		7.74E+01	4.9E-01	97.63%	4/30/02 05:41 a	1.0	GAMMA_GS
					Rec Limits:		.	.	0.0		SA	GER7\$1
CS-137	4.65E+01	9.5E+00	9.5E+00	4.64E+00	pCi/sample		5.01E+01	1.7E+00	92.79%	4/30/02 05:41 a	1.0	GAMMA_GS
					Rec Limits:		.	.	-0.1		SA	GER7\$1
EU-152	1.61E+02	2.3E+01	2.3E+01	1.01E+01	pCi/sample		1.54E+02	6.2E+00	104.79%	4/30/02 05:41 a	1.0	GAMMA_GS
					Rec Limits:		.	.	0.0		SA	GER7\$1

Number of Results: 3

Comments:

0013

SEVERN

TRENT

SERVICES

Data Review Checklist
 RADIOCHEMISTRY
 First Level Review

Lot Number: J2D230207

Client ID: BHE

Due Date: 5-14-02

QC Batch Number: 2113395

Method Test Parameter: GAMMA

Matrix: OTHER FILTER

SDG Number: W03760

Review Item	Yes (✓)	No (✓)	N/A (✓)
A. COC			
1. Is the ICOC page complete (includes all applicable analysts, dates, SOP numbers and revisions)?	✓		
B. QC Batch			
1. Do the Summary/Detailed Reports include a calculated result for each sample listed on the QC Batch Sheet?	✓		
2. Are the QC appropriate for the analysis included in the batch?	✓		
3. Is the Analytical Batch Worksheets complete (includes, as appropriate, volumes, count times, etc.)?	✓		
4. Does the Worksheets include a Tracer Vial label for each sample?			✓
C. QC & Samples			
1. Is the blank result, yield and MDA within contract limits?	✓	✓	
2. Is the LCS result, yield and MDA within contract limits?	✓	✓	
3. Are the MS/MSD results, yields and MDAs within contract limits?			✓
4. Are the duplicate results, yields and MDAs within contract limits?	✓	✓	
5. Are the sample yields and MDAs within contract limits?	✓	✓	
D. Raw Data			
1. Were results calculated in the correct units?	✓		
2. Were analysis volumes entered correctly?	✓		
3. Were yields entered correctly?			✓
4. Were spectra reviewed/meet contractual requirements?	✓		
5. Were raw counts reviewed for anomalies?			✓
E. Other			
1. Are all Nonconformances included and noted? J05417	✓		
2. Are all required forms filled out?	✓		
3. Was the correct methodology used?	✓		
4. Was transcription checked?	✓		
5. Were all calculations checked at a minimum frequency?	✓		

Comments on any "No" response:

First Level Review:

Date: 5-2-02

Data Review Checklist
RADIOCHEMISTRY
Second Level Review



QC Batch Number: 211 3395

Review Item	Yes (✓)	No (✓)	N/A (✓)
A. Calibration			
1. Is the calibration documentation included?			✓
B. Sample Analysis			
1. Are the sample yields within acceptance criteria?			✓
2. Is the sample Minimum Detectable Activity < the Contract Detection Limit?		✓	
C. QC Samples			
1. Is the blank yield within acceptance criteria?			✓
2. Is the Minimum Detectable Activity for the blank result ≤ the Contract Detection Limit?		✓	
3. Does the blank result meet the Contract criteria?		✓	
4. Is the blank result < the Contract Detection Limit?		✓	
5. Is the blank result > the Contract Detection Limit but the sample result < the Contract Detection Limit?		✓	
6. Is the LCS result within acceptance criteria?	✓		
7. Is the LCS yield within acceptance criteria?			✓
8. Is the LCS Minimum Detectable Activity ≤ the Contract Detection Limit?		✓	
9. Do the MS/MSD results and yields meet acceptance criteria?			✓
10. Do the duplicate sample results and yields meet acceptance criteria?	✓		
D. Other			
1. Are all Nonconformances included and noted?	✓		
2. Are all required forms filled out?	✓		
3. Was the correct methodology used?	✓		
4. Were units checked?	✓		

Comments on any "No" response: _____

Second Level Review: me mby Date: 5/2/02

Clouseau Nonconformance Memo

SEVERN
TRENT
SERVICES

NCM #: J05417	Classification: Anomaly
NCM Initiated By: Dale OConnell	Status: PMREVIEW
Date Opened: 05/02/02	Production Area: Environmental - Prep
Date Closed: N/A	Tests: Gamma by GER
	Lot #'s (Sample #'s): J2D230000 (395); J2D230207 (1)
	QC Batch: 2113395
Nonconformance: QC Result Out of Limits	
Subcategory: MDA exceeds RDL	

Problem Description / Root Cause

<u>Name</u>	<u>Date</u>	<u>Description</u>
Dale OConnell	05/02/02	Cause is limited sample quantity, 1 sample, and sequential analysis within a short turn time frame.
Dale OConnell	05/02/02	Insufficient sample to generate a duplicate.

Corrective Action

<u>Name</u>	<u>Date</u>	<u>Corrective Action</u>
Dale OConnell	05/02/02	Report results with MDAs achieved.
Dale OConnell	05/02/02	Precision determination achieved by recounting sample on a different detector.

Approval History

<u>Name</u>	<u>Date Approved:</u>	<u>Position</u>
Dale OConnell	05/02/02	

Data Review Checklist
RADIOCHEMISTRY
First Level Review

Lot Number: J2D230207 P
Client ID: BHI
Due Date: 5/14/82
QC Batch Number: 2113394
Method Test Parameter: TH-TSR
Matrix: air
SDG Number: W03760

Review Item	Yes (✓)	No (✓)	N/A (✓)
A. COC			
1. Is the ICOC page complete (includes all applicable analysts, dates, SOP numbers and revisions)?	✓		
B. QC Batch			
1. Do the Summary/Detailed Reports include a calculated result for each sample listed on the QC Batch Sheet?	✓		
2. Are the QC appropriate for the analysis included in the batch?	✓		
3. Is the Analytical Batch Worksheets complete (includes, as appropriate, volumes, count times, etc.)?	✓		
4. Does the Worksheets include a Tracer Vial label for each sample?	✓		
C. QC & Samples			
1. Is the blank result, yield and MDA within contract limits?	✓		
2. Is the LCS result, yield and MDA within contract limits?	✓		
3. Are the MS/MSD results, yields and MDAs within contract limits?			
4. Are the duplicate results, yields and MDAs within contract limits?			✓
5. Are the sample yields and MDAs within contract limits?	✓		
D. Raw Data			
1. Were results calculated in the correct units?	✓		
2. Were analysis volumes entered correctly?	✓		
3. Were yields entered correctly?	✓		
4. Were spectra reviewed/meet contractual requirements?			✓
5. Were raw counts reviewed for anomalies?	✓		
E. Other			
1. Are all Nonconformances included and noted?			✓
2. Are all required forms filled out?	✓		
3. Was the correct methodology used?	✓		
4. Was transcription checked?	✓		
5. Were all calculations checked at a minimum frequency?	✓		

Comments on any "No" response: _____

First Level Review: Pam Anderson Date: 5-13-82

**SEVERN
TRENT
SERVICES**

Review Item	Yes (✓)	No (✓)	N/A (✓)
A. Calibration			✓
1. Is the calibration documentation included?			✓
B. Sample Analysis	✓		X emb
1. Are the sample yields within acceptance criteria?	✓		
2. Is the sample Minimum Detectable Activity < the Contract Detection Limit?	✓		
C. QC Samples	✓		
1. Is the blank yield within acceptance criteria?	✓		
2. Is the Minimum Detectable Activity for the blank result ≤ the Contract Detection Limit?	✓		
3. Does the blank result meet the Contract criteria?	✓		
4. Is the blank result < the Contract Detection Limit?	✓		
5. Is the blank result > the Contract Detection Limit but the sample result < the Contract Detection Limit?			✓
6. Is the LCS result within acceptance criteria?	✓		
7. Is the LCS yield within acceptance criteria?	✓		
8. Is the LCS Minimum Detectable Activity ≤ the Contract Detection Limit?	✓		
9. Do the MS/MSD results and yields meet acceptance criteria?			✓
10. Do the duplicate sample results and yields meet acceptance criteria?			✓
D. Other			
1. Are all Nonconformances included and noted?			✓
2. Are all required forms filled out?	✓		
3. Was the correct methodology used?	✓		
4. Were units checked?	✓		

Second Level Review: B. M. B. Date: 5/13/02

CHAIN OF CUSTODY

U-21038

25423-01

Bechtel Hanford Inc.		CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST				B02-032-2		Page 1 of 2	
Collector R. Thoren		Company Contact I.D. Jacques		Telephone No. 372-9651		Project Coordinator TRENT, SJ		Price Code 9L Data Turnaround 21 Days	
Project Designation B Plant - Air Filter Analysis		Sampling Location B-Plant		SAF No. B02-032		Air Quality <input checked="" type="checkbox"/>			
Ice Chest No. Sml-001		Field Logbook No. EL-1562-1		COA B221BT2W14		Method of Shipment Gov. Vehicle			
Shipped To Severn Trent Incorporated, Richland		Offsite Property No. NA		Bill of Lading/Air Bill No. NA					
POSSIBLE SAMPLE HAZARDS/REMARKS Tie TO B14J82 J83 #B14J84 Special Handling and/or Storage None				Preservation None					
				Type of Container P					
				No. of Container(s) 1					
				Volume 60mL					
SDG W03760 SAMPLE ANALYSIS Due 5-13-02 J2D230207				See item (1) in Special Instructions.					
Sample No.	Matrix *	Sample Date	Sample Time						
B14JB2 E0AM9	OTHER SOLID	4/23/02	1030	X					
CHAIN OF POSSESSION				SPECIAL INSTRUCTIONS				Matrix *	
Relinquished By/Removed From R. Thoren		Date/Time 4/23/02		Received By/Stored In D. Thoren		Date/Time 4/23/02		(1) Gamma Spectroscopy (Cesium-137); Strontium-89,90 - Total Sr, Activity Scan; RCF GEA Shipping Screen	
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time		S=Soil SS=Sediment SO=Solid SI=Sludge W=Water O=Oil A=Air DS=Drum Solid DL=Drum Liquid T=Tissue WI=Wipe L=Liquid V=Vegetation X=Other	
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time		NOTE: This sample - original Sample # B14J82, B14J83, B14J84 sent to RCF. Samples picked up from RCF & Given Sample # B14JB2 on New COC with new label for STL	
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time		page 2 of 2 shows RCF TRANSFER	
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
LABORATORY SECTION		Received By		Title					
FINAL SAMPLE DISPOSITION		Disposal Method		Disposed By				Date/Time	

P.2/3

1200 APR 19 '02 08:46AM BHI S&D MANAGEMENT 509 372 9487

Bechtel Hanford Inc.		CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST				B02-032-01		Page 1 of 1	
Collector R. Nielson / T. Myers		Company Contact Jacques, ID		Telephone No. 372-9651		Project Coordinator TRENT, ST		Price Code 9L	
Project Designation B Plant - Air Filter Analysis		Sampling Location 200 East		SAF No. B02-032		Air Quality <input checked="" type="checkbox"/>		Data Turnaround 21 Days	
Ice Chest No.		Field Logbook No. BL-1562-1		COA B221BT2W14		Method of Shipment Government Vehicle			
Shipped To Seyern Trust Incorporated, Richland		Offsite Property No.		Bill of Lading/Air Bill No.					
POSSIBLE SAMPLE HAZARDS/REMARKS				Preservation		None			
Special Handling and/or Storage				Type of Container		D/P			
				No. of Container(s)		1			
				Volume		GAL			
SAMPLE ANALYSIS				See Para (1) in Special Instructions					
Sample No.		Matrix *		Sample Date		Sample Time			
B14J82		OTHER SOLID		4-8-02		1357		X	
B14J83		OTHER SOLID		4-10-02		1547		X	
B14J84		OTHER SOLID		4-11-02		1326		X	
CHAIN OF POSSESSION				SPECIAL INSTRUCTIONS				Matrix *	
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time		<p>** RCP is to store the samples for additional analysis at commercial laboratory.</p> <p>(1) Gamma Spectroscopy (Cesium-137); Strontium-90,50 -- Total Sr, Activity Scan; RCP QRA Shipping Screen</p> <div style="border: 2px solid black; padding: 5px; text-align: center; font-weight: bold; font-size: 1.2em;">AIR QUALITY DOCUMENT</div>	
Troy Myers		4/9/02 1130		R. Nielson		4/19/02 1130			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
R. Nielson		4/19/02 1300		T. Duffey		4-19-02 1300			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
T. Duffey		4-23-02							
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
Relinquished By/Removed From		Date/Time		Received By/Stored In		Date/Time			
LABORATORY SECTION		Received By		Title				Date/Time	
FINAL SAMPLE DISPOSITION		Disposal Method		Disposed By				Date/Time	

20F2
RF 423.02

NOTE: These samples combined, given new sample # B14J82 for TRANSFER TO STL for analysis. Page 1 of 2 shows custody transfer.

ERC Radiological Counting Facility Analysis Report

RCF Number RCF10203

Sample Date & Time 4/8/2002 1357

Project ID: 200 EAST

SAF Number: B02-032

Date Analyzed 4/22/2002 8:28

Sample ID: B14J82

Gamma Energy Analysis

Nuclide	Activity (pCi/g)	Error (pCi/g)	MDC (pCi/g)
K-40	< 1.0E+02		1.0E+02
Co-60	< 1.2E+01		1.2E+01
Cs-137	< 8.7E+00		8.7E+00
Eu-152	< 2.1E+01		2.1E+01
Eu-154	< 3.3E+01		3.3E+01
Eu-155	< 1.7E+01		1.7E+01
Tl-208	< 2.3E+01		2.3E+01
Pb-212	< 9.4E+01		9.4E+01
Bi-214	< 6.5E+01		6.5E+01
Pb-214	< 1.5E+01		1.5E+01
Ra-226	< 1.4E+02		1.4E+02
Ac-228	< 3.4E+01		3.4E+01
Pa-234	< 1.6E+01		1.6E+01
Th-234	< 6.6E+01		6.6E+01
U-235	< 4.0E+01		4.0E+01
Am-241	< 1.1E+01		1.1E+01

Tie to
B14JBZ

**QUALITATIVE
ONLY**

Total GEA (pCi/g) +/-

	Activity (pCi/g)	Error (pCi/g)
Gross Alpha**	N/R	+/- N/R
Gross Beta	N/R	+/- N/R

Definitions:

All errors reported at 2 standard deviations.

N/R = no result or analysis not requested. <MDA = Less than detection limit.

All GEA results reported as "<" list the Minimum Detectable Concentration (MDC) value for that radionuclide.

Rounding error may result in the reported total GEA activity differing from the sum of the > MDA GEA values in the second significant digit.

For soils and natural samples, the following applies:

The analysis of U-238 is based on the activity of Pa-234m.

The analysis of Np-237 is based on the activity of Pa-233.

U-238da is the activity of Pb-214 and Bi-214, short lived daughter products of U-238. Equilibrium between parent and daughter products probably does not exist in disturbed materials.


Th-232da is the activity of Ac-228, Pb-212, and Tl-208, short lived daughter products of Th-232. Equilibrium between parent and daughter products may not exist in disturbed materials.

Other samples, not containing natural materials, may have inapplicable results for the Th, U, transuranics and daughter products. The results must then be balanced for the gross alpha analysis.

**The gross alpha results are not corrected for mass absorption

No peaks for this radionuclide were visible above background in the spectrum. The result was reported as less than MDC.

Analyst


P.J. Snider

4/22/2002

Report To

Duane Jacques

SJ Trent

Joan Kessner

Fax

373-7711

372-9292

969-4823

Report Printed: Monday, April 22, 2002

0022

ERC Radiological Counting Facility Analysis Report

RCF Number RCF10204Project ID: 290 EAST

SAF Number: B02-032

Sample Date & Time 4/10/2002 1347

Sample ID: B14183

Date Analyzed 4/22/2002 9:31

Gamma Energy Analysis

Nuclide	Activity (pCi/g)	Error (pCi/g)	MDC (pCi/g)
K-40	< 8.8E+01		8.8E+01
Co-60	< 1.1E+01		1.1E+01
Cs-137	< 6.1E+00		6.1E+00
Ba-152	< 2.2E+01		2.2E+01
Eu-154	< 2.7E+01		2.7E+01
Ba-155	< 1.7E+01		1.7E+01
Tl-208	< 2.7E+01		2.7E+01
Pb-212	< 1.1E+01		1.1E+01
Bi-214	< 7.7E+01		7.7E+01
Pb-214	< 1.6E+01		1.6E+01
Ra-226	< 1.4E+02		1.4E+02
Ac-228	< 3.7E+01		3.7E+01
Pa-234	< 1.6E+01		1.6E+01
Th-234	< 5.5E+01		5.5E+01
U-235	< 4.5E+01		4.5E+01
Am-241	< 9.2E+00		9.2E+00

THE TO
B14JBZ

**QUALITATIVE
ONLY**

Total GEA (pCi/g) +/-

Activity (pCi/g)	Error (pCi/g)
Gross Alpha** N/R +/-	N/R
Gross Beta N/R +/-	N/R

Definitions:

All errors reported at 2 standard deviations.

N/R = no result or analysis not requested. <MDDA = Less than detection limit.

All GEA results reported as "<" but the Minimum Detectable Concentration (MDC) value for that radionuclide.

Rounding error may result in the reported total GEA activity differing from the sum of the > MDA GEA values in the second significant digit.

For soils and natural samples, the following applies:

The analysis of U-238 is based on the activity of Pa-234m.

The analysis of Np-237 is based on the activity of Pa-233.

U-238 decay is the activity of Pb-214 and Bi-214, short lived daughter products of U-238. Equilibrium between parent and daughter products probably does not exist in disturbed materials.

Th-232 decay is the activity of Ac-228, Pb-212, and Tl-208, short lived daughter products of Th-232. Equilibrium between parent and daughter products may not exist in disturbed materials.

Other samples, not containing natural materials, may have inapplicable results for the Th, U, transuramics and daughter products. The results must then be balanced for the gross alpha analysis.

**The gross alpha results are not corrected for mass absorption

No peaks for this radionuclide were visible above background in the spectrum. The result was reported as less than 1%DC.

Analyst

T.J. Snider

4/22/2002

Report To

Dunace Jacques

Fax

373-7711

SJ Trant

372-9292

Joan Kestner

909-4823

Report Printed: Monday, April 22, 2002

0023

ERC Radiological Counting Facility Analysis Report

RCF Number RCF10206Project ID: 200 EASTSAF Number: B02-032Sample Date & Time 4/11/2002

1326

Date Analyzed 4/22/2002 10:5Sample ID: B14J84

Gamma Energy Analysis

Nuclide	Activity (pCi/g)	Error (pCi/g)	MDC (pCi/g)
K-40	< 7.8E+01		7.8E+01
Co-60	< 1.0E+01		1.0E+01
Cs-137	< 9.7E+00		9.7E+00
Ba-152	< 2.1E+01		2.1E+01
Ba-154	< 3.3E+01		3.3E+01
Bu-155	< 1.6E+01		1.6E+01
Tl-208	< 2.0E+01		2.0E+01
Pb-212	< 8.3E+01		8.3E+01
Bi-214	< 7.0E+01		7.0E+01
Pb-214	< 1.5E+01		1.5E+01
Ra-226	< 1.5E+02		1.5E+02
Ac-228	< 3.1E+01		3.1E+01
Pa-234	< 1.6E+01		1.6E+01
Th-234	< 6.4E+01		6.4E+01
U-235	< 4.0E+01		4.0E+01
Am-241	< 9.6E+00		9.6E+00

TIC TO
B14JB2

**QUALITATIVE
ONLY**

Total GEA (pCi/g)	+/-
Activity (pCi/g)	Error (pCi/g)
Gross Alpha**	N/R +/- N/R
Gross Beta	N/R +/- N/R

Definitions:

All errors reported at 2 standard deviations.

N/R = no result or analysis not requested. <MDA = Less than detection limit.

All GEA results reported as "<" list the Minimum Detectable Concentration (MDC) value for that radionuclide.

Rounding error may result in the reported total GEA activity differing from the sum of the > MDA GEA values in the second significant digit.

For soils and natural samples, the following applies:

The analysis of U-238 is based on the activity of Pa-234m.

The analysis of U-235 is based on the activity of Th-231.

U-238 and U-235 are the activity of Pa-214 and Bi-214, short lived daughter products of U-238. Equilibrium between parent and daughter products probably does not exist in disturbed materials.

Th-232 and U-235 are the activity of Ac-228, Pa-212, and Tl-208, short lived daughter products of Th-232. Equilibrium between parent and daughter products may not exist in disturbed materials.

Other samples, not containing natural materials, may have inappropriate results for the Th, U, transuramics and daughter products. The results must then be balanced for the gross alpha analysis.

**The gross alpha results are not corrected for mass absorption.

No peaks for this radionuclide were visible above background in the spectrum. The result was reported as less than MDC.

Analyst



4/22/2002

Report Printed: Monday, April 22, 2002

Report To

Dunne Jacques

SU Trant

Joan Keaner

For

373-7711

372-9292

969-4823

0024



Sample Check-in List

- Date/Time Received: 4-23 1045
- Client: BHI SDG #: W03760 NA ☐ SAF #: B02-032 NA ☐
- Work Order Number: J2D230207 Chain of Custody # B02-032-2
- Shipping Container ID: cha Air Bill # cha
- Custody Seals on shipping container intact? NA ☐ Yes ☐ No ☒
 - Custody Seals dated and signed? NA ☐ Yes ☐ No ☒
 - Chain of Custody record present? Yes ☒ No ☐
 - Cooler temperature: _____ NA ☒ 5. Vermiculite/packing materials is NA ☒ Wet ☐ Dry ☐
 - Number of samples in shipping container: 3
 - Sample holding times exceeded? NA ☒ Yes ☐ No ☐
 - Samples have:
____ tape
____ custody seals
____ hazard labels
____ appropriate samples labels
 - Samples are:
____ in good condition
____ broken
____ leaking
____ have air bubbles
(Only for samples requiring head space)
 - Sample pH taken? NA ☒ pH < 2 ☐ pH > 2 ☐
 - Sample Location, Sample Collector Listed? * Yes ☐ No ☒
*For documentation only. No corrective action needed.
 - Were any anomalies identified in sample receipt? Yes ☐ No ☒
 - Description of anomalies (include sample numbers): _____

Sample Custodian: Hildeberg

Date: 4-23-02

Client Sample ID	Analysis Requested	Condition	Comments/Action

Client Informed on _____ by _____ Person contacted _____

☐ No action necessary; process as is.

Project Manager _____ Date _____

4/23/2002 3:00:40 PM

Sample Preparation/Analysis

PRIORITY

Balance Id: P8302-S

127642, BECHTEL HANFORD, INC.
Bechtel Hanford, Inc.

AW Gamma PrpRC5017

TA Gamma by HPGE

SI CLIENT: HANFORD

SEQUENTIAL

Pipet #: *N/A*Report Due: 05/14/2002 *W03700*

Sep1 DT/Tm Tech:

Batch: 2113395 FILTERS pCi/g

PM, Quote: BG1, 27038

Sep2 DT/Tm Tech:

SEQ Batch, Test: None

Prep Tech: *gr*

Work Order, Lot, Sample Date	Total Amt /Unit	Total Acidified/Unit	Initial Aliquot Amt/Unit	Adj Aliq Amt (Un-Acidified)	QC Tracer Prep Date	QC Vial 2 Prep Date	Count Time Min	Detector Id	Count On Off (24hr) Circle	CR Analyst, Init/Date
1 E0AM9-1-AA J2D230207-1-SAMP			25.49 / Sa.			25	200	63	0900	4/30/02
04/23/2002 10:30		AmtRec: 3XFILTERS		#Containers: 1		Scr Rst:		Alpha:		Beta:
2 E0AM9-1-AE-X J2D230207-1-DUP								6868	0903	5/1/02
04/23/2002 10:30		AmtRec: 3XFILTERS		#Containers: 1		Scr Rst:		Alpha:		Beta:
3 E0AW6-1-AA-B J2D230000-395-BLK			25.00					68	0900	4/30/02
04/23/2002 10:30		AmtRec:		#Containers: 1		Scr Rst:		Alpha:		Beta:
4 E0AW6-1-AC-C J2D230000-395-LCS			25.00		QCAG0766 PR.2/22/02 EX.8/27/02			67	0901	
04/23/2002 10:30		AmtRec:		#Containers: 1		Scr Rst:		Alpha:		Beta:

Comments: Don't count until 4-30-02. *gr 4-29-02*

All Clients for Batch:

127642, BECHTEL HANFORD, INC.

Bechtel Hanford, Inc.

, BG1, 27038

E0AM91AA-SAMP Constituent List:

Co-60	RDL:5.00E-02	pCi/g	LCL:	UCL:	RPD:	Cs-137	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:
Eu-152	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:	Eu-154	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:
Eu-155	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:						

E0AM61AA-BLK:

Co-60	RDL:5.00E-02	pCi/g	LCL:	UCL:	RPD:	Cs-137	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:
Eu-152	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:	Eu-154	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:
Eu-155	RDL:1.00E-01	pCi/g	LCL:	UCL:	RPD:						

E0AM61AC-LCS:

STL Richland

Key: In - Initial Amt, fi - Final Amt, di - Diluted Amt,

Richland Wa.

r - Reference date, ec-Enrichment Cell, ct-Cocktailed Added

Page 1

WO Cnt: 4

ICOC v4.5.3.2

4/23/2002 3:00:40 PM

Sample Preparation/Analysis

Balance Id:

AW Gamma PrpRC5017

TA Gamma by HPGE

PRIORITY

Pipet #:

Report Due: 05/14/2002

5I CLIENT: HANFORD

Sep1 DT/Tm Tech:

Batch: 2113395

pCi/g

Sep2 DT/Tm Tech:

SEQ Batch, Test: None

Prep Tech:

Work Order, Lot, Sample Date	Total Amt /Unit	Total Acidified/Unit	Initial Aliquot Amt/Unit	Adj Aliq Amt (Un-Acidified)	QC Tracer Prep Date	QC Vial 2 Prep Date	Count Time Min	Detector Id	Count On Off (24hr) Circle	CR Analyst, Init/Date
---------------------------------	--------------------	-------------------------	-----------------------------	--------------------------------	------------------------	------------------------	-------------------	----------------	---------------------------------	--------------------------

E0AM91AA-SAMP Calc Info:

Uncert Level (#s): 2 Decay to SaDt: Y Blk Subt.: N Sci.Not.: Y ODRs: B

E0AM61AA-BLK:

Uncert Level (#s): 2 Decay to SaDt: Y Blk Subt.: N Sci.Not.: Y ODRs: B

E0AM61AC-LCS:

Uncert Level (#s): 2 Decay to SaDt: Y Blk Subt.: N Sci.Not.: Y ODRs: B

5/2/02 1:56:10 PM

ICOC Fraction Transfer/Status Report

ByDate: 2/24/02, 5/3/02, Batch: '2113395', User: *All Order by BatchNbr,WorkOrderNbr,DateTimeAccepting

Q Batch	Work Ord	CurStatus	Accepting	Comments
2113395				
AC		InCnt1	WAGNERJ 4/29/02 2:55:58 PM	
SC		WagarR	IsBatched 4/23/02 3:00:31 PM	ICOC_RADCALC v4.5.3.2
SC		WAGNERJ	InPrep2 4/29/02 2:55:58 PM	RICH-RC-5017 REVISION 3
SC		DAWKINSO	InCnt1 4/29/02 4:09:32 PM	RICH-RD-0007 REVISION 2
AC		DAWKINSO	4/29/02 4:09:32 PM	

0028

AC: Accepting Entry; SC: Status Change

STL Richland
Richland Wa.

4/23/2002 3:00:39 PM

Sample Preparation/Analysis

Balance Id: PB302-S

127642, BECHTEL HANFORD, INC.
Bechtel Hanford, Inc.CI Sr-Total PrpRc5016, SepRC5006
TH Total Strontium by GPC
SI CLIENT: HANFORD**PRIORITY**

Pipet #:

Report Due: 05/14/2002

Sep1 DT/Tm Tech: 5-10-02 12:05

Batch: 2113394 FILTERS

pCi/g

PM, Quote: BG1, 27038

Sep2 DT/Tm Tech:

SEQ Batch, Test: None

Prep Tech: *zw*

Work Order, Lot, Sample Date/Time	Total Amt/Unit	Initial Aliquot Amt/Unit	QC Tracer Prep Date	QC Vial 2 Prep Date	Dish Size	Ppt or Geometry	Count Time Min	Detector Id	Count On Off (24hr) Circle	CR Analyst, Init/Date
1 E0AM9-1-AC J2D230207-1-SAMP		(Sa)				(mg)				
		SRTA7071 PR.3/15/02 EL.2/21/03			1.5"	86.9	100	28A	1956	5/10/2002-02
04/23/2002 10:30	AmtRec: 3XFILTERS	#Containers: 1				Scr Rst:		Alpha:		Beta:
2 E0AM9-1-AC J2D230207-1-DUP										
04/23/2002 10:30	AmtRec: 3XFILTERS	#Containers: 1				Scr Rst:		Alpha:		Beta:
3 E0AW4-1-AA-B J2D230000-394-BLK		(Sa) SRTA7072 PR.3/15/02 EL.2/21/03				88.6		28B	1956	5/10/2002-02
04/23/2002 10:30	AmtRec:	#Containers: 1				Scr Rst:		Alpha:		Beta:
4 E0AW4-1-AC-C J2D230000-394-LCS		STSB0579 PR.4/12/02 EL.2/21/03				83.9		28C	1956	5/10/2002-02
04/23/2002 10:30	AmtRec:	#Containers: 1				Scr Rst:		Alpha:		Beta:

Comments: Please remove dup. ²⁻³⁻⁰² 5-3-02

All Clients for Batch:

127642, BECHTEL HANFORD, INC.

Bechtel Hanford, Inc.

, BG1, 27038

E0AM91AC-SAMP Constituent List:

Sr-89/90	RDL:1.00E+00	pCi/g	LCL:	UCL:	RPD:
E0AM41AA-BLK:					
Sr-89/90	RDL:1.00E+00	pCi/g	LCL:	UCL:	RPD:

E0AM41AC-LCS:

E0AM91AC-SAMP Calc Info:

Uncert Level (#s): 2	Decay to SaDt: Y	Blk Subt.: N	Sci.Not.: Y	ODRs: B
----------------------	------------------	--------------	-------------	---------

E0AM41AA-BLK:

STL Richland
Richland Wa.Key: In - Initial Amt, fi - Final Amt, di - Diluted Amt,
r - Reference date, ec-Enrichment Cell, ct-Cocktailed Added

Page 1

WO Cnt: 4

ICOC v4.5.3.2

5/13/02 1:55:46 PM

ICOC Fraction Transfer/Status Report

ByDate: 4/13/02, 5/14/02, Batch: '2113394', User: 'All Order by BatchNbr,WorkOrderNbr,DateTimeAccepting

Q Batch	Work Ord	CurStatus	Accepting	Comments
2113394				
AC	InRev1	WAGNERJ	4/29/02 2:56:10 PM	
SC		WagarR	IsBatched	4/23/02 3:00:31 PM
SC		WAGNERJ	InPrep2	4/29/02 2:56:10 PM
SC		WAGNERJ	Prep2C	5/6/02 9:00:03 AM
SC		SteffensB	InSep1	5/10/02 11:02:57 AM
SC		SteffensB	Sep1C	5/10/02 1:54:06 PM
SC		DAWKINSO	Cnt1C	5/10/02 4:35:13 PM
SC		BlackCL	CalcC	5/11/02 8:15:33 AM
SC		KenitzerP	InRev1	5/13/02 10:33:47 AM
AC		WAGNERJ	5/6/02 9:00:03 AM	ICOC_RADCALC v4.5.3.2
AC		SteffensB	5/10/02 11:02:57 AM	RICH-RC-5016 REVISION 3
AC		SteffensB	5/10/02 1:54:06 PM	RICH-RC-5016 REVISION 3
AC		DAWKINSO	5/10/02 4:35:13 PM	RICH-RC-5006 REVISION 4
AC		BlackCL	5/11/02 8:15:33 AM	RICH-RC-5006 REVISION 4
AC		KenitzerP	5/13/02 10:33:47 AM	RICH-RD-0003 REVISION 2
				RICH-RD-0003 REVISION 2
				RICH-RC-0002 REVISION 5

0030

AC: Accepting Entry; SC: Status Change

STL Richland

Richland Wa.