

00000020

SAF-RC-001

Industrial Hygiene Sampling

FINAL DATA

NO DISTRIBUTION REQUIRED

COMMENTS:

SDG 051-4736-01 SAF-RC-001


Rad only ☒ Chem only Rad & Chem

☒ Complete Partial

300 Area 334A Bldg

RECEIVED
FEB 16 2006
EDMC

Report Identification Number: 05I-4736-01
 Subcontract Number: 0000X-BO-G0058-B-Mod#4
 Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby
 Laboratory Identification Number: DCHM
 SAF#: RC-001 / R33400 J451
 Payroll#: 72570



Sample Information

Sample Date	Customer Sample Number	Laboratory Sample Number	Method	Analytical Batch Identification	Sample Matrix
14 Nov 2005	J10JN9	05I44451	NMAM 7300M	G05BJ011	G WIPE
14 Nov 2005	J10JN8	05I44452	NMAM 7300M	G05BJ011	G WIPE
14 Nov 2005	J10JN7	05I44453	NMAM 7300M	G05BJ011	G WIPE
14 Nov 2005	J10JN6	05I44454	NMAM 7300M	G05BJ011	G WIPE

I certify that this electronic image and all hardcopies produced from this image accurately represent the data and are in compliance with the contract specific requirements, both technically and for completeness, other than the conditions detailed above or in the sample data package narrative. Release, by submission through email, the data contained in this electronic image and the computer-readable EDD (as applicable), has been authorized by the laboratory Manager or the Manager's designee.

Name: Lisa M. Reid
 Title: Chemist
 Date: November 18, 2005

Report Identification Number: 05I-4736-01
Subcontract Number: 0000X-BO-G0058-B-Mod#4
Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby
Laboratory Identification Number: DCHM
SAF#: RC-001 / R33400 J451
Payroll#: 72570

General Set Information: There are 4 samples in set 05I-4736-01 and 4 samples in set 05I-4737-01 for a total of 8 samples. The samples were analyzed for cadmium and beryllium on Ghost Wipe. No problems were encountered with the receipt of these samples and no contact with the CTR was required.

Method Summary: Samples were transferred to 50 ml centrifuge tubes and digested in the presence of 5 mL of nitric acid and 5 mL of ASTM Type II water. Samples were digested in a hot block set at 110°C for 60 minutes. Samples were then diluted to a 25 mL volume with ASTM Type II Water. Samples were shaken and delivered for ICP analysis.

Sample Preparation: All samples were prepared in accordance with DCL SOP "IH-AN-021" and NIOSH method NMAM 7300 modified for hot block digestion.

Holding Times: The holding times were met for both sample preparation and analysis.

Instrument Calibration: Instrument calibration was performed in accordance with NIOSH method NMAM 7300.

Initial and Continuing Calibration Verification Analysis: Beryllium and cadmium recoveries in all Initial Calibration Verification (ICV) and Continuing Calibration Verification (CCV) samples are within the quality control limits of $\pm 10\%$.

Initial and Continuing Calibration Blank Analysis: No beryllium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Contract Required Detection Limits (CRDL) of 0.01 ug/sample. No cadmium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Contract Required Detection Limits (CRDL) of 0.1 ug/sample.

Method Blank Analysis: No beryllium or cadmium was found in any of the media blank samples above the Contract Required Detection Limit (CRDL).

Dilution(s): None

Laboratory Control Sample and Duplicate Analysis: One Laboratory Control Sample (LCS) and one Laboratory Control Sample Duplicate (LCSD) were prepared and analyzed with the sample batch.

The LCS results were within the control limits of $\pm 20\%$. The Relative Percent Differences (RPDs) between the LCSs and the LCSDs were within the control limit of 20%.

Replicate Analysis: One sample in this batch was replicated. The RPDs between the samples and the replicates were within the control limit of 20%. If the result of the sample or replicate is below the CRDL, replicate analysis is negligible.

Flagging Codes: None

Nonconformance/Corrective Action Report (NC/CAR): N/A

Sample Calculation: The final results are calculated by the following equation:

Final result for aqueous samples ($\mu\text{g/sample}$) = (A) x (B) x (C)

Where:

A = Analyte concentration from instrument determination ($\mu\text{g/L}$)

B = Concentration factor from sample preparation

= $\frac{\text{Final Volume of Digestate (L)}}{\text{Sample}}$

C = Dilution performed at time of analysis

Example Calculation: $(1 \mu\text{g/L}) \times (0.025 \text{ L/sample}) \times (1) = 0.025 \mu\text{g/sample}$

Miscellaneous Comments: None.



Report Page

Page 4 of 7

Report Identification Number: 05I-4736-01
Subcontract Number: 0000X-BO-G0058-B-Mod#4
Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby
Laboratory Identification Number: DCHM
SAF#: RC-001 / R33400 J451
Payroll#: 72570

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Beryllium µg/sample		Cadmium µg/sample	
J10JN9	05I44451	17 Nov 2005	<0.01	U	<0.1	U
J10JN8	05I44452	17 Nov 2005	<0.01	U	<0.1	U
J10JN7	05I44453	17 Nov 2005	<0.01	U	<0.1	U
J10JN6	05I44454	17 Nov 2005	<0.01	U	<0.1	U
Limit of Detection (LOD)			0.01		0.1	
Required Detection Limit (RDL)						

U - Parameter not detected above LOD.

J - Parameter between LOD and RDL.

Report Identification Number: 05I-4736-01
Subcontract Number: 0000X-BO-G0058-B-Mod#4
Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby
Laboratory Identification Number: DCHM
SAF: RC-001 / R33400 J451
Payroll#: 72570

Batch ID: G05BJ011

QC Sample ID	QC Type	Analyte	Units	Result	Parent Result	Target	Percent Rec.	Relative Percent Diff.
BL-238337-1	MB	Beryllium	µg/sample	ND	NA	NA	NA	NA
BL-238337-1	MB	Cadmium	µg/sample	ND	NA	NA	NA	NA
QC-238337-1	LCS	Beryllium	µg/sample	11.0	NA	10.0	110.	NA
QC-238337-1	LCS	Cadmium	µg/sample	32.5	NA	30.0	108.	NA
QD-238337-1	LCSD	Beryllium	µg/sample	11.2	11.0	10.0	112.	2.36
QD-238337-1	LCSD	Cadmium	µg/sample	33.4	32.5	30.0	111.	2.50

MB - Method Blank

LCS - Laboratory Control Sample

LCSD - Laboratory Control Sample Duplicate

MS - Matrix Spike

MSD - Matrix Spike Duplicate

LD - Laboratory Duplicate

NA - Not Applicable

ND - Parameter not detected above LOD

LCS, LCSD Percent Rec. = (Result / Target) * 100.0

MS, MSD Percent Rec. = ((Result - Parent) / Target) * 100.0

LCS, LCSD Relative Percent Diff. = ((|LCS - LCSD|) / ((LCS + LCSD)/2.0)) * 100.

MS, MSD Relative Percent Diff. = ((|MS - MSD|) / ((MS + MSD)/2.0)) * 100.

LD Relative Percent Diff. = ((|Parent - LD|) / ((Parent + LD)/2.0)) * 100

WCH-SH-202 (08/29/2005)Page 1 of 2

Web Page: www.datachem.com
E-mail: lab@datachem.com

Enter on line below the first Sample Number from Page One:

J10JN9

CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST			
SIGN / PRINT NAMES / USE MILITARY TIME			
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
CJ Williams Cynthia Williams	11-14-05 / 1615	3746 Biolb RM 16 locked cabinet	11-14-05 / 1615
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
David St. John / Patrick Vidlet	11-16-05 1430	David St. John	11/16/05 1430
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
David St. John w/CH David St. John	11/16/05 1500	David St. John	11/16/05 1500
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
Teddy		Merick Spelman	11/16/05
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
Merick Spelman	11/16/05		
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
Relinquished By/To:	DATE/TIME	Received By/From:	DATE/TIME
LABORATORY SECTION	Received By	Title	DATE/TIME
	Merick Spelman		11/16/05

REVIEWED BY: _____ DATE: _____
PRINT/SIGN NAME