SAF-RC-001 Industrial Hygiene Sampling FINAL DATA

NO DISTRIBUTION REQUIRED

COMMENTS:

SDG 06I-0423-01 SAF-RC-001

Rad only X Chem only Rad & Chem

X Complete Partial

300 Area 333 Bldg







Cover Page

Report Identification Number:

06I-0423-01

Subcontract Number:

0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist:

Denise Pitts / Henry Ruby

Laboratory Identification Number:

DCHM

SAF#: RC-001 / r37120 2000

Payroll#: 72520

Web Page: www.datachem.com

E-mail: lab@datachem.com

Sample Information

Sample Date	Customer Sample Number	Laboratory Sample Number	Method	Analytical Batch Identification	Sample Matrix
01 Feb 2006	J113B4	06103387	NMAM 7300M	G061200Y	MCE
01 Feb 2006	J113B5	06103388	NMAM 7300M	G061200Y	MCE
01 Feb 2006	J113C6	06I03389	NMAM 7300M	G061200Y	MCE
01 Feb 2006	J113C7	06I03390	NMAM 7300M	G061200Y	MCE
01 Feb 2006	J113B6	06I03391	NMAM 7300M	G061200Y	MCE
01 Feb 2006	J113B7	06103392	NMAM 7300M	G061200Y	MCE

I certify that this electronic image and all hardcopies produced from this image accurately represent the data and are in compliance with the contract specific requirements, both technically and for completeness, other than the conditions detailed above or in the sample data package narrative. Release, by submission through email, the data contained in this electronic image and the computer-readable EDD (as applicable), has been authorized by the laboratory Manager or the Manager's designee.

Phone: (801) 266-7700

FAX: (801) 268-9992

Name: Joanna C. Sanchez

Title: Chemist

Date: February 06, 2006



Web Page: www.datachem.com

E-mail: lab@datachem.com



Case Narrative Page

Report Identification Number: 06I-0423-01

Subcontract Number: 0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise Pitts / Henry Ruby

Laboratory Identification Number: DCHM

SAF#: RC-001 / r37120 2000

Payroll#: 72520

There are 6 samples in set 06I-0423-01 which were analyzed for beryllium, General Set Information: lead and cadmium on MCE filter. No problems were encountered with the receipt of these samples and no contact with the CTR was required.

Method Summary: Samples were transferred to 50 ml centrifuge tubes and digested in the presence of 10 mL of 1:1 (v/v) nitric acid. Samples were digested in a hot block set at 110°C for 40 minutes. Samples were then diluted to a 25 mL volume with ASTM Type II Water. Samples were shaken and delivered for ICP analysis.

Sample Preparation: All samples were prepared in accordance with DCL SOP "IH-AN-021" and NIOSH method NMAM 7300 modified for hot block digestion.

Holding Times: The holding times were met for both sample preparation and analysis.

Instrument Calibration: Instrument calibration was performed in accordance with NIOSH method NMAM 7300.

Initial and Continuing Calibration Verification Analysis: Beryllium, cadmium and lead recoveries in all Initial Calibration Verification (ICV) and Continuing Calibration Verification (CCV) samples are within the quality control limits of +/- 10%.

Initial and Continuing Calibration Blank Analysis: No beryllium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Limit of Quantitation (LOQ) of 0.01 ug/sample. No cadmium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Limit of Quantitation (LOQ) of 0.08 ug/sample. No lead results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Limit of Quantitation (LOQ) of 1. ug/sample.

Method Blank Analysis: No beryllium, cadmium or lead was found in the media blank sample above the Contract Required Detection Limit (CRDL).

Phone: (801) 266-7700

FAX: (801) 268-9992

Dilution(s): NA.



Web Page: www.datachem.com

E-mail: lab@datachem.com



Case Narrative Page

Laboratory Control Sample and Duplicate Analysis: One Laboratory Control Sample (LCS) and one Laboratory Control Sample Duplicate (LCSD) were prepared and analyzed with the sample batch. The LCS result was within the control limit of +/- 20%. The Relative Percent Differences (RPD) between the LCS and the LCSD was within the control limit of 20%.

Replicate Analysis: One sample was replicated with this analysis run. The RPD between the sample and the replicate was within the control limit of 20%. If the result of the sample or replicate is below the CRDL, replicate analysis is negligible.

Flagging Codes: None

Nonconformance/Corrective Action Report (NC/CAR): N/A

Sample Calculation: The final results are calculated by the following equation:

Final result for aqueous samples (μ g/sample) = (A) x (B) x (C)

Where:

A = Analyte concentration from instrument determination ($\mu g/L$)

B = Concentration factor from sample preparation

= Final Volume of Digestate (L)

Sample

C = Dilution performed at time of analysis

Example Calculation: $(1 \mu g/L) \times (0.025 L/sample) \times (1) = 0.025 \mu g/sample$

Phone: (801) 266-7700

FAX: (801) 268-9992

Miscellaneous Comments: None.



Report Page

Report Identification Number:

06I-0423-01

Subcontract Number:

0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise Pitts / Henry Ruby

Laboratory Identification Number: DCHM

SAF#: RC-001 / r37120 2000

Payroll#: 72520

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Berylliun µg/sample		Berylliun μg/m³	1	Air Volume L
J113B4	06I03387	03 Feb 2006	< 0.01	U	< 0.047	U	211.
J113B5	06I03388	03 Feb 2006	< 0.01	U	< 0.024	U	410.
J113C6	06I03389	03 Feb 2006	< 0.01	U	< 0.027	U	367.
J113C7	06I03390	03 Feb 2006	< 0.01	U	< 0.026	U	378.
J113B6	06I03391	03 Feb 2006	< 0.01	U	**		**
J113B7	06I03392	03 Feb 2006	< 0.01	U	**		**
	Limit of I	· 0.	.01			•	
R	equired Detecti						

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Lead µg/sample	e	Lead μg/m³		Cadmiun µg/sampl	
J113B4	06I03387	03 Feb 2006	<1.	U	<4.7	U	< 0.08	U
J113B5	06I03388	03 Feb 2006	<1.	U	<2.4	U	< 0.08	U
J113C6	06I03389	03 Feb 2006	1.1		2.9		< 0.08	U
J113C7	06I03390	03 Feb 2006	<1.	U	<2.6	U	<0.08	U
J113B6	06I03391	03 Feb 2006	<1.	U	**		<0.08	U
J113B7	06I03392	03 Feb 2006	<1.	U	**		<0.08	U
	Limit of I		1.			0.	.08	
R	equired Detect	ion Limit (RDL)						

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Cadmium µg/m³
J113B4	06103387	03 Feb 2006	<0.38 U
J113B5	06103388	03 Feb 2006	<0.20 U
J113C6	06103389	03 Feb 2006	<0.22 U
J113C7	06I03390	03 Feb 2006	<0.21 U
J113B6	06I03391	03 Feb 2006	**
J113B7	06I03392	03 Feb 2006	**

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com

E-mail: lab@datachem.com



Web Page: www.datachem.com E-mail: lab@datachem.com



Report Page

Phone: (801) 266-7700

FAX: (801) 268-9992

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Cadmium µg/m³
R			

U - Parameter not detected above LOD.

J - Parameter between LOD and RDL.





QC Summary Page

Report Identification Number: 06I-0423-01

Subcontract Number: 0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise Pitts / Henry Ruby

Laboratory Identification Number: DCHM

SAF: RC-001 / r37120 2000

Payroll#: 72520

Batch ID: G061200Y

	QC				Parent		Percent	Relative Percent
QC Sample ID	Type	Analyte	Units	Result	Result	Target	Rec.	Diff.
BL-240883-1	MB	Beryllium	μg/sample	ND	NA	NA	NA	NA
BL-240883-1	MB	Lead	μg/sample	ND	NA	NA	NA	NA
BL-240883-1	MB	Cadmium	μg/sample	ND	NA	NA	NA	NA
QC-240883-1	LCS	Beryllium	μg/sample	10.8	NA	10.0	108.	NA
QC-240883-1	LCS	Lead	μg/sample	106.	NA	100.	106.	NA
QC-240883-1	LCS	Cadmium	μg/sample	33.8	NA	30.0	113.	NA
QD-240883-1	LCSD	Beryllium	μg/sample	10.9	10.8	10.0	109.	1.46
QD-240883-1	LCSD	Lead	μg/sample	107.	106.	100.	107.	0.845
QD-240883-1	LCSD	Cadmium	μg/sample	34.3	33.8	30.0	114.	1.48

MB - Method Blank

LCS - Laboratory Control Sample

LCSD - Laboratory Control Sample Duplicate

MS - Matrix Spike

MSD - Matrix Spike Duplicate

LD - Laboratory Duplicate

NA - Not Applicable

ND - Parameter not detected above LOD

LCS, LCSD Percent Rec. = (Result / Target) * 100.0

MS, MSD Percent Rec. = ((Result – Parent) / Target) * 100.0

LCS, LCSD Relative Percent Diff. = ((|LCS - LCSD|) / ((LCS + LCSD)/2.0)) * 100.

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com


E-mail: lab@datachem.com

MS, MSD Relative Percent Diff. = ((|MS - MSD|) / ((MS + MSD)/2.0)) * 100.

LD Relative Percent Diff. = ((|Parent - LD|) / ((Parent + LD)/2.0)) * 100



COC Page



Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com

E-mail: lab@datachem.com



COC Page

CHAIN OF CUSTODY	CHAIN OF CUSTODY SAMPLE ANALYSIS REOUEST	-
	SIGN / PRINT NAMES / USE MILITARY TIME	
Authorities By Newson		,
CLAN & ALPHUNG My CHATTAGE WILLIAMS . 21-56 / 1630	D 3746 BIOG am 16 locked cubinst 2-1-66	56 /1633
Joséed about blog 3746 Ruth		,
Coffee (Jeldie Melhon 02:02:06 150	12 Stoffly 12. Styll 2-2-06	0051/ 91
New Authority Annual Control of the	INITIAL PROPERTY OF THE PROPER	
Fed Ex	the Worldon 2131c	213106 1030
Models 3 by	Table 1 in the second s	is Si Si
E temporal by the J	Because By Branck	
Married by Stead	Income in Supers	
Potegrical by ward	Pared B.Smit Daily (1985	
Respond to South	Exercit Informat	
Richard Libywood	Burnel by Great	
LABORATORY ROSPING BY COLOUTED THE SECTION	1872	2/8/de (US)
REVIEWED BY:	DATE	
FRIN ISKR NAME		
WCH-SH-302 (IM2Y2UU)		Page 2 a Q_

DataChem Laboratories, Inc. 960 West Levoy Drive Salt Lake City, Utah 84123-2547

Enter on line below the first Sample Number from Page One:

Phone: (801) 266-7700 FAX: (801) 268-9992 Web Page: www.datachem.com E-mail: lab@datachem.com



77. H		16.		OF C	USTODY/S.			JIS KE				Data Tues		
Collector:	illiams	Company C Denise A. P	ontact itts and Henry W	. Ruby		531-1229 SPECIAL INSTRUCTIONS			Project Coordinator Joan H. Kessner SAF No.			Data Turriaround 24 hs		
Payroli #:	-	Sampling L												
70 S Type of Sample:	<u> </u>	300 ay	la/333.	B106	All relevant CO	As must be pr	pyided:	-	RC-001 Method of S	hinment		<u> </u>		
	onals	als							MICHIGA OF 2	Fe	d E	X		
Shipped To: Data U	0.000	Wipe Samp	le Media: IYes □ No		NIOSN 7300 Bill of Lading/Air Bill No.									
5.14101.5	1 1 17		Tie2 []140						RS	-44	9439	5 Y	1726	•
SVITTALL C	IAZARD/RI MARKS	MAT	RIX		<u> </u>		ļ	_ .	7	-/_ / -	<u> </u>		<u> </u>	<u> </u>
Be/co		A AD WI WI			n (i.e., cooling red, etc.)	No	No	No	No	No	No	No	No	Λ
Special Handling an) / O	x or		ıcılım	cu. cit./		"		""	140		14.		
		SAMPLE ANALY	SIS			Ashestos	Lead	Beryllius	n Berylliu		Lead	Cd	Ca	
SAMPLE NO.	MATRIX	SAMPLE DATE	VOLUME (L)		Comments	Airborne	Airborne	Airborne			Wipe	Wipe	Airbonse	
J113B4	A	2-1-06	211				X	Х	X	1 m			X	
511385	4	1	410		\		X	X		1			X	
5113cb			367		CIM		X	X			CIN		X	2,01
511307			378		2-1-06		X	X		2-1	06		Χ	
عا <u>8 3 11 3</u>	V		Nla				Х	X					X	
311387	A	2-1-06	Nla				X	X	\@\((1)	3 0			X	
									6					
					<u>CJ</u>	ى ك					ELD S	AMDI	COE	N
					2-1-	Ob								
VCH-SH-202 (08/29	/2005)			-1.			-	4			Page	100)	

Enter on line below	the first Sample Number from Page	: One:								
J113B	<u> </u>									
CHAIN OF CUSTODY/SAMPLE ANALYSIS REQUEST										
SIGN / PRINT NAMES / USE MILITARY TIME										
Relangueshed By/Sugerl		DATE / FINE	Record By Starat:	DATE (194L						
Cunham	Means/continuwilliam	5 2-1-86/1630	3746 BIOG AM 16 locked cubinet	2-1-06/1638						
Retrictinged By Spaced	Swet bldg 3746 R	4-HOLTE / TIME	Reus val By Seeral	DAR FRIED						
Costale	Man Cynthawillian Sout bldg 3746 R	DL-62-06 (580)	RZ Steffler P. Steffle	2-2-06 /1500						
RZ Steff	Ver R. 2. Stills	2-2-06 /1600	Fed EX	DATE / CIMI.						
Statement decid by America		DATE/TIME	Rounnel By State :	DAJE / HME.						
Reducipe deal by Securi		DATE / TIME.	Becaused By-Mount	PAILZIML						
R. langer-in, d by //www.il		DATE! FIME	Reacted By Steach	DATE CIM						
Reference the 18 y Stored		DATE / TIME	Rosavai 8y/Surrat	DAIL (1888)						
Retungarshed By Assect		DATE / 11ME	Received By Serret	DATE / CIME						
Retiniprobad By Sieri d		DATE / TIME	Recoved By/Stanot:	DATE / UMB						
Raimga h.d Bythre.d		DATE: TIME	Received by/Stored:	DATE / CIME						
LABORATORY SECTION	Received By	Title	<u> </u>	DATE / TIME						

DATE:

PRINT/SIGN NAME

REVIEWED BY: