# SAF-RC-001 Industrial Hygiene Sampling FINAL DATA

NO DISTRIBUTION REQUIRED

**COMMENTS:** 

SDG <u>06I-0618-01</u> SAF-RC-001

Rad only X Chem only Rad & Chem

X Complete

Partial

300 Area 303M Bldg





### Cover Page

Page 1 of 7

rovd 2/16/06

Web Page: www.datachem.com

E-mail: lab@datachem.com

Report Identification Number: 06I-0618-01

Subcontract Number:

0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist:

Denise A. Pitts / Henry W. Ruby

Laboratory Identification Number: DCHM

SAF#: RC-001 / R303MO J451

Payroll#: 73338

Sample Information

Sample Date	Customer Sample Number	Laboratory Sample Number	Method	Analytical Batch Identification	Sample Matrix
09 Feb 2006	J114R9	06104907	NMAM 7300M	G061F013	G WIPE
09 Feb 2006	J114T0	06104908	NMAM 7300M	G061F013	G WIPE
09 Feb 2006	J114T1	06104909	NMAM 7300M	G061F013	G WIPE
09 Feb 2006	J114T2	06I04910	NMAM 7300M	G061F013	G WIPE
09 Feb 2006	J114T3	06I04911	NMAM 7300M	G061F013	G WIPE

I certify that this electronic image and all hardcopies produced from this image accurately represent the data and are in compliance with the contract specific requirements, both technically and for completeness, other than the conditions detailed above or in the sample data package narrative. Release, by submission through email, the data contained in this electronic image and the computer-readable EDD (as applicable), has been authorized by the laboratory Manager or the Manager's designee.

Phone: (801) 266-7700

FAX: (801) 268-9992

Name: Joanna C. Sanchez

Title: Chemist

Date: February 15, 2006





#### Case Narrative Page

Report Identification Number: 06I-0618-01

Subcontract Number: 0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby

Laboratory Identification Number: DCHM

SAF#: RC-001 / R303MO J451

Payroll#: 73338

General Set Information: There are 4 samples in set 06I-0561-02, 7 samples in set 06I-0564-01, 4 samples in set 06I-0565-01, 4 samples in set 06I-0566-01, and 5 samples in set 06I-0618-01 which were analyzed for cadmium, lead and beryllium on Ghost Wipe. No problems were encountered with the receipt of these samples and no contact with the CTR was required.

Method Summary: Samples were transferred to 50 ml centrifuge tubes and digested in the presence of 5 mL of nitric acid and 5 mL of ASTM Type II water. Samples were digested in a hot block set at 110°C for 60 minutes. Samples were then diluted to a 25 mL volume with ASTM Type II Water. Samples were shaken and delivered for ICP analysis.

Sample Preparation: All samples were prepared in accordance with DCL SOP "IH-AN-021" and NIOSH method NMAM 7300 modified for hot block digestion.

**Holding Times:** The holding times were met for both sample preparation and analysis.

Instrument Calibration: Instrument calibration was performed in accordance with NIOSH method NMAM 7300.

Initial and Continuing Calibration Verification Analysis: Beryllium, cadmium and lead recoveries in all Initial Calibration Verification (ICV) and Continuing Calibration Verification (CCV) samples are within the quality control limits of +/- 10%.

Initial and Continuing Calibration Blank Analysis: No beryllium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Contract Required Detection Limits (CRDL) of 0.02 ug/sample. No cadmium results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Contract Required Detection Limits (CRDL) of 0.07 ug/sample. No lead results were found in the Initial Calibration Blank (ICB) or Continuing Calibration Blanks (CCB) at levels above the Contract Required Detection Limits (CRDL) of 2. ug/sample.

Method Blank Analysis: No beryllium, cadmium or lead was found in any of the media blank samples above the Contract Required Detection Limit (CRDL).

**Dilution(s):** None.

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com E-mail: lab@datachem.com





#### Case Narrative Page

Laboratory Control Sample and Duplicate Analysis: Two Laboratory Control Samples (LCSs) and two Laboratory Control Sample Duplicates (LCSDs) were prepared and analyzed with the sample batch. The LCS results were within the control limits of +/- 20%. The Relative Percent Difference (RPD) between the LCSs and the LCSDs were within the control limit of 20%.

Replicate Analysis: Four samples in this batch were replicated. The RPD between the samples and the replicates was within the control limit of 20%. If the result of the sample or replicate is below the CRDL, replicate analysis is negligible.

Flagging Codes: None

Nonconformance/Corrective Action Report (NC/CAR): N/A

<u>Sample Calculation:</u> The final results are calculated by the following equation: Final result for aqueous samples ( $\mu$ g/sample) = (A) x (B) x (C)

Where:

A = Analyte concentration from instrument determination ( $\mu$ g/L)

B = Concentration factor from sample preparation

= Final Volume of Digestate (L)

Sample

C = Dilution performed at time of analysis

Example Calculation:  $(1 \mu g/L) \times (0.025 L/sample) \times (1) = 0.025 \mu g/sample$ 

Miscellaneous Comments: The LOQ for sample 06I04583 for cadmium is 0.6 because of interference.

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com

E-mail: lab@datachem.com





### Report Page

Report Identification Number: 06I-0618-01

Subcontract Number: 0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby

Laboratory Identification Number: DCHM

SAF#: RC-001 / R303MO J451

Payroll#: 73338

Customer Sample Number	Laboratory Sample Number	Date Analyzed	Beryllium µg/sample		Cadmium µg/sample	1	Lead µg/sampl	e
J114R9	06104907	14 Feb 2006	< 0.02	U	< 0.07	U	<2.	U
J114T0	06104908	14 Feb 2006	< 0.02	U	< 0.07	U	<2.	U
J114T1	06104909	14 Feb 2006	< 0.02	Ū	< 0.07	Ü	<2.	U
J114T2	06104910	14 Feb 2006	< 0.02	U	< 0.07	U	<2.	U
J114T3	06104911	14 Feb 2006	< 0.02	U	< 0.07	U	<2.	U
	Limit of I	0.	02	0.	07		2.	
R	equired Detecti	ion Limit (RDL)						

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com

E-mail: lab@datachem.com

U - Parameter not detected above LOD.

J - Parameter between LOD and RDL.





### QC Summary Page

Report Identification Number: 06I-0618-01

Subcontract Number: 0000X-BO-G0058-B-Mod#4

Name of Industrial Hygienist: Denise A. Pitts / Henry W. Ruby

Laboratory Identification Number: DCHM

SAF: RC-001 / R303MO J451

Payroll#: 73338

Batch ID: G061F013

[								Relative
ĺ	QC				Parent		Percent	Percent
QC Sample ID	Type	Analyte	Units	Result	Result	Target	Rec.	Diff.
BL-241305-1	MB	Beryllium	μg/sample	ND	NA	NA	NA	NA
BL-241305-1	MB	Cadmium	μg/sample	ND	NA	NA	NA	NA
BL-241305-1	MB	Lead	μg/sample	ND	NA	NA	NA	NA
QC-241305-1	LCS	Beryllium	μg/sample	11.3	NA	10.0	113.	NA
QC-241305-1	LCS	Cadmium	μg/sample	32.9	NA	30.0	110.	NA
QC-241305-1	LCS	Lead	μg/sample	104.	NA	100.	104.	NA_
QD-241305-1	LCSD	Beryllium	μg/sample	11.0	11.3	10.0	110.	3.28
QD-241305-1	LCSD	Cadmium	μg/sample	31.8	32.9	30.0	106.	3.29
QD-241305-1	LCSD	Lead	μg/sample	102.	104	100.	102.	2.39

MB - Method Blank

LCS - Laboratory Control Sample

LCSD - Laboratory Control Sample Duplicate

MS - Matrix Spike

MSD - Matrix Spike Duplicate

LD - Laboratory Duplicate

NA - Not Applicable

ND - Parameter not detected above LOD

LCS, LCSD Percent Rec. = (Result / Target) \* 100.0

MS, MSD Percent Rec. = ((Result - Parent) / Target) \* 100.0

LCS, LCSD Relative Percent Diff. = ((|LCS - LCSD|) / ((LCS + LCSD)/2.0)) \* 100.

Phone: (801) 266-7700

FAX: (801) 268-9992

Web Page: www.datachem.com

E-mail: lab@datachem.com

MS, MSD Relative Percent Diff. = ((|MS - MSD|)/((MS + MSD)/2.0)) \* 100.

LD Relative Percent Diff. = ((|Parent - LD|) / ((Parent + LD)/2.0)) \* 100




# **COC** Page

					-v, ·															···		
ing.						• .	1	٠.				نا	2,04				_			_/		
		Town Town		ų.	. j	í	0///	2	_	V.	1			_				_		4		4
		Data Ternament		24 hr		, , ,		.2	!	ಶಕ್ಷ		×	7	×	×	×			90			
					ર્ફ	9.7	7435	- 5	2	i i		×	×	×	×	×	_	Ä	-	_	_	
\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\		ration			EX CACAS	5	8544	: 	2		<u>a</u>		2-9-0	ماد			$\overline{A}$	-			-	
~~ ~		Project Coordinator	Joan H. Kessner	SAF No.	Method of Sin	Bill of Lading/Air Bill No	9	;	£	Beryllian		×	×	×	×	×/	<u>/                                    </u>					
M. DONAN		A CIC	-		<u> </u>	1			ž	Berylline		Cathron-	80	98	2	7						
	7.1.1	ANALI	8		, s	in Section of the sec			2	_	Airbame	New Year			3.	<del>(-43</del> )						
		AMPLE	531-1229	TRUCTIONS As must be p	ISHC	7300 B			2		Airhope								i. ,			
	1	CHAIN OF CUSTODY/SAMPLE ANALTS SAMPLES		SPECIAL INSTRUCTIONS All referred COAs must be provided:	REGENOME	ANALYSIS METHOD (SPECLFIL): Niosh 7300 Be, Pb, cd		(i.e. cooling	required, etc.)	*	Comments	ہد	Blak	2	4	X/7	; ]					
		OFCU	V. Ruby	1,1	<u> </u>			Powervation	and the same		<u>.</u> 1	12   12   12   12   12   12   12   12	18	2			<del></del>			1		
			tect s and Henry A	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Cons wronk)	Medic Section	2 <u>4</u>	Ħ	品質	9	VOLUME (L) or Aven 100 cm²	7	7		8 5	9			-	<u>/_</u>	-	
			Company Contact Desire A. Pitts and Henry W. Roby	Sempling Location	\$ 000 1000 1000 1000 1000 1000 1000 100	Wipe Sample Media:	E	MATERIA	X WIFE	SAMPLE ANALYSIS	SAMPLE DATE	4				2-5-06			3/3			
				4	477	42		DREMARKS	gina		MATSOX		3	<del>                                     </del>	_	3					:	(500
Eg			Collector	Payroll #:	TYSS &	Shipped To	Dute Che	POSSBEE SAMPLE HAZARDREM	Special Hamiling and/or Storage	ď.	SANGTLE NO.	1.	•	OTPIN.		X 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	1 /					WCH-SH-202 (08/29/2005)



## COC Page



DataChem Laboratories, Inc. 960 West Levoy Drive Salt Lake City, Utah 84123-2547 Phone: (801) 266-7700 FAX: (801) 268-9992

Web Page: www.datachem.com E-mail: lab@datachem.com



			CHAIN C	OF CUSTODY/S	AMPLE	ANALY	SIS RI	EQUEST					
Chips Sch	uby	Project Coordin Joan H. Kessner		Data Turnaround									
Payroll#:	All relevant COAs must be provided:				SAF No. RC-001 24 h~								
Type of Sample: Be Pb Cd	- 1 / / / / / / / / / / / / / / / / / /			R303m		neoriio).		Method of Ships Federal					
Shipped To:		Wipe Samp	le Media:	ANALYSIS	7300 PS	•		Bill of Lading/	ir Bill No.				
Duter che Sult La		Other	Yes No	MOSH	7300 101	-, . 0,		85	T44	943	5	4770	770
POSSIBLE SAMPLE	HAZARD/REMARKS	MAI			1			<u> </u>	1			1	
Special Handling ar		A - AII WI - WI X - OT	PE	eservation (i.e., cooling required, etc.)	No	No	No	No	No	No	No	No	
		SAMPLE ANALY	YSIS			11	Ddiio	Parding.		Lead	Cd	ga	
SAMPLE NO.	MATRIX	SAMPLE DATE	VOLUME (L) or Ares 100 cm <sup>2</sup>	Comments	Asbertos Airbone	Lead Airborne	Berylliu Airborr		Noid	Wipe	Wipe	Airborne	
3114 R9	wi_	2-9-06	N/2	Black				X		Χ	X		
211470			N/A	Blak				X	الم	X	X	7	
311471			100	N/A				Х	Š	X	Χ		1
211472			100	N/A		n/c		X		X	X		
14473	w <sub>1</sub>	2-4-06	100_	N/A		\$		X		\_X_	X		
						\ 							
									ELD !	SAMP	LE C	OPY	
ſ		3.7.04								5%	2/06		
									<u>.</u>				

2114R9			14R9	711
--------	--	--	------	-----

	C	HAIN OF CUSTODY/SAN	IPLE ANALYSI	S REQUEST	
		SIGN / PRINT NAMES	<del></del>		
Retinquished By/Stored:		DATE/TIME	Received By/Staned:		DATE/TOGE
CL STORY	Jamis Sch. My 2.	-9-06 1615	Cocked Colema	1 Rullo 3746 Bldg 2-9-	06 1615
locked cale	anet bldg 3746 Ru#	C PALETIME	necessed by/station:	'00/ 1/ 1/	DATE/TIME
Caldy Malh	a Goldie Halhan	02-13-06 / 1430	KZ Ste	11/0 1.7. Stiff	2-13-06 / 1430
RZ Steffle	A. Still	2+13-06 / 1600	Fed E	- / // // · ·	DATE / TIME
Relinquished By/Stered:		DATE / TIME	Received By/Stared:		DATE/TIME
Relinquished By/Stored		DATE / TIME	Received Dy/Stores:		DATE / TIME
Relinquished By/Sqrpt:		DATE / TIME	Received By/Stead:		DATE / TIME
Relinquahed By/Stored:		DATE / TIME	Received By/Stored:	·	DATE / TIME
Relimquished By/Stored		DATE/TIME	Received By/Stand:		DATE/ITME
Retinquished By/Stored:	·	DATE/IME	Received By/Stared:		DATE/TIME
Relinquished By/Stored:		DATE/TIME	Received By/Stered:		DATE/TIME
LABORATORY SECTION	Received By	Title			DATE / TIME
REVIEWED BY:			· <del>-</del> -	DATE:	
		PRINT/SIGN NAME			